TRACHEO SET

Safe & Easy Tracheostomy

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TRACHEO SET

Safe & Easy Tracheostomy

The innovative forceps in Tracheostomy

The dilating forceps that allows the combination of the Griggs and the Ciaglia technique, gaining their full advantages and eliminating their disadvantages

FEATURES

- The forceps works through compression and not traction
- A valve allows the passage of the guide wire, without losing it
- The shape of the forceps is similar to a mono-dilator

ADVANTAGES

- Better control during the stoma formation
- Possibility to adjust the maximum opening of the forceps
- Fast and simple procedure for cannula insertion

TRACHEO SET

Istruzioni per l'uso (IT)	pag. 1
Instruction for use (EN)	pag. 5
Instrucciones de uso (ES)	pag. 9
تعليمات التشغيل (AR)	pag. 13
Инструкции по использованию (RU)	pag. 17
Kullanım için talimatlar (TR)	pag. 22
Gebrauchsanweisung (DE)	pag. 26
Mode d'emploi (FR)	pag. 30
Upute za korištenje (HR)	pag. 35
Hướng dẫn sử dụng (VIE)	pag. 39
Használati utasítás (HU)	pag. 43

INSTRUCTION FOR USE (EN)

Contents:

Product Description	5
Intended use	6
Warnings	6
Method of use	7
Preparing the components	7
Procedural technique	7
Other useful information	8

Please read the instructions for use carefully before using the medical device TRACHEO SET. X-MED S.r.l. guarantees the sterility and performance of the medical device until the expiry date stated on the label of each individual product. X-MED S.r.l. shall not, under any circumstances, be liable for property and/or people in the event of use not in compliance with the contents of these instructions for use or if used by unqualified medical personnel.

The device is sterilised using ethylene oxide. **STERILE EO**

CAUTION: The medical device TRACHEO SET must be used exclusively by qualified medical professionals with expertise in resuscitation and intensive care. Do not use the product if the packaging is damaged. **TRACHEO SET is a disposable medical device and it may not be reused in any way, in order to avoid the risk of infection for the patient.**

Product Description


The TRACHEO SET medical device is a disposable device for percutaneous dilated tracheostomy consisting of a series of accessories necessary for performing of the operation. TRACHEO SET is available in two versions: with or without dilatation forceps (Tracheo Forceps®).

The tracheostomy tube is available in 5 sizes and 2 different models in order to be able to adapt to the varying requirements of patients. The complete kit consists of the following components:

- 1) Tracheo Forceps® dilating forceps (if included in the pack);
- 2) Tracheostomy tube available with or without a subglottic aspiration line and in different sizes*, with a standard connector (15 mm);
- 3) Tracheostomy tube introducer;
- 4) Syringe (10 mL and 20 mL versions available);
- 5) Scalpel blade;
- 6) Cannula needle (14 G);
- 7) Stiff blue dilator;
- 8) Soft white dilator;
- 9) Seldinger guide wire;
- 10) Tracheostomy tube holder.

Sizes available for the internal diameter of tracheostomy tubes*

Tracheostomy tube (Int. Diam.)	Tracheostomy tube with subglottic suction line (Int. Diam.)
7 mm	7 mm
7.5 mm	7.5 mm
8 mm	8 mm
8.5 mm	8.5 mm
9 mm	9 mm

The above package drawing is for illustrative purposes only to show where the elements that compose the TRACHEO SET are positioned. The layout of the components may differ from the drawing depending on the REF code of TRACHEO SET ordered.

The **Tracheo Forceps®** dilating forceps was designed for a quicker and safer tracheostomy than that performed using the Ciaglia and Griggs techniques.

The special design of the forceps overcomes the drawbacks associated with the Griggs technique, because they open by compression rather than traction. This feature affords greater control during stoma formation. The forceps also features an opening lock, making it possible to establish the maximum width of the stoma to be created while avoiding tracheal damage.

The distal part of the forceps has a curved conical conformation that acts as a dilator during stoma formation. There is a hole on the right arm of the forceps for the guide wire, a feature that keeps it attached to the forceps during dilation.

Intended use

The TRACHEO SET medical device is intended for carrying out elective percutaneous dilational tracheostomies for respiratory tract management in adult patients only. The procedure described must be carried out in a controlled environment, such as an ICU or operating theatre, by medical personnel with expertise in resuscitation and intensive care.

Warnings

- Do not use the tracheostomy tube on patients with anatomical respiratory tract anomalies, as this could result in a partial or total obstruction of the airways;
- Do not perform the tracheostomy procedure in the presence of skin infections;
- Do not use the TRACHEO SET if the packaging is damaged;
- The components of the TRACHEO SET device are supplied sterile: do not reuse or resterilise the products and dispose of them in accordance with the applicable regulations, making sure cutting and sharp parts are disposed of safely.
- The TRACHEO SET must only be used by medical personnel trained to perform a percutaneous dilated tracheostomy with skills in resuscitation and intensive care;
- Apply strict aseptic techniques during the placement of the device;
- Before removing the tracheostomy tube, deflate the cuff completely to avoid damaging the trachea and/or stoma;

- Cuff pressure must be monitored constantly. Excessively high pressure could damage the trachea;
- Do not subject the patient to MRI because the cannula supplied in the package is not compatible with the MRI environment;
- The tracheostomy tubes must be replaced on a regular basis in order to meet the specific needs of the patient and in any case following the instructions of the medical director;
- During and after connection of the breathing system to the tracheostomy tube connector, avoid exerting excessive torsional or linear pressure on the cannula, in order to avoid the risk of causing accidental disconnection or obstruction;
- If the guide wire is damaged, it may be difficult to continue the procedure. In these cases, only proceed if the damaged section can be made to progress through the trachea and if the length of undamaged guide wire is sufficient for continuing the procedure. If the guide wire is not fit for use, a new pack or new guide wire must be used;
- The introducer must only be used to insert the tracheostomy tube, do not attempt to use it to dilate the stoma;
- Do not allow the tube and/or the cuff of the tracheostomy tube to come into contact with sharp edges, in order to avoid damaging the cannula cuff or tube;
- Do not attempt to reposition the tracheostomy tube once it is in place, if the cuff has been inflated;
- Clean the area around the tracheostomy tube regularly.

Method of use

Preparing the components

Carry out the following steps before using the TRACHEO SET device:

- Test the balloon cuff of the tracheostomy tube and the inflation system by introducing air using the syringe in the pack;
- Having checked the integrity of the balloon and the system, deflate the cuff completely to avoid the risk of failure during the introduction of the tube;
- Check that the guide wire slides freely through the dilator, forceps and tube introducer.

Procedural technique

It is recommended to carry out the procedure under endoscopic control.

1. Place the patient in a supine position, using a support to hyperextend the head and neck if necessary;
2. Prepare the sterile field;
3. Apply suction to the pharynx, withdraw the tracheostomy tube (where present) and position the cuff above the vocal cords to avoid damaging them; deflate the cuff of the endotracheal tube and withdraw it to an appropriate distance above the introduction site, but still inside the trachea. Re-inflate the cuff once the endotracheal tube is properly positioned.
4. Introduce the bronchoscope (**Fig. 1**);
5. Locate and mark the landmarks for stroma formation → usually between the first and second tracheal rings;
6. At the previously chosen site, use a scalpel blade to make a vertical incision large enough to introduce the tracheostomy tube, i.e. about 1.5-2 cm long. The size of the incision can be calculated using the equation $OD \cdot \pi / 2$;
7. Refill the syringe with a few cc of sterile saline and connect to the cannula needle;
8. Insert the cannula needle at the chosen point, perpendicular to the skin, until it reaches the trachea. When the cannula needle reaches the trachea air bubbles will appear in the solution in the syringe (**Fig. 2**);
9. Check that the needle is properly positioned (it must be perpendicular to the axis of the trachea) by means of continuous endoscope monitoring (**Fig. 2**);
10. Remove the needle and syringe, leaving the tube in place;
11. Check that it is positioned properly again by connecting the syringe to the tube and making sure air bubbles appear in the solution;
12. Connect the introducer of the guide wire dispenser to the tube and slide the guide wire into the trachea, in a caudal position, by about 10-15 cm, leaving the rest of the wire outside (**Fig. 3**);
13. Check that the guide wire slides freely inside the tube, leaving the guide wire in place;

CAUTION! It is important that the Seldinger guide wire is not damaged during the procedure. To avoid this, always hold the free tip of the guide wire between the index finger and thumb while introducing the components.

14. Remove the tube from the guide wire;

15. Introduce the dilator on to the guide wire and turn it gently to advance it through the soft tissues until the resistance of the anterior trachea wall is felt. By applying constant pressure, introduce the dilator into the trachea until the first 2 cm of it are visible on the endoscope. This confirms that the stoma has been properly dilated. Check that the guide wire slides easily (**Fig. 4**);

CAUTION! If the dilator comes into contact with a tracheal ring, do not force the introduction manoeuvre. Repeat the procedure by changing the cannula needle introduction site.

16. Leave the dilator in place for at least 30 seconds in order to hold the stoma open and facilitate the introduction of the forceps;

If there are no Tracheo Forceps® forceps inside the package, go to point 25, otherwise continue to point 17 below.

17. Adjust the maximum opening of the forceps by turning the pin between the two arms (**Fig. 5a**);

18. Remove the dilator from the guide wire;

19. Insert the free guide wire into the hole on the front of the right arm of the forceps (**Fig. 5 and Fig. 5b**);

20. Holding the free part of the guide wire, advance the forceps towards the anterior wall of the trachea, in a caudal direction. During this step, if the previously created stoma has remained open, no resistance should be encountered, and the tip of the forceps will enter the trachea with ease;

Otherwise, apply constant pressure on the trachea to overcome its resistance. **It is recommended for these manoeuvres to be carried out under endoscopic control.**

CAUTION! If insertion is nevertheless impossible, reinsert the dilator and repeat the procedure from point 15 of the instructions.

21. With the arms of the forceps closed (rest position), lift the forceps keeping them in line with the trachea, so that the tip of the forceps penetrates deeper (**Fig. 6**);

22. Check that the guide wire slides freely through the forceps at all times;

23. Proceed with the formation of the stoma by carrying out opening/closure movements and advancing the forceps. Repeat the sequence, leaving the guide wire in place, until the stoma opening is large enough to allow the placement of the tracheostomy tube (**Fig. 7**);

24. At the end of this step, check the integrity and proper positioning of the guide wire, which must slide freely through the arm of the forceps, then withdraw the forceps;

25. Place the introducer inside the tracheostomy tube (the holes on the introducer should be perpendicular to the axis of the tube). Insert the guide wire into the introducer and start positioning the tube in the trachea, (**Fig. 8**). A tube other than that provided may be used, provided it has the same inner diameter. Both armoured and standard extra-long tubes can be inserted on the introducer;

26. Remove the introducer and guide wire, leaving the tracheostomy tube in place (**Fig. 9**). Apply suction to remove any secretions in the trachea, in order to prevent any blood not removed by suction from clotting and obstructing the airways;

27. Inflate the cuff of the tracheostomy tube by connecting the syringe to the Luer connector of the indicator balloon to ensure a perfect seal;

28. Connect the artificial ventilation circuit to the tracheostomy tube;

29. Secure the tracheostomy tube with the tube holder provided;

30. Remove the endoscope and the tracheal tube from the patient's mouth (**Fig. 10**).

Other useful information

It is possible to watch a demonstration of the percutaneous dilational tracheostomy procedure using TRACHEO SET at:

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