OPTIVISTA



A unique combination of optical and digital enhancements, for improved in vivo diagnosis





A leap forward for in vivo histology with i-scan and i-scan OE, a unique combination of digital and optical enhancements



Sharp High-Definition image, combined with **i-scan imaging** and **excellent illumination** for a more detailed view of the mucosa and enhanced detection.

OPTIVISTA - Enlighten your perspectives


The new premium solution for advanced endoscopy

The success of a procedure can hinge on your ability to clearly visualize the epithelial surface pitpattern and vascular pattern. In addition to the already established i-scan modes, the OPTIVISTA now features i-scan OE (Optical Enhancement), creating a very unique platform where both digital and optical enhancements are available. This unique combination provides extra information for a more accurate *in vivo* diagnosis through improved vessel and mucosal pattern characterization.

The new PENTAX Medical OPTIVISTA EPK-i7010 video processor, together with the latest i10 HD+ endoscopes series, or with the MagniView, optical zoom endoscopes, provides our best image ever. Whether you are performing a routine procedure or a complex intervention, the EPK-i7010 can help you to achieve the best clinical outcomes.







ilvia Sanduleanu, Maastricht LIMC+, The Neth



Improved in vivo diagnosis

A unique combination of **digital enhancement** (i-scan) and **optical enhancement** (i-scan OE) provides extra information for a more accurate in vivo diagnosis through improved vessel and mucosal pattern characterization.



Powerful educational platform





The integrated audio and HD+ video recorder enables to collect and share the latest




i-scan for enhanced detection and characterization

i-scan technology is a digital, image enhanced endoscopy (IEE) technology by PENTAX Medical. As a virtual chromo-endoscopy tool, i-scan provides an enhanced view of the mucosal structures and vascular patterns, supporting early detection, demarcation and characterization.

i-scan 1 (SE) for detection

By enhancing the mucosal structure in a natural color tone, SE (Surface Enhancement) supports detection, especially of flat lesions by highlighting abnormalities.



With SE, the difference in luminance intensity between the pixels concerned and the surrounding pixels is analyzed and the edge components are enhanced. SE adjustment of the noise erasure function allows more evident enhancement of the edges, which corresponds to minor changes in structure.¹

i-scan 2 (TE) for pattern characterization

TE (Tone Enhancement) enhances the changes of vascular and mucosal structure with a color tone change, supporting pattern characterization.





R component



The RGB components of an ordinary endoscope image are disintegrated into each component (R, G and B), and each component thus isolated is converted independently along the tone curve, followed by a reconstruction of the three components in one single image.

¹ Kodashima et al: Novel image-enhanced endoscopy with i-scan technology. World J Gastroenterol 2010








i-scan OE for vessel characterization


PENTAX Medical has recently developed the i-scan OE (Optical Enhancement), an optical filter that produces bandwidth-limiting light. When combined with image enhancement processing technology it clearly displays the surface structures of blood vessels, glandular ducts and mucosal membranes in higher contrasts than white light.

i-scan OE technology



i-scan 3 (OE) for vessel characterization

i-scan OE displays the surface structures of the blood vessels, glandular ducts and mucosal membranes in a color tone and supporting vessel's characterization.



² K Mönkemüller et al.: Present and future perspectives of virtual chromoendoscopy with i-scan and optical enhancement technology. Digestive Endoscopy 2013; 26.



Improved *in vivo* diagnosis A unique combination of digital and optical enhancements



The main wavelengths of spectral transmission correspond to the peaks of the hemoglobin absorption spectrum and the light between these peaks is continuously slightly emitted by raising baseline transmittance.²

OPTIVISTA, a powerful educational platform

The OPTIVISTA EPK-i7010 video processor is a state-of-the-art educational tool with customizable functions operated easily by the touch screen, or via the endoscope buttons.




Teach and assess with the unique TwinMode

- TwinMode is useful in teaching the appropriate interpretation of image enhanced endoscopy, providing simultaneous comparison of side by side endoscopic images.
- Among experts, TwinMode is respected appreciated as an educational tool for

non-experts in "building the bridge" between HD+ white light images and the different i-scan modes and findings.

• The simultaneous comparison of enhanced clinical images is useful in teaching the appropriate characterization of lesions.




Colon adenoma



Collect and share with the unique integrated video and audio recording

- The video recording function enables capture of HD+ video files through a USB storage device, for fast and easy sharing of findings with peers. Audio recording for video is captured through external microphone.
- Contributes to cost savings in the endoscopy room by avoiding unnecessary external HD recording devices.
- For best image collection, the OPTIVISTA also incorporates freeze scan technology which automatically selects the sharpest picture.



Exchange and expand your knowledge of i-scan by visiting the i-scan website www.i-scanimaging.com



HD processors



A complete diagnostic tool, bringing together the best of digital and optical enhancement for complimentary use in pattern and vessel visualization



A High-Definition video processor offering standard functionality supported by state-of-the-art imaging technology.



The video processor to enter the High-Definition and i-scan world with an excellent price performance efficiency.

Endoscopes





Slim caliber, bigger working channel and Close Focus, all combined in an i10 HD+ endoscope, the high performance endoscopes series, for the day-to-day advanced endoscopy practice.



With the 90i endoscope series, you can experience HD+ images combined with exceptional illumination, during routine examinations.



90K endoscope series are the optimal choice for an exceptional high resolution image combined with outstanding illumination.

EC REP EMEA Headquarter Germany

PENTAX Europe GmbH Julius-Vosseler-Straße 104 22527 Hamburg Tel.: +49 40 / 5 61 92 - 0 Fax: +49 40 / 5 60 42 13 E-mail: info.emea@pentaxmedical.com

United Kingdom

PENTAX U.K. Limited PENTAX House Heron Drive, Langley Slough SL3 8PN Tel.: +44 17 53 / 79 27 33 Fax: +44 17 53 / 79 27 94 E-mail: medical.uk@pentaxmedical.com

France

PENTAX France Life Care S.A.S. 116 quai de Bezons B.P. 204 95106 ARGENTEUIL CEDEX Tel: +33 1 / 3025 7575 Fax: +33 1 / 3025 7445 E-mail: contact.medical@pentax.fr

Netherlands

PENTAX Nederland B.V. Edisonring 6 6669NB Dodewaard Tel.: +31 88 / 5 30 30 30 Fax: +31 88 / 5 30 30 40 E-mail: medical@pentax.nl

Italy

PENTAX Italia S.r.I. Via Dione Cassio, 15 20138 Milano Tel.: +39 / 02 50 99 58 1 Fax: +39 / 02 50 99 58 60 E-mail: marketing.lifecare@pentaxitalia.it

Spain

SIMMEDICA – Sistemas Integrales de Medicina, S.A. Avenida del Sistema Solar 25 28830 San Fernando de Henares · Madrid Tel.: +34 91 / 301 62 40 Fax: +34 91 / 751 31 15 E-mail: sim@simmedica.com

Russia

PENTAX Europe GmbH Representative office in Moscow Sadovnicheskaya str, 82, build 2, entrance 6 Regus, Business center "Aurora" Office 2012, 2013 Moscow, 115035 Tel.: +7 495 / 792 52 00 Fax: +7 495 / 792 35 66

Turkey

PENTAX Turkey Veko Giz Plaza, Meydan Sokak No:3/43 343396 Maslak – Istanbul Tel.: +90 212 / 705 05 26 Fax: +90 212 / 705 05 00



6-10-1 Nishi-shinjuku Shinjuku-ku Tokyo 160-00231

