Generic Name : Total Protein Kit (Biuret Method)

Abbreviated name : TP

Order Information

Cat. No.	Package size
TP0102	R 4×40 mL
TP1102	R 1×17 mL
TP0103	R 6×40 mL
TP0104	R 6×60 mL
TP0105	R 4×250 mL

Intended use

In vitro test for the quantitative determination of TP concentration in serum and plasma on photometric systems.

Summary¹

Serum total protein, including albumin, is take charge of transporting substances, including macromolecules, and maintaining the plasma osmotic pressure. Hypoproteinemia can be caused by anti-body deficiency syndrome, liver cirrhosis, impaired kidney function, diarrhea and nutritional loss. Hyperproteinemia is seldom discovered. Only serious chronic inflammation or self-immunity disease may lead to hyperproteinemia.

Method

Biuret method

Reaction Principle

$\text{Cu}^{2+} + \text{Protein} \xrightleftharpoons{\text{pH} > 12} \text{Copper-protein complex (blue-violet colour)}$

At an alkaline solution (pH > 12) copper ions combine with protein to produce a blue-violet colour complex. The absorbency increase is directly proportional to the concentration of protein.

Reagents

Components and concentrations

R:	Sodium-potassium tartrate	32 mmol/L
	Sodium hydroxide	200 mmol/L
	Potassium iodide	30 mmol/L
	Cupric sulfate	12 mmol/L

Warnings and precautions

1. For in vitro diagnostic use only.
2. Take the necessary precautions for the use of laboratory reagents.
3. Preservative contained. Do not swallow. Avoid contact with skin and mucous membranes.
4. Disposal of all waste material should be in accordance with local

guidelines.

5. Material safety data sheet is available for professional user on request.

Reagent Preparation

Single reagent is ready to use.

Storage and stability

Up to expiration date indicated on the label, when stored unopened at 2-8°C and protected from light.

Once opened, the reagents are stable for 28 days when refrigerated on the analyzer or refrigerator.

Contamination of the reagent must be avoided.

Do not freeze the reagent.

Reagent blank absorbency

The absorbance of reagent blank at 546 nm should be <0.3 A.

Materials required but not provided

1. Calibrator and controls as indicated below.

2. NaCl solution 9 g/L.

3. General laboratory equipments.

Specimen collection and preparation ²

1. Serum, heparin or EDTA plasma is suitable for samples. Whole blood, hemolysis is not recommended for use as a sample. Freshly drawn serum is the preferred specimen.

2. Use the suitable tubes or collection containers and follow the instruction of the manufacturer; avoid effect of the materials of the tubes or other collection containers.

3. Centrifuge samples containing precipitate before performing the assay.

4. Stability: 3 days at 2-8°C
6 months at -20°C

Assay procedure

	Blank	Sample
R	1000 µL	1000 µL
Dist. water	20 µL	—
Sample	—	20 µL
Mix thoroughly at 37 °C, and read the absorbance 10 min. later.		
$\Delta A = [\Delta A \text{ sample}] - [\Delta A \text{ blank}]$		

Application sheets for BS series analyzers are available in this document. Refer to the appropriate operator manual for the analyzer-specific assay instructions.

Calibration

1. It is recommended to use the Human multi-calibrator from Mindray and 9

g/L NaCl for two-point calibration. Traceability of the multi-calibrator can refer to the calibrator instructions for use of Mindray Company.

2. Calibration frequency:

After reagent lot changed.

As required following quality control procedures.

Quality control

At least two levels of control material should be analyzed with each batch of samples. In addition, these controls should be run with each new calibration, with each new reagent cartridge, and after specific maintenance or troubleshooting procedures as detailed in the appropriate system manual. We recommend using the Human Assayed Control made by Mindray to verify the performance of the measurement procedure; other suitable control material can be used in addition.

Each laboratory should establish its own internal quality control scheme and procedures for corrective action if controls do not recover within the acceptable tolerances.

Calculation

The analyzer calculates the TP concentration of each sample automatically after calibration.

Conversion factor: $\text{g/dL} \times 10^{-1} = \text{g/L}$

Or: $C_{\text{sample}} = (\Delta A_{\text{sample}} / \Delta A_{\text{calibration}}) \times C_{\text{calibration}}$

Reference Intervals¹

Each laboratory should establish its own reference intervals based upon its patient population. The reference intervals measured at 37°C listed below were taken from literature:

Sample Type		S.I. Units
Serum / Plasma	Adults	66-83 g/L
	Prematures	57-80 g/L
	Newborns	41-63 g/L

Performance Characteristics

Representative performance data obtained from Mindray system (Mindray BS series analyzers / Mindray TP Reagent) is given below. Results may vary if a different instrument, an individual laboratory or a manual procedure is used.

Limitations-interference

The following substances were tested for interference with this methodology. Criterion: Recovery within $\pm 10\%$ of initial value.

Substance	Level Tested	Observed Effect
Ascorbic acid	30 mg/dL	NSI*
Bilirubin	30 mg/dL	NSI

Hemoglobin	250 mg/dL	NSI
Lipemia	2000 mg/dL	NSI
Dextran	1000 mg/dL	NSI

* NSI: No Significant Interference (within $\pm 10\%$)

Linearity range

The Mindray System provides the following linearity range:

Sample Type	S.I. Units
Serum / Plasma	2 g/L-120 g/L

If the value of sample exceeds 120 g/L, the sample should be diluted with 9 g/L NaCl solution (e.g. 1+ 1) and rerun; the result should be multiplied by 2.

Analytic Sensitivity/Limit of Detection

The lowest measurable TP concentration that can be distinguished from zero is 2 g/L with 99.7% confidence.

Precision

Precision performance using the CLSI Approved Guideline EP5-A2 to assay serum control appears in the table below³. U: g/L

Type of Imprecision	Level II			Level III		
	Mean	SD	CV %	Mean	SD	CV %
Within-run	58.361	0.640	1.096	43.704	0.319	0.729
Between-run		0.267	0.457		0.312	0.715
Between-day		0.154	0.264		0.388	0.889
Within-device		0.710	1.217		0.592	1.354

Method Comparison

A comparison between Mindray System (Mindray BS series analyzers /Mindray TP Reagent) (y) and Hitachi/Roche System (Hitachi /Roche TP) (x) using 40 samples gave following correlation (g/L): $y=0.9468x+1.8543$, $R^2=0.9933$. Details of the comparison experiments are available on request.

References

- 1.Thomas L. Clinical Laboratory Diagnostics. 1st ed. Frankfurt: TH-Books Verlagsgesellschaft; 1998. p. 644-7.
- 2.Tietz NW, ed. Clinical Guide to Laboratory Tests, 3rd ed. Philadelphia, Pa: WB Saunders Company, 1995:518-522.
- 3.CLSI. Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline-Second Edition. CLSI document EP5-A2 [ISBN 1-56238-542-9. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087 USA, 2008.

Graphical symbols

In Vitro Diagnostic
medical device

Batch Code

European
Conformity

Authorized representative in
the European Community

Use By

Consult
Instructions for use

Temperature
Limit

Manufacturer

Catalogue
number

© 2024 Shenzhen Mindray Bio-Medical Electronics Co., Ltd. All rights Reserved

Manufacturer: Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

Address: Mindray building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057 P.R.China

E-mail Address: service@mindray.com

Website: www.mindray.com

Tel: +86-755-81888998

Fax: +86-755-26582680

EC-Representative: Shanghai International Holding Corp.
GmbH(Europe)

Address: Eiffestraße 80, Hamburg 20537, Germany

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Riferimenti

1. Thomas L. Clinical Laboratory Diagnostics. 1st ed. Frankfurt: TH-Books Verlagsgesellschaft; 1998. p. 644-7.
2. Tietz NW, ed. Clinical Guide to Laboratory Tests, 3rd ed. Philadelphia, Pa: WB Saunders Company, 1995:518-522.
3. CLSI. Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline-Second Edition. CLSI document EP5-A2 [ISBN 1-56238-542-9. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087 USA, 2008.

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Indirizzo: Mindray building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057 P.R.China

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© 2024 Shenzhen Mindray Bio-Medical Electronics Co., Ltd. Tüm hakları saklıdır
Üretici:Shenzhen Mindray Bio-Medical Electronics Co., Ltd.
AT Temsilcisi:Shanghai International Holding Corp.GmbH(Avrupa)