
Volume Viewer

GE HealthCare's premium 3D advanced visualization and image processing platform brings you powerful capabilities to help you visualize and interpret your multimodality imaging data with confidence and ease.

Imaging modalities provide a wealth of diagnostic information, but also present a challenge to your radiology workflow as image volumes continue to increase while your time for reading and reporting decreases.

Volume Viewer provides you with a rich 3D image processing toolset aimed at creating and displaying the views you need with little user input and streamlining interpretation and reporting by providing the visualization tools you need with minimum clicks.

With an intuitive, modern user interface, Volume Viewer helps a wide diversity of users to learn and master the expanding portfolio of sophisticated tools and applications that it provides. This customizable user interface maximizes the real estate allocated for displaying images and provides access to the user's favorite tools directly from the image viewport. In addition, a customizable toolbar provides one click access to commonly used tools to facilitate a productive review.

CONFIDENTIAL

SECRET DE AFACERI

Overview

Volume Viewer provides excellent 3D visualization and processing capabilities for reading and comparing CT, MR, 3D X-ray, PET, PET/MR and PET/CT datasets. Volume Viewer also features a broad portfolio of high-performance analysis tools, automating routine tasks and helping to make 3D image processing a stress-free component of your routine workflow.


Volume Viewer is available on VolumeShare 7, a multi-modality advanced visualization workflow solution that helps to enhance diagnostic precision and productivity.

Volume Viewer options*

- Volume Illumination
- Q.Freeze 2
- Autolaunch / Pre-processing

Compatible products**

- 3D Suite
- Integrated Registration
- Bone VCAR
- GSI Viewer

Features

- ✓ Unified **user interface**, with one click access to your favorite tools and more space for displaying images
- ✓ Customizable **protocol-driven workflows** designed to help enhance your clinical review productivity for a wide range of care areas, with 3D and 4D real-time interaction
- ✓ "Smart layout" protocol, which adapts the screen layout to the type of loaded exams
- ✓ High resolution, real-time rendering modes: MIP/Min IP/Average; standard volume rendering and enhanced volume rendering with Volume Illumination* on the fly; Navigator view for endoluminal fly-through; lumen and curved views
- ✓ Compare mode for reviewing multiple exams and modalities, simplified with **Dynamic Load**, to load and register new series into the current review session, with **Integrated Registration****
- ✓ Rich set of **2D/3D ROI tools**, including **AutoContour** for 3D semi-automatic contouring of structures of interest (CT, MRI, PET), and **ROI color coded** for tissue classification based on voxel values
- ✓ Advanced tools, to take advantage of full 3D capabilities: **AutoSelect**, for easy point-and-click segmentation; One or two clicks **Quick Vessel Trace** to analyze all vessels, in curved reformat, lumen, or MPR view; Advanced 3D processing tools
- ✓ **Summary Table** extends reporting capabilities, while collecting measurements as they are deposited on the images

CONFIDENTIAL
SECRET DE AFACERI

Feature details

Unified user interface

- The user interface provides one-click access to favorite tools and more space for clinical images
- The page layout is organized to facilitate navigation through the review steps. Guided protocols walk the unfamiliar user through a clear workflow of creating and manipulating volumetric data, helping to reduce the learning curve by providing a consistent workflow.
- The customizable toolbar can be docked at the top, bottom, left, or right, by a simple drag and drop. This preferred location is saved for future use for each user.
- An extensive set of advanced tools are grouped into categories on the toolbar to facilitate quick retrieval when needed
- Each category and tool are easily configurable by a simple drag and drop in the Toolbar so that they are only one click away for your future reviews
- To reduce long mouse trips, up to six favorite tools can be accessible from each image via a right mouse menu
- Customization of the toolbar and tool behaviors with a user's preferences is provided in a single intuitive Preferences menu
- Keyboard shortcuts and quick access tools for on-viewport controls

The following mouse modes are available to help you manipulate the images:

- Interactive VR adjustment lets you adjust VR opacity interactively with the mouse as Window Level/Window Width
- Direct Paging (continuous or contiguous)
- Free-hand rotation on 3D and oblique views
- Percentage fusion, to easily change the transparency of objects on fused images. There are several fusion modes available for your use.
- Mag glass displays zoomed-up imagery around the cursor

Protocol management and loading

- Customizable protocol-driven workflows for a wide range of clinical applications, including support for multiple review steps. You are guided through the review steps by a Navigation bar at the top of the User Interface.
- Review protocols created by user with Volume Viewer can be saved as favorite protocols and accessed easily in routine
- Progressive Load enhances image loading performance. This feature has two modes:
 - Images are loaded in sequential mode (e.g., Reformat)
 - Images are loaded in interlaced mode (e.g., 3D/VR)
- In either case, image review can begin as soon as the first image is displayed
- Navigation through series and exams can be performed without exiting the patient list

Smart layout

Volume Viewer General Review protocol presents an adaptive layout for single and multi-modality PACS-like reading of data. This protocol distributes and displays all loaded series of a single exam in an intelligent layout to maximize use of available monitor space. With multiple exams, the protocol automatically detects and distributes "like" series for optimal comparison.

Advanced visualization capabilities

Volume Viewer offers various capabilities to display advanced rendering modes:

- 3D Volume Rendering / Volume Illumination* images
- MPR views with different slice thickness. The following rendering modes are available for all Thick Slab: MIP, Min IP, Volume Rendering, Average.
- Navigator views which display interactive endoluminal views exportable as movies
- Curvilinear reformatting allowing display of curved, lumen and cross-section views for various structures (vessels, spine, etc.)
- Merge multiple volume rendered models into a single view

Multiple volume rendering and volume illumination* objects

- Volume Viewer allows you to merge up to eight Volume Rendered models from the same series into a single 3D view
- You can adjust independently the threshold, colors, and transparency of each Volume Rendering or Volume Illumination*

Predefined cut planes

- Volume Viewer lets the user define cut planes to isolate specific structures in the VR model, in 16 configurations, including: Left / Right / Inferior / Superior / Anterior / Posterior / Front
- You can then display the voxels values on each plane

Lumen view

Lumen View provides an unfolded 3D view around a user-defined centerline. The lumen view can be interactively adjusted (rotation around the centerline, width, field of view).

Auto-center fly through with smart cursor

- In fly-through studies (airways, colon, angiography) navigation along the centerline of the structure of interest can be recorded step by step and can be exported as a movie
- In addition to the regular 180° or less Navigator view, Volume Viewer enables a Fish-Eye View with any wide camera angle value from 180° to 360°. This provides a view of structures both in front of and behind the user's virtual location on the same image.

3D and 4D capabilities

- Volume Viewer allows 3D cursor synchronization in any orientation on the fly
- Real-time reformat in oblique planes is available for your exams, as well as simple and double oblique interactive modes
- Quick access to cross reference display gives you the ability to display the current slice location or all slices, or other planes
- Real-time interaction in 4D mode, by using Cine tool on all planes and 3D views for CT, MR, and PET multiphase data

Compare mode

- With Volume Viewer several volumes can be loaded either from the same exam or from different exams. This is useful for exam comparison (Compare Mode) or for Multiphase examinations (of the liver for example). The Compare Mode is also available for PET/CT examinations and takes advantage of a dual screen configuration.
- Multi-phase CT and MR data are listed as separate phases in the Series Selection screen, allowing selection of all phases or a subset of phase data for review. The Series Selection Panel displays the multi-phase data split into arterial and portal phases, for example.
- Dynamic Load, compatible with both Volume Viewer and Fusion protocols, lets you drag & drop 3D volumes from CT, MR, PET, and 3D XA modalities into a desired viewport. With Integrated Registration**, Dynamic Load allows new volumes to be registered and loaded on the fly. Save State series can be restored as a separate session. For dual monitor configurations, a Save State series can be displayed on the right monitor together with a current session on the left monitor. For single monitor configurations, a user may simply switch between the two sessions.
- The zoom and pan functions are propagated to all images displayed in the same orientation

2D/3D ROI for quantitative measurements

- Volume Viewer contains a set of basic 2D/3D tools: distance, angle, report cursor, arrow tool annotation, free hand ROI tool, that can be placed on fused images as well
- Easy deposition and labelling of measurements
- All 3D ROIs can be customized, and color coded to display statistics computed on different intensity range and help assessing tissue classification based on voxel values
- All measurements are considered bookmarks and can be accessed by navigating through the list of measurements

AutoContour tool

- AutoContour provides consistent tools for 3D contouring ROI, with the same workflow used for CT, MR, and PET images
- For CT and MR images, the contouring tool will suggest contour boundaries after defining one point or a diameter on a structure of interest, which can be edited prior to validation
- For PET images, the contouring tool has been adapted for standard SUV calculations

Advanced 3D tools

Segmentation tools

- AutoSelect segmentation tool allows you to add or remove any continuous structures of interest with a simple point and click
- You can also access advanced segmentation tools such as Scalpel, Threshold, Dilate, Erode, Subtraction, Intersection, Addition, Filter Floaters

Quick Vessel Trace (CT and MR images)

- You can perform vessel analysis with one or two points clicks from any CT or MR Angiography images
- Quick Vessel Trace automatically extracts the vessel after user clicks and launches vessel tracking for fast review in curved reformat, cross-sectional, best L-section, lumen, and MRP view
- Tracking can be extended proximally or distally for a full view of the vessel

Summary table

- The Summary Table collects measurements and organizes them according to finding and exam date
- The Summary Table allows interactive navigation through findings in the image viewports and can be docked or expanded
- The Summary Table can be copied into the client clipboard and then pasted in email, reports, office software, etc. (on AW Server only)

Export capabilities

Volume Viewer contains multiple standard options for exporting the results of a review session:

- Save images to the database, as new DICOM series
- Save the significant images as Key Image Note objects in the database. Key Image Notes and End Review allow you to flag images of interest as Key Image Notes (IHE profile) and push them to the archiving system when you exit the application.
- Save State object are used to save the work in a new series containing all the post processing data (3D Model, displays, measurements, annotations, etc), for future review
- Save STL, OBJ, VRML or 3MF file from 3D Volume Rendering / Volume Illumination* viewport for 3D printing purposes through 3D Suite**

Additional export tools are available to create specific sequences of images to be exported:

- The Batch tool allows creating a sequence of reformatted images or a sequence of rotating 3D views
- The Movie tool creates a comprehensive movie including different rotations, zooms, and pan of the image, which can be exported as DICOM series or mpeg file
- Quick Export: Exports in a single click a full batch of contiguous images at the displayed thickness for 2D images, or a batch of rotations of a 3D View
- Cardiac Review and Export: Processing and reviewing cardiac exams for CT, MR and PET with manual oblique reformatted protocols can be exported into a multi-phase Cine movie that allows the referring physician to review the exam in a dynamic mode
- Capture and send images from the viewport to your computer clipboard on the fly (Copy/Paste on AW Server only)

CONFIDENTIAL

SECRET DE AFACERI

Additional modality solutions

CT image quality

- CT filters to enhance or smooth CT images to help improve image quality on the fly
- The lung filter helps enhance contours of images reconstructed in standard mode for excellent visualization of lung structures
- A smoothing filter reduces noise while maintaining high image contrast

CT Dual energy protocols

Three protocols are available to enable review of dual energy images acquired on compatible GE HealthCare scanners and enhancing pixels with specific ratio: DE Calcium ($HU_{80}/HU_{140} \geq 1.25$), DE Uric Acid ($HU_{80}/HU_{140} \leq 1.25$) and DE Custom (custom threshold on HU_{Low}/HU_{High}).

GSI Viewer**

GSI Viewer provides dedicated capabilities for Gemstone Spectral Imaging (GSI) visualization. It is the primary method of viewing and manipulating spectral images. It allows for:

- Review of monochromatic energy images at user selectable energy levels
- Detailed analysis using material decomposed images (such as water-iodine, water calcium, liver fat, etc.) and complementary information using the Effective-Z images
- Creation of virtual unenhanced (VUE) images
- Gout protocol with automated workflow

Bone VCAR**

- Bone VCAR is available in Volume Viewer toolbar and can be launched on the fly
- Vertebrae are automatically recognized and labelled (deep learning algorithm)
- A dedicated tool provides a specific layout of the spine including curved views
- Combines spine labelling on the fly for all CT review and a dedicated layout to improve spine review

Review of MRI studies

Volume Viewer offers support for review of MR studies, through dedicated protocols.

- The MR-specific parameters (PSD name, image weighting, scan plane, etc.) are recognized at loading and are used to customize the layout display
- Anatomy-based protocols for standard review (e.g. spine, brain stroke...), with dedicated and guided review steps, to help review standard MR exams

- Dedicated MR Cardiac Viewer for single and multiple Cine Review enables comparison between cardiac series such as Time Course and Myocardial Delayed Enhancement (MDE) images
- Support of MR Multi Echo, Multiphase, Diffusion series
- Support of parametric series (e.g., ADC, MTT, and other parametric maps) created in READY View, allowing direct functional measurements
- 2D and 3D ROI propagation to other phases / series / exams allowing you to easily correlate information from multiple sources
- In combination with Integrated Registration**, enables direct access to MR image registration (inter/intra exam), image fusion and Whole Body MR Review protocols
- Ability to automatically bind several MR axial series corresponding to sequential axial locations into a single series. This is useful for consolidating multi-stage acquisitions for displaying Whole Body MR scanning.

SUV measurements on PET images

- Volume Viewer supports Standardized Uptake Values (SUV) for image display and measurements. Several SUV scales are available like SUVbw, SUV lbm, SUVbsa, as well as SUV Peak.
- Window/level presets may be user-defined for PET images based on SUV values
- The SUV values are available in all the basic 2D/3D ROI tools of Volume Viewer and in the AutoContour tool

PET image quality improvements

- Enhanced 3D visualization algorithm to fully support Q.Clear PET images
- Q.Freeze 2* provides a respiratory motion-corrected Gated PET volume, statistically similar to conventional static PET with significantly reduced or eliminated blurring effects due to patient respiration. It provides also a corrected Gated PET series where counts are summed back to each respiratory phase.

PET/CT and PET/MR dedicated protocols

- Volume Viewer contains a list of predefined protocols allowing the review of PET/CT and PET/MR data. The screen layouts contain fused views between the morphological and functional images, as well as 3DMIP rendering of PET images.
- Factory protocols are customizable for your own review
- PET 4D protocol allows users to load and display gated and dynamic PET series, as well as summing or reframing these series

CONFIDENTIAL
SECRET DE AFACERI

Platform features

Smart compression

- The Smart Compression technology automatically displays images with full fidelity when the image is still, then uses the selected compression level for increased interaction speed during user interaction. This allows for diagnostic reads on full fidelity static images with responsive dynamic display even at low bandwidth. On-image visual indicators notify the user when compression is in effect.

Pre-processing* features

- **Preprocessing*** is a Volume Viewer option that automatically performs routine processing tasks and saves the results, so they are ready for you
- When new exams are transferred to the database, Volume Viewer recognizes user-defined keywords in the Series Description and launches the appropriate pre-processing. The results are stored in a Save State object with the original study so you can load it when you are ready for review.
- **Preprocessing*** is compatible with the following purchasable medical devices as separate products: AutoBone™ Xpress, CardIQ Xpress Reveal, CardIQ Xpress Function, CardIQ Xpress Process, Advanced CTC Pro 3D EC, Colon VCAR (not available for sales in the US), and CT Perfusion 4D and FastStroke

AutoLaunch* features on AW workstation

- **AutoLaunch*** is a Volume Viewer feature that automatically preloads exams, so they are ready for your review
- When you are processing a study and a new exam is transferred to the AW Workstation, this feature automatically launches the new exam with an appropriate Volume Viewer protocol in the background. When you are ready, one click in the AutoLaunch* window instantly brings up the exam in the Volume Viewer application, eliminating waiting time and extra steps to load the exam into computer memory for processing.
- **AutoLaunch*** is compatible with CT, MR and PET single volume protocols of Volume Viewer
- When combined with optional applications (see "Pre-processing" section), this feature gives access to data already preprocessed, that you can review by just clicking on AutoLaunch* window

Summary of operation

Volumetric models are loaded by selecting the exam or series. The user can select a protocol category from an anatomical selector or go directly to a Review Layout. In either case, images are loaded progressively in the background; this gives control to the user in just a few seconds after selecting the images. Selecting a Review Layout launches a volumetric display protocol with predefined layout preferences. Review Layouts may be customized and combined to suit the user's workflow. Selecting a protocol category unlocks a variety of visual protocols that include the layout, threshold, rendering mode and filming formats. Some of these protocols direct the user through the process, providing capabilities to interactively view and manipulate the model, increasing productivity and consistency for all modalities.

Additional modality solutions

- **Dual energy images** supported within dedicated review protocols
- **MR review support:** Recognition of PSD name, and specific MR parameters at image loading, which simplifies protocol customization
- Dedicated review workflows for PET/CT and PET/MR images, including PET SUV measurements
- **Q Freeze 2*** combines the quantitative benefits of 4D-PET gated imaging into an image that provides both frozen patient motion and reduced image noise
- **Reformat XA** provides a dedicated review with MPR 3D and endoluminal navigation for 3D X-ray angiography (also known as Cone-Beam CT)

System requirements

- Refer to AW Applications platform compatibility document
- Recommended monitor resolution is up to dual 2MP (1600 x 1200) or a single 3MP (1536 x 2048)

Intended use/ Indications for use

Volume Viewer is a medical diagnostic software that allows the processing, review, analysis and communication of 3D reconstructed images and their relationship to originally acquired images from CT, MR, X-Ray Angiography and PET Scanning devices. The combination of acquired images, reconstructed images, annotations and measurements performed by the clinician are intended to provide to the referring physician clinically relevant information for diagnosis, and surgery and treatment planning.

CONFIDENTIAL

SECRET DE AFACERI

About GE HealthCare

GE HealthCare is a leading global medical technology, pharmaceutical diagnostics, and digital solutions innovator, dedicated to providing integrated solutions, services, and data analytics to make hospitals more efficient, clinicians more effective, therapies more precise, and patients healthier and happier. Serving patients and providers for more than 125 years, GE HealthCare is advancing personalized, connected, and compassionate care, while simplifying the patient's journey across the care pathway. Together our Imaging, Ultrasound, Patient Care Solutions, and Pharmaceutical Diagnostics businesses help improve patient care from diagnosis, to therapy, to monitoring. We are an \$19.6 billion business with approximately 51,000 colleagues working to create a world where healthcare has no limits.

Follow us on [LinkedIn](#), [X \(formerly Twitter\)](#), and [Insights](#) for the latest news, or visit our website [gehealthcare.com](https://www.gehealthcare.com) for more information.

CONFIDENTIAL
SECRET DE AFACERI

→ Scan the QR Code or visit us at
[gehealthcare.com/products/advanced-visualization/
all-applications/volume-viewer](https://www.gehealthcare.com/products/advanced-visualization/all-applications/volume-viewer)

Regulatory compliance:

This product complies with the Regulation (EU) 2017/45 of the European Parliament and the Council of medical device (MDR)

Rx Only

This product, its features, options, and compatible products may not be available in some other countries or regions. Please contact your sales associate.

References

* Volume Viewer option requiring license key.

** Separate medical device requiring specific license key – refer to dedicated product datasheet.

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→ Scanați codul QR sau vizitați-ne la
[gehealthcare.com/products/advanced-visualization/
all-applications/volume-viewer](https://www.gehealthcare.com/products/advanced-visualization/all-applications/volume-viewer)

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
GE Healthcare

CT Perfusion 4D

Comprehensive perfusion analysis available for
neurology, oncology and cardiology.

Easy-to-use, protocol-driven workflow for all organs.

Perfusion 4D has been expanded beyond stroke and oncology to now include dynamic myocardial perfusion. For any organ where blood perfusion is of interest Perfusion 4D is a complete package to quantitatively provide the necessary information to guide treatment decisions.

The protocol-driven design leads the user step-by-step through the process, reducing keystrokes and improving repeatability. Get the information you want quickly and reliably.

CONFIDENTIAL
SECRET DE AFACERI


Overview

GE's Perfusion 4D is a fast, easy-to-use automated software for analyzing CT Perfusion images related to stroke, tumor angiogenesis and dynamic myocardial perfusion. Its simple user interface and automated perfusion post-processing make it easy for you to diagnose quickly and accurately – and make treatment decisions more confidently.

GE leverages its deconvolution expertise, with its innovative delay correction algorithm. Perfusion 4D takes perfusion assessment to the next level, with the addition of CT Dynamic Myocardial Perfusion. Providing the classic color overlays with quantitative measurements to evaluate myocardial blood perfusion.

Highlights

- Improved IQ of functional maps in presence of noise
- Improved processing times
- Improved ROI mirroring for brain protocols
- Deep Learning algorithm for ventricle removal¹
- Automated midline detection for brain symmetry
- Absolute and relative values with customizable thresholds for tissue classification¹
- Permits injection rates of 4cc/sec
- Dynamic registration
- Dynamic myocardial perfusion maps

Visit us:

<https://www.gehealthcare.com/en/products/advanced-visualization/all-applications/ct-perfusion-4d-neuro>

Visit us:

<https://www.gehealthcare.com/en/products/advanced-visualization/all-applications/ct-perfusion-4d-multi>

Features

- Visualize all information in true volumetric form with the ability to employ all the volumetric-based image analysis tools in the AW Volume Viewer.
- Deep Learning algorithm for brain ventricles removal¹
- Automated midline detection for brain symmetry
- Absolute and relative values with customizable thresholds for tissue classification¹
- Background processing will send processed images and functional maps in email² format to stroke team²
- Whole heart perfusion with semi-automatic guided workflow to aid in assessment of myocardial ischemia
- Smart Maps using the intelligent 4D noise suppression algorithm improves functional map image quality in the presence of noise.
- Dynamic non-rigid registration protocols for dynamic acquisitions within the body.
- Productivity enhancements include:
 - **Ability to select automated or manual vessel detection for Brain Perfusion**
 - Faster processing times for functional map computations
 - Volumetric visualization of functional maps
 - Access from anywhere using AW Server.
- Time stamps on functional maps and Regions of interest (ROIs).
- Interactive large vessels exclusion.
- Many basic features, including the ability to:
 - Process double phase acquisition data.
 - Create 4D Regions of Interest (ROIs).
 - Display ROIs statistical information.
 - Display averaged time intensity information of the ROIs.
 - Simultaneously review any number of functional maps.
 - Customize protocols
 - Save parameters
- Includes CT perfusion protocols for stroke, brain & body tumor, myocardium as well as dynamic registration
- **Functional maps:**
 - Regional Blood Volume (BV)
 - **Regional Blood Flow (BF)**
 - Regional Mean Transit Time (MTT)
 - Contrast arrival delay (IRF T0)
 - Transit time to peak of impulse residue function (TMax)
 - Mean Slope of Increase
 - Capillary permeability surface area for lengthened acquisition protocols (PS)

- Average image
- Base image
- Tissue Classification

System Requirements³

- AW Workstations AW 4.7 Ext. 14 or higher
- AW Server 3.2 Ext. 3.2 or higher

Regulatory Compliance

This product complies with Regulation (EU) 2017/745 of the European Parliament and of the Council on medical devices (MDR). This product or its features may not be available in some other countries or regions. Please contact your sales associate.

Rx Only

Indications for Use

CT Perfusion 4D is an image analysis software package that allows the user to produce dynamic image data and to generate information with regard to changes in image intensity over time. **It supports the analysis of CT Perfusion Images (in the head and body)** after the intravenous injection of contrast, in calculation of the various perfusion-related parameters (i.e. regional blood flow, regional blood volume, mean transit time and capillary permeability). The results are displayed in a user-friendly graphic format as parametric images.

This software will aid in the assessment of the extent and type of perfusion, blood volume and capillary permeability changes, which may be related to stroke or tumor angiogenesis and the treatment thereof.

1. Requires CT Perfusion 4D *Neuro Enhanced* license
CT Perfusion 4D *Neuro Enhanced* may not be available in all countries or regions
Please contact your sales associate.
2. Requires *Send by Email* license and AW Server version with email service, not compatible with stand-alone AW systems
Send by Email may not be available in all countries or regions
Please contact your sales associate.
3. Send by email functionality requires an AW Server version with email service available

CONFIDENTIAL

SECRET DE AFACERI

ation at work

© 2021 General Electric Company.
All rights reserved. Data subject to change.
GE and GE Monogram are trademarks of General Electric Company.
* Trademark of General Electric Company.
DOC1996217

Features Detail

Basic Features

- Process double phase acquisition data (data with minimal temporal sampling and data with longer temporal sampling in one acquisition to cover entire functional perfusion mechanism).
- Create 4D ROIs (multi-locations ROI with time-varying position).
- Display ROI statistical information (Area, Min., Max, Std Dev.).
- Display input ROI from several locations on the same graph.
- Display the averaged time-intensity information from an ROI (VOI).
- Simultaneously review any number of functional maps.
- Customize any viewport to series, graph or functional view.
- View series or functional data in axial, sagittal, coronal or oblique view.
- View series or functional data in 3D or VR mode.
- Customize protocols (advanced settings, color ramp and/or W/L used for each generated functional map) and to save these user preferences.
- Save the state of a protocol including parameters, ROIs, functional maps and pre/post-processing results for further reloading in CT Perfusion.
- Save all images in the database.
- Save ROIs statistical information using screensave.
- Print to film or color paper printer.
- Save parameters, used to generate functional maps, as DICOM objects and to restore them including ROIs with the corresponding maps in CT Perfusion.

Protocols

CT Perfusion 4D has fully integrated protocols that assist the user in obtaining results efficiently for neurology, oncology and cardiology care areas.

CT Perfusion 4D package protocols are:

CT Brain Stroke Protocol. The CT brain stroke protocol will open a loading panel containing information about the selected perfusion series in a tabulated format. It is possible to remove unwanted time points from the loading sequence.


Once the exam is loaded the SW will perform the following algorithmic steps:

- **Image Registration** to minimize the effects of patient movement.
- **Bone Removal** to remove the cranium automatically from the dataset.
- **Arterial Input** used as the input for the deconvolution algorithm. Artery is automatically selected and can be easily adjusted if needed.
- **Venous output** for partial volume correction. Vein is automatically selected and can be easily adjusted if needed.
- **Last pre-enhancement / First post-enhancement** image are automatically selected to ascertain correct input for functional maps calculations.
- **Final Settings and Compute.** The software will compute the functional maps with the saved settings, but can be made to stop and recompute the functional maps if changes to the algorithm settings were implemented. The changes may include any of the following:
 - Algorithm Resolution, Smart Smoothing size.
 - Maximum Blood Flow, Brain Density, Hematocrit Ratio.

The software will use deconvolution to generate all functional maps, *Blood Flow, Blood Volume, Mean Transit Time, and Transit Time to IRF Peak (Tmax)*.

The Brain Stroke protocol has multiple review steps to guide the user through the functional information.


- **Review settings:** When loading the dynamic exam to process, this first review step is initiated. The general purpose of this step is to confirm the automatic selection of the input vessels for the brain stroke or tumor protocols, or to conduct a general review of the time resolved perfusion data and place ROI's to select vessels.


GE imagination at work

CONFIDENTIAL
SECRET DE AFACERI


- **Set Symmetry:** The software will automatically define the symmetry plane to be used for mirroring ROI's and relative thresholds. For the brain stroke protocol, the application provides a "Check view" image for a fast, intuitive way to indicate the plane and reference hemisphere identified.


If edition is needed a dedicated review step allows users to quickly adjust the symmetry axis.

- **Review Functional Maps** enables the viewing and assessment of functional maps such as Blood Flow, Blood Volume, Mean Transit Time, and Tmax. ROI mirroring functionality makes it easier to compare hemispheres.


- **Tissue Classification**, specific to the brain stroke protocol, enables the visualization of regions that are segmented from absolute or relative values, customizable thresholds and user selectable input maps.

The Brain Stroke protocol is typically used for the assessment of cerebral hemodynamic disturbances e.g. stroke and traumatic brain injury.

- **Email**, when working on AW Server and invoking the preprocessing for full automation of the Brain Stroke protocol, the processed images can be configured to auto send to the stroke team in email format.

CT Brain Tumor Protocol. Essentially identical to the brain stroke protocol but without the Tissue Classification step.

GE imagination at work

CONFIDENTIAL

SECRET DE AFACERI

CT Body Tumor Protocol. The Body tumor protocol produces functional maps automatically once the user selects an artery using a simple and intuitive workflow.

The Body Tumor protocol also leverages most of the advanced algorithm settings that are part of the Brain Stroke protocol. The protocol has the following algorithms

- **Noise Threshold.** Thresholds can be customized to exclude air and bone from the calculations.
- **Artery Input.** System allows for an intuitive selection of the artery.
- **Last pre-enhancement / Last post-enhancement image / Last Second Phase Image** (last image of the phase 2 data). These can be selected to ascertain correct input into the parameter calculations.
- **Final Settings and Compute.** The user can compute immediately or customize the following advanced settings:
 - Algorithm Resolution, Smart Smoothing.
 - Maximum Blood Flow, Brain Density, Hematocrit Ratio.

The Compute command then engages the exclusive de-convolution algorithm to produce the parametric images.

CT Pancreas Perfusion Protocol. Identical to the Body Tumor Protocol but optimized for pancreas tumors.

CT Prostate Perfusion Protocol. Identical to the Body Tumor Protocol but optimized for prostate tumors.

CT Kidney Perfusion Protocol. Identical to the Body Tumor Protocol but optimized for kidney tumors.

CT Soft Tissue Perfusion Protocol. Identical to the Body Tumor Protocol, with an additional step to define the vein output for partial volume correction (automatic or manual venous selection).


CT Spleen Perfusion Protocol. Identical to the Body Tumor Protocol but optimized for spleen tumors.

CT Bone Perfusion Protocol. Identical to the Body Tumor Protocol but optimized for bone tumors.


CT Non-Rigid Registration. Fully automated protocols used for dynamic body acquisitions to register data for use in perfusion analysis. Provides simultaneous review of original data with the registered results for full confidence in the newly created series. Saved registered results are to be loaded into the appropriate body or cardiac perfusion protocol for functional analysis

Two registration protocols are provided. CT Cardiac Dynamic Registration and CT Liver Dynamic Registration.

CT Liver Perfusion Protocol. Identical to the Body Tumor Protocol, with an additional step to define the portal vein input (automatic or manual venous selection).

Processing times have improved exponentially, now with an average processing time of one to three minutes depending on slice thickness. Hepatic Arterial Fraction maps have been improved and can be an efficient tool for tumor detection. A new map, hepatic arterial blood flow (HABF), has been added to assist in tumor assessment.

GE imagination at work


CONFIDENTIAL

SECRET DE AFACERI


Dynamic CT Myocardial Perfusion Protocol

Complete guided workflow for assessing dynamic myocardial perfusion. The dynamic myocardial perfusion protocol is laid out in a series of three easy to follow steps. Start by selecting the series that has been generated from the non-rigid registration protocol and then follow the guided workflow:

- 1) **Cardiac Reformat** – images are automatically orientated into short and long axis views,
- 2) **Segmentation** –define the valve plane and the apex of the heart. The software will then calculate the endocardial and epicardial borders and provide a 17 segmented map of the myocardium along with a graphed view for each segment

- 3) **Functional Maps** – the third and final step produces the quantitative data with the color overlay and bullseye maps for blood flow, blood volume, mean slope of increase and mean transit time.

Quantitative data can be exported as a .csv file for a complete segmental and regional quantitative evaluation.

Functional Maps Available

The key perfusion parameters that CT Perfusion 4D package offers (measured in their appropriate physiological units) are:

1. **Regional Blood Volume (BV; ml/100g).**
Characterizes the amount of blood present in a local region; this is affected by the size of the blood vessels and the number of open capillaries.
2. **Regional Blood Flow (BF; ml/min/100 g).**
It is a measure of the volume of blood that flows through a local tissue region within a period of time.
3. **Regional Mean Transit Time (MTT; sec).**
Mean Transit time is the time that blood requires to traverse the tissue circulation from arterial inlet(s) to the venous outlet(s). Mathematically:
 $rMTT = rCBV / rCBF$ (Cerebral), or
 $rMTT = rTBV / rTBF$ (Tissue)
4. **Contrast Arrival Delay (IRF T0)**
It is the delay of arrival of contrast to any given location from its arrival at the artery ROI.
5. **Transit Time to IRF Peak (Tmax; sec)**
It is the transit time to the peak of the Impulse Residue Function (IRF).
Mathematically:
 $Tmax = IRF\ T0 + MTT/2$
6. **Capillary Permeability Surface Area Product (PS; ml/min/100 g)**
It is the unidirectional rate of transfer of blood solutes (contrast media) through the capillary endothelium into the interstitial space. PS is given in the same units as blood flow, ml/min/100g. Permeability Surface maps are intended to be used with lengthened acquisition times, typically 60 seconds are greater.
7. **Base Image** corresponds to an average of all pre-enhancement images. It is used in baseline subtraction and may also be as a reference image.
8. **Average Image** corresponds to the as the average of all images at a given slice location.

CONFIDENTIAL


SECRET DE AFACERI

GE imagination at work

© 2021 General Electric Company.
All rights reserved. Data subject to change.
GE and GE Monogram are trademarks of General Electric Company.
* Trademark of General Electric Company
DOC1999217

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
VessellQ™ Xpress with AutoBone™ Xpress

Fast and Easy to Use Vessel Analysis with Automated Real-Time Tracking, Auto-Labeling and Simplified User Interface.

CONFIDENTIAL

SECRET DE AFACERI

Contents

1. Features.....	4
2. System Requirements.....	4
3. Indications for Use.....	4
4. Regulatory Compliance.....	4
5. Features Detail.....	4
5.01 Anatomy-based protocols.....	4
5.02 AutoBone™ Xpress.....	5
5.03 Auto Aorta.....	5
5.04 Fast Tracking.....	5
5.05 Vessel Label Data Base.....	5
5.06 Dynamic AVA.....	5
5.07 Advanced Vessel Analysis Tools.....	6
5.08 Stenosis Detection & ColorID.....	6
5.09 Quick AVA.....	6
5.10 Summary Table & Report.....	7
5.11 Multiple Study Comparison.....	7
5.12 Save State & Protocol.....	7
5.13 Summary.....	7

CONFIDENTIAL
SECRET DE AFACERI

Peripheral vascular disease (PVD) is a condition similar to coronary artery disease and carotid artery disease. In PVD, fatty deposits build up in the inner linings of the artery walls. These blockages restrict blood circulation, mainly in arteries leading to the kidneys, stomach, arms, legs and feet.

There are different clinical techniques used to diagnose PVD. Computed Tomography Angiography (CTA) can potentially provide sufficient information for therapeutic decision-making and reduced imaging costs for the management of PVD. One of the advantages of computed tomography (CT) is that it shows both vascular & nonvascular anatomy and also provides 3D information for best view. Also, CT may be more useful for patients with limited or no peripheral access or with stents, who may not be candidates for magnetic resonance angiography (MRA).

Complementary to the ability to perform CTA is to perform clinical assessment by a dedicated post-processing tool.

Highlights

- Automated detection of aorta and iliacs with auto-labeling of vessels
- Automated real-time fast tracking for all vasculature
- New contextual wheel menu including easy access to advanced options
- New centerline edition capabilities
- Improved predefined measurements and edit capabilities
- New reporting interface which includes measurements in summary table
- New user interface

Overview

VesselIQ™ Xpress with Autobone Xpress is a post-processing software package for the Advantage Workstation (AW) platform, AW Server, CT scanners and PACS reading stations. It is a tool to analyze CT angiographic data including stenosis analysis, thrombus, pre and post stent planning procedures, and directional vessel tortuosity visualization. Autobone Xpress provides zero-click bone segmentation for head, neck and other anatomy. The automation and ease of use streamlines workflow.

CONFIDENTIAL
SECRET DE AFACERI

Visit us:

<https://www.gehealthcare.com/products/advanced-visualization/all-applications/autobone-vesseliq-xpress>

Features

- Anatomy based protocols for guided workflow in the analysis for your peripheral vascular CTA acquisitions.
- Automatic bone removal for efficient review of vasculature.
- Automated abdominal aorta CTA protocol with vessel tracking and vessel labeling.
- Fast tracking detects entire length of vessel and displays the centerline based on a single start point.
- Advanced editing tools bringing new flexibility for faster editing based on targeted anatomy.
- Thrombus detection with measurements and comparison over time.
- Dynamic AVA provides access to vessel tracking at any time, from any protocol, on any vessel.
- Advanced measurement tools for stenosis detection and color identification of HU value for quantification of different densities.
- Create a queue of cases that are pre-loaded and pre-processed in the background while you read another exam.
- Measurements are organized in the summary table allowing for easy viewing and exporting.
- Fully customizable workflows with adjustable layouts, personalized parameter and settings, custom review steps.
- Save State lets you save and restore the state of the processed images at any stage

System Requirements

Refer to AW Applications platform compatibility document.

Indications for Use

VessellQ Xpress is an optional, non-invasive, optimized, post-processing application intended to provide images and tools to analyze vascular anatomy and pathology, aiding physicians in diagnosis and determination of treatment paths, from a set of Computed Tomography (CT) Angiographic images. VessellQ Xpress is an option for the Advantage Workstation (AW) platform, CT Scanner, and/or PACS, which can be used in the analysis of 2D and 3D CT Angiography images/data for the purpose of cardiovascular and vascular disease assessment. This software-only device is designed to support physician assessment for a wide variety of clinical uses such as stenosis analysis, pre/post stent planning, pre/post valve replacement planning, and directional vessel tortuosity visualization. VessellQ Xpress' automatic visualization tools provide users the capability to segment bony structures for accurate identification of the vessels.

Indications for Use (cont.)

Additional tools enable analysis of the vascular anatomy including the aorta, valves, and branching vessels for: anatomical sizing; density and volume analysis of segmented vasculature and calcified / non-calcified plaque; and measurements of abnormalities. The TAVI Analysis option for VessellQ Xpress is a planning tool used for Trans Aortic Valve Implantation (TAVI) procedures. It automatically segments the aorta and displays the aortic valve in multiple views for measurements of anatomic structures commonly needed for aortic annulus replacement planning. TAVI Analysis provides guided workflow and semi-automated tools to aid in evaluation of appropriate access pathways for interventional procedure planning.

AutoBone Xpress option is a software package that is intended to facilitate segmentation of bony structures and calcifications from abdominal and extremity CT Angiography data.

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Features Detail

Anatomy-based protocols

- The program offers you predefined protocols for fast, easy review of all vascular structures:
 - Intra-cerebral
 - Carotids
 - Thoracic Aorta
 - Abdominal Aorta
 - Runoff
 - Spinal
 - Extremity
- Protocols may be customized to meet your requirements and saved for quick access.

CONFIDENTIAL
SECRET DE AFACERI

Features Detail *(cont.)*

AutoBone™ Xpress

- Anatomical driven protocols with automatic bone removal.
- Improve vessel visualization by removing obstructive bony detail


- Automatic calcium removal for vessel lumen visualization

Auto Aorta

- With a single click, the program automatically removes bones and performs centerline tracking of the aorta and iliacs and labels vessels.



- Simply confirm the automated tracking or modify it by adding additional branch points, deleting a single point, or deleting multiple points from any view.

Fast Tracking

- Semi-automatic detection of complete vessel from a single starting point

- Advanced editing with variable interpolation tool for quick centerline edits adjusted for length of vessel.
 - Bulls eye adjustment capability on x-section views to perfect tracking of the smallest most tortuous vessels.
 - Volume measurements of the thrombus allow you to track change over time. Additional measurements include :
 - Minimum diameter
 - Maximum diameter
 - Mean diameter
 - Cross-sectional area
- 

Vessel Label Data Base

- Each protocol is associated with a database of vessel names
- Customize vessel label database to meet your requirements for every anatomical region.
- A right mouse click provides access to the list of vessel names per anatomy.

Dynamic AVA

- Real-time tracking of the vessel center line and instant visualization of an unfolded view for quick vessel inspection.
- Immediate access to vessel analysis tools.
- Edit contours with a selection of drawing tools with variable interpolation for fast editing.



CONFIDENTIAL

SECRET DE AFACERI

Features Detail (cont.)

Advanced Vessel Analysis Tools


- Semi-automated vessel analysis lets you track, extract, visualize, and size vessels.
- Track a vessel and its branches for quick review in curved, oblique longitudinal or cross-sectional views.
- Multi-point deposition capabilities for improved center line tracking of challenging cases.
- When you edit a vessel contour, the program automatically interpolates the edit to neighboring slices.
- Generate and save rotating movies
 - Curved planar reformation and cross-sectional (lumen) views
 - Best L-section reformatted images
- Select vessel segments for quantitative or qualitative analysis.

- Available measurements include:
 - Distance
 - Cross-sectional area
 - Single or dual reference point comparison
 - Relative percent stenosis / aneurysm
 - Color Identification
 - 3D Angle
 - Thrombus Extraction
 - Volume
 - Mean, min, max diameter
- Move freely between review steps and protocols without leaving the exam.

Stenosis Detection & ColorID

- With the customizable one- or two-point deposit technique, you can automatically quantify degree of stenosis and segment length.
- By assigning customizable color mapping to CT number ranges within an image, you can identify plaque on an axial image, curved reformat, or MPVR image. Up to ten levels of color coding are available. Colors transition smoothly from one density to another.
- Make volume and area measurements to track plaque size over time.
- Custom sizeable 3D ROI lets you apply color over all anatomy.
- Generate multiple ROIs on one vessel or multiple vessels.
- Automatically include each defined plaque, its volume and area into the report.

Quick AVA

- Dynamic AVA, provides access to vessel analysis from any protocol.
- Launch the Quick AVA feature from the main tool bar.
- Perform vessel analysis from any 3D or reformatted image with one- or two-point clicks.
- If you detect a lesion, you can analyze it by dropping a point above and below the section. The vessel analysis tool launches automatically for fast review in curved reformat, cross-sectional, lumen, and MRP view.
- Extend tracking proximally or distally for extended views at any time.

CONFIDENTIAL
SECRET DE AFACERI


Features Detail *(cont.)*

Summary Table & Report

- Measurements are automatically captured in summary table for easy review.


- Simply click on any measurement in the summary table and you will be taken to that measurement location.
- Prepare report with a single click or generate an advanced report with customized options
- Send the report to database, filmer or PDF.

Multiple Study Comparison

With VesselIQ Xpress, compare a previous study with the current study. You also have the flexibility to launch comparison within Quick AVA.

Save State & Protocol

- Save the current state of tracking and measurements for easy reloading with a single mouse click. Save state captures:
 - 3D images
 - Vessel tracking
 - Measurements
 - Vessel contours.
- Transfer the Save State to the PACS or other AW workstations for future review.
- Save custom protocols including anatomical labels and measurements.

Summary

VesselIQ Xpress provides user-friendly tools and an intuitive interface to help streamline your vessel analysis workflow and help you generate more comprehensive reports for referring physicians.

CONFIDENTIAL
SECRET DE AFACERI

About GE HealthCare

GE HealthCare is a leading global medical technology, pharmaceutical diagnostics, and digital solutions innovator, dedicated to providing integrated solutions, services, and data analytics to make hospitals more efficient, clinicians more effective, therapies more precise, and patients healthier and happier. Serving patients and providers for more than 100 years, GE HealthCare is advancing personalized, connected, and compassionate care, while simplifying the patient's journey across the care pathway. Together our Imaging, Ultrasound, Patient Care Solutions, and Pharmaceutical Diagnostics businesses help improve patient care from prevention and screening, to diagnosis, treatment, therapy, and monitoring. We are an \$18 billion business with 51,000 employees working to create a world where healthcare has no limits.

Follow us on [Facebook](#), [LinkedIn](#), [Twitter](#), [Instagram](#) and [Insights](#) for the latest news, or visit our website [gehealthcare.com](https://www.gehealthcare.com) for more information.

CONFIDENTIAL
SECRET DE AFACERI

Products mentioned in the material may be subject to government regulations and may not be available in all countries. Shipment and effective sale can only occur after approval from the regulator. Please check with local GE HealthCare representative for details.

©2024 GE HealthCare. GE is a trademark of General Electric Company used under trademark license.

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
GE HealthCare

Lung VCAR

For lung nodule analysis - automatic visualization,
measurement, reporting and follow-up

CONFIDENTIAL
SECRET DE AFACERI

Contents

1. Features	4
1.01 Review	4
1.02 Analysis	4
2. Indications for use	4
3. System requirements	4
4. Regulatory compliance	4
5. Features detail	5
5.01 Scanning	5
5.02 Review correlated workflow	5
5.03 Lung segmentation	5
5.04 Bookmark management	6
5.05 Layouts	6
5.06 Analysis	6
5.07 Patient report	7
5.08 Summary	7

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SECRET DE AFACERI

Clinical and technical background


CT scanning is utilized as the exam of choice to visualize and analyze complex lung pathology. The detection of pulmonary nodules and assessment of their evolution with CT are of major importance in chest imaging.

Overview

Lung VCAR brings efficient CT pulmonary nodule assessment and diagnosis. The innovative Digital Contrast Agent (DCA) feature automatically visualizes lung nodules to help you confirm the presence or absence of suspicious lesions from 2 to 12 mm in size. Lung VCAR allows automated follow-up for lesion matching by the registration of two or more datasets, automatic lesion classification.

Highlights

- Synchronized 2D, Digital Contrast Agent (DCA) and segmentation analysis.
- Automatic nodule visualization.
- Automatic nodule analysis (volume, doubling time, %growth).
- Automatic follow-up.
- Reporting workflow.

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Visit us:
<https://www.gehealthcare.com/products/advanced-visualization/all-applications/lung-vcar>

Features

Review

- Ability to synchronize multiple images for nodule comparison.
- Ability to review single or multiple exams and compare axial, sagittal, oblique, coronal, and volume-rendered images.
- Automatically propagates previous exam bookmarks to current exam.
- Automatically segments both right and left lung to reduce visual distraction.
- Digital Contrast Agent (DCA) automatically highlights spherical shapes to enhance visualization of suspicious nodules.

Analysis

- Performs automatic segmentation of all nodule types.
- Provides automatic nodule analysis, including % growth, doubling time.
- Automatic Nodule Contour to verify pixels within the volume.

Indications for use

Lung VCAR/AdvantageALA is intended to provide an optimized non-invasive application to measure abnormalities in the lung (for example, nodules, lesions, etc.) from a set of computed tomography (CT) images.

The software is designed to support the physician in confirming the presence or absence of physician identified lung lesions (e.g. nodules). The software allows measurement of volume over time using a consistent standardized measurement protocol, thus providing an estimation of the volume doubling time. The Lung VCAR/AdvantageALA software allows analysis and displays statistics for nodule characterization all the different nodule types. Lung VCAR optional Digital Contrast Agent (DCA) module is an automated highlight feature for the visual identification of possible lesions. Digital Contrast Agent (DCA) is a 3D filter that produces images that highlight spherical (S) or cylindrical (C) anatomical regions, such as nodules, cysts, scars, and vessels. Images are made available to the physician to aid in characterization of suspicious nodules and thus, the patient management care decision process. Lung VCAR/AdvantageALA provides the physician with additional information, meant to complement diagnosis based on classical techniques.

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
SECRET DE AFACERI

Colon VCAR EC

A comprehensive reading workflow solution for colonic lesion detection.


Colorectal cancer is a major health concern for both men and women. Recent studies conclude that early detection of colonic polyps can reduce the mortality rate of colon cancer. You need a protocol that lets you quickly and accurately detect and report results of CT colon scans; a program that optimizes your workflow and improves your reading productivity.

Overview

Colon VCAR EC gives you a comprehensive reading workflow solution for detecting colonic lesions with flexibility and efficient performance. With it you can visualize anatomy that would otherwise be obscured by tagged fluid. Plus digital contrast agent filter highlights potentially suspicious regions.

What's new

- The program allows you to read and problem solve in 2D, 3D, or 360 degree dissection views.
- Prone/Supine polyp linking provides quick navigation between datasets.
- Digital Contrast Agent (DCA) provides automatic visualization of shapes characteristic of polyps.
- Electronic Cleansing (EC) subtracts tagged stool and fluid, making it easier to identify lesions.

CONFIDENTIAL

SECRET DE AFACERI


Visit us:

www.gehealthcare.com/aw/applications/


Features

- Prone and supine synchronized reviewing for problem solving.
- You synchronize prone and supine navigation based on polyp bookmark location and distance from center line.
- Provides you with quick, easy to use polyp and center line measurement tools for complete reporting.
- You can leverage any combination of synchronized 2D, 3D, and 360 views.
- Dissection views make primary reading and problem solving more efficient.

- You can vary navigation speed.
- Lets you track entire colon from rectum to cecum for thorough, precise exams.
- You are provided with unrestricted viewing angles of all colon surfaces.
- You can apply color to bookmarked lesions for better visualization.
- The program allows you to assess luminal, mural and extracolonic anatomy and pelvic structures.

- A selection of application-specific tools ensures ease of use and thorough exams.
- Your workflow benefits from fast, interactive electronic patient reporting.

System Requirements

- AW Workstation
- AW Server

Indications for Use

Colon VCAR EC is a CT, non-invasive, image analysis software package that allows the visualization of 2D, 3D and dissected medical image data of the colon derived from DICOM 3.0 compliant CT scans. ColonVCAR is designed to aid the physician in evaluating the lumen and internal wall of the colon to confirm the presence or absence of colonic lesions (e.g. polyps). It provides functionality for 2D/3D rendering, bookmarking of suspected lesions, synchronized viewing of the 2D, 3D and 360 dissection views for data sets acquired in any position, and an object oriented endoluminal display. In comparison to Colonoscopy, this tool has an advantage of depth penetration due to its 3D presentation capability.

It is intended for use by clinicians to process, render, review, archive, print and distribute colon image studies.

The Colon VCAR DCA (Digital Contrast Agent) module is an automated highlight feature for the visual identification of spherical structures within the colon and is intended to be used as concurrent reading device. Digital Contrast Agent is a 3D filter that produces images that highlight spherical anatomical regions, such as polyps, and/or stool. Colon VCAR uses color to display these highlighted spheres. Images are made available to the physician to aid in characterization of potential polyps and thus, the patient management care decision process.

Regulatory Compliance

This product complies with the European CE marking regulation following Medical Devices Directive: Directive 93/42/EEC.

This product has not been approved for sale in the United States by the United States Food and Drug Administration (FDA).

GE imagination at work

SECRET DE
AFACERI

CONFIDENTIAL

© 2015 General Electric Company.
All rights reserved. Data subject to change.
GE and GE Monogram are trademarks of General Electric Company.
* Trademark of General Electric Company.

Colon VCAR EC

Features Detail

Electronic Cleansing


Electronic Cleansing, accessible in both 2D and 3D review, subtracts tagged stool and fluid from the generated data making it easier to identify lesions that may be visible in only one view.

Tagged material is identified based on Hounsfield units. Processing tasks such as Auto Segmentation, 3D flythrough, and Dissection also take advantage of tagged material subtraction.

Auto Segmentation shows the tagged material while 3D flythrough enables quick anatomy visualization without the tagged material. With this method, the tagged material is removed yet shown in a different color at an adjustable opacity level.

Similarly, the tagged material can be removed from the dissection view as seen in the images below.

User Interface

The Colon VCAR EC user interface has been optimized to facilitate fast, accurate examination of the colon's full extent. It provides seamless integration between automatic and manual interactions for a level of flexibility designed to reduce reading times.

Simultaneous, Indexed Viewing

The program can display both prone and supine

views, which can be manually synchronized for simultaneous and indexed viewing. During review, you can "Lock to Target" with a single button click to quickly problem solve during the primary read. And you are able to bookmark locations of suspected polyps.


360-degree Dissection

This powerful tool provides a 360-degree view of the entire internal large bowel simultaneously displayed in supine and prone positions. This ability speeds viewing and analysis.

This feature lets you leverage all available views: 2D, 3D, and 360-degree dissection during the primary read. For more traditional reading methods, a 2D, 3D, or Dissection view can be laid out and reserved for primary reading with reference to other views for problem solving.


Joystick Navigation

Colon VCAR EC lets you adjust the speed and direction of a flythrough simply by using the mouse. Flythroughs can also be set on autopilot for hands-off reviewing. Additionally, a single mouse click lets you lock onto a region of interest and fly around the target anomaly.

Computer Aided Reading

Colon VCAR EC's Contrast Agent (DCA), a 3D filter, highlights spherical anatomy such as polyps in color. The resulting images facilitate the characterization of potential polyps.

Virtual Biopsy

Problem solving while 3D reading is facilitated by the Virtual biopsy tool, which allows you to view areas in the navigator in color-coded or black/white scales. Different densities, determined by Hounsfield units, can be displayed using the color scale.


GE imagination at work

SECRET DE
AFACERI

CONFIDENTIAL

Colon VCAR EC -1

Auto-Segmentation


Colon VCAR EC quickly performs the time-intensive task of colon segmentation for both prone and supine data sets. All you do is select the protocol.

The program automatically defines a luminal centerline in both prone and supine positions. This can be applied to a high percentage of cases, depending on the level of insufflation and collapse.

Semi-automatic and manual segmentation modes are also available. Semi-automatic segmentation provides complementary, and often necessary, assistance in cases when the colon is collapsed and operator intervention is required. This mode is typically used to connect colon segments around a collapse to form a single centerline. Manual segmentation is typically used to extend the centerline to the end of the cecum for completeness. Manual editing is also available to refine automatic tracking if needed.

Dissection (centerline) Contouring Verification

Using a color mapping technique, Dissection Contouring Verification outlines colon anatomy included in the centerline trace. You can review verification in any orthogonal plane.

Automatic High Resolution Flythrough

Real time 512-matrix resolution automatic flythrough enhances navigational image quality, and provides operational flexibility, functionality, and simplicity. Anatomy remains sharp with minimal pixilation. Steering remains centered as you navigate through the colon. You can specify your navigation speed, and 3D navigation is correlated with 2D and 360-degree Dissection views for rapid synchronized referencing in primary or secondary reading.

512 matrix resolution

Save State

This feature gives you the ability to save any centerline edits, bookmarks, and colored polyps for later reloading. It is also convenient for quick recall during consultation or exam review.

Prone/Supine Polyp Linking


After bookmarking a polyp in prone and supine datasets, the program provides you with the capability to match and label them for easy tracking. This gives you the ability to read registered prone/supine datasets for a one pass read.

Synchronized Prone/Supine 3D Dissection View Navigation

Once prone/supine polyp linkage is made, the program co-registers and links the datasets and polyps, letting you visualize opposing anatomy for a synchronized review.

Small Bowel Extraction

This feature quickly segments the small bowel, giving you an unobstructed view of the large intestine.

Analysis Mode

The program's Analysis Mode provides a variety of tools to help streamline your workflow:

- Lock Mode lets you target and possible lesion and stay locked on it at all times.
- Synchronized display can be customized to show two oblique images, axial, sagittal, coronal, volume rendered and/or 3D views within the colon in a Lock mode. Leveraging multiple synchronized planes and views in a problem-solving mode may enhance your reading productivity.

GE imagination at work


**SECRET DE
AFACERI**

CONFIDENTIAL

Colon VCAR EC -2

Color Centerline Tracking

With Colon VCAR EC you can display a real-time 3D model of the colon with a defined centerline to quickly orient the location of the area you are viewing to the overall colon. This synchronization is always linked and available for quick reference.

Patient Report*

Customizable, interactive reporting capability gives greater structure and flexibility in generating your reports. You can choose how many images to display per sheet, place your facility's logo, apply physician signatures, list patient demographics, and detail your comments and impressions. Export your report to CD, HTTP, or print it to paper.

The Quick Report panel gives you instant access to report input fields and provides access to automatic bookmark annotation on the colon map.

Measurement Tools


With Colon VCAR EC, you can measure polyp distance from distal rectum to polyp location in millimeters, perform 2D measurements, and obtain statistics of ROI in Hounsfield units.

Archive

You can define the number of images per polyp you want to archive to either a networked location or to magnetic optical disk.


Customizable Review and Analysis Layouts

Create layouts that match your reading style and save them as protocols for future recall. Viewports can be shaped into a variety of square or rectangular configurations.


Five of the most popular protocol layouts are pre-built for quick loading. Choose from several dissection configurations:

- 360-dissection gives you convenience and productivity benefits.
- 120-dissection lets you read around distortions if necessary.

- Use the 2D read configuration if the exam has poor insufflation and multiple collapses. 2D can also be used for extracolonic review.

- Use the 3-series layout to load a prone, supine, and decubitus series in parallel to obtain a correlated view.

- You can also customize new layouts to suit your own reading style.
- If you work on a one-monitor Advantage Workstation, a single monitor display is also available.

Movie Generator

You can generate movie clips of the area of interest simply by defining the start and end locations with a mouse click. A cross-sectional movie will be generated. You can also generate lock-to-target movies. Movies can be saved as mpeg, avi, HTTP and saved to CD.

GE imagination at work

**SECRET DE
AFACERI**

CONFIDENTIAL

Colon VCAR EC -3

Entire Colon Movie

Generate a fly-through movie of the entire colon with one button click. Export the result as an mpeg or as screen captures to a PACS workstation for review or consultation.

Additional Features

- Automatic center line tracking
- SmartCursor™ for easier navigation
- Synchronized reformatted views
- Multiple correlations of anatomy with CT colonography
- Exploration of suspicious areas in primary 2D reformatted sections
- Primary 3D reformatted sections. Obliques, Navigator, Quick Virtual Dissection, and axial views are all correlated to the same anatomy of interest.
- Black/white or color display

Summary

Colon VCAR EC helps make colon evaluation a less traumatic procedure for patients. Preparation for the exam is less severe; the procedure is less time intensive and less expensive than optical colonoscopy. Moreover, the tools and features this program puts into your hands will streamline your workflow and potentially enable greater diagnostic confidence.

Indications for Use

Colon VCAR EC is a CT, non-invasive, image analysis software package that allows the visualization of 2D, 3D and dissected medical image data of the colon derived from DICOM 3.0 compliant CT scans. ColonVCAR is designed to aid the physician in evaluating the lumen and internal wall of the colon to confirm the presence or absence of colonic lesions (e.g. polyps). It provides functionality for 2D/3D rendering, bookmarking of suspected lesions, synchronized viewing of the 2D, 3D and 360 dissection views for data sets acquired in any position, and an object oriented endoluminal display. In comparison to Colonoscopy, this tool has an advantage of depth penetration due to its 3D presentation capability. It is intended for use by clinicians to process, render, review, archive, print and distribute colon image studies.

The Colon VCAR DCA (Digital Contrast Agent) module is an automated highlight feature for the visual identification of spherical structures within the colon and is intended to be used as concurrent reading device. Digital Contrast Agent is a 3D filter that produces images that highlight spherical anatomical regions, such as polyps, and/or stool. Colon VCAR uses color to display these highlighted spheres. Images are made available to the physician to aid in characterization of potential polyps and thus, the patient management care decision process.

*Report feature is only available on AW Workstation

**SECRET DE
AFACERI**


CONFIDENTIAL

GE imagination at work

© 2015 General Electric Company.
All rights reserved. Data subject to change.
GE and GE Monogram are trademarks of General Electric Company.
* Trademark of General Electric Company.

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