

Passive Middle Ear Implant

Stapesplasty Prosthesis

AW58482_1.0

hearLIFE

English	2
Deutsch	6
Français	11
Italiano.....	16
Español	21
Português	26
Nederlands	30
Dansk.....	35
Svenska	39
Suomi.....	43
Norsk	47
Polski	51
Türkçe	56
Magyar	60
Čeština.....	64
Slovenčina.....	68
Hrvatski	72
Slovenščina.....	76
Srpski jezik	80
Български език.....	84
Română.....	89
Русский.....	93
Ελληνικά.....	98
Lietuvių kalba	103
Latviešu valoda.....	107
Eesti keel	111
Қазақша.....	115
한국어	120
Tiếng Việt	124
中文.....	128
中文.....	131
137	العربية

Instructions for use

Passive Middle Ear Implant Stapesplasty Prosthesis

Device description

The MED-EL passive middle ear implant – stapesplasty prosthesis is an ossicular replacement prosthesis made of titanium. The ossicular chain can be reconstructed to restore the conductive hearing mechanism within the middle ear.

Intended use

The stapesplasty prosthesis is intended to be used as a replacement of the stapes arch or stapes arch and incus in case of a fixed stapes footplate. This passive implant restores the mechanical sound transmission from the tympanic membrane to the oval window. The stapesplasty prosthesis is intended to treat patients of all ages.

Indications

The stapesplasty prosthesis is indicated in case of congenital or acquired defects of the stapes due to e.g.

- otosclerosis
- congenital fixation of the stapes
- traumatic injury
- malformation of the ossicular chain/middle ear

The stapesplasty prosthesis is indicated to treat patients with inadequate conductive hearing from previous stapes surgery.

Contraindications

The stapesplasty prosthesis should not be used if:

- the patient is known to be intolerant to the implant materials
- an acute middle ear infection is present
- otitis externa and/or otological mycosis are present
- tympanic membrane perforation and/or chronic otitis media are present
- a balance problem, e.g. active Meniere's disease, is present

Caution:

The physician must fully assess the potential risks and benefits for the patient and his/her realistic expectations with the device prior to the decision to implant. The physician must exercise medical judgement and consider the patient's complete medical history.

Possible adverse effects

Surgery of the middle ear involves manipulation of the fragile ossicular bones and exposes the inner ear to the risk of surgical trauma. Serious complications may arise either during or after middle ear surgery that may result in irreparable damage to otologic structures causing irreversible partial or total loss of hearing. Subsequent surgical procedures may be required to correct these conditions, if possible.

Complications that may occur include: sensorineural hearing loss due to trauma during surgery; impaired acoustic signal transfer due to friction between prosthesis piston and the stapes footplate perforation; vertigo; tinnitus; granuloma; incus arrosion/erosion/necrosis; retraction or perforation of the tympanic membrane; post-surgery displacement of the prosthesis due to the development of scar tissue; post-surgical translocation of the prosthesis due to trauma or incorrect position; extrusion of the prosthesis; and post-surgical infection, including otitis media. It is also possible that the prosthesis could cause a significant loss in residual hearing. The responsible physician should ensure that the patient is fully apprised of the possible adverse effects of the implant procedure.

Additionally, patients are exposed to the normal risk of surgery and anaesthesia. Major ear surgery may result in numbness, swelling, or discomfort around the ear, the possibility of facial palsy, disturbance of balance or taste, or neck pain. If these occur, they are usually transient and resolve within a few weeks after surgery.

Warnings and precautions

Pre-operative

- The device is intended for single use only. Do not reuse or resterilize.
- Do not use if the sterile packaging is damaged, or if the use-by date (expiration date) has passed.
- Open the package only in a sterile environment to avoid contamination.
- Do not alter the device and only use it as intended.
- Do not use a damaged device.
- To eliminate unwanted electrostatic effects, water the prosthesis with a few drops of sodium chloride solution (0.9%) before you remove the prosthesis from its package.
- Do not implant if the patient is known to be intolerant to the implant material (titanium). Consider patient intolerances for surgical tool selection.
- Do not implant the device in direct contact with passive middle ear implants other than from MED-EL.

Intra-operative

Ossiculoplasty is a sophisticated surgical intervention and should be performed by a trained surgeon. Inadequate surgical technique can lead to irreparable damage.

MED-EL offers a variety of ossicular prostheses. The correct type and length of the prosthesis to be implanted is determined intra-operatively. Care must be taken when choosing the proper type and length of the prosthesis to be implanted to avoid post-operative complications (e.g. damage to the inner ear if the selected prosthesis is inserted too deep into the inner ear; insufficient contact between the prosthesis and the incus/malleus). A prosthesis insertion depth of 0.5 mm into the inner ear is recommended.

For prostheses with shaft it is possible to carefully bend the shaft max. three times, max. 45° in any direction.

Avoid accidental bending of the prosthesis in order to prevent functional damage.

The middle ear structures must be treated with care.

Assure the inner ear is sealed around the prosthesis piston with tissue to prevent perilymphatic fistula formation.

Post-operative

The patient shall be advised regarding proper hygiene of the operated ear.

After surgery, both surgeon and patient must be vigilant for any signs of complications (such as development of a granuloma or loss of benefit over time).

The physician must inform the patient that the patient should avoid situations such as sudden, loud noises or severe ambient pressure variations and that the patient should consult his/her physician with questions regarding post-operative activities.

The Implant Card shall be filled out and handed over to the patient.

MRI safety information

Non-clinical testing has demonstrated that the device is MR Conditional. A patient with this device can be safely scanned in an MR system meeting the following conditions:

- Static magnetic field of up to 7T
- Maximum spatial field gradient of 79T/m (7900G/cm)
- Maximum force product of 356T²/m (356000000G²/cm)

In non-clinical testing, the worst-case image artefact caused by the device extends approximately 4.59mm if imaged with a gradient echo pulse sequence in a 7T MR system.

Storage and disposal

The product must be stored in its sterile packaging. Devices must not be used after the use-by date specified on the package. Packaging should be disposed of in accordance with local requirements.

Warranty statement

Please contact your local MED-EL representative for information on our warranty provisions.

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© 2008 Pearson Education, Inc.
MED-EL Elektromedizinische Geräte GmbH