B BRAUN
SHARING EXPERTISE

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Vasco[®] Sensitive

NON STERILE EXAMINATION AND PROTECTIVE GLOVES | DATA SHEET

B. Braun Melsungen AG confirms that Vasco[®] Sensitive gloves comply with the following standards and regulations:

EC CERTIFICATES AND APPLIED STANDARDS

Medical Device Class I according to Medical Device Regulation (EU) 2017/745

EN 455 1-4, ASTM D3578

Personal Protective Equipment Category III according to Personal Protective Equipment Regulation (EU) 2016/425

EN 420, EN 374, EN 16523, ISO 16604, ASTM F1671

QUALITY CERTIFICATES

ISO 9001, ISO 13485

PERSONAL PROTECTIVE EQUIPMENT

Information and Declaration of Conformity according to PPER (EU) 2016/425:

www.bbraun.com/gloves-declarations-of-conformity

<https://www.sriranggloves.com/en/update/document>

 Sri Trang Gloves (Thailand), Public Company Limited
10 Soi 10, Phetkasem Road, Hat Yai, Songkhla 90110, Thailand
www.sriranggroup.com

B. Braun Melsungen AG

Dr. Hans-Ulrich Gaudin
Head of Global Regulatory Affairs OPM Germany

Vasco[®] Sensitive

NON STERILE EXAMINATION AND PROTECTIVE GLOVES | REGULATORY INFORMATION

MEDICAL DEVICE INFORMATION

MDR (EU) 2017/745 (CLASS I), EN 455

FOOD COMPLIANCE

Conformity for food contact according to 1935/2004/EEC

PERSONAL PROTECTIVE EQUIPMENT INFORMATION

2777 PPE Regulation (EU) 2016/425 (Cat. III); EN 420:2003+A1:2009

Tested in accordance with:

ISO 374-1/Type B

KPT

Code letter	Test chemical	EN 374-1:2016 Permeation level	EN 374-4:2013 Mean degradation
K	Sodium hydroxide 40 %	Level 6	-18,2 %
P	Hydrogen peroxide 30 %	Level 3	3,3 %
T	Formaldehyde 37 %	Level 5	-28,2 %

Tested acc. to EN 16523-1:2015

Performance levels acc. EN 374-1:2016 +A1:2018	1	2	3	4	5	6
Measured breakthrough times (mins)	> 10	> 30	> 60	> 120	> 240	> 480

Degradation levels indicate the change in puncture resistance of the gloves after exposure to the challenge chemical. NOTE: Where the test specimens gave an increased puncture force after chemical exposure, the result is reported as a negative degradation.

ISO 374-5:2016

VIRUS

AQL < 1.5

Resistance to bacteria and fungi	pass
Resistance to virus	pass

This information does not reflect the actual duration of protection in the workplace and the differentiation between mixtures and pure chemicals. The chemical and penetration resistance has been assessed under laboratory conditions from samples taken from the palm only and relates only to the chemical tested. It can be different if the chemical is used in a mixture. It is recommended to check that the gloves are suitable for the intended use because the conditions at the workplace may differ from the type test depending on temperature, abrasion and degradation. When used, protective gloves may provide less resistance to the dangerous chemical due to changes in physical properties. Movements, snagging, rubbing, degradation caused by the chemical contact etc. may reduce the actual use time significantly. For corrosive chemicals, degradation can be the most important factor to consider in selection of chemical resistant gloves. Before usage, inspect the gloves for any defect or imperfections.

Vasco[®] Sensitive

NON STERILE EXAMINATION AND PROTECTIVE GLOVES | TECHNICAL DATA

SIZE	REF 100/90* pcs.	GLOVE DIMENSIONS (EN 455)	
		Width of palm	Total length
XS	6067500	≤ 80 mm	
S	6067526	80 ± 10 mm	
M	6067549	95 ± 10 mm	≥ 240 mm
L	6067565	110 ± 10 mm	
XL*	6067590	≥ 110 mm	

PHYSICAL PROPERTIES

		Min. specification	Typical value
Wall thickness	Finger	0.08 mm	0.14 mm
	Palm	0.08 mm	0.12 mm
	Cuff		0.08 mm
Force at break	During shelf life	6 N	8.1 N after ageing
Elongation at break	Before ageing	650 %	816 %
	After ageing	500 %	916 %
Tensile strength	Before ageing	18 MPa	28 MPa
	After ageing	14 MPa	25 MPa

GLOVE DESIGN

Colour	natural white
Shape	straight fingers, ambidextrous fitting
Cuff	rolled rim, regular cuff
Surface finish	fingertip textured
Inner glove surface	polymer coated, powder-free
Outer glove surface	chlorinated

GLOVE MATERIAL

Natural rubber latex (NRL)	Protein content ≤ 50 µg/g lower claims are not considered to be reliable given the expected process variation in manufacture and inter-laboratory testing (EN 455-3:2020)
Latex allergy risk	containing natural rubber latex which may cause allergic reactions including anaphylactic reactions

ACCELERATORS

Zn-dithiocarbamate	
Free of thiurames and mercaptobenzothiazoles (MBT)	

LOGISTIC INFORMATION

Dispenser pack	100 / 90 pcs.	240 x 122 x 65 mm (L x W x H)
Transportation carton	10 dispenser packs	340 x 249 x 250 mm (L x W x H)
Shelf life	3 years	
Storage conditions	store at room temperature, protect from dust, humidity, sun light and ozone	

Packaging is made from recycled material

Vasco[®] Sensitive

NON STERILE EXAMINATION AND PROTECTIVE GLOVES | BARRIER PROPERTIES – CHEMICALS

Tested by SATRA, UK in accordance with

EN 16523-1: Determination of material resistance to permeation by chemicals.

CHEMICAL	CAS REGISTRY NO.	PERMEATION PERFORMANCE LEVEL	BREAKTHROUGH TIME
Acetone	67-64-1	not recommended	immediate
Acetonitrile	75-05-8	not recommended	immediate
Chloroform	67-66-3	not recommended	immediate
Dichloromethane	75-09-2	not recommended	immediate
Diethyl amine	109-89-7	not recommended	immediate
Diethyl ether	60-29-7	not recommended	immediate
Dimethylsulfoxide DMSO	67-68-5	not recommended	immediate
Ethanol 70 %	64-17-5	not recommended	immediate
Ethidium bromide 1 %	1239-45-8	level 6	> 480 min
Ethyl acetate	141-78-6	not recommended	immediate
Formaldehyde 37 %	50-00-0	level 5	> 240 min
Gasoline	8032-32-4	not recommended	immediate
Heptane-n	142-82-5	not recommended	immediate
Hexane-n	110-54-3	not recommended	immediate
Hydrogen peroxide 30 %	7722-84-1	level 3	> 60 min
Methanol p.a.	67-56-1	not recommended	immediate
Nitric acid 10 %	7697-37-2	level 1	> 10 min
Nitric acid 65 %	7697-37-2	level 1	> 10 min
Sodium hydroxide 40 %	1310-73-2	level 6	> 480 min
Sulphuric acid 47 %	7664-93-9	level 1	> 10 min
Sulphuric acid 96 %	7664-93-9	level 1	> 10 min
Toluene	108-88-3	not recommended	immediate
Trichlorethane	71-55-6	not recommended	immediate
Xylene	95-47-6	not recommended	immediate

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Wir

We

Sri Trang Gloves (Thailand) Public Company Limited,
(10 Soi 10, Phetkasem Road, Hat Yai, Songkhla 90110 Thailand),

Manufacturer according to Medical Device Regulation (EU) 2017/745

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hereby declare in our own responsibility
that the products

Vasco® Sensitive

Natural Rubber Latex examination and
protective gloves for single use. Non sterile.
Ambidextrous.

Basic UDI-DI: 88591306LC01TD
(article numbers see attachment I)

are in conformity with the requirements of the
Medical Device Regulation (EU) 2017/745

Conformity Assessment Procedure
according to article 52 section 7
of the Regulation named above

Classification
according to annex VIII of the Regulation named
above
Class I

Authorized Representative for the European
Union (Article 11)
MDSS GmbH
Schiffgraben 41
30175 Hannover
Germany

Furthermore, we declare in our own
responsibility that the above-mentioned
medical device/s meet/s the requirements of
the following EU Regulation

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Hat Yai, 2021-05-01

Nattawat.

Nattawat Promthong
Senior Assistant Technical Product Management

**Regulation (EU) 2016/425 on personal
protective equipment**

Classification
Category III (high risks)

Notified Body
SATRA Technology Europe Ltd
Bracetown Business Park
Clonee
Dublin
D15 YN2P
Ireland
Identification number 2777

Valid until 18.05.2023

according to our valid:
Certificate of EN ISO 13485: 2016
certificate number: Q5 099188 0004 Rev. 05
certification body: TÜV SÜD Product Service
GmbH
valid until: 31 October 2023

Certificate of ISO 9001: 2015
certificate number: QS5 099188 0010 Rev. 02
certification body: TÜV SÜD America Inc.
valid until: 11 October 2023

EU Type-Examination Certificate: 2777/10466-
04/E04-01
Issued on: 15.04.2020
valid until: 18.05.2023
Notified Body: SATRA Technology Europe Ltd.

Hat Yai, 2021-05-01

Sureerat

Sureerat Choosri
Product Manager (Glove)

Anlage I / Attachment I

Art.-Nr. / Art. No.	Produktname / Product name	Klasse / Class	Kategorie / Category
6067500	Vasco Sensitive	I	III
6067526	Vasco Sensitive	I	III
6067549	Vasco Sensitive	I	III
6067565	Vasco Sensitive	I	III
6067590	Vasco Sensitive	I	III

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Noi,

**Sri Trang Gloves (Thailand) Public Company Limited,
(10 Soi 10, Phetkasem Road, Hat Yai, Songkhla 90110 Thailand),**

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
Instructions for use are to be used in combination with the specific product related information on each product packaging. The gloves are sold as a bundled unit within the shipping carton. In case this bundled unit is dismantled and products are sold separately, the distributor must ensure that the instructions for use are accompanied with each separate unit.

The gloves are classified as Personal Protective Equipment (PPE) Category III according to PPE Regulation (EU) 2016/425 and have been shown to comply with this regulation through the applicable harmonised European standards. These gloves are designed to provide protection against specific chemicals tested, micro-organisms and particulate radioactive contamination (if applicable). The gloves meet the EN/ISO standards shown on each specific packaging. This PPE is single-use only and to be disposed after contamination.

EXPLANATION OF STANDARDS AND PICTOGRAMS

EN ISO 374-1	Permeation levels are based on breakthrough times (tested acc. EN 16523-1:2015) as follows:						
Type A / B / C	Permeation level acc. EN ISO 374-1:2016 +A1:2018	1	2	3	4	5	6
	Minimum breakthrough time in minutes	10	30	60	120	240	480
	Type A = chemical breakthrough time >30 minutes against at least 6 chemicals from the list Type B = chemical breakthrough time >30 minutes against at least 3 chemicals from the list Type C = chemical breakthrough time >10 minutes against at least 1 chemical from the list						
ABCDEFGHIJKLMNPST	TEST CHEMICALS: A = Methanol / B = Acetone / C = Acetonitrile / D = Dichloromethane / E = Carbon disulphide / F = Toluene / G = Diethylamine / H = Tetrahydrofuran / I = Ethyl acetate / J = n-Heptane / K = Sodium hydroxide 40% / L = Sulphuric acid 96% / M = Nitric acid 65% / N = Acetic acid 99% / O = Ammonium hydroxide 25% / P = Hydrogen peroxide 30% / S = Hydrofluoric acid 40% / T = Formaldehyde 37%						
EN 374-4:2013	The degradation (in %) indicates the change in puncture resistance of the gloves after exposure to the respective challenge chemical.						

This information does not reflect the actual duration of protection in the workplace and the differentiation between mixtures and pure chemicals. The chemical resistance has been assessed under laboratory conditions from samples taken from the palm only and relates only to the chemical tested. It can be different if the chemical is used in a mixture. It is recommended to check that the gloves are suitable for the intended use because the conditions at the workplace may differ from the type test depending on temperature, abrasion and degradation. When used, protective gloves may provide less resistance to the dangerous chemical due to changes in physical properties. Movements, snagging, rubbing, degradation caused by the chemical contact etc. may reduce the actual use time significantly. For corrosive chemicals, degradation can be the most important factor to consider in selection of chemical resistant gloves. Before usage, inspect the gloves for any defect or imperfections.

EN ISO 374-5: 2016	Tested for resistance to penetration according to EN 374-2:2014	EN 421:2010	Protection against particulate radioactive contamination.
	Tested for resistance to penetration by blood-borne pathogens according to EN ISO 374-5 / ASTM F1671 (virus resistance). Resistance to bacteria and fungi – pass Resistance to virus – pass The penetration resistance has been assessed under laboratory conditions and relates only to the tested specimen.		These gloves do not protect against mechanical risks.
			PPE is for single-use only and must not be reused.
	XXXX = Identification number of notified Body responsible for the EU type examination and supervising ongoing conformity.		Before usage read instructions for use carefully.
		EN 420:2003+A1:2009	

PRECAUTIONS FOR USE

Always check the gloves for possible mechanical damage, e.g. holes or tears, before use. Do not use damaged gloves. Glove length is appropriate to the end use where the risk to the wrist area is minimal.

INGREDIENTS / HAZARDOUS INGREDIENTS

Some gloves might contain ingredients which are known to possibly cause skin irritations or allergic reactions with sensitised persons. Check warning information on specific packaging carefully. Formulation available on request.

TEMPERATURE LIMIT

STORAGE INSTRUCTIONS

Keep storage area cool, dry and dust free, avoid ventilation and storage close to photocopy equipment. Protect gloves against ultraviolet light sources, sunlight, oxidizing agents and ozone. Store in original packaging according to the temperature limit, provided on the packaging.

DISPOSAL INSTRUCTIONS

Dispose of the gloves in accordance with the valid regulations for this material. Gloves contaminated with chemical substances must be disposed of in accordance with the regulations for the relevant chemicals.

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