

STERISEPT PLUS

Concentrated disinfection and cleaning agent for all water-tolerant surfaces

- free from phenols, aldehydes and phosphates
- broad spectrum of microbiological activity
- good cleaning properties and material compatibility



Sterisept Plus combines both cleaning and disinfection, which makes it a universal cleaning disinfectant for all surfaces.

Due to its specially selected active ingredients and surfactants, it offers a good material compatibility along with cleaning properties. Sterisept Plus is suitable for use in hospitals, food and pharmaceutical industries, veterinary area, laboratories, dental clinics and in other medical institutions.

Short contact times and low concentrations enable its economical use. Perfect product in higher hygiene requirements, where additional protection is needed.

Sterisept Plus has a broad spectrum of bactericidal (incl. MRSA, VRE), mycobactericidal (incl. *Mycobacterium avium* and *Mycobacterium terrae*), fungicidal and yeasticidal activity. The product is also active against viruses (incl. HBV, HIV, HCV, Vaccinia, Noro, Rota and Adeno).

USAGE:

Sterisept Plus is concentrated. To prepare a working solution, dilute the required amount of concentrate in water at room temperature. For example, mix 5 ml of Sterisept Plus with 995 ml of tap water to create a 0.5% working solution. Apply a sufficient amount of the prepared solution to the surfaces using a wipe or cloth, ensuring complete and even coverage. The contact time required depends on the targeted microorganisms. Rinsing the treated surfaces is not necessary.

High organic load, such as blood and secretions, should be removed before disinfection using a disposable cloth moistened with disinfectant solution.

The ready-made working solution is useable for 30 days.

Store it in airtight containers at room temperature.

Attention - Sterisept Plus cannot be used together with aldehydes-based preparations, therefore avoid using products which contain aldehydes before and after disinfecting with Sterisept Plus.



Sterisept Plus is suitable for use on the following surface materials:

- Metals such as stainless steel, aluminum, copper, and brass.
- Floor coverings such as PVC, linoleum, and ceramic tiles.
- Plastics, silicone, rubber.

Do not use on materials such as polycarbonate or varnished surfaces.

MICROBIOLOGICAL EFFICACY:

EFFICACY	EU STANDARDS	CONTACT TIME	CONCENTRATION
Bactericidal (incl. MRSA, VRE)	EN 13727; EN 14561; EN 13697	15 min	0,25%
Yeasticidal	EN 13624; EN 13697	15 min	0,25%
Mycobactericidal	EN 14348	15 min	0,5%
All enveloped viruses	EN 14476	5 min 15 min	0,5% 0,25%
Rotavirus	EN 14476	5 min 15 min	0,5% 0,25%

Additional test reports: EN 13624 (fungicidal), EN 14476 (adeno, noro); EN 17387 (bactericidal and fungicidal); EN 16615 (bactericidal and yeasticidal).

PRECAUTIONS:

Harmful if swallowed. Toxic to aquatic life with long lasting effects. Causes severe skin burns and eye damage. May cause damage to organs through prolonged or repeated exposure (Oral). May cause respiratory irritation. Do not breathe vapours. Avoid release to the environment.

Wear protective gloves/face protection/protective clothing/protective footwear.

IF SWALLOWED: Rinse mouth. Do NOT induce vomiting.

IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water or shower.

IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Immediately call a poison center/doctor.

Dispose of contents and / or containers in accordance with regulations on hazardous waste or packaging and packaging waste respectively.

INGREDIENTS: 100 g of product contains active ingredients: Didecyl-Dimethyl-Ammonium Chloride (DDAC) **12,5g** (CAS: 7173-51-5); N-(3-aminopropyl)-N-dodecylpropane-1,3-diamine **12,5g** (CAS: 2372-82-9).

PACKAGING:



1000 ml
dosing bottle



5000 ml
canister

