STERILE UV and Blue Light Filtering Hydrophobic Acrylic Foldable  
Diffractive Aspheric Posterior Chamber Intraocular Lens

#### DESCRIPTION

The Clareon™ PanOptix™ Trifocal Hydrophobic IOL Model CNWTT0 (henceforth referred to as Clareon™ PanOptix™ Trifocal IOL) and Clareon™ PanOptix™ Toric Trifocal Hydrophobic IOL Models CNWTT2, CNWTT3, CNWTT4, CNWTT5, and CNWTT6 (henceforth referred to as Clareon™ PanOptix™ Toric Trifocal IOL) are ultraviolet and blue light filtering foldable single-piece posterior chamber intraocular lenses. These lenses are composed of a high refractive index, soft hydrophobic acrylic material, which contains a covalently bonded blue light filtering chromophore. Alcon's proprietary chromophore filters blue light in a manner that approximates the human crystalline lens in the 400 to 475 nm wavelength range (Boettner and Wolter, 1962). In addition to UV absorption, the blue-light filtering chromophore reduces transmittance of blue light wavelengths.

The material is capable of being folded prior to insertion, allowing placement through an incision smaller than the optic diameter of the lens. After surgical insertion into the eye, the lens gently unfolds to a full-size lens body. Each lens has an optical portion and mechanical support elements (haptics). The haptics provide proper positioning of the lens optic within the capsular bag. The optic portion is biconvex and includes an aspheric surface. The optic diffractive structure is in the central 4.5 mm portion of the optic zone and divides the incoming light to create a +2.17 D intermediate and a +3.25 D near add power. The anterior surface of the Clareon™ PanOptix™ Trifocal and Clareon™ PanOptix™ Toric Trifocal Hydrophobic IOLs are designed with negative spherical aberration to compensate for the positive spherical aberration of the cornea. The effects of this aspheric design feature have not been clinically assessed on these Models. Additionally, the Clareon™ PanOptix™ Toric Trifocal Hydrophobic IOLs have a toric component on the posterior surface with axis marks to denote the flat meridian (plus cylinder axis). Alignment of the toric axis marks with the post-operative steep corneal meridian allows the lens to correct preexisting corneal astigmatism. The physical properties of this lens are described in Table 1 and Figures 1, 2 and 3.

Table 1: Physical Characteristics of Clareon™ PanOptix™ Trifocal IOL

| Physical Characteristics             | Description                                                                                                                                                                          |        |        |        |        |        |
|--------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------|--------|--------|--------|--------|
| IOL Model Number                     | CNWTT0                                                                                                                                                                               | CNWTT2 | CNWTT3 | CNWTT4 | CNWTT5 | CNWTT6 |
| Optic Type                           | Single-piece IOL with diffractive aspheric biconvex optic                                                                                                                            |        |        |        |        |        |
| Optics Material                      | Ultraviolet and Blue Light Filtering Hydrophobic Acrylate/Methacrylate Copolymer                                                                                                     |        |        |        |        |        |
| Spectral Transmission                | 10% transmittance at 403 nm (UV) for +20.0 diopter IOL                                                                                                                               |        |        |        |        |        |
| Index of Refraction                  | 1.55 at 35 °C                                                                                                                                                                        |        |        |        |        |        |
| Optic Powers                         | +6.0 to +30.0 diopters (in 0.5 diopter increments);<br>+31.0 through +34.0 diopter in 1.0 diopter increments<br>with a +2.17 diopter intermediate and a +3.25 diopter near add power |        |        |        |        |        |
| IOL Cylinder Powers (Diopters)       | NA                                                                                                                                                                                   | 1.00   | 1.50   | 2.25   | 3.00   | 3.75   |
| Haptic Configuration                 | STABLEFORCE™ Modified-L Haptics                                                                                                                                                      |        |        |        |        |        |
| Haptic Material                      | Ultraviolet and Blue Light Filtering Hydrophobic Acrylate/Methacrylate Copolymer                                                                                                     |        |        |        |        |        |
| Optic Diameter/ $\varnothing_B$ (mm) | 6.0                                                                                                                                                                                  |        |        |        |        |        |
| Overall Length/ $\varnothing_T$ (mm) | 13.0                                                                                                                                                                                 |        |        |        |        |        |
| Haptic Angle                         | 0°                                                                                                                                                                                   |        |        |        |        |        |

#### BIOCOMPATIBILITY TESTING

Potential patient safety risks to the material(s) of this device were evaluated through nonclinical physicochemical characterization and biocompatibility testing in accordance with international standards applicable to IOL devices. Nonclinical testing demonstrated no safety concerns for local or systemic toxicity, that the IOL material was physically and optically stable, and that there were no leachable substances arising from the manufacturing process (including sterilization) or device material(s) that posed a safety risk. The device possesses an acceptable patient safety profile when used in accordance with the Directions for Use (DFU) for its intended clinical purpose as an ocular implant device.

#### HOW SUPPLIED

The Clareon™ PanOptix™ Trifocal Hydrophobic IOL is supplied dry, in a package terminally sterilized with ethylene oxide, and must be opened only under aseptic conditions (see DIRECTIONS FOR USE).

Figure 1: Design of the Clareon™ PanOptix™ Trifocal IOL  
Model CNWTT0  
(All dimensions in millimeters)


Figure 2: Design of the Clareon™ PanOptix™ Toric Trifocal IOL  
Models CNWTT2-CNWTT6  
(All dimensions in millimeters)



Figure 3: Spectral Transmittance Curves

#### NOTES:

- The cutoff wavelength and the spectral transmittance curves presented here represent the range of transmittance values of IOLs made from hydrophobic acrylate/methacrylate copolymer with bonded UV-absorber and Alcon's proprietary blue light filtering chromophore.
- Measurements were by direct transmittance using Clareon™ IOLs with center thickness equivalent to the marketed range.
- Human crystalline lens data from Boettner and Wolter (1962).

#### INTENDED USE

The Clareon™ PanOptix™ Trifocal IOL is intended for use by a trained ophthalmic surgeon. The IOL is intended to be placed in the capsular bag in the posterior chamber of the eye, replacing the human crystalline lens. This position allows the IOL to function as a refractive medium in the correction of aphakia. The aspheric biconvex optic compensates for the positive spherical aberration of the cornea as compared to a standard spherical optic. The diffractive structure on the anterior surface divides incoming light to provide a range of vision from distance to intermediate to near. This IOL provides an option for clinicians to give patients with an intermediate add power of +2.17 D and a near add power of +3.25 D.

The Clareon™ PanOptix™ Toric Trifocal IOL is intended for use by a trained ophthalmic surgeon. The IOL is intended to be placed in the capsular bag in the posterior chamber of the eye, replacing the human crystalline lens. This position allows the IOL to function as a refractive medium in the correction of aphakia. The aspheric biconvex optic compensates for the positive spherical aberration of the cornea as compared to a standard spherical optic. The diffractive structure on the anterior surface divides incoming light to provide a range of vision from distance to intermediate to near. This IOL provides an option for clinicians to give patients with an intermediate add power of +2.17 D and a near add power of +3.25 D. Additionally, these IOLs have a toric component on the posterior surface with axis marks to denote the flat meridian (plus cylinder axis). Alignment of the toric axis marks with the post-operative steep corneal meridian allows the lens to correct preexisting corneal astigmatism. The performance of the blue light filtering chromophore in reducing the incidence of retinal disorders has not been established. The astigmatic correction at the corneal plane for each Clareon™ PanOptix™ Toric Trifocal IOL model is shown in Table 2.

**Table 2: Cylinder Power and Corneal Astigmatism Correction Range**

| Lens Model | Cylinder Power |                            |
|------------|----------------|----------------------------|
|            | IOL Plane      | Corneal Plane <sup>a</sup> |
| CNWT2      | 1.00           | 0.65                       |
| CNWT3      | 1.50           | 0.98                       |
| CNWT4      | 2.25           | 1.47                       |
| CNWT5      | 3.00           | 1.96                       |
| CNWT6      | 3.75           | 2.45                       |

<sup>a</sup>Based on an average pseudophakic human eye using an SRK/T optical A-constant of 119.1

#### INDICATIONS

The Clareon™ PanOptix™ Trifocal Hydrophobic IOL is indicated for primary implantation in the capsular bag in the posterior chamber of the eye for the visual correction of aphakia secondary to removal of a cataractous lens in adult patients with and without presbyopia, who desire near, intermediate and distance vision with increased spectacle independence.

The Clareon™ PanOptix™ Toric Trifocal IOL is intended for primary implantation in the capsular bag in the posterior chamber of the human eye for the visual correction of aphakia and pre-existing corneal astigmatism secondary to removal of a cataractous lens in adult patients with and without presbyopia, who desire near, intermediate, and distance vision with increased spectacle independence.

#### QUALIFIED COMBINATIONS FOR IOL IMPLANTATION

During implantation of the Clareon™ PanOptix™ Trifocal and Clareon™ PanOptix™ Toric Trifocal IOLs, an Alcon qualified delivery system and ophthalmic viscosurgical device (OVD) combination should be used. The use of an unqualified combination may cause damage to the lens and potential complications during the implantation process. The qualified combinations that can be used with this lens are listed in Table 3.

**Table 3: Qualified Combinations of Compatible Products**

| IOL Model                | Diopter Range | Cartridge                      | Handpiece                                             | OVD            |
|--------------------------|---------------|--------------------------------|-------------------------------------------------------|----------------|
| CNWT0<br>CNWT2-<br>CNWT6 | +6.0 to +25.0 | MONARCH™ III D<br>(8065977763) | CLAREON™ MONARCH™ IV<br>(silver/gray)<br>(8065977774) | VISCOAT™ OVD   |
|                          |               |                                |                                                       | PROVISC™ OVD   |
|                          |               |                                |                                                       | DISCOVISC™ OVD |
|                          | +6.0 to +30.0 | MONARCH™ III C<br>(8065977762) |                                                       | VISCOAT™ OVD   |
|                          |               |                                |                                                       | PROVISC™ OVD   |
|                          |               |                                |                                                       | DISCOVISC™ OVD |
|                          | +6.0 to +34.0 | MONARCH™ II B<br>(8065977758)  |                                                       | VISCOAT™ OVD   |
|                          |               |                                |                                                       | PROVISC™ OVD   |
|                          |               |                                |                                                       | DISCOVISC™ OVD |

Contact Alcon for any additional Alcon qualified ophthalmic viscosurgical devices, handpieces, and cartridges for use with this lens.

#### CALCULATION OF LENS POWER

Accurate biometry is essential for successful visual outcomes. Preoperative calculation of required lens power for the Clareon™ PanOptix™ Trifocal IOL and Clareon™ PanOptix™ Toric Trifocal IOL should be determined by the surgeon's experience and preference. A reference provisional SRK/T A-Constant value for optical biometry equipment such as IOLMaster<sup>2</sup> or LenStar<sup>2</sup> as well as contact ultrasound biometry is listed on the outer label. This reference A-Constant anticipates the use of both corneal power and axial length values from optical biometry or contact ultrasound equipment with standard settings for a typical patient population and a spectacle far point at 6 meters. In general, A-Constants must be "personalized" to compensate for such things as differences in instrumentation, surgical techniques, and IOL power calculation that may exist between clinical practices. IOL power calculation methods are often included with biometry equipment, and they are also described in the references (Hoffer 1993; Holladay 1997; Olsen 2007; Retzlaff, Sanders & Kraff 1990; Haignis 2014).

#### CONTRAINDICATIONS

There are no known contraindications to the use of Clareon™ PanOptix™ Trifocal IOL and Clareon™ PanOptix™ Toric Trifocal IOL when used as recommended.

#### WARNINGS

- The Clareon™ PanOptix™ Trifocal and Clareon™ PanOptix™ Toric Trifocal IOLs are intended for implantation in the capsular bag only. There are no clinical data to demonstrate its safety and performance for placement in the ciliary sulcus.
- DO NOT re-sterilize the IOL by any method.
- DO NOT reuse the IOL. This device is for single use only. Reuse of this single-use device may result in serious injury, such as endophthalmitis.
- DO NOT implant the IOL if the sterility has been compromised or if the sterile package has been unintentionally opened before use.
- Some visual effects may be expected due to the superposition of focused and unfocused multiple images. These may include some perceptions of halos, radial lines around point sources of light (starbursts) under nighttime conditions, or glare, as well as other visual symptoms. As with other multifocal IOLs, there is a possibility that visual symptoms may be significant enough that the patient will request explant of the multifocal IOL.
- A reduction in contrast sensitivity as compared to that expected with a monofocal IOL may be experienced by some patients and may be more prevalent in low lighting conditions. Therefore, patients implanted with multifocal IOLs should exercise caution when driving at night or in poor visibility conditions.
- The physician should consider the following points that are important to the use of this and any other multifocal IOL:
  - The surgeon must target emmetropia to achieve optimal visual performance.
  - Patients with significant preoperative (determined by keratometry) or expected postoperative astigmatism  $\geq 1.0$  D may not achieve optimal visual outcomes.
  - Care should be taken to achieve IOL centration as lens decentration may result in a patient experiencing visual disturbances under certain lighting conditions.
- Physicians considering lens implantation under any of the following circumstances should weigh the potential risk/benefit ratio: Patients in whom the posterior capsule is ruptured, zonules are damaged, or primary posterior capsulotomy is planned.

- Rotation of the Clareon™ PanOptix™ Toric Trifocal IOL away from its intended axis can reduce the astigmatic correction. Each degree of misalignment of a toric IOL may reduce the cylinder power effect by approximately 3.3% (Ma 2008). If IOL repositioning is necessary, it should occur as early as possible, prior to IOL encapsulation. Some clinical cases suggest encapsulation is complete within four weeks of implantation (Nishi 2002; Sacu 2005).
- Carefully remove all OVD from both the anterior and posterior sides of the lens at the close of surgery with emphasis on the space between the posterior capsule and IOL. Residual OVD may cause complications including lens rotation resulting in misalignment of this IOL with the intended axis of placement.

#### PRECAUTIONS

- A high level of surgical skill is required for intraocular lens implantation. Surgeons should have completed proper ophthalmic training where they observed and/or assisted in numerous implantations and successfully completed one or more courses on intraocular lens implantation before attempting to implant intraocular lenses.
- Prior to surgery, prospective patients should be informed of the possible risks and benefits associated with this IOL as well as the risks and benefits associated with cataract surgery. After surgery, physicians should provide patients with the implant card.
- As with all multifocal IOLs, the need to wear corrective lenses when reading small print or looking at small objects may vary for each patient after implantation of this IOL.
- Posterior capsule opacification (PCO) may significantly affect the vision of patients with multifocal IOLs sooner in its progression than patients with monofocal IOLs. This may be due to the reduced contrast sensitivity observed with multifocal IOLs.
- The safety and performance of this IOL have not been studied in patients with certain pre-existing conditions and/or intraoperative conditions (listed in Tables 4 and 5) as these patients were excluded from the clinical studies. Patients with pre-existing conditions may not achieve the visual acuity of patients without such conditions. As with the implantation of any IOL, careful preoperative evaluation and sound clinical judgment should be used by the surgeon to decide the benefit/risk ratio before implanting a lens in a patient with one or more of these conditions.
- DO NOT store intraocular lenses at temperatures over 30 °C (86 °F).
- Refer to directions for use (see below) for information on the maximum allowed time for the IOL to stay in the folded condition. Failure to adhere to manufacturer's recommendations may result in IOL damage.
- Care should be taken to remove all OVD from the eye prior to completing surgery.
- Accurate keratometry and biometry in addition to the use of the Alcon Online Toric Calculator (<http://www.mylcon-toriccalc.com>) are recommended to achieve optimal visual results. Variability in any of the preoperative measurements can influence patient outcomes, and the performance of treating eyes with lower amounts of preoperative corneal astigmatism.
- Anatomical and/or surgical factors may be related to the likelihood that a toric IOL could be placed incorrectly or rotate away from the intended position after placement. Some of these factors can be identified before or during the surgery, but others cannot. If a secondary surgical intervention is necessary to reposition the IOL, explanation should be considered as some subjects may have recurrent or persistent issues related to rotational instability and misalignment.
- As with any surgical procedure, there is risk involved. Thus, the patient should be followed postoperatively as medically appropriate. Potential complications accompanying cataract and/or IOL implantation surgery may include the following: lens epithelial cell on-growth, tissue damage (e.g., capsular tissue damage, iris prolapse), posterior capsule opacification, corneal cell endothelial damage, infection (endophthalmitis), inflammatory reactions [e.g., toxic anterior segment syndrome (TASS), hypopyon, vitritis, and cystic membrane], hyphema, retinal detachment, cystoid macular edema, corneal edema, corneal decompensation, pupillary block, pigment dispersion, color misperception, transient or persistent glaucoma, and secondary surgical interventions. Secondary surgical interventions include: lens repositioning, lens replacement, vitreous aspiration, iridectomy for pupillary block, wound leak repair, and retinal detachment repair.

**Table 4: Preexisting Conditions with No Safety and Performance Data**

|                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <ul style="list-style-type: none"> <li>Irregular corneal astigmatism</li> <li>Significant irregular corneal aberration</li> <li>Corneal irregularity (including irregularity due to dry eye syndrome)</li> <li>Retinal conditions or predisposition to retinal conditions, previous history of, or a predisposition to, retinal detachment or proliferative diabetic retinopathy, in which future treatment may be compromised by implanting this lens.</li> <li>Multifocal IOLs may decrease the level of retinal detail on exam or during treatment slightly and this could make laser and retinal surgeries and the diagnosis of some conditions more challenging (for example, early diabetic retinopathy when only 1 or 2 microaneurysms are present).</li> <li>Subjects with diagnosed degenerative visual disorders (e.g., macular degeneration or other retinal disorders) that are predicted (by subjective assessment of the retina) to cause future acuity losses to a level worse than 0.3 logMAR</li> <li>Amblyopia</li> <li>Clinically severe corneal dystrophy (e.g., epithelial, stromal, or endothelial dystrophy), keratitis, keratoconjunctivitis, keratouveitis, keratopathy, or kerectasia</li> </ul> | <ul style="list-style-type: none"> <li>Any Inflammation or edema (swelling) of the cornea</li> <li>Rubella, congenital, traumatic, or complicated cataracts</li> <li>Extremely shallow anterior chamber, not due to swollen cataract</li> <li>History of or current anterior or posterior segment inflammation of any etiology</li> <li>Aniridia</li> <li>Iris neovascularization</li> <li>Glaucoma (uncontrolled or controlled with medication)</li> <li>Microphthalmos</li> <li>Optic nerve atrophy</li> <li>Previous corneal transplant</li> <li>Pre-existing ocular conditions which may negatively impact stability of the implant</li> <li>Diabetic retinopathy</li> <li>Previous refractive surgery</li> <li>Color vision deficiencies</li> <li>Pregnancy</li> <li>Current or recent usage of alpha-1-selective adrenoceptor blocking agent or antagonist of alpha 1A adrenoceptor [e.g., Flomax<sup>®</sup> (tamsulosin HCL), Hytrin<sup>®</sup>, or Cardura<sup>®</sup>]</li> </ul> |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

**Table 5: Intraoperative Conditions with No Safety and Performance Data**

|                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <ul style="list-style-type: none"> <li>• Other planned ocular surgery procedures, including LASIK, astigmatic keratotomy, and limbal relaxing incisions</li> <li>• Excessive iris mobility</li> <li>• Mechanical or surgical manipulation required to enlarge the pupil to 4.5 mm or larger just prior to IOL implantation</li> <li>• Significant vitreous loss</li> <li>• Significant anterior chamber bleeding</li> <li>• Uncontrollable positive intraocular pressure</li> </ul> | <ul style="list-style-type: none"> <li>• Complications in which the IOL stability could be compromised, including: <ul style="list-style-type: none"> <li>◦ Zonular damage, separation, or rupture</li> <li>◦ Capsulotomy by any technique other than a circular tear or femtosecond laser</li> <li>◦ The presence of radial tears known or suspected at the time of surgery</li> <li>◦ Situations in which the integrity of the circular tear cannot be confirmed by direct visualization</li> <li>◦ Cataract extraction by techniques other than phacoemulsification or liquefaction</li> <li>◦ Situations where the need for a large capsulotomy can be anticipated (e.g., diabetics, retinal detachment in the fellow eye, peripheral retinal pathology, etc.)</li> <li>◦ Capsular rupture or capsulorhexis tear</li> <li>◦ Bag-sulcus, sulcus-sulcus or unknown placement of the haptics</li> </ul> </li> </ul> |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

#### SELECTION OF TORIC IOL

For selection of the appropriate toric IOL, Alcon provides a web-based tool ([www.myalcon-toricalc.com](http://www.myalcon-toricalc.com); Abulafia, Barrett, et al. 2015 and Abulafia, Hill, et al. 2015) that uses pre-operative biometry data, incision location, and surgically induced corneal astigmatism. The astigmatism to be corrected should be determined from biometry data rather than refractive data since the presence of lenticular astigmatism in the crystalline lens to be removed may influence results.

#### AXIS PLACEMENT OF TORIC IOL

For optimal results, the surgeon must ensure the correct placement and orientation of the lens within the capsular bag. The posterior surface of the IOL is marked with indentations (three at each end) at the haptic/optic junction that identify the flat meridian of the toric IOL optic. These indentations form an imaginary line representing the plus cylinder axis (note: IOL cylinder steep meridian is 90° away). The cylinder axis marks should be aligned with the post-incision steep corneal meridian (intended axis of placement) or as determined by the web-based IOL calculator.

Prior to surgery, the lens placement axis should be marked. One way used in clinical practice to mark the eye is as follows: With the patient sitting upright to prevent cyclotorsion, clearly and precisely mark two reference positions with a surgical skin marker or a marking pencil indicated for ophthalmic use. Using these marks as reference points, an axis marker can be used immediately prior to or during surgery to mark the optimal axis of lens placement identified by the web-based IOL calculator ([www.myalcon-toricalc.com](http://www.myalcon-toricalc.com)).

Precisely align the axis marking indentations on the toric IOL with the marked intended axis of lens placement. Remove all OVD from both the anterior and posterior sides of the lens since residual OVD may allow the lens to rotate causing misalignment of the toric IOL with the intended axis of placement. This may be accomplished by manipulating the IOL optic with the I/A tip and using standard irrigation/aspiration techniques to remove all OVD from the eye. Bimanual techniques may be used, if preferred, to ensure removal of OVD from behind the lens implant. Special care should be taken to ensure proper positioning of the toric IOL at the intended axis following OVD removal.

Misalignment of the axis of the lens versus the intended axis of placement may compromise its astigmatic correction. Such misalignment can result from inaccurate keratometry or marking of the cornea, inaccurate placement of the toric IOL axis during surgery, an unanticipated surgically induced change in the cornea, or physical rotation of the toric IOL after implantation. In order to minimize this effect, the surgeon should be careful to ensure that preoperative keratometry and biometry is accurate and that the IOL is properly oriented prior to the end of surgery.

#### DIRECTIONS FOR USE

1. Examine the label on the outer box for model, spherical equivalent powers (base power and add powers), proper configuration, and expiration date.
2. After the outer box is open, verify that the lens primary label information (e.g., model, powers, and serial number) is consistent with information on the outer box labeling.
3. Inspect the primary package carefully for tears, cuts, punctures, or other signs that the pouch has been opened or damaged. This device is sterile until the inner primary package is opened. DO NOT implant the IOL if the sterility has been compromised or if the sterile package has been unintentionally opened before use (see RETURNED GOODS POLICY).
4. Open the undamaged primary package and transfer the case to a sterile environment. Carefully open the case to expose the lens.
5. Use instrumentation that is scrupulously clean to minimize the occurrence of marks. Any forceps used for lens handling must have round edges and smooth surfaces.
6. Handle the IOL by the haptics only when removing the lens from the case. DO NOT grasp the optical area with forceps. Handle lenses carefully to avoid damage to lens surfaces or haptics. DO NOT attempt to reshape haptics in any way.
7. Examine the IOL carefully prior to insertion to ensure that particles have not adhered during handling.
8. Implant the IOL with the most appropriate surgical procedure for the patient. Surgeons should verify that appropriate instrumentation is available prior to surgery. Refer to **Table 2** for qualified combinations of compatible products.

**NOTE: During lens loading and insertion, DO NOT allow the IOL to remain in a folded condition within the selected IOL Delivery System for more than 3 minutes prior to completing insertion into the capsular bag.**

#### MAGNETIC RESONANCE COMPATIBILITY

The Clareon™ PanOptix™ Trifocal and Clareon™ PanOptix™ Toric Trifocal IOLs are magnetic resonance (MR) Safe. The IOL consists of acrylate/methacrylate copolymer material, which is a non-conducting, non-metallic, non-magnetic material that poses no known hazards in all magnetic resonance imaging environments.

#### LIFETIME OF THE IOL

Based upon material characterization of the Clareon™ IOL material, the IOL is expected to be stable indefinitely over the lifetime of the patient.

#### PATIENT IMPLANT CARD

The Patient Lens Implant Card included in the package is to be completed and given to the patient, together with instructions to keep the card as a permanent record to be shown to any eye care practitioner that the patient consults in the future.

To fill out the patient implant card:

1. Remove implant card sticker from the label set inside carton. The implant card sticker is found within the lower left quadrant of the label set and contains blank lines for fillable information.
2. Adhere sticker to back of implant card.
3. Fill out the following information on the card:
  - Date of surgery,
  - Eye implanted [mark left (L) or right (R)],
  - Patient name,
  - Surgeon name, and
  - Hospital or health institution name and address.

#### SERIOUS INCIDENT REPORTING

Any serious incident that may reasonably be regarded as device related should be reported to Alcon Canada, Inc.:

By Phone:

In Canada-(800) 613-2245

Website:

<http://www.alcon.ca>

Each IOL is identified by a serial number which provides traceability, and this information should be given to Alcon.

NOTE: In the EU, serious incidents must also be reported to the competent authority for medical devices of the appropriate Member State.

#### EXPIRATION DATE

Sterility is guaranteed until the use-by date unless the primary sterilization package is damaged or opened. The use-by date is clearly indicated on the outer box label of the carton. Any lens held after the use-by date should be returned to Alcon Laboratories, Inc. (see RETURNED GOODS POLICY).

#### RETURNED GOODS POLICY

In Canada, returned lenses will only be accepted for credit. All returns must be accompanied by an Alcon Canada, Inc. Returned Goods Number and should be shipped via traceable means. A Returned Goods Number is obtained by contacting Alcon Canada, Inc. Customer Service Department. Issuance of this number does not constitute final acceptance of the returned products. For detailed policy guidelines including exchange, please contact your Sales or Customer Service Representative.

#### CLINICAL BENEFITS OF Clareon™ PanOptix™ TRIFOCAL HYDROPHOBIC IOL and Clareon™ PanOptix™ TORIC TRIFOCAL IOL

Clinical evidence from published peer-reviewed clinical literature, clinical experience, and clinical investigations with parent IOLs establishes an acceptable safety and performance profile for Clareon™ PanOptix™ Trifocal IOL and Clareon™ PanOptix™ Toric Trifocal IOL. These IOLs were designed for visual correction of aphakia in adult patients with or without presbyopia and has the following potential benefits:

- Due to the proprietary hydrophobic copolymer formulation used in the Clareon™ PanOptix™ Trifocal IOL and Clareon™ PanOptix™ Toric Trifocal IOL, these IOLs demonstrate low levels of surface haze, subsurface nanoglistenings (SSNGs), and glistenings.
- The proprietary precision edge design of the Clareon™ PanOptix™ Trifocal IOL and Clareon™ PanOptix™ Toric Trifocal IOL is intended to minimize edge glare and maintain a low Nd:YAG rate.
- The Clareon™ PanOptix™ Trifocal IOL and Clareon™ PanOptix™ Toric Trifocal IOL mitigate the effects of presbyopia by providing near, intermediate, and distance vision following crystalline lens removal and IOL implantation.
- Further, the Clareon™ PanOptix™ Trifocal and Clareon™ PanOptix™ Toric Trifocal IOLs also correct ametropia.
- The Clareon™ PanOptix™ Trifocal and Clareon™ PanOptix™ Toric Trifocal IOLs provide consistent satisfactory visual acuity across multiple clinical investigations with parent IOLs.
- Finally, the Clareon™ PanOptix™ Toric Trifocal IOL provides correction of pre-existing corneal astigmatism.

#### CLINICAL STUDY INFORMATION

Several prospective clinical studies with parent IOLs were performed that support the safety and performance of the Clareon™ PanOptix™ Trifocal IOL (Model CNWTT0) and are described in this labeling.

1. A clinical study was conducted in the United States to describe visual outcomes and assess safety of the AcrySof™ IQ PanOptix™ IOL Model TFNT00 following bilateral implantation.
2. A clinical study was conducted in the United States to assess the safety and performance of the AcrySof™ IQ ReSTOR™ +3.0 D Multifocal Toric IOL Models SND1T3 to SND1T6.
3. A clinical study was conducted in the United States to evaluate the safety and performance of the Clareon™ IOL when implanted into the capsular bag following phacoemulsification. IOL rotational stability of the Clareon™ IOL was assessed as part of this study.

Summaries of these clinical studies are provided below to describe the performance characteristics of Clareon™ PanOptix™ Trifocal IOL. Please use caution when comparing these results with results from similar device studies due to potential differences in patient cohorts, test methods, etc. The summary of safety and clinical performance for Clareon™ PanOptix™ Trifocal and Clareon™ PanOptix™ Toric Trifocal IOLs will be made available via Eudamed (<https://ec.europa.eu/tools/eudamed>), which can be located using the model numbers [CNWTT0, CNWTT2-CNWTT6] for Clareon™ PanOptix™ Trifocal IOL per Regulation after Eudamed go-live.

**NOTE:** Data from clinical studies with parent IOLs support the safety and performance of the Clareon™ PanOptix™ Trifocal IOL (Model CNWTT0) and Clareon™ PanOptix™ Toric Trifocal IOL (CNWTT2-CNWTT6):

- The diffractive structure on the anterior side of the Clareon™ PanOptix™ Trifocal IOL and Clareon™ PanOptix™ Toric Trifocal IOL is identical to that of the clinically studied AcrySof™ IQ PanOptix™ IOL.
- The aspheric feature, which is also on the anterior side of the Clareon™ PanOptix™ Trifocal IOL and Clareon™ PanOptix™ Toric Trifocal IOL, is identical to the clinically studied AcrySof™ IQ PanOptix™ IOL.
- The cylinder power on the posterior side of the Clareon™ PanOptix™ Toric Trifocal IOL is the same design as the clinically studied AcrySof™ IQ ReSTOR™ +3.0 D Multifocal Toric IOL.

- The Clareon™ IOL material will not affect the performance of the diffractive or aspheric features as the optic design, mechanical properties, and optical properties are maintained between the Clareon™ PanOptix™ Trifocal IOL and the AcrySof™ IQ PanOptix™ IOL.
- The clinically studied Clareon™ IOL material will not affect the performance of the diffractive or aspheric features as the optic design, mechanical properties, and optical properties are maintained between the Clareon™ PanOptix™ Toric Trifocal IOL and the parent IOLs AcrySof™ IQ PanOptix™ Trifocal IOL and AcrySof™ IQ ReSTOR™ +3.0 D Multifocal Toric IOL.

# 1. CLINICAL STUDY OF AcrySof™ IQ PanOptix™ IOL

## Summary of Clinical Study

The clinical study was a prospective, nonrandomized, vision assessor-masked, parallel-group study and was designed for bilateral implantation of 250 subjects in total, with 125 subjects implanted with the investigational AcrySof™ IQ PanOptix™ IOL Model TFNT00 (referred to as the PanOptix™ IOL below), and 125 subjects implanted with the FDA approved AcrySof™ monofocal IOL Model SN60AT (referred to as the Monofocal IOL below), at 12 investigational sites in the United States.

The purpose of this clinical study was to compare the visual outcomes of the AcrySof™ IQ PanOptix™ IOL Model TFNT00 against that of a monofocal lens, the AcrySof™ Monofocal IOL Model SN60AT, in order to demonstrate comparable distance vision and superior near and intermediate vision.

All eyes with successful IOL implantation in at least one eye and at least one post-operative visit were considered evaluable for the All Implanted analyses. All eyes successfully implanted that had at least one postoperative visit and had no preoperative ocular pathology or macular degeneration at any time, and no major protocol deviations, were evaluable for Best Case analyses. The Best Case data set was the primary data set of analysis for contrast sensitivity and binocular defocus. All eyes with attempted IOL implantation (successful or aborted after contact with the eye) were considered evaluable for the safety analyses.

## Clinical Study Results

Tables presented below covering clinical results from this clinical study use the following conventions. In a column header, "N" is the number in the treatment group. The number of subjects with data "n" are indicated in the table body.

## Subject Population

A total of 243 subjects were implanted in this clinical study with 129 subjects receiving the PanOptix™ IOL and 114 subjects receiving the control monofocal IOL. The study consisted of 67.5% females and 32.5% males. Stratifying by race, the proportions were 86.0% White, 7.8% Black or African American, 3.3% Asian, 0.8% Native Hawaiian or Other Pacific Islander and 2.1% designated "Other". Ethnicity of the study population designated 4.5% as Hispanic or Latino. The mean (± SD) age for the study population was 67 ± 7 years. The Best Case cohort consisted of 129 PanOptix™ IOL subjects and 111 monofocal IOL subjects. Data are reported for the 6-month visit.

## Monocular Visual Acuity

Visual Acuity was assessed using a computerized test system (CTS, M&S Technologies, Niles, IOL). The PanOptix™ IOL met the clinical performance co-primary target (non-inferiority margin of 0.10 logMAR) for mean photopic monocular Best Corrected Distance VA (4 m) for the first operative eye and the co-primary target (superiority) for mean photopic monocular Distance Corrected Visual Acuity (40 cm) for the first operative eye. In addition, the PanOptix™ IOL clinical performance for Distance Corrected Intermediate Visual Acuity (66 cm) for the first operative eye was superior to that of the monofocal IOL.

## Binocular Visual Acuity

There were clinically relevant differences in mean photopic binocular Distance Corrected Visual Acuity (DCVA) at 40 cm and 66 cm for subjects implanted with the PanOptix™ IOL compared with subjects implanted with the control monofocal IOL.

The following is a summary of binocular visual acuity (VA) results for subjects who completed the Form 4A (6 months after second eye implantation) visit. The data are presented in **Tables 6-8** below.

**Table 6: Cumulative Binocular Near (40 cm) Visual Acuity by Lens Model, All Implanted**

|                             |               | N   | Total | 0.04 logMAR [20/20] or better | 0.14 logMAR [20/25] or better | 0.24 logMAR [20/32] or better | 0.34 logMAR [20/40] or better | Worse than 0.34 logMAR [20/40] |
|-----------------------------|---------------|-----|-------|-------------------------------|-------------------------------|-------------------------------|-------------------------------|--------------------------------|
| Uncorrected Photopic        | PanOptix™ IOL | 129 | 127   | 63 (49.6)                     | 117 (92.1)                    | 125 (98.4)                    | 127 (100.0)                   | 0 (0.0)                        |
|                             | Monofocal IOL | 114 | 111   | 1 (0.9)                       | 10 (9.0)                      | 38 (34.2)                     | 67 (60.4)                     | 44 (39.6)                      |
| Distance Corrected Photopic | PanOptix™ IOL | 129 | 127   | 60 (47.2)                     | 122 (96.1)                    | 127 (100.0)                   | 127 (100.0)                   | 0 (0.0)                        |
|                             | Monofocal IOL | 114 | 111   | 0 (0.0)                       | 0 (0.0)                       | 18 (16.2)                     | 44 (39.6)                     | 67 (60.4)                      |
| Distance Corrected Mesopic  | PanOptix™ IOL | 129 | 127   | 5 (3.9)                       | 33 (26.0)                     | 85 (66.9)                     | 119 (93.7)                    | 8 (6.3)                        |
|                             | Monofocal IOL | 114 | 111   | 0 (0.0)                       | 2 (1.8)                       | 6 (5.4)                       | 12 (10.8)                     | 99 (89.2)                      |

Percentage calculated as (n / Total) \* 100

**Table 7: Cumulative Binocular Photopic Intermediate (66 cm) Visual Acuity by Lens Model, All Implanted**

|                    |               | N   | Total | 0.04 logMAR [20/20] or better | 0.14 logMAR [20/25] or better | 0.24 logMAR [20/32] or better | 0.34 logMAR [20/40] or better | Worse than 0.34 logMAR [20/40] |
|--------------------|---------------|-----|-------|-------------------------------|-------------------------------|-------------------------------|-------------------------------|--------------------------------|
| Uncorrected        | PanOptix™ IOL | 129 | 127   | 93 (73.2)                     | 119 (93.7)                    | 124 (97.6)                    | 127 (100.0)                   | 0 (0.0)                        |
|                    | Monofocal IOL | 114 | 111   | 25 (22.5)                     | 56 (50.5)                     | 85 (76.6)                     | 102 (91.9)                    | 9 (8.1)                        |
| Distance Corrected | PanOptix™ IOL | 129 | 127   | 104 (81.9)                    | 124 (97.6)                    | 127 (100.0)                   | 127 (100.0)                   | 0 (0.0)                        |
|                    | Monofocal IOL | 114 | 111   | 6 (5.4)                       | 29 (26.1)                     | 71 (64.0)                     | 92 (82.9)                     | 19 (17.1)                      |

Percentage calculated as (n / Total) \* 100

**Table 8: Cumulative Binocular Photopic Distance (4 m) Visual Acuity by Lens Model, All Implanted**

|                |               | N   | Total | 0.04 logMAR [20/20] or better | 0.14 logMAR [20/25] or better | 0.24 logMAR [20/32] or better | 0.34 logMAR [20/40] or better | Worse than 0.34 logMAR [20/40] |
|----------------|---------------|-----|-------|-------------------------------|-------------------------------|-------------------------------|-------------------------------|--------------------------------|
| Uncorrected    | PanOptix™ IOL | 129 | 127   | 93 (73.2)                     | 117 (92.1)                    | 126 (99.2)                    | 127 (100.0)                   | 0 (0.0)                        |
|                | Monofocal IOL | 114 | 111   | 87 (78.4)                     | 105 (94.6)                    | 110 (99.1)                    | 111 (100.0)                   | 0 (0.0)                        |
| Best Corrected | PanOptix™ IOL | 129 | 127   | 123 (96.9)                    | 127 (100.0)                   | 127 (100.0)                   | 127 (100.0)                   | 0 (0.0)                        |
|                | Monofocal IOL | 114 | 111   | 109 (98.2)                    | 111 (100.0)                   | 111 (100.0)                   | 111 (100.0)                   | 0 (0.0)                        |

Percentage calculated as (n / Total) \* 100

**Table 9** shows the proportion of subjects achieving each Snellen level or better uncorrected binocular visual acuity for all distances (distance – 4 m, intermediate – 66 cm, near – 40 cm). 95.3% of the PanOptix™ IOL subjects achieved 20/32 at all distances.

**Table 9: Proportion of Subjects Achieving VA Thresholds for the Near, Intermediate, and Distance Uncorrected Photopic Binocular Visual Acuity, All Implanted**


| Snellen Category               | PanOptix™ IOL (N = 129) n (%) | Monofocal IOL (N = 114) n (%) |
|--------------------------------|-------------------------------|-------------------------------|
| Total                          | 127                           | 111                           |
| 0.04 logMAR [20/20] or better  | 50 (39.4)                     | 1 (0.9)                       |
| 0.14 logMAR [20/25] or better  | 106 (83.5)                    | 9 (8.1)                       |
| 0.24 logMAR [20/32] or better  | 121 (95.3)                    | 37 (33.3)                     |
| 0.34 logMAR [20/40] or better  | 127 (100.0)                   | 66 (59.5)                     |
| Worse than 0.34 logMAR [20/40] | 0 (0.0)                       | 45 (40.5)                     |

Percentage calculated as (n / Total) \* 100

## Binocular Defocus Curves

Binocular defocus curves were obtained at 6 months for the PanOptix™ IOL and the Monofocal IOL and are shown in **Figure 4**. Data were obtained from best-case subjects in each arm using a computerized visual acuity test system (CTS, M&S Technologies, Niles, IL). The curves display two peaks and one peak respectively that demonstrate the PanOptix™ IOL versus monofocal performance. The main peak, or single peak for the monofocal IOL, is at the zero defocus baseline position, which corresponds to optical infinity. For the PanOptix™ IOL, an additional peak demonstrates the improved performance compared to a monofocal IOL. The PanOptix™ IOL provided mean performance of 0.14 logMAR [20/25] or better vision (depth of focus) from -2.5 D to 0.00 D, corresponding to a range of distances from approximately 40 cm to infinity.

**Figure 4: Mean Binocular Defocus Curves with 95% Confidence Limits by Lens Model, Best Case**

#### Need for Eyeglasses/Contact Lenses

A validated questionnaire consisting of sponsor-developed questions was used in the study to assess need for eyeglass/contact lens following implantation with the IOL. **Table 10** provides the proportions for each response to the question "Overall, in the past 7 days, how often did you need to wear eyeglasses to see?"

**Table 10: Proportion of Subjects Who Respond "Never" to Q1 of the IOLSAT Questionnaire at 6 Months All-Implanted**

|       | PanOptix™ IOL<br>(N = 129)<br>n (%) | Monofocal IOL<br>(N = 114)<br>n (%) | Difference |                |
|-------|-------------------------------------|-------------------------------------|------------|----------------|
|       |                                     |                                     | %          | (95% CI)       |
| Total | 123                                 | 110                                 |            |                |
| Never | 99 (80.5)                           | 9 (8.2)                             | 71.2       | (61.87, 80.46) |

Percentage calculated as (n / Total) \* 100  
Difference = PanOptix™ IOL – Monofocal IOL  
Treatment 1 - Treatment 2 estimate based on Mantel-Haenszel common difference in proportions stratified by site  
CI = Confidence Interval for the common difference, Response scored per user manual

#### Adverse Events

The incidences of cumulative adverse events for the PanOptix™ IOL and the control Monofocal IOL as compared to the ISO 11979-7:2014 historical grid rates are provided in **Tables 11** and **12**. If the same event occurred multiple times in an eye, only the first occurrence is counted in the table below. The rate of secondary surgical interventions (SSIs) did not exceed the SPE grid rate for the PanOptix™ IOL or monofocal IOL group.

**Table 11: Cumulative and Persistent Serious Adverse Events and SPE Rates, First Eye, Safety Analysis Set**

|                                          | PanOptix™ IOL    |                   |                         |          |
|------------------------------------------|------------------|-------------------|-------------------------|----------|
|                                          | (N = 129)<br>n % | 2-sided<br>95% CI | 1-sided 95%<br>Lower CL | SPE<br>% |
| <b>Cumulative Serious Adverse Events</b> |                  |                   |                         |          |
| Cystoid macular edema                    | 0 (0.0)          | (0.00, 2.82)      | 0.00                    | 3.0      |
| Hypopyon                                 | 0 (0.0)          | (0.00, 2.82)      | 0.00                    | 0.3      |
| Endophthalmitis                          | 0 (0.0)          | (0.00, 2.82)      | 0.00                    | 0.1      |
| Lens dislocated from posterior chamber   | 0 (0.0)          | (0.00, 2.82)      | 0.00                    | 0.1      |
| Pupillary block                          | 0 (0.0)          | (0.00, 2.82)      | 0.00                    | 0.1      |
| Retinal detachment                       | 0 (0.0)          | (0.00, 2.82)      | 0.00                    | 0.3      |
| Secondary surgical intervention          | 1 (0.8)          | (0.02, 4.24)      | 0.04                    | 0.8      |
| Other                                    |                  |                   |                         |          |
| Retinal tear                             | 1 (0.8)          | (0.02, 4.24)      | 0.04                    | N/A      |
| <b>Persistent Serious Adverse Events</b> |                  |                   |                         |          |
| Corneal stroma edema                     | 0 (0.0)          | (0.00, 2.82)      | 0.00                    | 0.3      |
| Cystoid macular edema                    | 0 (0.0)          | (0.00, 2.82)      | 0.00                    | 0.5      |
| Iritis                                   | 0 (0.0)          | (0.00, 2.82)      | 0.00                    | 0.3      |
| Raised IOP requiring treatment           | 0 (0.0)          | (0.00, 2.82)      | 0.00                    | 0.4      |

CI = Confidence Interval, CL = Confidence Limit, SPE = Safety and Performance Endpoints  
Persistent = present or ongoing at the final scheduled visit,  
IOP = Intraocular Pressure  
If an eye has multiple occurrences of an AE, the eye is presented only once in the respective eye count column (n) for the corresponding AE. Percentages are calculated as (n/N) \* 100  
The SPE rate is considered not exceeded if the 1-sided 95% lower CL for an AE is less than the SPE%  
"Other" includes the MedDRA Preferred Term for ocular SAEs that do not belong to any predefined SPE categories

The single secondary surgical intervention that occurred with the first eye for PanOptix™ IOL was an explant of the IOL due to subjective complaints of dissatisfaction with the level of vision. This SSI was determined to be related to the optical properties of the IOL.

**Table 12: Cumulative and Persistent Serious Adverse Events and SPE Rates, Second Eye, Safety Analysis Set**

|                                          | PanOptix™ IOL    |                   |                         |          |
|------------------------------------------|------------------|-------------------|-------------------------|----------|
|                                          | (N = 127)<br>n % | 2-sided<br>95% CI | 1-sided 95%<br>Lower CL | SPE<br>% |
| <b>Cumulative Serious Adverse Events</b> |                  |                   |                         |          |
| Cystoid macular edema                    | 1 (0.8)          | (0.02, 4.31)      | 0.04                    | 3.0      |
| Hypopyon                                 | 0 (0.0)          | (0.00, 2.86)      | 0.00                    | 0.3      |
| Endophthalmitis                          | 0 (0.0)          | (0.00, 2.86)      | 0.00                    | 0.1      |
| Lens dislocated from posterior chamber   | 0 (0.0)          | (0.00, 2.86)      | 0.00                    | 0.1      |
| Pupillary block                          | 0 (0.0)          | (0.00, 2.86)      | 0.00                    | 0.1      |
| Retinal detachment                       | 0 (0.0)          | (0.00, 2.86)      | 0.00                    | 0.3      |
| Secondary surgical intervention          | 2 (1.6)          | (0.19, 5.57)      | 0.28                    | 0.8      |
| Other                                    |                  |                   |                         |          |
| Device dislocation                       | 1 (0.8)          | (0.02, 4.31)      | 0.04                    | N/A      |
| Vitreous prolapse                        | 1 (0.8)          | (0.02, 4.31)      | 0.04                    | N/A      |
| <b>Persistent Serious Adverse Events</b> |                  |                   |                         |          |
| Corneal stroma edema                     | 0 (0.0)          | (0.00, 2.86)      | 0.00                    | 0.3      |
| Cystoid macular edema                    | 0 (0.0)          | (0.00, 2.86)      | 0.00                    | 0.5      |
| Iritis                                   | 0 (0.0)          | (0.00, 2.86)      | 0.00                    | 0.3      |
| Raised IOP requiring treatment           | 0 (0.0)          | (0.00, 2.86)      | 0.00                    | 0.4      |

CI = Confidence Interval, CL = Confidence Limit, SPE = Safety and Performance Endpoints  
Persistent = present or ongoing at the final scheduled visit,  
IOP = Intraocular Pressure  
If an eye has multiple occurrences of an AE, the eye is presented only once in the respective eye count column (n) for the corresponding AE. Percentages are calculated as (n/N) \* 100  
The SPE rate is considered not exceeded if the 1-sided 95% lower CL for an AE is less than the SPE%  
"Other" includes the MedDRA Preferred Term for ocular SAEs that do not belong to any predefined SPE categories

The first secondary surgical intervention that occurred with the second eye for PanOptix™ IOL was a vitrectomy performed due to a vitreous prolapse. The second secondary surgical intervention that occurred with the second eye for PanOptix™ IOL was a lens repositioning procedure due to a tilted/displaced IOL. These SSIs occurred in different subjects and neither were determined to be related to the optical properties of the IOL.

#### Contrast Sensitivity

Binocular best corrected distance contrast sensitivity was performed using a backlit sine wave grating chart system (CSV1000, VectorVision, Greenville, OH) at 6 months under four conditions: photopic without glare, photopic with glare, mesopic without glare, and mesopic with glare. The mean contrast sensitivity was similar between the 2 groups, across all tested conditions, and irrespective of glare source, with an observed maximum difference in means of  $\leq 0.11$  log units between the 2 groups, falling below the level of clinical significance from ISO11979-7: 2014.

#### Visual Disturbances

A validated Patient Reported Outcome Measure instrument was developed and used in this clinical study to assess visual disturbances. Subjects were first asked if they experienced a particular visual disturbance. If the subject responded affirmatively, he or she was asked to rate the severity, frequency, and bothersomeness. As demonstrated in **Table 13**, reports of visual disturbances were similar between the PanOptix™ IOL and the control Monofocal IOL groups at 6 months. The highest rate of most bothersome reports ("Bothered Very Much") of visual disturbances/distortions at 6 months was for starbursts at 4.8% for the PanOptix™ IOL and 0.9% for the control monofocal IOL.

**Table 13: Visual Disturbance Bothersomeness, Safety Analysis Set**

|                    | PanOptix™ IOL<br>N=129 |                                             |                |             |               |             | Monofocal IOL<br>N=114 |                                             |                |             |               |             |
|--------------------|------------------------|---------------------------------------------|----------------|-------------|---------------|-------------|------------------------|---------------------------------------------|----------------|-------------|---------------|-------------|
| Visual Disturbance | n                      | Bothered                                    |                |             |               |             | n                      | Bothered                                    |                |             |               |             |
|                    |                        | Not experienced or<br>Not bothered at all % | A Little bit % | Some-what % | Quite a bit % | Very much % |                        | Not experienced or<br>Not bothered at all % | A Little bit % | Some-what % | Quite a bit % | Very much % |
|                    |                        |                                             |                |             |               |             |                        |                                             |                |             |               |             |
| Glare              | 126                    | 54.8                                        | 18.3           | 18.3        | 7.1           | 1.6         | 111                    | 69.4                                        | 15.3           | 8.1         | 6.3           | 0.9         |
| Halos              | 127                    | 51.2                                        | 21.3           | 16.5        | 8.7           | 2.4         | 110                    | 83.6                                        | 10.9           | 3.6         | 0.9           | 0.9         |
| Starbursts         | 125                    | 55.2                                        | 16.8           | 16.0        | 7.2           | 4.8         | 109                    | 79.8                                        | 10.1           | 8.3         | 0.9           | 0.9         |
| Hazy vision        | 125                    | 86.4                                        | 6.4            | 6.4         | 0.8           | 0.0         | 110                    | 89.1                                        | 5.5            | 3.6         | 0.9           | 0.9         |
| Blurred vision     | 127                    | 81.1                                        | 10.2           | 6.3         | 2.4           | 0.0         | 111                    | 86.5                                        | 4.5            | 3.6         | 3.6           | 1.8         |
| Double vision      | 125                    | 96.0                                        | 2.4            | 1.6         | 0.0           | 0.0         | 110                    | 98.2                                        | 0.0            | 1.8         | 0.0           | 0.0         |
| Dark Area*         | 127                    | 89.8                                        | 7.1            | 3.1         | 0.0           | 0.0         | 111                    | 92.8                                        | 3.6            | 2.7         | 0.9           | 0.0         |

Percentage calculated as (n / N) \* 100  
\*Dark Area corresponds to negative dysphotopsia

#### Patient Satisfaction

A validated questionnaire consisting of sponsor-developed questions was used in the study to evaluate subject satisfaction following implantation with the IOL. **Table 14** provides the results.

**Table 14: IOLSAT: Satisfaction with Your Vision (Collected at 6 Months), All-Implanted**

| Question                                                                                               | Response                           | PanOptix™ IOL<br>(N = 129) | Monofocal IOL<br>(N = 114) |
|--------------------------------------------------------------------------------------------------------|------------------------------------|----------------------------|----------------------------|
|                                                                                                        |                                    | n (%)                      | n (%)                      |
| In the past 7 days, how satisfied were you with your vision?                                           | Total                              | 127                        | 110                        |
|                                                                                                        | Very Dissatisfied                  | 2 (1.6)                    | 0 (0.0)                    |
|                                                                                                        | Dissatisfied                       | 2 (1.6)                    | 3 (2.7)                    |
|                                                                                                        | Neither Satisfied nor Dissatisfied | 2 (1.6)                    | 7 (6.4)                    |
|                                                                                                        | Satisfied                          | 27 (21.3)                  | 34 (30.9)                  |
|                                                                                                        | Very Satisfied                     | 94 (74.0)                  | 66 (60.0)                  |
| Given your vision today, if you had to do it all over, would you have the same lenses implanted again? | Total                              | 127                        | 111                        |
|                                                                                                        | No                                 | 1 (0.8)                    | 14 (12.6)                  |
|                                                                                                        | Yes                                | 126 (99.2)                 | 97 (87.4)                  |
| Given your vision today, would you recommend the lenses you had implanted to your family or friends?   | Total                              | 127                        | 110                        |
|                                                                                                        | No                                 | 2 (1.6)                    | 5 (4.5)                    |
|                                                                                                        | Yes                                | 125 (98.4)                 | 105 (95.5)                 |
| Percentage calculated as (n / Total) * 100                                                             |                                    |                            |                            |

## 2. CLINICAL STUDY with AcrySof™ IQ ReSTOR™ +3.0 D MULTIFOCAL TORIC IOL (Models SND1T2- SND1T6)

### Summary of Clinical Study

Data relevant to Clareon™ PanOptix™ Toric Trifocal IOL from a clinical study of the AcrySof™ IQ ReSTOR™ +3.0 D Multifocal Toric IOL Models SND1T3-SND1T6 (referred to as ReSTOR™ Toric +3.0 D IOL below) in comparison to AcrySof™ ReSTOR™ (+4.0 D Add) Multifocal IOL Model SA60D3 are described below. The clinical study was a prospective, nonrandomized, unmasked, parallel-group study with a follow up of 12 months (defined as 330 to 420 days postoperative after second eye implant). The investigational ReSTOR™ Toric +3.0 D IOL was designed with a near reading distance of 40 cm, and the control ReSTOR™ +4.0 D IOL was designed with a near reading distance of 33 cm. All eyes with successful IOL implantation in at least one eye were considered evaluable for the All Implanted analyses.

### Subject Population

A total of 574 subjects were bilaterally implanted with 386 subjects receiving the ReSTOR™ Toric +3.0 and 188 subjects receiving the control ReSTOR™ +4.0 D IOL. The study consisted of 65.5% females and 34.5% males. Stratifying by race, there were 93.7% White, 4.5% Black or African American, 0.9% Asian and 0.9% designated "Other". The mean age for the study population was 67 ± 9 years.

## RESULTS RELEVANT TO THE CLAREON™ PanOptix™ TORIC TRIFOCAL IOL

### Orientation of Lens Axis

The results for lens axis orientation at all postoperative visits were compared to those at surgery to determine lens axis rotation. The difference between the achieved lens axis orientation at month 12 and the achieved axis placement at surgery was 2.7° ± 5.8 in the first operative eyes and 2.2° ± 2.7 in the second operative eyes (Table 15). Eight subjects had lens axis rotation of twenty degrees or more at month 12, two of whom had incorrect lens axis orientation measurements and three of whom underwent lens repositioning and have improved outcomes with the lens implanted (post repositioning rotation was less than 6 degrees). All eight subjects had improved visual performance at month 12.

**Table 15: Descriptive Statistics for the Absolute Difference Between Lens Axis Orientation at the Post-operative Visit and Achieved Axis Placement (Degrees) at the Operative Visit (All Implanted Set)**

|                                                                                                           |            | Absolute Rotation   |                      |
|-----------------------------------------------------------------------------------------------------------|------------|---------------------|----------------------|
|                                                                                                           |            | First Implanted Eye | Second Implanted Eye |
| 1 month                                                                                                   | n          | 367                 | 368                  |
|                                                                                                           | Mean (SD)  | 2.2 (5.1)           | 2.1 (2.7)            |
|                                                                                                           | (Min, Max) | (0, 85)             | (0, 24)              |
|                                                                                                           | 95% CI     | (1.6, 2.7)          | (1.8, 2.4)           |
| 6 months                                                                                                  | n          | 363                 | 364                  |
|                                                                                                           | Mean (SD)  | 2.3 (5.2)           | 2.3 (3.0)            |
|                                                                                                           | (Min, Max) | (0, 85)             | (0, 27)              |
|                                                                                                           | 95% CI     | (1.7, 2.8)          | (2.0, 2.6)           |
| 12 months                                                                                                 | n          | 356                 | 357                  |
|                                                                                                           | Mean (SD)  | 2.7 (5.8)           | 2.2 (2.7)            |
|                                                                                                           | (Min, Max) | (0, 84)             | (0, 24)              |
|                                                                                                           | 95% CI     | (2.1, 3.3)          | (1.9, 2.5)           |
| For subjects with missing Operative Visit axis placement data, Day 1 (Visit 1) data were used as baseline |            |                     |                      |

Furthermore, the rotational stability of the ReSTOR™ Toric +3.0 D IOL was maintained between 2 consecutive visits at least 3 months apart (between 1 month and 6 months). As recommended by the 2010 ANSI standard for toric intraocular lenses, the data demonstrate that at least 90% of ReSTOR™ Toric +3.0 D IOL subjects achieved a rotational stability of 5 degrees or less between 2 consecutive visits, at least 3 months apart (Table 16).

**Table 16: Number and Percentage of Subjects by Lens Axis Rotation Between 1 Month and 6 Months (All Implanted Set)**

|                      |                           | ReSTOR™ Toric +3.0 D IOL |        |
|----------------------|---------------------------|--------------------------|--------|
|                      |                           | n                        | (%)    |
| First Implanted Eye  | Total                     | 359                      |        |
|                      | Lens Movement ≤ 5 degrees | 338                      | (94.2) |
|                      | Lens Movement >5 degrees  | 21                       | (5.8)  |
| Second Implanted Eye | Total                     | 361                      |        |
|                      | Lens Movement ≤ 5 degrees | 339                      | (93.9) |
|                      | Lens Movement >5 degrees  | 22                       | (6.1)  |

ReSTOR™ Toric +3.0 D IOL = AcrySof™ IQ ReSTOR™ +3.0 D Multifocal Toric IOL Models SND1T3/ SND1T4/SND1T5/SND1T6

Subjects with missing observations at either 1 month or 6 months were excluded

### Reduction of Cylinder

The ReSTOR™ Toric +3.0 D IOLs are effective in the reduction of corneal astigmatism in the range of 0.75 D to 2.82 D. As demonstrated in Table 17, the percent reduction in cylinder with respect to target cylinder was calculated and descriptive statistics were computed at the 12 month postoperative visit. Target cylinder was defined as the amount of anticipated residual astigmatism as calculated by the study specific web-based IOL Clinical Calculator.

**Table 17: Number and Percentage of Subjects With Reduction of Cylinder Within the Target Cylinder at 12 months for ReSTOR™ Toric +3.0 D IOL (All Implanted)**

|             | First Implanted Eye<br>(N = 373) |        | Second Implanted Eye<br>(N = 371) |        |
|-------------|----------------------------------|--------|-----------------------------------|--------|
|             | n                                | (%)    | n                                 | (%)    |
| Within 0.5D | 278                              | (74.5) | 295                               | (79.5) |
| Within 1.0D | 351                              | (94.1) | 362                               | (97.6) |
| > 1.0D      | 22                               | (5.9)  | 9                                 | (2.4)  |

## 3. CLINICAL STUDY OVERVIEW FOR CLAREON™ ASPHERIC HYDROPHOBIC ACRYLIC IOL

A prospective, multicenter, open-label clinical study in adult subjects requiring cataract surgery with IOL implantation was performed on the Clareon™ IOL in the United States and completed in 2019. The objective of this study was to demonstrate favorable visual acuity and adverse event outcomes for the Clareon™ IOL compared to historical safety and performance endpoint (SPE) rates as reported in EN ISO 11979-7:2014. The results achieved by the subjects successfully followed for 12 Months postoperatively (defined as 330 to 420 days postoperative) provide reasonable assurance that the Clareon™ IOL is a safe and effective device for the visual correction of aphakia following cataract surgery.

The primary performance and safety endpoints for the clinical study are provided below to describe the performance characteristics of the IOL. Please use caution when comparing these results with results from similar device studies due to potential differences in patient cohorts, test methods, etc.

### STUDY RESULTS FOR CLAREON™ ASPHERIC HYDROPHOBIC ACRYLIC IOL

The subject population in the clinical study consisted of 350 subjects in need of cataract surgery that were unilaterally implanted, the percentage of females (60.9%) exceeded the percentage of males (39.1%) enrolled. The mean age for the total population was 69.7 years.

### Visual Acuity

The best corrected distance visual acuity (BCDVA) achieved at one year postoperatively by the Best Case Analysis Set is presented in Table 18 and visual acuity achieved by All-Implanted Analysis Set is presented in Table 19. Similar data for BCDVA expressed as Snellen values are presented in Tables 20 and 21. The Best Case analysis included all eyes successfully implanted with the test article, had at least one postoperative visit, did not have any ocular pathology at Visit 0 (preoperative), had no macular degeneration at any visit, and did not have any previous surgery for the correction of refractive errors. All-Implanted Analysis set includes all eyes with successful test article implantation.

The one-sided exact 95% upper confidence limit for the percentage of subjects with monocular BCDVA of 0.3 logMAR or better at 12 months postoperative (Visit 5) is greater than or equal to the SPE rates as reported in EN ISO 11979-7:2014.

**Table 18: Best Corrected Visual Acuity in the Best Case Analysis Set at Month 12, Clareon™ IOL**

| Visual Acuity                                                                                                                                                                                            | n     | %    |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------|------|
| 0.0 logMAR or better                                                                                                                                                                                     | 271   | 83.1 |
| 0.1 logMAR or better                                                                                                                                                                                     | 307   | 94.2 |
| 0.2 logMAR or better                                                                                                                                                                                     | 323   | 99.1 |
| 0.3 logMAR or better                                                                                                                                                                                     | 325   | 99.7 |
| Worse than 0.3 logMAR                                                                                                                                                                                    | 1     | 0.3  |
| ISO SPE Rate for % of 0.3 logMAR or better                                                                                                                                                               | 96.7% |      |
| N                                                                                                                                                                                                        | 326   |      |
| ISO SPE Rate = BS EN ISO 11979-7:2014, Ophthalmic Implants - Intraocular Lenses - Part 7: Clinical Investigations, Table B.4 - Posterior Chamber IOL Adverse Event Safety and Performance Endpoint Rates |       |      |

**Table 19: Best Corrected Visual Acuity in the All-Implanted Analysis Set at Month 12, Clareon™ IOL**

| Visual Acuity                                                                                                                                                                                            | n     | %    |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------|------|
| 0.0 logMAR or better                                                                                                                                                                                     | 277   | 81.0 |
| 0.1 logMAR or better                                                                                                                                                                                     | 318   | 93.0 |
| 0.2 logMAR or better                                                                                                                                                                                     | 339   | 99.1 |
| 0.3 logMAR or better                                                                                                                                                                                     | 341   | 99.7 |
| Worse than 0.3 logMAR                                                                                                                                                                                    | 1     | 0.3  |
| ISO SPE Rate for % of 0.3 logMAR or better                                                                                                                                                               | 92.5% |      |
| N                                                                                                                                                                                                        | 342   |      |
| ISO SPE Rate = BS EN ISO 11979-7:2014, Ophthalmic Implants - Intraocular Lenses - Part 7: Clinical Investigations, Table B.3 - Posterior Chamber IOL Adverse Event Safety and Performance Endpoint Rates |       |      |

**Table 20: Best Corrected Visual Acuity in the Best Case Analysis Set at Month 12 (Snellen), Clareon™ IOL**

| Visual Acuity             | n                   | %                   |
|---------------------------|---------------------|---------------------|
| 20/20 or better           | 287                 | 88.0                |
| 20/25 or better           | 316                 | 96.9                |
| 20/32 or better           | 324                 | 99.4                |
| 20/40 or better           | 325                 | 99.7                |
| Worse than 20/40          | 1                   | 0.3                 |
| N                         | 326                 |                     |
| Note: 20/20 = 0.04 logMAR | 20/25 = 0.14 logMAR | 20/32 = 0.24 logMAR |
| 20/40 = 0.34 logMAR       |                     |                     |

**Table 21: Best Corrected Visual Acuity in the All-Implanted Analysis Set at Month 12 (Snellen), Clareon™ IOL**

| Visual Acuity             | n                   | %                   |
|---------------------------|---------------------|---------------------|
| 20/20 or better           | 297                 | 86.8                |
| 20/25 or better           | 331                 | 96.8                |
| 20/32 or better           | 340                 | 99.4                |
| 20/40 or better           | 341                 | 99.7                |
| Worse than 20/40          | 1                   | 0.3                 |
| N                         | 342                 |                     |
| Note: 20/20 = 0.04 logMAR | 20/25 = 0.14 logMAR | 20/32 = 0.24 logMAR |
| 20/40 = 0.34 logMAR       |                     |                     |

#### IOL Rotation

The rotational stability of the Clareon™ IOL was evaluated in a subset of the total study population, which included 141 subjects at 6 sites. Results were compared to performance targets in EN ISO 11979-7:2014, defined as the difference between IOL axis orientation on day of surgery and at the Month 6 visit. A total of 98.4% of subjects had IOL rotation <10 degrees and 99.2% had rotation <20 degrees or <30 degrees at Month 6 compared to the day of surgery, which met the pre-specified performance targets. Results shown in **Table 22**.

**Table 22: Categorical Statistics for Absolute IOL Rotation (10 Degrees Increment) at Month 6 from Day of Surgery, Rotation Analysis Set, Clareon™ IOL**

|                                                                                                                                                                         | (N=141) |      |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------|------|
| Category                                                                                                                                                                | n       | %    |
| Total                                                                                                                                                                   | 124     | -    |
| Less than 10 degrees                                                                                                                                                    | 122     | 98.4 |
| Less than 20 degrees                                                                                                                                                    | 123     | 99.2 |
| Less than 30 degrees                                                                                                                                                    | 123     | 99.2 |
| Greater than 30 degrees                                                                                                                                                 | 1       | 0.8  |
| N = Number of eyes in the analysis set<br>n = Number of eyes in specified category<br>Total = Number of eyes with data<br>Percentages are calculated as (n/Total) * 100 |         |      |

In an additional analysis between Month 1 and Month 6, 98.4% of subjects implanted with the Clareon™ IOL had an absolute rotation of less than or equal to 5 degrees. The remaining subjects (1.6%) had rotation between 5 to 10 degrees.

Absolute rotation of Clareon™ IOL was also described at all post-operative visits (1 Day, 1 Week, 1 Month, and 6 Months) compared to the day of surgery. Descriptive statistics results at each postoperative visit are shown in **Table 23**.

**Table 23: Descriptive Statistics for Absolute IOL Rotation (Degrees), Rotation Analysis Set, Clareon™ IOL**

| Visit                                                                                                                                                                                                                                                                                             | Statistic  | (N = 141)*       |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------|------------------|
| 1 Day                                                                                                                                                                                                                                                                                             | n          | 127              |
|                                                                                                                                                                                                                                                                                                   | Mean (SD)  | 1.8483 (3.72254) |
|                                                                                                                                                                                                                                                                                                   | Median     | 0.942            |
|                                                                                                                                                                                                                                                                                                   | (Min, Max) | (0.002, 38.232)  |
| 1 Week                                                                                                                                                                                                                                                                                            | 95% CI     | (1.1946, 2.5020) |
|                                                                                                                                                                                                                                                                                                   | n          | 126              |
|                                                                                                                                                                                                                                                                                                   | Mean (SD)  | 1.9942 (3.86948) |
|                                                                                                                                                                                                                                                                                                   | Median     | 1.213            |
| 1 Month                                                                                                                                                                                                                                                                                           | (Min, Max) | (0.020, 40.547)  |
|                                                                                                                                                                                                                                                                                                   | 95% CI     | (1.3119, 2.6764) |
|                                                                                                                                                                                                                                                                                                   | n          | 127              |
|                                                                                                                                                                                                                                                                                                   | Mean (SD)  | 2.2335 (3.95863) |
| 6 Months                                                                                                                                                                                                                                                                                          | Median     | 1.414            |
|                                                                                                                                                                                                                                                                                                   | (Min, Max) | (0.008, 41.110)  |
|                                                                                                                                                                                                                                                                                                   | 95% CI     | (1.5384, 2.9287) |
|                                                                                                                                                                                                                                                                                                   | n          | 124              |
|                                                                                                                                                                                                                                                                                                   | Mean (SD)  | 2.2696 (3.87231) |
|                                                                                                                                                                                                                                                                                                   | Median     | 1.449            |
|                                                                                                                                                                                                                                                                                                   | (Min, Max) | (0.014, 40.033)  |
|                                                                                                                                                                                                                                                                                                   | 95% CI     | (1.5812, 2.9579) |
| N = Number of eyes in the analysis set<br>n = Number of eyes at visit with data<br>SD = Standard Deviation, CI = Confidence Interval<br>* Of 141, 10 had images at baseline which were not readable and at each visit, individual subjects may have had images which were missing or not readable |            |                  |

#### Serious Adverse Events

The cumulative and persistent adverse events rates for the Safety Analysis Set are presented in **Table 19**. Cumulative adverse events are those that occur at any point in the study while persistent adverse events are those that occur or are ongoing at the 12 Month Visit. The Safety Analysis Set included all eyes with attempted implantation with the test article (successful or aborted after contact with the eye). In this study, the number of eyes in the Safety Analysis Set was equivalent to the All-Implanted Analysis Set.

For both cumulative and persistent adverse events, the exact one-sided 95% lower confidence limit of cumulative and persistent adverse events were within limits compared to the SPE rates reported in EN ISO 11979-7:2014.

**Table 24: Cumulative and Persistent Adverse Events in the Safety Analysis Set, Clareon™ IOL**

|                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  | (N = 350)<br>n (%) | 2-sided<br>95% CI | 1-sided 95%<br>Lower CL | SPE<br>% |
|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------|-------------------|-------------------------|----------|
| <b>Cumulative Serious Adverse Events</b>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |                    |                   |                         |          |
| Cystoid macular edema                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            | 3 (0.9)            | (0.18, 2.48)      | 0.23                    | 3.0      |
| Hypopyon                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         | 0 (0.0)            | (0.00, 1.05)      | 0.00                    | 0.3      |
| Endophthalmitis                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  | 0 (0.0)            | (0.00, 1.05)      | 0.00                    | 0.1      |
| Lens dislocated from posterior chamber                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           | 0 (0.0)            | (0.00, 1.05)      | 0.00                    | 0.1      |
| Pupillary block                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  | 0 (0.0)            | (0.00, 1.05)      | 0.00                    | 0.1      |
| Retinal detachment                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               | 0 (0.0)            | (0.00, 1.05)      | 0.00                    | 0.3      |
| Secondary surgical intervention                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  | 6 (1.7)            | (0.63, 3.69)      | 0.75                    | 0.8      |
| Other                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |                    |                   |                         |          |
| Herpes virus infection                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           | 1 (0.3)            | (0.01, 1.58)      | 0.01                    | NA       |
| Macular fibrosis                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 | 1 (0.3)            | (0.01, 1.58)      | 0.01                    | NA       |
| Macular hole                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     | 1 (0.3)            | (0.01, 1.58)      | 0.01                    | NA       |
| Punctate keratitis                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               | 1 (0.3)            | (0.01, 1.58)      | 0.01                    | NA       |
| Retinal tear                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     | 2 (0.6)            | (0.07, 2.05)      | 0.10                    | NA       |
| <b>Persistent Serious Adverse Events</b>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |                    |                   |                         |          |
| Corneal stroma edema                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             | 0 (0.0)            | (0.00, 1.05)      | 0.00                    | 0.3      |
| Cystoid macular edema                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            | 1 (0.3)            | (0.01, 1.58)      | 0.01                    | 0.5      |
| Iritis                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           | 0 (0.0)            | (0.00, 1.05)      | 0.00                    | 0.3      |
| Raised IOP requiring treatment                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   | 0 (0.0)            | (0.00, 1.05)      | 0.00                    | 0.4      |
| CI = Confidence Interval<br>CL = Confidence Limit<br>SPE = Safety and Performance Endpoints<br>If an eye has multiple occurrences of an AE, the eye is presented only once in the respective eye count column (n) for the corresponding AE. Percentages are calculated as (n/N) * 100.<br>The SPE rate is considered not exceeded if the 1-sided 95% lower CL for an AE is less than the SPE%.<br>"Other" includes the MedDRA Preferred Term for ocular SAEs that do not belong to any predefined SPE categories.<br>Persistent = Present or ongoing at the final scheduled visit.<br>IOP = Intraocular Pressure |                    |                   |                         |          |

#### IOL Observations

At all five scheduled postoperative visits and at unscheduled visits, slit-lamp examinations were performed. The list of pre-specified slit-lamp findings included IOL observations, which were described as IOL glistenings, scratches/cracks on the IOL, and surface haze on the IOL. Across 1,852 postoperative (including unscheduled) visits, no IOL observations (IOL glistenings, scratches/cracks on the IOL, and surface haze on the IOL) were noted for the Clareon™ IOL.

## REFERENCES

Boettner, E.A. and Wolter, J.R. Transmission of the ocular media. *Invest. Ophthalmol.* 1962;1:776-83.  
 Haigis W. ULIB - User group for Laser Interference Biometry. Updated 2014 Jul 17; cited 2015 Mar 13.  
 Available at: <http://ocusoft.de/ulib/>.  
 Hoffer KJ. The Hoffer Q formula: a comparison of theoretic and regression formulas. *J Cataract Refract Surg.* 1993;19(6):700-12.  
 Holladay JT. Standardizing constants for ultrasonic biometry, keratometry, and intraocular lens power calculations. *J Cataract Refract Surg.* 1997;23(9):1356-70.  
 Olsen T. Calculation of intraocular lens power: a review. *Acta Ophthalmol Scand.* 2007;85(5):472-85.  
 Retzlaff JA, Sanders DR, Kraff M. *Lens Implant Power Calculation*. 3rd ed. Thorofare (NJ): Slack, Inc.; 1990.

### STAND ALONE SYMBOLS FROM ISO 7000/ISO 7001\* USED ON LABELING (ISO 7000/IEC 60417 Title: Graphical Symbols for Use on Equipment) (\*ISO 7001 Title: Graphical Symbols – Public Information Symbols)

| Symbol | Reference Number from ISO 7000 / IEC 60417* / ISO 7001* | Symbol Title / Explanatory Text                                         |
|--------|---------------------------------------------------------|-------------------------------------------------------------------------|
|        | 1051                                                    | Do not re-use                                                           |
|        | 2608                                                    | Do not re-sterilize                                                     |
|        | 2607                                                    | Use-by date                                                             |
|        | 2501                                                    | Sterilized using ethylene oxide                                         |
|        | 2498                                                    | Serial number                                                           |
|        | 2493                                                    | Catalogue number                                                        |
|        | 0434A                                                   | Caution                                                                 |
|        | 3082                                                    | Manufacturer                                                            |
|        | 0533                                                    | Upper limit of temperature                                              |
|        | 1641                                                    | Consult instructions for use or consult electronic instructions for use |
|        | 2606                                                    | Do not use if package has been damaged and consult instructions for use |
|        | 2497                                                    | Date of Manufacture                                                     |
|        | 3079                                                    | Open Here                                                               |
|        | 3010                                                    | RFID tag, general                                                       |
|        | 6050*                                                   | Model number                                                            |
|        | 6049*                                                   | Country of manufacture                                                  |
|        | 5662*                                                   | Date                                                                    |
|        | 5664*                                                   | Patient identification                                                  |
|        | PI PF 044*                                              | Health care centre or doctor                                            |
|        | PI PF 002'                                              | Hospital                                                                |
|        | 3707                                                    | Single sterile barrier system                                           |

\*This symbol is from ISO 7001.

\*This symbol is from IEC 60417.

### STAND ALONE SYMBOL FROM ASTM F2503-13 USED ON LABELING (Title: Standard Practice for Marking Medical Devices and Other Items for Safety in the Magnetic Resonance Environment)

| Symbol | Reference Number from ASTM F2503-13 | Symbol Title / Explanatory Text |
|--------|-------------------------------------|---------------------------------|
|        | N/A                                 | MR (magnetic resonance) safe    |

## ABBREVIATIONS or OTHER SYMBOLS USED ON LABELING

| Abbreviation or Symbol | Symbol Title / Explanatory Text                                                            |
|------------------------|--------------------------------------------------------------------------------------------|
|                        | Medical device                                                                             |
|                        | Patient information website                                                                |
|                        | Unique device identifier                                                                   |
| IOL                    | Intraocular lens                                                                           |
| OVD                    | Ophthalmic viscosurgical device                                                            |
|                        | UV and Blue Light Filter                                                                   |
|                        | Eye                                                                                        |
|                        | Posterior chamber IOL                                                                      |
| UV                     | Ultraviolet                                                                                |
| D                      | Diopter                                                                                    |
| Ø <sub>B</sub>         | Body diameter (Optic diameter)                                                             |
| Ø <sub>T</sub>         | Overall diameter (Overall length)                                                          |
| L                      | Left                                                                                       |
| R                      | Right                                                                                      |
| ADD                    | Add power                                                                                  |
| PWR                    | Spherical Equivalent Power                                                                 |
|                        | D-size nozzle of MONARCH™ cartridge*                                                       |
|                        | C-size nozzle of MONARCH™ cartridge*                                                       |
|                        | B-size nozzle of MONARCH™ cartridge*                                                       |
|                        | CLAREON™ MONARCH™ IV                                                                       |
|                        | Not made with natural rubber latex                                                         |
|                        | Does not contain PHT (phthalates)                                                          |
|                        | Caution: Federal (USA) Law restricts this device to sale by or on the order of a physician |
|                        | Authorized Representative in the European Community                                        |

\*The recommendation shown on the labeling is for the smallest qualified cartridge nozzle size per diopter. For all qualified combinations of products, please review **Table 2**.

Alcon Laboratories, Inc.  
 6201 South Freeway  
 Fort Worth, Texas 76134-2099 USA

† Trademarks are the property of their respective owners.

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© 2022 Alcon Inc.