Instruction for use

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«HSV 1,2 IgM-ИФА»

**A SOLID-PHASE ENZYME IMMUNOASSAY
FOR THE QUALITATIVE DETERMINATION
OF IgM ANTIBODIES TO HSV 1/2
IN HUMAN SERUM OR PLASMA**

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
Authorized Representative in EU:
Polmed.de
Steinacker 20, D-73773
Aichwald, Germany
e-mail: info@polmed.de

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CONTENT

1. INTENDED USE	11
2. SUMMARY AND EXPLANATION	11
3. PRINCIPLE OF THE TEST	11
4. WARNINGS AND PRECAUTIONS	12
5. KIT COMPONENTS	13
6. SPECIMEN COLLECTION AND STORAGE	14
7. TEST PROCEDURE	14
8. QUALITY CONTROL	15
9. CALCULATION OF RESULTS	15
10. EXPECTED VALUES	15
11. PERFORMANCE CHARACTERISTICS	16
12. LITERATURE	16

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Instruction for use

A SOLID-PHASE ENZYME IMMUNOASSAY FOR THE QUALITATIVE DETERMINATION OF IGM ANTIBODIES TO HSV 1/2 IN HUMAN SERUM OR PLASMA

1. INTENDED USE

A solid-phase enzyme immunoassay for the qualitative determination of IgM antibodies to HSV 1/2 in serum or plasma.

This kit is designed for measurement of IgM antibodies to HSV 1/2 in serum or plasma. For possibility of use with other sample types, please, refer to Application Notes (on request). The kit contains reagents sufficient for 96 determinations and allows to analyze 46 unknown samples in duplicates.

2. SUMMARY AND EXPLANATION

Herpes simplex virus (HSV)infection is a chronic recurrent infectious disease caused by Herpes simplex viruses, types 1 and 2. The disease may be manifested as a local skin and/or mucosal vesicular rash or as a generalized form with multi-organ damage.

During HSV infection, a sequential synthesis of HSV-specific IgM, IgG and IgA occurs. IgM antibodies appear 10 or more days following infection, their level continuously rising during two weeks. HSV-specific IgG appear 7–10 days following appearance of IgM and persist lifelong.

A primary HSV infection is diagnosed by appearance of HSV-specific IgM and/or a 4-fold increase of HSV-specific IgG in sequential samples taken with 14-20 days interval. Appearance of HSV-specific IgM in an HSV-infected person indicates exacerbation of the disease.

3. PRINCIPLE OF THE TEST

This test is based on two-site sandwich enzyme immunoassay principle. Tested specimen is placed into the microwells coated by murine monoclonal antibodies to IgM. Antibodies from the specimen bind coated murine monoclonal antibodies to IgM on the microwell surface. Unbound material is removed by washing procedure. HSV 1,2, labelled with peroxidase enzyme, is then added into the microwells. After subsequent washing procedure, the remaining enzymatic activity bound to the microwell surface is detected and quantified by addition of chromogen-substrate mixture, stop solution and photometry at 450 nm. Positivity index (PI, %) is calculated by the formula (see Calculations).Optical density in the microwell is directly related to the quantity of the measured analyte in the specimen.

4. WARNINGS AND PRECAUTIONS

4.1. For professional use only.

4.2. This kit is intended for in vitro diagnostic use only.

4.3. INFECTION HAZARD: There is no available test methods that can absolutely assure that Hepatitis B and C viruses, HIV-1/2, or other infectious agents are not present in the reagents of this kit. All human products, including patient samples, should be considered potentially infectious. Handling and disposal should be in accordance with the procedures defined by an appropriate national biohazard safety guidelines or regulations.

4.4. Avoid contact with stop solution containing 5.0% H₂SO₄. It may cause skin irritation and burns.

4.5. Wear disposable latex gloves when handling specimens and reagents. Microbial contamination of reagents may give false results.

4.6. Do not use the kit beyond the expiration date.

4.7. All indicated volumes have to be performed according to the protocol. Optimal test results are only obtained when using calibrated pipettes and microplate readers.

4.8. Do not smoke, eat, drink or apply cosmetics in areas where specimens or kit reagents are handled.

4.9. Chemicals and prepared or used reagents have to be treated as hazardous waste according to the national biohazard safety guidelines or regulations.

4.10. Do not mix reagents from different lots.

4.11. Replace caps on reagents immediately. Do not swap caps.

4.12. Do not pipette reagents by mouth.

4.13. Specimens must not contain any AZIDE compounds – they inhibit activity of peroxidase.

4.14. Material Safety Data Sheet for this product is available upon request directly from XEMA Co., Ltd.

4.15. The Material Safety Data Sheet fit the requirements of EU Guideline 91/155 EC.

5.1. Contents of the Kit

5. KIT COMPONENTS

Symbol	Description	Qty	Units	Colour code	Stability of opened/ diluted components
1	SORB MTP HSV 1, 2 IgM EIA strips, 8x12 wells	1	pcs		until exp.date
2	CONTROL -, CONTROL + Control sera (0.5 ml and 0.2 ml, resp.)	2	pcs	colour- less; red	2 months
3	CONJ 11X Conjugate concentrate, 1.25 ml	1	pcs	blue	Concentrate - until exp.date Diluted - 1 day at 2-8°C
4	DIL CONJ Conjugate dilution buffer, 12 ml	1	pcs	bright yellow	until exp.date
5	DIL EIA buffer 14 ml	1	pcs	blue	until exp.date
6	SUBS TMB Substrate solution, 14 ml	1	pcs	colour- less	until exp.date
7	BUF WASH 26X Washing solution concentrate 26X, 22 ml	1	pcs	colour- less	Concentrate - until exp.date Diluted washing solution - 45 days at 2-8 °C or 15 days at RT
8	STOP Stop solution, 14 ml	1	pcs	colour- less	until exp.date
9	N003 Plate sealing tape	2	pcs		N/A
10	K104MI Instruction HSV 1, 2 IgM EIA	1	pcs		N/A
11	K104MQ QC data sheet HSV 1, 2 IgM EIA	1	pcs		N/A

5.2. Equipment and material required but not provided

- Distilled or deionized water;
- Automatic or semiautomatic multichannel micropipettes, 90–250 µl, is useful but not essential;
- Calibrated micropipettes with variable volume, range volume 10–250 µl;
- Dry thermostat or thermostat shaker for 37 °C ±2 °C;
- Calibrated microplate photometer with 450 nm wavelength and OD measuring range 0–3.0

5.3. Storage and stability of the Kit

Store the whole kit at +2...+8 °C upon receipt until the expiration date.

After opening the pouch keep unused microtiter wells TIGHTLY SEALED BY ADHESIVE TAPE (INCLUDED) to minimize exposure to moisture.

6. SPECIMEN COLLECTION AND STORAGE

This kit is intended for use with serum or plasma (ACD- or heparinized). Grossly hemolytic, lipemic, or turbid samples should be avoided.

Specimens may be stored for up to 48 hours at +2...+8 °C before testing. For a longer storage, the specimens should be frozen at -20 °C or lower. Repeated freezing/thawing should be avoided.

7. TEST PROCEDURE**7.1.** Reagent Preparation

- All reagents (including unsealed microstrips) should be allowed to reach room temperature (+18...+25 °C) before use.
- All reagents should be mixed by gentle inversion or vortexing prior to use. Avoid foam formation.
- It is recommended to spin down shortly the tubes with calibrators on low speed centrifuge.
- Prepare washing solution from the concentrate BUF WASH 26X by 26 dilutions in distilled water.

7.2. Procedural Note:

It is recommended that pipetting of all calibrators and samples should be completed within 3 minutes.

7.3. Assay flowchart

See the example of calibration graphic in Quality Control data sheet.

7.4. Assay procedure

1	Put the desired number of microstrips into the frame; allocate 4 wells for control samples CONTROL - CONTROL + (3 and 1 wells resp.) and two wells for each unknown sample. DO NOT REMOVE ADHESIVE SEALING TAPE FROM UNUSED STRIPS.
2	Pipet 90 µl of EIA buffer into each well.
3	Pipet 10 µl of control samples CONTROL - CONTROL + and unknown samples into the wells. Cover the wells by plate adhesive tape (included into the kit).
4	Incubate 30 minutes at +37 °C and continuous shaking at 500-600 rpm.
5	Prepare washing solution by 26X dilution of washing solution concentrate BUF WASH 26X by distilled water. Minimal quantity of washing solution should be 250 µl per well. Wash strips 3 times.
6	Dispense 100 µl of working conjugate solution into the wells. Cover the wells by plate adhesive tape.
7	Incubate 30 minutes at +37 °C and continuous shaking at 500-600 rpm.
8	Wash the strips 5 times.
9	Dispense 100 µl of SUBS TMB into the wells.
10	Incubate 10–20 minutes at +18...+25 °C.
11	Dispense 100 µl of STOP into the wells.
12	Measure OD (optical density) at 450 nm.
13	Set photometer blank on air.

Alternative incubation:

1	Put the desired number of microstrips into the frame; allocate 4 wells for control samples CONTROL - CONTROL + (3 and 1 wells resp.) and two wells for each unknown sample. DO NOT REMOVE ADHESIVE SEALING TAPE FROM UNUSED STRIPS.
2	Pipet 90 µl of EIA buffer into each well.
3	Pipet 10 µl of control samples CONTROL - CONTROL + and unknown samples into the wells. Cover the wells by plate adhesive tape (included into the kit).
4	Incubate 60 minutes at +37 °C.
5	Prepare washing solution by 26X dilution of washing solution concentrate BUF WASH 26X by distilled water. Minimal quantity of washing solution should be 250 µl per well. Wash strips 3 times.
6	Dispense 100 µl of working conjugate solution into the wells. Cover the wells by plate adhesive tape.
7	Incubate 30 minutes at +37 °C.

(Continuation see page 15)

8	Wash the strips 5 times.
9	Dispense 100 µl of SUBS TMB into the wells.
10	Incubate 10–20 minutes at +18...+25 °C.
11	Dispense 100 µl of STOP into the wells.
12	Measure OD (optical density) at 450 nm.
13	Set photometer blank on air.

8. QUALITY CONTROL

It is recommended to use control samples according to state and federal regulations. The use of control samples is advised to assure the day to day validity of results.

The test must be performed exactly as per the manufacturer's instructions for use. Moreover the user must strictly adhere to the rules of GLP (Good Laboratory Practice) or other applicable federal, state, and local standards and/or laws. This is especially relevant for the use of control reagents. It is important to always include, within the test procedure, a sufficient number of controls for validating the accuracy and precision of the test.

For the assay to be valid, the following requirements should be met:

1. OD450 for CONTROL+ should be nlt 1 AU.
2. OD450 for CONTROL- should not be more than 0.15 AU for all replicates.
3. OD450 for any CONTROL- replicate should be within 50%-150% of the mean OD450 value for CONTROL-. If any value lies outside this range (although meets requirement #2), it should be discarded and not used for calculation of the mean OD450 value for CONTROL-.

9. CALCULATION OF RESULTS

9.1. Calculate the mean absorbance values (OD450) for CONTROL- in triplicates and each pair of samples.

9.2. Calculate the cut-off value: (mean OD450 for CONTROL-) +0.2.

9.3. Calculate Positivity Index (PI) for each sample:

PI = mean OD450(sample) / Cut-off

10. EXPECTED VALUES AND ASSAY LIMITATIONS

If PI value is greater than 1.1, the result is POSITIVE.

If PI value is less than 0.9, the result is NEGATIVE.

If PI value is between 0.9 and 1.1, the result is EQUIVOCAL. Such samples should be retested. If the result is equivocal again, a new sample should be obtained 2-4 weeks later and tested again. If the result remains equivocal, the sample should be considered negative.

Therapeutic consequences should not be based on results of IVD methods alone – all available clinical and laboratory findings should be used by a physician to elaborate therapeutic measures.

NOTE: the patients that have received murine monoclonal antibodies for radioimaging or immunotherapy develop high titered anti-mouse antibodies (HAMA). The presence of these antibodies may cause false results in the present assay. Sera from HAMA positive patients should be treated with depleting adsorbents before assaying.

11. PERFORMANCE CHARACTERISTICS

11.1. Analytical specificity

Specificity of the test was evaluated on 47 serum specimens found negative in DiaSorin (Italy) and NovaTec (Germany). All tested specimens were found negative. Based on these data, specificity of the test is 100%.

11.2. Analytical sensitivity

The kit HSV1,2 IgM EIA was evaluated using panel of naturally occurring serum and plasma specimens from Boston Biomedica, Inc (Anti-HSV1/2 Mixed Titer Performance Panel PTH201).

11.3. Precision

Intra-assay precision is shown below:

Serum, no	duplicated	value, K	CV1, %	CV2, %
1	32	0.547	7.4	6.2
2	32	1.548	5.1	4.8

Inter-assay precision is shown below:

Serum, no	duplicated	value, K	CV, %
1	8	0.785	8.2
2	8	1.95	6.4

12. LITERATURE

1. Isakov V.A., Borisova V.V., Isakov D.V. Herpes infection: pathogenesis and laboratory diagnostics. A guideline for physicians – St.-Petersburg, Lanj - 1999 - p.192
2. Cowan F.M., Johnson A.M., Ashley R. et al. Relation ship between antibodies to herpes simplex virus (HSV) and symptoms of HSV infection // J. Infect. Dis. - 1996 - v.174-p.470-475
3. Munday P.E., Vuddamalay J., Slomka M.J. et al. Role of type specific herpes simplex virus serology in the diagnosis and management of genital herpes // Sex Transm. Inf - 1998-v.74 - p.175-178
4. Roest R.W., van der Meijden W.I., van Dijk G. et al. Prevalence and association between herpes simplex viruses types 1- and 2-specific antibodies in attendees at a sexually transmitted disease clinic. // Int Epidemiol. - 2001 - v.30 - p.580-588
5. Whitley R/I/ Herpes simplex virusinfections. In: Infection Diseases of the Fetus and Newborn Infants (Remington JS, Klein JO eds.) - Philadelphia: WB Saunders, 1990-p.282-305

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Russian Diagnostic
Manufacturers Association

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Russian Association
of Medical Laboratory
Diagnosticians

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