Aesculap® Acculan 4

en Instructions for use/Technical description
Dermatome GA340 / Dermatome 0.1 mm GA341

USA Note for U.S. users

This Instructions for Use is NOT intended for United States users. Please discard. The Instructions for Use for United States users can be obtained by visiting our website at www.aesculapusa.com. If you wish to obtain a paper copy of the Instructions for Use, you may request one by contacting your local Aesculap representative or Aesculap's customer service at 1-800-282-9000. A paper copy will be provided to you upon request at no additional cost.

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





Dermatome GA340 / Dermatome 0.1 mm GA341

Legend

- 1 Product (Dermatome)
- 2 Trigger (for oscillation frequency control)
- 3 Safety catch
- 4 Battery bay
- 5 Locking pin
- 6 Sterile funnel
- 7 Battery
- 8 Lid
- 9 Lid release element
- 10 Battery removal device
- 11 Dermatome blade
- 12 Sliding surfaces
- 13 Blade cover
- 14 Roll pin
- 15 Clamping lever
- 16 Blade guide peg
- 17 Pusher dog
- 18 Setting lever
- 19 Lock (of setting lever)
- 20 Dial wheel
- 21 Surface
- 22 Guide slots
- 23 Flap rod
- 24 Flaps
- 25 Flap tab
- 26 Nut
- 27 Symbol for alignment of flaps (on the flap rod)

The depictions are schematic only.

Symbols on product and packages

	Caution Observe important safety information such as warnings and precautions in the instructions for use.
 YYYY-MM	Maintenance label Indication of the next maintenance appointment (Date: Year-Month)
	Two-dimensional machine readable code The code contains a unique serial number that can be used for electronically tracking single instruments. The serial number is based on the worldwide standard sGTIN (GS1).
	Manufacturer
	Date of manufacture
	Follow the instructions for use

	Labeling of electrical and electronic devices pursuant to directive 2002/96/EG (WEEE)
	Classification Type BF
	Manufacturer's batch designation
	Manufacturer serial number
	Manufacturer order number
	Temperature limits during transport and storage
	Air humidity limits during transport and storage
	Atmospheric pressure limits during transport and storage
S9	Operating mode
Rx only	According to US federal law, this product may only be sold by a physician or by the order of a physician
	CE mark according to Regulation (EU) 2017/745
	Medical device
	Motor speed control
	Rotational direction for loosening the nut
	Direction of rotation for tightening the nut
	Alignment of flaps on the flap rod

Contents

1.	About this document	3
1.1	Scope	3
1.2	Warnings	3
2.	General information	3
2.1	Intended use	3
2.2	Main functions and design characteristics	3
2.2.1	Operating mode	3
2.3	Indications	4

2.4	Absolute contraindications	4
2.5	Relative contraindications	4
3.	Safe handling	4
4.	Product description	4
4.1	Scope of supply	4
4.2	Components required for operation	4
4.3	Operating principle	4
5.	Preparation	4
6.	Working with the device	4
6.1	System set-up	4
6.1.1	Inserting the rechargeable battery	5
6.1.2	Intraoperative battery change	5
6.1.3	Removing the rechargeable battery	5
6.1.4	Protection against inadvertent activation	5
6.1.5	To mount the flap rod and flaps	5
6.1.6	Inserting the dermatome blade	6
6.1.7	Removing the dermatome blade	6
6.1.8	Intraoperative storage	6
6.2	Function checks	6
6.3	Safe operation	6
6.3.1	Adjusting the cutting thickness	6
6.3.2	Adjusting the cutting width	6
6.3.3	Operating the product	6
6.3.4	Taking skin grafts	7
7.	Validated reprocessing procedure	7
7.1	General safety notes	7
7.2	General information	7
7.3	Reusable products	7
7.4	Preparations at the place of use	7
7.5	Disassembling the product before carrying out the reprocessing procedure	7
7.5.1	Removing the flap rod of the dermatome	7
7.6	Preparation before cleaning	7
7.7	Single-use products	7
7.8	Product-specific safety instructions for the processing procedure	8
7.9	Manual cleaning with wipe disinfection	8
7.10	Automatic cleaning/disinfection with manual pre-cleaning	9
7.10.1	Manual pre-cleaning with a brush	9
7.10.2	Mechanical alkaline cleaning and thermal disinfection	9
7.11	Inspection, maintenance and checks	10
7.12	Packaging	10
7.13	Steam sterilization	10
7.14	Storage	10
8.	Maintenance	10
9.	Troubleshooting list	10
10.	Technical Service	12
11.	Accessories/Spare parts	12
12.	Technical data	12
12.1	Classification according to Regulation (EU) 2017/745	12
12.2	Performance data, information about standards	12
12.3	Nominal operating mode	12
12.4	Ambient conditions	12
13.	Disposal	13

1. About this document

1.1 Scope

These instructions for use apply to the products listed below:

Art. no.	Designation
GA340	Dermatome
GA341	Dermatome 0.1 mm

► For instructions for use of specific items and information on material compatibility, see B. Braun eIFU at eifu.bbraun.com

1.2 Warnings

Warnings indicate risks to patients, users, and/or the product that could develop while using the product. Warnings are marked as follows:

DANGER

Indicates a potential danger. If it is not prevented, death or severe injuries may result.

WARNING

Indicates a potential danger. Minor or moderate injuries may result if this is not prevented.

CAUTION

Indicates a potential damage to property. The product may become damaged if it is not prevented.

2. General information

2.1 Intended use

Task/Function

When combined with a dermatome blade, the Dermatome GA340 / Dermatome 0.1 mm GA341, is used to obtain split skin grafts with a thickness that can be set.

Application Environment

The product fulfills the requirements for type BF pursuant to IEC/DIN EN 60601-1 and is used in operating rooms in sterile environments of explosion risk areas (such as areas with pure oxygen or anaesthesia gases).

2.2 Main functions and design characteristics

Oscillation frequency	min. 0 rpm to max. 6 500 rpm
-----------------------	------------------------------

2.2.1 Operating mode

Operation with non-periodic load and speed changes (type S9 pursuant to IEC EN 60034-1)

- 60 second application, 60 second pause
- 10 repetitions
- 30 min cooling time
- Max. temperature 48 °C

Electrical systems generally heat up during continual operation. It is advised to give the system a break after use to cool down, as listed in the table on operating mode.

Heating depends on the tool used and the load. After a certain number of repetitions, the system should cool down. This procedure prevents the system overheating as well as possible injury to the patient or user.

The user is responsible for the use and adherence to the pause sequence described.

2.3 Indications

Type and area of application depend on the selected tool.

2.4 Absolute contraindications

The product is not approved for use in the central nervous system or central circulatory system.

2.5 Relative contraindications

The safe and effective use of the product depends greatly on factors that only operators themselves can control. The information provided therefore only represents general conditions.

The clinically successful use of the product depends on the knowledge and experience of the surgeon. They must decide which structures can be sensibly treated and must take the safety information and warnings in the instruction for use into account.

3. Safe handling

WARNING

Risk of injury and material damage if this product is not used as intended!

- Use the product only for its intended purpose.

WARNING

Risk of injury and damage to property due to improper handling of the product!

- Follow the instructions for use of all products used.

- General risk factors associated with surgical procedures are not described in this documentation.
- It is the operating surgeon's responsibility to ensure that the surgical procedure is performed correctly.
- The operating surgeon must have a thorough understanding of both the hands-on and conceptual aspects of the established operating techniques.
- Remove the transport packaging and clean the new product, either manually or mechanically, prior to its initial sterilization.
- Check the functionality and proper condition before using the product.
- Observe "Guidance on Electromagnetic Compatibility (EMC) for Acculan Components" TA022450, see B. Braun eIFU at eifu.bbraun.com
- To prevent damage caused by improper setup or operation, and in order not to compromise warranty and manufacturer liability:
 - Use the product only according to these instructions for use.
 - Follow the safety and maintenance instructions.
 - Only combine Aesculap products with each other.
- Ensure that the product and its accessories are operated and used only by persons with the requisite training, knowledge, or experience.
- Keep the instructions for use accessible for the user.
- Comply with valid standards.
- Ensure that the electrical installation of the room is consistent with the requirements of IEC/DIN EN.
- Do not operate the product in explosion-hazard areas.
- Sterilize product before use.
- When using ECCOS holder systems, observe relevant instructions for use TA009721, see B. Braun eIFU at eifu.bbraun.com

Note

The user must report all serious incidents occurring in connection with the product to the manufacturer and to the competent authority of the country in which the user is registered.

4. Product description

4.1 Scope of supply

Art. no.	Name
GA340	Dermatome
GA341	– or – Dermatome 0.1 mm
GA678	Sterile funnel
TA014549	Instructions for use of Dermatome GA340/GA341 (flyer)

4.2 Components required for operation

- Li-Ion battery short GA346 (charged)
- Sterile funnel GA678
- Lid GA675
- Dermatome blade GB228R

4.3 Operating principle

The product **1** contains an electric motor, which is powered by a replaceable battery **7**.

The charged non-sterile battery **7** is inserted into the product **1** via the sterile funnel **6** and sealed with the lid **8** so that it is sterile.

The product moves the dermatome blade **11** in oscillation. The oscillation frequency is controlled electronically and can be continuously regulated with the trigger **2**.

5. Preparation

Aesculap assumes no liability if the following rules are not followed:

- Do not use products from open or damaged sterile packaging.
- Prior to use, inspect the product and its accessories for any visible damage.
- Use the products and their accessories only if they are in perfect technical condition.

6. Working with the device

6.1 System set-up

WARNING

Risk of infection and contamination!

Product is delivered unsterilized!

- Sterilize the product before use pursuant to the operating instructions.

WARNING

Risk of injury and material damage due to accidental activation of the product!

- Products which are not being actively used must be secured against accidental activation (position OFF).

WARNING

Risk of injury and material damage due to inappropriate use of tools!

- Always follow the safety advice and information given in the instructions for use.
- When coupling/uncoupling, handle tools with cutting edges with care.

⚠ WARNING

Damage to the product if dropped!

- ▶ Use the products only if they are in perfect technical condition, see Function check

⚠ WARNING

Risk of burns to skin and tissue through blunt tools or if product has not been maintained properly!

- ▶ Use tools only if they are in perfect condition.
- ▶ Replace blunt tools.
- ▶ Maintain the product properly, see maintenance guide.

⚠ WARNING

Risk of injury and/or malfunction due to unrecognized change to cutting setting!

The basic calibration of the cutting setting can change if the product is dropped, for instance.

- ▶ Do not use product.
- ▶ Have the manufacturer inspect the product.

6.1.1 Inserting the rechargeable battery

- ▶ Product 1 with the battery shaft 4 must be turned upwards and the sterile funnel 6 (sterile) attached, see Fig. A.
- ▶ Battery 7 (non-sterile) is inserted into the battery shaft 4 by a second (non-sterile) person, see Fig. A.

Note

After the battery has been inserted, multiple signals will sound, which indicate that the product is ready to use.

- ▶ After insertion of the battery, have the sterile funnel 6 (non-sterile) removed by a second person.
- ▶ Attach the lid 8 (sterile) in such a manner that it clicks into place with both release catches 9.

Note

The sterility of the product is only guaranteed with a correctly placed lid.

6.1.2 Intraoperative battery change**⚠ WARNING**

Risk of injury from hot battery!

Battery may be hot after use in the machine.

- ▶ Remove the battery with the battery removal device and let it cool.

The battery removal device is used for replacing the battery during surgery while ensuring that sterile conditions are maintained.

- ▶ Turn the product 1 with the battery shaft 4 upwards.
- ▶ Press both safety catches 9 on the lid 8 at the same time and remove the lid 8.
- ▶ Attach the sterile battery removal device 10, see Fig. B.
- ▶ Gently shake the product 1 with the attached battery removal device 10 with battery compartment 4 facing downwards.
The battery 7 slides easily into the battery removal device 10.
- ▶ Hand battery removal device 10 including dead battery 7 to a non-sterile person.
- ▶ Insert a charged battery 7, see Inserting the rechargeable battery.

6.1.3 Removing the rechargeable battery**⚠ WARNING**

Risk of injury from hot battery!

Battery may be hot after use in the machine.

- ▶ Let the battery cool in the machine and then remove.
– or –
- ▶ Remove the battery with the battery removal device and let it cool.

⚠ CAUTION

Damage to the battery caused by knocking it against hard objects!

- ▶ Only remove battery by tapping the product on the palm of the hand.

⚠ CAUTION

Damage to, or destruction of the batteries due to processing!

- ▶ Do not sterilize batteries.

Rechargeable battery must be removed after each surgical procedure and prior to processing of the device.

Note

The battery removal device can be used for easier battery removal, see Fig. B.

- ▶ Turn the product 1 with the battery shaft 4 upwards.
- ▶ Press both safety catches 9 on the lid 8 at the same time and remove the lid 8.
- ▶ Tightly grasp the battery shaft 4 at the lower end.
- ▶ Tap the battery shaft 4 on the palm of the hand until the battery 7 slides out of the battery shaft 4 and can be removed.

6.1.4 Protection against inadvertent activation

To prevent inadvertent activation of the product during blade change, for example, the trigger can be locked.

To lock trigger:

- ▶ Twist the safety catch 3 to position OFF.
The trigger 2 is blocked and the product 1 cannot be operated.

To unlock trigger:

- ▶ Twist the safety catch 3 to position ON.
The trigger 2 is unlocked and the product 1 can be operated.

6.1.5 To mount the flap rod and flaps

- ▶ Install the flaps 24 on flap rod 23 to the edge of the rear surface. Observe the symbol 27 on the face of the thread, see Fig. C.
- ▶ Screw the nut 26 counterclockwise onto the thread of flap rod 23 (left-handed thread).
- ▶ Turn the nut 26 up to the edge of the visible surface.
- ▶ Insert and twist the complete flap rod 23 in its guide slots 22 on the side surfaces.
- ▶ Push the flap rod 23 laterally all the way so that the cross pin of the flap rod 23 comes to rest in the guide slot 22.
- ▶ Tighten nut 26 by turning it counterclockwise.

6.1.6 Inserting the dermatome blade

DANGER

Risk of infection and contamination due to torn or cut surgical gloves!

- Avoid contact with the cutting edge when inserting the dermatome blade.

DANGER

Risk of injury and material damage due to inadvertent activation of the dermatome while inserting the dermatome blade!

- Prior to inserting the dermatome blade, engage the safety catch, see Protection against inadvertent activation.
- Engage safety catch 3.
- Push both clamping levers 15 in the direction of the arrows, see Fig. C. Blade cover 13 is now unlocked.
- Open blade cover 13 in the direction of the arrow.
- Securely insert dermatome blade 11 in tappet 17 and blade guide pin 16, see Fig. D.
- Close blade cover 13.
- Tighten both clamping levers 15 so that roll pin 14 engages in the recess of clamping lever 15.

6.1.7 Removing the dermatome blade

- Release both clamping levers 15.
- Open blade cover 13.
- Remove the dermatome blade 11 from tappet 17 and blade guide pin 16.

6.1.8 Intraoperative storage

Note

The ECCOS holder can be used for depositing during surgery. However, the dermatome must not be processed in this position. For the correct processing position, refer to see Validated reprocessing procedure.

- Engage safety catch 3.
- Fold back the clamp of the ECCOS holder.
- Insert product 1 in the ECCOS holder, see Fig. E

6.2 Function checks

The function checks must be carried out prior to each use and after each intraoperative battery change.

- Make certain there is a battery in the battery compartment.
- Ensure that the lid has completely clicked into place.
- Ensure that the dermatome blades are not mechanically damaged.
- Ensure that a dermatome blade is correctly seated in the dermatome.
- Check that the dermatome blade is securely seated.
- Check that the roll pin is securely seated.
- Release product for use (ON position).
- Briefly operate product at maximum oscillation frequency.
- Make certain that the rotational direction is correct in each case.
- Do not use the product if it is damaged or defective.
- Set aside the product if it is damaged.

6.3 Safe operation

WARNING

Risk of infection from aerosol formation!

Risk of injuries caused by particles scaling off the tool!

- Use suitable protection (such as waterproof protective clothing, face mask, safety gases, suction).

WARNING

Risk of injury and/or malfunction!

- Always carry out a function check prior to using the product.

WARNING

Risk of injury when using the product beyond the field of view!

- Apply the product only under visual control.

Note

The drive motor of the product is operated by a magnetic sensor system. In order to prevent inadvertent activation of the motor, the product may not be exposed to any magnetic fields (such as magnetic instrument pads).

Note

The slight whistling noise when the product is starting up is due to the construction.

6.3.1 Adjusting the cutting thickness

Note

The skin condition (e.g. due to age) of the patient must be considered when adjusting the cutting depth (graft thickness).

- Set the cutting depth with setting lever 18. The lock 19 of the setting lever 18 prevents any unintentional change to the setting of the cutting depth, see Fig. C.
- Read the cutting depth on dial wheel 20 (1/10 mm gradation).
- Follow the troubleshooting information, see Troubleshooting list.

6.3.2 Adjusting the cutting width

WARNING

Risk of injury/cuts by the dermatome blade if the flaps are opened incorrectly!

- Always use the flap tabs for lifting the flaps.

Note

The width of the skin graft will be approx. 2mm less than the cutting width set by the flaps.

- Set the required cutting width by opening the appropriate number of flaps 24.
- Safely lift flaps 24 by pressing flap tabs 25, see Fig. F.

6.3.3 Operating the product

Note

The slight whistling noise when the product is starting up is due to the construction.

- Activate trigger 2.
The saw runs at the continuously controlled oscillation frequency.

6.3.4 Taking skin grafts

Note

To familiarize with the operation of the Dermatome prior to its first surgical application, the user should carry out some test cuts on a specimen.

- ▶ Stretch flabby skin.
- ▶ With the sliding surface in level contact with the skin, push product 1 forward in a smooth movement, applying gentle pressure. Be careful not to cant the product 1.
- ▶ For very thin grafts, use a forceps, from time to time, to lay up the resected portion of the graft on surface 21, see Fig. C.
- ▶ Turn off the motor and pull the freshly cut skin graft out of the device – or –
- ▶ Lower the handle and cut off the skin graft with the product 1.

7. Validated reprocessing procedure

7.1 General safety notes

Note

Adhere to national statutory regulations, national and international standards and directives, and local, clinical hygiene instructions for reprocessing.

Note

For patients with Creutzfeldt-Jakob disease (CJD), suspected CJD, or possible variants of CJD, observe the relevant national regulations concerning the reprocessing of products.

Note

Mechanical reprocessing should be favored over manual cleaning as it gives better and more reliable results.

Note

It should be noted that successful reprocessing of this medical device can only be guaranteed following prior validation of the reprocessing method. The operator/reprocessing technician is responsible for this.

Note

If there is no final sterilization, then a virucidal disinfectant must be used.

Note

See also B.Braun eIFU at eifu.bbraun.com for current information on processing and compatibility of material

The validated steam sterilization procedure was carried out in the Aesculap sterile container system.

7.2 General information

Dried or affixed surgical residues can make cleaning more difficult or ineffective and lead to corrosion. Therefore the time interval between application and processing should not exceed 6 h; also, neither fixating pre-cleaning temperatures >45 °C nor fixating disinfecting agents (active ingredient: aldehydes/alcohols) should be used.

Excessive neutralizing agents or basic cleaners may result in a chemical attack and/or fading and the laser marking becoming unreadable either visually or by machine.

On stainless steel, residues containing chlorine or chloride (such as surgical residues, drugs, saline solutions in water for cleaning, disinfection and sterilization) may lead to corrosion (pitting corrosion, tensile corrosion) and thus to the destruction of the product. These must be removed by rinsing thoroughly with demineralized water and then drying.

Perform additional drying, if necessary.

Only process chemicals that have been tested and approved (e.g. VAH or FDA approval or CE mark) and which are compatible with the product's materials according to the chemical manufacturers' recommendations may be used for processing the product. All the chemical manufacturer's application specifications must be strictly observed. Failure to do so can result in the following problems:

- Material damage (such as corrosion, cracks, breaks, premature aging or swelling).
- ▶ Do not use metal cleaning brushes or other abrasives that would damage the product surface and could cause corrosion
- ▶ Further detailed advice on hygienically safe and material-/value-preserving reprocessing can be found at www.a-k-i.org link to "AKI-Brochures", "Red brochure".

7.3 Reusable products

The service life of the product is limited by damage, normal wear, type and duration of the application, handling, storage and transportation of the product.

Careful visual and functional testing prior to next use is the best way to identify a malfunctioning product.

7.4 Preparations at the place of use

- ▶ Remove all attached components from the product (dermatome blade, battery, accessories).
- ▶ Remove any visible surgical residues as much as possible with a damp, lint-free cloth.
- ▶ Place the dry product in a sealed waste container and forward it on for cleaning and disinfection within 6 hours.

7.5 Disassembling the product before carrying out the reprocessing procedure

- ▶ Remove battery 7 from the product 1, see Removing the rechargeable battery.
- ▶ Open the clamping levers 15 to remove the dermatome blade 11, see Removing the dermatome blade.
- ▶ Disassemble the product 1 immediately after use, according to instructions.

7.5.1 Removing the flap rod of the dermatome

- ▶ Loosen nut 26 by turning it clockwise (left-handed thread).
 - ▶ Screw the nut 26 back to the end of the visible surface.
 - ▶ Press on nut 26 and push flap rod 23 approx. 4 mm to one side.
 - ▶ Twist flap rod 23 until it can be removed.
 - ▶ Remove flap rod 23.
 - ▶ Slide off the flaps 24 from flap rod 23.
- The flap rod 23 has been dismantled.

7.6 Preparation before cleaning

- ▶ Prior to first mechanical cleaning/disinfection: Mount the ECCOS holders in a suitable basket (e.g. JC254R) or use the ECCOS basket GB256R equipped with holders.
- ▶ Insert the products in the correct position into the ECCOS holders, see Fig. G.

7.7 Single-use products

Dermatome blade GB228R

⚠ WARNING

Risk of infection for patients and/or users and impairment of product functionality due to reuse. Risk of injury, illness or death due to contamination and/or impaired functionality of the product!

- ▶ Do not reprocess the dermatome blade GB228R.

7.8 Product-specific safety instructions for the processing procedure

⚠ CAUTION

Damage to the product due to inappropriate cleaning/disinfecting agents and/or excessive temperatures!

- ▶ Use cleaning and disinfecting agents according to the manufacturer's instructions which
 - are approved for plastic material and high-grade steel,
 - and which do not affect softeners (e.g., in silicone).
- ▶ Do not use cleaning agents that contain acetone.
- ▶ Observe information concerning concentration, temperature and exposure time.
- ▶ Do not exceed the maximum temperature of 60 °C in chemical cleaning and/or disinfection.

- ▶ Do not exceed maximum temperature during thermal disinfection using 96 °C DI water.
- ▶ Dry the product for at least 10 minutes at a maximum of 120 °C.

⚠ CAUTION

Damage or destruction of the batteries due to processing!

- ▶ Protect battery from moisture.

Note

The drying time listed is only indicative. It has to be checked taking into account the specific conditions (e.g. load) and adjusted if necessary.

7.9 Manual cleaning with wipe disinfection

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemicals
I	Pre-cleaning	RT (cold)	≥2	–	D–W	Until visually clean
II	Cleaning with enzyme solution	RT (cold)	≥2	0.8	D–W	pH-neutral*
III	Intermediate rinse	RT	≥5	–	D–W	–
IV	Drying	RT	–	–	–	–
V	Wipe disinfection	–	>1	–	–	Meliseptol HBV wipes 50 % Propan-1-ol
VI	Final rinse	RT (cold)	0.5	–	FD–W	–
VII	Drying	RT	–	–	–	–

D–W: Drinking water

FD–W: Fully de-ionized water (de-mineralized, microbiological, with minimum quality of drinking water)

RT: Room temperature

* Suitable enzyme solution: Helizyme, Cidezyme (the latter for validation purposes)

- ▶ Do not clean the product in a ultrasonic bath and do not immerse the product in any fluids. Let any fluid incursions drain out immediately, otherwise there is a danger of corrosion and loss of function.

Phase I

- ▶ Mobilize non-rigid components during cleaning.
- ▶ Clean the product under running water, using a suitable cleaning brush until all visible residues have been removed from the surfaces.
- ▶ Brush difficult to access surfaces with a suitable plastic cleaning brush for at least 1 min.

Note

For details on difficult to access surfaces, see Acculan pre-cleaning and care information TA016000 (available in Aesculap Extranet at <https://extranet.bbraun.com>).

Phase II

- ▶ Follow the operating instructions of the enzyme cleaner with regard to correct concentration, dilution, temperature and water quality.
- ▶ Spray products with a pH neutral enzyme solution, let soak in for at least 2 minutes and then wipe off.
- ▶ Contamination should be removed with a lint-free cloth or soft brush moistened with enzyme cleaner.

- ▶ Flexible components should be rinsed for 20 seconds with the water pistol (cold water, at least 2.5 bar).
- ▶ After manual cleaning, check visible surfaces and areas of flexible components for residues.

Phase III

- ▶ Mobilize non-rigid components during cleaning.
- ▶ Rinse product under running tap water for at least 5 minutes.
- ▶ If necessary, repeat the cleaning process (phase 1 to 3).

Phase IV

- ▶ Dry the product in the drying phase with suitable equipment (such as lint-free cloths, pressurized air).

Phase V

- ▶ Wipe all surfaces of the product with a single-use disinfecting wipe.

Phase VI

- ▶ Rinse disinfected surfaces after the prescribed reaction time for at least 1 minute under running demineralized water.
- ▶ Drain any remaining water fully.

Stage VII

- ▶ Dry the product in the drying phase with suitable equipment (such as lint-free cloths, pressurized air).

7.10 Automatic cleaning/disinfection with manual pre-cleaning

Note

The cleaning and disinfection device must fundamentally have a tested efficacy (such as FDA approval or CE label pursuant to DIN EN ISO 15883).

Note

The cleaning and disinfection machine used for processing must be serviced and checked at regular intervals.

7.10.1 Manual pre-cleaning with a brush

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemistry/Note
I	Rinse	RT (cold)	–	–	D–W	Until visually clean
II	Brush	RT (cold)	–	–	D–W	Until visually clean

D–W: Drinking water

RT: Room temperature

- ▶ Do not clean the product in a ultrasonic bath and do not immerse the product in any fluids. Let any fluid incursions drain out immediately, otherwise there is a danger of corrosion and loss of function.

Phase I

- ▶ Dismantle the flap rod for cleaning. Store flap rods, flaps and nut separately in tray.
- ▶ Mobilize non-rigid components during cleaning.
- ▶ Thoroughly clean the product under running water.

Phase II

- ▶ Mobilize non-rigid components during cleaning.
- ▶ Brush difficult to access surfaces with a suitable plastic cleaning brush for at least 1 min.
- ▶ After manual preparation, check visible surfaces for residue and repeat the pre-cleaning process as needed.

Note

For details on difficult to access surfaces, see Acculan pre-cleaning and care information TA016000 (available in Aesculap Extranet at <https://extranet.bbBraun.com>).

7.10.2 Mechanical alkaline cleaning and thermal disinfection

Machine type: single-chamber cleaning/disinfection device without ultrasound

Phase	Step	T [°C/°F]	t [min]	Water qual- ity	Chemical
I	Pre-rinse	<25/77	3	D–W	–
II	Cleaning	55/131	10	FD–W	■ Concentrate, alkaline: – pH ~ 13 – <5% anionic surfactants ■ Working solution 0.5% – pH ~ 11*
III	Intermediate rinse	>10/50	1	FD–W	–
IV	Thermal disinfection	90/194	5	FD–W	–
V	Drying	max. 120/248	min. 10 min	–	–

D–W: Drinking water

FD–W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

*Recommended: BBraun Helimatic Cleaner alkaline

- ▶ Insert the product in its proper position in the ECCOS holder. Attach the battery compartment over the clamp and swivel device into the ECCOS holder until the knob 2 points upwards.
- ▶ Ensure that the blade cover 13 is kept open through the entire processing procedure.
- ▶ Dismantle the flap rod 23 for cleaning. Store flap rod 23, flaps 24 and nut 26 separately in tray.
- ▶ After automatic cleaning/disinfection, check visible surfaces for residues and repeat the cleaning/disinfection process as needed.

7.11 Inspection, maintenance and checks

- ▶ Allow the product to cool down to room temperature.

Note

Aesculap recommends lubricating the flap rod, flaps, sliding surfaces, and other moving parts, (e.g. knob, cover lids) using STERILIT Power Systems oil spray GB600 prior to every sterilization. Remove excess oil using a lint-free wipe.

- ▶ Inspect the product after each cleaning and disinfecting cycle for: cleanliness, damage, function, abnormal operation noise, excessive heat, or heavy vibration.
- ▶ Set aside the product if it is damaged.

7.12 Packaging

- ▶ Observe the instructions for use of the respective packaging and holders (e.g. Instructions for use TA009721 for Aesculap ECCOS holder systems).
- ▶ Insert the products in the correct position into the ECCOS holders, see Fig. G.
- ▶ Pack trays appropriately for the sterilization process (e.g. in Aesculap sterile containers).
- ▶ Ensure that the packaging will prevent a recontamination of the product.

7.13 Steam sterilization

⚠ CAUTION

Damage or destruction of the batteries due to processing!

- ▶ Do not sterilize the rechargeable battery.

Note

The product may only be sterilized with the blade, battery and lid as well as opened blade cover all demounted.

Note

The product can be sterilized with a mounted flap tab and flaps in the open position.

- ▶ Check to ensure that the sterilizing agent will come into contact with all external and internal surfaces (e.g. by opening any valves and faucets).
- ▶ Use a validated sterilization method:
 - Steam sterilization using fractional vacuum process
 - Steam sterilizer DIN EN 285 and validated pursuant to DIN EN ISO 17665
 - Sterilization using fractionated vacuum process at 134 °C/holding time 5 min

When sterilizing multiple products in one steam sterilizer:

- ▶ Ensure that the maximum permitted load specified by the manufacturer for the steam sterilizer is not exceeded.

7.14 Storage

- ▶ Store sterile products in germ-proof packaging, protected from dust, in a dry, dark, temperature-controlled area.
- ▶ Store sterile single-use GB228R product in germ-proof packaging in a dry, dark and temperature-controlled room.

8. Maintenance

To ensure reliable operation, the product must be maintained in accordance with the maintenance labeling or at least once a year.

YYYY-MM

For technical service, please contact your national B. Braun/Aesculap agency, see Technical Service.

9. Troubleshooting list

- ▶ Have defective products repaired by Aesculap Technical Service, see Technical Service.

Issue	Cause	Detection	Remedy
Product not running	No battery	No battery in the battery shaft	Insert battery.
	Battery not charged	No signal tone when inserting the battery	Charge battery in charger.
	Battery defective	No signal tone when inserting the battery	Have the manufacturer repair the battery.
	Product in locked OFF position	Safety catch is in OFF position	Turn the safety catch to the ON position.
	Product defective	Product not running	Have the manufacturer repair the product.
Product becomes too hot	Excessive use	Product becomes hot	Observe operating instructions (nominal operating mode).
	Processing/care not performed properly	Product heating	Follow operating instructions (preparation, care). Preventive measure: Oil product before every sterilization.
	Fall damage, product defective	Product becomes hot	Return the product to the manufacturer for repairs.
	Blunt tool	Heating of tool and product	Change tool.

Issue	Cause	Detection	Remedy
Insufficient power	Product defective	Insufficient product performance	Observe instructions for use (processing, care). Preventative: lubricate product before every sterilization. Have the manufacturer repair the product.
		Steep rise in temperature after a short period of time	Follow operating instructions (operating mode). Return the product to the manufacturer for repairs
	Blunt tool	Tool blades worn out	Change tool.
Loud running noise	Defective product drive/ball bearings	Loud, perceptible noise during operation	Follow operating instructions (preparation, care). Preventive measure: Oil product before every sterilization. Return the product to the manufacturer for repairs
Lid cannot be fitted/removed	Lid not compatible	Lid does not click into place	Use the compatible lid for GA340/GA341.
	Lid deformed/broken	Lid is hard to fit/remove or cannot be fitted/removed	Have the manufacturer repair the lid.
	Connection to product deformed/defective	Lid is difficult or cannot be mounted/removed	Have the manufacturer repair the product.
Battery cannot be inserted/removed	Battery not compatible	Battery cannot be completely inserted into the battery shaft	Use the compatible battery for GA340/GA341.
	Battery deformed/defective	Battery is hard to fit/remove or cannot be fitted/removed	Have the manufacturer repair the battery.
	Battery compartment on product deformed/broken	Battery is difficult or cannot be mounted/removed	Have the manufacturer repair the product
Sterile funnel or battery removal device cannot be attached	Sterile funnel or battery removal device not compatible	Sterile funnel or battery removal device cannot be attached to the battery compartment	Use compatible sterile funnel or battery removal device for GA340/GA341.
	Sterile funnel or battery removal device deformed/broken	Sterile funnel or battery removal device difficult to attach or cannot be attached	Change sterile funnel or battery device.
	Battery shaft to product deformed/defective	Sterile funnel or battery removal device difficult or not attachable	Have the manufacturer repair the product.
Speed control knob not working	Product in locked OFF position	Safety catch is in OFF position	Turn the safety catch into the ON position.
	Speed control knob sticks/is defective	Speed control knob not working	Have the manufacturer repair the product.
Dermatome blade not moving	Defective gearings	High operating noise levels of product	Return the product to the manufacturer for repairs
Insufficient cutting performance of dermatome blade	Dermatome blade blunted	Worn cutting edge of dermatome blade	Replace dermatome blade
	Dermatome maladjusted	Performance of product is inadequate	Have the manufacturer repair the product
Cutting depth uneven/deviating from setting	Skin grafts uneven/too thick/too thin	Mechanical effect (e.g. dropped)	Return the product to the manufacturer for repairs

10. Technical Service

DANGER

Danger to life of patients and users in case of malfunctions and/or failure of protective measures!

- ▶ Do not perform any servicing or maintenance work under any circumstances while the product is being used on a patient.

CAUTION

Modifications of medical devices may result in a loss of potential guarantee/warranty claims and forfeiture of applicable licenses.

- ▶ Do not modify the product.
- ▶ Contact national B. Braun/Aesculap representative for service and repair.

Service addresses

Aesculap Technischer Service

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78532 Tuttlingen / Germany

Phone: +49 (7461) 95 -1601

Fax: +49 (7461) 14 -939

E-Mail: ats@aesculap.de

Other service addresses can be obtained from the address indicated above.

11. Accessories/Spare parts

Art. no.	Name
GA346	Li-Ion battery short
GA340345	Nut
GA643817	Flap rod
GA670210	Flap
GA675	Lid
GA678	Sterile funnel
GA679	Battery removal device
GB228R	Dermatome blade
GB256R	ECCOS set for dermatome
GB487R	ECCOS fixation for lid
GB488R	ECCOS fixation for sterile funnel
GB489R	ECCOS holder for battery removal device
GB498R	ECCOS fixation for straight machine
GB600	STERILIT Power Systems oil spray
TA014548	Instructions for use of Dermatome GA340/GA341 (A4 for ring binder)
TA014549	Instructions for use of Dermatome GA340/GA341 (flyer)

12. Technical data

12.1 Classification according to Regulation (EU) 2017/745

Art. no.	Name	Class
GA340	Dermatome	Ila
GA341	0.1 mm Dermatome	Ila

12.2 Performance data, information about standards

Max. power	approx. 250 W
Max. oscillation frequency	6 500 rpm
Blade travel	approx. 3.1 mm
Cutting depth setting	0.2 mm to 1.2 mm (with GA340) 0.1 mm to 1.1 mm (with GA341) Scale gradations 1/10 mm
Cutting width	max. approx. 78 mm, set by 10 flaps of width 8 mm
Weight (ready for operation)	1.65 kg ±10 %
Dimensions (L x W x H, ready for operation)	277 mm x 113 mm x 83 mm ±5 %
Application part	Type BF
EMC	IEC/DIN EN 60601-1-2
Compliance with standards	IEC/EN 60601-1

The product was subject to 500 preparation cycles as a test by the manufacturer and has passed.

12.3 Nominal operating mode

Operation with non-periodic load and speed changes (type S9 pursuant to IEC EN 60034-1)

- 60 second application, 60 second pause
- 10 repetitions
- 30 min cooling time
- Max. Temperature 48 °C

12.4 Ambient conditions

	Operation	Transport and storage
Temperature	10 °C to 27 °C	-10 °C to 50 °C
Relative humidity	30 % to 75 %	10 % to 90 %
Atmospheric pressure	700 hPa to 1 060 hPa	500 hPa to 1 060 hPa

13. Disposal

WARNING

Risk of infection from contaminated products!

- Observe national regulations when disposing of or recycling the product, its components, and their packagings.

Note

The operator must process the product before disposal, see Validated reprocessing procedure.

The recycling pass can be downloaded from the Extranet as a PDF document under the respective article number. (The recycling pass includes disassembling instructions for the product, as well as information for proper disposal of components harmful to the environment.)

Products carrying this symbol are subject to separate collection of electrical and electronic devices. Within the European Union, disposal is taken care of by the manufacturer as a free-of-charge service.

- For questions regarding the disposal of the product, please contact your local B. Braun/Aesculap representative, see Technical Service.

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


Skin mesh dermatome BA720R

Legend

- 1 Carrying handle
- 2 Cutting cylinder
- 3 Counter cylinder
- 4 Cutting device
- 5 Carrier plate for skin graft
- 6 Setting wheel for cutting depth
- 7 Gearwheel transmission
- 8 Cutting cylinder lock
- 9 Hinge
- 10 Shaft of cutting cylinder 2
- 11 Fastening screws for upper part 17
- 12 Hand ratchet for driving the cutting cylinder
- 13 Screwdriver end for screws 11 and shaft 10
- 14 Lower part
- 15 Guide strips for carrier plates
- 16 Insertion tray for carrier plates 5
- 17 Upper part, hinged
- 18 Scale of the setting wheel
- 19 Removal tray for carrier plates

Symbols on product and packages

	Caution, general warning symbol Caution, see documentation supplied with the product
	Scale symbol for "Skin graft cutting"
	Date of manufacture

Contents

1.	Applicable to	2
2.	Safe handling	3
3.	Product description	3
3.1	Scope of supply	3
3.2	Intended use.	3
3.3	Operating principle	3
4.	Preparation.	3
5.	Working with the skin mesh dermatome BA720R	3
5.1	System set-up	3
5.2	Function checks.	3
5.3	Safe operation	3
6.	Validated reprocessing procedure	5
6.1	General safety instructions.	5
6.2	Single-use products.	5
6.3	General information	5
6.4	Preparations at the place of use.	6
6.5	Preparation before cleaning	6
6.6	Cleaning/disinfection	6
6.7	Manual cleaning/disinfecting	6
6.8	Mechanical cleaning/disinfecting	7
6.9	Inspection, maintenance and checks.	8
6.10	Packaging	8
6.11	Steam sterilization.	8
6.12	Sterilization for the US market.	8
6.13	Storage	8
7.	Maintenance	8
8.	Troubleshooting list.	8
9.	Technical Service	9
10.	Accessories/Spare parts.	9
11.	Technical data	9
11.1	Ambient conditions.	9
12.	Disposal	9
13.	Distributor in the US/Contact in Canada for product information and complaints	9

1. Applicable to

- For item-specific instructions for use and information on material compatibility, see also the Aesculap Extranet at <https://extranet.bbraun.com>

2. Safe handling

CAUTION

Federal law restricts this device to sale by, or on order of a physician!

The operation and reprocessing instructions for the skin mesh dermatome may differ from these instructions for use!

These instructions for use describe only the operation and reprocessing of the skin mesh dermatome as of series number 4000.

- For upgrades, please contact your national B. Braun/Aesculap agency.

- Remove the transport packaging and clean the new product, either manually or mechanically, prior to its initial sterilization.
- Prior to use, check that the product is in good working order.
- To prevent damage caused by improper setup or operation, and in order not to compromise warranty and manufacturer liability:
 - Use the product only according to these instructions for use.
 - Follow the safety and maintenance instructions.
 - Only combine Aesculap products with each other.
- Ensure that the product and its accessories are operated and used only by persons with the requisite training, knowledge, or experience.
- Keep the instructions for use accessible for the user.

3. Product description

3.1 Scope of supply

Designation	Art. no.
Skin mesh dermatome, complete with hand ratchet	BA720R
Carrier plate for expansion factor 1.5; sterile	BA721
Carrier plate for expansion factor 3; sterile	BA722
Carrier plate for expansion factor 6; sterile	BA723

3.2 Intended use

The skin mesh dermatome BA720R is used for the cutting of skin grafts in general, burns and plastic surgery.

3.3 Operating principle

The cutting cylinder **2** is turned anticlockwise with hand ratchet **12**. The carrier plate **5** transports the skin graft between cutting cylinder **2** and counter cylinder **3**. During this transport, the skin graft is cut.

The cutting depth is set with setting wheel **6**, using scale **18**. The staggered slits in the skin graft allow expanding the graft to rhombic skin mesh.

The skin mesh can be used for covering defects of up to six times the size of the grafting site.

4. Preparation

Non-compliance with the following instructions will preclude all responsibility and liability in this respect on the part of Aesculap.

- Prior to use, inspect the skin mesh dermatome BA720R and its accessories for visible damage.
- Only use the skin mesh dermatome BA720R and its accessories if they are in perfect condition.
- Be certain that the skin mesh dermatome is set up on a sufficiently stable support (e.g. table, equipment cart, etc.).

5. Working with the skin mesh dermatome BA720R

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The device and its accessories are delivered in unsterile condition.

- Sterilize the device and accessories before use.

- Slide hand ratchet **12** on shaft **10** so that it clicks into position.

5.2 Function checks

Carry out a trial run of the skin mesh dermatome BA720R prior to each use:

- Verify that gearwheel transmission **7** is moved when turning hand ratchet **12**.
- Make certain that setting wheel **6**, once engaged, can be easily rotated through the range of scale **18**.
- Only use the device if it is in perfect condition.

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Setting the cutting cylinder

Damage to the cutting cylinder! The cutting cylinder can only be set when running freely.

- Make certain that the carrier plate is not engaged under the cutting cylinder.

- Slightly pull out setting wheel **6** until it disengages.
- Turn setting wheel **6** to the required position.

Skin graft cutting	
Possible settings at the setting wheel	Scale symbols
extremely thick skin graft	
normal setting for cutting	
reserve setting for cutting	

- Slightly push in setting wheel **6** until it engages.

Processing the skin graft

Note

Only place individual carrier plates on the insertion tray, one after the other.
For long skin grafts, sequentially insert several carrier plates.

CAUTION

Coiling of the skin graft! The skin graft can be lifted off from carrier plate as soon as it has been cut into by the cutting cylinder!

- ▶ Press the skin graft onto the carrier plate, using fine tweezers, directly behind the cut-in line.
- ▶ If the cutting cylinder has been turned too far already, turn back the cutting cylinder, using the screwdriver end of the ratchet.

- ▶ Position the required carrier plate 5, with the grooved surface facing up, on insertion tray 16.
- ▶ Position the skin graft on carrier plate 5.
- ▶ Push carrier plate 5 between guide strips 15 under cutting cylinder 2 until the cutting cylinder engages.
- ▶ Turn hand ratchet 12 anticlockwise.
The carrier plate is transported.
- ▶ Keep turning hand ratchet 12 until the complete carrier plate 5 is in position on removal tray 9.
This completes the processing of the skin graft.

Changing the cutting cylinder

Note

The cutting cylinder cannot be sharpened.

Dismounting the cutting cylinder

WARNING

Cuts on hands or other body parts, caused by sharp edges on the cutting cylinder!

- ▶ Touch the cutting cylinder only at the shaft.

- ▶ Turn cutting device 4 upside down so that the four screws 11 are visible, see Fig. 1.


Fig. 1

- ▶ Loosen screws 11 by hand or, if they are difficult to turn, with the screwdriver end 13 of hand ratchet 12.
- ▶ Return cutting device 4 to its upright position.

- ▶ Open upper part 17, see Fig. 2.


Fig. 2

- ▶ Hold down lock 8 and, at the same time, pull at shaft 10 to decouple cutting cylinder 2, see Fig. 3.


Fig. 3

- Turn cutting cylinder **2** while pulling it out of bearings **20**, see Fig. 4. This completes the dismounting of cutting cylinder **2**.


Fig. 4

Mounting the cutting cylinder

Note

To prevent transport damage, cutting cylinder **2** is protected by a piece of cardboard on delivery. Remove the cardboard protection only after cutting cylinder **2** has been mounted completely.

- Insert the cutting cylinder **2** into the two bearings **20**, see Fig. 4.
- Hold cutting cylinder **2** at shaft **10** and press it into bearings **20** until it clicks into position, see Fig. 2.
- Turn shaft **10** and check cutting cylinder **2** for correct seating.
- Close upper part **17** on lower part **14**, see Fig. 2.
- Turn cutting device **4** upside down so that the four screws **11** are visible, see Fig. 1.
- Tighten screws **11** by hand.
- Return cutting device **4** to its upright position.
- Remove the piece of cardboard from cutting cylinder **2**. This completes the mounting of cutting cylinder **2**.

6. Validated reprocessing procedure

WARNING

The reprocessing instructions for the skin mesh dermatome may differ from these instructions for use!

These instructions for use describe only the operation and reprocessing of the skin mesh dermatome as of series number 4000.

- For upgrades, please contact your national B. Braun/Aesculap agency.

6.1 General safety instructions

Note

Adhere to national statutory regulations, national and international standards and directives, and local, clinical hygiene instructions for sterile processing.

Note

For patients with Creutzfeldt-Jakob disease (CJD), suspected CJD or possible variants of CJD, observe the relevant national regulations concerning the reprocessing of products.

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Note

Successful processing of this medical device can only be ensured if the processing method is first validated. The operator/sterile processing technician is responsible for this.

The recommended chemistry was used for validation.

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Note

For up-to-date information about reprocessing and material compatibility, see also the Aesculap Extranet at <https://extranet.bbraun.com>

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6.2 Single-use products

WARNING

Infection hazard for patients and/or users and impairment of product functionality due to reuse. Risk of injury, illness or death due to contamination and/or impaired functionality of the product!

- Do not reprocess the product!

Designation	Article no.
Carrier plate	BA721
	BA722
	BA723

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Excessive measures of neutralizing agents or basic cleaners may result in a chemical attack and/or to fading and the laser marking becoming unreadable visually or by machine for stainless steel.

Residues containing chlorine or chlorides e.g. in surgical residues, medicines, saline solutions and in the service water used for cleaning, disinfection and sterilization will cause corrosion damage (pitting, stress corrosion) and result in the destruction of stainless steel products. These must be removed by rinsing thoroughly with demineralized water and then drying.

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- Optical changes of materials, e.g. fading or discoloration of titanium or aluminum. For aluminum, the application/process solution only needs to be of pH >8 to cause visible surface changes.
- Material damage such as corrosion, cracks, fracturing, premature aging or swelling.

- ▶ Do not use metal cleaning brushes or other abrasives that would damage the product surfaces and could cause corrosion.
- ▶ Further detailed advice on hygienically safe and material-/value-preserving reprocessing can be found at www.a-k-i.org, link to Publications, Red Brochure – Proper maintenance of instruments.

6.4 Preparations at the place of use

- ▶ Remove any visible surgical residues to the extent possible with a damp, lint-free cloth.
- ▶ Transport the dry product in a sealed waste container for cleaning and disinfection within 6 hours.

6.5 Preparation before cleaning

- ▶ Carry out non-fixating/NaCl-free pre-cleaning immediately after use.

6.6 Cleaning/disinfection

Product-specific safety instructions for the reprocessing procedure

Damage to the product due to inappropriate cleaning/disinfecting agents!

- ▶ Only use cleaning/disinfecting agents approved for surface cleaning. Follow the manufacturer's instructions for the respective cleaning/disinfecting agent.

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- ▶ Use cleaning and disinfecting agents according to the manufacturer's instructions which
 - are approved for plastic material and high-grade steel,
 - do not attack softeners (e.g. in silicone).
- ▶ Observe specifications regarding concentration, temperature and exposure time.
- ▶ Do not exceed the maximum permitted cleaning temperature of 60 °C.

- ▶ Dry the product for at least 10 minutes at a maximum of 120 °C.

Note

The indicated drying temperature is a guide temperature only. It must be checked taking into account the specific conditions (e.g. load) and if applicable adjusted.

Validated cleaning and disinfection procedure

Validated procedure	Special features	Reference
Manual cleaning and wipe disinfection	■ When cleaning products with movable hinges, ensure that these are in an open position and, if applicable, move the joint while cleaning. ■ Drying phase: Use a lint-free cloth or compressed air for medical purposes ■ Ensure that the product is positioned in such a way that water will not enter the product e.g. through coupling interfaces. (Immediately remove any fluid that entered the product inadvertently.) ■ To avoid damage to cutting cylinder 2, clean the cutting cylinder with a soft brush.	Chapter Manual cleaning/disinfecting and subsection: ■ Chapter Manual cleaning and wipe disinfecting
Mechanical alkaline cleaning and thermal disinfecting	■ Install Eccos holder GB688R in a suitable wire basket (e.g. JF214R). ■ Insert the product in its proper position in Eccos holder GB682R.	Chapter Mechanical cleaning/disinfecting and subsection: ■ Chapter Mechanical alkaline cleaning and thermal disinfecting

6.7 Manual cleaning/disinfecting

- ▶ Prior to manual disinfecting, allow water to drip off for a sufficient length of time to prevent dilution of the disinfecting solution.
- ▶ After manual cleaning/disinfection, check visible surfaces visually for residues.
- ▶ Repeat the cleaning/disinfection process if necessary.

Manual cleaning and wipe disinfecting

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Cleaning	RT (cold)	-	-	D-W	-
II	Drying	RT	-	-	-	-
III	Wipe disinfection	-	>1	-	-	Meliseptol HBV wipes 50 % Propan-1-ol
IV	Final rinse	RT (cold)	0.5	-	FD-W	-
V	Drying	RT	-	-	-	-

D-W: Drinking water

FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

RT: Room temperature

Phase I

- Clean the product under running faucet water, using a suitable cleaning brush until all visible residues have been removed from the surfaces.
- Mobilize non-rigid components, such as set screws, links, etc. during cleaning.

Phase II

- Dry the product in the drying phase with suitable equipment (e.g. cloth, compressed air), see Validated cleaning and disinfection procedure.

Phase III

- Wipe all surfaces of the product with a single-use disinfectant wipe.

Phase IV

- After the specified exposure time (at least 1 min), rinse the disinfected surfaces under running FD water.
- Drain any remaining water fully.

Phase V

- Dry the product in the drying phase with suitable equipment (e.g. cloth, compressed air), see Validated cleaning and disinfection procedure.

6.8 Mechanical cleaning/disinfecting

Note

The cleaning and disinfection device must be of tested and approved effectiveness (e.g. FDA approval or CE mark according to DIN EN ISO 15883).

Note

The cleaning and disinfection device used for processing must be serviced and checked at regular intervals.

Mechanical alkaline cleaning and thermal disinfecting

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Phase	Step	T [°C/°F]	t [min]	Water quality	Chemical/Note
I	Prerinse	<25/77	3	D-W	-
II	Cleaning	55/131	10	FD-W	■ Concentrate, alkaline: - pH = 13 - <5 % anionic surfactant ■ 0.5 % working solution - pH = 11*
III	Intermediate rinse	>10/50	1	FD-W	-
IV	Thermal disinfecting	90/194	5	FD-W	-
V	Drying	-	-	-	According to the program for cleaning and disinfection device

D-W: Drinking water

FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

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- Check visible surfaces for residues after mechanical cleaning/disinfecting.

6.9 Inspection, maintenance and checks

- ▶ Allow the product to cool down to room temperature.
- ▶ Inspect the product after each cleaning and disinfecting cycle to be sure it is: clean, functional, and undamaged.
- ▶ Set aside the product if it is damaged.

6.10 Packaging

- ▶ Follow the instructions for use for the applied packaging and storage systems (e.g. instructions for use TA009721 for Aesculap Eccos storage system).
- ▶ Insert the product in its proper position in the Eccos holder, or put it on a tray in such a way that the product is protected against damage. Ensure that all cutting edges are protected.
- ▶ Pack trays appropriately for the sterilization process (e.g. in Aesculap sterile containers).
- ▶ Ensure that the packaging will prevent a recontamination of the product.

6.11 Steam sterilization

Note

The carrier plates BA721, BA722 and BA723 are gamma-sterilized single-use products. These products must not be steam-sterilized and are intended for single use only.

Note

The product may only be sterilized with the upper part 17 open and the cutting cylinder 2 removed.

- ▶ Sterilize the skin mesh dermatome immediately after cleaning.
- ▶ Store the skin mesh dermatome in a suitable wire basket (e.g. JF214R).
- ▶ Check to ensure that the sterilizing agent will come into contact with all external and internal surfaces (e.g. by opening any valves and faucets).
- ▶ Validated sterilization process
 - Disassemble the product
 - Steam sterilization through fractionated vacuum process
 - Steam sterilizer according to DIN EN 285 and validated according to DIN EN ISO 17665
 - Sterilization using fractionated vacuum process at 134 °C/holding time 5 min

- ▶ When sterilizing several products at the same time in a steam sterilizer, ensure that the maximum load capacity of the steam sterilizer specified by the manufacturer is not exceeded.

6.12 Sterilization for the US market

- Aesculap advises against sterilizing the device by flash sterilization or chemical sterilization.
- Sterilization may be accomplished by a standard prevacuum cycle in a steam autoclave.

To achieve a sterility assurance level of 10^{-6} , Aesculap recommends the following parameters:

Aesculap Orga Tray/Sterile container (perforated bottom) Minimum cycle parameters*			
Sterilization method	Temp.	Time	Minimum drying time
Prevacuum	270 °F/275 °F	4 min	20 min

*Aesculap has validated the above sterilization cycle and has the data on file. The validation was accomplished in an Aesculap sterile container cleared by FDA for the sterilization and storage of these products. Other sterilization cycles may also be suitable, however individuals or hospitals not using the recommended method are advised to validate any alternative method using appropriate laboratory techniques. Use an FDA cleared accessory to maintain sterility after processing, such as a wrap, pouch, etc.

6.13 Storage

- ▶ Store sterile products in germ-proof packaging, protected from dust, in a dry, dark, temperature-controlled area.

7. Maintenance

In order to ensure reliable operation, the product must be maintained after 300 reprocessing cycles or at least once a year.

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8. Troubleshooting list

Malfunction	Detection	Cause	Remedy
Cutting cylinder 2 fails to rotate	Foreign body in device	Cutting cylinder jammed	Remove foreign body
	–	Hand ratchet 12 not engaged	Engage hand ratchet 12
	Excessive axial slack of cutting cylinder	Cutting cylinder loose	Engage cutting cylinder and secure with lock 8
Cutting cylinder 2 stiff	Strong force required on hand ratchet 12	Defective bearing	Repair by manufacturer
Skin not cut through completely	Incorrect symbol on fine adjustment scale 18	Incorrect setting	Select the required fine adjustment, see Safe operation
	–	Screws 11 loose	Tighten screws by hand or with screwdriver 13
	Cutting tips shiny, hand ratchet 12 stiff	Cutting cylinder 2 worn	Replace cutting cylinder 2, see Safe operation
	Cutting tips broken off	Cutting cylinder 2 defective	Replace cutting cylinder 2, see Safe operation

9. Technical Service

Risk of injury and/or malfunction!

► **Do not modify the product.**

► For service and repairs, please contact your national B. Braun/Aesculap agency.

Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.

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10. Accessories/Spare parts

Art. no.	Designation
BA725R	Replacement cutting cylinder
BA721	Carrier plate factor 1.5
BA722	Carrier plate factor 3
BA723	Carrier plate factor 6
BA726R	Hand ratchet
BA727R	Basket storage aid
GB688R	Eccos holder
TA008023	Instructions for use BA720R

11. Technical data

Classification acc. to Directive 93/42/EEC

Art. no.	Designation	Class
BA720R	Skin mesh dermatome	Ila

Dimensions (H x W x D)	150 mm x 195 mm x 105 mm
Weight	4 800 g
Length of carrier plate	220 mm
Width of carrier plate	75 mm

11.1 Ambient conditions

	Operation	Storage and transport
Temperature		
Relative humidity		
Atmospheric pressure		

12. Disposal

Note

The user institution is obliged to process the product before its disposal, see Validated reprocessing procedure.

The skin mesh dermatome BA720R is made of stainless steel.

► Dispose of the device with other recyclable metals.

► Detailed information concerning the disposal of the product is available through your national B. Braun/Aesculap agency, see Technical Service.

13. Distributor in the US/Contact in Canada for product information and complaints

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

without base perforation = JK


with base perforation = JN

Lids for individual color combination

Inner lid

Outer lid

		Outside dimensions Length x Width x Height (mm)
JK440	JN440	592 x 274 x 90
JK441	JN441	592 x 274 x 120
JK442	JN442	592 x 274 x 135
JK444	JN444	592 x 274 x 187
JK446	JN446	592 x 274 x 247

	Color		Color
JK485	red	JK475A	red
JK486	blue	JK476A	blue
JK487	green	JK477A	green
JK488	gold	JK478A	gold
JK489	silver	JK479A	silver

Dimensions complete - Length x Width x Height (mm)





	Outside dimensions Inside dimensions	Outside dimensions Inside dimensions		Outside dimensions Inside dimensions	Outside dimensions Inside dimensions
Body 1/1	with inner lid	with outer lid	Body 1/1	with inner lid	with outer lid
JK440	592 x 285 x 108 544 x 258 x 75	592 x 291 x 115 544 x 258 x 75	JN440	592 x 285 x 108 544 x 258 x 63	592 x 291 x 115 544 x 258 x 63
JK441	592 x 285 x 138 544 x 258 x 105	592 x 291 x 145 544 x 258 x 105	JN441	592 x 285 x 138 544 x 258 x 93	592 x 291 x 145 544 x 258 x 93
JK442	592 x 285 x 153 544 x 258 x 120	592 x 291 x 160 544 x 258 x 120	JN442	592 x 285 x 153 544 x 258 x 108	592 x 291 x 160 544 x 258 x 108
JK444	592 x 285 x 205 544 x 258 x 172	592 x 291 x 212 544 x 258 x 172	JN444	592 x 285 x 205 544 x 258 x 160	592 x 291 x 212 544 x 258 x 160
JK446	592 x 285 x 265 544 x 258 x 232	592 x 291 x 272 544 x 258 x 232	JN446	592 x 285 x 265 544 x 258 x 220	592 x 291 x 272 544 x 258 x 220

All dimensions: Length x Width x Height (mm)


Basket **with** feet





		Outside dim. Inside dim.		Outside dim. Inside dim.		Outside dim. Inside dim.		Outside dim. Inside dim.
 36	JF111R	243 x 253 x 36 241 x 251 x 22	JF251R	406 x 253 x 36 404 x 251 x 22	JF211R	485 x 253 x 36 483 x 251 x 22	JF221R	540 x 253 x 36 538 x 251 x 22
 56	JF112R	243 x 253 x 56 241 x 251 x 42	JF252R	406 x 253 x 56 404 x 251 x 42	JF212R	485 x 253 x 56 483 x 251 x 42	JF222R	540 x 253 x 56 538 x 251 x 42
 76	JF113R	243 x 253 x 76 241 x 251 x 62	JF253R	406 x 253 x 76 404 x 251 x 62	JF213R	485 x 253 x 76 483 x 251 x 62	JF223R	540 x 253 x 76 538 x 251 x 62
 106	JF114R	243 x 253 x 106 241 x 251 x 92	JF254R	406 x 253 x 106 404 x 251 x 92	JF214R	485 x 253 x 106 483 x 251 x 92	JF224R	540 x 253 x 106 538 x 251 x 92

Instrument baskets

Basket **without** feet

		Outside dim. Inside dim.		Outside dim. Inside dim.		Outside dim. Inside dim.		Outside dim. Inside dim.
 24	JG111R	243 x 253 x 24 241 x 251 x 22	JG251R	406 x 253 x 24 404 x 251 x 22	JG211R	485 x 253 x 24 483 x 251 x 22	JG221R	540 x 253 x 24 538 x 251 x 22
 44	JG112R	243 x 253 x 44 241 x 251 x 42	JG252R	406 x 253 x 44 404 x 251 x 42	JG212R	485 x 253 x 44 483 x 251 x 42	JG222R	540 x 253 x 44 538 x 251 x 42
 64	JG113R	243 x 253 x 64 241 x 251 x 62	JG253R	406 x 253 x 64 404 x 251 x 62	JG213R	485 x 253 x 64 483 x 251 x 62	JG223R	540 x 253 x 64 538 x 251 x 62
 94	JG114R	243 x 253 x 94 241 x 251 x 92	JG254R	406 x 253 x 94 404 x 251 x 92	JG214R	485 x 253 x 94 483 x 251 x 92	JG224R	540 x 253 x 94 538 x 251 x 92

All dimensions: Length x Width x Height (mm)

Reusable filters

Standard	Outside dim.
JK090	Ø 190

Aesculap reusable filter, round for use in up to 1,000 cleaning and sterilization cycles.

Order unit (PAK): Pack of 10

Standard	Outside dim.
JK089	225 x 225

Aesculap reusable filter, square for legacy Standard containers, service life 1,000 cleaning and sterilization cycles.

Order unit (PAK): Pack of 10

Mini	Outside dim.
Dental	
JK091	95 x 215

Aesculap reusable filter for Mini and Dental containers with filter holder JK098, service life 1,000 cleaning and sterilization cycles.

Order unit (PAK): Pack of 2

Dental	Outside dim.
JK092	171 x 231

Aesculap reusable filter for Dental containers without filter holder JK098, service life 1,000 cleaning and sterilization cycles.

Order unit (PAK): Pack of 2

Marking pens

JK086	
-------	--

Marking pen for Aesculap reusable filters

Order unit (PAK): Pack of 3