ULTRAVERSE® 035

PTA Dilatation Catheter

ENGLISH

INSTRUCTIONS FOR USE

CAUTION: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

Device Description

The ULTRAVERSE® 035 PTA Dilatation Catheter is a semi-compliant balloon catheter consisting of an over the wire (OTW) catheter with an angioplasty balloon fixed at the distal tip. Two radiopaque markers delineate the working length of the balloon and aid in balloon placement. The catheter also includes a radiopaque, atraumatic tip. GEOALIGN™ Marker Bands are designated on the catheter shaft by 1cm increment bands. Each 10cm increment is labeled with the distance from the distal balloon tip. Thicker bands denote the midway point (5cm) between the labeled distances. GEOALIGN™ Marker Bands are designed to be used as a location reference tool. GEOALIGN™ Marker Bands are also designed to be used as a guide to assist with geographic alignment when used with an adjunctive therapy that utilizes the same GEOALIGN™ Marker Bands. The ULTRAVERSE® 035 PTA Dilatation Catheter is compatible with .035" guidewires. The proximal portion of the catheter includes a female luer lock hub connected to the catheter with a guidewire lumen and an inflation lumen. Packaged with every product is a protective sheath that is positioned over the balloon and should be removed prior to use. A stylet is placed into the tip of the catheter to aid in rewrap/refolding of the balloon. A re-wrapping tool is also provided on the catheter shaft. These products are not made with natural rubber latex.

Figure 1



Figure 1a
(Enhanced Graphic of Catheter Shaft from Figure 1)

NOTE: GEOALIGN™ Marker Bands are non-radiopaque and designed to be utilized outside the introducer sheath/guide catheter.

Indications For Use

The ULTRAVERSE® 035 PTA Dilatation Catheter is intended to dilate stenoses in the peripheral arteries, to treat obstructive lesions of native or synthetic AV fistulae and/or re-expand endoluminal stent graft elements in the iliac arteries. This device is also recommended for post-dilatation of balloon expandable and self-expanding stents in the peripheral vasculature. This catheter is not for use in coronary arteries.

Contraindications

None known.

Warnings

1. Contents supplied **STERILE** using ethylene oxide (EO). Non-Pyrogenic. Do not use if sterile barrier is opened or damaged. Do not reuse, reprocess or resterilize.
2. This device has been designed for single use only. Reusing this medical device bears the risk of cross-patient contamination as medical devices – particularly those with long and small lumina, joints, and/or crevices between components – are difficult or impossible to clean once body fluids or tissues with potential pyrogenic or microbial contamination have had contact with the medical device for an indeterminable period of time. The residue of biological material can promote the contamination of the device with pyrogens or microorganisms which may lead to infectious complications.
3. Do not resterilize. After resterilization, the sterility of the product is not guaranteed because of an indeterminable degree of potential pyrogenic or microbial contamination, which may lead to infectious complications. Cleaning, reprocessing, and/or resterilization of the present medical device increases the probability that the device will malfunction due to potential adverse effects on components that are influenced by thermal and/or mechanical changes.
4. To reduce the potential for vessel damage, the inflated diameter and length of the balloon should approximate the diameter and length of the vessel just proximal and distal to the stenosis.
5. When the catheter is exposed to the vascular system, the location of the balloon should be confirmed while under high quality fluoroscopic observation. Do not advance or retract the catheter unless the balloon is fully deflated. If

resistance is met during manipulation, determine the cause of the resistance before proceeding. Applying excessive force to the catheter can result in tip or catheter breakage, catheter kink, or balloon separation.

6. Do not exceed the RBP recommended for this device. Balloon rupture may occur if the RBP rating is exceeded. To prevent over pressurization, use of a pressure monitoring device is recommended.
7. After use, this product may be a potential biohazard. Handle and dispose of in accordance with acceptable medical practices and applicable local, state, and federal laws and regulations.

Precautions

1. Carefully inspect the catheter prior to use to verify that catheter has not been damaged during shipment and that its size, shape, and condition are suitable for the procedure for which it is to be used. Do not use if product damage is evident.
2. ULTRAVERSE® 035 PTA Dilatation Catheter shall only be used by physicians experienced in the performance of percutaneous transluminal angioplasty.
3. The minimal acceptable introducer sheath French size is printed on the package label. Do not attempt to pass the PTA catheter through a smaller size introducer sheath than indicated on the label.
4. Use the recommended balloon inflation medium (25% contrast medium/75% sterile saline solution). It has been shown that a 25%/75% contrast/saline ratio has yielded faster balloon inflation/deflation times. Never use air or other gaseous medium to inflate the balloon.
5. If resistance is felt during post procedure withdrawal of the catheter through the introducer sheath/guide catheter, determine if contrast is trapped in the balloon with fluoroscopy. If contrast is present, push the balloon out of the introducer sheath/guide catheter and then completely evacuate the contrast before proceeding to withdraw the balloon.
6. If resistance is still felt during post procedure withdrawal of the catheter, it is recommended to remove the balloon catheter and introducer sheath/guide catheter as a single unit.
7. Do not continue to use the balloon catheter if the catheter shaft has been bent or kinked.
8. Prior to re-insertion through the introducer sheath/guide catheter, the balloon should be wiped clean with gauze and rinsed with sterile normal saline.
9. Balloon re-wrapping should only occur while the balloon catheter is supported with a guidewire or stylet.
10. GEOALIGN™ Marker Bands are designed to be used only as an additional reference tool to accompany the interventionalist standard operation procedure.

Potential Adverse Reactions

The complications that may result from a peripheral balloon dilatation procedure include:

- Additional intervention
- Allergic reaction to drugs or contrast medium
- Aneurysm or pseudoaneurysm
- Arrhythmias
- Embolization
- Hematoma
- Hemorrhage, including bleeding at the puncture site
- Hypotension/hypertension
- Inflammation
- Occlusion
- Pain or tenderness
- Pneumothorax or hemothorax
- Sepsis/infection
- Shock
- Short term hemodynamic deterioration
- Stroke
- Thrombosis
- Vessel dissection, perforation, rupture, or spasm

DIRECTIONS FOR USE

Handling & Storage

Store in a cool, dry, dark place. Do not store near radiation or ultra-violet light sources.

Rotate inventory so that the catheters and other dated products are used prior to the "Use By" date. Do not use if packaging is damaged or opened.

Equipment For Use

- Contrast medium
- Sterile saline solution
- Luer lock syringe/inflation device with manometer (10 ml or larger)
- Appropriate introducer sheath/guide catheter and dilator set

- .035" guidewire
- Optional: Catheter stabilization device

Dilatation Catheter Preparation

Note: Remove the stylet prior to removing the protective sheath as part of the preparation for use.

1. Remove the balloon catheter from package. Verify the balloon size is suitable for the procedure and the selected accessories accommodate the catheter as labeled.
2. Remove the stylet and protective sheath by grasping the balloon catheter just proximal to the balloon and with the other hand, gently grasp the mandrel and sheath protector and slide distally off of the balloon catheter.
3. Prior to use, the air in the balloon catheter should be removed. To facilitate purging, select a syringe or inflation device with a 10 ml or larger capacity and fill approximately half of it with the appropriate balloon inflation medium (25% contrast medium/75% sterile saline solution). Do not use air or any gaseous medium to inflate the balloon.
4. Connect a stopcock to the balloon inflation female luer hub on the dilatation catheter.
5. Connect the syringe to the stopcock.
6. Hold the syringe with the nozzle pointing downward, open the stopcock and aspirate for approximately 15 seconds. Release the plunger.
7. Repeat step #6 two more times or until bubbles no longer appear during aspiration (negative pressure). Once completed, evacuate all air from the barrel of the syringe/inflation device.
8. Prepare the wire lumen of the catheter by attaching a syringe to the wire lumen hub and flushing with sterile saline solution.

Use of the ULTRAVERSE® 035 PTA Dilatation Catheter

1. Backload the distal tip of the ULTRAVERSE® 035 PTA Dilatation Catheter over the pre-positioned guidewire and advance the tip of the balloon to the introduction site.
2. Advance the catheter through the introducer sheath/guide catheter and over the wire to the site of inflation. If the stenosis cannot be crossed with the desired dilatation catheter, use a smaller diameter catheter to pre-dilate the lesion to facilitate passage of a more appropriately sized dilatation catheter.


Note: If using a hydrophilic guidewire, ensure that it is kept hydrated with sterile normal saline at all times.

3. Position the balloon relative to the lesion to be dilated, ensure the guidewire is in place, and inflate the balloon to the appropriate pressure.
4. Apply negative pressure to fully evacuate fluid from the balloon. Confirm that the balloon is fully deflated under fluoroscopy and that no contrast is left in the balloon.
5. While maintaining negative pressure and the position of the guidewire, grasp the balloon catheter just outside the introducer sheath/guide catheter and withdraw the deflated dilatation catheter over the wire through the introducer sheath/guide catheter. Use of a gentle counterclockwise motion may be used to help facilitate catheter removal through the introducer sheath/guide catheter.

Optional: To ensure that the introducer sheath/guide catheter does not move during the procedure, use a catheter stabilization device to stabilize the introducer sheath/guide catheter at the access site.

Optional: When used as a reference tool, utilize the GEOALIGN™ Marker Bands to assist with measuring how far the PTA Dilatation Catheter has advanced.

Figure 2

Balloon Reinsertion

Precaution: Do not continue to use the balloon catheter if the shaft has been bent or kinked.

Precaution: Prior to re-insertion through the introducer sheath/guide catheter, the balloon should be wiped clean with gauze, rinsed with sterile normal saline, and carefully re-wrapped. Balloon re-wrapping should only occur while the balloon catheter is supported with a guidewire or stylet.

1. Load the balloon catheter onto a guidewire.
2. Advance the balloon catheter over the pre-positioned guidewire to the introduction site and through the introducer sheath/guide catheter. If resistance is encountered, replace the previously used balloon catheter with a new balloon catheter.
3. Continue the procedure according to the "Use of the ULTRAVERSE® 035 PTA Dilatation Catheter" section herein.

Adjunctive Balloon Use Utilizing the ULTRAVERSE® 035 PTA Dilatation Catheter with GEOALIGN™ Marker Bands

When using the ULTRAVERSE® 035 PTA Dilatation Catheter with an adjunctive therapy:

1. Load the new balloon catheter onto a guidewire.
2. Advance the balloon catheter over the pre-positioned guidewire to the introduction site and through the introducer sheath. If resistance is encountered, replace the previously used balloon catheter with a new balloon.
3. To ensure geographic alignment with the pre-dilatation PTA Dilatation Catheter, advance the catheter through the introducer sheath/guide catheter until the GEOALIGN™ Marker Bands are in the same location at the introducer sheath/guide catheter hub as the initial PTA Dilatation Catheter. (See Figure 3 and 3a)
4. Continue the procedure according to the "Use of the ULTRAVERSE® 035 PTA Dilatation Catheter" section herein.

Figure 3



Figure 3a

(Enhanced Graphic of Catheter Shaft and Introducer From Figure 3)

Warning: After use, this product may be a potential biohazard. Handle and dispose of in accordance with acceptable medical practices and applicable local, state and federal laws and regulations.

WARRANTY

Bard Peripheral Vascular warrants to the first purchaser of this product, that this product will be free from defects in materials and workmanship for a period of one year from the date of first purchase and liability under this limited product warranty will be limited, to repair or replacement of the defective product, in Bard Peripheral Vascular's sole discretion, or refunding your net price paid. Wear and tear from normal use or defects resulting from misuse of this product is not covered by this limited warranty.

TO THE EXTENT ALLOWABLE BY APPLICABLE LAW, THIS LIMITED PRODUCT WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, WHETHER EXPRESS OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. IN NO EVENT WILL BARD PERIPHERAL VASCULAR BE LIABLE TO YOU FOR ANY INDIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES RESULTING FROM YOUR HANDLING OR USE OF THIS PRODUCT.

Some countries do not allow an exclusion of implied warranties, incidental or consequential damages. You may be entitled to additional remedies under the laws of your country.

An issue or revision date and revision number for these instructions are included for the user's information on the last page of this booklet.

In the event 36 months have elapsed between this date and product use, the user should contact Bard Peripheral Vascular to see if additional product information is available.

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Copyright © 2014 C. R. Bard, Inc. All rights reserved. Printed in the USA.
Copyright © 2014 C. R. Bard, Inc. Tous droits réservés. Imprimé aux États-Unis.
Copyright © 2014 C. R. Bard, Inc. Alle Rechte vorbehalten. Gedruckt in den USA.
Copyright © 2014 C. R. Bard, Inc. Tutti i diritti riservati. Stampato negli Stati Uniti.
Copyright © 2014 C. R. Bard, Inc. Todos los derechos reservados. Impreso en EE.UU.
Copyright © 2014 C. R. Bard, Inc. Alle rechten voorbehouden. Gedrukt in de VS.