




You Set the Limits



A Healthier Future for Women







Voluson Expert

You Set the Limits

Voluson[™] A Healthier Future for Women



You Set the Limits

You set out to make a difference. You became a leader.

Medicine is who you are. Women's health—your life, your legacy.

Your eye is focused on 'right now', but you are envisioning 'what's next', finding answers to questions that have yet to be asked.

And the more you know, the more you want to know.

That's why we created the Voluson Expert 22 - to help you:

REVEAL THE INVISIBLE
DETECT THE UNDETECTABLE
ACHIEVE THE UNACHIEVABLE
EXPECT THE UNEXPECTED

It's a milestone in the making.

Not the end of a path, but the continuation of a journey.

Because like you, we don't just hope for a Healthier Future for Women.

We insist on creating it.



Voluson Expert 22 – You Set the Limits

It's more than you ever thought possible. And in your hands, there's endless potential.



Imaging - Reveal the Invisible

Experience next-level imaging capabilities to help see details

Early Detection - Detect the Undetectable

Uncover critical answers sooner for faster assessment, detection, and diagnosis

Workflow - Achieve the Unachievable


Remove obstacles and drive productivity

Support Excellence - Expect the Unexpected

More services. More support. More futurefocused solutions. Finally, a partnership that won't leave you wanting more.


Voluson Expert 22


- Voluson's first ultra premium ultrasound system
- The most powerful and adaptive platform to date to help clinicians establish new standards
- Lyric Architecture generates new levels of clarity, penetration, and resolution to reveal fine anatomy in 2D/3D/4D with ease
- Specialized imaging tools help users detect and diagnose uncovering answers faster to complex and routine exams
- Modernized iconic Voluson workflow, with new personalization options, help ease your workday
- Advanced efficiency in daily exams with Artificial Intelligence and automation to reduce repetitive actions and inconsistent results


















Reveal the Invisible Imaging Excellence

Reveal the Invisible – Imaging Excellence



The Voluson Expert 22 is prepared to show you more than you ever thought possible. We are taking innovation further, making it easier to collect more information at record speeds for even faster, sharper, and clearer imaging. The Voluson Expert 22 will help you see critical details required for a confident diagnosis.

Lyric Architecture

Unlock new imaging and processing power to achieve high resolution, detailed

images

2D Imaging

281471 1860 Recline(08) BEHFI 10.20 - 2.50 Control of Fish? SER 4/CB3


Uncover answers by elevating your 2D imaging with exceptional spatial and contrast resolution for distinct border and tissue differentiation with UltraHD, Augment, and Radiant imaging

Volume Imaging



Visualizing anatomy from another view can be the difference between delaying or making a confident diagnosis. Discover the next generation of 3D/4D imaging featuring our new $HDlive^{TM}$ Studio+.

Color Imaging



Make blood flow assessment faster and easier with advanced Voluson technologies.
Increased resolution and sensitivity in color and pulsed wave Doppler for unprecedented clarity

Lyric Architecture

(ge)

How do you find a needle in a haystack? Build a powerful engine with the most progressive and adaptive capabilities yet.

The Lyric Architecture unlocks new imaging and processing power to achieve high resolution, detailed images – independent of body habitus and other difficult scanning conditions.

 Graphic-based beamforming that generates new levels of penetration, resolution, and frame rates to reveal fine anatomy in 2D/3D/4D with ease

Delivers uniformity throughout the image with increased spatial and contrast resolution

 Works in harmony with our unique probe technology offering advanced personalized imaging to address patient challenges

Opens doors to unique Voluson imaging capabilities

Generations ahead of other ultrasounds, the Lyric Architecture allows you to set new standards for years to come





The Ultimate Foundation – 2D Imaging



To uncover answers, you depend on details. The Voluson Expert 22 elevates 2D imaging with exceptional spatial and contrast resolution to enhance border distinction and tissue differentiation








Radiant - Reach New Heights




Enhance border visibility by changing elevation levels for a 3D-like appearance, which aids in crucial border differentiation in anatomy such as brain and fetal heart

Change the elevation visualization – Off – Min - Mid - Max








Off/Max Comparison

UltraHD - Visualize Super Fine Details



Obtain highly detailed images of the fetal brain, first trimester development, and other complex anatomy with increased axial and lateral resolution







On/Off Comparison

12

Augment - Simplify Challenging Exams



Utilize the Augment feature to reduce noise and increase penetration for a robust, cleaner image even in difficult to scan situations like high body mass index (BMI)



Off/On Comparison BMI >50 Off/On Comparison

Shadow Reduction



Reduces shadowing in B-mode imaging, even in difficult to scan patients



Without Shadow Reduction



With Shadow Reduction

A New Perspective – Volume Imaging



Visualizing anatomy from another view can be the difference between delaying or making a confident diagnosis. The Voluson Expert 22 builds on our world class volume imaging to deliver the next generation of 3D/4D imaging.







15

HDlive Studio+


(gg)

An essential problem-solving feature for a deeper understanding of relational anatomy and developing structures. The suite of features offer added opportunities for realistic views and assessment of anatomy for greater confidence.









HDlive

Illuminate anatomy and surrounding fluid with up to 3 independent light sources with various intensity and hues.

HDlive Silhouette

Volume imaging is taken to the next level with the ability to control what tissue type is displayed internally vs. externally.

Additional transparency tuning controls allow for more details, providing diagnostic assessment and confidence.

Perspective Rendering

"Camera Zoom Technology" type representation to enhance spatial and optical visualization for a realistic view and increased depth display.

HDlive Flow

Vascular structures are displayed with greater depth perception and dimension


HDlive Flow Silhouette

Visualize vascular structures, along with surrounding tissues, with various levels of border enhancement and through 16 transparency.

Advanced VCI



Adjusts slice thickness on 3D or 4D images to help enhance contrast resolution with use of render techniques such as bone and tissue renderings. Can be applied in the acquisition plane (VCI-A), static 3D volumes, color, or OmniView






OmniView



Obtain any plane from a 3D or 4D volume by simply drawing a line, curve, poly-line, or trace through a structure. This valuable technology enables views of even irregularly shaped structures not attainable in 2D imaging






Beyond Basic Blood Flow - Doppler Imaging



Make blood flow assessment faster and easier with advanced Voluson technologies. Increased resolution and sensitivity in color and pulsed wave Doppler for unprecedented clarity







19


Radiant flow TM



Achieve a new standard of color Doppler – Delivers easy, fast visualization of blood flow using the amplitude of the color Doppler signal to enhance robustness and create a 3D-like appearance









Slowflow



Observe low velocity vessels in 2D and 3D images with Slow*flow*HD and Slow*flow*3D – extending color capabilities for first trimester, placenta, and gynecological exams for enhanced diagnostic confidence







Slow*flow*HD Slow*flow*3D


Pulsed Wave Doppler



New levels of unprecedented sensitivity and clarity for confident assessment







Unique Probe Technologies



With the Lyric Architecture, we have taken our probe technology to the next level. Extraordinary image quality helps meet your clinical and ergonomic needs to simplify exams and your day.




eM6C - More Information in Less Time





Cutting edge technology to assist with the most complex cases while having the flexibility and comfort to integrate into daily workflow


- First commercially available curved matrix electronic 4D technology designed for women's health
- Amazing image quality based on Lyric Architecture even faster frame rates, processing power, and resolution
- Workflow enhancement features include Bi-Plane, SnapShot, eSTIC, and VCI-A




Bi-Plane Imaging



Real-Time 4D




VCI-A (Volume Contrast Imaging)



 $e {\sf STIC}$

eM6C Imaging Technologies











Bi-Plane imaging — Provides simultaneous display of high resolution, high frame rate images in two perpendicular planes. Technology may be used in 2D and color Doppler


Clinical use: Spine, palate, heart, brain, limbs, and procedures including amniocentesis, percutaneous umbilical cord blood sampling (PUBS), chorionic villus sampling (CVS); Confirmation of perpendicular (i.e., NT, fetal heart) – View potential obstructions during interventional procedures

eSTIC (Spatio-Temporal Image Correlation) enhances fetal cardiac exams with up to 75% reduction in acquisition time over traditional STIC. Reduces artifact from fetal movement. Improves B & C plane. Acquisition time reduced from 12 secs to 1.3 secs*

Clinical use: Fetal heart - Evaluate anatomy, pathology, and physiology in 3 dimensions

e4D Imaging Technologies







VCI (Volume Contrast Imaging) — Delivers excellent contrast resolution through thick slice volume of grey scale and color Doppler images.

Clinical use: Fetal brain sutures, extremities, spine, palate, soft tissue, heart in *any* plane - Added contrast resolution in soft tissue, brain, heart, etc.



Real-Time 4D - Presents exceptional visibility with an 86% increase in 4D volume rates for real-time display of motion*

Clinical use: Save time and be confident when evaluating anatomy, pathology, and/or physiology in 3 dimensions with high success rate even with moving fetus – excellent for fetal heart arrhythmias.

e4D SnapShot

(gg)

Optimize exam time with speed and ease.

- Seamlessly and quickly capture high quality 3D or eSTIC volumes without interruption in workflow
- SnapShot optimizes exam time with one button access from Real-Time 4D to acquire high resolution images







RIC6-12 — Earlier than Ever Before





Confidently explore fetal images in the earliest stages with the RIC6-12 probe

- High frequency, high resolution probe
- Delivers extremely detailed images required for visualizing and detecting anomalies in the first trimester



7-week gestational sac



10-week fetal heart



10-week embryo



Early cord visualization

RM7C – Seeing is Believing






From the first trimester to the third, our RM7C volume matrix probe helps you set new standards in performance and image quality

- XDclear[™] technology provides exceptional 2D, 3D, and 4D resolution and color sensitivity.
- High resolution matrix convex array probe
- Remarkable performance offering sharp imaging, even with difficult to scan patients

Increased user comfort - reduced weight by 80g or 17% over its RM6C predecessor




16-week fetal abdomen



18-week fetal spine



Fetal heart



SonoPelvicFloor


C2-9/C1-6 — Simply Stunning





Achieve exceptional tissue and detail resolution with our XDclear™ probes: C2-9 and C1-6. The combination of single crystal, acoustic amplification, and cool stack technology offers extraordinary 2D imaging, even with difficult to scan patients

- C1-6 Penetration and robustness covering wide range of patients, even the more challenging scans
- C2-9 High Frequency imaging for fine details and 1st trimester exams




Voluson Expert 22 |May 21, 2023 30

Shared Service Gallery





Breast cyst with RSP6-16 probe




Adult kidney shown with RM7C probe and Slow *flow*3D




Pediatric kidney with C2-9 probe



High detail pediatric brain with Radiant acquired with 6S probe



Carotid artery with 9L probe



Radiant and UltraHD offering clear tissue differentiation in adult gallbladder



Liver with Power Doppler



Abdomen with HD *live* Studio+ and Slow *flow 3*D



Detect the Undetectable Early Detection

Detect the Undetectable – Early Detection



Complex cases come with enough uncertainty. That is why the Voluson Expert 22 delivers extraordinary image quality and unique tools to help you focus on early detection and intervention.

1st Trimester




Examine developing structures earlier to uncover anomalies with high resolution probe & imaging technologies.

Fetal Heart



Increase detection, diagnose confidently, and effectively monitor the fetal heart.

Pelvic Health



Simplify reproductive health with advanced analysis tools and easy 3D technologies.

First Trimester

(gg)

Small Anatomy. Fine Details. Timely Diagnosis.

Improving detection rates of chromosomal and structural anomalies in the earliest stages of development is critical for better decision making. With small anatomy, exposing anomalies requires clear, high-resolution imaging and simplified assessment tools



UltraHD

Obtain highly detailed images of early pregnancy development with increased axial and lateral resolution for fine details



Radiant

Utilize Radiant for greater definition and clarity. Change the levels to enhance the 3D-like elevation effect for greater border visibility and sharpness.



Color Imaging


Increase resolution and sensitivity in color for unprecedented clarity with Radiant flow, Slow flow HD, and Slow flow 3D

First Trimester

(gg)

Small Anatomy. Fine Details. Timely Diagnosis.

Improving detection rates of chromosomal and structural anomalies in the earliest stages of development is critical for better decision making. With small anatomy, exposing anomalies requires clear, high-resolution imaging, and simplified assessment tools



HDlive Studio+

Essential problem-solving tools for a deeper understanding of relational anatomy and developing structures. Easily obtain volume images with unprecedented depth and clarity.



SonoNT/SonoIT

(Sonography-based Nuchal /Intracranial Translucency) – semi-automatic, standardized measurements of nuchal and intracranial translucencies



RIC6-12 Probe

Delivers extremely detailed images required for visualizing and detecting anomalies in the first trimester


Fetal Heart Evaluation


(gg)

Get to the heart of the matter

Identifying fetal cardiac abnormalities earlier means you can intervene sooner, plan for delivery, and potentially improve outcomes. The Voluson Expert 22 provides a full suite of progressive tools to help distinguish the tiniest structures with stunning clarity to detect patient answers faster.







Detection Diagnosis Monitor

Fetal Heart Detection



Utilize easy automation to help obtain and visualize the recommended fetal heart view for a complete exam.



SonoFHR

Automation tool that automatically places calipers in M-mode or pulsed wave Doppler to quickly obtain the fetal heart rate



fetalHS

Al-guided workflow to support views of the 4-chamber, 3 vessel view, and heart angle measurement (cardiac axis) to help detect fetal heart anomalies

fetalHS



An Al-guided workflow with step-by-step instructions to help identify normal anatomy with a focus on Fetal Situs, 4-Chamber Heart, 3-Vessels and Trachea View, and Cardiac Axis. Includes instructional messages, reference images, and diagrams that represent normal anatomy. 4CH, 3VV/3VT, and Cardiac Axis automatically generated from user acquired cine loops.





Fetal presentation - helps confirm proper situs of stomach and heart



Fetal spine - helps to determine fetal position and situs



Fetal stomach - helps confirm stomach location is on fetal left.



Stomach to 4-chamber sweep - From this sweep, *fetal*HS automatically locates and documents the 4-chamber view to confirm heart situs and measures a cardiac axis.



4-chamber to neck sweep - -from this sweep, *fetal*HS automatically finds and documents the 3 vessels and trachea view to confirm the size, alignment, and arrangement of the Pulmonary Artery, Aorta, SVC and Trachea to rule out anomalies

Fetal Heart Diagnosis




Find answers to your challenging exams with cutting edge tools that help provide more clarity, more speed, and more flexibility



Radiant

Utilize Radiant for greater definition and clarity. Change the levels to enhance the 3D-like elevation effect for greater border visibility and sharpness.



Radiant *flow*

Delivers easy, fast visualization of blood flow using the amplitude of the color Doppler signal to enhance robustness and create a 3D-like appearance



Slowflow

Increase resolution and sensitivity in color for unprecedented clarity with Radiant flow, Slow flow HD, and Slow flow 3D

Fetal Heart Diagnosis




Find answers to your challenging exams with cutting edge tools that help provide more clarity, more speed, and more flexibility




eM6C Bi-Plane Imaging

Provides simultaneous display of high resolution, high frame rate images in two perpendicular planes. Technology available in 2D and color Doppler modes



eM6C eSTIC

Enhances fetal cardiac exams with up to 75% reduction in acquisition time over traditional STIC and delivers improved resolution in the B and C planes*



Specialized Tools

Continuous Wave on curved probes, phased array probe for fetal and maternal cardiac imaging, and many more specialized tools

Voluson Expert 22 | May 21, 2023

^{*} As compared to volume rates on e4D vs mechanical probes (RAB6) on Voluson E10 BT20

Fetal Heart Monitoring



Not only anomalies can impact the fetal heart. Monitoring the Size, Shape, and Contractility of the fetal heart can provide additional data to help you make more informed decisions about maternal/fetal well-being and delivery planning.



Qualitative

- B-Mode
- Color Doppler
- Tissue Doppler

Quantitative

- Spectral Tissue Doppler
- TEI Index
- PW and CW Doppler
 - Valves
- Ejection Fraction



Tissue Doppler

Conduct a more comprehensive assessment of fetal cardiac function with Tissue Doppler for wall and valve motion analysis - use with pulsed wave for quantitative movement analysis

Fetal Echo Measurement Package

Dedicated measurements to include Zscores, Cardiovascular Profile Scores, and more

fetalHQ

Conduct an easy and comprehensive evaluation of the size, shape and contractility of the fetal heart from the 4chamber view using measurements based on 2D imaging and speckle tracking

fetalHQ



Conduct an easy and comprehensive evaluation of the size, shape, and contractility of the fetal heart from the 4-chamber view using measurements based on 2D imaging and speckle tracking. *fetal*HQ contains an in-depth report including Z-scores and percentiles for each of the cardiac measurements. AutoFlip and Quiver features help simplify the fetal heart orientation and identification of the endocardial border and valve insertions

Report package supports both qualitative and quantitative results to offer insights based on your preference

Qualitative

Visual representation of ventricular contractility

Quantitative

GSI for assessment of fetal heart size and shape


- Ventricular shape and contractility (assessed with speckle tracking)
- Comprehensive measurements extracted such as ejection fraction, cardiac output, TAPSE/MAPSE, etc.



Global Sphericity Index (GSI) helps assess size and shape of the fetal heart




Fetal heart contractility outlined in an in-depth report









Symptoms such as pelvic pain, post-menopausal bleeding, genitourinary dysfunction, and infertility can be confusing and concerning for your patients. Exceptional imaging and advanced analysis tools can help provide clinical insights into gynecological health.



Gynecology



Pelvic Floor



Reproductive Medicine

Uterine Trace



Taking the effort out of 3D with Uterine Trace – In 3 simple steps, obtain the coronal plane of the uterus. Links directly to Uterine Classification pictograms for easy documentation of uterine shape.

Select Uterine Trace Preset in 3D pre-mode


Trace endometrium on touch panel.
See guide area on main screen
3D sweep started automatically

Adjust OmniView line if needed

Classify uterine cavity according to ESHRE or ASRM
















Exceptional imaging and advanced analysis tools can help provide clinical insights into gynecological health



UltraHD

Obtain highly detailed images with increased axial and lateral resolution for fine details



Radiant

Utilize Radiant for greater definition and clarity. Change the levels to enhance the 3D-like elevation effect for greater border visibility and sharpness.




Color Imaging

Increase resolution and sensitivity in color for unprecedented clarity with Radiant flow, Slow flow HD, and Slow flow 3D




Exceptional imaging and advanced analysis tools can help provide clinical insights into gynecological health



Tomographic Ultrasound Imaging (TUI)

Simplifies analysis and documentation of dynamic studies with a simultaneous view of multiple slices of a volume data set



Advanced VCI with OmniView

Obtain any plane from a 3D or 4D volume by simply drawing a line, curve, poly-line, or trace through a structure. This valuable technology enables views of even irregularly shaped structures not attainable in 2D imaging




HDlive Studio+

Illuminate anatomy and surrounding fluid with up to 3 independent light sources with various intensity and hues – Perspective Rendering and additional silhouette controls to aid with making a confident diagnosis



Exceptional imaging and advanced analysis tools can help provide clinical insights into gynecological health



Uterine Classifications

Simplify identification and reporting of uterine abnormalities with **Uterine Classifications** based on ESHRE, ESGE, or ASRM guidelines

IETA/IDEA

Evaluate risk and assign scoring of the endometrial lining based on IETA (International Endometrial Tumor Analysis)

IDEA assessment tool helps guide through the protocol for quick and comprehensive evaluation of endometriosis and deep infiltrating endometriosis (DIE). Can be used with and without Scan Assistant



IOTA

Ovarian cancer risk estimate based on IOTA (International Ovarian Tumor Analysis)

Simple Rules IOTA LR2 IOTA ADNEX

First ultrasound supplier to offer all 3 IOTA models validated and approved by the IOTA group




Exceptional imaging and advanced analysis tools can help provide clinical insights into gynecological health

SonoPelvicFloor


- Al-based feature to help obtain pelvic floor measurements for more consistent and reproducible results
- Automates plane alignment of valsalva
- Automated measurement of levator hiatus (a-p, lateral and area/circumference)
- Comprehensive pelvic floor measurement package








Exceptional imaging and advanced analysis tools can help provide clinical insights into gynecological health



SonoAVC™ follicle

(Sonography-based Automated Volume Count follicle) – Automatically calculates the number, dimensions, and volume of hypoechoic structures in a volume sweep to help monitor patient follicles faster



SonoAVCantral^{2.0}

More accurately evaluate ovarian reserve and increase success for Assisted Reproductive Medicine treatments



Achieve the Unachievable Workflow Efficiency

Achieve the Unachievable - Workflow Efficiency

Modern, yet familiar, the Voluson Expert 22 is our most intuitive and customizable ultra-premium ultrasound. Designed with a new generation of Artificial Intelligence and automation tools, it helps increase exam speed and accuracy to streamline your day. Just think of the potential.

Design



Sleek and modern with advanced ergonomics and familiar user interface

Artificial Intelligence & Automation



Designed with a new generation of Artificial Intelligence to help increase exam speed and accuracy to streamline your day.

Personalization & Customization



Personalize the user interface for your unique preferences and select system colors to suit your mood

A Modern Classic Design

Incorporating the same intuitive functionality, Voluson's **familiar workflow** is simple and seamless

Work comfortably and see finer details with ease – 23.8" **High-definition ultrasound (HDU)** display offer 3 imaging sizes

Probe ports conveniently located for less reach and easier access

Gel warmer keeping gel at a comfortable temperature




Achieve the Unachievable



Artificial Intelligence (AI) has become part of our everyday lives – reducing tasks and simplifying our daily routines. It is also beginning to play an important role in women's health – helping to improve the delivery of care and efficiency of ultrasound exams.







SonoLyst

SonoPelvicFloor

fetalHS

The Voluson™ Expert 22 contains AI tools to help you be more efficient, allowing you to reduce repetitive tasks and focus on your patients and practice needs.



Optimized Workflow

Respond probe activation that automatically selects used probe and goes to live – **New & unique in the industry!**

Flow Profiles - predefined and optimized preset tool for color and pulsed wave Doppler modes for more efficient imaging

RFID – effortlessly incorporating security into daily workflow to log on







Personalization & Customization



Customize the user interface >23M options

- Large touch panel (individual button arrangements)
- Tailored hard keys

Personalize colors & lights >4000 combinations

- Trackball colors
- User interface illumination
- Ambient trolley lighting
- Touch panel color scheme





SonoLyst



SonoLyst is your virtual, on-board expert utilizing the power of AI to identify fetal anatomy seen on standard views and enhancing efficiency by adding annotations and measurements. For quality assurance, SonoLyst can also be used to compare the image or view acquired to standard criteria, ensuring exam quality and consistency.

SonoLystlR: Simply scan, then freeze and SonoLystlR (Image Recognition) does the rest. Using ISUOG practice guidelines for performance of the routine mid-trimester fetal ultrasound scan, SonoLystlR automatically detects the 21 recommended views. It identifies anatomy then selects all applicable annotations and measurements. Confirm, and data is entered into the Scan Assistant checklist and report, enhancing workflow and reducing variability between operators for improved consistency.

SonoLystX: Build and refine your skills with SonoLystX, your virtual, on-board ultrasound expert. Using AI, the system compares the image or view acquired to standard criteria accepted by experts to ensure it meets clinical standards. SonoLystX can help enhance accuracy and quality with anatomy diagrams plus the ability to insert image examples. Ideal for teaching and training, progress can be monitored for quality assurance to ensure the highest quality imaging standards and consistency.








SonoPelvicFloor



- Al based feature to help obtain pelvic floor measurements for more consistent and reproducible results
- Automates plane alignment in Valsalva
- Automated measurement of levator hiatus (a-p, lateral and area/circumference)
- Comprehensive pelvic floor measurement package







Voluson Expert 22 |

fetalHS



An Al-based guided workflow with step-by-step instructions to help identify normal anatomy with a focus on Fetal Situs, 4-Chamber Heart, 3-Vessels and Trachea View, and Cardiac Axis. Includes instructional messages, reference images, and diagrams that represent normal anatomy. 4CH, 3VV/3VT and Cardiac Axis automatically generated from user acquired cine loops.






SonoCNS



Al tool that supports consistent measurements and helps drive workflow efficiency by helping align and display recommended views and measurements of the fetal brain from a 3D volume





Achieve the Unachievable



Remove obstacles and drive productivity to transform your day with the Voluson Expert 22 Sono-Automation tools - so you can keep pushing what is possible.



SonoBiometry

Helps reduce keystrokes and improves exam time with automatic caliper placement and fetal measurements of bi-parietal diameter (BPD), head circumference (HC), abdominal circumference (AC), femur length (FL), humerus length (HL), lateral ventricle, cerebellum, and cisterna magna.



SonoNT/SonoIT

(Sonography-based Nuchal /Intracranial Translucency) Semi-automatic, standardized measurements of nuchal and intracranial translucencies

Achieve the Unachievable



Remove obstacles and drive productivity to transform your day with the Voluson Expert 22 Sono-Automation tools - so you can keep pushing what is possible.



SonoFHR

Automation tool that automatically places calipers in M-mode or pulsed wave Doppler to quickly obtain the fetal heart rate



SonoRender live

Simplifies volume rendering by automating render line placement in 3D and 4D imaging

Flow Profiles



Predefined and optimized preset tool for color and pulsed wave Doppler modes for more efficient imaging



Scan Assistant




Balance imaging schedules while increasing patient satisfaction with complete exams and less callbacks

These flexible, customizable exam protocols help reduce time required to conduct and document results by guiding you through an exam more efficiently, aiding in annotation, measuring, reporting, and transferring data to an image management system or PACS based system on your order sequence and output requirements

- · Status monitor onscreen
- Anatomy page completion during exam
- Comparison Image option

- Improved configuration page
- Image shuffling
- Link probes & presets to items

Improved efficiency with Scan Assistant offline configurator



Patient scan time reduced by as much as 45% with Scan Assistant 2.0 and ViewPoint, saving an average of 7 min x exam*



Expect the UnexpectedSupport Excellence



Leverage our entire ecosystem of support for peace of mind – now and in the future. Combined with a flexible service program,. Experience more. More services. More support. More future-focused solutions. Finally, a partnership that won't leave you wanting more.

Education & Training



Take advantage of clinical and technical education, along with training support, to help you provide the highest quality of care and achieve operational excellence.

Collaboration



Digital Expert Connect offers real-time collaboration opportunities where clinicians can easily consult with collaborators within their network anytime, anywhere

Maintenance, Repair & Remote Support Services



Manage your total cost of ownership while improving system uptime with maintenance contracts and features like OnWatch proactive monitoring with InSite™ connectivity.



Leverage our entire ecosystem of support for peace of mind – now and in the future. Combined with a flexible service program,. Experience more. More services. More support. More future-focused solutions. Finally, a partnership that won't leave you wanting more.

Updates, Upgrades, and Device Protection



Download critical Windows® patches, security updates and Voluson extension release software at your convenience with **eDelivery***. Protect your ultrasound system from cyber-attacks with Voluson's comprehensive **SonoDefense** solution.

Probe Repair and Hygiene Solutions



A continuum of support for your probes, including disinfection and hygiene solutions. Stay a step ahead with **Probe Check**, which assesses the health of your probes and can alert you to potential issues.

Performance Optimization Solutions



Maximize the performance, workflow and utilization of your Voluson ultrasound system using Imaging Insights, iCenter Performance Management tools and other digital solutions.



Leverage our entire ecosystem of support for peace of mind – now and in the future. Combined with a flexible service program,. Experience more. More services. More support. More future-focused solutions. Finally, a partnership that won't leave you wanting more.

Centralized Device Management



AVURI is a cloud-based solution that lets you manage devices remotely from one central location. Monitor devices across your fleet with the ability to view, back up, and deploy device configurations, all without leaving your desk.*

Voluson Club



Like family, we are here for you. The only ultrasound community dedicated to the education and collaboration of women's health practitioners to help expand your expertise and your practice.



Leverage our entire ecosystem of support for peace of mind – now and in the future. Combined with a flexible service program. Experience more. More services. More support. More future-focused solutions. Finally, a partnership that won't leave you wanting more.

ViewPoint™ 6



ViewPoint 6 is a dedicated ultrasound reporting and image management solution with enhanced features to help improve workflow.

Tricefy ™



Tricefy inside – Unique cloud- based solution to share, collaborate and archive ultrasound images and reports securely with colleagues and patients directly from your Voluson.

4D View



Optimize, manipulate and analyze volume ultrasound data offline – when and where you need it. 4DView extends the capabilities of your Voluson and gives you the flexibility to free up your ultrasound machine for additional patients.

Experience more. More services. More support. More future-focused solutions. Finally, a partnership that won't leave you wanting more





© GE, 2022

GE Healthcare reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your GE Healthcare representative for the most current information. GE, the GE Monogram, Voluson, Radiant flow, HD live, InSite and ViewPoint are trademarks of GE. SonoLyst incorporates the AI technology of Intelligent Ultrasound. Tricefy™ is a registered trademark of Trice Imaging, Inc. Windows is a registered trademarks of Microsoft Corporation. GE Healthcare, a division of GE. GE Medical Systems, Inc., doing business as GE Healthcare.

May 2022 JB18772XX