# Total T3

For use on IMMULITE® 2000 systems

**SIEMENS**

# IMMULITE® 2000 Total T3

---

## English

---

**Intended Use:** For *in vitro* diagnostic use with the IMMULITE® 2000 Systems Analyzers — for the quantitative measurement of total circulating triiodothyronin (T3) in serum. Measurements of T3 are used in the diagnosis and treatment of thyroid disease.

Catalog Numbers:  
**L2KT32** (200 tests)    **L2KT36** (600 tests)

Test Code: **T3**    Color: **Violet**

### Summary and Explanation

Under normal physiological conditions, T3 represents approximately 5 percent of the thyroid hormones in serum. Although present in lower concentration, T3 has a greater intrinsic metabolic activity, faster turnover and larger volume of distribution than circulating T4.<sup>1</sup> Reports that thyrotoxicosis may be caused by abnormally high concentrations of T3 — rather than T4 — have reinforced the importance of T3 measurements.<sup>2-5</sup> In addition, T3 determination is an important tool for monitoring hypothyroid patients receiving sodium liothyronine therapy. Unlike "T3 Uptake" tests, which estimate the saturation of thyroid hormone binding proteins, total T3 analysis actually measures circulating levels of triiodothyronine. Most reports indicate that T3 levels distinguish clearly between euthyroid and hyperthyroid subjects, but provide a less clear-cut separation between hypothyroid and euthyroid subjects.<sup>1</sup>

Numerous conditions unrelated to thyroid disease may cause abnormal T3 values.<sup>1,6-9</sup> Consequently, total T3 values should not be used on their own in establishing the thyroid status of an individual. The levels of serum T4, TBG (thyroid binding globulin), TSH (thyroid stimulating hormone) and other clinical findings must be considered as well.

### Principle of the Procedure

IMMULITE 2000 Total T3 is a solid-phase, competitive chemiluminescent enzyme immunoassay.

**Incubation Cycles:** 1 × 30 minutes

### Specimen Collection

The use of an ultracentrifuge is recommended to clear lipemic samples.

Hemolyzed samples may indicate mistreatment of a specimen before receipt by the laboratory; hence the results should be interpreted with caution.

EDTA plasma is not recommended for use.

Centrifuging serum samples before a complete clot forms may result in the presence of fibrin. To prevent erroneous results due to the presence of fibrin, ensure that complete clot formation has taken place prior to centrifugation of samples. Some samples, particularly those from patients receiving anticoagulant therapy, may require increased clotting time.

Blood collection tubes from different manufacturers may yield differing values, depending on materials and additives, including gel or physical barriers, clot activators and/or anticoagulants.

IMMULITE 2000 Total T3 has not been tested with all possible variations of tube types. Consult the section on Alternate Sample Types for details on tubes that have been tested.

**Volume Required:** 25 µL serum

**Storage:** 7 days at 2–8°C, or 2 months at –20°C.

## Warnings and Precautions

For *in vitro* diagnostic use.

### CAUTION! POTENTIAL BIOHAZARD

Contains human source material. Each donation of human blood or blood component was tested by FDA-approved methods for the presence of antibodies to human immunodeficiency virus type 1 (HIV-1) and type 2 (HIV-2) as well as for hepatitis B surface antigen (HBsAg) and antibody to hepatitis C virus (HCV). The test results were negative (not repeatedly reactive). No test offers complete assurance that these or other infectious agents are absent; this material should be handled using good laboratory practices and universal precautions.<sup>8-10</sup>

**CAUTION:** This device contains material of animal origin and should be handled as a potential carrier and transmitter of disease.

**Reagents:** Store at 2–8°C. Dispose of in accordance with applicable laws.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing.

**Chemiluminescent Substrate:** Avoid contamination and exposure to direct sunlight. (See insert.)

**Water:** Use distilled or deionized water.

## Materials Supplied

Components are a matched set. Labels on the inside box are needed for the assay.

### Total T3 Bead Pack (L2T312)

With barcode. 200 beads, coated with monoclonal murine anti-T3. Stable at 2–8°C until expiration date.

**L2KT32:** 1 pack

**L2KT36:** 3 packs

### Total T3 Reagent Wedge (L2T3A2)

With barcode. 11.5 mL alkaline phosphatase (bovine calf intestine) conjugated to T3, with preservative. Stable at 2–8°C until expiration date.

**L2KT32:** 1 wedge

**L2KT36:** 3 wedges

Before use, tear off the top of the label at the perforations, without damaging the barcode. Remove the foil seal from the top of wedge; snap the sliding cover down into the ramps on the reagent lid.

### Total T3 Adjustors (LT3L, LT3H)

Two vials (Low and High), 2.0 mL each, of T3 in processed human serum, with preservative. Stable at 2–8°C for 30 days after opening, or for 6 months (aliquotted) at –20°C.

**L2KT32:** 1 set

**L2KT36:** 2 sets

Before making an adjustment, place the appropriate Aliquot Labels (supplied with the kit) on test tubes so that the barcodes can be read by the on-board reader.

## Kit Components Supplied Separately

### Multi-Diluent 1

For the on-board dilution of high samples. One vial of concentrated (ready-to-use) processed, normal human serum, containing undetectable to low levels of total T3 with preservative. Stable at 2–8°C for 30 days after opening, or for 6 months (aliquotted) at –20°C.

**L2M1Z:** 25 mL

Barcode labels are provided for use with the diluent. Before use, place an appropriate label on a 16 × 100 mm test tube, so that the barcodes can be read by the on-board reader.

**L2M1Z:** 3 labels

**L2SUBM:** Chemiluminescent Substrate

**L2PWSM:** Probe Wash

**L2KPM:** Probe Cleaning Kit

**LRXT:** Reaction Tubes (disposable)

**L2ZT:** 250 Sample Diluent Test Tubes (16 × 100 mm)

**L2ZC:** 250 Sample Diluent Tube Caps

Also Required

Distilled or deionized water; test tubes; controls.

## Assay Procedure

Note that for optimal performance, it is important to perform all routine maintenance procedures as defined in the IMMULITE 2000 Systems Operator's Manual.

See the IMMULITE 2000 Systems Operator's Manual for preparation, setup, dilutions, adjustment, assay and quality control procedures.

**Recommended Adjustment Interval:**  
2 weeks

**Quality Control Samples:** Follow government regulations or accreditation requirements for quality control frequency.

Use controls or serum pools with at least two levels (low and high) of T3.

Siemens Healthcare Diagnostics recommends the use of commercially available quality control materials with at least 2 levels (low and high). A satisfactory level of performance is achieved when the analyte values obtained are within the Acceptable Control Range for the system, or within an established range determined by an appropriate internal laboratory quality control scheme.

## Expected Values

**Adult:** A total of 169 serum samples from apparently healthy subjects were processed by the IMMULITE 2000 Total T3 assay. Non-parametric analysis yielded a median value of 119 ng/dL and a central 95% range of

84–172 ng/dL (1.3–2.6 nmol/L)

Reference ranges were generated using serum samples collected in Becton Dickinson plastic vacutainer tubes without gel barriers.

**Pediatric:** Reference intervals for the pediatric population (children and adolescents) were established for the IMMULITE Total T3 assay in accordance with CLSI guideline EP28-A3C.<sup>11</sup>

For analysis of data, the population was divided into three age subgroups:

- Infants: subjects aged 1–23 months
- Children: subjects aged 2–12 years
- Adolescents: subjects aged 13–20 years

A non-parametric approach was used to establish the reference intervals for children and adolescents where the 2.5 and 97.5 percentiles of the distribution of values were calculated. For the infant population, a robust measure of location and spread, as developed by Horn and Pesce, was used to estimate the 2.5 and the 97.5 percentile reference intervals, accommodating the smaller sample size.<sup>11–13</sup>

The reference intervals detailed by age group and number of samples are presented in the Reference Intervals table.

### IMMULITE 2000 Total T3 Pediatric Reference Intervals

Age Group	<i>n</i>	Conventional ng/dL	SI nmol/L
Infants (1–23 Months)	90	116–241	1.8–3.7
Children (2–12 Years)	197	109–206	1.7–3.2
Adolescents (13–20 Years)	148	93–170	1.4–2.6

Consider these limits as *guidelines* only. Each laboratory should establish its own reference ranges.

## Limitations

Many physiological, pharmacological, pathological, and genetic factors affect the interpretation of total T3 results. (Refer to the Summary and Explanation section.)

More than 20% of patients with non-thyroidal disease who are critically ill have low serum total T3 levels.

Heterophilic antibodies in human serum can react with the immunoglobulins included in the assay components causing interference with *in vitro* immunoassays. [See Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. *Clin Chem* 1988;34: 27-33.] Samples from patients routinely exposed to animals or animal serum products can demonstrate this type of interference potentially causing an anomalous result. These reagents have been formulated to minimize the risk of interference; however, potential interactions between rare sera and test components can occur. For diagnostic purposes, the results obtained from this assay should always be used in

combination with the clinical examination, patient medical history, and other findings.

## Performance Data

See Tables and Graphs for data *representative* of the assay's performance. Results are expressed in ng/dL. (Unless otherwise noted, all were generated on serum samples collected in tubes without gel barriers or clot-promoting additives.)

### Conversion Factor:

ng/dL  $\times$  0.01536  $\rightarrow$  nmol/L

**Calibration Range:** 40 to 600 ng/dL  
(0.61 to 9.2 nmol/L)

The assay is traceable to an internal standard manufactured using qualified materials and measurement procedures.

**Analytical Sensitivity:** 19 ng/dL  
(0.29 nmol/L)

**Precision:** Samples were processed in duplicate over the course of 20 days, two runs per day, for a total of 40 runs and 80 replicates. (See "Precision" table.)

**Linearity:** Samples were assayed under various dilutions. (See "Linearity" table for representative data.)

**Recovery:** Samples spiked 1 to 19 with three T3 solutions (1660, 2980 and 3840 ng/dL), were assayed. (See "Recovery" table for representative data.)

**Specificity:** The antibody is highly specific for T3. (See "Specificity" table.)

**Bilirubin:** Presence of bilirubin in concentrations up to 200 mg/L has no effect on results, within the precision of the assay.

**Biotin:** Specimens that contain biotin at a concentration of 1500 ng/mL demonstrate a less than or equal to 10% change in results. Biotin concentrations greater than this may lead to incorrect results for patient samples.

**Hemolysis:** Presence of packed red blood cells in concentrations up to 30  $\mu$ L/mL has no effect on results, within the precision of the assay.

**Lipemia:** Presence of triglycerides in concentrations up to 5000 mg/dL has no effect on results, within the precision of the assay.

**Alternate Sample Type:** To assess the effect of alternate sample types, blood

was collected from 25 volunteers into plain, heparinized and Becton Dickinson SST<sup>®</sup> vacutainer tubes. Equal volumes of the matched samples were spiked with various concentrations of total T3, to obtain values throughout the calibration range of the assay, and then assayed by the IMMULITE 2000 Total T3 procedure.

(Heparin) = 1.06 (Serum) – 6.66 ng/dL  
 $r = 0.993$

(SST) = 1.06 (Plain Tubes) – 4.09 ng/dL  
 $r = 0.994$

Means:

186 ng/dL (Heparin)

183 ng/dL (SST)

179 ng/dL (Serum)

**Method Comparison:** The assay was compared to IMMULITE Total T3 on 340 samples. (Concentration range: approximately 70 to 540 ng/dL. See graph.) By linear regression:

(IML 2000) = 0.91 (IML) + 13.6 ng/dL  
 $r = 0.987$

Means:

135 ng/dL (IMMULITE 2000)

133 ng/dL (IMMULITE)

## References

- 1) Hollander CS, Shenkman L. Radioimmunoassays for triiodothyronine and thyroxine. In: Rothfeld B, editor. Nuclear medicine *in vitro*. Philadelphia: Lippincott, 1974:136-49.
- 2) Ahmed M, Doe RP, Nuttall FQ. Triiodothyronine thyrotoxicosis following iodide ingestion: a case report. *J Clin Endocrinol Metab* 1974;38:574-6.
- 3) Hollander CS, Nihei N, Burday SZ, Mitsuma T, Shenkman L, Blum M. Clinical and laboratory observations in cases of triiodothyronine toxicosis confirmed by radioimmunoassay. *Lancet* 1972;i:609-11.
- 4) Mitsuma T, Owens R, Shenkman L, Reiter E, Hollander CS. T3 toxicosis in childhood: hyperthyroidism due to isolated hypersecretion of triiodothyronine. *J Pediatr* 1972;81:982-4.
- 5) Sterling K, Refetoff S, Selenkow HA. Thyrotoxicosis due to elevated serum triiodothyronine levels. *JAMA* 1970;213:571-5.
- 6) Bates HM. *Clin Lab Prod* 1974;3:16.
- 7) Utiger RD. Serum triiodothyronine in man. *Annu Rev Med* 1974;25:289-302.
- 8) Larson PR. Triiodothyronine: review of recent studies of its physiology and pathophysiology in man. *Metabolism* 1972;21:1073-92.
- 9) Oppenheimer JH. Role of plasma proteins in the binding, distribution and metabolism of the thyroid hormones. *N Engl J Med* 1968;278:1153-62.
- 10) Centers for Disease Control. Update: Universal precautions for prevention of transmission of human immunodeficiency virus, hepatitis B virus and other bloodborne

pathogens in healthcare settings. MMWR, 1988;37:377–82, 387–8. 9) Clinical and Laboratory Standards Institute (formerly NCCLS). Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline - Third Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2005. NCCLS Document M29-A3. 10) Federal Occupational Safety and Health Administration, Bloodborne Pathogens Standard, 29 CFR 1910.1030. 11) Clinical and Laboratory Standards Institute. *Defining, Establishing, and Verifying Reference Intervals in the Clinical Laboratory; Approved Guideline—Third Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2010. CLSI Guideline EP28-A3C. 12) Horn PS, Pesce AJ, Reference Intervals. A User's Guide, Washington, DC:AACC Press; 2005. 13) Reed AH, Henry RJ, Mason WB. Influence of statistical method used on the resulting estimate of normal range. 1971; 17:275-284.

## Technical Assistance

In the United States, contact Siemens Healthcare Diagnostics Technical Services department. Tel: 877.229.3711. Outside the United States, contact your National Distributor.

[www.siemens.com/diagnostics](http://www.siemens.com/diagnostics)

The Quality System of Siemens Healthcare Diagnostics Products Ltd. is certified to ISO 13485.

## Tables and Graphs

### Precision (ng/dL)

	Mean <sup>3</sup>	Within-Run <sup>1</sup>		Total <sup>2</sup>	
		SD <sup>4</sup>	CV <sup>5</sup>	SD	CV
1	73	8.8	12%	11	15%
2	88	6.7	7.6%	7.4	8.6%
3	96	6.7	7.0%	9.8	10%
4	105	7.4	7.1%	12	11%
5	171	9.4	5.5%	12	7.0%
6	362	16	4.4%	19	5.3%

### Linearity (ng/dL)

	Dilution <sup>1</sup>	Observed <sup>2</sup>	Expected <sup>3</sup>	%O/E <sup>4</sup>
1	8 in 8 <sup>5</sup>	340	—	—
	4 in 8	187	170	110%
	2 in 8	86	85	101%
	1 in 8	41	43	95%
2	8 in 8	380	—	—
	4 in 8	166	190	87%
	2 in 8	84	95	88%
	1 in 8	40	48	83%
3	8 in 8	529	—	—
	4 in 8	277	265	105%
	2 in 8	131	132	99%
	1 in 8	75	66	114%

### Recovery (ng/dL)


	Solution <sup>1</sup>	Observed <sup>2</sup>	Expected <sup>3</sup>	%O/E <sup>4</sup>
1	—	127	—	—
	A	213	204	104%
	B	274	270	102%
	C	325	313	104%
2	—	264	—	—
	A	335	334	100%
	B	395	400	99%
	C	456	443	103%

### Specificity

Compound <sup>1</sup>	Amount Added <sup>2</sup> ng/dL	% Cross-reactivity <sup>3</sup>
Triiodo-L-thyronine (T3)	—	100%
Reverse T3	500,000	0.02%
L-Thyroxine (T4)	24,000	0.60%
D-Thyroxine	5000	1.3%
Triiodothyroacetic acid	10,000	0.7%
	1000	0.6%
Tetraiodothyroacetic acid	20,000	0.37%
Monoiodotyrosine	1,000,000	0.005%
Diiodo-L-tyrosine	1,000,000	ND
Methimazole	1,000,000	0.008%
5,5'-Diphenylhydantoin	1,000,000	0.004%
Phenylbutazone	1,000,000	ND
6- <i>n</i> -Propyl-2-thiouracil	1,000,000	ND

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IMMULITE is a trademark of Siemens Healthcare Diagnostics.
© 2018 Siemens Healthcare Diagnostics. All rights reserved.