Medtronic

Amphirion™ Deep

PTA Catheter

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
English

1 Description

The Amphirion Deep PTA catheter is dedicated to percutaneous transluminal angioplasty (PTA) of small peripheral arteries.

The Amphirion Deep PTA catheter is an over-the-wire catheter consisting of a dual lumen shaft ending proximally with a Y connector (hub) and distally with a balloon close to the catheter tip. The straight port of the Y connector is the guide wire entrance and the side port is used to inflate and deflate the balloon. Both the two lumens run all through the entire shaft length. The guide wire lumen permits the use of guide wires to facilitate advancement of the catheter to and through the stenosis to be dilated and it ends at the tip of the catheter. Maximum guide wire diameter is 0.014 in (0.36 mm).

A needle with a luer port to facilitate the insertion of the 0.014 in (0.36 mm) guidewire and a clip to aid catheter handling and holding on the sterile field are provided in the sterile packaging. The usable catheter length is 120 cm and 150 cm. The balloon is designed to reach specific diameters at specific pressures (see compliance chart included in the package). In order to correctly position the balloon under fluoroscopy, radiopaque markers are placed on the shaft under the balloon itself, defining its cylindrical area. The catheter includes a smooth, soft and atraumatic tip to facilitate advancement of the catheter through the stenosis. In order to facilitate catheter advancing through the vascular district and the vessel stenosis, a LFC hydrophilic coating is present on the distal shaft (23 cm to 25 cm) and the balloon.

The Amphirion Deep PTA catheter is available in different balloon sizes. Nominal balloon diameters and lengths are printed on the hub.

2 Intended purpose

The intended purpose of the Amphirion Deep PTA catheter is to restore lumen patency and blood flow to peripheral arteries, such as the femoral artery and below the knee vessels.

2.1 Intended patient population

The Amphirion Deep PTA catheter is intended for patients in need of percutaneous transluminal angioplasty to improve and maintain luminal diameter in peripheral arteries, such as the femoral artery and below the knee vessels.

2.2 Indications for use

The Amphirion Deep PTA catheter is specially designed for the percutaneous transluminal angioplasty of vascular districts as femoral artery and below the knee vessels with a reference diameter equal or larger than selected balloon sizes.

2.3 Contraindications

Inability to cross lesion with a guidewire. Amphirion Deep PTA catheter is NOT indicated for use in coronary arteries.

2.4 Clinical benefits

The clinical benefits of restoring lumen patency and blood flow with the Amphirion Deep PTA catheter in the treatment of occlusions or lesions in peripheral arteries include: resolution of arterial disease related symptoms, improvement in quality of life, and increased mobility. The Amphirion Deep PTA catheter has been found to be safe and effective with clinical benefits outweighing the risks in the intended population.

Table 1. Amphirion Deep PTA catheter clinical benefits

Data endpoints and published studies	
Freedom from amputation (limb salvage) at 12 months	Rate (n/N) [95% CI]
Gandini 2009 ^a	100.0% (31/31) [88.8, 100.0]
Jia 2021 ^b	98.2% (55/56) [90.4, 100.0]
Liistro 2013 ^c	98.4% (63/64) [91.6, 100.0]
Wang 2009 ^d	94.1% (32/34) [80.3, 99.3]
	Kaplan-Meier Estimate (N) [SE]
Deloose 2009 ^e	91.5% (N=35) [4.6%]
Troisi 2016 ^f	91.1% (N=24) [4.9%]
Freedom from TLR at 12 months	Rate (n/N) [95% CI]
Gandini 2009	93.5% (29/31) [78.6, 99.2]
Jia 2021	76.8% (43/56) [65.3, 87.7]
Pre-procedure vs post-procedure maximum walking distance	Mean ± SD (N) [95% CI]
Wei 2014 ^g (pre-procedure)	168.7 ± 87.3 m (N=30) [137.5, 199.9]
Wei 2014 (post-procedure)	654.3 ± 278.3 m (N=30) [554.7, 753.9]
Pre-procedure vs post-procedure Rutherford classification	Mean ± SD (N) [95% CI]
Wei 2014 (pre-procedure)	3.1 ± 2.0 (N=71) [2.6, 3.6]
Wei 2014 (post-procedure)	0.3 ± 1.1 (N=71) [0.04, 0.56]
Primary patency of peripheral artery at 12 months	Kaplan-Meier Estimate (N) [SE]
Deloose 2009	52.9% (N=24) [7.4%]
Resolution of pain symptoms post-procedure	Rate (n/N) [95% CI]
Troisi 2016	77.1% (27/35) [59.9, 89.6]
Wang 2009	73.5% (25/34) [55.6, 87.1]
Complete wound healing at 12 months	Rate (n/N) [95% CI]
Jia 2021	75.0% (24/32) [55.6, 88.5]
Liistro 2013	67.1% (43/64) [54.3, 78.4]
Troisi 2016	77.1% (27/35) [59.9, 89.6]
Wei 2014	100.0% (12/12) [73.5, 100.0]
Time to complete wound healing	Mean ± SD (N) [95% CI]
Liistro 2013	5.2 ± 1.6 months (N=43) [4.7, 5.7]
Wei 2014	3.2 ± 1.2 months (N=12) [2.5, 3.9]

^a Gandini R, Volpi T, Pampana E, Uccioli L, Versaci F, Simonetti G. Applicability and clinical results of percutaneous transluminal angioplasty with a novel, long, conically shaped balloon dedicated for below-the-knee interventions. *J Cardiovasc Surg (Torino)*. 2009;50(3): 365-371.

^b Jia X, Zhuang B, Wang F, et al. Drug-Coated Balloon Angioplasty Compared With Uncoated Balloons in the Treatment of Infrapopliteal

Artery Lesions (AcoArt II-BTK). *J Endovasc Ther*. 2021;28(2): 215-221.

^c Liistro F, Porto I, Angioli P, et al. Drug-eluting balloon in peripheral intervention for below the knee angioplasty evaluation (DEBATE-BTK): a randomized trial in diabetic patients with critical limb ischemia. *Circulation*. 2013;128(6):615-621.

^d Wang J, Zhu YQ, Zhao JG, et al. Infrapopliteal angioplasty with a long over-the-wire (OTW) balloon in the treatment of severe limb ischemia in diabetic patients: a retrospective study. *Acta Radiol*. 2009;50(4): 360-367.

^e Deloose K, Bosiers M, Peeters P. One year outcome after primary stenting of infrapopliteal lesions with the CHROMIS DEEP stent in the management of critical limb ischaemia. *EuroIntervention*. 2009;5(3): 318-324.

^f Troisi N, Ercolini L, Chisci E, et al. Use of Tapered Balloons to Recanalize Occluded Below-the-Knee Arteries in Diabetic Patients with Critical Limb Ischemia. *Ann Vasc Surg*. 2016;31:105-110.

^g Wei LM, Zhu YQ, Zhao JG, Wang J, Lu HT, Zhang PL. Retrograde transplatar arch angioplasty of below-the-knee arterial occlusions: outcomes compared to antegrade recanalization. *Acad Radiol*. 2014;21(11):1475-1482.

2.5 Intended users

Interventional procedures using the Amphirion Deep PTA catheter should only be performed by physicians who have experience with interventional techniques in the vascular system.

This device should only be used in a clinical setting in a sterile environment.

2.6 Device performance characteristics

The Amphirion Deep PTA catheter is an over the wire catheter with a balloon bonded onto the distal tip. The overall design and material selection allow for treatment of occlusions or lesions in the peripheral arteries. The materials provide flexibility and conformability to move through tortuous anatomy and track to distal lesions. The balloon is inflated to its predetermined diameter and exerts a constant, gentle outward force to establish patency.

3 Warnings

- Inspect the device prior to the procedure to verify the device is intact and functional. Do not use the device if the outer or inner packaging is damaged or opened.
- When the catheter is in the body, it should be manipulated under sufficient high quality fluoroscopy. Prior to withdrawing the balloon catheter from the lesion, the balloon must be completely deflated under vacuum. If resistance is met during manipulation, determine the cause of the resistance before proceeding.
- Do not use air or any gaseous medium to inflate the balloon. Use only the recommended inflation medium.
- Do not expose the device to organic solvents, for example, alcohol.
- Do not exceed the Rated Burst Pressure (RBP). The RBP is based on the results of in vitro testing. Use of a pressure-monitoring device is recommended to prevent over-pressurization.
- To reduce the risk of vessel damage, the inflated diameter of the balloon should approximate the diameter of the vessel just proximal and distal to the stenosis.
- Use the catheter prior to the Use-by date specified on the package.
- Do not use with Lipiodol™ or Ethiodol™ contrast media (or other such contrast media which incorporate the components of these agents).

4 Precautions

- This device was designed for single use only. Do not reuse, reprocess, or resterilize this device. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device or create a risk of contamination, which could result in patient injury, illness, or death.
- Appropriate drug therapy (anticoagulant, vasodilator, etc.) should be administered to the patient according to standard protocols for PTA before insertion of the dilatation catheter.
- Precautions to prevent or reduce clotting should be taken when any catheter is used. Flush or rinse all products entering the vascular system with sterile isotonic saline or a similar solution via the guidewire access port prior to use. Consider the use of systemic heparinization.

- Exercise care during handling in order to avoid possible damage to the catheter. Do not use a catheter that has been damaged.
- To minimize the possible introduction of air into the system, it is imperative that before proceeding careful attention is paid to the maintenance of tight catheter connections and through aspiration and flushing of the system.
- Amphirion Deep PTA catheter should be used with caution for procedures involving calcified lesions due to the abrasive nature of these lesions.
- Never advance the PTA catheter without the guidewire extending from the tip.
- Never attempt to move the guidewire when the balloon is inflated.
- Do not advance the catheter against significant resistance. The cause of resistance should be determined via fluoroscopy and remedial action taken.
- In order to activate the LFC hydrophilic coating, it is recommended to wet the Amphirion Deep PTA catheter with a sterile saline solution immediately prior to its insertion in the body.

5 Adverse effects

Complications associated with the use of Amphirion Deep PTA catheter are similar to the ones associated with standard PTA procedures.

Possible complications may include, but are not limited to the following:

- Death
- Hemorrhage or hematoma
- Restenosis of the dilated vessel
- Vessel dissection, perforation, rupture, spasm or injury
- Drug reactions, allergic reaction to contrast medium
- Hypotension
- Infection
- Arteriovenous fistula
- Pain and tenderness
- Arrhythmias
- Pseudo-aneurysm
- Thromboembolic episodes
- Thrombosis and/or embolization
- Sepsis/infection
- Endocarditis

If a serious incident related to the device occurs, immediately report the incident to Medtronic and the applicable competent authority or regulatory body.

6 Instructions for use

Prior to angioplasty, carefully examine all equipment to be used during the procedure, including the dilatation catheter, to verify proper function. Verify that the catheter and sterile packaging have not been damaged in shipment and that the catheter size is suitable for the specific procedure for which it is intended.

1. Inflation device preparation

- Prepare the inflation device according to the manufacturer's instructions.
- Prior to use, carefully examine the unit to verify that the catheter and sterile packaging have not been damaged during shipment. After moistening the balloon's protective sheath with a physiological saline solution, gently remove the protective sheath, taking care not to damage the balloon in the process. In case of resistance, twist the protective sheath with one hand while holding the shaft steady with the other.

2. Dilatation catheter selection

The nominal balloon size must be chosen equal to or smaller than the inner diameter of the artery proximal and distal to the lesion. If the stenosis cannot be crossed with the desired dilatation catheter, use a smaller diameter catheter to pre-dilate the lesion to facilitate the passage of a more appropriate-sized dilatation catheter.

3. Dilatation catheter preparation

- The catheter is packaged in a protective hoop; carefully remove the catheter from the package.
- Remove the stylet and the protective sheath from the balloon.
- The balloon catheter in deflated position contains tiny air bubbles that should be purged prior to inserting the balloon catheter. To do this, connect a three way stopcock to the inflation port fitting on the dilatation catheter. Flush through

the stopcock. Connect a Luer-lock syringe, partly filled with sterile normal saline and contrast medium, to the stopcock. Never use air or any gaseous medium to inflate the balloon. Orient the dilatation catheter with the distal tip and the balloon pointing in a downward vertical position. Pull back the plunger and aspirate for 15 seconds until the air is completely evacuated. Release the plunger. Disconnect the syringe and evacuate the collected air. Reconnect the syringe and repeat this operation a couple of times until the balloon is completely free of air bubbles.

- Flush the wire lumen properly through the connector luer hub.

4. Inflation device connection to catheter

- To remove any air lodged in the distal Luer fitting of the inflation device, purge approximately 1 mL (cm³) of contrast medium.
- With the stopcock in the closed position, disconnect the syringe used in preparation applying a slight positive pressure. A meniscus of contrast medium will appear in the balloon port when the syringe is removed. Verify that a meniscus of contrast medium is evident in both the dilatation catheter balloon port (hub) and the inflation device connection. Securely couple the inflation device to the balloon port of the balloon dilatation catheter.

5. Use of Amphirion Deep PTA catheter

- Insert a guidewire through the hemostatic valve following the manufacturer's instructions or standard practice. Advance the guidewire carefully into the guiding catheter/ introducer sheath. When complete, withdraw the guidewire/ introducer, if used.
- Attach a torque device to the wire, if desired. Under fluoroscopy, advance the guidewire to the desired vessel, then across the stenosis.
- In order to activate the coating, wet the balloon catheter with a sterile saline immediately prior to its insertion in the introducer sheath and thoroughly aspirate and flush the guiding catheter in preparation (when guiding catheter is used) for introduction of the dilatation catheter. Load the distal tip of the dilatation catheter onto the guidewire.
Note: To avoid kinking, advance the dilatation catheter slowly, in small increments until the proximal end of the guidewire emerges from the catheter.
- Advance the catheter through the hemostatic valve slowly, while the balloon is fully deflated. It should be observed that the hemostatic valve is only closed as much as to prevent blood return yet permitting easy movements of the dilatation catheter. If resistance is encountered, do not advance the catheter through the adapter.
- Under fluoroscopy, use the balloon radiopaque markers to position the balloon within the lesion to be dilated and inflate the balloon to the appropriate pressure (refer to balloon compliance table). Maintain negative pressure on the balloon between inflations.
- Completely deflate the balloon catheter. Withdraw the deflated dilatation catheter and guidewire from the guiding catheter / introducer sheath, through the hemostatic valve. Tighten the knurled knob on the hemostatic valve.
- If necessary the balloon catheter can be exchanged over the guidewire which remains in the vessel, for different balloon types or sizes.
- After use, this product may be a potential biohazard. Handle and dispose of all such devices in accordance with accepted medical practice and applicable local, regional, and national laws and regulations.

7 How supplied

The Amphirion Deep PTA catheter is supplied sterile and intended for single use only. It is sterilized by ethylene oxide gas and it will remain sterile as long as the packaging remains unopened and undamaged. Use product prior to Use-by date.

Caution: Do not use if the inner package is open or damaged.

8 Storage

Store the device in a dry area away from sunlight.

Do not store catheters where they are directly exposed to organic solvents, ionizing radiation, or ultraviolet light. Rotate inventory so that the device is used before the Use-by date on the package label.

9 Device disposal

Dispose of the device in accordance with applicable laws, regulations, and hospital procedures, including those regarding biohazards, microbial hazards, and infectious substances.

10 Disclaimer of warranty

Note: This disclaimer of warranty does not apply in Australia.

The warnings contained in the product labeling provide more detailed information and are considered an integral part of this disclaimer of warranty. Although the product has been manufactured under carefully controlled conditions, Medtronic has no control over the conditions under which this product is used. Medtronic, therefore, disclaims all warranties, both express and implied, with respect to the product, including, but not limited to, any implied warranty of merchantability or fitness for a particular purpose. Medtronic shall not be liable to any person or entity for any medical expenses or any direct, incidental, or consequential damages caused by any use, defect, failure, or malfunction of the product, whether a claim for such damages is based upon warranty, contract, tort, or otherwise. No person has any authority to bind Medtronic to any representation or warranty with respect to the product.

The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this disclaimer of warranty is held to be illegal, unenforceable, or in conflict with applicable law by a court of competent jurisdiction, the validity of the remaining portions of this disclaimer of warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this disclaimer of warranty did not contain the particular part or term held to be invalid.

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Artery Lesions (AcoArt II-BTK). *J Endovasc Ther.* 2021;28(2): 215-221.

- ^c Liistro F, Porto I, Angioli P, et al. Drug-eluting balloon in peripheral intervention for below the knee angioplasty evaluation (DEBATE-BTK): a randomized trial in diabetic patients with critical limb ischemia. *Circulation.* 2013;128(6):615-621.
- ^d Wang J, Zhu YQ, Zhao JG, et al. Infrapopliteal angioplasty with a long over-the-wire (OTW) balloon in the treatment of severe limb ischemia in diabetic patients: a retrospective study. *Acta Radiol.* 2009;50(4): 360-367.
- ^e Deloose K, Bosiers M, Peeters P. One year outcome after primary stenting of infrapopliteal lesions with the CHROMIS DEEP stent in the management of critical limb ischaemia. *EuroIntervention.* 2009;5(3): 318-324.
- ^f Troisi N, Ercolini L, Chisci E, et al. Use of Tapered Balloons to Recanalize Occluded Below-the-Knee Arteries in Diabetic Patients with Critical Limb Ischemia. *Ann Vasc Surg.* 2016;31:105-110.
- ^g Wei LM, Zhu YQ, Zhao JG, Wang J, Lu HT, Zhang PL. Retrograde transplatar arch angioplasty of below-the-knee arterial occlusions: outcomes compared to antegrade recanalization. *Acad Radiol.* 2014;21(11):1475-1482.

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^b Jia X, Zhuang B, Wang F et al. Drug-Coated Balloon Angioplasty Compared With Uncoated Balloons in the Treatment of Infrapopliteal Artery Lesions (AcoArt II-BTK). *J Endovasc Ther*. 2021;28(2):215-221.

^c Liistro F, Porto I, Angioli P et al. Drug-eluting balloon in peripheral intervention for below the knee angioplasty evaluation (DEBATE-BTK): a randomized trial in diabetic patients with critical limb ischemia. *Circulation*. 2013;128(6):615-621.

^d Wang J, Zhu YQ, Zhao JG et al. Infrapopliteal angioplasty with a long over-the-wire (OTW) balloon in the treatment of severe limb ischemia in diabetic patients: a retrospective study. *Acta Radiol*. 2009;50(4):360-367.

^e Deloose K, Bosiers M, Peeters P. One year outcome after primary stenting of infrapopliteal lesions with the CHROMIS DEEP stent in the management of critical limb ischaemia. *EuroIntervention*. 2009;5(3):318-324.

^f Troisi N, Ercolini L, Chisci E et al. Use of Tapered Balloons to Recanalize Occluded Below-the-Knee Arteries in Diabetic Patients with Critical Limb Ischemia. *Ann Vasc Surg*. 2016;31:105-110.

^g Wei LM, Zhu YQ, Zhao JG, Wang J, Lu HT, Zhang PL. Retrograde transplatar arch angioplasty of below-the-knee arterial occlusions: outcomes compared to antegrade recanalization. *Acad Radiol*. 2014;21(11):1475-1482.

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^a Gandini R, Volpi T, Pampana E, Uccioli L, Versaci F, Simonetti G. Applicability and clinical results of percutaneous transluminal angioplasty with a novel, long, conically shaped balloon dedicated for below-the knee interventions. *J Cardiovasc Surg (Torino)*. 2009; 50(3): 365-371.

^b Jia X, Zhuang B, Wang F, et al. Drug-Coated Balloon Angioplasty Compared With Uncoated Balloons in the Treatment of Infrapopliteal Artery Lesions (AcoArt II-BTK). *J Endovasc Ther*. 2021; 28(2): 215-221.

^c Liistro F, Porto I, Angioli P, et al. Drug-eluting balloon in peripheral intervention for below the knee angioplasty evaluation (DEBATE-BTK): a randomized trial in diabetic patients with critical limb ischemia. *Circulation*. 2013; 128(6): 615-621.

^d Wang J, Zhu YQ, Zhao JG, et al. Infrapopliteal angioplasty with a long over-the-wire (OTW) balloon in the treatment of severe limb ischemia in diabetic patients: a retrospective study. *Acta Radiol*. 2009; 50(4): 360-367.

^e Deloose K, Bosiers M, Peeters P. One year outcome after primary stenting of infrapopliteal lesions with the CHROMIS DEEP stent in the management of critical limb ischaemia. *EuroIntervention*. 2009; 5(3): 318-324.

^f Troisi N, Ercolini L, Chisci E, et al. Use of Tapered Balloons to Recanalize Occluded Below-the-Knee Arteries in Diabetic Patients with Critical Limb Ischemia. *Ann Vasc Surg*. 2016; 31: 105-110.

^g Wei LM, Zhu YQ, Zhao JG, Wang J, Lu HT, Zhang PL. Retrograde transplanter arch angioplasty of below-the-knee arterial occlusions: outcomes compared to antegrade recanalization. *Acad Radiol*. 2014; 21(11): 1475-1482.

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^a Gandini, R.; Volpi, T.; Pampana, E.; Uccioli, L.; Versaci, F.; Simonetti, G. Applicability and clinical results of percutaneous transluminal angioplasty with a novel, long, conically shaped balloon dedicated for below-the-knee interventions. *J Cardiovasc Surg (Torino)*. 2009; 50 (3): 365–371.

^b Jia, X.; Zhuang, B.; Wang, F. et al. Drug-Coated Balloon Angioplasty Compared With Uncoated Balloons in the Treatment of Infrapopliteal Artery Lesions (AcoArt II-BTK). *J Endovasc Ther*. 2021; 28 (2): 215–221.

^c Liistro, F.; Porto, I.; Angioli, P. et al. Drug-eluting balloon in peripheral intervention for below the knee angioplasty evaluation (DEBATE-BTK): a randomized trial in diabetic patients with critical limb ischemia. *Circulation*. 2013; 128 (6): 615–621.

^d Wang, J.; Zhu, Y. Q.; Zhao, J. G. et al. Infrapopliteal angioplasty with a long over-the-wire (OTW) balloon in the treatment of severe limb ischemia in diabetic patients: a retrospective study. *Acta Radiol*. 2009; 50 (4): 360–367.

^e Deloose, K.; Bosiers, M.; Peeters, P. One year outcome after primary stenting of infrapopliteal lesions with the CHROMIS DEEP stent in the management of critical limb ischaemia. *EuroIntervention*. 2009; 5 (3): 318–324.

^f Troisi, N.; Ercolini, L.; Chisci, E. et al. Use of Tapered Balloons to Recanalize Occluded Below-the-Knee Arteries in Diabetic Patients with Critical Limb Ischemia. *Ann Vasc Surg*. 2016; 31: 105–110.

^g Wei, L. M.; Zhu, Y. Q.; Zhao, J. G.; Wang, J.; Lu, H. T.; Zhang, P. L. Retrograde transplatar arch angioplasty of below-the-knee arterial occlusions: outcomes compared to anterograde recanalization. *Acad Radiol*. 2014; 21 (11): 1475–1482.

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^a Gandini R, Volpi T, Pampana E, Uccioli L, Versaci F, Simonetti G. Applicability and clinical results of percutaneous transluminal angioplasty with a novel, long, conically shaped balloon dedicated for below-the-knee interventions. *J Cardiovasc Surg (Torino)*. 2009;50(3):365-371.

^b Jia X, Zhuang B, Wang F, et al. Drug-Coated Balloon Angioplasty Compared With Uncoated Balloons in the Treatment of Infrapopliteal Artery Lesions (AcoArt II-BTK). *J Endovasc Ther*. 2021;28(2):215-221.

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^d Wang J, Zhu YQ, Zhao JG, et al. Infrapopliteal angioplasty with a long over-the-wire (OTW) balloon in the treatment of severe limb ischemia in diabetic patients: a retrospective study. *Acta Radiol*. 2009;50(4):360-367.

^e Deloose K, Bosiers M, Peeters P. One year outcome after primary stenting of infrapopliteal lesions with the CHROMIS DEEP stent in the management of critical limb ischaemia. *EuroIntervention*. 2009;5(3):318-324.

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^g Wei LM, Zhu YQ, Zhao JG, Wang J, Lu HT, Zhang PL. Retrograde transplanter arch angioplasty of below-the-knee arterial occlusions: outcomes compared to antegrade recanalization. *Acad Radiol*. 2014;21(11):1475-1482.q

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^a Gandini R, Volpi T, Pampana E, Uccioli L, Versaci F, Simonetti G. Applicability and clinical results of percutaneous transluminal angioplasty with a novel, long, conically shaped balloon dedicated for

below-the-knee interventions. *J Cardiovasc Surg (Torino)*. 2009;50(3): 365-371.

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^g Wei LM, Zhu YQ, Zhao JG, Wang J, Lu HT, Zhang PL. Retrograde transplanter arch angioplasty of below-the-knee arterial occlusions: outcomes compared to anterograde recanalization. *Acad Radiol*. 2014;21(11):1475-1482.

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- ^a Gandini R, Volpi T, Pampana E, Uccioli L, Versaci F, Simonetti G. Applicability and clinical results of percutaneous transluminal angioplasty with a novel, long, conically shaped balloon dedicated for below the knee interventions. *J Cardiovasc Surg (Torino)*. 2009;50(3): 365-371.
- ^b Jia X, Zhuang B, Wang F, et al. Drug-Coated Balloon Angioplasty Compared With Uncoated Balloons in the Treatment of Infrapopliteal Artery Lesions (AcoArt II-BTK). *J Endovasc Ther*. 2021;28(2): 215-221.
- ^c Liistro F, Porto I, Angioli P, et al. Drug-eluting balloon in peripheral intervention for below the knee angioplasty evaluation (DEBATE-BTK): a randomized trial in diabetic patients with critical limb ischemia. *Circulation*. 2013;128(6):615-621.
- ^d Wang J, Zhu YQ, Zhao JG, et al. Infrapopliteal angioplasty with a long over-the-wire (OTW) balloon in the treatment of severe limb ischemia in diabetic patients: a retrospective study. *Acta Radiol*. 2009;50(4): 360-367.
- ^e Delooste K, Bosiers M, Peeters P. One year outcome after primary stenting of infrapopliteal lesions with the CHROMIS DEEP stent in the management of critical limb ischaemia. *EuroIntervention*. 2009;5(3): 318-324.
- ^f Troisi N, Ercolini L, Chisci E, et al. Use of Tapered Balloons to Recanalize Occluded Below-the-Knee Arteries in Diabetic Patients with Critical Limb Ischemia. *Ann Vasc Surg*. 2016;31:105-110.
- ^g Wei LM, Zhu YQ, Zhao JG, Wang J, Lu HT, Zhang PL. Retrograde transplatar arch angioplasty of below-the-knee arterial occlusions: outcomes compared to antegrade recanalization. *Acad Radiol*. 2014;21(11):1475-1482.

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- ^aGandini R, Volpi T, Pampana E, Uccioli L, Versaci F, Simonetti G. Applicability and clinical results of percutaneous transluminal angioplasty with a novel, long, conically shaped balloon dedicated for below-the knee interventions. J Cardiovasc Surg (Torino). 2009;50(3):365-371.
- ^bJia X, Zhuang B, Wang F, et al. Drug-Coated Balloon Angioplasty Compared With Uncoated Balloons in the Treatment of Infrapopliteal Artery Lesions (AcoArt II-BTK). J Endovasc Ther. 2021;28(2):215-221.
- ^cLiistro F, Porto I, Angioli P, et al. Drug-eluting balloon in peripheral intervention for below the knee angioplasty evaluation (DEBATE-BTK): a randomized trial in diabetic patients with critical limb ischemia. Circulation. 2013;128(6):615-621.
- ^dWang J, Zhu YQ, Zhao JG, et al. Infrapopliteal angioplasty with a long over-the-wire (OTW) balloon in the treatment of severe limb ischemia in diabetic patients: a retrospective study. Acta Radiol. 2009;50(4):360-367.
- ^eDeloose K, Bosiers M, Peeters P. One year outcome after primary stenting of infrapopliteal lesions with the CHROMIS DEEP stent in the management of critical limb ischaemia. EuroIntervention. 2009;5(3):318-324.
- ^fTroisi N, Ercolini L, Chisci E, et al. Use of Tapered Balloons to Recanalize Occluded Below-the-Knee Arteries in Diabetic Patients with Critical Limb Ischemia. Ann Vasc Surg. 2016;31:105-110.
- ^gWei LM, Zhu YQ, Zhao JG, Wang J, Lu HT, Zhang PL. Retrograde transplanter arch angioplasty of below-the-knee arterial occlusions: outcomes compared to antegrade recanalization. Acad Radiol. 2014;21(11):1475-1482.

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^a Gandini R, Volpi T, Pampana E, Uccioli L, Versaci F, Simonetti G. Applicability and clinical results of percutaneous transluminal angioplasty with a novel, long, conically shaped balloon dedicated for below-the knee interventions. *J Cardiovasc Surg (Torino)*. 2009;50(3): 365-371.

^b Jia X, Zhuang B, Wang F, et al. Drug-Coated Balloon Angioplasty Compared With Uncoated Balloons in the Treatment of Infrapopliteal Artery Lesions (AcoArt II-BTK). *J Endovasc Ther*. 2021;28(2): 215-221.

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^g Wei LM, Zhu YQ, Zhao JG, Wang J, Lu HT, Zhang PL. Retrograde transplanter arch angioplasty of below-the-knee arterial occlusions: outcomes compared to antegrade recanalization. *Acad Radiol*. 2014;21(11):1475-1482.

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- ^a Gandini R, Volpi T, Pampana E, Uccioli L, Versaci F, Simonetti G. Applicability and clinical results of percutaneous transluminal angioplasty with a novel, long, conically shaped balloon dedicated for below-the-knee interventions. *J Cardiovasc Surg (Torino)*. 2009;50(3):365-371.
- ^b Jia X, Zhuang B, Wang F, et al. Drug-Coated Balloon Angioplasty Compared With Uncoated Balloons in the Treatment of Infrapopliteal Artery Lesions (AcoArt II-BTK). *J Endovasc Ther*. 2021;28(2):215-221.
- ^c Liistro F, Porto I, Angioli P, et al. Drug-eluting balloon in peripheral intervention for below the knee angioplasty evaluation (DEBATE-BTK): a randomized trial in diabetic patients with critical limb ischemia. *Circulation*. 2013;128(6):615-621.
- ^d Wang J, Zhu YQ, Zhao JG, et al. Infrapopliteal angioplasty with a long over-the-wire (OTW) balloon in the treatment of severe limb ischemia in diabetic patients: a retrospective study. *Acta Radiol*. 2009;50(4):360-367.
- ^e Deloose K, Bosiers M, Peeters P. One year outcome after primary stenting of infrapopliteal lesions with the CHROMIS DEEP stent in the management of critical limb ischaemia. *EuroIntervention*. 2009;5(3):318-324.
- ^f Troisi N, Ercolini L, Chisci E, et al. Use of Tapered Balloons to Recanalize Occluded Below-the-Knee Arteries in Diabetic Patients with Critical Limb Ischemia. *Ann Vasc Surg*. 2016;31:105-110.
- ^g Wei LM, Zhu YQ, Zhao JG, Wang J, Lu HT, Zhang PL. Retrograde transplatar arch angioplasty of below-the-knee arterial occlusions: outcomes compared to antegrade recanalization. *Acad Radiol*. 2014;21(11):1475-1482.

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^b Jia X, Zhuang B, Wang F a kol. Drug-Coated Balloon Angioplasty Compared With Uncoated Balloons in the Treatment of Infrapopliteal Artery Lesions (AcoArt II-BTK). *J Endovasc Ther*. 2021;28(2): 215-221.

^c Liistro F, Porto I, Angioli P a kol. Drug-eluting balloon in peripheral intervention for below the knee angioplasty evaluation (DEBATE-BTK): a randomized trial in diabetic patients with critical limb ischemia. *Circulation*. 2013;128(6):615-621.

^d Wang J, Zhu YQ, Zhao JG a kol. Infrapopliteal angioplasty with a long over-the-wire (OTW) balloon in the treatment of severe limb ischemia in diabetic patients: a retrospective study. *Acta Radiol*. 2009;50(4): 360-367.

^e Deloose K, Bosiers M, Peeters P. One year outcome after primary stenting of infrapopliteal lesions with the CHROMIS DEEP stent in the management of critical limb ischaemia. *EuroIntervention*. 2009;5(3): 318-324.

^f Troisi N, Ercolini L, Chisci E a kol. Use of Tapered Balloons to Recanalize Occluded Below-the-Knee Arteries in Diabetic Patients with Critical Limb Ischemia. *Ann Vasc Surg*. 2016;31:105-110.

^g Wei LM, Zhu YQ, Zhao JG, Wang J, Lu HT, Zhang PL. Retrograde transplatar arch angioplasty of below-the-knee arterial occlusions: outcomes compared to antegrade recanalization. *Acad Radiol*. 2014;21(11):1475-1482.

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- ^a Gandini R, Volpi T, Pampana E, Uccioli L, Versaci F, Simonetti G. Applicability and clinical results of percutaneous transluminal angioplasty with a novel, long, conically shaped balloon dedicated for below-the-knee interventions. *J Cardiovasc Surg (Torino)*. 2009;50(3):365-371.
- ^b Jia X, Zhuang B, Wang F và cộng sự. Drug-Coated Balloon Angioplasty Compared With Uncoated Balloons in the Treatment of Infrapopliteal Artery Lesions (AcoArt II-BTK). *J Endovasc Ther*. 2021;28(2):215-221.
- ^c Liistro F, Porto I, Angioli P và cộng sự. Drug-eluting balloon in peripheral intervention for below the knee angioplasty evaluation (DEBATE-BTK): a randomized trial in diabetic patients with critical limb ischemia. *Circulation*. 2013;128(6):615-621.
- ^d Wang J, Zhu YQ, Zhao JG và cộng sự. Infrapopliteal angioplasty with a long over-the-wire (OTW) balloon in the treatment of severe limb ischemia in diabetic patients: a retrospective study. *Acta Radiol*. 2009;50(4):360-367.
- ^e Deloosse K, Bosiers M, Peeters P. One year outcome after primary stenting of infrapopliteal lesions with the CHROMIS DEEP stent in the management of critical limb ischaemia. *EuroIntervention*. 2009;5(3):318-324.
- ^f Troisi N, Ercolini L, Chisci E và cộng sự. Use of Tapered Balloons to Recanalize Occluded Below-the-Knee Arteries in Diabetic Patients with Critical Limb Ischemia. *Ann Vasc Surg*. 2016;31:105-110.
- ^g Wei LM, Zhu YQ, Zhao JG, Wang J, Lu HT, Zhang PL. Retrograde transplatar arch angioplasty of below-the-knee arterial occlusions: outcomes compared to antegrade recanalization. *Acad Radiol*. 2014;21(11):1475-1482.

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© 2021 Medtronic
Medtronic, Medtronic with rising man logo, and Medtronic logo are trademarks of Medtronic. Third-party trademarks ("TM") belong to their respective owners. The following list includes trademarks or registered trademarks of a Medtronic entity in the United States and/or in other countries.
Medtronic, Inc.
Medtronic B.V.