Medtronic

Symplicity G3™

Renal Denervation RF Generator, RDNG3A

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Contents

1 Introduction	6
1.1 Overview	6
1.2 Contraindications	6
1.3 Warnings and precautions	6
2 Product overview	7
2.1 Product description	7
2.2 System components	8
2.3 Compatible components	8
2.4 System notes	8
3 Installation, setup, and operation	9
3.1 Symplicity G3 generator installation	9
3.2 Symplicity G3 generator setup	9
3.3 Attaching the dispersive electrode	9
3.4 Treatment procedure	10
4 Accessories	10
4.1 Cart (optional)	10
4.2 Remote control	11
5 Product details	12
5.1 Front panel	12
5.2 Touchscreen interface	12
5.3 Side panel connections	25
5.4 Back panel	27
5.5 Detection of previously used catheter	27
6 System software	28
7 Audio and display information	28
8 Check status conditions	29
9 Fault conditions	31
10 Troubleshooting	32
11 Cleaning	33
11.1 Symplicity G3 renal denervation RF generator	33
11.2 Remote control	33
11.3 Cart	33
12 Maintenance	34
12.1 Recommended annual electrical safety inspection: IEC62353	34
12.2 End of life disposition	34
13 Product specifications	34
14 Technical information	35
14.1 Power output diagrams	35
14.2 Electromagnetic compatibility (EMC)	35
14.3 Electromagnetic emissions per IEC 60601-1-2	36
14.4 Electromagnetic immunity per IEC 60601-1-2	36
15 Security	37
15.1 Data security	37
15.2 Cybersecurity events	37
16 Limited warranty	37

1 Introduction

1.1 Overview

The Symplicity G3 renal denervation RF generator is designed to deliver low-level radiofrequency (RF) energy, via the Symplicity Spyral multi-electrode renal denervation catheter, through the wall of the renal artery to achieve renal denervation. Components of the system include the disposable multichannel Symplicity Spyral catheter, the Symplicity G3 generator, and remote control. The catheter is used to deliver RF energy to a target site. The generator provides an automatically controlled treatment of RF energy delivered by the catheter to the target site. The remote control provides a method for channel selection and RF activation from the sterile field.

The single-use, Symplicity Spyral catheter is 6 Fr guide catheter compatible. Using standard interventional techniques, the Symplicity Spyral catheter is placed in the renal artery, so RF energy can be delivered through the catheter's electrodes to the arterial wall and surrounding tissue where sympathetic nerves reside. Refer to the Symplicity Spyral catheter's Instructions for Use for further information about using the catheter.

The unit is specifically configured for use with dispersive electrodes (also known as a neutral electrodes, return electrodes, or grounding pads) conforming to the IEC 60601-2-2 standard (see *Section 2.3* for compatibility information). The Symplicity G3 generator has been tested and found to comply with the electromagnetic emissions limits for medical devices given in IEC 60601-1-2. These limits are designed to provide reasonable protection against harmful interference in a typical medical installation.

The date of manufacture of the unit can be found by checking the unit's back panel.

1.2 Contraindications

There are no known absolute contraindications to the use of radiofrequency surgery. Please refer to the warnings and precautions section (*Section 1.3*) for more details.

1.3 Warnings and precautions

Note: Refer to the Symplicity Spyral catheter's Instructions for Use for additional warnings and precautions, including those related to patients.

Related to the use of radiofrequency in catheterization laboratories

- Carefully read all instructions before use.
- Radiofrequency surgery uses high-frequency output. Do not perform procedures if flammable or explosive media are present, such as flammable anesthetics or skin preparation agents.
- Interference produced by the operation of the Symplicity G3 generator may adversely influence the operation of other electronic equipment.
- The Symplicity G3 generator may produce hazardous electrical output and is designed to be used only by qualified medical personnel trained in the use of this equipment.
- Safe use of monopolar radiofrequency energy demands proper attachment of a dispersive electrode to the patient. Follow all of the manufacturer's directions for skin preparation and the placement of the dispersive electrode and for proper insulation between the patient and any metal surfaces.
- The dispersive electrode should be placed on the thigh or other nonbony area of the body and should be outside of the angiographic field of view. Shave the placement area for good contact between the dispersive electrode and the skin. Failure to achieve good skin contact by the entire adhesive surface of the dispersive electrode may result in a burn or high impedance measurements. Do not apply the dispersive electrode where fluid may pool.
- The entire area of the dispersive electrode should be reliably attached to a suitably prepared and appropriate area of the patient's body per the above instructions.
- Skin-to-skin contact (for example, between the arms and body of the patient) should be avoided. Use dry gauze to prevent skin-to-skin contact.
- The patient should not come into contact with metal parts that are grounded or have an appreciable capacitance to ground (for example, operating table supports). The use of antistatic sheeting is recommended.
- Monitoring electrodes of any physiological monitoring equipment used on the patient at the same time as RF energy delivery should be placed as far away as possible from the lower back region. Needle monitoring electrodes are not recommended. Physiological monitoring equipment incorporating high-frequency current-limiting devices is recommended. In all cases, monitoring systems incorporating high-frequency current-limiting devices are recommended.

Related to RF treatment

- Remove any guidewires that are not contained within the Symplicity Spyral catheter (such as a buddy wire) from the treatment site before activating the RF output.
- Implantable pacemakers (IPGs) and implantable cardioverter defibrillators (ICDs) or other active implants may be adversely affected by RF ablation. Refer to the implantable device's Instructions for Use for more information.

- Do not allow the Symplicity Spyral catheter electrode to come into contact with a metal instrument or surface during energy delivery. Superficial skin burns could occur.
- Loss of contact between the dispersive electrode and the patient will only result in an auditory alert when using the dispersive electrodes specified in *Section 2.3*.
- Do not touch the Symplicity Spyral catheter electrodes and the dispersive electrode at the same time during energy delivery. Superficial skin burns could occur.
- Initiate RF energy delivery only when the system is displaying the READY screen, as defined in *Chapter 5*.

Related to Symplicity G3 generator and ablation catheter

- The generator should be powered on and allowed to complete the system self-tests before introducing the catheter to the vasculature.
- The Symplicity G3 generator should be placed at least 30 cm (12 in) away from walls or corners.
- The Symplicity G3 generator should sit for a minimum of 10 min with power ON before initiating a treatment.
- Use the Symplicity G3 generator in a well-vented area. Do not drape the Symplicity G3 generator because proper ventilation is required throughout the procedure.
- Check all external surfaces of the Symplicity G3 generator, as well as components such as the power cord, foot switch, and remote for mechanical damage before each use. The equipment should not be used if any damage is found. If inspection reveals a defect, the Symplicity G3 generator should not be used until it has been properly repaired/serviced. The operator must immediately notify Medtronic of the defect.
- A failure of the Symplicity G3 generator could result in an unintended increase of output power.
- To avoid the risk of electric shock, this equipment must only be connected to AC supply mains with protective earth.
- If the Symplicity G3 generator malfunctions, stop using it and contact a Medtronic representative.
- The Symplicity G3 generator is intended to be used only with the Symplicity Spyral catheter.
- Use the Symplicity G3 generator only with the accessories in *Section 2.3*. Other active accessories or equipment may not have the appropriate ratings for the maximum output voltage of the Symplicity G3 generator.
- Do not connect any cables or equipment to the Symplicity G3 generator that are not specified by Medtronic. Doing so could potentially damage the system or cause injury.
- Do not remove the cover of the Symplicity G3 generator as electric shock may occur. Refer servicing to qualified personnel. There are no user-serviceable parts inside the Symplicity G3 generator.
- Do not disassemble the Symplicity G3 generator. No modification of this equipment is allowed. This will void the warranty and may cause serious injury to the operator or damage to the unit. If any problems are not resolved by the directions in *Chapter 10*, contact Medtronic for further instructions.
- The covered connectors are for Medtronic use only.
- When connecting to an auxiliary secondary monitor (DVI-D), the Symplicity G3 generator should only be connected to a monitor that is certified to IEC60950.

2 Product overview

2.1 Product description

The Symplicity G3 generator is designed to deliver precisely controlled RF energy treatments using the Symplicity Spyral catheter. The Symplicity G3 generator is an automatically-controlled electrosurgical generator that produces RF energy according to a programmed algorithm. During energy delivery, power is continuously monitored and controlled based on temperature and impedance measured at the electrode-tissue interface to ensure proper operation. The Symplicity G3 generator is compact for ease of use within the interventional suite. It is activated by using an optional foot switch, the hand-held remote control, or through the interface on the Symplicity G3 generator. The remote is nonsterile and needs to be draped if operated within the sterile field.

Warning: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 in) to any part of the Symplicity G3 generator, including cables specified by Medtronic. Otherwise, degradation of the performance of this equipment could result.

Figure 1. Symplicity G3 generator

2.2 System components

Components provided with the Symplicity G3 generator:

- Hospital-grade AC power cord
- Remote control
- DVI-D cable

Components available separately:

- Symplicity Spyril catheter
- Symplicity G3 generator cart
- Dispersive electrode (see *Section 2.3* for compatible dispersive electrodes)
- Foot switch (see *Section 2.3* for foot switch compatibility)

2.3 Compatible components

Catheter compatibility

The Symplicity G3 generator is designed to work specifically with the Symplicity Spyril catheter, and should only be used with Medtronic manufactured or specified compatible devices. Closely follow the Instructions for Use provided with the Symplicity Spyril catheter.

Dispersive electrode compatibility

The Symplicity G3 generator is designed to work specifically with the dispersive electrodes listed in *Table 1*.

Table 1. Compatible dispersive electrodes

Compatible dispersive electrodes	Covidien Valleylab REM Polyhesive Adult Patient Return Electrode, Model E7507 Covidien Valleylab REM Polyhesive Adult Patient Return Electrode, Model E7507DB
----------------------------------	--

Optional foot switch compatibility

The Symplicity G3 generator is designed to work specifically with the optional foot switch listed in *Table 2*.

Table 2. Compatible foot switch

Compatible foot switch	Herga Technology Foot Switch, Model 6210-0058
------------------------	---

2.4 System notes

Use only the Symplicity Spyril catheter with the Symplicity G3 generator.

The Symplicity Spyril catheter and dispersive electrodes are intended for single use only. Do not attempt to sterilize and reuse them. Once used, catheters and dispersive electrodes should be treated as biohazards and disposed of in accordance with hospital policies, local governing ordinances, and recycling plans.

3 Installation, setup, and operation

Note: Review all warnings and precautions before using the Symplicity G3 generator.

3.1 Symplicity G3 generator installation

Before initial use of the Symplicity G3 generator, complete the performance tests listed in this section. This ensures that the Symplicity G3 generator is operating properly when a new Symplicity G3 generator is initially put into service, or when the Symplicity G3 generator has been moved between facilities. There is no special equipment required for the performance tests outlined in this section.

Caution: Do not use the Symplicity G3 generator if it does not meet the inspection criteria detailed below. Contact your local Medtronic representative if the system does not behave as indicated for the following performance tests:

1. Visual inspection

Complete a visual inspection according to the following instructions:

- Verify that the labels on the Symplicity G3 generator are legible.
- Inspect the Symplicity G3 generator, connectors, AC power cord, remote control, and any additional components or accessories to be used with the generator for any evidence of damage.
- Inspect the catheter and catheter cable for breaks, cracks, nicks, and other damage before using the catheter. Do not use the catheter if any damage is found.

2. System self-test

Power on the Symplicity G3 generator and verify the following:

- Verify that the display illuminates with the appropriate start-up screen.
- Verify that all indicator lights temporarily illuminate.
- Verify that the Speaker Test audio alert sounds to indicate that the speaker is functioning.
- Verify that the system self-test was successful and that no faults occur. See *Chapter 9* for more information on faults.

3. Cooling fan operation

While the Symplicity G3 generator is powered on, verify the following:

- Verify that the fan at the back of the Symplicity G3 generator is functioning.
- Verify that the fan is free of dust or any foreign particles. If not, take a clean brush and dust the fan opening to provide air flow.

3.2 Symplicity G3 generator setup

Note: Perform the following steps before accessing the patient's vasculature or opening/connecting a catheter:

1. If applicable and desired, attach the Symplicity G3 generator to the cart by placing it into the latching mechanism on the top of the cart. Ensure the Symplicity G3 generator latches into place and is secure within the cart.
2. If applicable and desired, connect the remote control to the back of the Symplicity G3 generator. While in use, ensure the remote control is properly draped in order to be accessed within the sterile field.
3. If applicable and desired, insert the foot switch connector into the receptacle on the back panel. Keep the foot switch in an area where it is easy to recognize and operate without any confusion between the foot switch for the Symplicity G3 generator and any other foot pedal or foot switch present in the operating area.
4. Ensure the Symplicity G3 generator is at least 30 cm (12 in) away from any corner or wall so that it can properly ventilate during a procedure.
5. Power on the Symplicity G3 generator:
 - Verify the system self-test was successful (see the Symplicity G3 generator installation section).
 - If the Symplicity G3 generator does not complete the system self-test, do not use the Symplicity G3 generator (see *Chapter 10*). If the problem cannot be addressed, contact your Medtronic representative.
 - Following system self-test, the system is in STANDBY. No measurements or power controls are enabled.
6. Ensure the appropriate date and time settings are correct on the Symplicity G3 generator screen. If necessary, adjust the date and time in the **Settings** menu.

Note: The Symplicity G3 generator will not automatically update for local and annual time changes.

3.3 Attaching the dispersive electrode

Prepare the patient using standard techniques for electrosurgery. Ensure that the patient's entire body, including extremities, is insulated from contact with grounded metal parts. The Symplicity G3 generator is only compatible with the dispersive electrodes listed in *Section 2.3*. Closely follow the instructions provided by the manufacturer of the dispersive electrode. Failure to achieve good skin contact by the entire adhesive surface of the dispersive electrode may result in a burn or measured high impedance. The dispersive

electrode should be placed on the thigh or another nonbony area out of the angiogram field (if necessary, shave the placement area for good contact). Do not apply the dispersive electrode where fluid may pool.

3.4 Treatment procedure

Before activating the RF output, verify that the following steps are taken:

1. Connect the dispersive electrode to the appropriate receptacle on the Symplicity G3 generator.
2. Prepare the catheter per its Instructions for Use.
3. Connect the catheter to the Symplicity G3 generator.
4. Select **New Patient** or **Same Patient**, as appropriate, on the touchscreen interface on the Symplicity G3 generator.
5. Position the catheter to the targeted renal anatomy per the appropriate interventional procedure (see the catheter Instructions for Use).
6. During catheter positioning, instructions to guide the workflow of the procedure will appear on the screen of the Symplicity G3 generator.
7. When the electrodes of the catheter are not fully exposed to the vasculature (that is, if any electrodes are present within the guide catheter), the messaging screens will persist.
8. Once all catheter electrodes are positioned in the vasculature, distal to the guide catheter, screen messages will be replaced with impedance measurements. This is part of the software's READY state:
 - a. The impedance measurements should be used along with fluoroscopy to assess catheter positioning. The impedance plot for each channel provides the last 20 s of measurements for the respective channel.
 - b. Electrode stability is portrayed by the stability of the plot; electrode movement may be reflected by fluctuations of the plot.
Note: Electrode stability does not necessarily reflect adequate electrode contact against the arterial wall. Fluoroscopic imaging is required to ensure proper electrode contact.
 - c. If the screen continues to show messages and does not show impedance graphs and values after the electrodes have been positioned in the vasculature, then ensure that the catheter and the dispersive electrode are properly connected to the patient and generator. Proper and full contact of the entire adhesive surface of the dispersive electrode is required for operation of the generator and catheter; repositioning of the dispersive electrode may be required.
9. Electrodes may be deselected and selected using the respective channel numbers on the touch screen or remote control. Channel 1 reflects the most distal electrode of the catheter, while channel 4 reflects the most proximal electrode. When an electrode is deselected, RF energy is not applied to the deselected electrode.
10. Electrode contact to the renal artery needs to be established so that the therapy is delivered to the site of the electrodes. Refer to the catheter Instructions for Use for techniques to improve electrode contact with arterial wall.
11. Before initiation of RF energy, ensure that the RF channels and kidney demarcation are appropriately selected.
12. When it is appropriate to apply RF energy, press the foot switch, remote control **RF** button, or the **RF** button on the front of the Symplicity G3 generator to begin RF energy delivery.
13. During the procedure, the display shows temperature, percent impedance drop, the starting impedance, and the remaining treatment time. The RF generator delivers power using an automated algorithm. When the treatment is complete, power delivery automatically stops. At any point in the procedure, treatment on all active channels may be stopped by pressing the foot switch, remote control **RF** button, or the **RF** button on the front of the Symplicity G3 generator.
Note: If the impedance rises above the appropriate limits or if the temperature exceeds threshold levels, the Symplicity G3 generator stops delivery of RF energy on the affected channels (see *Chapter 8*).
14. From the Report screen, the treatment summary reports can be downloaded to a compatible flash drive via the USB port. A compatible flash drive is a USB 2.0 flash drive, FAT32 formatted, with at least 100 MB of free space. In addition, from the Report screen, treatment summary data can be viewed.
Note: Generator storage is limited to a maximum of 200 treatment summary reports. When storage limits are met, the oldest treatment summary reports will be deleted to create storage space for new treatment summary reports. Medtronic recommends downloading treatment summary reports after each patient or on a regular schedule.
15. See the Symplicity Spyril catheter IFU for specific instructions regarding catheter use.
16. At the end of the procedure, disconnect the catheter, power off the Symplicity G3 generator, disconnect the power cord from the power source, and fasten the power cord on the straps located on the back of the generator. Disconnect all remaining cables before storage. Store the Symplicity G3 generator in the appropriate storage location.

4 Accessories

4.1 Cart (optional)

An optional wheeled cart is shown in *Figure 2*. If available, the cart can be used for keeping the Symplicity G3 generator stable during therapy delivery and also to move the Symplicity G3 generator from one location to another. When kept on the cart, the Symplicity G3

generator automatically latches to the cart using its own weight. If the Symplicity G3 generator has to be detached from the cart, a latch lever (located on the back of the cart) releases the Symplicity G3 generator.

Figure 2. Symplicity G3 renal denervation RF generator on a wheeled cart

4.2 Remote control

As shown in *Figure 3*, a remote control is available and can be used by connecting it on the back panel port provided for the remote control. Refer to the Back panel section (*Section 5.4*) for details. The remote control can be stored in a detachable holster, which can be connected to the side panel of Symplicity G3 generator as shown in *Figure 4*.

Figure 3. Remote control

Figure 4. Remote control stored in holster on the side panel of the Symplicity G3 renal denervation RF generator

The remote control duplicates some of the control functionalities of the touch screen as mentioned.


1. The **RF** button on the remote control provides the same functionality as the **RF** button on the front panel of the Symplicity G3 generator.
2. The buttons numbered 1, 2, 3, and 4 provide the same functionality as electrode selection buttons on the touch screen.
3. The **L** and **R** buttons (with kidney shadows) allow the user to toggle between the **L** and **R** buttons on the touch screen. This functionality is available only during the READY state, and provides an annotation in the ablation tab header indicating whether

the ablation was performed in a renal artery leading to the left or right kidney. Labeling tabs is an optional step, and can be completed only before performing an ablation (it is not possible to label a tab once an ablation has been completed).

4. The arrow keys duplicate the function of the left and right tab buttons on the touch screen, which can be used to move between different treatment screens. Pressing the left arrow displays the previous treatment and pressing the right arrow displays the latest treatment.

Note: The remote control is not a sterile device. If the remote control is being used in the sterile area, then draping is required, as shown in *Figure 5*.

Figure 5. Remote control in draping for use in the sterile field


5 Product details

5.1 Front panel

Figure 6 shows the front panel display, **RF** button, and indicator lights.

Figure 6.

- 1 Display: The display is the primary means of providing information to the user. During start-up and use, the status of the unit and energy delivery is displayed on these screens. The display has a touchscreen function that enables the user to interact with the system.
- 2 **RF** button: Can be used to activate RF energy output while the system is in READY or deactivate RF energy output when the system is in RF ON (see *Chapter 6*). The **RF** button functions the same as the remote control or foot switch (see *Section 5.4* for information about the foot switch connection on the back panel).
- 3 Indicator Lights: blue=RF ON; red=FAULT. See *Chapter 7* for a functional description of indicator lights.


5.2 Touchscreen interface

The display of the Symplicity G3 generator contains a touchscreen interface that enables the user to interact with the system throughout the various states of the software, with the exception of the RF ON state. The touchscreen can be used to set the date and time, select a new patient or the same patient when a catheter is connected, deselect and select channels for RF activation, scroll

through the ablation history plots of a given patient, select which kidney is being treated, access and download treatment reports and summaries, and other activities (see *Chapter 6* for a description of the referenced states).

5.2.1 Standby State - Connect Catheter screen

Figure 7. Connect Catheter screen


The **Connect Catheter** screen is displayed in the STANDBY state of the system when no catheter is connected to the system.






Icon	Description
	The Connect Catheter label points to the side panel where the Symplicity Spyrax catheter is attached.
	The Medtronic logo is displayed on the Connect Catheter screen.
	From the Connect Catheter screen, selecting the Settings icon will enable the user to adjust the date, time, and language. (English is the only language available.)
	From the Connect Catheter screen, selecting the Reports icon will enable the user to select, view, and download treatment summaries.

5.2.2 Standby State - Settings screen

Figure 8. Settings screen


The Settings screen can be accessed from the **Connect Catheter** screen when no catheter is connected. This screen is accessible in the STANDBY state.

Icon	Description
	In the Settings screen, making selections in this field will modify the default language. Note: In order to change the setting, the Save selection must be made.
	In the Settings screen, making selections in these fields will modify the date. Note: In order to change the setting, the Save selection must be made.
	In the Settings screen, making selections in these fields will modify the time. Note: In order to change the setting, the Save selection must be made.
	In the Settings screen, these icons are used to scroll up or down through the values of date, time, and language selections. (English is the only language available).
	In the Settings screen, pressing the Save icon will save the currently selected date, time, and language. (English is the only language available).
	In the Settings screen, pressing the Exit icon exits out of the Settings screen.

5.2.3 Standby State - Patient Selection screen

Figure 9. Patient Selection screen


The Patient Selection screen is displayed once a catheter is connected in the STANDBY state.

Icon	Description
	In the Patient Selection screen, selecting the New Patient icon associates the connected catheter with a new patient which resets the ablation counter.
	In the Patient Selection screen, selecting the Same Patient icon associates the connected catheter with the previously connected catheter's patient history and maintains the ablation counter from the last connected catheter.

5.2.4 Ready State - Message screen

Figure 10. Message screen


In the READY state, the Message screen is used to display troubleshooting messages.

Icon	Description
	In the Message screen, this symbol is followed by troubleshooting and guidance messages.
	In the READY state, the left and right arrows can be used to move between different treatment screens.
 OR 	In the READY state, the tab is used to view previously completed ablations. The indicator L or R within the tab indicates that Left or Right selection was made. The number represents the ablation count. If the tab is marked by an orange box, it indicates that 1 or more channels ended in a treatment out of range or check status condition.
	From the READY state, selecting the Reports icon will enable the user to select, view, and download treatment summaries.

5.2.5 Ready State - Ready screen

Figure 11. Ready screen


In the READY state, the Ready screen is used to monitor the impedance of all active channels.

Icon	Description
 OR 	In the READY state, the icon indicates the most distal channel on the four-channel catheter. The icon is shaded green when selected or shaded gray when deselected.
 OR 	In the READY state, the icon indicates the second most distal channel on the four-channel catheter. The icon is shaded green when selected or shaded gray when deselected.
 OR 	In the READY state, the icon indicates the second most proximal channel on the four-channel catheter. The icon is shaded green when selected or shaded gray when deselected.

Icon	Description
 <p>OR</p>	<p>In the READY state, the icon indicates the most proximal channel on the four-channel catheter. The icon is shaded green when selected or shaded gray when deselected.</p>
 <p>OR</p>	<p>In the Ready screen, selecting the L or R icons associates the current ablation counter with the renal artery leading to the left or right kidney.</p>
	<p>In the READY state, the left and right arrows appearing on the screen can be used to move between different treatment screens. The left arrow is pressed to display the previous treatment and right arrow is pressed to display the latest treatment.</p>
	<p>In the READY state, the text displays the current status of the RF delivery.</p>
	<p>This is the symbol of ohm associated with the impedance value displayed.</p>
 <p>OR</p>	<p>In the READY state, the tab is used to view previously completed ablations. The indicator L or R within the tab indicates that Left or Right selection was made. The number represents the ablation count. If the tab is marked by an orange box, it represents that one or more channels ended in a treatment out of range or check status condition.</p>
	<p>From the READY state, selecting the Reports icon will enable the user to select, view, and download treatment summaries.</p>

5.2.6 RF ON state - RF ON screen

Figure 12. RF ON screen


In the RF ON state, the **RF ON** screen provides feedback on the progress of the ablation.

Icon	Description
	In the RF ON screen, this icon indicates the most distal channel on the four-channel catheter.
	In the RF ON screen, this icon indicates the second most distal channel on the four-channel catheter.
	In the RF ON screen, this icon indicates the second most proximal channel on the four-channel catheter.
	In the RF ON screen, this icon indicates the most proximal channel on the four-channel catheter.
	In the RF ON state, the text displays the current status of the RF delivery (for example, RF ON).
	This icon appears next to the temperature value displayed in degrees Celsius.
	This is the percent icon associated with impedance drop from the baseline impedance value in RF ON condition.

Icon	Description
 OR  OR 	In the RF ON screen, these icons show if the impedance has decreased, stayed the same or increased from the baseline impedance.


5.2.7 RF OFF state - Check Status screen

Figure 13. Check Status screen


The Check Status screen is displayed when the ablation has been shut down due to a check status condition.

Icon	Description
	This icon is displayed to show that RF delivery has been shut off due to a check status condition.
	The channel that had the check status condition is displayed in orange color.

Icon	Description
 OR  OR 	These icons, which show if the impedance decreased, stayed the same, or increased from the baseline impedance, are grayed out to illustrate that the impedance will not change as RF delivery has stopped.
	The stop icon is displayed on channels that are shut off due to a check status condition on a different channel.
	From the Check Status screen, selecting the Reports icon will enable the user to select, view, and download treatment summaries.

5.2.8 Check Dispersive Electrode screen

Figure 14. Check Dispersive Electrode screen


The **Check Dispersive Electrode** screen is displayed when the ablation has been shut down due to a check status condition related to the dispersive electrode.

Icon	Description
	This icon is displayed to show that RF delivery has been shut off due to a check status condition related to the dispersive electrode. Refer to <i>Chapter 10</i> for troubleshooting instructions.

5.2.9 RF OFF state - User-Ended screen

Figure 15. User-Ended screen


The User-Ended screen is displayed when a treatment is stopped by the user.

Icon	Description
	The stop icon is displayed on channels due to a user-ended treatment.
	From the User-Ended screen, selecting the Reports icon will enable the user to select, view, and download treatment summaries.



5.2.10 Report State - Report screen

Figure 16. Report screen


The Report screen is used to display, select, and download treatment summaries.

Icon	Description
	In the Report screen, the left and right arrows can be used to move between different pages with sets of up to 5 reports per page.
	In the Report screen, the box icon indicates if the report associated with that data set is selected for download. When selected, the box icon will be filled with a check symbol. Data column: Selecting box icons in the Data column is inadvisable for non-Medtronic personnel. The data is encrypted and only usable by Medtronic representatives for research purposes. The device will lock until encrypted data is exported. Summary column: Selecting box icons in the Summary column flags individual treatment summary reports for export.
	In the Report screen, selecting the Select All Summary icon will select all available treatment summary reports for download. It will not include data from the Data column.
	In the Report screen, selecting the Deselect All icon will deselect all available treatment summary reports and detailed reports.
	In the Report screen, selecting the Export icon will download all selected treatment summaries to the USB memory stick connected to the USB-A port.

Icon	Description
	From the Report screen, selecting the Exit icon will exit out of the Report state.
	In the Report screen, selecting the Summary icon will bring up the Treatment Summary screen.

5.2.11 Report State - Treatment Summary screen

Figure 17. Treatment Summary screen


The **Treatment Summary** screen is used to view treatment summaries of each ablation of the selected patient.

Icon	Description
	In the Treatment Summary screen, the left and right arrows can be used to move between different treatment summary pages within the same selected patient.
Electrode #	In the Treatment Summary screen, the label Electrode # refers to the electrodes in the four-channel catheter with 1 being the most distal and 4 being the most proximal.
Baseline Imp	In the Treatment Summary screen, the label Baseline Imp refers to the impedance value of the channel when the treatment was initiated.
Max Imp Drop	In the Treatment Summary screen, the label Max Imp Drop refers to the maximum percent reduction in impedance seen during the treatment.
Max Temperature	In the Treatment Summary screen, the label Max Temperature refers to the maximum temperature reached during the treatment.
Max Power	In the Treatment Summary screen, the label Max Power refers to the maximum power delivered during the treatment.
End Power	In the Treatment Summary screen, the label End Power refers to the power level at the end of the treatment.

Icon	Description
RF Time	In the Treatment Summary screen, the label RF Time refers to the time (indicated in seconds) that the therapy was delivered.
Kidney selection	In the Treatment Summary screen, the label Kidney selection displays R or L if the kidney selection button was pressed in the Ready screen for the associated treatment.
Code	In the Treatment Summary screen, the label Code will display Normal for a normal treatment or an associated Check Status Code per <i>Chapter 8</i> .

5.2.12 Fault State - FAULT screen

Figure 18. FAULT screen

The **FAULT** screen is displayed when a Fault condition develops.

Icon	Description
	This icon indicates that a fault has taken place (see <i>Chapter 9</i> for more information).

5.3 Side panel connections

The side panel connections are shown in *Figure 19*. The Symplicity Spyral catheter and dispersive electrode connect to the Symplicity G3 generator on the side panel.


Figure 19. Side panel layout

Side panel graphics		
1		USB-A: This port is for downloading treatment summary data to a memory stick.
2		Four-channel catheter connector: The four-channel catheter is connected on this port.
3		Dispersive electrode (floating): The dispersive electrode connects to the Symplicity G3 generator via this port. The dispersive electrode provides a path for the electrical current through the patient and back to the Symplicity G3 generator. It is important to properly attach the dispersive electrode to the patient (see Attaching the dispersive electrode in <i>Chapter 3</i>).

5.4 Back panel

Figure 20. Back panel layout

The back panel is represented in *Figure 20*.

Back panel graphics		
1		Allows for projecting the display of the Symplicity G3 generator onto an external monitor using a DVI-D cable.
2		Remote control connection: Allows for the connection of the remote control.
3		Volume adjust: Allows for adjustment of the output level of the volume-adjustable audio alerts.
4		Potential equalization terminal: Provides common electrical ground for other electronic devices used during the procedure.
5		Foot switch connection: Allows for connection of an optional foot switch. The pneumatic foot switch is connected via the connector adjacent to this symbol. Depress the foot switch to begin delivering RF energy. To stop the energy delivery during operation, depress the foot switch again.
6		AC power connector: Connector for the AC line power cable.
7		Power switch: Turns the Symplicity G3 generator on and off. Position the Symplicity G3 generator so that the power switch is easy to operate to disconnect power from the Symplicity G3 generator.
8		Fuse door: For use by trained Medtronic service personnel only.
9		Straps: Secures any cord slack, or is used for storage of the cord on the back of the generator when the generator is not in use.

5.5 Detection of previously used catheter

The system will recognize a previously used catheter and display a message to the user to replace the catheter with a new catheter. Refer to *Chapter 7*.

6 System software

The Symplicity G3 generator has the following system states that are controlled by software and user inputs. See *Chapter 7* for audio alerts and information provided to the user in each state.

State	Description
STANDBY	System state after the Symplicity G3 generator is turned on and has completed a successful system self-test. Entered only after self-test. Medtronic logo and software version are displayed.
READY	System state following STANDBY. Short pulses of nonharmful low-level RF power are delivered to allow for impedance measurements, which may aid in positioning the electrode. This state is accompanied by impedance plots showing on the display when the Symplicity G3 generator is ready to initiate power delivery. An informational message may be displayed if certain conditions are detected. See <i>Chapter 7</i> for more information.
RF ON	System state during RF power delivery that executes the programmed control algorithm. RF ON will be displayed along with the RF progress bar. If the programmed algorithm detects a measured treatment condition, for example, temperature or impedance out of its acceptable range, RF output power is automatically discontinued to the respective electrodes within this state, or the system transitions to the RF OFF state.
RF OFF	System state when RF power delivery is automatically discontinued after treatment time, or delivery is manually stopped by the user (through use of the RF button, optional foot switch, or remote control), or stopped by the control algorithm. This system state is also reached when a combination of low temperature rise and impedance decrease is detected at the end of the treatment period (see <i>Chapter 8</i>).
REPORT	System state when treatment summary data is being accessed and viewed on the display or downloaded onto an external storage device.
FAULT	System state when the Symplicity G3 generator detects an incorrect connection or setting, a self-test failure, or a failure in internal circuitry (see <i>Chapter 9</i>).

7 Audio and display information

The Symplicity G3 generator uses audio alerts, indicator lights, and the display screen to communicate information to the user. The information conveyed to the user during a procedure by these methods is summarized below.

Audio alert and display information			
Display message	Indicator light	Audio alert	Event
N/A [READY state]	RF ON pulsed	None	System in READY. Low-level RF output pulses for impedance measurements.
RF ON	RF ON lit (steady)	A repeating tone, followed by louder, repeating tones	System in RF ON and RF output power is being delivered.
RF ON (treatment nearing completion)	RF ON lit (steady)	A repeating tone faster than the standard RF ON tone to indicate that the treatment is about to complete	System in RF ON and RF output power is being delivered.
RF OFF (standard)	RF ON pulsed	A single, short tone	System transition to RF OFF after treatment has completed.
RF ON (with check status code during treatment)	RF ON lit (steady)	Three tones in short succession	System in RF ON and RF output power is no longer being delivered on the channel or channels with the Check Status condition.
RF OFF with check status code	RF ON pulsed	Three tones in short succession	System transition to RF OFF due to Check Status condition. See <i>Chapter 8</i> .
FAULT ##	FAULT lit	One continuous beep (lasting ~5 s.)	Fault condition exists. See <i>Chapter 9</i> . Turn off power and contact Medtronic.
Check catheter position	RF ON pulsed	None	Normally indicates that the catheter is partially in the guide catheter.
Insert catheter	RF ON pulsed	None	Indicates that the catheter has been detected but has not been inserted, or has been partially inserted, into the patient.

Audio alert and display information			
Display message	Indicator light	Audio alert	Event
High impedance. Check catheter position and/or connections.	RF ON pulsed	None	Indicates that the catheter is inserted into the patient but could still be in the guide catheter.
Replace catheter	RF ON pulsed	None	Indicates a potentially faulty connection in either the catheter or the cable. The catheter should be replaced.
Previously used catheter connected. Please insert a new catheter	RF ON pulsed	None	If enabled in software, this message indicates that a catheter used in a previous procedure has been attached to the Symplicity G3 generator.
Medtronic catheter not recognized. Please unplug and fully re-insert	RF ON pulsed	None	Indicates that the catheter has not been fully inserted into the receptacle. Remove and reinsert the catheter into the Symplicity G3 generator.
Report download in progress	N/A	None	Indicates that the selected reports are being downloaded to an external storage device.
Selected reports have been successfully exported	N/A	One tone	Indicates that the selected reports have finished downloading onto the external storage device.
Memory stick not usable. Please insert a different memory stick	N/A	None	Indicates that the selected reports cannot be downloaded onto the external storage device. Consider using a different external storage device.
Check dispersive electrode	RF ON not lit	Alternating tone (last ~5 s)	Indicates that a REM Alert condition has occurred and RF delivery to all channels has stopped.
Clock battery is dead. Time will not update when system is powered off. Contact Medtronic to have the battery replaced.	N/A	None	Indicates that the clock battery has died. Device is still capable of administering therapy. Press OK to continue and set the current date and time before administering treatment. Note: Date and time will not update while system is powered off until battery is replaced and will need to be manually updated before completing subsequent treatments. Contact your Medtronic representative to have the battery replaced.

8 Check status conditions

A Check Status condition occurs when the Symplicity G3 generator detects a treatment condition outside of its acceptable range. When a Check Status condition is detected, the Symplicity G3 generator will go into the RF OFF state or remain in the RF ON state, depending on the Check Status condition, and:

- Prevent or stop RF energy delivery on the respective electrodes
- Sound the Check Status audio alert
- Display a Check Status code and message

The following table provides the list of Check Status codes and conditions that could occur during a treatment. Suggested actions for addressing the condition are provided. The Check Status or user-ended treatment message can be cleared by pressing the > key on the display or remote control, or by pressing the foot switch 2 times within 3 s. If the condition persists after performing the suggested action, contact your Medtronic representative.

Symplicity G3 generator Check Status conditions			
Check Status code	Description	Action	Corresponding condition on the other channels
20	Treatment impedance out of range: High	If Check Status code 20 occurs within 15 s of starting treatment: Ensure no excessive tissue/catheter electrode contact, that is, vessel distension and that catheter is not positioned in guide, a branch, or small artery. Reposition the catheter electrode if necessary and attempt to restart treatment. If Check Status code 20 occurs after 15 s of starting treatment: Check the catheter/dispersive electrode connections and/or replace	The other channels stop power delivery, and system transitions to the check status state.

Symplicity G3 generator Check Status conditions

Check Status code	Description	Action	Corresponding condition on the other channels
		catheter/dispersive electrode. Proceed as appropriate.	
21	Treatment impedance out of range: Low	Remove or disconnect the catheter. Turn unit off, then on again. Ensure that other devices are not interfering with the Symplicity G3 generator. Proceed as appropriate.	The other channels stop power delivery, and system transitions to the check status state.
22a	Temperature too high	Image treatment site; wait approximately 1 min and image prior treatment site as increased vessel reactivity, such as spasm, may be encountered. Remove and inspect the catheter electrode. Proceed to the next treatment site as appropriate.	The other channels stop power delivery, and system transitions to the check status state.
22b	Temperature too high with indications of electrode movement	Image treatment site; wait approximately 1 min and image prior treatment site as increased vessel reactivity, such as spasm, may be encountered. Electrodes may have moved. Remove and inspect the catheter electrode. Proceed to the next treatment site as appropriate.	The other channels stop power delivery, and system transitions to the check status state.
23	Insufficient temperature rise	Check the catheter electrode position and reposition as necessary to improve contact with vessel wall.	The other channels continue power delivery. After RF activation is completed, check status message for the appropriate channel is shown.
24a	Unexpected impedance rise with indications of electrode movement	Image treatment site; wait approximately 1 min and image prior treatment site as increased vessel reactivity, such as spasm, may be encountered. Electrodes may have moved. Remove and inspect the catheter electrode. Proceed to the next treatment site as appropriate.	The other channels stop power delivery, and system transitions to the check status state.
24b	Insufficient impedance decrease	Electrodes may have moved. If not, reposition as necessary to improve contact with vessel wall. Proceed to the next treatment site as appropriate.	The other channels continue power delivery. After RF activation is completed, check status message for the appropriate channel is shown.
24c	Insufficient impedance decrease	Check the catheter electrode position and reposition as necessary to improve contact with vessel wall. Proceed to the next treatment site as appropriate.	N/A
26	Initial impedance out of range: High	Check that the catheter electrode is not within the guide catheter and at the desired treatment location. Ensure no excessive tissue/catheter electrode contact, that is, vessel distension and catheter is not positioned in a small branch. Reposition the catheter electrode to achieve lower starting impedance. Check catheter/dispersive electrode connections. Remove and inspect catheter electrode. Replace catheter/dispersive electrode if necessary. Proceed as appropriate.	The other channels stop power delivery, and system transitions to the check status state.
27	Initial impedance out of range: Low	Remove or disconnect catheter. Turn unit off, then on again. Ensure that other devices are not interfering with the Symplicity G3 generator.	The other channels stop power delivery, and system transitions to the check status state.

Simplicity G3 generator Check Status conditions

Check Status code	Description	Action	Corresponding condition on the other channels
		Proceed as appropriate.	
28	Starting temperature too high	Ensure that other devices are not interfering with the Simplicity G3 generator. Replace catheter.	The other channels stop power delivery, and system transitions to the check status state.
29	Starting temperature too low	Ensure that the catheter electrodes are not within the guide catheter and are within the desired treatment location. Also ensure that other devices are not interfering with the Simplicity G3 generator.	The other channels stop power delivery, and system transitions to the check status state.
30	Power shutoff due to continual temperature rise	Image treatment site; wait approximately 1 min and image prior treatment site as increased vessel reactivity, such as spasm, may be encountered. Remove and inspect catheter electrode. Proceed to the next treatment site as appropriate.	The other channels continue power delivery. After RF activation is completed, check status message for the appropriate channel is shown.
32	Unexpected temperature decrease	Electrodes may have moved. Reposition as necessary.	The other channels continue power delivery. After RF activation is completed, check status message for the appropriate channel is shown.
33	Low temp during RF ON	Replace catheter.	The other channels stop power delivery, and system transitions to the check status state.
34	Unexpected temperature rise on an electrode which has been turned OFF (not delivering RF)	Image treatment site; wait approximately 1 min and image prior treatment site as increased vessel reactivity, such as spasm, may be encountered. Remove and inspect catheter electrode. Proceed to the next treatment site as appropriate.	The other channels stop power delivery, and system transitions to the check status state.
35	Catheter movement/torquing	Entire catheter may have moved. Image treatment site to determine current position of electrodes. Proceed to next treatment site.	The other channels stop power delivery.
36	Shorting of the electrodes (in the Simplicity Spyrax catheter) during RF ON	The catheter may have moved. Image treatment site; wait approximately 1 min and image prior treatment site as increased vessel reactivity, such as spasm, may be encountered. Remove and inspect catheter electrode. Proceed to the next treatment site as appropriate.	The other channels stop power delivery, and system transitions to the check status state.
90	Malfunctioning catheter	Replace the catheter.	All channels stop power delivery, and system transitions to the check status state.

9 Fault conditions

A Fault condition occurs when the Simplicity G3 generator detects an incorrect connection or setting, a self-test failure, or a failure in internal circuitry. When a Fault condition is detected, the Simplicity G3 generator will go into the FAULT state and one or more of the following Fault state events can occur:

- Stop RF energy delivery
- Turn on the Fault indicator light
- Attempt to display a Fault code
- Sound the Fault audio alert
- No longer allow user input

The following table provides Fault conditions and codes that could occur during a treatment. Suggested actions for addressing the condition are provided. Fault conditions can only be reset by turning off power. Before turning off power, record the Fault code number.

Simplicity G3 generator Fault conditions		
Fault code	Description	Action
Undefined	Internal software or hardware fault	Turn unit off and contact your Medtronic representative.
11	Button failure	A button on the front panel or the foot switch may have been pressed while the unit was powered on. Turn unit off, check buttons/foot switch, and then turn unit on. If condition persists, contact your Medtronic representative.
46	Simplicity G3 generator internal temperature fault (onboard sensor)	If the internal Simplicity G3 generator temperature is outside of the allowed operating temperature range, this fault may occur. The unit must be operated in an environment where the temperature is within the bounds noted in the Operating environmental conditions under <i>Chapter 13</i> . If the unit is transferred from an environment where the temperature is outside these bounds, it will take some time for the internal temperature to come to the temperature of its surroundings. If this fault occurs, let the unit sit for 5 to 10 min with the power ON (the hardware fault beep will stop after 10 s). Power the unit off and then back on. If the fault occurs again, let the unit sit for several more minutes and power the unit off and then back on again. If the fault occurs yet again, turn unit off and contact your Medtronic representative.
48	Simplicity G3 generator internal temperature fault (CJC sensor)	If the internal Simplicity G3 generator temperature is outside of the allowed operating temperature range, this fault may occur. Follow the same course of action given for fault code 46.
76	Catheter connected at start-up	If a catheter is connected when the Simplicity G3 generator is powered on, the unit should be powered off. The catheter needs to be disconnected before the unit is powered back on. If condition persists, contact your Medtronic representative.
All other codes	Internal software or hardware fault	Turn unit off and contact your Medtronic representative.

10 Troubleshooting

The following table provides a list of probable symptoms, which may occur during routine operation. If a problem is encountered that is not listed here or cannot be addressed with the suggested actions, contact your Medtronic representative.

Troubleshooting	
Symptom	Action
The display or indicator lights do not illuminate when the power switch is turned on	Ensure that the Simplicity G3 generator is plugged into a working electrical outlet. Check the power switch on the rear of the unit. If problem persists, contact your Medtronic representative
Fault indicator light comes on and a Fault code appears in message window	See <i>Chapter 9</i> .
Check Dispersive Electrode appears in a message window	If this message appears, take the following actions: <ul style="list-style-type: none"> Verify that the dispersive electrode cord is correctly connected to the generator. Inspect the dispersive electrode plug and cord. If there is evidence of cracks, breaks, or other visible damage, replace the dispersive electrode or the cord. Verify that the compatible dispersive electrode is securely attached to the patient. (See <i>Table 1</i> for dispersive electrode compatibility information.)
Simplicity G3 generator appears frozen and will not accept user input	In certain fault conditions the Simplicity G3 generator may not be able to turn the Fault indicator on, beep, or display a Fault code. With the Simplicity G3 generator in this state, RF output cannot be turned ON and user input is disabled. Disconnect the catheter from the Simplicity G3 generator (if connected), turn unit off, and then back on again. If problem persists, contact your Medtronic representative.
RF power does not turn on when foot switch is depressed	Check that the Simplicity G3 generator is displaying the READY screen, as defined in <i>Chapter 5</i> . Check the foot switch—unplug the foot switch from the Simplicity G3 generator, place your thumb over the foot switch hose and depress the foot switch. If the foot switch is functioning properly, you should feel air pressure. Alternatively, RF power can be activated from the front panel of the generator.
RF interferes with ultrasound and other equipment	If the Simplicity Spyrax catheter along with the Simplicity G3 generator, causes interference with other equipment, which can be isolated by turning the Simplicity G3 generator off and on, try to correct the interference by one or more of the following measures:

Troubleshooting	
Symptom	Action
	<ul style="list-style-type: none"> • Reorient or relocate the receiving equipment. • Increase the separation between the equipment. • Connect the Symplicity G3 generator to an outlet on a circuit different from that to which the other equipment is connected. • Contact your Medtronic representative for assistance.
Impedance is less than 100 Ω	Remove catheter and inspect electrode and shaft; replace if necessary. Turn unit off, then on again. Ensure that other devices are not interfering with the Symplicity G3 generator. If the problem persists, contact your Medtronic representative.
Unable to clear informational message, Insert catheter	Confirm that the catheter is inserted into the guide catheter and into the patient's vasculature. Replace the catheter if message persists. If the problem persists, contact your Medtronic representative.
Unable to clear informational message, Check catheter position or dispersive electrode	Indicates that the catheter is inserted into the patient but could still be in the guide catheter, that the dispersive electrode is not properly connected, or that the system has encountered a high impedance condition on one or more of the channels. Check that all catheter electrodes are within the patient's vasculature, outside the guide catheter. Check that the dispersive electrode is properly connected to the Symplicity G3 generator and properly coupled to the patient. Consider changing the location of the dispersive electrode pad on the patient. If the problem persists, contact your Medtronic representative.
Replace catheter	Indicates a potentially faulty connection in the catheter. Before replacing the catheter, consider checking the dispersive electrode connection and changing its location on the patient. Replace the catheter if message persists. If the problem persists, contact your Medtronic representative.
Check catheter position	Indicates that the catheter is inserted into the patient but some electrodes may still be in the guide catheter. Verify that all catheter electrodes are within the patient's vasculature, outside the guide catheter. Check that the dispersive electrode is properly connected to the Symplicity G3 generator and properly coupled to the patient. Consider changing the location of the dispersive electrode pad on the patient. If the problem persists, contact your Medtronic representative.
Touchscreen icon press is missed	Press the desired icon longer or slower. Press touchscreen icons with a slow, methodical cadence to ensure optimal results.

11 Cleaning

11.1 Symplicity G3 renal denervation RF generator

- Do not use steam or heat sterilization. Do not soak in any disinfectants or fluids. Do not allow liquid to enter into any of the electrical connections or the interior of the unit.
- Disconnect the unit from either the AC or power inlet before cleaning.
- To clean, wipe the exterior surface of the unit with 70 to 90% isopropyl alcohol only.
- Allow all surfaces and connections to dry before reconnecting the system.

11.2 Remote control

- Do not use steam or heat sterilization. Do not soak in any disinfectants or fluids. Do not allow liquid to enter into any of the electrical connections or the interior of the remote control.
- Disconnect the remote control from the Symplicity G3 generator before cleaning.
- To clean, wipe the exterior surface of the remote control and cable with 70 to 90% isopropyl alcohol only.
- Allow all surfaces and connections to dry before reconnecting the remote control to the generator.

11.3 Cart

- Do not use steam or heat sterilization. Do not soak in any disinfectants or fluids. Do not allow liquid to enter the interior of the cart.
- Disconnect the Symplicity G3 generator from the cart before cleaning.
- To clean, wipe the exterior surface of the cart with 70 to 90% isopropyl alcohol only.
- Allow all surfaces and connections to dry before reattaching the Symplicity G3 generator.

12 Maintenance

12.1 Recommended annual electrical safety inspection: IEC62353

Follow test instructions in the appropriate area of the standard. The acceptable limits (normal condition) for the Symplicity G3 generator are as follows:

- Protective earth resistance (resistance between Mains Plug protective earth connector and protectively earthed accessible conductive parts): 0.3 Ω
- Earth (or Equipment) leakage current (all other AC operating voltages): 0.5 mA (500 μ A)
- Enclosure (or Touch) leakage current: 0.1 mA (100 μ A)
- Patient (or Applied Part) leakage current: 0.01 mA (10 μ A)

12.2 End of life disposition

Do not dispose of the Symplicity G3 generator or its accessories (cables, remote control, and cart) in the unsorted municipal waste stream. Follow local regulations for proper disposal.

Please refer to Medtronic website <http://recycling.Medtronic.com> for proper disposal guidance.

13 Product specifications

RF output power	6.5 W maximum per electrode $\pm 20\%$ accuracy with impedance of 175 Ω to 200 Ω $\pm 10\%$ accuracy with impedance of 201 Ω to 1200 Ω
Impedance measurement	175 Ω to 250 Ω , $\pm 20\%$ 251 Ω to 700 Ω , $\pm 10\%$ 701 Ω to 1200 Ω , $\pm 15\%$
Fuse	2 Slow-Blow (T3.15A L, 250 V), 20 x 5 mm
Temperature range	Measurement 37°C to 65°C, $\pm 3^\circ\text{C}$
Input power	100 to 240 V ~, 50/60 Hz universal power supply, 200 VA input power rating Disconnecting the AC power cord removes AC mains power from the unit.
Dimensions	29.9 cm x 36.7 cm x 18.9 cm (11.8 in x 14.5 in x 7.45 in) (height x width x depth)
Weight	9.07 kg (20 lb)
Front controls	Touchscreen buttons, RF button
Rear controls	Power On/Off, Volume adjust
Display	1 touchscreen display on front panel, 1024x768 XGA monitor displaying graphics, messages, impedance, temperature, and time
Connections	AC power, isolated USB, foot switch, four-channel catheter connector, remote control connector, DVI-D connector, and dispersive electrode connector
Protection	Class I, Defibrillator Proof—Type CF, intermittent operation with a duty cycle of 60 s on, 10 s off. The Symplicity Spyral catheter and Symplicity G3 generator are not suitable for use in the presence of a flammable anesthetic mixture with air, oxygen, or nitrous oxide. IEC 60601-1 Class I, IEC 60601-2-2 Safety of high-frequency surgical equipment; IEC 60601-1-2 Electromagnetic compatibility (EMC).
Rated accessory voltage	150 V _p
Maximum output voltage	150 V _p
Maximum current rating	700 mA averaged over 60 s
Maximum permitted length of accessory cords	DVI-D cable: 3 m (118 in) Power cord: 3 m (118 in) Remote control cord: Not applicable. The Symplicity G3 generator is only compatible with the provided remote control. The nominal remote control cord length is 2.8 m (109 in). Symplicity Spyral catheter cord: Not applicable. The Symplicity G3 generator is only compatible with the Symplicity Spyral catheter. The nominal catheter cord length is 3 m (118 in).

Transit conditions

Temperature	-35°C to +57°C (-31°F to +135°F)
Humidity	30% to 95% relative humidity, noncondensing
Pressure	595 hPa to 1060 hPa [-0.595 to 1.05 ATM]

Storage conditions

Temperature	15°C to 40°C (59°F to 104°F)
Humidity	10% to 90% relative humidity, noncondensing
Pressure	595 hPa to 1060 hPa [-0.595 to 1.05 ATM]

Operating environmental conditions


Temperature	15°C to 40°C (59°F to 104°F)
Humidity	30% to 70% relative humidity, noncondensing
Pressure	700 hPa to 1060 hPa [-0.7 to 1.05 ATM]

14 Technical information

14.1 Power output diagrams

Figure 21 shows the maximum power from the catheter as a sum of wattage from the active channels. Each individual electrode of the Symplicity Spyril catheter can deliver a maximum of 6.5 W.

Figure 21. Symplicity Spyril catheter total power delivered per active channel

- | | |
|-------------------|-----------------|
| 1 Impedance (Ω) | 4 3 channels on |
| 2 Total power (W) | 5 2 channels on |
| 3 4 channels on | 6 1 channel on |

Note: Total power delivered by the Symplicity Spyril catheter is presented as a function of the number of active electrodes.

14.2 Electromagnetic compatibility (EMC)

Interference produced by the operation of high-frequency surgical equipment, such as the Symplicity G3 generator, may adversely influence the operation of other electronic medical equipment such as monitors and imaging systems. Special precautions should be taken. The Symplicity Spyril catheter and Symplicity G3 generator should be installed and put into service according to the EMC information provided in this section.

Some degradation in RFID functionality may be experienced if the generator is in the proximity of other devices generating frequencies in the range of 12.5 MHz to 14.8 MHz. In such cases, the screen may display a message indicating that the catheter is not detected. To resolve this issue, the generator must be moved away from devices generating frequencies in this range or moved into another area.

Warning: The Symplicity Spyril catheter and Symplicity G3 generator are intended for use by healthcare professionals only. The Symplicity G3 generator may cause radio interference and disrupt the operation of nearby equipment. It may be necessary to take mitigation measures, such as reorienting or relocating the Symplicity G3 generator or shielding the location. Interference with

electronic medical equipment, such as monitors and imaging systems, is usually resolved or minimized by rearranging the cables such that the unit cables do not overlap with the cables from the monitoring equipment.

Warning: The Symplicity G3 generator should not be located immediately adjacent to or stacked with other equipment. If adjacent or stack use is necessary, the Symplicity G3 generator should be observed to verify normal operation in the configuration in which it will be used. Portable and mobile RF communications equipment can affect the Symplicity G3 generator.

Warning: The AC power cord set (AC power cable) used with the Symplicity G3 generator should be non-shielded, have a maximum length of less than 3 m, be at least 18 AWG (0.825 mm²) in USA/Canada, or 0.75 mm² in other countries and be approved for use in the country for which it is intended.

Warning: Use of accessories, transducers, and cables other than those specified or provided by Medtronic could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

Note: The Emissions characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 class A). If it is used in a residential environment (for which CISPR 11 class B is normally required) this equipment might not offer adequate protection to radio-frequency communication services. The user might need to take mitigation measures, such as relocating or re-orienting the equipment.

Note: The unit uses a 13.56 MHz RFID reader with ISO 14443/15693 protocols supporting OOK and ASK modulation. Effective Radiated Power (E.R.P.) = 0.0015 mW.

14.3 Electromagnetic emissions per IEC 60601-1-2

The Symplicity Spyril catheter and Symplicity G3 generator are suitable for use in the electromagnetic environment specified below. The user should ensure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment — guidance
RF emissions CISPR 11 (EN55011)	Group 1	When the Symplicity G3 generator is in the STANDBY state, RF energy is only used for its internal function. Therefore, its RF emissions are very low and are not likely to cause interference in nearby electronic equipment. During the RF ON state, the Symplicity G3 generator must emit electromagnetic energy to perform its intended function. Nearby electronic equipment may be affected.
RF emissions CISPR 11 (EN55011)	Class A	
Harmonic emissions IEC 61000-3-2	Class A	The Symplicity Spyril catheter and Symplicity G3 generator are suitable for use in all establishments other than domestic, that is, residential and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes, provided the warning following this table is heeded.
Voltage fluctuations/flicker emissions IEC 61000-3-3	Complies	

Warning: This equipment/system is intended for use by healthcare professionals only. This equipment/system may cause radio interference or may disrupt the operation of nearby equipment. It may be necessary to take mitigation measures, such as re-orienting or relocating the Symplicity G3 generator or shielding the location.

14.4 Electromagnetic immunity per IEC 60601-1-2

The Symplicity Spyril catheter and Symplicity G3 generator are suitable for use in the electromagnetic environment specified below. The user should ensure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment — guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±15 kV air	±8 kV contact ±15 kV air	Floors should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-5	±2 kV for power supply lines ±1 kV for input/output lines	±2 kV for power supply lines ±1 kV for input/output lines	AC mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±0.5 kV, ±1 kV line to line and ±0.5 kV, ±1 kV, 2 kV line to ground	±0.5 kV, ±1 kV line to line and ±0.5 kV, ±1 kV, 2 kV line to ground	AC mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions, and voltage variations on power supply input lines IEC 61000-4-11 U _T =230 Vac	0% U _T (100% dip in U _T) for 0.5 cycle 0% U _T (100% dip in U _T) for 1 cycle 70% U _T (30% dip in U _T) for 25/30 cycles	0% U _T (100% dip in U _T) for 0.5 cycle 0% U _T (100% dip in U _T) for 1 cycle 70% U _T (30% dip in U _T) for 25/30 cycles 0% U _T (100% dip in U _T) for 5 s	AC mains power quality should be that of a typical commercial or hospital environment. If the user of the Symplicity G3 generator requires continued operation during AC power mains interruptions, it is recommended that the

The Symplicity Spyral catheter and Symplicity G3 generator are suitable for use in the electromagnetic environment specified below. The user should ensure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment — guidance
	0% U_T (100% dip in U_T) for 5 s		Symplicity G3 generator be powered from an uninterruptible power supply or a battery.
Power frequency (50 Hz) magnetic field IEC 61000-4-8 U_T is the AC mains voltage before application of the test level.	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Conducted RF IEC 61000-4-6	3 V_{RMS} 150 kHz to 80 MHz 6 V_{RMS} ISM/Amateur Radio bands inside 150 kHz to 80 MHz	3 V_{RMS} 150 kHz to 80 MHz 6 V_{RMS} ISM/Amateur Radio bands inside 150 kHz to 80 MHz	Warning: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 in) to any part of the Symplicity G3 generator, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.7 GHz	3 V/m 80 MHz to 2.7 GHz	Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ^a should be less than the compliance level in each frequency range ^b .

Note: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflecting from structures, objects, and people.

Note: This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference in which case the user will be required to correct the interference at the user's own expense.

^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Symplicity G3 generator is used exceeds the applicable RF compliance level above, the Symplicity G3 generator should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the Symplicity G3 generator.

^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

15 Security

15.1 Data security

The Symplicity G3 generator uses and stores treatment data. The system does not protect exported data. Exported data should be handled in accordance with your facility's security policy for data handling and storage.

Medtronic recommends that you always export data to an encrypted mass storage device.

15.2 Cybersecurity events

If you suspect a cybersecurity event has occurred (such as strange or unexpected behavior, even if a fault or check status condition is not generated), stop using the generator (if possible). Contact your IT department or Medtronic support for information on how to confirm and respond to the suspected incident.

If you have further questions related to cybersecurity, contact your IT department or Medtronic support.

16 Limited warranty

Medtronic warrants that the Symplicity G3 RF generator, hereafter referred to as "generator," shall be free from material defects and defects in workmanship for a period of 12 months commencing on the date of delivery to the customer, provided that the customer provides Medtronic with written notice of the defect prior to the expiration of the warranty period and returns the generator to Medtronic no later than 10 days following the expiration of the warranty period.

Medtronic will, at its discretion, either:

- A. Without charge to the customer, repair the generator and bring it into full conformance with its specifications, including retest/reinspection by Medtronic; or,
- B. Without charge to the customer, replace the generator with a new device that is functionally comparable to the generator.

To qualify for this limited warranty, all of these conditions must be met:

- A. The generator must be used in accordance with the labeling, including but not limited to the instructions for use, recommendations, and manuals;
- B. The generator must not be subjected to unauthorized alteration, modification, or repair;
- C. The generator must not be subjected to misuse, abuse, accident, or improper handling;
- D. The generator must be properly maintained in accordance with the maintenance schedule and used solely in conjunction with products approved by Medtronic for use with the generator;
- E. The generator must be returned to Medtronic in accordance with the return guidelines stipulated in the user manual; the generator must not be moved outside the customer's site without prior notice and approval from Medtronic.

This limited warranty is limited to its express terms. In particular:

- A. Except as expressly provided by this limited warranty, Medtronic is not responsible for any direct, incidental, or consequential damages based on any defect, failure, or malfunction of the generator, whether the claim is based on warranty, contract, tort, or otherwise.
- B. This limited warranty is made only to the customer to whom the generator was originally sold. As to all others, Medtronic makes no warranty, express or implied, including but not limited to, any implied warranty of merchantability, or fitness for a particular purpose whether arising from statute, common law, custom, or otherwise. No such express or implied warranty to the customer shall extend beyond the 12 month period commencing on the date of delivery to the customer. This limited warranty shall be the exclusive remedy available to any person.
- C. The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene any mandatory provisions of applicable law. If any part or term of this limited warranty is held to be illegal, unenforceable, or in conflict with applicable law by a court of competent jurisdiction, the validity of the remaining portions of this limited warranty shall not be affected, and all rights and obligations shall be enforced as if this limited warranty did not contain the particular part or term held to be invalid.
- D. No person has any authority to bind Medtronic to any representation, condition, or warranty with respect to the generator that deviates in any respect from this limited warranty.
- E. This limited warranty is not applicable to catheters and accessories used with the generator.

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Medtronic, Medtronic with rising man logo, and Medtronic logo are trademarks of Medtronic. Third-party trademarks (“TM”) belong to their respective owners. The following list includes trademarks or registered trademarks of a Medtronic entity in the United States and/or in other countries.
Medtronic, Inc.