
VERTEBRAL BIOPSY NEEDLE

Intended Use: the device is intended to perform a vertebral biopsy during the kyphoplasty/vertebroplasty procedure in order to obtain a tissue sample to be analysed for diagnostic purposes.

CHARACTERISTICS AND PERFORMANCE:

The biopsy needle consists of an external cannula with a plastic handle, an internal stylet with a cutting wire at the distal end and a dedicated extractor to be used at the end of the biopsy procedure. The needle is introduced inside a working channel pre-inserted in the vertebral body during a kyphoplasty procedure. Turning the needle inside the working channel, as explained below, it's possible to extract a biopsy sample that remains inside the biopsy needle cannula. In order to collect the sample a dedicated extractor is used.

VERTEBRAL BIOPSY NEEDLE

This device has to be used only in combination with kyphoplasty medical devices "RENOVA SPINE" of Biopsybell.

INDICATIONS:

The device is part of the Kyphoplasty/vertebroplasty systems manufactured by Biopsybell used to treat pathological fractures of the vertebral body due to osteoporosis, cancer, or benign lesions. Cancer includes multiple myeloma and metastatic lesions, including those arising from breast or lung cancer, or lymphoma. Benign lesions include hemangioma and giant cell tumor.

CONTRA - INDICATIONS:

- Infections
- Patients with clotting disorders
- Patients with severe cardiac and / or pulmonary insufficiency
- Patients with known hypersensitivity or allergy to any of the components used in the procedure
- Vertebra plana or circumstances where safe percutaneous access to the vertebra cannot be guaranteed
- Unstable vertebral fractures due to posterior involvement
- Previous damage to the pedicle wall (transpedicular access)
- Lesions that feature narrowing of the spinal canal (more than 20%), including fracture or neoplasm, with or without myelopathy
- Retropulsing vertebral fragments with myelopathy
- A satisfactory response to conservative treatment
- Asymptomatic stable vertebral fractures

Target Patient:

The device is intended for the following patient population:

SEX: Male – Female

MATURATION OF THE SKELETAL SYSTEM: Skeletally mature patients

EXCLUSION: Newborns (<6 months) - Babies (>6 months to 2 years) - Skeletally immature patients

The device is used to treat patients with pathological fractures of the vertebral body due to osteoporosis, cancer, or benign lesions. Cancer includes multiple myeloma and metastatic lesions, including those arising from breast or lung cancer, or lymphoma. Benign lesions include hemangioma and giant cell tumor.

The patients to be excluded are the followings:

- Patients with infections
- Patients with clotting disorders
- Patients with severe cardiac and / or pulmonary insufficiency
- Patients with known hypersensitivity or allergy to any of the components used in the procedure
- Patients with vertebra plana or circumstances where safe percutaneous access to the vertebra cannot be guaranteed
- Patients with unstable vertebral fractures due to posterior involvement
- Patients with previous damage to the pedicle wall (transpedicular access)
- Patients with lesions that feature narrowing of the spinal canal (more than 20%), including fracture or neoplasm, with or without myelopathy
- Patients with retropulsing vertebral fragments with myelopathy
- Patients with a satisfactory response to conservative treatment
- Patients with asymptomatic stable vertebral fractures

INTENDED USERS AND ENVIRONMENT USE CONDITIONS:

The device is for professional use and, in particular, must be used by qualified and trained radiologists, neuro-radiologists, interventional radiologists, orthopaedic surgeons and neuro surgeons, pain therapists. The professional must have documented basic skills about biopsy and kyphoplasty. The device is to be used inside a sterile surgical field such as operating and angiographic rooms.

WARNINGS, PRECAUTIONS AND SIDE-EFFECTS:

- Do not use the instrument if the package is damaged.
- Do not use the instrument beyond the expiry date.
- The device is single-use and must be destroyed after use. Reuse is strictly prohibited. Reuse or resterilization may compromise the structural integrity of the device and/or cause damage to the device which may cause injury, disease or death of the patient. Reuse or resterilization may also create a risk of device contamination and/or cause infection of the patient or cross-infection, including but not limited to transmission of infectious diseases from one patient to the next. Device contamination may cause harm, disease or death of the patient.
- Do not reuse a component more than once.
- Only qualified and appropriately trained physicians should perform the vertebroplasty/kyphoplasty procedure.
- It is also recommended not to use the device if particularly risky manoeuvres are necessary where it will be subjected to high torsional or bending forces which might compromise the integrity and proper functioning of the system.
- The procedures must be performed under fluoroscopic guidance with appropriate X-ray equipment.
- It is strictly prohibited to modify, restructure, repair or resterilize the device for reuse.


INSTRUCTION FOR USE

Preliminary operations


- Before use check that the package is not damaged and that the set and its components are not damaged or kinked; if so, inform the manufacturer or your area representative. Take the devices out of the package operating in aseptic conditions.
- Sedate and position the patient. The entire procedure must be performed in real-time under combined CAT and fluoroscopic guidance. Digital angiography may also prove useful. Perform local or general anaesthesia.
- Perform the necessary skin disinfection procedures.
- Use in conjunction with a kyphoplasty instrument kit. Prepare the various instruments to be used for the operation. The inflation device, the balloon catheter and the entire kyphoplasty kit may only be used by medical staff.

Instructions

- 1) Create a working channel use the tools of the RENOVA SPINE kit of Biopsybell, following the specific instructions for use
- 2) Insert the needle (the complete system composed by the external cannula and the internal stylet with cutting wire) in the working channel
- 3) Then turn the complete needle VERTEBRAL BIOPSY NEEDLE of 360° for two complete rotations in such a way as to internally cut the fragment of tissue without damaging it and without the need to dislocate.


- 4) Finally, remove the complete system from the working channel, extract the internal **stylet with cutting wire** (A) and to deposit the sample taken onto the slide, insert the **finger-protection** device and the **Extractor** from the opposite side of the grip (B).

A

B

PACKAGING, STERILIZATION, SINGLE USE

The system is individually packed in rigid PVC blister and medical-grade paper. Do not use the device if the package is open or damaged or contaminated. It is sold in packs of 5 pieces. The needle is sterilized by Ethylene Oxide (EO) in compliance with ISO 11135 requirements. The device is not reusable. The device is single-use and must be destroyed after use.

DISPOSAL:

After use, the device may constitute a risk of contamination and/or infection and puncture or injuries caused by sharp parts of the device. Therefore, handle with care and dispose of the product in accordance with the regulations in force and medical practice.

NOTICE FOR USER AND/OR PATIENT: ANY SERIOUS INCIDENT THAT HAS OCCURRED IN RELATION TO THE DEVICE SHOULD BE REPORTED TO THE MANUFACTURER AND THE COMPETENT AUTHORITY OF THE MEMBER STATE IN WHICH THE USER AND/OR PATIENT IS ESTABLISHED.

Biopsybell declines all responsibility for any damage caused by improper use or different from what is indicated in this instruction leaflet.

Symbol	Description
	Manufacturer
	Catalogue number
	Lot number
	Quantity
	Date of manufacture
	Consult instructions for use
	Caution
	Use by date
	Sterilized using ethylene oxide
	Do not re-use
	Do not re-sterilize
	Compliant to MDD 93/42/EEC
	Caution: Federal law restricts this device to sale by or on the order of a physician.
	Do not use if package is damaged
	Keep dry

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