Thromboplastin L

REF 5265HL

REF 5265L

REF 5267L

Helena Biosciences Europe, Queensway South, Team Valley Trading Estate, Gateshead, Tyne and Wear, NE11 0SD, United Kingdom
Tel: +44 (0)191 482 8440
Fax: +44 (0)191 482 8442
Email: info@helena-biosciences.com
Web: www.helena-biosciences.com

HL-2-3035P 2015/10 (1)

Thromboplastin L

Instructions for use

en

INTENDED PURPOSE

The Thromboplastin L kit is intended for carrying out clot based haemostasis assays.

The first standardised one-stage prothrombin time test was developed by Dr. Armand Quick in 1935. It has now become the basic coagulation screening test for the diagnosis of congenital and acquired deficiencies of clotting factors from the extrinsic pathway (factors II, V, VII and X)^{1,2}. It is also used for the induction and monitoring of oral anticoagulant therapy^{3,4} and can be used to assess the protein synthesis capability of the liver in chronic or acute hepatic disorders. Thromboplastin L is of rabbit brain origin but resembles human preparations in its low International Sensitivity Index (ISI). The ISI of Thromboplastin L is approximately 1.1 and is calibrated against the WHO international reference preparation⁵. Thromboplastin L is particularly suited to the monitoring of oral anticoagulant therapy and, in conjunction with the appropriate factor deficient plasma, the measurement of factor activity in the extrinsic pathway. Tissue thromboplastin, in the presence of calcium ions, is an activator which initiates the extrinsic pathway of coagulation. When a mixture of tissue thromboplastin and calcium ions is added to normal citrated plasma, the clotting mechanism is activated, leading to a fibrin clot. If a deficiency exists within the extrinsic pathway, the time required for clot formation will be prolonged depending on the severity of the deficiency.

WARNINGS AND PRECAUTIONS

The reagents contained in this kit are for *in vitro* diagnostic use only – DO NOT INGEST. Wear appropriate personal protective equipment when handling all kit components. Refer to the product safety declaration for the link to appropriate hazard and precautionary statements where applicable. Dispose of components in accordance with local regulations.

COMPOSITION

Composition	Content	Description	Preparation
Thromboplastin L	2 x 5 mL (REF 5265HL) 8 x 5 mL (REF 5265L) 10 x 10 mL (REF 5267L)	Liquid Rabbit Brain Thromboplastin containing Calcium Chloride, stabilisers and preservatives.	The liquid, calcified thromboplastin is ready-for-use. No further calcium is required to carry out standard PT Assays. The contents of the vial should be mixed well before use. (5 minutes on roller).
Each kit contains Instructions For Use.			
Each kit contains lot specific reference values insert.			

ITEMS REQUIRED BUT NOT PROVIDED

The below products can be used in conjunction with Thromboplastin L:

REF 5519	ISI Calibrant Plasma Set
REF 5490	INR Reference Set

STORAGE, SHELF-LIFE AND STABILITY

Unopened reagents are stable until the given expiry date when stored under conditions indicated on the vial or kit label.

Thromboplastin L: Opened vials are stable for 2 months at *2–8°C, 5 days at *15°C (on-board Sysmex CA-1500) and 6 hours at *37°C (on-board AC-4 including reagent container and cap). A shift-use stability of 7 days (Sysmex CA-1500) can be achieved.

DO NOT FREEZE: Large clumps of particles or changes in expected values may indicate product deterioration.

SAMPLE COLLECTION AND PREPARATION

Plastic or siliconised glass should be used throughout. Blood (9 parts) should be collected into 3.2% or 3.8% sodium citrate anticoagulant (1 part). Separate plasma after centrifugation at 1500 x g for 15 minutes. Plasma should be kept at *18–24°C. Testing should be completed within 4 hours of sample collection, or plasma can be stored frozen at -20°C for 2 weeks or -70°C for 6 months. Thaw quickly at *37°C prior to testing. Do not keep at *37°C for more than 5 minutes⁶.

PROCEDURE

For accurate INR reporting, it is recommended to determine the laboratory specific ISI of the reagent with the testing system in use. The Helena Biosciences Europe ISI Calibrant Plasma Set (REF 5519) is recommended for this purpose^{7,8}. This should be performed for each new reagent batch. The Helena Biosciences Europe INR Reference Set (REF 5490) should be used to check for shifts in the local system ISI which have been noted with changes in laboratory temperature and post instrument servicing, amongst other local variances.

Manual Method

- Mix sufficient Thromboplastin L to complete the anticipated testing for the day and incubate at *37°C for no more than 4 hours.
- Prewarm 0.1 mL of the test plasma at *37°C for 2 minutes.
- Add 0.2 mL of freshly mixed thromboplastin reagent to the plasma while simultaneously starting a stopwatch.
- Note the time for clot formation to the nearest 0.1 seconds.

Automated Method

Refer to the appropriate instrument operator manual for detailed instructions or contact Helena Biosciences Europe for instrument specific application guides.

INTERPRETATION OF RESULTS

Results should be reported to the nearest 0.1 seconds and duplicates should agree within 5% of each other. %PT values can be interpolated from the calibration graph (%PT of PT Calibration plasmas versus measured clot time), which should be a straight line when plotted on log-log graph paper.

INR values can be calculated using the following formula: INR = (PT Time Patient / Mean Normal PT Time)^{ISI}

For clear guidance on the indications for and management of patients on warfarin, please refer to The British Society for Haematology, for their most current edition of 'Guidelines on oral anticoagulation with warfarin'. At time of printing this is the 2011 fourth edition⁹.

LIMITATIONS

The use of serial dilutions of a reference plasma for the %PT curve is not recommended as this can lead to discrepancies caused by the low fibrinogen in the reference plasma dilutions which are not reflected in patient samples having predominantly normal fibrinogen levels. Helena Biosciences Europe advise use of the 5504R %PT/Direct INR kit for this purpose.

QUALITY CONTROL

Each laboratory should establish a quality control program. Normal and abnormal control plasmas should be tested prior to each batch of patient samples, to ensure satisfactory instrument and operator performance. If controls do not perform as expected, patient results should be considered invalid.

Helena Biosciences Europe supplies the following controls available for use with this product:

REF 5186	Routine Control N
REF 5187	Routine Control A
REF 5183	Routine Control SA
REF 5490	INR Reference Set

REFERENCE VALUES

Reference values can vary between laboratories depending on the techniques and systems in use. For this reason each laboratory should establish its own reference ranges. This is particularly important for local ISI calibration. Using the Sysmex series of instruments, normal values ranging from 11.50 - 14.60 seconds; 0.930 - 1.160 INR; 79.10 - 112.80 %PT are typical.

PERFORMANCE CHARACTERISTICS

The following performance characteristics have been determined by Helena Biosciences Europe or their representatives using a Sysmex CA-1500 coagulation instrument. Each laboratory should establish its own performance data.

Reproducibility

Reproducibility						
	SD	CV (%)	SD	CV (%)	SD	CV (%)
Repeatability	0.07	0.59	0.24	1.09	0.45	1.11
Between-run	0.10	0.83	0.16	0.75	0.49	1.20
Between-day	0.04	0.32	0.06	0.27	0.25	0.62
Within-device / Laboratory	0.12	1.07	0.29	1.35	0.72	1.75

Interferences

Helena Thromboplastin L is insensitive to Heparin levels of up to 2 U/mL. Using a 5% interference threshold, there is no significant interference from Haemoglobin at concentrations up to 10 g/L. Using a 5% interference threshold, there is no significant interference from Bilirubin at concentrations up to 0.5 g/L for Thromboplastin L. Lipid interference testing demonstrates that lipid levels do not directly affect the clot time of the reagent up to 3.75g/L. Lipid concentrations in excess of this prevent clot detection.

Method Comparison

Comparison of clot time in seconds and INR values were determined using Thromboplastin L and Thromboplastin LI on 268 samples. The following correlations were obtained:

Thromboplastin L (Seconds) = 0.9911x + 0.1038 $r^2 = 0.9941$ n = 268

Thromboplastin L (INR) = 0.9853x + 0.0261 $r^2 = 0.9500$ n = 268

BIBLIOGRAPHY

- Quick AJ (1935) A Study of the Coagulation Defect in Hemophilia and Jaundice, *Am. J. Med. Sci.* 190: 501.
- Biggs R (1976) Human Blood Coagulation, Haemostasis and Thrombosis, 2nd Edition, Blackwell Scientific Publications, London.
- Hirsh J, Poller L, Deykin D, Levine J, Dalen JE (1989) Optimal Therapeutic Range for Oral Anticoagulants, *Chest*, 95: 55-115.
- Poller L (1986) Laboratory Control of Anticoagulant Therapy, *Sem. Thromb. Haemostasis*, 12: 13-19.
- World Health Organisation (1984) Expert Committee on Biological Standards, *Technical Series*, 700: 19.
- Clinical and Laboratory Standards Institute (2008) Collection, Transport and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Haemostasis Assays: Approved Guideline, 5th edn. CLSI: H21-A5.
- Poller L, Triplett DA, Hirsh J, Carroll J, Clarke K (1995) The value of plasma calibrants in correcting coagulometer effects on International Normalised Ratios (INR): An international multicentre study, *Amer. J. Clin. Pathol.*, 103: 358-365.
- Poller L, Triplett DA, Hirsh J, Carroll J, Clarke K (1995) A comparison of lyophilised artificially depleted plasmas and lyophilised plasmas from warfarin treated patients in correcting for coagulometer effects on International Normalised Ratios, *Amer. J. Clin. Pathol.*, 103: 366-371.
- Keeling D (2011) Guidelines on Oral Anticoagulation with warfarin: Forth Edition, *British Journal of Haematology*, 154(3): 311-324.

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BIOGRAFIA

- Quick AJ (1935) A Study of the Coagulation Defect in Hemophilia and Jaundice, *Am. J. Med. Sci.* **190**: 501.
- Biggs R (1976) Human Blood Coagulation, Haemostasis and Thrombosis, 2nd Edition, Blackwell Scientific Publications, London.
- Hirsh J, Poller L, Deykin D, Levine J, Dalen JE (1989) Optimal Therapeutic Range for Oral Anticoagulants, *Chest*, **95**: 55-115.
- Poller L (1986) Laboratory Control of Anticoagulant Therapy, *Sem. Thromb. Haemostasis*, **12**: 13-19.
- World Health Organisation (1984) Expert Committee on Biological Standards, *Technical Series*, **700**: 19.
- Clinical and Laboratory Standards Institute (2008) Collection, Transport and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Haemostasis Assays: Approved Guideline, 5th edn. CLSI: H21-A5.
- Poller L, Triplett DA, Hirsh J, Carroll J, Clarke K (1995) The value of plasma calibrants in correcting coagulometer effects on International Normalised Ratios (INR): An international multicentre study, *Amer. J. Clin. Pathol.*, **103**: 358-365.
- Poller L, Triplett DA, Hirsh J, Carroll J, Clarke K (1995) A comparison of lyophilised artificially depleted plasmas and lyophilised plasmas from warfarin treated patients in correcting for coagulometer effects on International Normalised Ratios, *Amer. J. Clin. Pathol.*, **103**: 366-371.
- Keeling D (2011) Guidelines on Oral Anticoagulation with warfarin: Forth Edition, *British Journal of Haematology*, **154**(3): 311-324.

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Declaration of Conformity

HL-7-0135DC DOI 2015/07 (7)

In Application of the Council Directive 98/79/EC on the approximation of the laws of the Member States relating to *In Vitro* Diagnostic Medical Devices & CE marking.

Declaration of conformance to applicable sections of Annex I - Essential Requirements and Annex III (EC Declaration of Conformity) imposed by sections 2 to 5. The below listed products are not classified under Annex II Lists A or B. Access to the appropriate technical files will be made available to the appropriate body in the event this is required.

Product Code	Description	GMDN Classification Code
5183	Routine Control SA	30590

I, the undersigned declare that the devices registered against the above GMDN Classification Code conforms to the said Directives.

Full Name: M.J. Stephenson

Title: Managing Director

Signed:

Date: 28 Jul 2015

Tel +44 (0)191 482 8440
Fax +44 (0)191 482 8442
info@helena-biosciences.com
www.helena-biosciences.com

Helena Biosciences Europe
Queensway South, Team Valley Trading Estate,
Gateshead, Tyne and Wear, NE11 0SD,
United Kingdom

Declaration of Conformity

HL-7-0137DC DOI 2015/07 (7)

In Application of the Council Directive 98/79/EC on the approximation of the laws of the Member States relating to *In Vitro* Diagnostic Medical Devices & CE marking.

Declaration of conformance to applicable sections of Annex I - Essential Requirements and Annex III (EC Declaration of Conformity) imposed by sections 2 to 5. The below listed products are not classified under Annex II Lists A or B. Access to the appropriate technical files will be made available to the appropriate body in the event this is required.

Product Code	Description	GMDN Classification Code
5186	Routine Control N	30590

I, the undersigned declare that the devices registered against the above GMDN Classification Code conforms to the said Directives.

Full Name: M.J. Stephenson

Title: Managing Director

Signed:

Date: 28 Jul 2015

Tel +44 (0)191 482 8440
Fax +44 (0)191 482 8442
info@helena-biosciences.com
www.helena-biosciences.com

Helena Biosciences Europe
Queensway South, Team Valley Trading Estate,
Gateshead, Tyne and Wear, NE11 0SD,
United Kingdom

Declaration of Conformity

HL-7-0512DC DOI 2015/08 (5)

In Application of the Council Directive 98/79/EC on the approximation of the laws of the Member States relating to *In Vitro* Diagnostic Medical Devices & CE marking.

Declaration of conformance to applicable sections of Annex I - Essential Requirements and Annex III (EC Declaration of Conformity) imposed by sections 2 to 5. The below listed products are not classified under Annex II Lists A or B. Access to the appropriate technical files will be made available to the appropriate body in the event this is required.

Product Code	Description	GMDN Classification Code
5556	Clauss Fibrinogen 50	55997

I, the undersigned declare that the devices registered against the above GMDN Classification Code conforms to the said Directives.

Full Name: M.J. Stephenson

Title: Managing Director

Signed:

Date: 12 Aug 2015

Tel +44 (0)191 482 8440

Fax +44 (0)191 482 8442

info@helena-biosciences.com

www.helena-biosciences.com

Helena Biosciences Europe

Queensway South, Team Valley Trading Estate,

Gateshead, Tyne and Wear, NE11 0SD,

United Kingdom

Declaration of Conformity

HL-7- 0640DC DOI 2015/07 (1)

In Application of the Council Directive 98/79/EC on the approximation of the laws of the Member States relating to *In Vitro* Diagnostic Medical Devices & CE marking.

Declaration of conformance to applicable sections of Annex I - Essential Requirements and Annex III (EC Declaration of Conformity) imposed by sections 2 to 5. The below listed products are not classified under Annex II Lists A or B. Access to the appropriate technical files will be made available to the appropriate body in the event this is required.

Product Code	Description	GMDN Classification Code
5504R	Calibration Plasma	55995

I, the undersigned declare that the devices registered against the above GMDN Classification Code conforms to the said Directives.

Full Name: M.J. Stephenson

Title: Managing Director

Signed:

Date: 30 Jul 2015

Tel +44 (0)191 482 8440
Fax +44 (0)191 482 8442
info@helena-biosciences.com
www.helena-biosciences.com

Helena Biosciences Europe
Queensway South, Team Valley Trading Estate,
Gateshead, Tyne and Wear, NE11 0SD,
United Kingdom

Declaration of Conformity

HL-7-0664DC DOI 2015/08 (1)

In Application of the Council Directive 98/79/EC on the approximation of the laws of the Member States relating to *In Vitro* Diagnostic Medical Devices & CE marking.

Declaration of conformance to applicable sections of Annex I - Essential Requirements and Annex III (EC Declaration of Conformity) imposed by sections 2 to 5. The below listed products are not classified under Annex II Lists A or B. Access to the appropriate technical files will be made available to the appropriate body in the event this is required.

Product Code	Description	GMDN Classification Code
5267L	Thromboplastin L	55983

I, the undersigned declare that the devices registered against the above GMDN Classification Code conforms to the said Directives.

Full Name: M.J. Stephenson

Title: Managing Director

Signed:

Date: 06 Aug 2015

Tel +44 (0)191 482 8440
Fax +44 (0)191 482 8442
info@helena-biosciences.com
www.helena-biosciences.com

Helena Biosciences Europe
Queensway South, Team Valley Trading Estate,
Gateshead, Tyne and Wear, NE11 0SD,
United Kingdom

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QUALITY MANAGEMENT SYSTEM - ISO 13485:2016 & EN ISO 13485:2016

This is to certify that:


Helena Laboratories (UK) Ltd
trading as Helena Biosciences Europe
Queensway South
Team Valley Trading Estate
Gateshead
Tyne and Wear
NE11 0SD
United Kingdom

Holds Certificate Number:

MD 69326

and operates a Quality Management System which complies with the requirements of ISO 13485:2016 & EN ISO 13485:2016 for the following scope:

The design, manufacture, supply, servicing and repair of in-vitro diagnostic devices, molecular biology products, immunochemistry products and medical laboratory equipment and consumables.

Gary E Slack

For and on behalf of BSI:

Gary E Slack, Senior Vice President - Medical Devices

Original Registration Date: 2002-10-25

Effective Date: 2021-04-14

Latest Revision Date: 2021-04-13

Expiry Date: 2024-04-13

Page: 1 of 2

...making excellence a habit.TM

This certificate was issued electronically and remains the property of BSI and is bound by the conditions of contract.
An electronic certificate can be authenticated [online](#).
Printed copies can be validated at www.bsigroup.com/ClientDirectory

Information and Contact: BSI, Kitemark Court, Davy Avenue, Knowlhill, Milton Keynes MK5 8PP. Tel: + 44 345 080 9000
BSI Assurance UK Limited, registered in England under number 7805321 at 389 Chiswick High Road, London W4 4AL, UK.
A Member of the BSI Group of Companies.

Certificate No: **MD 69326**

Location

Registered Activities

Helena Laboratories (UK) Ltd
trading as Helena Biosciences Europe
Sunderland Enterprise Park
Colima Avenue
Sunderland
SR5 3XB
United Kingdom

The design, manufacture, supply, servicing and repair of
in-vitro diagnostic devices, molecular biology products,
immunochemistry products and medical laboratory equipment
and consumables.

Helena Laboratories (UK) Ltd
trading as Helena Biosciences Europe
Queensway South
Team Valley Trading Estate
Gateshead
Tyne and Wear
NE11 0SD
United Kingdom

The design, manufacture, supply, servicing and repair of
in-vitro diagnostic devices, molecular biology products,
immunochemistry products and medical laboratory equipment
and consumables.

Original Registration Date: 2002-10-25

Effective Date: 2021-04-14

Latest Revision Date: 2021-04-13

Expiry Date: 2024-04-13

Page: 2 of 2

This certificate was issued electronically and remains the property of BSI and is bound by the conditions of contract.
An electronic certificate can be authenticated [online](#).

Printed copies can be validated at www.bsigroup.com/ClientDirectory

Information and Contact: BSI, Kitemark Court, Davy Avenue, Knowlhill, Milton Keynes MK5 8PP. Tel: + 44 345 080 9000
BSI Assurance UK Limited, registered in England under number 7805321 at 389 Chiswick High Road, London W4 4AL, UK.
A Member of the BSI Group of Companies.

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**"VACUTEST KIMA" vacuum blood and biological liquids
collection tubes in plastic.**

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Hereby we declare

*under our sole responsibility that the above mentioned devices meet the applicable provisions of the Directive 98/79/EC
as amended on "In Vitro Diagnostic Medical Devices".*

All the supporting documents, as required by Annex III, in order to prove conformity to the Essential
Requirements as listed in Annex I, are retained under the premises of the Manufacturer

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firma
signature

Instrument Training

Vital Scientific BV hereby declares that the participant has attended a four days seminar for service engineers and the participant is now a certified engineer for the declared instruments.

Participant: Mr. A. Legun

Company: Global Biomarketing Group-Moldova SRL
Moldova

Instrument: Vitalab: XL Series
E Series
Junior Series
Dry ISE
Micro Series
ProXS

Date of training: April 20th – April 23rd, 2010

System Support Manager:

A blue ink signature of the name "Jan Oostendorp".

Jan Oostendorp

System Support Engineer:

A blue ink signature of the name "Frank v.d. Korput".

Frank v.d. Korput

EC DECLARATION OF CONFORMITY

Lorne Laboratories Ltd declares that the following in vitro diagnostic reagent:

Product Name	Catalogue Number
Anti-A Monoclonal	600010

has been classified as List A (Directive 98/79/EC, Annex II) and complies with the essential requirements and provisions of Directive 98/79/EC of the European Parliament and of the Council (also SI 2002 No.618 which transposes the requirements of Directive 98/79/EC) and the Commission Decision on Common Technical Specifications 2009/108/EC.

and is in conformity with the national standards transposing harmonised standards:

- BS EN ISO 13485:2012
- BS EN 13612:2002
- BS EN 13641:2002
- BS EN ISO 14971:2012
- BS EN ISO 15223-1:2016
- BS EN ISO 18113-2:2011
- BS EN ISO 23640:2015

The conformity assessment procedure performed was in accordance with Annex IV of Directive 98/79/EC and was carried out by UL International (UK) Ltd, Wonersh House, The Guildway, Old Portsmouth Road, Guildford, Surrey GU3 1LR, United Kingdom, Notified Body Number 0843.

The certificates issued by UL-UK Ltd to show compliance are numbers 354.170425 and 355.130523.

This declaration of conformity is issued under the sole responsibility of Lorne Laboratories Ltd and is valid from 23 May 2017.

Eddy Velthuis
Technical Director

File No A12241;
ISO 13485:2003; ISO 9001:2008

Lorne Laboratories Limited
Unit 1 Cutbush Park Industrial Estate
Danehill, Lower Earley
Berkshire RG6 4UT United Kingdom

Tel: +44 (0) 118 921 2264
Fax: +44 (0) 118 986 4518
Email: info@lornelabs.com
www.lornelabs.com

Registered office as above. Registered in England No. 04540797. VAT No. 800 3655 66

EC DECLARATION OF CONFORMITY

Lorne Laboratories Ltd declares that the following in vitro diagnostic reagent:

Product Name	Catalogue Number
Anti-B Monoclonal	610010

has been classified as List A (Directive 98/79/EC, Annex II) and complies with the essential requirements and provisions of Directive 98/79/EC of the European Parliament and of the Council (also SI 2002 No.618 which transposes the requirements of Directive 98/79/EC) and the Commission Decision on Common Technical Specifications 2009/108/EC.

and is in conformity with the national standards transposing harmonised standards:

- BS EN ISO 13485:2012
- BS EN 13612:2002
- BS EN 13641:2002
- BS EN ISO 14971:2012
- BS EN ISO 15223-1:2016
- BS EN ISO 18113-2:2011
- BS EN ISO 23640:2015

The conformity assessment procedure performed was in accordance with Annex IV of Directive 98/79/EC and was carried out by UL International (UK) Ltd, Wonersh House, The Guildway, Old Portsmouth Road, Guildford, Surrey GU3 1LR, United Kingdom, Notified Body Number 0843.

The certificates issued by UL-UK Ltd to show compliance are numbers 354.170425 and 355.130523.

This declaration of conformity is issued under the sole responsibility of Lorne Laboratories Ltd and is valid from 23 May 2017.

A handwritten signature in black ink, appearing to read "Eddy Velthuis".

Eddy Velthuis
Technical Director

File No A12241;
ISO 13485:2003; ISO 9001:2008

Lorne Laboratories Limited
Unit 1 Cutbush Park Industrial Estate
Danehill, Lower Earley
Berkshire RG6 4UT United Kingdom

Tel: +44 (0) 118 921 2264
Fax: +44 (0) 118 986 4518
Email: info@lornelabs.com
www.lornelabs.com

Registered office as above. Registered in England No. 04540797. VAT No. 800 3655 66

EC DECLARATION OF CONFORMITY

Lorne Laboratories Ltd declares that the following in vitro diagnostic reagent:

Product Name	Catalogue Number
Anti-D Duoclone Monoclonal	740010

has been classified as List A (Directive 98/79/EC, Annex II) and complies with the essential requirements and provisions of Directive 98/79/EC of the European Parliament and of the Council (also SI 2002 No.618 which transposes the requirements of Directive 98/79/EC) and the Commission Decision on Common Technical Specifications 2009/108/EC.

and is in conformity with the national standards transposing harmonised standards:

- BS EN ISO 13485:2012
- BS EN 13612:2002
- BS EN 13641:2002
- BS EN ISO 14971:2012
- BS EN ISO 15223-1:2016
- BS EN ISO 18113-2:2011
- BS EN ISO 23640:2015

The conformity assessment procedure performed was in accordance with Annex IV of Directive 98/79/EC and was carried out by UL International (UK) Ltd, Wonersh House, The Guildway, Old Portsmouth Road, Guildford, Surrey GU3 1LR, United Kingdom, Notified Body Number 0843.

The certificates issued by UL-UK Ltd to show compliance are numbers 354.170425 and 355.130523.

This declaration of conformity is issued under the sole responsibility of Lorne Laboratories Ltd and is valid from 23 May 2017.

Eddy Velthuis
Technical Director

File No A12241;
ISO 13485:2003; ISO 9001:2008

Lorne Laboratories Limited
Unit 1 Cutbush Park Industrial Estate
Danehill, Lower Earley
Berkshire RG6 4UT United Kingdom

Tel: +44 (0) 118 921 2264
Fax: +44 (0) 118 986 4518
Email: info@lornelabs.com
www.lornelabs.com

Registered office as above. Registered in England No. 04540797. VAT No. 800 3655 66

EC CERTIFICATE

Lorne Laboratories Ltd

Unit 1 Cutbush Park Industrial Estate,
Danehill, Lower Earley, Berkshire RG6 4UT, UK

EC Certificate - Full Quality Assurance System Approval Certificate

Annex IV, (excluding sections 4 and 6) of Council Directive 98/79/EC on In Vitro Diagnostic Medical Devices

Scope of Certificate:

The design and manufacture of in vitro diagnostic reagents for identification of blood groups

Device Classification:

Annex II, List A and B

Device Descriptions:

Please refer to Attachment 1

Model:

Please refer to Attachment 1

File Number A12241

Cycle Start Date 23 May 2017

Certificate No. 354.170425

Effective Date 23 May 2017

Expiry Date 22 May 2022

Authorised by

B. Rodgers
Certification Manager
For and on Behalf of UL International (UK) Ltd

We hereby declare that an examination of the full quality assurance system has been carried out per report 11640248, following the requirements of the national legislation to which the undersigned is subject, transposing Annex IV (with the exemption of sections 4 and 6) of Council Directive 98/79/EC on In Vitro Diagnostic Medical Devices. We certify that the full quality assurance system conforms with the relevant provisions of the aforementioned directive and is subject to periodic surveillance as required by 98/79/EC, Annex IV, Section 5. For Annex II, List A devices where they are covered by this certificate, an EC Design Examination certificate according to 98/79/EC, Annex IV, Section 4 is required. This certificate is issued with 1 attachment listing model numbers.

Notified Body

0843

EC CERTIFICATE

Lorne Laboratories Ltd

Unit 1 Cutbush Park Industrial Estate,
Danehill, Lower Earley, Berkshire RG6 4UT, UK

Attachment 1 of 1

The products detailed below are covered under the scope of this certificate

Device Description	Model	Classification
Anti-A Monoclonal	600005/600010/600000	Annex II List A
Anti-B Monoclonal	610005/610010/610000	Annex II List A
Anti-A,B Monoclonal	620005/620010/620000	Annex II List A
Anti-C Monoclonal	690005	Annex II List A
Anti-E Monoclonal	691005	Annex II List A
Anti-c Monoclonal	692005	Annex II List A
Anti-e Monoclonal	693005	Annex II List A
Anti-K Monoclonal	760005/760010	Annex II List A
Anti-D Clone 2 Monoclonal	710010/710000	Annex II List A
Anti-D Clone 1 Monoclonal	730010/730000	Annex II List A
Anti-D Duoclone Monoclonal	740010/740000	Annex II List A
Anti-Jka Polyclonal	323002/323000	Annex II List B
Anti-Jkb Polyclonal	324002/324000	Annex II List B
Anti-Fyb Polyclonal	317002/317000	Annex II List B
AHG Elite Clear	415010/415100/415000	Annex II List B
AHG Elite Green	435010/435100/435000	Annex II List B
Anti-Fya Monoclonal	774000/774002	Annex II List B
Anti-C+D+E Monoclonal	700005/700010/700000	Annex II List A
Anti-Human IgG Clear	401010/401000	Annex II List B
Anti-Human IgG Green	402010/402000	Annex II List B
Monoclonal Rh Control	640010	Annex II List A
Monoclonal D Negative Control	650010	Annex II List A

File Number A12241
Certificate No. 354.170425

Cycle Start Date 23 May 2017
Effective Date 23 May 2017
Expiry Date 22 May 2022

Authorised by

B. Rodgers
Certification Manager
For and on Behalf of UL International (UK) Ltd

Notified Body

0843

IVDD A4 S3 FQ 00-NB-F0051 Issue: 6.0

UL International (UK) Limited
Womersh House, The Guildway, Old Portsmouth Road,
Guildford, Surrey, GU3 1LR, United Kingdom

Training certificate

helena
Biosciences Europe

This is to certify that

Sergiu Sorocovici

from

IM Global Biomarketing Group

has received training on the following:

Electrophoresis products: SAS-1/2,V8

Haemostasis products: C-series, AC-4, AggRAM and reagents

Service training: AC-4

Signed:

Date: 31st October - 4th November 2011

Tel +44 (0)191 482 8440
Fax +44 (0)191 482 8442

info@helena-biosciences.com
techsupport-hs@helena-biosciences.com

Queensway South, Team Valley Trading Estate,
Gateshead, Tyne and Wear, NE11 0SD, United Kingdom