MENDEC[®] SpINE Kit

HIGH BARIUM
CONTENT

for OUS market

TECRES[®]

ADVANCING HIGH TECHNOLOGY

TECRES S.p.A.

Via A. Doria, 6 - 37066 Sommacampagna - Verona - ITALY

Telefono +39.045.9217311 - Fax +39.045.9217330

info@tecres.it - www.tecres.it

for Tecres
internal use

CE
2797

MENDEC® Spine Kit

IT

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MENDEC® Spine Kit

EN

INSTRUCTIONS FOR USE To the medical Staff's attention

INSTRUCTIONS FOR USE

MENDEC SPINE KIT is a disposable, sterile system for the filling of pathologic vertebral bodies. The kit components are sterilized with ethylene oxide. The liquid in the vial is sterilized using filtration. The system allows acrylic resin to flow into the pathological vertebral body using an injection system. Mendec Spine Kit device must only be used with needles specific for the selected procedures with Luer-lok type connection.

DESCRIPTION OF THE DEVICE

MENDEC SPINE KIT consists of:

- MENDEC SPINE radio-opaque, high-fluidity acrylic resin
 - 20 g container of powder
 - Polymethylmethacrylate 67.5% w/w
 - Barium sulphate 30.0% w/w
 - Benzoyl peroxide 2.5% w/w
 - 9.4 g phial of liquid
 - Methylmethacrylate 99.1% w/w
 - N,N-dimethyl-p-toluidine 0.9% w/w
 - Hydroquinone 75 ppm
- Injection system:
 - Spine Gun
 - Extension tube

The injection system uses Luer-lok type connections

CONTRAINDICATIONS

Active infections

Emorrhagic diatesis

Extended vertebral destruction, with collapse superior to 2/3 of the standard thickness, destruction of the posterior wall with epidural extension of the pathologic tissue and clinical signs of medullar compression.

The use of MENDEC SPINE KIT must be carefully considered in presence of hypersensitivity to monomer or to any of the other components of the resin.

USE DURING PREGNANCY, BREASTFEEDING AND IN CHILDREN

There are no tests which demonstrate the utilisation safety of MENDEC SPINE KIT during pregnancy or breast-feeding. MENDEC SPINE KIT should not be used in the first three months of pregnancy; for the remaining gestation period, should only be used in life-endangering situations.

The use in children is advised only when it is believed impossible to operate through other forms of intervention.

PRECAUTIONS FOR USE

- Read this instruction booklet very carefully.
- Make sure that the inner package is undamaged and that the components are integral, i.e. the powder should not present yellow or brown discolouring and the liquid should not be syrupy.
- Mendec Spine Kit must only be used with needles specific for the selected procedures. It is advised the use of Tecres' Spine Needles. If needles of other manufacturers are used, do check the presence of Luer-Lok type connections and be sure that they can be tightly assembled to the safety locknut present on the Extension tube of Mendec Spine kit.
- Temperature has a major effect on the preparation characteristics of any acrylic resin. Temperatures of more than 23°C for the product, the preparation accessories and the environment accelerate the various stages in the preparation procedure. Lower temperatures retard the preparation stages. **Before using MENDEC SPINE KIT it is strongly advised to make sure that the package is stored at a temperature of 23°C ± 1°C for the previous 24 hours.**
- Do not open the phial of liquid over the mixing bowl to avoid the risk of glass fragments entering the dough.
- NEVER add other substances or foreign bodies to the acrylic resin.

To ensure safe application of MENDEC SPINE KIT, surgeons must be aware of the properties of the product, its characteristics, applicational limitations and correct application technique. It is strongly recommended that the surgical team carry out practical trials prior to use in patients under the same instrumental and environmental conditions.

Clinical data demonstrate the need to maintain strictly aseptic surgical techniques. It is important to keep in mind that the possible deep infection of a surgical wound is a serious risk for the successful outcome of the surgical operation.

Prolonged respiration of the vapours may cause drowsiness; excessive exposure to monomer vapours may irritate the respiratory apparatus and the eyes. Avoid monomer contact with the skin and mucous membranes. Cases of

contact dermatitis have been observed in susceptible subjects. It is therefore advisable to wear a second pair of surgical gloves and scrupulously observe the instructions for mixing the components in order to reduce the possibility of reactions caused by hypersensitivity. Do not allow the liquid component to come into direct contact with accessories in rubber or elastomers. The vapours of the liquid component may affect soft contact lenses. The liquid component is flammable and volatile and for this reason the operating theatre must be correctly ventilated. The liquid component and/or its vapours must never be directly exposed to naked flames or heated materials.

Caution: Acrylic resins reach temperatures higher than physiological temperatures during the exothermic polymerisation reaction.

Following surgery, if any form of infection should arise, patients must immediately consult their doctors to reduce the risk of infection.

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Times indicated in the table were obtained in the laboratory under temperature controlled conditions of 23°C±1°C using the injection kit system (including the extension tube) and with a bevelled needle (gauge 11, length 100 mm).

Any **increase in temperature** of the working environment, the product and the preparation accessories over 23°C **reduces the waiting and hardening times** of the acrylic resin. Equally, **lower temperatures increase such times**.

Operation	Definition	Phase duration
Mixing	Mixing of the components	60"
Delivery device filling	Transfer into the Spine Gun	60"
Waiting	The cement cannot be used	180"
Working	The cement can be delivered	390"
Hardening/Polimerization	The cement hardens increasing its viscosity and cannot be delivered anymore. Then the exothermic reaction takes place	690"

During the application of the resin, the viscosity of the resin can be fitted to the implant site by **anticipating or delaying the moment of application**. Do not exceed the indicated times since the resin's fluidity reduces with time and it becomes more difficult to eject.

INSTRUCTIONS FOR USE

- **Positioning the Spine Needle:** the needle and the stylet should be positioned inside the vertebral body using the access route that is best suited to the pathology and the condition of the vertebra. Correct positioning is checked using a fluoroscopic test and/or a CAT scan. When positioning is complete, leave the stylet inserted in the needle to keep the access route to the vertebral body closed.
- **Resin preparation:** Set the container of powder on a flat surface and open the cap. Open the vial and pour all of the liquid into the powder. Close the container and manually shake it until a liquid and homogenous mass is achieved.

WARNING: never arbitrarily change the ratio between liquid and solid components. Do not re-sterilize and/or re-use the device. It is designed for single-use on a single patient. Avoid dividing the product into two or more portions, in order to use it other times. This would be a re-utilization and may lead to an error in the proportion of the powder and liquid components. It can also cause a sterility loss. Re-sterilization should not be carried out since it can cause infection risks for the patient. Re-sterilization can also alter the morphology and the mechanical features of the device, causing a malfunction of the same with serious risks for the patient's health. All residues must be considered waste surgical material and therefore should be eliminated at the end of the operation.

- **Loading the Spine Gun:** Open the luer cap. Make sure the Spine Gun's threaded staff is completely inserted. Connect the Gun with the luer lock connection on the cap of the resin container, tightening the ring nut on the tip of the Gun onto the luer connection. Keep the Spine Gun grip pressed (which circumvents the safety catch). Turn the Spine Gun and the connected container upside down and pull the threaded staff until the Spine Gun is completely full. To facilitate the release of any air bubbles, after having released the Spine Gun grip, rotate the handle of the threaded staff keeping the tip of the Spine Gun upwards.
- **Removal of the powder container:** Turn the device upside down again and turn the ring nut on the tip of the Spine Gun to free the Luer-loc, and then pull out the powder container.
- **Application of the extension tube:** Bring the Luer-loc on the extension tube and the tip of the Spine Gun together. Tighten the ring nut on the tip of the Spine Gun onto the extension tube connector tight enough to create a solid grip.
- **Filling the extension tube:** Once the extension tube is attached to the Spine Gun, turn the knob of the Spine Gun's piston until the extension tube fills completely with resin. Let a few drops of resin come out of from the extension tube before connecting it to the injection cannula. This first leakage needs to be eliminated.
- **Connecting the needle to the extension tube:** Remove the stylet from the needle and connect the extension tube to the needle, tightening the tube's ring nut to the mouth of the needle.

- **Supplying resin:** Resin must be supplied to the vertebral body by turning the handle of the threaded system clockwise and continuously monitoring the filling of the vertebra using fluoroscopy. Dispense regularly by slow rotation of the knob of the threaded system. Max injection speed 2ml/min. Rotate the Spine Gun handle the opposite way to slow down or block the release of resin.
- **Removing the needle:** When all the resin has been introduced, remove the extension tube from the needle, reinsert the stylet completely and turn the needle a few times to free it from the resin. Then wait approximately one minute before removing the needle.

WARNINGS

- Store at a temperature below 25°C and away from light.
- Sterility is assured only if the unit container is not damaged or opened.
- Do not re-sterilise any of the components.
- Do not use the product if the powder has a yellowish or brownish colour or if the liquid is syrupy. These two conditions indicate that the product has not been stored correctly.

NEGATIVE EFFECTS

The following negative reactions may arise with the use of the acrylic resin, but they are not directly due to the acrylic resin itself.

The surgeon must be aware of these possible events and be ready to treat them if they appear.

- cardiac arrest
- myocardial infarct
- pulmonary embolism
- temporary drop in blood pressure
- haemorrhage and haematoma
- infection of surface/deep surgical lesion
- short-term irregularity of the cardiac rhythm
- sudden death
- diffusion of the resin outside the vertebral body: in the perivertebral veins (pulmonary embolism), in the epidural plex (mielopathy, radiculopathy), in the intervertebral disc.

IMPORTANT INFORMATION FOR THE PHYSICIAN

The use of bone cement requires extremely close cooperation between the surgeon and the anaesthetist.

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In some cases, events defined as bone cement implantation syndrome (BCIS) may occur, which are characterised by a series of clinical characteristics that include hypoxia, hypotension, cardiac arrhythmia, increased pulmonary vascular resistance and cardiac arrest, which must be controlled with the methods in use in modern anaesthesiology.

Even though these phenomena are rare, they are described in literature: percutaneous injection of bone cement can increase the intramedullary pressure and hence force both the bone marrow and the bone cement into the perivertebral venous system and these may migrate to the right side of the heart and into pulmonary circulation (pulmonary embolism) resulting in BCIS. BCIS has a very wide spectrum of severity and in many cases patients with pulmonary embolism detected during a chest X-ray remain asymptomatic.

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CONTENTS STERILE UNLESS PACKAGE IS DAMAGED OR OPENED

CAUTION: EXTERNAL SURFACE OF FIRST CONTAINER INSIDE IS NOT STERILE

FLAMMABLE LIQUID

CONSULT INSTRUCTIONS FOR USE

STERILIZED BY ETHYLENE OXIDE

MANUFACTURER

STERILE USING ASEPTIC TECHNIQUES

KEEP AWAY FROM SUNLIGHT

STERILE

DO NOT RESTERILISE

STORE AT A TEMPERATURE BELOW 25°C

DO NOT RE-USE

CAUTION

BATCH NUMBER

CATALOGUE NUMBER

USE BY

SYMBOLS


For further product information, please contact Customer Service:

Tecres S.p.A.
Via A. Doria, 6
37066 Sommacampagna
Verona Italy
Phone: +39 0459217311
Fax: +39 0459217330

INFORMATION

- The handling characteristics of bone cements are affected by operating room conditions, including the room temperature, temperature of the cement components prior to mixing, humidity, the geometry of the mixing apparatus, time spent mixing and the geometry of the delivery device. Any change in one or more of these conditions can alter the handling characteristics of bone cement.
- If the resin is injected early, it will result in a highly viscous material.
- If the resin is injected late, it will result in a highly viscous material.

CAUTION

The cement's preparation and application is strongly influenced by the temperature in the operating room. The temperature's effect on cement setting time was evaluated with a laboratory test. A graph on setting time according to temperature is provided for ease of reference (data was obtained in controlled environmental and storage conditions subjected to standard deviation). In addition to temperature and humidity, different factors can influence the cement's setting time. These include the mixing technique (speed, use of mixer), the thoroughness of mixing, the utilization of the entire liquid and solid component, the inclusion of external substances inside the cement such as blood, saline solution, etc.

THE EFFECT OF TEMPERATURE ON PREPARATION AND APPLICATION OF MENDEC SPINE

During the application of the resin, the viscosity of the resin can be fitted to the implant site by anticipating or delaying the moment of application. Do not exceed the time period indicated, since the fluidity of the resin decreases over time, and it becomes more difficult to dispense.

Operation	Definition	Phase duration
Mixing	Mixing of the components	60"
Delivery Device Filling	Transfer into Spine Gun	60"
Waiting	The cement cannot be used	180"
Working	The cement can be delivered	390"
Hardening/Setting	The cement hardens increasing its viscosity and cannot be delivered anymore. Then the exothermic reaction takes place	690"

Times indicated in the instructions were obtained by performing tests under temperature controlled conditions of $23^{\circ}\text{C} \pm 1^{\circ}\text{C}$ with the injection system of the kit (extension tube included) and a needle beveled tip 11 gauge 100 mm. Any increase in temperature of the working environment, the product and the preparation accessories over 23°C reduces the waiting and hardening times of the acrylic resin. Equally, lower temperatures increase such times.

TIMING FOR THE PREPARATION AND APPLICATION

- **Removal of the powder container:** Turn the device upside down again and turn the ring nut on the tip of the Spine Gun to free the Luer-lock, and then pull out the powder container.
 - **Application of the extension tube:** Bring the Luer-lock on the extension tube and the tip of the Spine Gun together. Tighten the ring nut on the tip of the Spine Gun onto the extension tube connector tight enough to create a solid grip.
 - **Filling the extension tube:** Once the extension tube is attached to the Spine Gun, turn the knob of the Spine Gun's piston until the extension tube fills completely with resin. Let a few drops of resin come out of from the extension tube before connecting it to the injection cannula. This first leakage needs to be eliminated.
 - **Connecting the needle to the extension tube:** Remove the stylet from the needle and connect the extension tube to the needle, tightening the tube's ring nut to the mouth of the needle.
 - **Supplying resin:** Resin must be supplied to the vertebral body by turning the handle of the threaded system clockwise and continuously monitoring the filling of the vertebra using fluoroscopy.
 - **Dispense regularly** by slow rotation of the knob of the threaded system. Max injection speed 2ml/min. Rotate the Spine Gun handle the opposite way to slow down or block the release of resin.
 - **Caution:** never deliver acrylic resin into the vertebral body without the use of fluoroscopic guidance.
 - **Removing the needle:** When all the resin has been introduced, remove the extension tube from the needle, reinsert the stylet completely and turn the needle a few times to free it from the resin. Then wait approximately one minute before removing the needle.
 - **Loading the Spine Gun:** Open the luer cap. Make sure the Spine Gun's threaded staff is completely inserted. Connect the Gun with the luer lock connection on the cap of the resin container, tightening the ring nut on the tip of the Gun onto the luer connection. Keep the Spine Gun grip pressed (which circumvents the safety catch). Turn the Spine Gun and the connected container upside down and pull the threaded staff until the Spine Gun is completely full. To facilitate the release of any air bubbles, after having released the Spine Gun grip, rotate the handle of the threaded staff keeping the tip of the Spine Gun upwards.
 - **Removal of the powder container:** Turn the device upside down again and turn the ring nut on the tip of the Spine Gun to free the Luer-lock, and then pull out the powder container.
 - **Application of the extension tube:** Bring the Luer-lock on the extension tube and the tip of the Spine Gun together. Tighten the ring nut on the tip of the Spine Gun onto the extension tube connector tight enough to create a solid grip.
 - **Filling the extension tube:** Once the extension tube is attached to the Spine Gun, turn the knob of the Spine Gun's piston until the extension tube fills completely with resin. Let a few drops of resin come out of from the extension tube before connecting it to the injection cannula. This first leakage needs to be eliminated.
 - **Connecting the needle to the extension tube:** Remove the stylet from the needle and connect the extension tube to the needle, tightening the tube's ring nut to the mouth of the needle.
 - **Supplying resin:** Resin must be supplied to the vertebral body by turning the handle of the threaded system clockwise and continuously monitoring the filling of the vertebra using fluoroscopy.
 - **Dispense regularly** by slow rotation of the knob of the threaded system. Max injection speed 2ml/min. Rotate the Spine Gun handle the opposite way to slow down or block the release of resin.
 - **Caution:** never deliver acrylic resin into the vertebral body without the use of fluoroscopic guidance.
 - **Removing the needle:** When all the resin has been introduced, remove the extension tube from the needle, reinsert the stylet completely and turn the needle a few times to free it from the resin. Then wait approximately one minute before removing the needle.
- All residues must be considered waste surgical material and therefore should be eliminated at the end of the operation.
- serious risks for the patient's health.
- Re-sterilization should not be carried out since it can cause infection risks for the patient. Re-sterilization can also alter the morphology and the mechanical features of the device, causing a malfunction of the same with cause a sterility loss.
- be a re-sterilization and may lead to an error in the proportion of the powder and liquid components. It can also avoid dividing the product into two or more portions, in order to use it other times. This would not re-sterilize and/or re-use the device. It is designed for single-use on a single patient.

INSTRUCTIONS FOR USE

- **Positioning the Spine Needle:** the needle and the stylet should be positioned inside the vertebral body using the access route that is best suited to the pathology and the condition of the vertebra. Correct positioning is checked using a fluoroscopic test and/or a CAT scan. When positioning is complete, leave the stylet inserted in the needle to keep the access route to the vertebral body closed.
- **Resin preparation and mixing:** Set the container of the powder on a flat surface and open the cap. Break open the vial and pour all the liquid into the powder container. Do not open the vial of liquid over the powder container to avoid the risk of glass fragments entering the dough. Close the container and manually shake it for 1 minute until a liquid and homogeneous mass is achieved. **Never add other substances or foreign bodies to the acrylic resin. Moreover, never arbitrarily modify the ratios between the liquid and solid components.**

pulmonary embolism detected during a chest X-ray remain asymptomatic.

(pulmonary embolism) resulting in BCIS. BCIS has a very wide spectrum of severity and in many cases patients with perivertebral venous system and these may migrate to the right side of the heart and into pulmonary circulation increase the intramedullary pressure and hence force both the bone marrow and the bone cement into the vertebral body. Even though these phenomena are rare, they are described in literature: percutaneous injection of bone cement can increase the intramedullary pressure and hence force both the bone marrow and the bone cement into the vertebral body. Even though these phenomena are rare, they are described in literature: percutaneous injection of bone cement can increase the intramedullary pressure and hence force both the bone marrow and the bone cement into the vertebral body. Even though these phenomena are rare, they are described in literature: percutaneous injection of bone cement can increase the intramedullary pressure and hence force both the bone marrow and the bone cement into the vertebral body.

The anaesthetist must be told when the bone cement is implanted during the operation.

In some cases, events defined as bone cement implantation syndrome (BCIS) may occur, which are characterised by a series of clinical characteristics that include hypoxia, hypotension, cardiac arrhythmia, increased pulmonary vascular resistance and cardiac arrest, which must be controlled with the methods in use in modern anaesthesiology. In some cases, events defined as bone cement implantation syndrome (BCIS) may occur, which are characterised by a series of clinical characteristics that include hypoxia, hypotension, cardiac arrhythmia, increased pulmonary vascular resistance and cardiac arrest, which must be controlled with the methods in use in modern anaesthesiology.

The use of bone cement requires extremely close cooperation between the surgeon and the anaesthetist. The use of bone cement requires extremely close cooperation between the surgeon and the anaesthetist.

arbitrarily modify the ratios between the liquid and solid components.

- Never add other substances (including antibiotics) or foreign bodies to the acrylic resin. Moreover, never arbitrarily modify the ratios between the liquid and solid components.
- these circulating concentrations have not been established.
- concentrations of methyl methacrylate/methacrylic acid although changes in blood pressure in response to methacrylic acid, and that a significant fraction of the circulating methacrylic acid is in the form of circulating poly(methyl methacrylate) bone cements. Recent data indicate that the monomer undergoes rapid hydrolysis to methacrylic acid, and that a significant fraction of the circulating methacrylic acid is in the form of circulating poly(methyl methacrylate) bone cements. Recent data indicate that the monomer undergoes rapid hydrolysis to methacrylic acid, and that a significant fraction of the circulating methacrylic acid is in the form of circulating poly(methyl methacrylate) bone cements.
- Adverse reactions affecting the cardiovascular system have been reported following the application of safety instructions that appear in the directions for use.
- The physician is responsible for any complication or harmful consequences which may result from an erroneous indication or operating technique, from inappropriate use of the material, or non-observation of the instructions for use.
- Percutaneous vertebroplasty procedures should only be performed in medical settings in which emergency decompressive surgery is available.

IMPORTANT PHYSICIAN INFORMATION

Although the majority of these adverse events present early within the post-operative period, there have been some reports of diagnoses beyond a year or more after the procedure.

Other reported adverse events for acrylic bone cements intended for vertebroplasty include:

Leakage of the bone cement beyond the site of its intended application with introduction into the vascular system resulting in embolism of the lung and/or heart or other clinical sequelae.

- sudden death
- short-term irregularity of the cardiac rhythm
- the significant downward force exerted during needle insertion
- no fracture in patients with diffuse osteopenia, especially during thoracic vertebroplasty procedures, due to
- pulmonary embolism
- pneumothorax
- pneumonia
- myocardial infarction
- intercostal neuralgia
- infection of surface/deep surgical lesion
- haemorrhage and haematoma
- fracture of a pedicle
- epidural plex (mielopathy, radiculopathy) or the intervertebral disc
- diffusion of the resin outside the vertebral body, such as the perivertebral veins (pulmonary embolism), the
- compression of the spinal cord with paralysis or loss of feeling
- collapse of a vertebra adjacent to the one injected due to an osteoporotic disease
- cerebrovascular accident
- cement extrusion into soft tissue
- cardiac embolism
- cardiac arrest

vertebroplasty include:

Serious adverse events, some with fatal outcome, associated with the use of acrylic bone cements for

ADVERSE EVENTS

devices in surgical sites near freshly implanted bone cements has been reported.

- Adequately ventilate the operating room to eliminate as much monomer vapor as possible. The liquid monomer is highly volatile and flammable. Ignition of monomer fumes caused by the use of electrocautery symptoms of infection occur.
- Following surgery, patients should be instructed to contact their physicians if any signs or symptoms of infection occur.

- Store at a temperature below 25°C and away from light.
- Do not use the product if the powder has a yellowish or brownish colour or if the liquid is syrupy.
- These two conditions indicate that the product has not been stored correctly.
- The users must have specific training and experience to be thoroughly familiar with the properties, handling characteristics and application of the product and percutaneous cement delivery. It is strongly recommended that the surgical team carry out practical trials prior to use in patients under the same environmental conditions using the same instruments for clinical use.
- Adverse patient reactions affecting the cardiovascular system have been associated with acrylic resins (bone cements); hypotensive reactions have occurred between 10 and 165 seconds. They have lasted monitored carefully for any change in blood pressure during and immediately following the application of bone cement.
- The liquid component is flammable and volatile: the operating theatre must be correctly ventilated. The liquid component and/or its vapours must never be directly exposed to naked flames or heated materials.
- MMA can cause hypersensitivity reactions in susceptible persons, which may result in an anaphylactic response.
- The use of this product is not recommended in patients that do not exhibit a pathological condition, such as primary or secondary osteoporosis or a tumor, which would impair the ability of the patient to heal using conservative treatment methods.
- Adhere strictly to good surgical principles and techniques. Deep wound infection is a serious postoperative complication and may require total removal of the embedded cement. Deep wound infection may be latent and not manifest itself even for several years postoperatively.
- Follow handling, mixing and preparation instructions carefully.
- Do not resterilize. Single patient use only. Sterile only if package is unopened and undamaged.
- Caution should be used to prevent excessive exposure to the concentrated vapors of the monomer, which may produce irritation of the respiratory tract, eyes and possibly liver damage.
- Do not use monomer if the liquid component shows any sign of thickening or premature polymerization.
- Do not allow the liquid component to contact rubber or latex gloves. The liquid component is a powerful liquid solvent. Should contact occur, the gloves may dissolve and tissue damage may occur. Wearing a second pair of gloves may diminish the possibility of hypersensitivity reactions.
- Do not allow personnel wearing contact lenses to be near or involved in mixing of the bone cement.
- Use appropriate imaging techniques to verify correct needle placement, absence of damage to surrounding structures, and appropriate location of the injected bone cement. Use imaging, such as fluoroscopy, to assess the ability of the vertebra to contain the injected bone cement.
- Cement leakage may cause tissue damage, nerve or circulatory problems, and other serious adverse events. Avoid over pressurization of the bone cement because this may lead to leakage of the bone cement beyond the site of its intended application. Cement leakage may cause tissue damage, and nerve or circulatory problems.
- Leaks can also occur when injecting if the needle is in a vein or if unseen microfractures are present.
- If bone cement is seen outside of the vertebral body or in the circulatory system during the procedure, immediately stop the injection.
- Inadequate fixation or unanticipated postoperative events may affect the cement-bone interface and lead to micro-motion of the cement against bone surface. A fibrous tissue layer may develop between the cement and the bone. Long term follow-up is advised for all patients on a regular scheduled basis.
- The completion of polymerization occurs in the patient and is an exothermic reaction with considerable liberation of heat. According to ISO 5833 standard, the temperature can be as high as 90°C ± 5°C. The heat produced in situ has a long term effect and positioning of the patient should be maintained securely to obtain proper fixation. For proper fixation, 1-2 hours or longer may be required, as determined by the patient's medical condition and the attending physician.
- The long term effects of the bone cement in the spine have not yet been established.

PRECAUTIONS

- The product should not be used after the expiration date
- Kit Mendec Spine and Delivery System must only be used with needles specific for vertebroplasty. It is advised to use the Tecres Spine Needles. If needles of other manufacturers are used, confirm the presence of Luer-Lock type connections and be sure that they can be tightly assembled to the safety locknut present on the Extension tube of delivery system of the kit.
- Avoid monomer contact with the skin and mucous membranes. Cases of contact dermatitis have been observed in susceptible subjects. It is therefore advisable to wear a second pair of surgical gloves and conspicuously observe the instructions for mixing the components in order to reduce the possibility of reactions caused by hypersensitivity.
- Temperature has a major effect on the preparation characteristics of any acrylic resin.
- Temperatures of more than 23°C for the product, the preparation accessories and the environment accelerate the various stages in the preparation procedure. Lower temperatures retard the preparation stages. Before using Kit Mendec Spine and Delivery System it is strongly advised to make sure that the package has been stored at a temperature of 23°C ± 1°C for the previous 24 hours.
- Acrylic resins reach temperatures higher than physiological temperatures during the exothermic polymerization reaction. The released heat may damage bone or tissue surrounding the implant.
- Extrusion of the resin beyond the region of its intended application may occur.
- The polymer component should be disposed of in an authorized waste facility. The liquid component can be evaporated under a ventilated hood or absorbed by an inert material and transferred in a suitable container for disposal.

MENDEC[®] Spine Kit

INSTRUCTIONS FOR USE

For the attention of the surgeon

CAUTION

Federal law restricts this device to sale by or on the order of a physician. Users should be trained in use of percutaneous vertebral body injection procedures using bone cement for the treatment of vertebral pathologies and in the management of complications resulting from procedures. Read this document carefully. Failure to follow instructions may cause malfunction and/or cause patient injury.

DESCRIPTION

The Kit Mendec Spine and delivery system consists of the acrylic resin Mendec Spine and of a delivery system. Mendec Spine is a highly radiopaque acrylic resin of medium viscosity. The delivery system consists of an injection system and it may include an injection needle (spine needle) depending on the variant of the kit. It allows acrylic resin to flow into the pathological vertebral body. The injection system consists of Spine Gun and Extension tube. It is to be used only with needles specific for vertebroplasty with Luer-lock type connections.

The kit is sterile and single-use.

The injection system, the powder container and the liquid ampoule are contained in a hermetically sealed blister pack sterilized with ethylene oxide. The liquid is sterilized by means of filtration and the powder by means of ethylene oxide treatment. The needle is sterilized with ethylene oxide.

Composition of MENDEC SPINE:

POWDER	20 g container of powder	67.5% w/w
	Polymethylmethacrylate	30.0% w/w
	Barium sulphate	2.5% w/w
	Benzoyl peroxide	99.1% w/w
LIQUID	9.4 g ampoule of liquid	0.9% w/w
	Methylmethacrylate	75 ppm
	N,N-dimethyl-p-toluidine	
	Hydroquinone	

INDICATIONS

The Mendec Spine and delivery system kit is indicated for the treatment of pathological fractures of the vertebrae using a vertebroplasty procedure. Painful vertebral compression fractures may result from osteoporosis, benign lesions (hemangioma), and malignant lesions (metastatic cancers, myeloma).

CONTRAINDICATIONS

- Active or incompletely treated infections at the site where the resin is to be applied
- Hemorrhagic diathesis
- Extended vertebral destruction, with collapse superior to 2/3 of the standard thickness, destruction of the posterior wall with epidural extension of the pathologic tissue and clinical signs of medullar compression
- Coagulation disorders or severe cardiopulmonary disease
- Spinal stenosis (>20% by retro-pulsed fragments)
- Compromise of the vertebral fractures due to posterior involvement
- Patient clearly improving on medical therapy
- Prophylaxis in metastatic or osteoporotic patients with no evidence of acute fracture
- Non-pathological, acute traumatic fractures of the vertebra
- Vertebra plana (collapse>90%)
- Compromise of the vertebral body or the walls of the pedicles
- Unstable vertebral fractures due to posterior involvement
- In patients previously sensitized to the monomer or any of the other components of the resin

USE DURING PREGNANCY, BREAST-FEEDING AND CHILDREN

The safety and effectiveness of the device has not been studied in pregnant women, breast-feeding women or children.

WARNINGS

- Read and understand instructions before using this product
- This product is intended for use only by clinicians with adequate training and experience in the percutaneous vertebral body injection procedures. The physician is advised to consult the medical literature regarding techniques, complications and hazards associated with the intended procedures.
- Resin as a manufacturer does not recommend a surgical technique. It is the responsibility of the physician to determine the appropriateness of this product and specific technique for each patient

TECRES S.p.A.
 Via A. Doria, 6 - 37066 Sommacampagna - Verona - ITALY
 Telefono +39.045.9217311 - Fax +39.045.9217330
 info@tecres.it - www.tecres.it

for US market

HIGH BARIUM
 CONTENT

MENDEC® Spine Kit