CERTIFIED COMPANY UNI EN ISO 9001:2008 & UNI EN ISO 13485:2012

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SPESORE 40 MY

POLYPROPYLENE (PP) BAGS FOR AUTOCLAVE STERILIZATION, RESISTANT UP TO 134°C. BAGS SINGLE FOLD TYPE.

THESE BAGS ARE PARTICULARLY USEFUL IN ORDER TO REMOVE, BY AUTOCLAVE STERILIZATION PROCESS, PATHOGENIC SUBSTANCES NORMALLY PRESENT IN CONTAMINATED CONTAINERS (PIPETTES, TEST TUBES, MICROTITER PLATES, ECC.).

“BIOHAZARD” SYMBOL AND OPERATING INSTRUCTIONS DIRECTLY PRINTED ON THE BAG.

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“BIOHAZARD” SYMBOL AND OPERATING INSTRUCTIONS DIRECTLY PRINTED ON THE BAG. WRITTEN IN YELLOW COLOUR CONFORM TO UNI 7545/1-10 (PANTONE 109C / RAL 1023) PRINTED ON TRANSPARENT BAG.

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Intended purpose: "GENERAL LABORATORY USE". Product suitable for autoclave sterilization. PRODUCT NOT SUBJECT TO CE MARKING.

For use in professional test laboratory only

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- > Containers filled with liquid, must not be sealed or tapped;
- > Do not introduce sharp objects, such as broken glassware, into the autoclave bags;
- > Add a little high of water into the bags containing solid residues. The water will vaporize channeling outside the residual air once reached the sterilization temperature internally to the bags;
- > Do not close the bag as this will avoid to the air to move out during the sterilization process;
- > Do not overload the autoclave and leave enough space for the steam circulation;
- > For the decontamination and deactivation of biological waste, particularly resistant, to do autoclave process at +134 °C.

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Keep out of flame or heat sources which might damage the product

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Do not use after expiry date or if packing is opened

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Do not re-use: Disposable device

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Do not vary the intended purpose of the product

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Store in dry place, Temperature range: min -10°C max +50°C

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Disposal: use appropriate personal protective equipment and act according to applicable regulations

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Before use with particular substances, check the resistance / compatibility chart on our catalogue

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