

User Manual for Sterilization Roll

Product Name: Heat Sealing Sterilization Flat Roll

Specification: 50mmx200m, 55mmx200m, 75mmx200m, 100mmx200m, 150mmx200m, 200mmx200m, 250mmx200m, 300mmx200m, 350mmx200m, 400mmx200m

Intended Use: it is widely used in sterilization and preservation of medical instruments in hospital supply room, operating room and stomatology department, etc.

Structure: Composed of a layer of medical dialysis paper resistant to high temperature and high-pressure steam and a layer of high temperature resistant transparent composite film. On the paper, it is with EO and steam sterilization indicator.

Instruction for Use:

1. Cut the appropriate size according to the medical device.
2. Put the clean and dry medical instrument into the sterilization bag.
3. Seal the bag with heat sealing machine.
4. Put the sterilization bag into STEAM sterilization machine or ETO sterilization machine, Sterilize according to conventional methods.
5. Mark the sterilization time and other information on the bag, confirm that the sterilization indicator has changed color.
6. When needed, tear open the sterilization bag and take it out for use.
7. DO NOT USE A SECOND TIME.

Suitable for Sterilization

1. The sterilization condition of the lower row pressure steam sterilizer by the sterilization bag is 102.9 kPa, the temperature reaches 121 °C, the sterilization time is 30 minutes, and biological monitoring is carried out regularly.
2. The sterilization condition of the sterilization bag to the pre vacuum pressure steam sterilizer is 205.8kpa, the temperature reaches 134 °C, the sterilization time is 4 minutes, and biological monitoring is carried out regularly.
3. The sterilization bag is also suitable for ETO ethylene oxide sterilization. Or packaging materials with other sterilization methods can be selected.
4. The sterilization bag is not suitable for dry heat sterilization.

Use Suggestions

1. After standard sterilization, the sterilized instruments and articles stored in the sterilization bag are stored in a qualified environment. When the temperature, humidity and cleanliness are appropriate, it is recommended to use them within 12 months after steam or ETO ethylene oxide sterilization.
2. Before use, check the integrity of the sterilization bag. Do not use it if it is damaged.
3. The shelf life of the sterilization bag is 5 years.