
GE Healthcare

NM 830

Data Sheet

DOC2121937 Rev. 3

SECRET DE AFACERI
CONFIDENTIAL

Overview

NM830 is a premium, all-purpose, dual detector free-geometry integrated nuclear imaging system, featuring an advanced, all-digital Elite NXT detector technology, slim gantry, cantilevered patient table, acquisition station, and the Xeleris® processing & review workstation.

Elite NXT detectors feature either 3/8" or 5/8" crystal for all-purpose nuclear imaging.

Key features of NM830 design include:

- Slim-profile, wide-bore, robotic gantry design
- 180° and 90° orientations of the detectors for high SPECT and WB scanning efficiency
- Rapid simultaneous multi-axis gantry motions
- Upright and horizontal detector orientations for exceptional clinical versatility, including patients in hospital bed, standing or sitting during scan
- Multi-functional, dual-axis imaging table
- Automatic "home" positioning enables easy setup of the gantry and the table using pre-programmed detectors geometries and imaging modes
- Real-time automatic body contouring
- User friendly, intuitive Linux-based user interface
- Ignite one-click, integrated workflow with Xeleris processing and review workstation designed to help enhance departmental productivity

Primary Benefits

- Excellent image quality based on advanced Elite NXT detector technology and SPECT optimized design
- NM dose reduction without compromising image quality, enabled by Evolution™ family and SwiftScan option
- Exceptional productivity enabled through half-time SPECT scans, fast and flexible robotic gantry motions for maximal clinical versatility and Ignite streamlined workflow
- Advanced integrated Xeleris workstation clinical applications with remote PC and PACS processing and review solutions

System Components

Elite NXT detectors

NM830 image quality is based on two Elite NXT slim, large rectangular field-of-view digital detectors, featuring five corrections performed on each detected event in real time even at high count rates:

- Uniformity
- Linearity
- Energy
- Isotope decay
- Center of Rotation (COR)

3/8" NM Detector

- 59 circular PMTs - 53 x 3" (76 mm) and 6 x 1.5" (38 mm)
- Crystal thickness: 3/8" (9.5 mm)
- One ADC per PMT, 30.0 MHz sampling rate
- UFOV: 21.25" x 15.75" (54 x 40 cm)
- Energy range: 40 - 620 keV

5/8" NM Detector

- 59 circular PMTs - 53 x 3" (76 mm) and 6 x 1.5" (38 mm)
- Crystal thickness: 5/8" (15.9 mm)
- One ADC per PMT, 30.0 MHz sampling rate
- UFOV: 21.25" x 15.75" (54 cm x 40 cm)
- Energy range: 40 - 620 keV

Gantry

NM830 features a wide 70 cm bore and slim gantry, enabling SPECT, WB and planar imaging in various geometries for a wide patient population.

Key features include:

- Automated detector radial motion (in/out), rotation around the ring, transitions between 180° and 90° geometries or other orientations.
- Flexible design enables a variety of scanning orientations including upright seated or standing patients and imaging patients on stretchers.
- The stationary gantry is secured to the floor, for tomographic center-of-rotation precision.
- Camera setup is performed interactively by the remote control handset and via user-definable pre-programmed acquisition-specific "home" positions. A gantry display unit displays status of the gantry's moving parts and the patient table.
- Real-time, infrared-based Automatic Body Contouring (ABC) designed to enhance scanning efficiency and resolution in 90° & 180° SPECT, and whole body scanning procedures and help minimize patient-detector distance for maximal image quality.

Patient table

NM830 features dual-axis cantilevered table, used for planar, whole body, SPECT and other scanning procedures performed in horizontal patient orientation.

Key features include:

- Maximum patient weight - 500 lb. (227 kg)
- Whole Body scan range 78.7" (200 cm) or 88.6" (225 cm) with optional table extender
- Minimum table height 23.2" (59 cm) - facilitates patient loading and unloading from wheelchair or stretcher
- A dual-axis cantilevered table is used for planar, whole body, and SPECT applications.
- Low attenuation carbon fiber tabletop includes mattress pad/straps for maximum patient comfort.
- Automated positioning helps simplify setup
- Optional interactive bedside touch-ruler for easy setup of scan range
- Manual emergency patient egress

SECRET DE AFACERI
CONFIDENTIAL

- Mobile design enables easy swiveling of the table away from the gantry around a pivot point at the rear, facilitating dual-collimator exchange and imaging and scanning of patients who are seated or on hospital bed / stretcher
- Firm anchoring to the floor with locking pins and floor plates. Free access from both sides for patient loading / unloading with IV, EKG or other devices.
- Optional EKG trigger integrated in the patient-table-integrated to help enable easier patient setup
- Optional accessories include a head holder, table extender, arm support, leg support, infant scan support and table pads/straps

Acquisition station

The NM830 acquisition station employs a Graphic User Interface for exam scheduling, scan acquisition, and utilities for protocol editing, routine quality control and analysis and networking, including:

- Universal connectivity via DICOM 3.0 (as per DICOM conformance statement) and Interfile 3.3 TCP/IP based protocols
- HIS-integrated workflow including DICOM Worklist

Display

- 1280 x 1024 true-color display
- Threshold and windowing control in multiple window settings
- Cinematic display of dynamic and all multi-frame datasets
- Online, live display of acquired data and imaging parameters

Main Features

Acquisition

The NM830 acquisition station is based on Linux operating system with an icon-based graphical user interface, common with the Xeleris workstation.

NM Data acquisition may be performed in any of the following imaging modes: Static, Dynamic, Multi-Gated, Whole Body Scanning, SPECT, and Gated SPECT.

Operation of the NM830 is done with the interactive, graphical GE Common User Interface (CUI). The CUI is designed to help maximize flexibility and productivity and it includes the following features:

- Pre-defined or user-configurable protocols for rapid recall and setup
- On-the-fly digital linearity, energy and uniformity corrections
- Energy spectrum histogram (PHA) display with up to 64 independent windows per detector. Ensures acquisition into correct energy window for given isotope(s)
- Simultaneous acquisition of up to 64 independent windows (e.g.: four frames with 16 independent windows per frame, or sixteen frames with 4 independent windows per frame)
- Acquisition termination by preset time, preset count or manual stop
- Pan/zoom and rotate modes
- Ability to resume paused acquisitions for whole body, SPECT, and gated SPECT
- Ability to create NM radiopharmaceuticals dose structured report according to Dicom standard.

Utilities

Acquisition software includes control of camera maintenance activities including:

- Disk space management

- Pulse Height Analysis (PHA)
- Center of Rotation (COR)
- Single isotope uniformity correction map for all isotopes
- Single isotope Energy, sensitivity, and linearity maps for all isotopes
- Daily/periodic QC including gantry calibrations
- Customizable system parameters
- Definition and setup of acquisition sequences
- Use of preset acquisition protocols

NM830 is capable of acquiring single or multi-isotope scans in the following acquisition modes:

Static Acquisition

- Single frame image (word mode), Frame sizes: 32², 64², 128², 256², 512², 1024²
- Optional bi-lateral detector motion with Pinhole collimator, motion range: 13.8" (35 cm)

Dynamic Acquisition

- Frame sizes: 32², 64², 128², 256²
- Maximum frame rate: 50 frames/sec

Whole Body Scanning

- Infrared-based real-time Automatic Body Contouring (ABC)
- Sequential multi-spot ("step-and-shoot") or continuous scanning modes
- Maximum scan range 78.7" (200 cm) or up to 88.6" (225 cm) with optional table extender
- Minimum WB continuous scan speed: 1.97"/min (5 cm/min)
- Maximum WB continuous scan speed: 11.8"/min (30 cm/min)
- Simultaneous anterior/posterior dual view scan
- Image matrix: 256 x 1024

SwiftScan Planar

Static and Whole-Body scanning using "Clarity 2D" post processing bundled with LEHRS collimator for improved small lesion detectability² and up to 25% time or dose reduction³

Multi-Gated Acquisition

- Equi-time and equi-phase gating modes with real-time irregular beat rejection
- ECG display during acquisition
- On-line R-to-R histogram display
- Multi-gated acquisition capacity:
- 8, 16, 24, or 32 frames per cycle

SPECT Acquisition

- Infrared-based on-line Automatic Body Contouring system
- Step-and-shoot mode with: (2°, 3°, 5°, 6°, 90, 10°, 12°, 15°, 18°, 30°, 36°, 45°, 60° and 90° angular sampling
- Continuous mode with 30, 40, 50, or 60 angular sampling
- Maximum Scan range – 78.7" 200 cm + non-continuous 40 cm (15.7") with extender
- Consecutive multiple SPECT scans combined with table translation enable Volumetrix pasting of up to 7 FOV's (covering most body lengths)
- Either 90° ("L-mode") or 180° ("H-mode") detector geometries
- Matrix size: 64 x 64, 128 x 128, 256 x 256

Gated SPECT acquisition

- Gated tomography with real-time irregular beat rejection

SECRET DE AFACERI

- User-modified R-R acceptance window
- An extra bin sums all data, both accepted and rejected, securing non-gated study completion in case of exceptionally irregular heart-rate.
- Number of frames per R-R interval may vary between 4,8,12,16,24

Dynamic SPECT acquisition (option)

- SPECT scan time:
 - Minimum - 1 min/3600
 - Maximum - 15 min/1800
- Up to 40 consecutive SPECT scans in H-mode
- Continuous mode acquisition, angular sampling: 30, 40, 50, 60
- Supports Pan and zoom
- Supports multi-isotope and multi-peak energy settings

SwiftScan SPECT

SPECT and Gated SPECT scanning using the "Step and Shoot Continuous" mode bundled with LEHRS collimator for improved small lesion detectability⁴ and up to 25% time or dose reduction⁵

Ignite

The Ignite feature for NM830 and Xeleris is designed to help streamline clinician workflow via three simple steps:

- Select the patient in workload
- Set up patient and utilize auto-home positioning
- Click once to begin acquisition of scan and automatic transfer and processing of results on Xeleris.

NM data processing and review

The NM830 acquisition station allows networking to local area networks.

Data acquired on NM830 is transferred under Ignite workflow, a one-click productivity solution, to Xeleris processing & review workstations via DICOM 3.0 standard for processing, archiving, reporting and hard copying

SPECT processing & review

Xeleris offers Volumetrix[®] MI, a hybrid imaging tomography review package to superimpose and display matching pairs of nuclear and CT tomograms.

Volumetrix MI is an interactive display allowing users to navigate in 3D throughout the entire tomographic study, including comparative display of corrected and uncorrected images.

Options

High Precision Collimators

NM830 features a choice of optional collimators which cover a comprehensive selection of applications in a wide range of energies including fan beam and pinhole collimators. For a detailed list of collimators refer to the collimators table within this datasheet.

Pinhole & LEHR Collimator Cart

An optional collimator cart is available, including:

- Pinhole collimator (for Detector 1)
- LEHR collimator (for Detector 2)

Enhances workflow by enabling successive Thyroid / Parathyroid acquisitions without the need to change collimators between scans.

Volumetrix Suite

Optional software for Volumetrix MI that expands the review capabilities and offers extended display and reporting flexibility.

Volumetrix MI

A tomographic data viewing and processing package - a consolidated application for SPECT and PET data processing with or without anatomical data (CT or MR), capable of Multi-SPECT sessions pasting, follow up on multiple studies and input of anatomical oblique slices.

Its Integrated SPECT reconstruction and correction algorithm includes interactive selection of optimal image reconstruction parameters.

- Inputs anatomical oblique slices and features on-the-fly oblique creation/manipulation.
- Facilitates dual isotope and study follow-up reading and accepting multiple SPECT, PET, CT and MR Datasets.
- Supports standard uptake value (SUV) calculation and display mode for PET (LBM, BSA, & BW methods).
- Flexible and user customized layouts.
- Auto Condense mode to automatically condense slices to fit on one screen for printing or filming.
- Multiple fusion modes, including weighted fusion, which allows the user to adjust the contribution of the NM and CT images to the fused results.
- Both the anatomical and the functional images retain resolution
- Fast CT scroll

Volumetrix 3D

Optional software which integrates 3D Fusion into nuclear medicine workflow, bringing together the benefits of advanced 3D visualization with the productivity of traditional 2D image analysis.

- 3D VR of fused SPECT-CT, PET-CT or MR images
- NM segmentation tools to include or exclude portions of either volume in the 3D rendered images, including removal of the table from the CT image and threshold based semi-automatic segmentation.
- Clip & Cut Planes to integrate traditional Axial, Sagittal and Coronal slices simultaneously in to the 3D rendered objects

Volumetrix IR

Optional software providing a choice of DICOM 3.0 CT data sets to be registered to NM or PET data sets in the NM workflow:

- Rigid registration of NM to CT or MR
- Manual adjustment via shift/rotate
- Full or Regional registration (fusion & collect supported)
- Auto accept or user intervention are allowed

MDC - Motion Detection and Correction (option)

- Automated SPECT motion correction

SECRET DE AFACER!

CONFIDENTIAL

- Sinogram and linogram images for QC analysis
- Tools for manual adjustment and correction

SPECT reconstruction

- Adjustment of any of the reconstructions settings
- On the fly image correction setup including resolution recovery (option), attenuation and scatter.
- Filter selection with filter gallery review to determine the most effective filter and filter parameters. 2D Hanning, Metz, Butterworth, Hamming, and Wiener filters are available for FBP reconstruction. 3D Butterworth, Hanning and Gauss filters are available for OSEM reconstruction

Multi-FOV Pasting

- Automatic pasting (user overridden) up to 7 FOVs (covering the whole body)
- Preview pasted volume on SPECT/CT (NM & fused available - up to two FOVs)

Registration

- Manual adjustment via shift/rotate
- Single landmark registration
- Regional registration
- Fusion and spots collection are supported on registered images
- Auto accept or user intervention are allowed at any one of the above packages

Evolution for Bone¹

Optional iterative reconstruction with resolution recovery well-suited for Bone SPECT studies. The Evolution for Bone algorithm, which was developed at John Hopkins University & UNC Chapel Hill, models the collimator-detector response, improves Bone SPECT resolution, signal to noise ratios and reduces noise variability.

Use of the option enables:

1. Improved resolution of Bone SPECT studies acquired over standard acquisition time; or
2. Non-inferior image quality with up to 50% Reduction in count density, enabling:
 - a. Imaging at ½ acquisition time compared to standard Bone SPECT protocols; or
 - b. Imaging at ½ the dose compared to standard Bone SPECT imaging protocols; or
 - c. Various combinations of (a) and (b).

Evolution toolkit¹

Optional package enabling improved resolution and reduced noise for SPECT studies of Tc99m, I123, In111, Ga67, I131 by the use of the Evolution reconstruction technique with resolution-recovery. Compared to standard FBP or iterative reconstruction, Evolution Toolkit can enable improved visual clarity.

Evolution for Bone Planar¹

Optional add-on feature to Whole Body and Spots Bone Review application. Evolution for Planar Bone is an Adaptive Structure Matching Non-Local Filter which preserves fine structures in the image. It employs effective noise-reduction algorithm, well-suited pixel size and optimal energy window settings. Evolution for Planar Bone enables up to 50% reduction in count density compared to

standard protocol, delivering any of the following:

1. Improved WB planar image quality*** for the same scan time
2. Shorter WB planar scan time while preserving image quality
3. Reduced injected dose of the same scan time while preserving image quality

Note: Evolution for Bone Planar should not be used in conjunction with SwiftScan Planar.

Evolution for Cardiac¹

An optional OSEM resolution recovery reconstruction algorithm which models the collimator-detector response developed at John Hopkins University & UNC Chapel Hill. It improves cardiac SPECT resolution, signal to noise ratios, and noise variability. The package provides non-inferior image quality with up to 50% reduction in count density of Cardiac SPECT studies, enabling:

1. Imaging at ½ acquisition time compared to standard MPI protocols, or
2. Imaging at ½ the dose compared to standard MPI protocols, Or
3. Various combinations of 1 and 2.

Cardiac Morphing

Optional software providing the following performance:

- Delivers enhanced perfusion image quality for Gated SPECT data
- Overcomes non-uniform blurring of the myocardium due to cardiac motion, enhancing the visual clarity of the images
- Based on elastic summation of tomograms representing all bins in a cardiac cycle - rescaled to match end-diastolic LV size

CT attenuation correction (CTAC)

Optional feature for Myovation and Xeleris application that enables the use of imported CT AC scan of same patient for accurate attenuation correction of multiple SPECT scans, assuming patient position, orientation, (including all movable organs) and body size are similar, yielding:

1. Help enable enhanced image quality, compared to non-corrected image.
2. Help enable elimination of the need for acquisition of separate CT scan for attenuation correction of each SPECT study, leading to simplification and potentially higher productivity
3. Help enable reduced exposure to CT radiation by eliminating the need for separate CT scan for attenuation correction of each SPECT study
4. Note: this feature was validated using data from GE Healthcare NM and CT scanners.

SECRET DE AFACERI
CONFIDENTIAL

Other Options

- Carbon fiber axial head holder for brain SPECT scans
- Internal NORAV ECG trigger
- External IVY R-Wave Trigger, with & without ECG chart recorder, with or without stand
- Table Leg Extender
- Interactive bedside touch-ruler
- QA bar phantom
- Rectangular Co-57 flood (Site license required)
- Butterfly arm support – optimized for cardiac SPECT studies
- Leg support
- Infant Scan Positioner
- Uninterruptible Power Supply (UPS)
- Flat Floor Plate
- Mobile PC cart for acquisition station

Siting Requirements

Recommended Suite Size

Exam Room Minimum Size
16'9" x 12'3" (5.12 m x 3.74 m)

These values meet minimum clearance requirements under U.S. Federal Regulations and National Standards: 29 CFR 1910 (OSHA), NFPA 70E (Standard for electrical Safety in the Workplace), and NFPA 101 (Life Safety Code). Specific room requirements may also need to comply with local and regulatory requirements.

Power Requirements

Power conditioning is incorporated into the primary power supply of the system. The system can operate on line voltage based on local conditions and codes.

208-240 VAC nominal, 1 PH + Neutral 50/60Hz.

Power consumption is 6 kVA, 2.5 kVA average.

Operating Conditions

Parameter	Maximum	Minimum	Maximum Change Rate
Temperature	26° C (79° F)	18° C (64° F)	3° C/hr (5° F/hr)
Humidity	60% non-condensing relative	30% non-condensing relative	5% / hr

Regulatory Compliance

This product is designed to comply with applicable safety standards.

The NM830 system complies with the requirements of IEC-60601-1-2 and relevant collaterals, UL 60601-1, Safety Standards for Medical Electrical Equipment

This product is a CE-compliant device that satisfies regulations regarding Electro-Magnetic Compatibility (EMC) and Electro-Magnetic Interference (EMI), pursuant to IEC-60601.

GE Healthcare has been certified to be ISO-13485 compliant.

Warranty

Please contact your sales representative for information.

Mechanical specifications

Table	Specification
Table Weight	860 lb. (390 kg)
Maximum patient load capacity	500 lb (227 kg)
Maximum NM WB Scan Length	78.7" (200 cm) + 9.9" (25 cm) with table extender
Maximum SPECT scan range	78.7" (200 cm) + non continuous 15.7" (40 cm) with table extender
Table Width	20.9" (53 cm)
Tabletop width	15.7" (40 cm)
Table length	9'2" (280 cm)
Height	Minimal: 23.2" (59 cm) Maximal: 38.6" (98 cm)
Attenuation	<10% for 140 keV gamma rays
Maximum pallet Height at center	34.1" (86.7 cm)
Vertical travel time (full range)	Slow: 34 sec Fast: 26 sec
Horizontal Speed (Manual, handheld controller)	Slow: 0.98"/s (24.8 mm/s) Fast: 3.94 "/s (100 mm/s)

SECRET DE AFACERI

CONFIDENTIAL

Gantry	Specification
Depth (length)	53.5" (136 cm)
Width	60.2" (153 cm)
Height	83.1" (211 cm)
Gantry bore size (Diameter)	27.6" (70.0 cm)
Weight (with LEHR collimators)	4,523 lb. (2,052 kg)
Rotational (Axial) Motion Speed	0.0333 to 3.0 rpm automatically; 1.0 or 3.0 rpm manually
Independent Radial Motion Speed	19.7" or 29.5"/min (50 or 75 cm/min)
Emergency Stops	3
Rotation	540°
NM Gantry Orientations SPECT Planar	90° and 180° opposing 90° and 180° opposing, Horizontal 0° Planar, Vertical 180° Planar
Swivel Tilt Range	0° to 180° about detector axis
Caudal/Cephalic Tilt Range	Each detector can swivel independently $\pm 45^\circ$ (available only on hospital bed/gurney)
Detector Scan Radius (With LEHR Collimators)	Minimum: 3.9" (10 cm) radius/7.8" (20 cm) diameter Maximum: 13.8" (35 cm) radius/27.6" (70 cm) diameter
Lateral Motion Speed	19.7" or 29.5"/min (50 or 75 cm/min)

Main programmed orientations of NM830 gantry

NM830 features summary

SUBSYSTEM	KEY FEATURE	USER/PATIENT BENEFIT
Elite NXT DETECTORS	2 Elite NXT rectangular digital detectors with real-time corrections for sensitivity, linearity, energy, isotope decay, and COR	High resolution large FOV detectors with excellent image quality, and stability
	Slim detectors	Low gantry load for excellent mechanical precision and wide detectors clearance for wide referral population
	Real-time Automatic Body Contouring	Automatically follows the contour of the patient for both SPECT and whole body imaging, maintaining minimal patient to detector distance for high resolution images; minimizes time to position patients
	3/8" or 5/8" crystal - Shielded for 40 - 620 keV range	Optimized for wide energy range and multiple isotope studies
COLLIMATION	Collision sensitive pads on collimators surface and on body contouring device	Protect against collision with patient body; automatically halt detector/table motions for patient safety
	High Precision Collimators	Maintain excellent image quality in all applications
	Collimator storage/exchange device combined	Productivity, ease of use and safety; both collimators removed/ replaced during one exchange session

SUBSYSTEM	KEY FEATURE	USER/PATIENT BENEFIT
DUAL AXIS IMAGING TABLE	Single table for all study types	Ease of use and productivity for both whole body and SPECT procedures without changing table positions
	Patient weight load of up to 227 kg (500 lb.)	Accommodates greater portion of general population compared with previous GE NM camera models
	Tabletop longitudinal travel accommodates 6'7" (200 cm) patient	Accommodates greater portion of general population compared with previous GE NM camera models
	Table vertical travel range: Minimal: 23.2" (59 cm) Maximal: 38.6" (98 cm)	Ease of patient transfer from wheelchair or stretcher; well suited for geriatric and pediatric patients
	Mobile design pivots on rear floor pin	Study flexibility for seated and stretcher patients, as well as collimator exchange; reduces tripping hazards
	Low attenuation carbon-fiber tabletop	Patient comfort on wide table can help reduce patient movement and, therefore, improve image quality
	Hand grips on sides and cradle release on rear of table	Ease of table positioning and fast patient egress (cradle release) in case of emergency
	Accommodates head holder	Fixed patient positioning for high quality brain SPECT studies
	Accommodates wide arm support straps	Patient comfort
	Optional bedside touch ruler	Simplify and speed up scan range setup
	Accommodates tabletop extender	Ability to image taller patients in legs-in position and extend WB scan range from standard 200 cm to 225 cm
	Accommodates Infant Scan support	Ability to image infants with improved detector-patient proximity scanning due to narrow accessory pallet.
HANDHELD CONTROLLER	Icon-based design	Ease of use for patient setup
	Pole-mounted handheld controller on top of gantry	User can operate controller from either side of the gantry for flexibility; user can start exam without leaving the patient's side
GANTRY	Slim gantry system fits 16'9" x 12'3" (5.12 m x 3.74 m) room layout	Minimal room layout requirements
	Detectors shape optimized for 90° geometry (valid for LEHR and ELEGP)	Shaped detector edges fit together tightly during 90° cardiac mode, minimizing dead space for high sensitivity and resolution
	Secured position on floor	Gantry does not move on rails, providing for mechanical stability and reliability; reduces tripping hazards
	Automatic study "Home" positions	Ease of use and quick patient setup with factory home positions of: collimator exchange, 180° SPECT, 180° Brain SPECT, 180° Whole Body, Vertical orientation facing out, 90° Cardiac SPECT, 90° SPECT Spine, and Stretcher
	Externally mounted dual detectors	Patient friendly; ease and speed of patient positioning
	Emergency Stop buttons	Stops all system motions upon emergency for patient safety
	Motorized radial detector movement in 90° position	Highly flexible patient positioning for cardiac imaging; COR stability and reproducibility
	Free geometry capability of various 0°, 90°, and 180° configurations	High throughput configurations for cardiac SPECT and whole body imaging, as well as providing flexibility of a single-head camera
	Flexible detectors positioning including upright standing or seated patients and hospital bed scanning	High clinical utility including "single-head positioning flexibility"
	Ultra-fast transition between the various gantry geometries	Technologist productivity; ease and speed of study set-up

SECRET DE AFACERI

CONFIDENTIAL

SUBSYSTEM	KEY FEATURE	USER/PATIENT BENEFIT
ACQUISITION STATION	Graphical user interface with Xeleris workstation	Facilitates ease of use
	Linux operating system	Parallel operations allow simultaneous acquisition and display to help enhance productivity
	Fully-integrated workflow solution with Ignite	Workflow productivity for streamlined imaging procedures including camera /patient setup, acquisition, processing and archiving.
	Networks to Xeleris Processing & Review	Non-proprietary communication protocols; open system
DICOM CONNECTIVITY	Networks to DICOM-compliant systems	Communicates with 3 rd party DICOM-compliant workstations
	DICOM Modality Worklist	Streamlines patient admittance and scheduling

Elite NXT NM Detector NEMA Performance Specifications[†] Summary

SPECIFICATION	PARAMETER		RANGE	DATA (3/8" Crystal)	DATA (5/8" Crystal)
Field-of-View	UFOV	Length	=	54 cm	54 cm
		Width		40 cm	40 cm
Intrinsic Energy Resolution (Tc-99m @ 20 kcps)	UFOV	FWHM	Ⓔ	9.5%	9.5%
Intrinsic Spatial Resolution	CFOV	FWHM	Ⓔ	3.7 mm	4.5 mm
		FWTM	Ⓔ	6.8 mm	8.5 mm
	UFOV	FWHM	Ⓔ	3.8 mm	4.6 mm
		FWTM	Ⓔ	6.9 mm	8.7 mm
Intrinsic Uniformity	CFOV	Differential	Ⓔ	2.1%	2.1%
		Integral	Ⓔ	3.0%	3.0%
	UFOV	Differential	Ⓔ	2.3%	2.3%
		Integral	Ⓔ	3.6%	3.6%
Intrinsic Spatial Linearity	CFOV	Differential	Ⓔ	0.2 mm	0.2 mm
		Absolute	Ⓔ	0.4 mm	0.5 mm
	UFOV	Differential	Ⓔ	0.2 mm	0.2 mm
		Absolute	Ⓔ	0.4 mm	0.5 mm
Multiple Window Spatial Registration			≤	0.5 mm	1.0 mm
Intrinsic Count Rate	Maximum Count Rate		≥	460 kcps	460 kcps
	Maximum @ 20% Window		≥	400 kcps	400 kcps
	20% Loss @ 20% Window		≥	250 kcps	250 kcps
SPECT Reconstructed Spatial Resolution with Scatter (LEHR Collimators)	Central		Ⓔ	9.9 mm	10.3 mm
	Radial		Ⓔ	9.9 mm	10.3 mm
	Tangential		Ⓔ	7.5 mm	7.9 mm

SECRET DE AFACERI

CONFIDENTIAL

Evolution for Bone SPECT performance specifications

SPECIFICATION	PARAMETER	RANGE	Full time	Half time
Evolution for Bone SPECT Reconstructed Spatial Resolution with Scatter ¹ (LEHR Collimators)	Central	≤	6.4 mm	7.0 mm
	Radial	≤	5.7 mm	6.0 mm
	Tangential	≤	5.1 mm	5.4 mm

¹per NEMA NU-1 2007

²reconstructed without post-filtering

NM830 Parallel Hole Collimators

DESCRIPTION	NAME	CATALOG NUMBER (a)	RECOMMENDED ISOTOPE	FIELD OF VIEW (cm) (b)	CALCULATED PENETRATION (%)	SYSTEM SENSITIVITY (cpm/μCi) @100 mm 3/8"/5/8" Per Detector (c)	SYSTEM SENSITIVITY (cps/MBq) @100 mm 3/8"/5/8" Per Detector (c)	SYSTEM RESOLUTION FWHM (mm) @100mm 3/8"/5/8" (d)	TYPE OF HOLE	HOLE DIAMETER (mm)	SEPTAL THICKNESS (mm)	HOLE LENGTH (mm)	WEIGHT (kg/lb) 1 pcs
Low Energy** High Resolution Sensitivity	LEHRS	H3909AD	Tc99m	54 x 40	2.3 (Tc-99m)	206/N.A. (Tc-99m)	92/N.A. (Tc-99m)	7.4/N.A. 6.1 ^{ml}	hex	1.43	0.13	32	43/95
Low Energy** High Sensitivity	LEHS	H2506TV	Tc99m	54 x 40	5 (Tc-99m)	470/N.A. (Tc-99m)	211/N.A. (Tc-99m)	10.8/N.A.	hex	2.31	0.152	34	60/132
Low Energy** Ultra-High Resolution	LEUHR	H2506TH	Tl-201/ Tc99m Studies	54 x 40	0.3 (Tc-99m)	83/N.A. (Tc-99m)	38/N.A. (Tc-99m)	6.1/N.A.	hex	1.22	0.15	38	80/176
Low Energy** High Resolution	LEHR	H2506TB	Tl-201/ Tc99m Studies	54 x 40	0.3 (Tc-99m)	160/165 (Tc-99m)	72/74 (Tc-99m)	7.4/7.7	hex	1.5	0.2	35	60/132
Extended** Low Energy General Purpose	ELEGP	H2506TD	I-123/ Kr-81 Studies	54x 40	0.3 (I-123) 2.3 (Kr-81)	320/330 (Tc-99m) 224/245 (I-123)	144/148 (Tc-99m) 101/110 (I-123)	10.3/10.6	hex	2.5	0.4	40	62/136
Medium Energy General Purpose	MEGP	H2506TC	Ga-67/ In-111 studies	54 x 40	2.0 (Ga-67)	144/150 (Ga-67)	65/67 (Ga-67)	9.4/9.8	hex	3.0	1.05	58	103/227
High Energy General Purpose	HEGP	H2506TE	I-131 studies	54 x 40	2.0 (I-131)	97/165 (I-131)	43/73 (I-131)	12.0/12.5	hex	4.0	1.8	66	131/289

- a. Each commercial item includes 2 collimators mounted on a cart
b. Collimator field-of-view
c. Measured with 20% window using relevant isotope for each collimator, with ±10% tolerance
d. Measured at 100 mm distance from collimator face with ±4% tolerance
e. Measured on 3/8" configuration with Clarity2D 40% blending

**L-mode SPECT (90° detectors geometry) is effective with low energy collimators (LEHS, LEUHR, LEHR, LEHRS and ELEGP)

CONFIDENTIAL

SECRET DE AFACERI

NM830 Pinhole Collimator

DESCRIPTION	NAME	CATALOG NUMBER	RECOMMENDED APPLICATION/ ISOTOPE	FIELD OF VIEW (mm)	Weight (kg/lb)	Insert hole diameter (mm)	SYSTEM SENSITIVITY (cpm/μCi) @100 mm 3/8"/5/8" Per Detector (b)	SYSTEM SENSITIVITY (cps/MBq) @100 mm 3/8"/5/8" Per Detector (b)	SYSTEM RESOLUTION FWHM (mm) @100mm 3/8"/5/8" (c)
General Purpose Pin Hole (3 inserts)	GPPH(a)	H2506TF	Thyroid/Tc99m, I123, I131	200 diameter	98/216 – 1 pcs. 169/370 – 2 pcs.	2 4.45 8	43/47 200/216 570/600	19/20 90/95 258/268	3.8/4.2 6.5/7.1 11.4/12.1

- a. Commercial item includes a single pinhole collimator mounted on a cart or single pinhole collimator with LEHR parallel hole collimator, mounted on the same cart.
b. Sensitivity measured with Co57 (point source at 100 mm from insert center with PSD cover, 20% window) extrapolated to Tc-99m with ±10% tolerance
c. Resolution measured with ±4% tolerance

NM830 Fan Beam Collimator

The optional Fan Beam Collimator enables higher efficiency brain SPECT studies compared to LEHR collimator.

DESCRIPTION	NAME	CATALOG NUMBER (a)	RECOMMENDED APPLICATION/ ISOTOPE (b)	FIELD OF VIEW (mm) (b)	CALCULATED PENETRATION (%)	TYPE OF HOLE	HOLE DIAMETER (mm)	SEPTAL THICKNESS (mm)	HOLE LENGTH (mm)	Volume sensitivity (c) 3/8"/5/8"	SPECT resolution 3/8"/5/8" (d)	Weight (kg/lb) 1 pcs
Fan Beam	Fan Beam	H2506TG	Brain/ Tc99m	540 x >250	0.1 (Tc99m)	hex	1.5	0.2	40	39,960 [cpm/ μ Ci/cm ²] 18,000 [cps/MBq/cm ²]	Central 8.0/8.2 mm Peripheral Radial 8.5/8.7 mm Peripheral Tangential 6.7/7.0 mm	64/140

a. Commercial item includes 2 collimators mounted on a cart

b. Collimator field-of-view


c. Sensitivity measured on two detectors, Tc99m with 15% window, with $\pm 10\%$ tolerance

d. SPECT resolution measured with $\pm 4\%$ tolerance

Notes:

- Recommended brain scan radius is 14-16 cm
- Focal distance from collimator surface is ≥ 350 mm
- Effective collimator thickness is approximately ~ 57 mm

NM830 Minimal Room Layout: 16'9" x 12'3" (5.12 m x 3.74 m)

The sample floor plan of minimum room layout (16'9" x 12'3" / 5.12 m x 3.74 m) meets egress and service clearance requirements under U.S. Federal Regulations and National Standards: 29 CFR 1910 (OSHA), NFPA 70E (Standard for Electrical Safety in the Workspace), and NFPA 101 (Life Safety Code). Specific room layouts may be subject to additional local and regulatory requirements.

For non-U.S. installations, specific installations are subject to country, local and regulatory requirements. For more information, please contact your local GE representative.

Evolution Disclaimer:

- (1) In clinical practice, Evolution options (Evolution for Bone, Evolution for Cardiac, Evolution for Bone Planar) and Evolution Toolkit are recommended for use following consultation of a NM physician, physicist and/or application specialist to determine the appropriate dose or scan time reduction to obtain diagnostic image quality for a particular clinical task, depending on the protocol adopted by the clinical site.
- Evolution claims are supported by simulation of count statistics using default factory protocols and imaging of Tc-99m based radiotracers with LEHR collimator on anthropomorphic phantom or realistic NCAT + SIMSET phantom followed by quantitative and qualitative images comparison.
 - Evolution Toolkit claims are supported by simulation of full count statistics using lesion simulation phantom images based on various radiotracers and collimators and by showing that SPECT image quality reconstructed with Evolution Toolkit provide equivalent clinical information but have better signal-to-noise, contrast, and lesion resolution compare to the images reconstructed with FBP / OSEM.
 - Evolution is supported with LEHRs, LEHR, MEGP, HEGP, ELEGP and LEHS collimators.

SwiftScan Planar Disclaimer:

- (2) As demonstrated in phantom testing using a model observer. Compared to using the LEHR collimator without Clarity 2D.
- (3) Compared to using the LEHR collimator without Planar Clarity 2D. As demonstrated in phantom testing using a bone scan protocol and a model observer. Because model observer results may not always match those from a human reader, the actual time/dose reduction depends on the clinical task, patient size, anatomical location and clinical practice. A radiologist should determine the appropriate scan time/dose for the particular clinical task.

SwiftScan SPECT Disclaimer:

- (4) As demonstrated in phantom testing using a model observer. Compared to using the LEHR Collimator and a SPECT Step and Shoot acquisition.
- (5) Compared to using the LEHR collimator with SPECT Step & Shoot scan mode. As demonstrated in phantom testing using a bone scan protocol, Evolution processing, and a model observer. Because model observer may not always match those from a human reader, the actual time/dose reduction depends on the clinical task, patient size, anatomical location, and clinical practice. A radiologist should determine the appropriate scan time/dose for the particular clinical task.

SECRET DE AFACERI
CONFIDENTIAL

About GE Healthcare

GE Healthcare provides transformational medical technologies and services to meet the demand for increased access, enhanced quality and more affordable healthcare around the world. GE (NYSE: GE) works on things that matter - great people and technologies taking on tough challenges. From medical imaging, software & IT, patient monitoring and diagnostics to drug, biopharmaceutical manufacturing technologies and performance improvement solutions, GE Healthcare helps medical professionals deliver great healthcare to their patients.

GE Healthcare
3000 North Grandview Blvd
Waukesha, WI 53188
U.S.A
www.gehealthcare.com

Imagination at work

©2019 General Electric Company - All rights reserved.

General Electric Company reserves the right to make changes in specification and features shown herein, or discontinue the product described at any time without notice or obligation.

GE, GE Monogram and Imagination at work are trademarks of General Electric Company.

GE Healthcare, a division of General Electric Company.

*Trademark of General Electric Company

**Only LEHS, LEHRS, LEHR, LEUHR, and ELEGP collimators can be used for SPECT at 90° geometry

***See BrightSpeed Elite Product Data Sheet for details

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
GE Healthcare
3000 North Grandview
Blvd Waukesha, WI
53188
U.S.A
www.gehealthcare.com

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GE, GE Monogram and Imagination at work are trademarks of General Electric Company. GE Healthcare, a division of General Electric Company.

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PRODUCT DETAILS

I. NM 830

Cty.	Catalog Ref	Description
1	H3910AR	NM 830
1	H2506TB	LEHR Collimators for D670/D630 (2 with Cart)
1	H3909AX	NM800 NM600 Russian Operating Documents
1	H3909BF	NM800 NM600 Romanian Operating Documents
1	H2506KK	17 Inch LCD Acquisition Display Monitor
1	B77322CA	CT chair (no armrest)
1	H3100PE	D670 QC Point Source HLDR
1	H3100PF	D670 QC Flood Scldr Kit
1	H3100XK	QC Bar Phantom (without CE mark)
1	H3602SL	Point source holder for gamma-camera quality assurance
1	H2506KR	Norav ECG Gating for D630
1	H3100NW	D670 Axial Head Holder
1	H3100PS	NM800/NM600 Touch Ruler
1	H3100NP	Straps and Pad Kit
1	S8390AJ	X4 DR WS SPECT
1	H3901TV	Xeleris DICOM Viewer
1	H3901NJ	X3.1 DVD/CD Archiving Lic
1	H3901RH	X3.1 Cedars Suite 1st N 2
1	H3900PD	Manchester Benetogram Analysis Package
1	H3901TW	Xeleris Wide Monitor
1	H3901TW	Xeleris Wide Monitor
1	H3903DR	Xeleris Second Monitor and Dual Monitor License
1	H3900JZ	Non-US PC & Monitor Power Cables
1	H3904BV	X4 DR Romanian Lang Kit
1	E46011BR	CE PDR CT370-0/tight
1	E45971ND	ABB UPS Power Value 11T B2 Q2 CE UPS, 6 kVA, 230V, 50Hz, 3 min runtime, plug-and-play
1	E85151AA	Rectangular refillable flood phantom for cameras with 40 x 54 cm detectors
1	E8500NB	Patient Arm Support System for Nuclear, PET/CT, MRI, width 48cm, coated, velcro attachment
1	E8500NC	Patient Leg Rest for Nuclear, PET/CT, MRI, coated
1	E6900AE	48In Wide Workstation Table
1	R12023AC	Mandatory A3 SSA Warranty package
1	A82000NM	10 days training education for customers
1	A33331NM	Digital Academy in Nuclear Medicine

SECRET DE AFACER!

CONFIDENTIAL

LONG DESCRIPTIONS

All pictures featured are for illustration purposes only, and do not necessarily relate to products or services mentioned in the text.

H3910AR - NM 830

NM830 is a premium, all-purpose, dual detector free-geometry integrated nuclear imaging system, featuring an advanced, all-digital Elite NXT detector technology, slim gantry, cantilevered patient table, acquisition station, and the Xeleris* processing & review workstation.

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H2506TB - LEHR Collimators for D670/D630 (2 with Cart)

Set of two low energy high resolution collimators with cart for Discovery NM/CT 670 or Discovery NM 630 gamma-camera.

Includes:

- Two LEHR Collimators
- One mobile cart for storage

H2506KK - 17 Inch LCD Acquisition Display Monitor

The 17 inch acquisition display monitor is offered for the operator display of the gamma-camera acquisition station;

It Includes a compact LCD display monitor to be integrated with the acquisition station.

H3100PE - D670 QC Point Source HLDR

An L-shaped metal plate attachable to the wall with an opening for a syringe in order to acquire point source-based flood acquisition at a few meters distance from vertically positioned detector for QA purposes.

H3100PF - D670 QC Flood Srchldr Kit

Quality Control Flood Source Holder Kit

A large plate mounted at a small distance above the NM detector on which the flood source is positioned in order to perform acquisition of flood studies for QA/QC purposes.

SECRET DE AFACER
CONFIDENTIAL

H3100XK - QC Bar Phantom (without CE mark)

Bar phantom for spatial resolution and linearity tests of gamma cameras. The phantom consists of four quadrants with different bar specification.

For each of the quadrant, bar spacing is 2.5mm, 3.2mm, 3.5mm 4.0mm (without CE mark)

H3602SL - Point source holder for gamma-camera quality assurance

Center of rotation source holder for Quality assurance, easily attached to Infinia or Ventri table.

H2506KR - Norav ECG Gating for D630

Norav ECG Gating for D630

Package Contents

The PC-ECG 1200 package contains the following elements:

- Acquisition box
- Patient cable

INDICATIONS FOR USE

Trigger testing intended use: The ECG Trigger provides ECG and R-Wave pulse for display and triggering applications. It is intended for use when TTL level pulse should be triggered during medical tests.

E.g. MRI, CT, Ultra Sound, Nuclear Medicine etc.

H3100NW - D670 Axial Head Holder

The Axial Head holder is ergonomically designed to position patient's head outside of the patient tabletop pallet, enabling brain SPECT orbiting as close as possible to the patient's skull with maximal coverage of the target tissue.

H3100PS - NM800/NM600 Touch Ruler

An interactive touch-sensitive device mounted at one side of the patient table, used to define nuclear imaging scan range (start and stop points), saving the need to enter these values manually from the operator console. This saves necessary manual operation and greatly aids clinical efficiency.

H3100NP - Straps and Pad Kit

A set of tabletop pad and straps used to secure the position and reduce motion of patients scanned on Discovery NM 600 series

SECRET DE AFACE
CONFIDENTIAL

S8390AJ - X4 DR WS SPECT

Xeleris* 4 DR SPECT molecular imaging workstation is a Nuclear Medicine, PET, NM/CT, and PET/CT processing, analysis, and review system. Designed to leverage the latest SPECT quantitative applications for routine clinical use, it accelerates workflow and improves diagnostic confidence. The Xeleris 4 DR opens the doors to the new era of digital healthcare delivery through the enablement of Healthcare Cloud potential and advanced applications to help solve some of the most complex clinical presentations.

Combining streamlined workflow with a comprehensive clinical library and extensive networking capabilities on a molecular imaging workstation, Xeleris 4 DR is at the nucleus of productivity in the clinical imaging department along with enhanced security features. Utilizing the GE Healthcare-wide graphical user interface, Xeleris 4 DR is the processing and review platform of the Discovery*, Optima* and Brivo* NM and NM/CT series, Infinia* Hawkeye* 4, Ventri, Discovery PET/CT 600 series, and all other molecular imaging cameras in GE Healthcare's current offering.

Xeleris 4 DR provides the automated processing and connectivity necessary in today's demanding environment.

Xeleris* 4 DR SPECT includes Motion detection & correction software.

H3901TV - Xeleris DICOM Viewer

The Xeleris DICOM Viewer enhances access to molecular imaging data by providing a DICOM viewing application for general and hybrid Nuclear Medicine studies in a common Xeleris user interface. The application may be distributed on CD, DVD or network destination via Xeleris interface. The Xeleris DICOM Viewer has two components:

Essential Viewer which includes:

- Planar,
- Dynamic,
- DICOM SC,
- DICOM MFSC support.

Basic Workspace review tools including:

- Scroll,
- Pan and Zoom,
- Color map,
- Window level Measurement.

Operates on:

- Windows XP,
- Windows Vista,
- Windows 7,
- Windows 8.

Lumetrix Viewer which includes:

- 3D MIP display and triangulation to 2D images
- Synchronized triangulation between SPECT-CT and PET-CT images for optimum localization
- SPECT-CT and PET-CT fusion with blending tools
- Predefined layouts for optimum SPECT-CT and PET-CT display
- Oblique reformatting

SECRET DE AFAC-F
CONFIDENTIAL

H3901NJ - X3.1 DVD-RAM Archiving Lic

Archiving of patient data on DVD-RAM without change in archived data characteristics, functionality or intended use. The DVD archiving is performed through standard Windows tools, increasing data storage versatility.

H3901RH - X3.1 Cedars Suite 1st N 2

Cedars Sinai Cardiac Packages (option)

A comprehensive set of nuclear cardiology protocols for advanced cardiac analysis, including:

- Cedars Sinai Quantitative Perfusion SPECT(option)
- Automatic 3-Dimensional software approach to quantitative Perfusion SPECT.
- Cedars Sinai Quantitative Gated SPECT(option)
- An application calculating the ejection fraction of the left ventricle and a 3D surface display is generated.
- Cedars Sinai Companion (option)
- Optional module for QGS and QPS applications features
 - 17 segment scores and templates in QPS
 - Diastolic filling parameters in QGS
 - Eccentricity ratio in QGS

H3900PD - Manchester Renogram Analysis Package

Manchester Renogram Analysis Package

The Manchester Renogram Analysis Package is an Aladdin program, which has been developed by the Nuclear Medicine Department at Manchester Royal Infirmary (UK) for analysis of renogram data on GE nuclear medicine workstations.

The program incorporates several features which have been established as clinically useful over the authors many years of experience in renogram analysis at Manchester Royal Infirmary. These include:

- Background subtraction using either the Rutland/Patlak plot, two perirenal regions of interest or a single manually drawn background region.
- Ability to deal with duplex kidneys and solitary kidneys.
- Verification of regions of interest against a dynamic display of the whole renogram study to ensure correct region placement.
- Allowance for attenuation due to kidney depth by a choice of three methods.
- Renogram curves scaled to percentage of administered activity using predetermined sensitivity factors for each camera/collimator used.
- Fixed axis ranges for curve display to aid 'at a glance' visualisation of whether curves represent good or poor function.
- Representative images with regions superimposed are displayed along with curves to show where regions were placed.
- Results of multiple processing of the same data may be saved for redisplay at a later date.
- Separate summed image display to give better visualisation of images alone.
- Ability to quantify activity remaining in kidneys on any late static images.

H3901TW - XELERIS WIDE SCREEN LCD

Xeleris 24" Wide Screen Display

24" wide-screen LCD color monitor

SECRET DE AFACET
CONFIDENTIAL

GE Healthcare

H3903DR - Xeleris Second Monitor and Dual Monitor License

LCD Monitor & License for a single Xeleris Workstation. This item contains: • One 24" WideScreen format monitor for Xeleris MI workstation. • One Dual monitor license Dual Monitor license enables the option on Xeleris Workstation and all its XFL clients

H3900JZ - Xeleris Power Cable

Xeleris Power Cable

E8500NB - Patient arm support for nuclear PET/CT, MRI

Features/Benefits

- Ergonomically designed positioning product developed to prevent arm and shoulder pain and to also reduce related motion artifact that commonly occurs during supine imaging
- This safe, clinically proven Patient Arm Support combines total arm support and passive restraint, thereby increasing patient comfort during extended procedures
- Designed to accommodate virtually all patients, this product not only enhances patient comfort it improves patient throughput

Specifications

- Velcro-mounted; easily cleaned
- Sold per Each

Compatibility

- Compatible with most Nuclear Imaging systems and can also be used in MRI, CT and PET applications

E8500NC - Contoured patient leg rest for nuclear, PET/CT, MRI

Features/Benefits

- The Contoured Leg Rest is a unique support product created to prevent the cycle of low back stress and pain that commonly occurs during supine imaging and treatment
- Positioned under the knees, the Contoured Leg Rest produces mild lumbar spine extension and hip flexion, causing the pelvis to rotate forward creating greater low back support and torso stability. The unique leg recesses create additional lateral support thereby minimizing both rotational and translational body motion

Specifications

- Constructed of full support polyfoam with a seamless coated finish
- Velcro-mounted; easily cleaned
- 18 cm tall from table surface to top of support
- Sold per Each

Compatibility

- Compatible with most Nuclear Imaging systems and can also be used in MRI, CT and PET applications

SECRET DE AFACET
CONFIDENTIAL

E6900AE - 48In Wide Workstation Table

Two monitor workstation table with ergonomic work surface.

Features:

- Work surface that is individually height adjustable
- Cold rolled steel frame with epoxy power coated paint finish
- Particle board core top with high pressure laminate finish
- Tilt, height adjustable, and retractable keyboard tray
- Grommet outlets, CPU cradle, mobile file pedestals, wire management

Double casters (4), 2 locking, fully assembled

R12023AC - Standard Service Package/License

The Standard License provides access to service tools used to perform basic level service on the Equipment and is included at no charge for the warranty period.

A33331NM - Digital Academy in Nuclear Medicine

Our **Level Up** program for a new equipment installation is powered by the **Digital Academy** with an unlimited access to your user community GECARES.com, Live Expert webinars and a unique access to an expert for live assistance. Your team will also benefit from our Learning Management system (LMS), giving them access to training curricula and E-learning.

GE CARES.COM	LIVE EXPERT WEBINARS	E-LEARNING MANAGEMENT*	REMOTE SUPPORT
> CONNECT with your peers and GE experts > SHARE your experience with your peers and publish articles > LEARN new techniques with video tutorials, webinar replays, protocols and clinical cases > CONTACT GE remote experts for all your application needs	> PARTICIPATE to interactive webinars delivered by experts and Key Opinion Leaders on clinical imaging topics > ASK your experience with your peers and publish articles	> INCREASE your technical and clinical knowledge and learnings and download certificates of completion > ASSIGN courses to your team, follow their progress and report activity > RECEIVE learning recommendations based on your needs and profile * Includes 10 licenses and 1 administrator license	> BENEFIT from our screenshare capabilities, directly on your console > RECEIVE remote support during the lifecycle of your equipment, with quick and precise answers to your questions > BOOK an appointment with our online Application Specialists

Target attendees: Radiographers, Technologists, Radiologists, Physicians, Cardiologists

Certification: GE Healthcare Education Certificate

Language: English or local language when available

Access to the Learning Management System includes: 10 licenses and one administrator license

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