Operation Manual

Auto Hematology Analyzer Reagent (M-52D Diluent)

English (3 - 5)
Spanish (6 - 8)
French (9 - 11)
Portuguese (12 - 14)
Russian (15 - 17)
Turkish (18 - 20)

© 2012-2022 Shenzhen Mindray Bio-Medical Electronics Co., Ltd. All rights Reserved.

Intellectual Property Statement

SHENZHEN MINDRAY BIO-MEDICAL ELECTRONICS CO., LTD. (hereinafter called Mindray) owns the intellectual property rights to this Mindray product and this manual. This manual may refer to information protected by copyright or patents and does not convey any license under the patent rights or copyright of Mindray, or of others.

Mindray intends to maintain the contents of this manual as confidential information. Disclosure of the information in this manual in any manner whatsoever without the written permission of Mindray is strictly forbidden.

Release, amendment, reproduction, distribution, rental, adaptation, translation or any other derivative work of this manual in any manner whatsoever without the written permission of Mindray is strictly forbidden.

mindray

MINDRAY

are the trademarks, registered or otherwise, of Mindray in China and other countries. All other trademarks that appear in this manual are used only for informational or editorial purposes. They are the property of their respective owners.

Responsibility on the Manufacturer Party

Contents of this manual are subject to change without prior notice.

All information contained in this manual is believed to be correct. Mindray shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance, or use of this manual.

Mindray is responsible for the effects on safety, reliability and performance of this product, only if:

- 1) Use of this product is conducted as instructed in this operation manual by Mindray authorized personnel or skilled/trained clinical professionals.;
- 2) The product is stored and transported as instructed in the operation manual and the labeling;
- 3) The product is used only on the specified models.

Warranty

Free service will be provided to:

- Products that are covered by the warranty terms of Mindray.

Pay service will be provided to:

- 1) Expired products;
- 2) Products within the warranty period, but are damaged by improper use or human error, stored in improper environment, or damaged by force majeure or other reason not caused by the product itself.

Return Procedure

Follow the procedure below to return products to Mindray if needed.

Contact Mindray customer service department, inform them the product name and lot No. (the name and lot No. are labeled on the outer package of the product). If the name and lot No. are not legible, your return request will not be accepted. Please specify the product name, model, lot No. and the return reason.

M-52D Diluent

1 Product name

M-52D Diluent

2 Model

M-52D

3 SPECIFICATIONS

20L×1

4 Intended use

M-52D diluent is used in sample dilution, blood cell counting, volume measurement and hemoglobin measurement.

5 Principle

M-52D diluent provides stable solution environment for blood cells, dilute blood sample, generates sheath fluid and provide conductance in blood cell analysis process.

6 Ingredients

Tris Buffer.....≤0.5%

Sodium Chloride.....≤0.1%

7 Storage condition and shelf life

The product must be stored in the temperature of 2°C~30°C, and the relative humidity must not exceed 90%. The storage room must be well ventilated and without corrosive gas. The shelf life is two years.

The working temperature range of the product is consistent with that of its applicable instruments. The open-vial validity is 60 days.

8 Applicable analyzer models

The product applies to BC-5000, BC-5150, BC-5120, BC-5130 and BC-5140 Auto Hematology Analyzers manufactured by Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

9 SAMPLE REQUIREMENT

Fresh human whole blood anticoagulated. Do not use contaminated samples.

10 MATERIALS REQUIRED BUT NOT PROVIDED

The following materials are required but not provided with the product: Mindray-manufactured measurement instruments and matched reagents.

11 Instructions for use

- 1) Restore the product to usage temperature.
- 2) Refer to Mindray Operation Manual to obtain instrument-specific instructions for appropriate use of Auto Hematology Analyzer Reagent.

12 CUT-OFF VALUE/REFERENCE INTERVALS

Not applicable.

13 RESULT ELABORATION

Not applicable.

14 LIMITATIONS

Not applicable.

15 PRODUCT SPECIFICATIONS

- 1) Appearance: transparent liquid, no sediments, grains or flocculi.
- 2) Blank count results: the blank count results of M-52D diluent tested on Mindray Analyzers meet the requirements in Table 1.

Table 1 Requirements of M-52D diluent blank count

Test parameter	Blank count requirements
WBC	$\leq 0.20 \times 10^9/L$
RBC	$\leq 0.02 \times 10^{12}/L$
HGB	$\leq 1g/L$
PLT	$\leq 10 \times 10^9/L$

16 PRECAUTIONS

1. For professional use in vitro diagnosis only.
2. Read the package insert carefully before using this product. It shall be used before the expiration date and disposed of properly when expired.
3. Do not use the reagent if it is frozen.
4. If the blank count is abnormal after the reagent is transported, place it still for 24 hours at room temperature before using it.
5. If the reagent is polluted or affected by other factors and becomes abnormal, stop using it and replace it with a normal one.
6. Dispose of the waste, residual and contaminated packing based on local regulations.
7. The following factors can affect the sample analysis: expired or ineffective reagent; reagent polluted by dust in the air; improper disposal of the sample; mixed with or used with reagents produced by other company; mixed use of residual from the old container and the newly-opened; used in condition other than specified.
8. Take the necessary precautions for the use of the product. Do not swallow. Avoid contact with skin and mucous membranes. If you accidentally take the reagent into your mouth, or the reagents accidentally spill on your skin or into your eyes, wash them off with plenty of water and go seek medical treatment if necessary.
9. Disposal of waste liquid and materials should be in accordance with local guidelines.
10. The Material Safety Data Sheet (SDS) is available upon request
11. All identified risks have been reduced as far as possible by generally acknowledged state of art, and the overall residual risk is acceptable.
12. Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

17 GRAPHICAL SYMBOLS

		
Batch code	Use-by date	Temperature limit
		
HUMIDITY LIMITATION	Manufacturer	<i>In vitro</i> diagnostic medical device

CATALOGUE NUMBER

Authorized representative in
the European Community

Consult instructions for use

European Conformity

Unique Device Identifier

18 REFERENCES

Not applicable.

19 COMPANY CONTACT

Manufacturer: Shenzhen Mindray Bio-Medical Electronics Co., Ltd.
Address: Mindray Building, Keji 12th Road South, High-tech industrial park,
Nanshan, Shenzhen 518057, P.R.China
E-mail Address: service@mindray.com
Tel: +86 755 81888998
Fax: +86 755 26582680

EC-Representative: Shanghai International Holding Corp. GmbH(Europe)
Address: Eiffeustraße 80, 20537 Hamburg, Germany
Tel: 0049-40-2513175
Fax: 0049-40-255726

20 Approval Date of the Operator's Manual

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