# INSTRUCTIONS FOR USE

## CONCERTO™ Detachable Coil and I.D. (Instant Detacher)

This page is intentionally left blank.

## TABLE OF CONTENTS

### ***CONCERTO™ Detachable Coil and I.D. (Instant Detacher)***

---

English . . . . .	5
Français. . . . .	7
Deutsch. . . . .	10
Italiano . . . . .	12
Español . . . . .	15
Svenska. . . . .	17
Nederlands. . . . .	20
Português . . . . .	22
Suomi. . . . .	25
Dansk . . . . .	27
Ελληνικά . . . . .	30
Čeština . . . . .	32
Magyar . . . . .	35
Русский. . . . .	37
Polski . . . . .	40
Türkçe. . . . .	42
Norsk . . . . .	45
Slovenčina . . . . .	47
Română. . . . .	50
Български . . . . .	52
한국어. . . . .	55
عربي . . . . .	59
<b>Symbol Glossary . . . . .</b>	<b>60</b>

This page is intentionally left blank.

# English

## Instructions for Use

EN

### CONCERTO™ Detachable Coil and I.D. (Instant Detacher)

#### CAUTION

Federal (USA) law restricts this device to sale, distribution, and use by or on the order of a physician.

This device should be used only by physicians with a thorough understanding of angiography and percutaneous interventional procedures.

#### DESCRIPTION

The CONCERTO™ Detachable Coil consists of a platinum embolization coil attached to a composite implant delivery pusher with a radiopaque positioning marker and a hand-held I.D. (Instant Detacher) which when activated detaches the coil from the delivery pusher tip. Some of the CONCERTO™ Detachable Coils are enlaced with PGLA fibers or Nylon fibers. I.D. (Instant Detacher) is sold separately.

#### DEVICE COMPATIBILITY

The following devices are required for use with the CONCERTO™ Detachable Coil:

	Type	Diameter (mm)	Minimum micro catheter inside diameter (in)	Detacher
<b>Concerto™ Detachable Coil System</b>	Nylon	2 to 4	0.0165	I.D. (Instant Detacher)
<b>Concerto™ Detachable Coil System</b>	Nylon	5 to 10	0.021	I.D. (Instant Detacher)
<b>Concerto™ Detachable Coil System</b>	PGLA	2 to 10	0.0165	I.D. (Instant Detacher)
<b>Concerto™ Detachable Coil System</b>	PGLA	12 and above	0.021	I.D. (Instant Detacher)

Other Accessory Products (required to perform a procedure)

- 6-8F Guide Catheter\*
- Micro catheter with two marker bands (see above)\*
- Guidewires compatible with micro catheter\*
- Continuous saline/heparin saline flush set\*
- Rotating haemostatic valves (RHV)\*
- 3-Way stopcock\*
- 1-Way stopcock\*
- IV pole\*
- Femoral Sheath\*

\*Not provided as part of the system; chosen based upon physician experience and preference.

#### INDICATIONS FOR USE

The CONCERTO™ Detachable Coil System is indicated for arterial and venous embolizations in the peripheral vasculature.

#### POTENTIAL COMPLICATIONS

Potential complications include, but are not limited to:

- Puncture site hematoma
- Thromboembolic episodes
- Vessel perforation
- Neurological deficits including stroke and death
- Vasospasms
- Vascular thrombosis
- Hemorrhage
- Ischemia

#### WARNINGS

- The CONCERTO™ Detachable Coil, the dispenser track, and the introducer sheath are supplied in a sterile and non-pyrogenic, unopened and undamaged package. The package should be checked for potential damage. Damaged CONCERTO™ Detachable Coils must not be used, as they may result in patient injury.
- CONCERTO™ Detachable Coils are intended for one use only. I.D. (Instant Detacher) is supplied sterile and intended for single patient use. After use do not re-sterilize and/or reuse. Reprocessing or re-sterilization may compromise the structural integrity of the device and/or lead to device failure which in turn, may result in patient injury, illness or death.
- Do not use if sterile packaging has been compromised or damaged.
- Damaged implant delivery pusher and/or coils may affect coil delivery to, and stability inside, the vessel or aneurysm, possibly resulting in coil migration or stretching.
- Do not rotate the implant delivery pusher during or after delivery of the coil into the aneurysm. Rotating the delivery pusher during or after coil delivery into the aneurysm may result in a stretched coil or premature detachment of the coil from the implant delivery pusher, which could result in coil migration.
- Do not use hemostats in an attempt to advance delivery pusher. This may result in a kinked pusher which may lead to premature detachment.
- Verify that the distal shaft of the micro catheter is not under stress before the CONCERTO™ Detachable Coil detachment. Axial compression or tensile forces could be stored in the micro catheter causing the tip to move during the CONCERTO™ Detachable Coil delivery. Micro catheter tip movement could cause the aneurysm or vessel to rupture.
- Advancing the delivery pusher beyond the micro catheter tip once the coil has been deployed and detached involves risk of aneurysm or vessel perforation.

#### WARNINGS

- If undesirable movement of the CONCERTO™ Detachable Coil can be seen under fluoroscopy following coil placement and prior to detachment, remove the coil and replace with another more appropriately sized CONCERTO™ Detachable Coil. Movement of the coil may indicate the coil could migrate once it is detached. Angiographic controls should also be performed prior to detachment to ensure that the coil mass is not protruding into the parent vessel.
- High quality, digital subtraction fluoroscopic road mapping is mandatory to achieve safe catheterization of the aneurysm or vessel and correct placement of the first coil. With smaller aneurysms this is a particularly important step.
- If CONCERTO™ Detachable Coil repositioning is necessary, take special care to retract coil under fluoroscopy in a one-to-one motion with the implant pusher. If the coil does not move with a one-to-one motion, or repositioning is difficult, the coil has been stretched and could possibly break. Gently remove and discard both the catheter and coil.
- Due to the delicate nature of the CONCERTO™ Detachable Coil, the tortuous vascular pathways that lead to certain aneurysms and vessels, and the varying morphologies of intracranial aneurysms, a coil may occasionally stretch while being maneuvered. Stretching is a precursor to potential malfunctions such as coil breakage and migration.
- If resistance is encountered while withdrawing an CONCERTO™ Detachable Coil, which is at an acute angle relative to the catheter tip, it is possible to avoid coil stretching or breaking by carefully repositioning the distal tip of the catheter at the ostium of the aneurysm, or just slightly inside the parent artery.
- Take care not to puncture gloves or sterile drape while handling implant delivery pusher.
- Multiple placements of CONCERTO™ Detachable Coils may be required to achieve the desired occlusion of some aneurysms or vessels.
- The long term effect of this product on extravascular tissues has not been established so care should be taken to retain this device in the intravascular space.
- This device is supplied STERILE for single use only. Do not reprocess or re-sterilize. Reprocessing and re-sterilization increase the risks of patient infection and compromised device performance.

#### PRECAUTIONS

- Handle the CONCERTO™ Detachable Coil with care to avoid damage before or during treatment.
- Do not advance the CONCERTO™ Detachable Coil against a noted resistance until the cause of the resistance is cleared by fluoroscopy. This may lead to the destruction of the coil and/or catheter or perforation of the vessel.
- It is essential to confirm catheter compatibility with the CONCERTO™ Detachable Coil. The outer diameter of the CONCERTO™ Detachable Coil should be checked to ensure that the coil will not block the catheter.
- The dispenser track, introducer sheath and the I.D. (Instant Detacher) are not intended to make contact within the patient.
- Do not use the CONCERTO™ Detachable Coil and I.D. (Instant Detacher) after the expiration date printed on the product label.
- In order to achieve optimal performance of the CONCERTO™ Detachable Coil and to reduce the risk of thromboembolic complications, it is critical that a continuous infusion of appropriate flush solution be maintained.
- Advance and retract CONCERTO™ Detachable Coils slowly and smoothly, especially in tortuous anatomy. Remove the coil if unusual friction or "scratching" is noted. If friction is noted in a second coil, carefully examine both the coil and the catheter for possible damage such as catheter shaft buckling or kinking, or an improperly fused joint.
- If delivery pusher buckling or kinking is noted, grasp the distal most portion of the delivery pusher, distal to the kink, buckle or break and remove from micro catheter.
- Do not advance the coil with force if the coil becomes lodged within or outside the micro catheter. Determine the cause of resistance and remove the system when necessary.
- If resistance is encountered when withdrawing the implant delivery pusher, draw back on the infusion catheter simultaneously until the delivery pusher can be removed without resistance.
- If resistance is noted during coil delivery, remove the system and check for possible damage to the catheter.

#### STORAGE

Store CONCERTO™ Detachable Coils in a cool dry place.

Some of the CONCERTO™ Detachable Coils are enlaced with PGLA fibers. Store these configurations in a cool, dry place with the maximum storage temperature not exceeding 50°C (122°F). A temperature indicator is located on each unit box and pouch for the PGLA configurations of the CONCERTO™ Detachable Coil. If the product has been exposed to temperatures greater than 50°C (122°F), the temperature indicator on the box and/or pouch will change to red.

#### WARNING

- Do not use the product if the temperature indicator on either the pouch or the carton box is red. A red indicator means the product has been exposed to temperature greater than 50°C (122°F). Use of the product which has been exposed to temperature greater than 50°C (122°F), may compromise patient safety.

#### PREPARATIONS FOR USE

- In order to achieve optimal performance of the CONCERTO™ Detachable Coil and to reduce the risk of thromboembolic complication, it is advised that a continuous saline flush be maintained between
    - the femoral sheath and the guiding catheter,
    - the micro catheter and guiding catheter and
    - the micro catheter and the implant delivery pusher and the CONCERTO™ Detachable Coil.
  - Place the appropriate guiding catheter following recommended procedures. Connect a rotating haemostatic valve (RHV) to the hub of the guiding catheter. Attach a 3-way stopcock to the side arm of the RHV, then connect a line for the continuous flush.
  - Attach a second RHV to the hub of the micro catheter. Attach a 1-way stopcock to the side arm of the RHV, then connect a line for continuous flush.
- For CONCERTO™ Detachable Coils:** One drop from the pressure bag every 3-5 seconds is suggested
- For PGLA or Nylon Fiber CONCERTO™ Detachable Coils:** One drop from the pressure bag every 1-3 seconds is suggested.
- Check all fittings so that air is not introduced into guiding catheter or micro catheter during continuous flush.

#### DIAGNOSTIC MR IMAGING


Non-clinical testing has demonstrated that CONCERTO™ Detachable Coils are MR Conditional. CONCERTO™ Detachable Coils can be scanned safely under the following conditions:

- Static magnetic field of 3-Tesla or less
  - Spatial gradient field of 720 Gauss/cm or less
  - Maximum averaged specific absorption rate (SAR) of 3 W/kg for 15 minutes of scanning.
- In non-clinical testing, CONCERTO™ Detachable Coils produced a temperature rise of less than 0.3°C at a maximum averaged specific absorption rate (SAR) of 3 W/kg for 15 minutes of MR scanning in a 3-Tesla MR scanner system (Excite, Software G3.0-052B, General Electric Healthcare, Milwaukee, WI).

## DIRECTIONS FOR USE

1. Slowly and simultaneously remove the CONCERTO™ Detachable Coil and introducer sheath from the dispenser track. Inspect proximal implant delivery pusher for irregularities. If irregularities exist, replace with a new CONCERTO™ Detachable Coil.
2. Slowly advance the CONCERTO™ Detachable Coil out of the introducer sheath into the palm of your gloved hand and inspect for irregularities of the coil or the detachment zone. Due to potential risks of irregularities, a visual proof should be performed. If irregularities exist, replace with a new CONCERTO™ Detachable Coil.

Figure 1

- |                          |                                   |
|--------------------------|-----------------------------------|
| 1. Implant               | 4. Delivery Pusher                |
| 2. Detachment Zone       | 5. Hypotube Break Indicator (HBI) |
| 3. Coil Alignment Marker | 6. Positive Load Indicator        |


3. Gently immerse the CONCERTO™ Detachable Coil and its detachment zone in heparinized saline. Take care not to stretch the coil during this procedure, in order to preserve the coil memory. While still immersed in the heparinized saline, point introducer sheath vertically into saline and gently retract the distal tip of the coil into the introducer sheath.
4. Insert the distal end of the introducer sheath through the rotating hemostatic valve (RHV) and into the hub of the micro catheter until the sheath is firmly seated. Tighten the RHV around the introducer sheath to prevent back flow of blood, but not so tight as to damage the coil during its introduction into the catheter.
5. Transfer the CONCERTO™ Detachable Coil into the micro catheter by advancing the implant pusher in a smooth, continuous manner (1-2cm strokes). Once the flexible portion of the implant pusher has entered the catheter shaft, loosen the RHV and remove the introducer sheath over the implant pusher's proximal end. Do not attempt to advance the entire implant pusher into the introducer sheath this may cause kinking in the proximal pusher. Stop approximately 15 cm from the distal end of implant pusher. Once completed, tighten the RHV around the implant pusher. Leaving the introducer sheath in place will interrupt normal infusion of flushing solution and allow back flow of blood into the micro catheter.
6. Visually verify that the flushing solution is infusing normally. Once confirmed, loosen the RHV enough to advance the implant pusher, but not enough to allow back flow of blood into the implant pusher.
7. Advance the CONCERTO™ Detachable Coil under fluoroscopy and position carefully at the desired site. If coil placement is unsatisfactory, slowly withdraw by pulling on the implant pusher, then slowly advance again to reposition the coil. If the coil size is inappropriate, remove and replace with an appropriately sized coil.
  - a. If delivery pusher buckling or kinking is noted, grasp the distal most portion of the delivery pusher, distal to the kink, buckle or break and remove from micro catheter.

## WARNING

- Do not use hemostats in an attempt to advance delivery pusher. This may result in a kinked pusher which may lead to premature detachment.

8. Continue to advance the CONCERTO™ Detachable Coil until the coil alignment marker of the implant pusher is just distal of the proximal marker of the micro catheter (See Figure 2).
  - a. Advance coil alignment marker just beyond proximal catheter marker band, then retract delivery pusher until coil alignment marker creates a "T" with the proximal catheter marker band under fluoroscopy. This relieves forward stress that may lead to false positive detachment.

Figure 2


- |                                  |                                |
|----------------------------------|--------------------------------|
| 1. Coil Alignment Marker         | 3. Distal Catheter Marker Band |
| 2. Proximal Catheter Marker Band |                                |

9. Tighten RHV to prevent movement of the implant pusher.

10. Remove the I.D. (Instant Detacher) from its protective packaging and place it within the sterile field. I.D. (Instant Detacher) is packaged separately as a sterile device for single patient use only.
11. Confirm again under fluoroscopy that the coil alignment marker of the implant pusher creates a "T" with the proximal marker of the micro catheter.
12. Verify that the RHV is firmly locked around the implant pusher before attaching the I.D. (Instant Detacher) to ensure that the coil does not move during the connection process. Ensure that the implant pusher is straight between the RHV and the I.D. (Instant Detacher). Straightening this section of the implant pusher optimizes alignment to the I.D. (Instant Detacher).
13. Hold the proximal end of the implant pusher at the distal end of the load indicator. Advance the I.D. (Instant Detacher) over the proximal end of the implant pusher until the load indicator fully enters the funnel and the pusher is firmly seated in the actuator (See Figure 3).


Figure 3

- |  |   |
|--|---|
| 1. Proximal end of Load Indicator  | 4. Ready to Detach  |
| 2. Grip pusher here when advancing AXIUM™ I.D. (Instant Detacher) over implant pusher. | 5. RHV and Implant pusher are in a straight line and distal end of Load Indicator is fully inserted in the AXIUM™ I.D. (Instant Detacher) funnel. |
| 3. Advance AXIUM™ I.D. (Instant Detacher) over implant pusher.                         |   |

**Note:** If indicator band is still showing as in Figure 4 below, the I.D. (Instant Detacher) should be advanced further until the implant pusher is fully seated in the funnel as in Figure 3 above.


Figure 4

1. Not Fully Loaded
14. To detach coil, place the I.D. (Instant Detacher) into palm and retract the Thumb-slide back until it stops and clicks, and slowly allow the Thumb-slide to return to its original position. Remove the I.D. (Instant Detacher).
 

**Note:** The I.D. (Instant Detacher) can also be removed at the end of stroke if desired. To remove the I.D. (Instant Detacher) at the end of stroke, hold the thumb slide in its most rearward position and remove the I.D. (Instant Detacher) (see Figure 5).

Figure 5


- |                    |  |
|--------------------|--|
| 1. Ready to Detach | 2. Distal end of Load indicator fully in I.D. (Instant Detacher) funnel. |
|--------------------|--|

15. Successful coil detachment must be verified by fluoroscopic monitoring to ensure that the coil has detached. Slowly pull back the implant pusher while watching the fluoroscopy to make sure that the coil does not move. In the unlikely event the coil moves, repeat steps 12-14. If necessary, advance the implant pusher to re-establish the coil and catheter marker alignment. Verify coil detachment as above.
16. If you wish to confirm detachment, grasp the positive load indicator between your thumb and forefinger of your left hand and the proximal end of the implant delivery pusher with your thumb and forefinger of your right hand. Gently pull on the proximal end of the implant delivery pusher. If it moves freely from the hypo-tube, the system is detached properly. If it does not, repeat steps 13-15.
 


**Note:** If the Coil does not detach after 3 attempts, discard the I.D. (Instant Detacher) and replace with a new I.D. (Instant Detacher).
17. In the rare event that the coil does not detach and cannot be removed from the implant delivery pusher, use the following steps for detachment.
  - a. Grip the hypotube approximately 5 cm distal of the positive load indicator at the hypotube break indicator and bend the implant delivery pusher just distal to the HBI 180 degrees.
  - b. Next straighten the pusher back, continue bending and straightening until the pusher tubing opens exposing the release element (Figure 6).

Figure 6

- c. Gently separate the proximal and distal ends of the open pusher. Then, under fluoroscopy, pull the proximal portion of the implant delivery pusher approximately 2-3 cm to confirm implant detachment per IFU (Figure 7).

Figure 7

18. Once coil detachment had been detected and fluoroscopically confirmed slowly withdraw the implant pusher from the micro catheter.

### WARNING

In case of:

- a. **A False Positive** (Attempted Detachment failed), remove coil from treatment area and micro catheter and replace with a new CONCERTO™ Detachable Coil.
  - b. **A False Negative** (Coil becomes prematurely detached), remove implant pusher and:
    - i. Advance next coil to push remaining tail of prematurely detached coil into treatment area
    - ii. Remove prematurely detach coil with the appropriate retrieval device.
19. Repeat steps 1-18, if additional coil placements are required.
20. Once the procedure is complete, discard the I.D. (Instant Detacher).

### WARNING

- The I.D. (Instant Detacher) is intended for a maximum of 25 cycles.

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	Sterilized using ethylene oxide		Keep away from sunlight
	Do not re-use		Keep dry
	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician		Upper limit of temperature
	Do not resterilize		Authorized representative in the European Community
	Caution		Catalogue number
	<p>MR Conditional Compatible avec l'IRM sous certaines conditions Bedingt MR-tauglich Compatibile con la RM Condicional RM MR-villkorlig Voorwaardelijk MRI-veilig Condicionado ao ambiente de ressonância magnética (RM) Turvallinen magneettikuvauksessa tietyn edellytyksin MR-betinget Συμβατό με τη μαγνητική τομογραφία υπό προϋποθέσεις</p> <p>Podmíněně kompatibilní s magnetickou rezonancí MR-feltételes MP-совместимо Warunkowe stosowanie w środowisku MR MR Koşullu MR-betinget Podmienené bezpečné v prostredí MR Nu prezintă probleme de rezonanță magnetică în condiții specifice de utilizare Условно зависимо от MP MR 조건부 آمن للاستخدام في التصوير بالرنين المغناطيسي</p>		Manufacturer
			Use-by date
			Batch code
	<p>Non-pyrogenic Apyrogène Nichtpyrogen Il prodotto non è pirogeno No pirogénico Ikke-pyrogen Niet-pyrogeen Apirogénico Pyrogeeniton Ikke-pyrogen Μη πυρογόνο</p> <p>Apyrogenni Nem pirogén Апирогенно Niepirogenny Pirojenik deçildir Ikke-pyrogen Nepyrogeenne Nepirogenic Апирогенно 비발열성 لا يسبب الإصابة بالحمى</p>		<p>Contents of Package Contenu de l'emballage Packungsinhalt Contenuto della confezione Contenido del embalaje Förpackningens innehåll Inhoud van de verpakking Conteúdo da embalagem Pakkauksen sisältö Pakkens indhold Περιεχόμενα της συσκευασίας</p> <p>Obsah balení A csomag tartalma Содержимое упаковки Zawartość opakowania Ambalajın İçindekiler Pakkens innhold Obsah balenia Conținutul ambalajului Съдържание на опаковката 패키지 내용물 محتويات العبوة</p>

Manufacturer  
Micro Therapeutics, Inc.  
d/b/a ev3 Neurovascular  
9775 Toledo Way  
Irvine, CA 92618  
USA  
Tel: 1.949.837.3700  
[www.ev3.net](http://www.ev3.net)