

virotype[®] CSFV 2.0 RT-PCR Kit Validation Report

For the detection of RNA from the *Classical Swine
Fever Virus* (CSFV)



Contents

- 1 Introduction 4
 - 1.1 Intended use 4
 - 1.2 General information 4
 - 1.3 Description of the test principle 4
 - 1.4 Kit contents 5
 - 1.5 Storage 5
 - 1.6 Equipment and reagents to be supplied by user 6
 - 1.7 RNA extraction 6
 - 1.8 Important notes 8
- 2 Procedure 9
 - 2.1 Important points before starting 9
 - 2.2 Things to do before starting 9
 - 2.3 Test procedure 9
- 3 Data interpretation 11
- 4 Characteristics of the test 13
 - 4.1 Analytical sensitivity 13
 - 4.1.1 Analytical sensitivity using the Bio-Rad CFX96 instrument 13
 - 4.1.2 Analytical sensitivity – Limit of detection 16
 - 4.1.3 Analytical sensitivity of pooled samples 18
 - 4.2 Specificity 19
 - 4.2.1 Comparative analysis of the virotype CSFV 2.0 RT-PCR Kit and other commercially available RT-PCR Kits using reference samples (in-house testing) 19
 - 4.2.2 Detection of CSFV and other genetically related *pestiviruses* 20
 - 4.2.3 Discrimination of pathogens for differential diagnosis 24
 - 4.3 Diagnostic sensitivity, specificity and efficiency 26
 - 4.3.1 Definition diagnostic sensitivity 26
 - 4.3.2 Definition diagnostic specificity 26
 - 4.3.3 Definition diagnostic efficiency 26
 - 4.3.4 Validation of the virotype CSFV 2.0 RT-PCR Kit 26
 - 4.4 Repeatability 29
 - 4.4.1 Intra-assay variance 30
 - 4.4.2 Inter-assay variance 33

4.4.3	Batch-to-batch comparison	36
4.4.4	Comparison of real-time PCR thermocyclers	39
4.5	Stability testing	42
4.5.1	Freeze-thaw-cycles	42
4.5.2	Heparin inhibition.....	44
4.5.3	EDTA inhibition	45

1 Introduction

1.1 Intended use

The virotype CSFV 2.0 RT-PCR Kit is intended for the detection of RNA from *Classical Swine Fever Virus* (CSFV) in blood, tissue, swabs, and samples in stabilizing transport media from pigs and wild boar.

The kit is approved by the Friedrich-Loeffler-Institut and licensed in accordance with § 11 (2) of the German Animal Health Act (FLI-C 113) for use in Germany for veterinary diagnostic procedures.

For veterinary use only.

1.2 General information

The virotype CSFV 2.0 RT-PCR Kit is a highly sensitive and specific solution for the detection of RNA from *Classical Swine Fever Virus* (CSFV) in samples from pigs and wild boar. Classical Swine Fever (CSF) is economically one of the most important viral infectious diseases of swine. CSF is widespread in domestic pig and wild boar populations. CSF is an internationally notifiable animal disease. The causative agent, *Classical Swine Fever Virus*, is a single-stranded RNA virus and a member of the genus *Pestivirus* which belongs to the *Flaviviridae* family like *Bovine Viral Diarrhea Virus* (BVDV) in cattle and *Border Disease Virus* (BDV) in sheep.

The high sensitivity of the virotype CSFV 2.0 RT-PCR Kit allows the early detection of the pathogen in individual or pooled blood samples, in individual or pooled tissue samples, as well as from swabs, and samples in stabilizing transport media from pigs and wild boar. The virotype CSFV 2.0 RT-PCR Kit detects all known CSFV-strains of various genotypes. In rare cases the test may detect CSF vaccine virus. Positive results of animals from areas where CSFV vaccination is performed should therefore be verified.

1.3 Description of the test principle

Polymerase chain reaction (PCR) is based on the amplification of specific regions of the pathogen genome. In real-time PCR, the amplified product is identified using fluorescent dyes. These are usually linked to oligonucleotide probes that bind specifically to the amplified product. Monitoring the fluorescence intensities during the PCR run (i.e., in real-time) allows detection of the accumulating product without the need to re-open the reaction tubes afterward.

The virotype CSFV 2.0 RT-PCR Kit contains all the necessary reagents for the detection of CSFV RNA, including a Positive and Negative Control. With this kit, both reverse transcription and PCR are performed in one reaction tube, reducing the risk of contamination.

An Internal Control excludes the possibility of false-negative results. The kit uses two specific primer/probe combinations:

- FAM™ fluorescence for RNA of CSFV
- JOE™ fluorescence for the endogenous Internal Control (EC; β -actin present within the sample)

A Positive Control serves to verify the functionality of the reaction mix for the amplification of one of the CSFV RNA target.

1.4 Kit contents

virotype CSFV 2.0 RT-PCR Kit	(96)
Cat. no.	VT281825
Number of reactions	96
Master Mix (tube with orange cap), includes primers, probes and enzymes	2 x 980 μ l
Positive Control (tube with red cap)	1 x 150 μ l
Negative Control (tube with blue cap)	1 x 150 μ l
Handbook	1

1.5 Storage

The components of the virotype CSFV 2.0 RT-PCR Kit should be stored at -30°C to -15°C and are stable until the expiration date stated on the label. Avoid repeated thawing and freezing (> 3x), as this may reduce assay sensitivity. Freeze the components in aliquots if they will only be used intermittently.

1.6 Equipment and reagents to be supplied by user

When working with chemicals, always wear a suitable lab coat, disposable gloves and protective goggles. For more information, consult the appropriate safety data sheets (SDSs), available from the product supplier.

- Pipets
- Nuclease-free, aerosol-resistant pipet tips with filters
- Sterile 1.5 ml Eppendorf® tubes
- Nuclease-free (RNase/DNase-free) consumables. Special care should be taken to avoid nuclease contamination of all reagents and consumables used to set up PCR for sensitive identification of viral nucleic acids
- Cooling device or ice
- Benchtop centrifuge with rotor for 1.5 ml tubes
- Real-time cycler with appropriate fluorescent channels
- Appropriate software for chosen real-time cycler
- Appropriate strip tubes and caps or 96-well optical microplate with optical sealing film or cover for chosen real-time cycler

1.7 RNA extraction

The virotype CSFV RT-PCR Kit can be used for the detection of CSFV RNA from blood, tissue, swabs, and samples in stabilizing transport media from pigs and wild boar.

Due to the high sensitivity of the test individual or pooled samples can be tested. For domestic pigs, up to 20 individual serum, plasma, or EDTA-blood samples or up to 10 tissue samples can be used. Furthermore, for wild boar, pools can consist of up to 10 serum, plasma, EDTA-blood, or tissue samples. Also swab samples and samples in stabilizing transport media from pigs and wild boar can be tested.

Note: For use in Germany the specifications described in the „*Amtliche Methodensammlung*“ apply.

Prior to real-time PCR, viral RNA must be extracted from the starting material. INDICAL offers a range of validated kits for the extraction of RNA from animal samples.

Extraction based on magnetic beads:

- IndiMag® Pathogen Kit* (SP947457)
- IndiMag Pathogen Kit* w/o plastics (SP947257)
- IndiMag Pathogen IM48 Cartridge (SP947654P608, SP947654P224)
- IndiMag Pathogen KF96 Cartridge (SP947855P196)

Extraction based on spin columns:

- IndiSpin® Pathogen Kit* (SP54104, SP54106)
- IndiSpin QIAcube® HT Pathogen Kit (SP54161)

* suitable for simultaneous extraction of CSFV RNA und ASFV DNA

Note: When using difficult sample material, it is recommended to use INDICAL's „Pretreatment T4 (phenol extraction)“.

If real-time RT-PCR is not performed immediately after extraction, store the RNA at -20°C or at -80°C for longer storage.

For further information on automated and manual extraction of CSFV RNA from different sample types, refer to the respective handbook or contact INDICAL Support at support@indical.com.

1.8 Important notes

General precautions

The user should always pay attention to the following:

- Use nuclease-free pipet tips with filters.
- Store and extract positive materials (specimens, positive controls and amplicons) separately from all other reagents and add them to the reaction mix in a spatially separated facility.
- Thaw all components on ice before starting the assay.
- When thawed, mix the components by inverting and centrifuge briefly.
- Do not use components of the test kit past the expiration date.
- Keep samples and controls on ice or in a cooling block during the setup of reactions.

Negative control

At least one negative control reaction should be included in each PCR run, containing all the components of the reaction except for the pathogen template. This enables assessment of contamination in the reaction.

Positive control

When performing PCR on unknown samples, it is recommended to include a positive control reaction in the PCR run, containing a sample that is known to include the targeted viral RNA. A positive control serves to prove the functionality of the pathogen assay, for example, the correct setup of the reaction mix. Use 5 µl of the Positive Control provided with the virotype CSFV 2.0 RT-PCR Kit to test for successful amplification of the target.

Extraction and amplification control

For increased process safety and convenience, one extraction and amplification control assay is included in the test kit.

An endogenous internal control (EC) detects the β -actin gene present within the sample. This allows extraction and amplification to be monitored.

2 Procedure

2.1 Important points before starting

- Please read „Important notes“ before starting.
- Include at least one positive control (Positive Control) and one negative control (Negative Control) per PCR run.
- Before beginning the procedure, read through the protocol and ensure that you are familiar with the operation of the chosen real-time PCR cyclers.
- RNA is unstable. Perform the protocol without interruption.

2.2 Things to do before starting

- Thaw all reagents on ice and protect from light.
- Maintain reagents on ice during PCR setup.
- Maintain reagents on ice or in a cooling block during PCR setup.

2.3 Test procedure

1. Before use, mix the Master Mix by inverting 5 times or until mixed thoroughly, then centrifuge briefly to collect the fluids.
2. Pipet 20 µl of the Master Mix into each reaction tube. Then add 5 µl of the sample RNA (Table 1).

Include positive and negative control reactions.

Positive Control: Use 5 µl of the positive control (Positive Control) instead of sample RNA.

Negative Control: Use 5 µl of the negative control (Negative Control) instead of sample RNA.

Table 1. Preparation of reaction mix

Component	Volume
Master Mix	20 µl
Sample	5 µl
Total volume	25 µl

3. Close the reaction tubes or seal the plate and invert 5 times or until mixed thoroughly. Then centrifuge briefly to collect the fluids.

4. Set the filters for the reporter dyes in the software of your thermal cycler according to Table 2.

Table 2. Filter settings for the reporter

Pathogen/ Internal Control	Reporter
CSFV	FAM
Endogenous Internal Control	HEX/ JOE™ ¹
Passive reference ²	Texas Red/ ROX™

¹ Use the option appropriate for your thermal cycler.

² Internal reference for use with Applied Biosystems® ABI PRISM® Sequence Detection Systems

5. Run the real-time PCR protocol according to Table 3.

Table 3. Real-time RT-PCR protocol for CSFV 2.0.

Step	Temperature	Time	Number of cycles
Reverse Transcription	50°C	10 min	1
Initial Activation	95°C	2 min	1
2-step cycling			
Denaturation	95°C	5 s	40
Annealing/Extension*	60°C	30 s	

* Fluorescence data collection, approximate run time 67 min (Mx3005P®, Agilent Technologies, Inc.)

Note: The above protocol can also be used in combination with the virotype ASFV 2.0 PCR Kit.

3 Data interpretation

Interpretation of results

For the assay to be valid, the Positive Control must give a signal in the FAM and HEX/ JOE channels with a $C_T < 35$. The Negative Control must give no signal.

The following results are possible if working with unknown samples. The possible sample results are also summarized in Table 4.

The sample is positive for CSFV, and the assay is valid, if the following criteria are met:

- The sample yields a signal in the FAM channel (regardless of any signal in the HEX/ JOE channel).
- The Positive Control yields a signal in the FAM and HEX/ JOE channel.
- The Negative Control does not yield a signal in the FAM and HEX/ JOE channel.

Note that very high concentrations of CSFV RNA in the sample may lead to a reduced signal or no signal for the endogenous Internal Control (EC; HEX/ JOE channel) due to competition with the internal control.

The sample is negative for CSFV, and the assay is valid, if the following criteria are met:

- The sample **does not yield** any signal in the FAM channel.
- The sample yields a signal in only the HEX/ JOE channel.
- The Positive Control yields a signal in the FAM and HEX/ JOE channels.
- The Negative Control does not yield a signal in the FAM and HEX/ JOE channels.

The sample results are inconclusive, and the assay is invalid, if the following occurs:

- The sample yields no signal in the FAM and HEX/ JOE channels.

If no signal is detected in the FAM (CSFV) and the HEX/ JOE (endogenous Internal Control, EC) channel, the result is inconclusive. The absence of a signal for the housekeeping gene indicates strong PCR inhibition and/or other malfunctions, e.g., during extraction.

To check for inhibition, we recommend 1:5 dilution of the sample RNA in nuclease free water, to repeat the RNA extraction, or repeat the whole test procedure starting with new sample material.

Check that there is a fluorescence signal in the FAM channel for the positive control reaction (Positive Control). Absence of a signal for the Positive Control indicates an error, which could be due to incorrect setup of the reaction mix or incorrect cycling conditions.

Table 4. Results interpretation table*

Sample result	FAM (CSFV)	HEX/JOE (EC)
CSFV positive	X	X
CSFV strong positive	X	
CSFV negative		X
Inconclusive		

* Interpretation of sample results can be determined provided positive and negative control reactions are performed. The positive control must yield a signal in the FAM and HEX/ JOE channels. The negative control must yield no signal in any channel. For a complete explanation of possible sample results, please refer to "Data interpretation".

4 Characteristics of the test

4.1 Analytical sensitivity

4.1.1 Analytical sensitivity using the Bio-Rad CFX96 instrument

The high analytical sensitivity of the virotype CSFV 2.0 RT-PCR Kit was verified by a titration series of *Classical Swine Fever Virus* (CSFV) *in vitro* RNA [$10^6 - 1$ copies/well], performed in triplicates of relevant dilutions using the virotype CSFV 2.0 RT-PCR protocol on the Bio-Rad CFX96 instrument.

Results / Conclusion

The virotype CSFV 2.0 RT-PCR Kit is able to detect up to ten CSFV virus copies per sample (Table 5, Figure 1 and Figure 2). There is a high correlation between RNA copy number and amplification results. A correlation coefficient of 0.985 with an efficiency of 108.1 % for the *in vitro* RNA was calculated when using the virotype CSFV 2.0 RT-PCR Kit on the Bio-Rad CFX96 instrument (Figure 2).

Table 5. Individual and mean C_T values of **CSFV (FAM)** *in vitro* RNA titration series in triplicates. The test was performed on the Bio-Rad CFX96 instrument using the virotype CSFV 2.0 real-time RT-PCR protocol.

Type	Copy number	C _T (FAM)	C _T mean	SD	Result
Standard	10 ⁶	22.51			+
Standard	10 ⁶	22.60	22.55	0.05	+
Standard	10 ⁶	22.55			+
Standard	10 ⁵	25.87			+
Standard	10 ⁵	25.73	25.70	0.19	+
Standard	10 ⁵	25.49			+
Standard	10 ⁴	28.89			+
Standard	10 ⁴	28.39	28.75	0.32	+
Standard	10 ⁴	28.97			+
Standard	10 ³	31.76			+
Standard	10 ³	31.84	31.77	0.07	+
Standard	10 ³	31.70			+
Standard	100	35.10			+
Standard	100	34.39	34.78	0.36	+
Standard	100	34.85			+
Standard	50	36.62			+
Standard	50	34.93	35.83	0.85	+
Standard	50	35.95			+
Standard	25	38.41			+
Standard	25	36.54	37.08	1.16	+
Standard	25	36.30			+
Standard	10	39.83			+
Standard	10	-	39.83	-	-
Standard	10	-			-
Standard	1	-			-
Standard	1	-	-	-	-
Standard	1	-			-

SD = standard deviation, - = no C_T




Figure 1. Individual values of a titration series of **CSFV (FAM) *in vitro*** RNA in triplicates. The test was performed on the Bio-Rad CFX96 instrument using the virotype CSFV 2.0 real-time RT-PCR protocol.




Figure 2. Standard curve of obtained C_T values for a titration series of **CSFV (FAM) *in vitro*** RNA. The test was performed on the Bio-Rad CFX96 instrument using the virotype CSFV 2.0 real-time RT-PCR protocol.

4.1.2 Analytical sensitivity – Limit of detection

The limit of detection (LOD) for the target sequence of the *Classical Swine Fever Virus* was determined by testing individual titration series of *in vitro* RNA of this sequence in octuplicates. The limit of detection with 95 % confidence interval (LOD_{95 %}: mean number of copies yielding a probability of detection of 0.95) was determined using the web tool <https://quodata.de/content/validation-qualitative-pcr-methods-single-laboratory>.

Results / Conclusion

Results are summarized in Table 6 and Figure 3. Using the virotype CSFV 2.0 RT-PCR Kit, a high correlation between RNA copy number and the amount of amplified product was demonstrated for the CSFV targeted sequence.

The LOD_{95%} is 32.9 copies/reaction with a 95 % confidence interval of [18.632, 58.582] (Figure 3).

Table 6. Limit of detection for CSFV *in vitro* RNA tested in octuplicates on the Bio-Rad CFX96 instrument using the virotype CSFV 2.0 real-time RT-PCR protocol.

Copies/test	Total number of replicates	Number of replicates positive	Number of replicates negative
1000,000	8	8	0
100,000	8	8	0
10,000	8	8	0
1,000	8	8	0
100	8	8	0
50	8	8	0
25	8	8	0
10	8	3	5
1	8	1	7
0.1	8	0	8




Figure 3. POD (probability of detection) curve and LOD_{95%} for CSFV. The blue diamonds characterize the laboratory-specific rates of detection. The blue curve denotes the mean POD curve along with the corresponding 95 % confidence range highlighted as the grey band. The POD curve under ideal conditions is displayed as the black dashed curve.

4.1.3 Analytical sensitivity of pooled samples

Two positive CSFV samples were diluted in a CSFV-negative RNA sample from domestic pig blood to simulate pools of 5, 10 and 20 samples. The resulting pools were tested with the virotype CSFV 2.0 RT-PCR Kit on the Bio-Rad CFX96 instrument using the virotype CSFV 2.0 RT-PCR protocol.

Results / Conclusion

The C_T values of the CSFV (FAM) and the endogenous Internal Control (JOE) are shown in Table 7. Pooled samples up to 20 could be detected with the virotype CSFV 2.0 RT-PCR Kit.

Table 7. Analysis of the **CSFV (FAM)** and the **endogenous Internal Control (JOE)** signals for simulated pool samples with the virotype CSFV 2.0 RT-PCR Kit tested on the Bio-Rad CFX96 instrument using the virotype CSFV 2.0 real-time RT-PCR protocol.

Pool size	Sample	Species	virotype CSFV 2.0 RT-PCR Kit	
			C_T CSFV	C_T Internal Control
0	FLI-1	Wild boar	29.51	22.35
5			33.26	22.46
10			34.51	22.19
20			35.33	22.54
0	CSFV 1047	Pig	27.67	22.72
5			31.76	22.94
10			32.69	22.57
20			33.09	22.86

4.2 Specificity

4.2.1 Comparative analysis of the virotype CSFV 2.0 RT-PCR Kit and other commercially available RT-PCR Kits using reference samples (in-house testing)

Twenty-one CSFV positive RNA samples (Table 8) from cell culture, wild boars and pigs and dilutions thereof were provided by PIWet (Polish National Reference Laboratory for CSFV / PAŃSTWOWY INSTYTUT WETERYNARYJNY), the European Reference Laboratory for CSFV (TiHo Hannover) as well as the German National Reference Laboratory for CSFV (Friedrich-Loeffler-Institut, Riems, Germany). The RNA samples were tested using the virotype CSFV 2.0 RT-PCR Kit on the Bio-Rad CFX96 instrument.

Results / Conclusion

Results obtained with the virotype CSFV 2.0 RT-PCR Kit show better or equal sensitivity compared to some competitor RT-PCR kits, especially kits A and C (Table 9).

Table 8. Samples for comparative testing of the virotype CSFV 2.0 RT-PCR Kit and other commercially available RT-PCR Kits

Sample	CSFV status	Origin	Material	Strain (genotype)	Dilution
1	positive	PIWet	Cell culture	Alfort 187 (gt 1.1)	1:10
2	positive	PIWet	Cell culture	Brescia (gt 1.1)	1:10
3	positive	EURL	Pig blood	CSF385 (gt 2.3)	1:100
4	negative	INDICAL	Wild boar blood	-	-
5	positive	EURL	Pig serum	CSF1060 (gt 2.2)	1:10
6	positive	EURL	Pig plasma	CSF309 (gt 3.4)	1:10
7	NC	INDICAL		Negative Control (virotype CSFV 2.0 RT-PCR Kit, alpha test)	
8	PC	INDICAL		Positive Control (virotype CSFV 2.0 RT-PCR Kit, alpha test)	
9 - 16		FLI (FLI-1)	Pig blood	unknown	10 ⁻¹ - 10 ⁻⁶
17 - 24		FLI (FLI-2)	Pig tonsil	Israel (gt 2.1)	10 ⁻¹ - 10 ⁻⁵

NC = Negative Control, PC = Positive Control, PIWet = PAŃSTWOWY INSTYTUT WETERYNARYJNY, EURL = European Reference Laboratory, FLI = Friedrich-Loeffler-Institut

Table 9. Comparative analysis of the virotype CSFV 2.0 RT-PCR Kit and competitor kits for CSFV (FAM signal).

Sample	Material	Status	C _T (CSFV)					
			virotype CSFV 2.0 RT-PCR Kit	Kit A	Kit B	Kit C	Kit D	
#1	Cell culture	pos	30.26	31.54	30.05	33.14	31.21	
#2	Cell culture	pos	27.09	27.20	27.45	30.36	28.17	
#3	Blood	pos	32.77	-	35.92	-	36.47	
#4	Blood	neg	-	-	-	-	-	
#5	Blood	pos	29.77	29.99	30.97	32.58	30.99	
#6	Blood	pos	30.10	30.23	31.74	33.40	31.18	
#7	NC	neg	-	-	-	-	-	
#8	PC	pos	29.47	31.21	-	-	-	
#9	Blood	FLI 1 (10 ⁻¹)	20.38	20.80	21.72	25.23	22.65	
#10		FLI 1 (10 ⁻²)	23.50	23.75	25.22	29.01	26.75	
#11		FLI 1 (10 ⁻³)	26.42	27.15	28.94	32.18	30.12	
#12		FLI 1 (10 ⁻⁴)	29.51	30.39	31.83	35.68	32.80	
#13		FLI 1 (10 ⁻⁵)	33.26	35.26	35.28	39.70	37.29	
#14		FLI 1 (10 ^{-5.5})	34.51	35.97	35.61	-	36.33	
#15		FLI 1 (10 ^{-5.8})	35.33	-	36.34	-	36.39	
#16		FLI 1 (10 ⁻⁶)	36.72	-	37.58	-	38.23	
#17		Tonsil	FLI 2 (10 ⁻¹)	24.74	26.31	26.72	28.75	28.52
#18			FLI 2 (10 ⁻²)	27.67	30.14	29.72	32.56	31.44
#19			FLI 2 (10 ⁻³)	31.76	34.47	33.22	38.39	35.42
#20			FLI 2 (10 ^{-3.5})	32.69	35.56	33.79	36.38	34.96
#21			FLI 2 (10 ^{-3.8})	33.09	-	34.87	37.73	37.39
#22			FLI 2 (10 ⁻⁴)	34.74	-	36.57	38.72	37.82
#23			FLI 2 (10 ^{-4.5})	36.15	-	36.85	-	38.37
#24			FLI 2 (10 ⁻⁵)	-	-	-	-	-

4.2.2 Detection of CSFV and other genetically related *pestiviruses*

RNA samples from the Epizone panel of different relevant CSFV genotypes (Friedrich-Loeffler-Institute) were tested with the virotype CSFV 2.0 RT-PCR Kit on the Bio-Rad CFX96 instrument using the virotype CSFV 2.0 RT-PCR protocol.

Results/Conclusion

All available CSFV genotypes were correctly detected and no cross-reactivity to other relevant pestiviruses was observed when using the virotype CSFV 2.0 RT-PCR Kit (Table 10).

Table 10. Testing cell culture samples of the Epizone panel with of the virotype CSFV 2.0 RT-PCR on the Bio-Rad CFX96 instrument using the virotype CSFV 2.0 real-time RT-PCR protocol.

	Sample	Dilution tested	virotype CSFV 2.0 RT-PCR Kit C_T (CSFV)
<i>Classical Swine Fever Virus (CSFV)</i>	CSFV-C strain_gt1.1	1:10	33.38
	CSFV-Eystrup91_gt1.1	1:100	38.52
	CSFV-Alfort187_gt1.1	1:100	38.52
	CSFV-Koslov1128_gt1.2	1:100	37.73
	CSFV-Brescia_gt1.2	1:10	33.37
	CSFV-Schweiz II_gt2.1	1:100	37.80
	CSFV-Pader_gt2.1	1:100	39.62
	CSFV-Bergen_gt2.2	1:10	33.62
	CSFV-D4886/82/Ro_gt2.2	-	31.31
	CSFV-Uelzen_gt2.3	-	29.69
	CSFV-Spante_gt2.3	1:10	31.81
	CSFV-Congenital Tremor_gt3.1	-	32.01
	CSFV-Kanagawa_gt3.4	-	31.15
	CSFV-CSF 1027	1:100	24.10
	CSFV-CSF 1024	1:100	28.08
	CSFV-CSF 0867	1:100	24.49
	CSFV-CSF 0866	1:100	28.34
CSFV-CSF 0864	1:100	23.66	

	CSFV-CSF1015	1:100	24.09
	CSFV-CSF 0852	1:100	27.62
	CSFV-CSF 0854	1:100	23.69
	CSFV-CSF 848	1:100	30.23
	CSFV-CSF 847	1:100	25.22
	CSFV-CSF 0840	1:100	26.49
	CSFV-CSF 0838	1:100	29.81
	CSFV-CSF 0822	1:100	25.01
Border Disease Virus (BDV)	BDV-Moredun_gt1	-	-
	BDV-Rudolph_gt2	-	-
	BDV-Gifhorn_gt3	-	-
	BDV-Isard_gt4	-	-
	BDV-2 ST 1507	-	-
	BDV-1 SF6/87	-	-
	BDV-1 137/4	-	-
Bovine Viral Diarrhea Virus (BVDV)	BVDV-1-WUS 5708	-	-
	BVDV-1-BO806-17	-	-
	BVDV-1-Arnsby 1599	-	-
	BVDV-1b-Grub	-	-
	BVDV-1-17R507	-	-
	BVDV-1-NADL_gt1a	-	-
	BVDV-1-Paplitz_gt1b	-	-
	BVDV-1-PI809_gt1d	-	-
	BVDV-1-NC3807-1251/1_gt1e	-	-
	BVDV-1-Egbert_gt1f	-	-
	BVDV-1-BO806-17_gt1g	-	-
	BVDV-1-BO807-3_gt1h	-	-
	BVDV-1-NC3807-8757_gt1x	-	-
	BVDV-2-Bure	-	-
	BVDV-2-10/01	-	-
	BVDV-2-P01600	-	-

	BVDV-2-Vepfe	-	-
	BVDV-2-8644_gt2a G	-	-
	BVDV-2-Bure_gt2a US	-	-
	BVDV-2-Walter_gt2b	-	-
	BVDV-2-PO1600_gt2c	-	-
	BVDV-Bayonce	-	-
<i>Pestivirus</i>	Hobi_gtatypical	-	-
	Giraffe H138_gtatypical	-	-

- = no Ct

4.2.3 Discrimination of pathogens for differential diagnosis

Cross-reactivity was tested with samples positive for *African Swine Fever Virus* (ASFV), *Porcine Circovirus-2* (PCV-2), *Swine Influenza Virus* (SIV) und *Porcine Reproductive and Respiratory Syndrome Virus* (PRRSV). The samples were kindly provided by the Friedrich-Loeffler-Institut and other State Veterinary Laboratories.

Results/Conclusion

No cross-reactivity to other relevant porcine viral pathogens was detected using the virotype CSFV 2.0 RT-PCR Kit (Table 11).

Table 11. Cross-reactivity of the virotype CSFV 2.0 RT-PCR Kit to other swine-related pathogens on the Bio-Rad CFX96 instrument using the virotype CSFV 2.0 real-time RT-PCR protocol.

	Sample	Sample material	virotype CSFV 2.0 RT-PCR Kit	Reference assay*
			C_T (CSFV)	C_T (Pathogen)
African Swine Fever Virus (CSFV)	Arm07 (gt II)	Lymph node	-	23.77
	Sardinia-ws12-4	Blood	-	17.08
	Kenia05-hs07-7	Blood	-	25.82
Porcine Circovirus-2 (PCV-2)	PCV-2_1	Serum	-	29.73
	PCV-2_2	Serum	-	30.76
	PCV-2_3	Serum	-	38.87
	PCV-2_4	Serum	-	24.63
	PCV-2_5	Serum	-	29.54
Swine Influenza Virus (SIV)	SIV-01	Blood	-	22.81
	SIV-02	Blood	-	26.62
	SIV-03	Blood	-	29.77
Porcine Reproductive and Respiratory Syndrome Virus (PRRSV)	Intervet 10-3	Culture	-	22.38
	Stendal V852-10-3	Culture	-	22.34
	USA 18-18-10-4	Culture	-	26.09
	PRRSV-22	Serum	-	31.59
	PRRSV-23	Serum	-	24.42
	PRRSV-24	Serum	-	25.33

* Reference assay were as follows: virotype ASFV 2.0 PCR Kit (ASFV samples), virotype PCV2/PCV3 Reagent (PCV-2 samples), virotype Influenza A RT-PCR Kit (SIV samples), virotype PRRSV RT-PCR Kit (PRRSV samples)

4.3 Diagnostic sensitivity, specificity and efficiency

4.3.1 Definition diagnostic sensitivity

Percentage of positive samples in the new test of a population of true positive samples. True positive samples giving negative results in the new test are termed false negative.

Calculation: $[\text{true positives} / (\text{true positives} + \text{false negatives})] * 100$

4.3.2 Definition diagnostic specificity

Percentage of negative samples in the new test of a population of true negative samples. True negative samples giving positive results in the new test are termed false positive.

Calculation: $[\text{true negatives} / (\text{false positives} + \text{true negatives})] * 100$

4.3.3 Definition diagnostic efficiency

Diagnostic efficiency refers to the amount of agreement between the results from the new test and those from the reference test. It is expressed as a proportion of correctly identified samples among all samples.

Calculation: $[(\text{true positives} + \text{true negatives}) / (\text{true positives} + \text{true negatives} + \text{false positives} + \text{false negatives})] * 100$

4.3.4 Validation of the virotype CSFV 2.0 RT-PCR Kit

For validation of the virotype CSFV 2.0 RT-PCR Kit, $n = 472$ samples (blood/ blood swab ($n = 156$), tissue ($n = 107$), plasma ($n = 85$), serum ($n = 55$) and cell culture ($n = 69$) were tested. Tissue samples comprised of undefined tissues samples ($n = 50$), tonsils ($n = 38$), pancreas ($n = 12$) and lymph nodes ($n = 7$).

Positive reference samples were kindly provided by the Tiermedizinische Hochschule Hannover (TiHo Hannover, EURL for CSFV), the Friedrich-Loeffler-Institut (German NRL for CSFV) and PIWet (Polish NRL for CSFV). The CSFV-positive samples comprised of eight different genotypes (1.1, 1.2, 1.3, 2.1, 2.2, 2.3, 3.1, and 3.4).

Negative reference samples from domestic pigs and wild boars were kindly provided by the Friedrich-Loeffler-Institut and PIWet. Additionally, CSFV-negative blood and blood swab samples were acquired from two German state veterinary laboratories or were collected in-house during different hunting seasons. The latter samples were processed using the IndiMag Pathogen Kit (INDICAL Bioscience) according to the manufacturer's instructions and subsequently tested using the virotype CSFV 2.0 RT-PCR Kit as well as the reference Kit A.

Results/ Conclusion

The summary is shown in Table 12. Individual results are shown in Figure 4.

All positive CSFV samples, even the very weak positives, were detected correctly with the virotype CSFV 2.0 RT-PCR Kit. The virotype CSFV 2.0 RT-PCR Kit demonstrated a diagnostic sensitivity of 100 %, a diagnostic specificity of 100 % and a diagnostic efficiency of 100 %. In this study the virotype CSFV 2.0 RT-PCR Kit demonstrated an overall higher sensitivity (with lower C_T values) compared to in-house reference methods.

Table 12. Diagnostic sensitivity, specificity and efficiency of the virotype CSFV 2.0 RT-PCR Kit.

virotype CSFV 2.0 RT-PCR Kit (CSFV; FAM)		Comparative data			
Total	472	Reference positive	296	Reference negative	176
positive	296	true-positive	296	false-positive	0
negative	176	false-negative	0	true-negative	176

Diagnostic sensitivity: 100 %

Diagnostic specificity: 100 %

Diagnostic efficiency: 100 %

**Comparison virotype CSFV 2.0 RT-PCR Kit vs.
virotype CSFV RT-PCR Kit
(CSFV; FAM; $n = 431$)**




Figure 4. Comparison of C_T values from CSFV-positive and negative samples tested with the virotype CSFV 2.0 RT-PCR Kit compared to the virotype CSFV RT-PCR Kit. All samples situated above the black diagonal line showed lower C_T values with the virotype CSFV 2.0 RT-PCR Kit than tested with the virotype CSFV RT-PCR Kit. Please note that samples that were tested in a 1:10/ 1:100 dilution (virotype CSFV 2.0 RT-PCR Kit) and values for the virotype CSFV RT-PCR Kit (undiluted samples) were omitted from this graph.

4.4 Repeatability

The same sample panel comprising five CSFV-positive RNA samples (samples 1-3, 5 and 6), one CSFV-negative sample (sample 3) and the controls (NC, PC) of the test kit was used for assessment of intra-assay variance, inter-assay variance, batch-to-batch variance, stability testing and for comparison of test results obtained with the virotype CSFV 2.0 RT-PCR Kit tested on different real-time PCR thermocyclers (Table 13).

Table 13. Sample panel for assessment of intra-assay variance, inter-assay variance, batch-to-batch variance, stability and for comparison of test results obtained with the virotype CSFV 2.0 RT-PCR Kit on different real-time PCR thermocyclers. Samples 1-3 and 5-6 were diluted in extracted pig blood RNA to obtain signals for the endogenous Internal Control (JOE).

Sample	CSFV status	Origin	Material	Strain (genotype)	Dilution
1	positive	PIWet	Cell culture	Alfort 187 (gt 1.1)	1:10
2	positive	PIWet	Cell culture	Brescia (gt 1.1)	1:10
3	positive	EURL	Pig blood	CSF385 (gt 2.3)	1:100
4	negative	INDICAL	Wild boar blood	-	-
5	positive	EURL	Pig serum	CSF1060 (gt 2.2)	1:10
6	positive	EURL	Pig plasma	CSF309 (gt 3.4)	1:10
NC	NC	INDICAL		Negative Control (virotype CSFV 2.0 RT-PCR Kit, alpha test)	
PC	PC	INDICAL		Positive Control (virotype CSFV 2.0 RT-PCR Kit, alpha test)	

NC = Negative Control, PC = Positive Control, PIWet = PAŃSTWOWY INSTYTUT WETERYNARYJNY, EURL = European Reference Laboratory, FLI = Friedrich-Loeffler-Institut

4.4.1 Intra-assay variance

The sample panel listed in Table 13 was tested in a sevenfold setup in the same PCR run with the virotype CSFV 2.0 RT-PCR Kit (batch Valid-5) on the Bio-Rad CFX96 instrument using the virotype CSFV 2.0 RT-PCR protocol.

Results / Conclusion

The intra-assay variance is on average 0.72 % for CSFV (FAM) and 0.89 % for the endogenous Internal Control (JOE) (Table 14, Table 15, and Figure 5).

Table 14. Intra-assay variance for **CSFV** (FAM) for the virotype CSFV 2.0 RT-PCR Kit using the Bio-Rad CFX96 instrument.

Intra-assay variance for CSFV (FAM)											
Sample	CSFV status	Reactions (C_T values)							Mean	SD	CV%
		1	2	3	4	5	6	7			
1	pos	30.52	30.64	31.02	30.75	30.43	30.47	30.72	30.65	0.205	0.669
2	pos	27.06	27.20	27.50	27.21	27.10	27.07	27.24	27.20	0.153	0.563
3	pos	33.78	33.24	33.24	33.10	33.65	33.26	33.36	33.38	0.244	0.731
4	neg	-	-	-	-	-	-	-	-	-	-
5	pos	30.31	30.02	30.27	30.34	30.44	30.37	30.36	30.30	0.134	0.443
6	pos	30.76	30.41	31.03	30.37	30.63	30.85	30.78	30.69	0.237	0.773
NC	neg	-	-	-	-	-	-	-	-	-	-
PC	pos	30.31	29.77	29.74	29.48	30.17	30.27	29.63	29.91	0.333	1.112
Mean									0.72		

NC = Negative Control, PC = Positive Control, neg = negative, pos = positive, SD = standard deviation, CV = coefficient of variation, - = no C_T

Table 15. Intra-assay variance for the **endogenous Internal Control (JOE)** for the virotype CSFV 2.0 RT-PCR Kit using the Bio-Rad CFX96 instrument.

Intra-assay variance for the endogenous Internal Control (JOE)											
Sample	CSFV status	Reactions (C_T values)							Mean	SD	CV%
		1	2	3	4	5	6	7			
1	pos	22.19	22.06	22.14	22.17	22.46	22.17	22.25	22.20	0.127	0.572
2	pos	21.86	21.94	22.28	22.28	22.00	21.95	22.03	22.05	0.167	0.757
3	pos	22.78	22.76	23.04	22.62	22.58	22.97	22.75	22.79	0.169	0.741
4	neg	20.88	20.59	20.54	20.47	20.34	20.44	19.99	20.47	0.268	1.307
5	pos	22.88	22.56	23.03	23.04	22.82	23.06	22.45	22.83	0.244	1.068
6	pos	22.59	22.70	23.02	22.68	22.87	22.97	23.08	22.85	0.188	0.824
NC	neg	-	-	-	-	-	-	-	-	-	-
PC	pos	24.68	24.67	24.26	24.05	24.45	24.57	24.48	24.45	0.230	0.941
Mean										0.89	

NC = Negative Control, PC = Positive Control, neg = negative, pos = positive, SD = standard deviation, CV = coefficient of variation, - = no C_T

Intra-assay variance virotype CSFV 2.0 RT-PCR Kit




Figure 5. Boxplots of intra-assay variance for **CSFV (FAM)** and the **endogenous Internal Control (JOE)** for the virotype CSFV 2.0 RT-PCR Kit tested using the Bio-Rad CFX96 instrument.

4.4.2 Inter-assay variance

The sample panel listed in Table 13 was tested in seven independent PCR runs using the virotype CSFV 2.0 RT-PCR Kit (batch Valid-5) on the Bio-Rad CFX96 instrument using the virotype CSFV 2.0 RT-PCR protocol.

Results / Conclusion

The inter-assay variance is on average 0.82 % for CSFV (FAM) and 0.92 % for the endogenous Internal Control (JOE) (Table 16, Table 17, and Figure 6).

Table 16. Inter-assay variance for **CSFV** (FAM) for the virotype CSFV 2.0 RT-PCR Kit using the Bio-Rad CFX96 instrument.

Inter-assay variance for CSFV (FAM)											
Sample	CSFV status	RT-PCR runs (C_T values)							Mean	SD	CV%
		1	2	3	4	5	6	7			
1	pos	30.72	30.59	31.09	30.58	30.63	30.72	30.98	30.76	0.198	0.644
2	pos	27.24	27.05	27.51	27.25	27.34	27.16	27.51	27.29	0.172	0.630
3	pos	33.36	33.41	34.24	33.24	34.03	33.65	34.01	33.71	0.389	1.154
4	neg	-	-	-	-	-	-	-	-	-	-
5	pos	30.36	29.81	30.40	30.15	30.38	30.06	30.48	30.23	0.240	0.794
6	pos	30.78	30.47	30.64	30.42	30.43	30.29	30.85	30.55	0.207	0.677
NC	neg	-	-	-	-	-	-	-	-	-	-
PC	pos	29.63	29.70	29.34	29.62	29.61	29.08	30.03	29.57	0.296	1.002
Mean									0.82		

NC = Negative Control, PC = Positive Control, neg = negative, pos = positive, SD = standard deviation, CV = coefficient of variation, - = no C_T

Table 17. Inter-assay variance for the **endogenous Internal Control (JOE)** for the virotype CSFV 2.0 RT-PCR Kit using the Bio-Rad CFX96 instrument.

Inter-assay variance for the endogenous Internal Control (JOE)											
Sample	CSFV status	RT-PCR runs (C_T values)							Mean	SD	CV%
		1	2	3	4	5	6	7			
1	pos	22.25	22.36	22.35	21.86	22.21	22.65	22.41	22.30	0.240	1.076
2	pos	22.03	22.09	22.46	22.25	22.67	22.72	22.59	22.40	0.280	1.250
3	pos	22.75	23.10	23.00	23.09	22.89	22.94	23.12	22.98	0.136	0.590
4	neg	19.99	20.15	19.85	20.19	19.98	19.94	20.03	20.02	0.118	0.589
5	pos	22.45	22.28	22.19	22.25	22.91	22.57	22.28	22.42	0.253	1.130
6	pos	23.08	23.30	23.22	22.97	22.71	22.86	23.02	23.02	0.203	0.881
NC	neg	-	-	-	-	-	-	-	-	-	-
PC	pos	24.48	24.21	24.87	24.67	24.74	24.80	24.74	24.64	0.229	0.927
Mean										0.92	

NC = Negative Control, PC = Positive Control, neg = negative, pos = positive, SD = standard deviation, CV = coefficient of variation, - = no C_T

Inter-assay variance virotype CSFV 2.0 RT-PCR Kit




Figure 6. Boxplots of inter-assay variance for **CSFV (FAM)** and the **endogenous Internal Control (JOE)** for the virotype CSFV 2.0 RT-PCR Kit tested using the Bio-Rad CFX96 instrument.

4.4.3 Batch-to-batch comparison

The sample panel listed in Table 13 was tested in the same PCR run using three different batches of the CSFV 2.0 RT-PCR Kit (batch 1 = Valid-4; batch 2 = Valid- 5; batch 3 = Valid-6) on the Bio-Rad CFX96 instrument using the virotype CSFV 2.0 RT-PCR protocol.

Results / Conclusion

The batch-to-batch performance showed on average variance of 0.82 % for CSFV (FAM) and 1.00 % for the endogenous Internal Control (JOE) (Table 18, Table 19, and Figure 7).

Table 18. Batch-to-batch variance for **CSFV** (FAM) for the virotype CSFV 2.0 RT-PCR Kit using the Bio-Rad CFX96 instrument.

Sample	CSFV status	Batch number (C _T values)			Mean	SD	CV%
		1	2	3			
		1	pos	30.97			
2	pos	27.42	27.27	27.41	27.36	0.086	0.314
3	pos	33.60	33.07	33.59	33.42	0.300	0.898
4	neg	-	-	-	-	-	-
5	pos	30.82	30.46	30.73	30.67	0.188	0.613
6	pos	30.91	30.81	30.75	30.82	0.082	0.266
NC	neg	-	-	-	-	-	-
PC	pos	29.58	29.31	29.06	29.32	0.263	0.896
Mean					0.82		

NC = Negative Control, PC = Positive Control, neg = negative, pos = positive, SD = standard deviation, CV = coefficient of variation, - = no C_T

Table 19. Batch-to-batch variance for the **endogenous Internal Control** (JOE) for the virotype CSFV 2.0 RT-PCR Kit using the Bio-Rad CFX96 instrument.

Batch-to-batch variance for the endogenous Internal Control (JOE)							
Sample	CSFV status	Batch number (C_T values)			Mean	SD	CV%
		1	2	3			
1	pos	23.40	23.17	22.90	23.15	0.250	1.079
2	pos	22.31	22.43	22.01	22.25	0.215	0.966
3	pos	23.23	23.10	22.33	22.89	0.487	2.126
4	neg	20.85	20.63	20.83	20.77	0.124	0.597
5	pos	23.34	23.49	23.74	23.52	0.198	0.843
6	pos	23.79	23.69	23.56	23.68	0.118	0.496
NC	neg	-	-	-	-	-	-
PC	pos	25.31	25.18	25.63	25.37	0.232	0.914
Mean							1.00

NC = Negative Control, PC = Positive Control, neg = negative, pos = positive, SD = standard deviation, CV = coefficient of variation, - = no C_T

Batch-to-batch variance virotype CSFV 2.0 RT-PCR Kit




Figure 7. Boxplots of batch-to-batch variance for CSFV (FAM) and the endogenous Internal Control (JOE) for the virotype CSFV 2.0 RT-PCR Kit tested using the Bio-Rad CFX96 instrument.

4.4.4 Comparison of real-time PCR thermocyclers

The virotype CSFV 2.0 RT-PCR Kit can be used on different standard real-time PCR cyclers. Table 20 gives an overview of selected PCR cyclers and their approximate run times, using the virotype CSFV 2.0 RT-PCR protocol.

Note: The use of the RT-PCR Kit is not limited to the mentioned instruments.

Table 20. Selected overview of real-time thermocyclers and their approximate run times for the virotype CSFV 2.0 RT-PCR protocol.

Thermocycler				
	Manufacturer	Model	Filters	Run time [minutes]
A	Bio-Rad Laboratories, Inc., Hercules, California, USA	CFX96	FAM, HEX	72
B	Agilent Technologies, Santa Clara, California, USA	AriaMx	FAM, HEX	59
C	Agilent Technologies, Santa Clara, California, USA	Stratagene Mx3005P	FAM, HEX	66
D	Thermo Fisher Scientific Inc., Waltham, Massachusetts, USA	Applied Biosystems™ 7500 Fast ¹	FAM, JOE	75
E	Thermo Fisher Scientific Inc., Waltham, Massachusetts, USA	QuantStudio® 5	FAM, VIC	65
F	QIAGEN® GmbH, Hilden, Germany	Rotor-Gene® Q 5	Green, Yellow	85
<i>Theoretical run time of the virotype CSFV 2.0 RT-PCR Kit:</i>				36

¹ Standard Mode setting

The sample panel listed in Table 13 was tested with the virotype CSFV 2.0 RT-PCR Kit (batch Valid-5) on six different real-time PCR thermocycler instruments (named A-F, see Table 20) using the virotype CSFV 2.0 RT-PCR protocol. For the Thermo Fisher ABI 7500 Fast instrument, the Standard Mode with ROX as passive reference dye was performed.

Results / Conclusion

The results are summarized in Table 21 (CSFV pathogen/ FAM channel) and Table 22 (endogenous Internal Control/ JOE channel). All samples tested on different real-time PCR thermocycler instruments showed comparable results.

Based on this data set, INDICAL can recommend using the virotype CSFV 2.0 RT-PCR Kit on the Bio-Rad CFX96, Agilent Technologies AriaMx, Agilent Technologies Stratagene Mx3005P, Thermo Fisher ABI 7500 Fast, Thermo Fisher QuantStudio 5 and QIAGEN Rotor-Gene instruments.

Table 21. Inter-thermocycler variance for the pathogen **CSFV** (FAM).

Inter-thermocycler variance for CSFV (FAM)								
Sample	Material	CSFV status	Thermocycler (C _T values)					
			A	B	C	D	E	F
1	Cell culture	pos	30.26	30.79	31.41	31.72	31.65	29.10
2	Cell culture	pos	27.09	26.91	27.92	27.15	27.72	25.85
3	Blood	pos	32.77	34.14	34.49	34.05	34.33	33.29
4	Blood	neg	-	-	-	-	-	-
5	Serum	pos	29.77	29.99	30.79	31.97	30.44	28.63
6	Plasma	pos	30.10	30.03	30.62	30.74	30.33	28.79
NC	NC	neg	-	-	-	-	-	-
PC	PC	pos	29.47	28.97	29.50	29.01	28.30	27.93

NC = Negative Control, PC = Positive Control, neg = negative, pos = positive, SD = standard deviation, CV = coefficient of variation, - = no C_T

Table 22. Inter-thermocycler variance for the **endogenous Internal Control (JOE)**.

Inter-thermocycler variance for the endogenous Internal Control (JOE)								
Sample	Material	CSFV status	Thermocycler (C_T values)					
			A	B	C	D	E	F
1	Cell culture	pos	22.63	23.26	23.36	23.87	23.46	22.25
2	Cell culture	pos	22.79	22.91	23.13	24.00	23.91	22.60
3	Blood	pos	23.20	22.96	22.66	23.85	23.40	22.19
4	Blood	neg	21.15	19.76	20.15	20.54	19.59	19.36
5	Serum	pos	23.28	23.17	23.86	24.36	24.10	22.79
6	Plasma	pos	23.30	23.00	23.11	24.05	23.62	22.56
NC	NC	neg	-	-	-	-	-	-
PC	PC	pos	25.21	25.27	24.31	25.40	23.01	21.74

NC = Negative Control, PC = Positive Control, neg = negative, pos = positive, SD = standard deviation, CV = coefficient of variation, - = no C_T

4.5 Stability testing

4.5.1 Freeze-thaw-cycles

The sample panel listed in Table 13 was used for validating the stability of the virotype CSFV 2.0 RT-PCR Kit. One kit batch (batch: Valid-6) was tested at the time of production and after six freeze/thaw cycles. The mean value (Mean), standard deviation (SD) and coefficient of variation (CV) were calculated.

Results / Conclusion

The virotype CSFV 2.0 RT-PCR Kit shows excellent stability with an average variance of 0.82 % for CSFV (FAM) and 0.68 % for the endogenous Internal Control (JOE) (Table 23, Table 24).

Table 23. Stability testing for **CSFV** (FAM) of the virotype CSFV 2.0 RT-PCR Kit using the Bio-Rad CFX96 instrument.

Sample	CSFV status	Stability for CSFV (FAM)		Mean	SD	CV%
		Freeze-thaw-cycle (C _T values)				
		1	6			
1	pos	30.95	30.41	30.68	0.378	1.232
2	pos	27.27	27.87	27.57	0.429	1.554
3	pos	33.07	33.21	33.14	0.095	0.288
4	neg	-	-	-	-	-
5	pos	30.46	30.09	30.27	0.264	0.872
6	pos	30.81	31.18	30.99	0.261	0.843
NC	neg	-	-	-	-	-
PC	pos	29.31	29.26	29.29	0.034	0.116
Mean						0.82

NC = Negative Control, PC = Positive Control, neg = negative, pos = positive, SD = standard deviation, CV = coefficient of variation, - = no C_T

Table 24. Stability testing for the **endogenous Internal Control (JOE)** of the virotype CSFV 2.0 RT-PCR Kit using the Bio-Rad CFX96 instrument.

Sample	CSFV status	Stability for Internal Control (JOE)		Mean	SD	CV%
		Freeze-thaw-cycle (C _T values)				
		1	6			
1	pos	23.17	23.01	23.09	0.109	0.471
2	pos	22.43	22.36	22.39	0.050	0.224
3	pos	23.10	22.99	23.04	0.083	0.360
4	neg	20.63	21.05	20.84	0.301	1.444
5	pos	23.49	23.24	23.37	0.174	0.743
6	pos	23.69	23.96	23.82	0.188	0.790
NC	neg	-	-	-	-	-
PC	pos	25.18	25.43	25.31	0.178	0.704
Mean						0.68

NC = Negative Control, PC = Positive Control, neg = negative, pos = positive, SD = standard deviation, CV = coefficient of variation, - = no C_T

4.5.2 Heparin inhibition

To test the stability of the virotype CSFV 2.0 RT-PCR Kit, sample inhibition was simulated by treating the sample FLI-2 (Israel, genotype 2.1, 1: 100 dilution) with an increasing concentration of heparin (0.17 – 3.40 U/reaction). The sample was tested in duplicate using the virotype CSFV 2.0 RT-PCR Kit (batch: Valid-6) and the Kit D.

Results / Conclusion

Whilst complete inhibition of the sample by heparin was observed at 1.70 U/reaction with the competitor Kit D, the CSFV signal (FAM) was fully inhibited at a higher heparin concentration of 3.40 U/reaction when using the virotype CSFV 2.0 PCR Kit (Table 25). The inhibition study using heparin showed that the virotype CSFV 2.0 RT-PCR Kit shows better resilience to inhibition by heparin than the competitor Kit D.

Table 25. Stability testing of the virotype CSFV 2.0 RT-PCR Kit for inhibition by heparin.

Stability (heparin inhibition)		
Heparin [U/reaction]	CSFV (FAM) (C_T value)	
	virotype CSFV 2.0 RT-PCR Kit	Kit D
-	31.12	31.32
-	30.44	31.27
0.17	30.44	32.03
0.17	30.43	32.04
0.34	30.30	34.21
0.34	30.50	34.24
1.70	35.04	-
1.70	34.57	-
3.40	-	-
3.40	-	-

- = no C_T

4.5.3 EDTA inhibition

To test the stability of the virotype CSFV 2.0 RT-PCR Kit, sample inhibition was simulated by treating the sample FLI-2 (Israel, genotype 2.1, 1: 100 dilution) with an increasing concentration of Ethylenediaminetetraacetic acid (EDTA; 0.5 – 6.0 mM final concentration in reaction mix). The sample was tested in duplicate using the virotype CSFV 2.0 RT-PCR Kit (batch: Valid-6) and the Kit D.

Results / Conclusion

The inhibition study using EDTA showed that complete sample inhibition (CSFV, FAM) occurs at high concentrations only (> 5.0 mM) when using the virotype CSFV 2.0 RT-PCR Kit (Table 26).

However, using the competitor Kit D, inhibition of the pathogen signal (CSFV, FAM) already occurred at concentration > 2.5 mM EDTA. Thus, the virotype CSFV 2.0 RT-PCR Kit shows better resilience to inhibition by EDTA than the competitor Kit D.

Table 26. Stability testing of the virotype CSFV 2.0 RT-PCR Kit for inhibition by EDTA.

Stability (EDTA inhibition)			
EDTA [mM]	CSFV (FAM) (C_T value)		
	virotype CSFV 2.0 RT-PCR Kit		Kit D
-	30.08		31.68
-	30.02		31.37
0.5	30.05		n.d.
0.5	30.11		n.d.
1.0	30.17		n.d.
1.0	30.41		n.d.
1.5	30.60		n.d.
1.5	30.51		n.d.
2.0	31.30		35.49
2.0	30.85		-
2.5	31.28		38.48
2.5	30.70		-
3.0	32.11		-
3.0	31.97		-
3.5	32.43		-
3.5	32.84		-
4.0	33.41		-
4.0	33.58		-
5.0	34.18		-
5.0	-		-
6.0	-		-
6.0	-		-

n.d. = not done, - = no C_T