# SPACER TRIALS

TRIALS FOR SPACER FOR HIP,  
TRIALS FOR SPACER FOR KNEE,  
TRIALS FOR SPACER FOR KNEE ATS  
COMPONENT, TRIALS FOR SPACER  
FOR SHOULDER

# SPACER TRIALS

## TRIALS FOR SPACER FOR HIP, TRIALS FOR SPACER FOR KNEE, TRIALS FOR SPACER FOR KNEE ATS COMPONENT, TRIALS FOR SPACER FOR SHOULDER

Spacer Trials are surgically invasive medical devices for transient use to be utilised in orthopaedic surgery. Spacer Trials are trial devices for the choice of the right size of the corresponding hip, knee and shoulder Tecres Spacers.

The devices are re-usable; prior to use, the devices must be cleaned and sterilized following a validated cycle.

### INTENDED USERS

The devices are intended for use by health care professionals:

- Nurses and surgical team members: for the preparation of the medical devices (cleaning, sterilization) and their disposal;
- Orthopaedic Surgeons: for the use of the medical devices (insertion and explantation in the joint).

### PATIENT TARGET GROUP(S)

The devices are intended for use in skeletally mature patients suffering from a peri-prosthetic joint infection (PJI) of the joint.

### CLINICAL BENEFIT & PERFORMANCES

Spacer Trials help to:

- select the proper spacer size
- prepare the host site for spacer implantation
- evaluate spacer stability

### WARNINGS

The devices must be used only in the short term in order to determine the right size of the related Spacer to be implanted. The trial devices must not be implanted.

Do not use if the device appears to be damaged (deformations, colour changes, presence of bubbles). The trial must be replaced if it shows damage following a visual inspection.

### CLEANING & STERILIZATION

The devices are supplied non-sterile.

The user must perform cleaning (manual or automated) and sterilization of the devices prior to their use. Only the processes described in this document have been validated by the manufacturer.

Equivalent cleaning and sterilisation programs may be used, provided that the users have checked that they have been verified / validated.

#### Manual Cleaning

##### Pre-cleaning steps

1. Rinse the device using running tap water for 3 minutes.
2. Fully immerse the device in 300 mL of a freshly prepared solution of Neodisher MediClean Forte (or equivalent) at a concentration of 0.5%.
3. Leave the device to soak for 10 minutes at room temperature.
4. During soaking time, use a soft bristled brush to gently clean the device for 1 minute. Scrub the device with small circular motions.
5. Remove the device from the detergent solution and rinse all the surfaces with Ultra Purified Water (UPW) for 3 minutes.
6. Inspect the device under normal lighting conditions to determine if all visible adherent contaminants have been removed from the surfaces.

##### Manual cleaning steps

1. After the device has been pre-cleaned and visually inspected, place the device in an ultrasonic bath with 300 mL of a freshly prepared solution of Neodisher MediClean Forte (or equivalent) at a concentration of 0.5%, for 15 minutes at room temperature.
2. During sonication, use a soft-bristled brush to gently clean the device for 1 minute. Scrub the device with small circular motions.
3. Remove the device from the ultrasonic bath and rinse all surfaces with UPW for 3 minutes.
4. Dry all the surfaces by placing the device at 105°C for 15 minutes.
5. Inspect the device under normal lighting conditions to determine if all visible adherent contaminants have been removed from the surfaces. Check that the device marking is still readable.

### Automated Cleaning

#### Pre-cleaning steps

1. Rinse the device using running tap water for 3 minutes
2. Fully immerse the device in 300 mL of a freshly prepared solution of Neodisher MediClean Forte (or equivalent) at a concentration of 0.5%.
3. Leave the device to soak for 10 minutes at room temperature.
4. During soaking time, use a soft-bristled brush to gently clean the device for 1 minute. Scrub the device with small circular motions.
5. Remove the device from the detergent solution and rinse all the surfaces with Ultra Purified Water (UPW) for 3 minutes.
6. Inspect the device under normal lighting conditions to determine if all visible adherent contaminants have been removed from the surfaces.

#### Automated cleaning steps

1. Place the device in a metallic basket.
2. Place the metallic basket in the washer-disinfector.
3. Clean the device using the program illustrated in table 1.
4. Inspect the device under normal lighting conditions to determine if all adherent visible contaminants have been removed from the surfaces. Check that the device marking is still readable.

PHASE	TEMPERATURE	TIME	DETERGENT
Prewash	Cold tap water	2 min	-
Wash 1	Warm tap water (50°C)	3 min	Neodisher mediclean Forte (0.2%)
Wash 2	Warm tap water (50°C)	2 min	Neodisher mediclean Forte (0.2%)
Rinse	Warm PW (50°C)	2 min	-
Drying	105°C	15 min	-

Table 1. Automated Cleaning Program

### Sterilization

1. Pack the device in a double steam sterilization pouch or in a packaging which performs equivalently when subjected to a steam sterilization procedure.
2. Perform the steam sterilization using one of the following cycles.

	CYCLE 1	CYCLE 2
Sterilizer Type:	Pre-vacuum	Pre-vacuum
Preconditioning pulse:	4	4
Temperature:	132°C	121°C
Exposure time:	4 minutes	15 minutes
Dry time:	30 minutes	30 minutes

Table 2. Sterilization Cycle Parameters

Sterile, packaged devices should be stored in a designated, limited access area that is well ventilated and provides protection from dust, moisture, insects, vermin, and temperature/humidity extremes.

### PREOPERATIVE WARNING

Before use, check for possible visible damage on the sterile packaging: if damage is observed, the devices must be packaged and sterilized again.

Before use, check for possible visible damage to the devices: if damage is observed, the devices must be replaced.

### SPACER TRIAL PLACEMENT

Spacer Trials act like a probe as their geometry and shape mimic the implantable Spacers.

The probe device is implanted in the joint so that the correct spacer size, which better fits and which provides the best stability, can be chosen. At the same time it aids in guiding the preparation of the host site prior

to implant of the spacer device.

Before performing any surgical operation, the surgeon must know the surgical procedure, the devices and their characteristics.

- Choose the appropriate trial on the basis of the dimensions of the removed prosthesis.
- Position the trial in the joint space to determine if it is anatomically correct and stable.
- Reduce the joint.

**Note:** when using the Spacer-K Trial, be aware that the knee joint should not be too tight during trialing as it will tighten further upon cement fixation. Tightness may be relieved through downsizing and/or recontouring the bone to achieve a satisfactory fit.

- Once the trial has been tested and verified, remove it from the patient, and use it as a reference when selecting the implantable Spacer.

**Caution only for the Spacer-K ATS Trial Kit:** Spacer-K ATS Trial Kit must be used in combination with the Spacer-K Trial and Spacer-K XL Trial. The allowed combinations between Spacer-K ATS Trial Kit and Spacer-K Trial/Spacer-K XL Trial are reported in the table below:

<b>SPACER-K ATS TRIAL KIT (REF CODE: SPK90Z0)</b>	<b>COMBINATION WITH SPACER-K TRIAL (REF CODE: SPK03) AND SPACER-K XL TRIAL (REF CODE: SPK03Z0)</b>
Spacer-K ATS Trial 60/07	Spacer-K Trial (54mm-60mm) or Spacer-K Trial (64mm-70mm)
Spacer-K ATS Trial 60/12	
Spacer-K ATS Trial 80/07	Spacer-K Trial (74mm-80mm) or Spacer-K XL Trial (84mm-90mm)
Spacer-K ATS Trial 80/12	

Table 3. Combination allowed for Spacer-K ATS Trial Kit

#### DISPOSAL

Repeated re-processing operations of the Spacer Trials according to the instructions included in this document have a minimal impact in terms of useful lifetime. Internal laboratory testing has shown that the Spacer Trials are able to withstand at least 50 consecutive sterilization cycles without suffering damage and without any dimensional change.

The useful lifetime of the Spacer Trials is mainly affected by wear and damage caused by surgical use. Spacer Trials need to be replaced if they appear damaged after a visual inspection.

The devices must be disposed of as biological waste in accordance with local waste regulations, as the devices might have been contaminated with potentially infectious substances of human origin.

#### SERIOUS INCIDENT REPORTING

Any serious incident potentially attributable to Spacer Trials should be reported to the manufacturer, TECRES S.p.A., and to the competent authority in which the user and/or the patient is established.

The manufacturer can be reached at the following contacts:

- email: [info@tecres.it](mailto:info@tecres.it)
- tel.: +39 045 9217311

The contact information of the competent authority located in the European Member State in which the patient is established is reported on European Commission website, under the section "medical devices" at the following link: [https://ec.europa.eu/health/md\\_sector/contact\\_en](https://ec.europa.eu/health/md_sector/contact_en)

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TRIALS FOR SPACER FOR HIP, TRIALS FOR SPACER FOR KNEE, TRIALS FOR SPACER FOR KNEE ATS COMPONENT, TRIALS FOR SPACER FOR SHOULDER

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