Integrated Registration

Multi-modality image management

In radiology practice today, multiple complementary imaging modalities are frequently relied on to image the same patient, as each modality provides unique anatomical and physiological information. Three dimensional (3D) modalities such as computed tomography (CT) and magnetic resonance imaging (MRI) provide information about the relationship of anomalous findings to surrounding structures and tissues.

The Integrated Registration application provides you with a streamlined workflow to combine and display images from any two of the five major modalities (CT, MRI, positron emissions tomography (PET), single-photon emission computed tomography (SPECT) and X-ray angiography (XA)—aimed at improving your ability to provide more clinically relevant diagnostic information and improving surgical and treatment planning.


CONFIDENTIAL

SECRET DE
AFACERI

Overview

- Integrated Registration provides you with the capability to align and fuse two volumetric acquisitions from either the same or different acquisition modalities. With it, you can easily compare 3D anatomical images from CT, MRI with PET, SPECT and XA¹ for a comprehensive analysis.
- Integrated Registration is available on VolumeShare 7, a multi-modality advanced visualization workflow solution that helps to enhance diagnostic precision and productivity.

Highlights

- Load and register on the fly with drag and drop
- Easy registration using automatic, manual and landmark methods
- Multiple 2D and 3D fusion capabilities
- Automatic registration propagation across series and one-click access to pre-defined protocols
- Rigid and deformable registration
- Drawing and saving of contours as RTSS DICOM objects
- Access to regional registration on the fly for an easy review of challenging areas
- Compatible with AdvantageSim™ MD and other standard treatment planning software

Features

- ✓ Allows you to load exams and series from CT, MRI, PET, SPECT or XA together for registration in the same session
- ✓ Offers you a choice of four automatic registration methods
- ✓ Allows you to customize screen layouts to suit your review needs and provides comparison of a current exam with previous exams to let you evaluate disease progression and treatment efficacy
- ✓ Combines functional and anatomical images from different modalities, providing you with enhanced perspective of the area of interest
- ✓ Allows you to define one or more contours around relevant anatomy and save them as RTSS DICOM objects
- ✓ Provides you the ability to save registered data as a new DICOM series or as a registration DICOM object (does not include SPECT)
- ✓ Define volume of interest in one model for automatic report to another model for subsequent use in Radiation Therapy (RT) planning²

CONFIDENTIAL

SECRET DE
AFACERI

Feature details

Loading

Integrated Registration allows you to load multiple exams and series from CT, MRI, PET, SPECT or XA for registration in the same session. Within an open registration protocol session, you can easily add additional series using a simple drag and drop method. You can also choose to have the registration protocol run automatically at loading or dynamic loading. Reference and Moving series are preselected; however, you can easily control and adjust these selections.

Registration

With Integrated Registration, you can employ four rigid registration methods:

- Automatic registration based on “mutual information” and “matching boundaries” for CT, MR, PET or XA. (SPECT can only be registered using manual and landmark registration)
- Automatic regional registration, which can be run by defining a region of interest (ROI) upon which the registration methods (listed above) are applied only within the designated boundaries
- Interactive manual registration, which, in real time, conveys translation and rotation made to the alignment of one data set to the second data set. This method can be used with all modality pair combinations.
- Landmark registration, which is based on manual definition of common points between the two data sets. This method can be used with all modality pair combinations.

You can combine these four methods to achieve optimal results. The methods allow registration of volumes with linear FOV rescaling and 6 degrees of freedom (translation and rotation).

Non-rigid automatic registration can be employed for CT to CT data sets only. The non-rigid transformation of the CT series can then be propagated to any other MR, PET, SPECT, or XA series.

Global registrations are possible for both rigid and non-rigid registrations. Rigid registration also benefits from regionally optimized registration using an ROI.

Multiple series acquired in the same exam (i.e., with the same frame of reference) are grouped together and registration that is applied to one series of the group is automatically propagated to other series of the same group by default. You can manually group or ungroup series if desired. Propagation will then be applied to manually grouped series while ungrouped series may be registered independent of their original group. Use manual ungrouping specifically to independently register series of one identical exam during which the patient has moved.

The Center, Reset and Undo/Redo options improve productivity during the registration process.


Quality control

Two methods are available to validate the accuracy of the registration:

- Qualitatively by visually inspecting registration results using all available display modes (see Visualization below)
- Quantitatively by defining pairs of common landmarks on the two exams. This includes calculation of the Maximal Error on a pair and the Root Mean Square (RMS) error on all pairs

CONFIDENTIAL

SECRET DE
AFACERI

Visualization

Optimized screen layouts are provided by default and you can also customize them to suit your needs.

The common cursor, active by default, provides a triangulated display in which the cursor is positioned at the same geometrical coordinates in all views.

You can apply dedicated display modes to fused views during your review:

Windows	An adjustable viewing window of the registered exam that can be moved anywhere over the reference exam and shows the corresponding part of the moving exam.
Tiles	The selected view is subdivided into four tiles that meet at the position of the 3D cursor. Two tiles display regions of the reference exam; two tiles display corresponding regions of the moving exam.
Fused	The moving exam is interlaced with the reference exam with an adjustable fusion percentage.
Inversed fused	The moving exam is displayed over the reference exam in inverse linear grey levels with an adjustable fusion percentage.
Threshold	The overlay exam is displayed with lower threshold on top of reference exam, all with an adjustable fusion percentage.
Swipe	Interactively splits the image vertically in 2 parts: the overlaid exam on one side and the fused exams on the other side.

You can create 2D, 3D fused, and hybrid 2D/3D fused views with a simple drag and drop method. Save the layouts you define for personalized display of registration results.

Contours

With Integrated Registration you can define one or multiple contours around relevant anatomical structures on the exam.

To help define the area, the contour appears on all exams as you draw it. You can draw the contour slice by slice or on two non-contiguous slices with interpolation on slices in between.

You can modify or remove contours from a single slice or from the entire exam. You can also save them in the patient database, linked with the corresponding exams in DICOM RT Structure Set format for analysis in AdvantageSim MD* or other DICOM RT-compatible planning systems. On AW Workstation, you have direct access to AdvantageSim MD* during the contouring and saving processes.

Integrated Registration protocols also allow you to define, save, and load RTSS objects.

Saving

You can save registered results in several ways:

- As a new reformatted series in the Reference or Moving exam geometry
- As Save State after the registration is completed; display, segmentations, and measurements are also saved
- As a DICOM spatial registration object in the same exam with the reference series in the Patient List

Operation

Registration is possible in two contexts:

- Loading a dedicated registration protocol
 - Load a DICOM 3 CT, MR, PET, SPECT or XA data set and another data set from one of these modalities into the Integrated Registration protocol. Registration is performed based on reference and moving series selection
 - Review the registration quality with visualization tools and validate the result
 - Optionally, you can define and save the contours of structures of interest
 - Save registration results
- Loading multiple volumes into a generic Volume Viewer protocol, after which registration can be performed on demand using Dynamic Load feature. The registration features are available in Visualization toolbar.

CONFIDENTIAL

SECRET DE
AFACERI

System requirements

- Recommended monitor resolution is up to dual 2MP (1600 x 1200) or a single 3MP (1536 x 2048)
- Refer to AW Applications platform compatibility document

Image requirements

Acquisition position and plane, slice resolution, number, thickness and spacing may be different in the acquisition of the two exams to be registered.

Each image set should meet basic requirements, however:

- Field of view, matrix size and display center should be the same for all images inside one series
- Orientation should be the same for all images in the series
- Series should include more than one image
- Tilted acquisitions are not supported by mutual information based automatic algorithms; they can be registered with manual, landmark, or automatic “matching boundaries” registration methods. Datasets used as reference and moving should have some matching anatomical location.

Indications for use

Integrated Registration provides easy means for comparison of 3D images from CT, MRI, PET or SPECT and XA. To help physicians in diagnostic radiology or therapy planning, Integrated Registration allows 3D registration between volumetric acquisitions that may come from the same acquisition modality or from different acquisition modalities.

CONFIDENTIAL

SECRET DE
AFACERI

Summary

Integrated Registration gives you the capability to combine findings from different modality exams into a fused model demonstrating relationships between anatomical structures and physiological processes. The program makes registration and fusion easy and fast with features such as:

- Robust automatic registration that can be run over a matter of only a few seconds
- The ability to perform automatic registration upon loading the exam from the patient list or from Dynamic Load, reducing the need for you to manually launch the registration
- Automatic registration propagation to all series acquired in the same exam (i.e., same frame of reference), speeding the registration process (specifically helpful for MRI exams where multiple series are acquired)
- Ability to load additional series into an active session. Registration applied to already loaded series of the same exam (i.e., same frame of reference) automatically propagates to the additional series
- Interactive manual registration provides flexible registration optimization based on real time visual feedback

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GE HealthCare is a leading global medical technology, pharmaceutical diagnostics, and digital solutions innovator, dedicated to providing integrated solutions, services, and data analytics to make hospitals more efficient, clinicians more effective, therapies more precise, and patients healthier and happier. Serving patients and providers for more than 125 years, GE HealthCare is advancing personalized, connected, and compassionate care, while simplifying the patient's journey across the care pathway. Together our Imaging, Ultrasound, Patient Care Solutions, and Pharmaceutical Diagnostics businesses help improve patient care from diagnosis, to therapy, to monitoring. We are an \$19.6 billion business with approximately 51,000 colleagues working to create a world where healthcare has no limits.

Follow us on [LinkedIn](#), [X](#) (formerly Twitter), and [Insights](#) for the latest news, or visit our website gehealthcare.com for more information.

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SECRET DE
AFACERI

→ Scan the QR Code or visit us at
gehealthcare.ca/en-ca/products/advanced-visualization/all-applications/integrated-registration

Regulatory compliance:

This product complies with the Regulation (EU) 2017/45 of the European Parliament and the Council of medical device (MDR)

Rx Only

This product and its feature and compatible products may not be available in some other countries or regions. Please contact your sales associate.

*Separate medical device requiring an additional license.

References:

1. For XA modality series, Integrated Registration currently supports only 3D X-ray Angiography images (stored as CT Image Storage DICOM objects) acquired with GE HealthCare Innova™ equipment and reconstructed with the Innova3DXR application.
2. Requires DICOM RT capability on RT planning system.

© 2024 GE HealthCare. AdvantageSIM and Innova are trademarks of GE HealthCare. GE is a trademark of General Electric Company used under trademark license.

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Despre GE HealthCare Technologies Inc.

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SECRET DE
AFACERI

→ Scanați codul QR sau vizitați-ne la [gehealthcare.ca/en-ca/products/advanced-visualization/all-applications/integrated-registration](https://www.gehealthcare.ca/en-ca/products/advanced-visualization/all-applications/integrated-registration)

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SECRET DE
AFACERI

GE HealthCare

VesseliQ™ Xpress with AutoBone™ Xpress

CONFIDENTIAL

SECRET DE
AFACERI

Fast and Easy to Use Vessel Analysis with Automated Real-Time Tracking, Auto-Labeling and Simplified User Interface.

Contents

1. Features..... 4

2. System Requirements..... 4

3. Indications for Use..... 4

4. Regulatory Compliance..... 4

5. Features Detail..... 4

5.01 Anatomy-based protocols..... 4

5.02 AutoBone™ Xpress..... 5

5.03 Auto Aorta..... 5

5.04 Fast Tracking..... 5

5.05 Vessel Label Data Base..... 5

5.06 Dynamic AVA..... 5

5.07 Advanced Vessel Analysis Tools..... 6

5.08 Stenosis Detection & ColorID..... 6

5.09 Quick AVA..... 6

5.10 Summary Table & Report..... 7

5.11 Multiple Study Comparison..... 7

5.12 Save State & Protocol..... 7

5.13 Summary..... 7

Peripheral vascular disease (PVD) is a condition similar to coronary artery disease and carotid artery disease. In PVD, fatty deposits build up in the inner linings of the artery walls. These blockages restrict blood circulation, mainly in arteries leading to the kidneys, stomach, arms, legs and feet.

There are different clinical techniques used to diagnose PVD. Computed Tomography Angiography (CTA) can potentially provide sufficient information for therapeutic decision-making and reduced imaging costs for the management of PVD. One of the advantages of computed tomography (CT) is that it shows both vascular & nonvascular anatomy and also provides 3D information for best view. Also, CT may be more useful for patients with limited or no peripheral access or with stents, who may not be candidates for magnetic resonance angiography (MRA). Complementary to the ability to perform CTA is to perform clinical assessment by a dedicated post-processing tool.

Overview

VessellQ™ Xpress with Autobone Xpress is a post-processing software package for the Advantage Workstation (AW) platform, AW Server, CT scanners and PACS reading stations. It is a tool to analyze CT angiographic data including stenosis analysis, thrombus, pre and post stent planning procedures, and directional vessel tortuosity visualization. Autobone Xpress provides zero-click bone segmentation for head, neck and other anatomy. The automation and ease of use streamlines workflow.

Highlights

- Automated detection of aorta and iliacs with auto-labeling of vessels
- Automated real-time fast tracking for all vasculature
- New contextual wheel menu including easy access to advanced options
- New centerline edition capabilities
- Improved predefined measurements and edit capabilities
- New reporting interface which includes measurements in summary table
- New user interface

CONFIDENTIAL

SECRET DE
AFACERI

Visit us:

<https://www.gehealthcare.com/products/advanced-visualization/all-applications/autobone-vesseliq-xpress>

Features

- **Anatomy based protocols for guided workflow in the analysis for your peripheral vascular CTA acquisitions.**
- Automatic bone removal for efficient review of vasculature.
- Automated abdominal aorta CTA protocol with vessel tracking and vessel labeling.
- Fast tracking detects entire length of vessel and displays the centerline based on a single start point.
- Advanced editing tools bringing new flexibility for faster editing based on targeted anatomy.
- Thrombus detection with measurements and comparison over time.
- Dynamic AVA provides access to vessel tracking at any time, from any protocol, on any vessel.
- **Advanced measurement tools for stenosis detection and color identification of HU value for quantification of different densities.**
- Create a queue of cases that are pre-loaded and pre-processed in the background while you read another exam.
- Measurements are organized in the summary table allowing for easy viewing and exporting.
- Fully customizable workflows with adjustable layouts, personalized parameter and settings, custom review steps.
- Save State lets you save and restore the state of the processed images at any stage

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Indications for Use

VessellIQ Xpress is an optional, non-invasive, optimized, post-processing application intended to provide images and tools to analyze vascular anatomy and pathology, aiding physicians in diagnosis and determination of treatment paths, from a set of Computed Tomography (CT) Angiographic images. **VessellIQ Xpress is an option for the Advantage Workstation (AW) platform, CT Scanner, and/or PACS, which can be used in the analysis of 2D and 3D CT Angiography images/data for the purpose of cardiovascular and vascular disease assessment. This software-only device is designed to support physician assessment for a wide variety of clinical uses such as stenosis analysis, pre/post stent planning, pre/post valve replacement planning, and directional vessel tortuosity visualization. VessellIQ Xpress' automatic visualization tools provide users the capability to segment bony structures for accurate identification of the vessels.**

Indications for Use *(cont.)*

Additional tools enable analysis of the vascular anatomy including the aorta, valves, and branching vessels for: anatomical sizing; density and volume analysis of segmented vasculature and calcified / non-calcified plaque; and measurements of abnormalities. The TAVI Analysis option for VessellIQ Xpress is a planning tool used for Trans Aortic Valve Implantation (TAVI) procedures. It automatically segments the aorta and displays the aortic valve in multiple views for measurements of anatomic structures commonly needed for aortic annulus replacement planning. TAVI Analysis provides guided workflow and semi-automated tools to aid in evaluation of appropriate access pathways for interventional procedure planning.

AutoBone Xpress option is a software package that is intended to facilitate segmentation of bony structures and calcifications from abdominal and extremity CT Angiography data.

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Features Detail

Anatomy-based protocols

- **The program offers you predefined protocols for fast, easy review of all vascular structures:**
 - Intra-cerebral
 - Carotids
 - Thoracic Aorta
 - Abdominal Aorta
 - Runoff
 - Spinal
 - Extremity
- Protocols may be customized to meet your requirements and saved for quick access.

CONFIDENTIAL

SECRET DE
AFACERI


Features Detail *(cont.)*

AutoBone™ Xpress

- Anatomical driven protocols with automatic bone removal.
- Improve vessel visualization by removing obstructive bony detail

- Automatic calcium removal for vessel lumen visualization

Auto Aorta


- With a single click, the program automatically removes bones and performs centerline tracking of the aorta and iliacs and labels vessels.

- Simply confirm the automated tracking or modify it by adding additional branch points, deleting a single point, or deleting multiple points from any view.

Fast Tracking


- Semi-automatic detection of complete vessel from a single starting point

- Advanced editing with variable interpolation tool for quick centerline edits adjusted for length of vessel.
- Bulls eye adjustment capability on x-section views to perfect tracking of the smallest most tortuous vessels.

- Volume measurements of the thrombus allow you to track change over time. Additional measurements include :

- Minimum diameter
- Maximum diameter
- Mean diameter
- Cross-sectional area

Vessel Label Data Base

- Each protocol is associated with a database of vessel names
- Customize vessel label database to meet your requirements for every anatomical region.
- A right mouse click provides access to the list of vessel names per anatomy.

Dynamic AVA

- Real-time tracking of the vessel center line and instant visualization of an unfolded view for quick vessel inspection.
- Immediate access to vessel analysis tools.
- Edit contours with a selection of drawing tools with variable interpolation for fast editing.

CONFIDENTIAL

SECRET DE
AFACERI

Features Detail *(cont.)*

Advanced Vessel Analysis Tools


- Semi-automated vessel analysis lets you track, extract, visualize, and size vessels.
- Track a vessel and its branches for quick review in curved, oblique longitudinal or cross-sectional views.
- Multi-point deposition capabilities for improved center line tracking of challenging cases.
- When you edit a vessel contour, the program automatically interpolates the edit to neighboring slices.
- Generate and save rotating movies
 - Curved planar reformation and cross-sectional (lumen) views
 - Best L-section reformatted images
- Select vessel segments for quantitative or qualitative analysis.

- Available measurements include:
 - Distance
 - Cross-sectional area
 - Single or dual reference point comparison
 - Relative percent stenosis / aneurysm
 - Color Identification
 - 3D Angle
 - Thrombus Extraction
 - Volume
 - Mean, min, max diameter
- Move freely between review steps and protocols without leaving the exam.

Stenosis Detection & ColorID

- With the customizable one- or two-point deposit technique, you can automatically quantify degree of stenosis and segment length.
- By assigning customizable color mapping to CT number ranges within an image, you can identify plaque on an axial image, curved reformat, or MPVR image. Up to ten levels of color coding are available. Colors transition smoothly from one density to another.
- Make volume and area measurements to track plaque size over time.
- Custom sizeable 3D ROI lets you apply color over all anatomy.
- Generate multiple ROIs on one vessel or multiple vessels.
- Automatically include each defined plaque, its volume and area into the report.

Quick AVA

- Dynamic AVA, provides access to vessel analysis from any protocol.
- Launch the Quick AVA feature from the main tool bar.
- Perform vessel analysis from any 3D or reformatted image with one- or two-point clicks.
- If you detect a lesion, you can analyze it by dropping a point above and below the section. The vessel analysis tool launches automatically for fast review in curved reformat, cross-sectional, lumen, and MRP view.
- Extend tracking proximally or distally for extended views at any time.


CONFIDENTIAL

SECRETE DE
AFACERI


Features Detail *(cont.)*

Summary Table & Report

- Measurements are automatically captured in summary table for easy review.

- Simply click on any measurement in the summary table and you will be taken to that measurement location.
- Prepare report with a single click or generate an advanced report with customized options
- Send the report to database, filmer or PDF.

Multiple Study Comparison

With VessellQ Xpress, compare a previous study with the current study. You also have the flexibility to launch comparison within Quick AVA.

Save State & Protocol

- Save the current state of tracking and measurements for easy reloading with a single mouse click. Save state captures:
 - 3D images
 - Vessel tracking
 - Measurements
 - Vessel contours.
- Transfer the Save State to the PACS or other AW workstations for future review.
- Save custom protocols including anatomical labels and measurements.

Summary

VessellQ Xpress provides user- friendly tools an intuitive interface to help streamline your vessel analysis workflow and help you generate more comprehensive reports for referring physicians.

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GE HealthCare

CardIQ Xpress 2.0

Bringing integration and automation to your CT Cardiac post processing needs.

CONFIDENTIAL

SECRET DE
AFACERI

Contents

1. Features.....	4
2. System Requirements.....	4
3. Regulatory Compliance.....	4
4. Indications for Use.....	4
5. Features Detail.....	4
5.01 Pre-processing.....	4
5.02 Auto Launch.....	4
5.03 Automated Coronary Vessel Analysis (CVA).....	5
5.04 2D Review.....	5
5.05 Quick/Dynamic AVA.....	6
5.06 Vessel Data Base.....	6
5.07 Relative Perfusion.....	6
5.08 IVUS Views.....	6
5.09 PlaqID.....	6
5.10 VR Heart.....	6
5.11 Multi-Phase Image Review.....	6
5.12 Transparency View.....	6
5.13 3D Angiographic View.....	7
5.14 One Touch Cath Views.....	7
5.15 3D Ejection Fraction.....	7
5.16 Layouts.....	7
5.17 4D Movie Mode.....	7
5.18 SnapShot Freeze (SSF).....	7
5.19 Protocols.....	7

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SECRET DE
AFACERI

Cardiac disease is one of the leading health concerns worldwide. Successful treatment of the many conditions that cause and perpetuate heart disease requires that physician's approach specific cardiac problems with as much information as today's technology can provide. To supply your referring physicians with this information you need software that gives insight into cardiac anatomy and tools for thorough evaluation into the extent of disease.

Overview

CardIQ Xpress 2.0 is an integrated post processing imaging analysis application dedicated to cardiovascular imaging on GE HealthCare Advantage Workstation (AW) and AW Server Enterprise System (AWS). The CardIQ Xpress 2.0 software option can be used to display, reformat, and analyze 2D or 3D cardiac CT images for qualitative or quantitative assessment of heart anatomy and coronary artery vessels from a single or multiple cardiac phase image data sets. Cardiac motion is a very real challenge that can occur at any heart rate. CardIQ Xpress 2.0 is designed to work with SnapShot Freeze* images to automatically process and display images generated with reduced motion blur artifact.

Highlights

- Automatically segment coronary tree across phases.
- Automatically tracks and labels coronary arteries.
- Improved centerline editing tools for faster edits.
- Right mouse wheel menu for quick access to renaming, deleting, and editing centerlines.
- **Plaque ID tool assists in visualizing and quantifying plaque burden**
- Relative perfusion highlights and quantifies hypo-dense areas of myocardium

CONFIDENTIAL

SECRET DE
AFACERI

Visit us:

<https://www.gehealthcare.com/products/advanced-visualization/all-applications/cardiq-xpress-20-reveal>

Features

- Pre-processing automatically recognizes cardiac datasets and performs all segmentations as data arrives on the system saving processing time
- Auto Coronary analysis automatically segments the coronary tree, tracks and labels the coronary arteries
- Three pre-defined orientations provide quick access to planes for best visualization of the coronaries in 2D
- Quick AVA allows access to vessel tracking at any time
- Relative perfusion color codes and quantifies percent of hypo-dense areas of myocardium with four selectable color maps and fusion overlays with the coronary tree
- **Plaque ID provides volume measurements for four distinct Hounsfield ranges to aid with identification and visualization of coronary plaque**
- Automatically display SnapShot Freeze¹ processed images for reduced motion blur
- Robust, automatic calculation of ejection fraction and stroke volumes from the 3D endocardium volumes
- **Stenosis measurements, IVUS views and preset volume rendering models assist in communication of specific findings**
- Measure ES and ED for ejection fraction & volume with automatic extraction of the left ventricle
- Create multiphase beating hearts
- Select oblique reformat views in the standard cath angles for easy analysis of the coronary vessels
- Display 4D valve views with a single click

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Rx Only

Indications for Use

CardIQ Xpress 2.0 is intended to provide an optimized non-invasive application to analyze cardiovascular anatomy and pathology and aid in determining treatment paths from a set of Computed Tomography (CT) Angiographic images. CardIQ Xpress 2.0 is a CT, non-invasive, image analysis software package, which aids in diagnosing of cardiovascular disease to include, coronary artery disease, functional parameters of the heart, heart structures and follow-up for stent placement, bypasses and plaque imaging. CardIQ Xpress 2.0 offers unique tools such as automatic tracking, which will pre-process the CT data into multiple viewing ports to allow for an expedited read time improving workflow. With CardIQ Xpress 2.0, the user can color code the myocardial tissue to show hypo/hyper-dense areas in the myocardial tissue of the heart. With the IVUS-like view the user can color code the HU units of the plaque to better visualize the difference between calcified and non-calcified plaque in the wall of the vessel and the lumen to determine the amount of atherosclerosis. The user can see the different valve planes along with a variety of new layouts to align the heart. The IVUS-like view is created by applying GE HealthCare's Volume Rendering on a cross-section perpendicular to the detected centerline. This view merely displays a cross section as in IVUS imaging and color codes like IVUS images. No new or additional diagnostic information is added. CardIQ Xpress 2.0 is for use on the Advantage Workstation (AW) platform, CT Scanner, PAC or Centricity stations, which can be used in the analysis of 2D or 3D CT angiography images/data derived from DICOM 3.0 CT scans.

Features Detail

Pre-processing

- Automated processing of single or multiphase data
- Customizable menu to choose which protocols are pre-processed
- Pre-processed protocols include: Auto Coronary Analysis, Labeling of vessels, VR heart, Angiographic View, and Auto EF

Auto Launch

- Multiple exams can be loaded into the auto launch to ready them for review (AW only)
- Color-code display of ready-to-read exam listing (AW only)
- The ability to switch between exams without having to quit out of the application and reload


CONFIDENTIAL

SECRET DE
AFACERI

Features Detail *(cont.)*

Automated Coronary Vessel Analysis (CVA)


- CVA allows the users to track, extract, visualize and measure coronary arteries from either single or multiple cardiac phase images
- Auto launch case selection and display of automatically tracked and labeled coronary vessels
- Simple, corrections to vessel branch tracking
- Curved, oblique, longitudinal, and cross-sectional reformatted views are automatically generated in the software for quick review

- Current-state tracking points within vessel analysis may be saved for future review and/or manipulation
- Color-coded plaque analysis with volume measurements for evaluating the change in plaque size over time
- Generate and save rotational movies from curved planar reformation, best section and cross-sectional (lumen) views
- Quantitative or qualitative coronary vessel assessment on user-selected vessel segments.
- Vessel measurements including:
 - Distance and volume
 - Cross-sectional area and mean diameter
 - Single or Dual reference point comparison
 - Relative percent stenosis
- Single-screen filming capabilities with multi-views within the screen to show the entire picture of the vessel

Automated Coronary Vessel Analysis (CVA) *(cont.)*

- Automatically display of the coronary vessel tree
- Using pre-set protocols, 2D or 3D coronary vessel tree models can be generated and displayed in an automatic fashion to give users a qualitative overview of coronary vessel structures

- New right mouse wheel menu for quick access to renaming, deleting and editing centerlines
- Improved centerline editing tools for faster edits
- The vessel tree models can be in the form of either 3D volume rendered or Maximum Intensity Projection (MIP)

2D Review

- Three pre-defined orientation protocols for easy review of the coronary vessels and chambers
- Dual reformat review allows automatic review of the coronaries from an axial image while linking to oblique views in longitudinal and cross-sectional planes
- Reformat review predefines workflow steps to automatically present thick, multi-planer views of cardiac anatomy
- One-click access to quick vessel analysis

CONFIDENTIAL

SECRET DE
AFACERI


Features Detail *(cont.)*

Quick/Dynamic AVA

- Single or dual click vessel analysis from any protocol; 3D, reformat or vessel tree
- Real-time tracking of the center line with instant visualization of a vessel for quick inspection
- Ability to extend tracking proximally or distally for full view of the vessel
- Ability to select a point above and below potential lesions for automated vessel analysis

Vessel Data Base

- A vessel name data base for quick access to vessel labels
- Anatomically relevant listings
- Ability to add new vessel names into the data base
- Automated processing of single or multiphase data

Relative Perfusion

- Hypo-dense areas of myocardium can be high-lighted with color to aid in the visualization of ischemic heart disease
- Pre-defined layouts for density enhancement
- A hybrid display view to show the relationship of the vessel tree with the perfusion defect
- Quantification of a hypo-perfused area as it relates to myocardial defects
- Four selectable color maps to help display the hypo-dense areas
- Exportable statistics via one touch

IVUS Views

- Interactive volume rendered images to better define the calcium, lumen and non-calcified plaque related to lumen narrowing or wall abnormalities
- Can be applied to any 2D MIP image to include; cross-section, MPR, best L-section images
- Display cross-sections perpendicular to the vessel centerline to create an IVUS-like image. No new diagnostic information is added

PlaqID

- Customizable color mapping to HU ranges for easy Identification of plaque in axial, curved reformat or MPVR images
- Four distinct colors to aid in differentiating vessel lumen, non-calcified plaque and calcified plaque
- Easy modification of the color and transparency of plaque
- Volume, area to track size of plaque over time
- Smooth transition of color from one density to another
- Plaque volume and area automatically captured in summary table.

VR Heart

- One-touch automatic segmentation algorithm designed to extract the cardiac anatomy from within the chest image
- Optimize algorithms for the easy review of heart and bypass grafts
- Multiple optimized VR curve settings to enhance different structures within the heart.
- Automatic segmentation of both single or multi-phase cardiac image data sets
- 4D beating heart images with the ability to page through the phases, rotate the image with real-time functionality

Multi-Phase Image Review

- Multiphase images can be reviewed in any protocol
- 10 phase images show the heart throughout the complete cycle
- Quickly edit phase images to keep only the phase or phases needed for analysis of the coronary vessels once multi-phase review is complete

Transparency View


3D vessel tree overlaid on the heart chambers with the ability to adjust opacity of chambers of heart, myocardium vessels tree and bones

CONFIDENTIAL

SECRET DE
AFACERI

Features Detail *(cont.)*

3D Angiographic View

- One-touch automatic visualization of the coronary arteries in an x-ray angiographic view
- Single click cardiac cath views to orientate the angiographic view into RAO/LAO, or CRANIAL/CAUDAL orientation
- Added or removed vessels are easily archived with the Auto Select feature
- Toggle between gray scale or inverse video views

One Touch Cath Views

- Default cardiac cath views are provided with CardIQ Xpress 2.0 software to provide 3D VR heart, 3D vessel tree or reformatted image in the orientation of a standard cardiac catheterization
- Single click on the angulation area to change between one-touch views
- Interactive RAO/LAO & CRANIAL/CAUDAL orientations allow for manually positioning the view
- Customizable orientations can be easily generated and saved for future use

3D Ejection Fraction

- Automatic extraction of the left ventricle across all phases and the automatic detection of systole and diastole delivers calculated ejection fraction and stroke volumes
- Robust, automatic calculation of ejection fraction and stroke volumes from the 3D endocardium volumes

Layouts

- Configurable layouts
- Dual monitor support
- Configurable default layout for use at AVA launch
- Easily capture the screen layout for identical visualization on PACs or film

4D Movie Mode

- Movie mode allows the user to view and save multiple views of the heart and display as a movie sequence
- Load multiple phases from systole to diastole to create a beating heart movie
- Save movies as a DICOM image set or export in JPEG/MPEG format

SnapShot Freeze (SSF)

- Automatically display motion blur-reduced SnapShot Freeze processed images
- Reprocess SnapShot Freeze images for artifacts caused by motion blurring after manual vessel editing
- Motion blurring correction requires the use images produced by a CT using the optional SnapShot Freeze feature

Protocols

- CardIQ Xpress 2.0 is supplied with a set of pre-defined protocols that are easily adaptable to customized protocols
- The protocols have comprehensive instructions and tools which minimize the need to memorize the procedures or to refer continuously to user documentation.

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Reference:

* Trademark of General Electric Company.

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Referință:

* Marcă comercială a General Electric Company.

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SECRET DE
AFACERI


CardIQ Function Xpress

Real time cardiac review the instant you're ready to read.

CONFIDENTIAL

SECRET DE
AFACERI

Contents

1. Features..... 4

2. System Requirements..... 4

3. Image Requirements..... 4

4. Indications for Use..... 4

5. Regulatory Compliance..... 4

6. Features Detail..... 4

6.01 Features..... 4

6.02 Pre-processing..... 5

6.03 Auto Launch (AW only)..... 5

6.04 Automatic Chamber Detection..... 5

6.05 Automatic Ejection Fraction..... 5

6.06 Display Table..... 5

6.07 Myocardial Analysis..... 6

6.08 Simplified Workflow..... 6

CONFIDENTIAL

SECRET DE
AFACERI

Cardiac disease is one of the leading health concerns worldwide. Successful treatment of the many conditions that cause and perpetuate heart disease requires that physicians approach specific cardiac problems with as much information as today's technology can provide. To supply your referring physicians with this information you need software that gives you insight into cardiac anatomy and function supported by quantifiable data.


Overview

CardIQ Function Xpress post-processing software helps you evaluate cardiac function and aids in the diagnosis of cardiovascular disease with a high degree of confidence. Providing accurate and reproducible quantification of left and right ventricular volumes, ejection fractions, and myocardial mass, CardIQ Function Xpress is optimized to assess cardiac function using multiphase, multi-slice cardiac CT images.

The CardIQ Function Xpress option provides an easy-to-use and time-effective way for you to perform cardiovascular functional analysis.

Highlights

- Automatically pre-processes and loads exams (AW only).
- Automatically detects all chambers in all phases.
- Automatically performs ejection fraction and chamber volume analysis with 91% reliability on LV and RV automatic ejection fraction.
- Provides single click myocardial analysis with bull's eye images.
- Performs simultaneous case review in real time.
- Load multiphase data with no load time.

CONFIDENTIAL

SECRET DE
AFACERI

Visit us:

<https://www.gehealthcare.com/products/advanced-visualization/all-applications/cardiq-function-xpress>

Features

- Performs behind-the-scene function data process for real time ejection fraction review, volume analysis, and myocardial analysis.
- Enables you to extract, render, and display 3D volumetric models of the endocardium for ejection fraction calculation.
- Automatic detection of the endocardial and epicardial walls for wall motion, wall thickness, wall thickening and myocardial mass analysis.
- Automatic volume analysis of all chambers.
- Automatically calculates left atrium volume while excluding the pulmonary vein.
- Activate visual wall motion with short axis images in basal, mid, and distal orientations along with a two-chamber long axis view in a single click.
- Perform myocardial analysis with wall motion, wall thickness, and mass calculations.

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Image Requirements

CardIQ Function Xpress option can be applied to gated CT images.

These images can be acquired on multi-slice CT scanners. The GE HealthCare family product line must have the CardIQ SnapShot option to acquire gated images.

Indications for Use

CardIQ Function Xpress is intended to provide an optimized noninvasive application to analyze cardiovascular anatomy and pathology and aid in determining treatment paths from a set of Computed Tomography (CT) Angiographic images. CardIQ Function Xpress in conjunction with CT cardiac images to automatically calculate and display various left ventricular and right ventricular functional parameters as ejection fraction, end systolic and end diastolic volumes, stroke volumes, wall motion, wall thickening, cardiac output, myocardial mass, systemic and pulmonary vascular resistance. Volume measurement of each chamber of the heart is also available. With CardIQ Function Xpress atrium volumes may be used to determine volume assessment of atrial disease to include but not limited to atrial fibrillation. CardIQ Function Xpress is a CT, non-invasive image analysis software package, which aids in the assessment of cardiac function and in determination of cardiovascular disease diagnosis and management. CardIQ Function Xpress is for use on the Advantage Workstation (AW) platform, CT Scanner, PAC or Centricity stations, which can be used in the analysis of 2D or 3D CT angiography images/data derived from DICOM 3.0 CT scans.

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
Features Detail

Features

CardIQ Function xpress allows you to accurately and reproducibly quantify myocardial mass, left and right ventricular volumes and ejection fractions. The software automatically detects left ventricular endocardial and epicardial contours so that you can more easily assess left ventricular (LV) and right ventricular (RV) functional parameters. The program is optimized to perform assessment of cardiac function using GE HealthCare multiphase, multi-slice cardiac CT images.

LV Segmentation

RV Segmentation

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**SECRET DE
AFACERI**

Features Detail (cont.)

Pre-processing

- The program automatically processes multi-phase data for functional analysis.
- You can select configurations to process only the chambers you require.
- For more streamlined reading, ejection fraction, volume analysis, and myocardial analysis are completed before you view the case.

Auto Launch (AW only)

- Load multiple cases into auto launch so they are ready to review when you are.
- Switch between exams without having to quit the application and reload.
- Exams ready for review display in green.

Automatic Chamber Detection

- The program loads single or multi-phase data sets and selects the heart's chambers for volume analysis.
- Analyze individual chambers with automatic contouring of the heart's four chambers.
- The system automatically detects the left/right ventricles and atriums and contours the endocardium for ejection fraction calculation.
- Works with single or multi-phase datasets


Automatic Ejection Fraction

- The system features a protocol that, when chosen, automatically detects the end diastolic and end systolic phases and calculates an accurate ejection fraction of the left ventricle, right ventricle, and left atrium.
- Ejection fraction measurements are available for the right and left ventricles.
- Stroke volume, end diastolic and end systolic volumes are automatically calculated and shown in a graphical representation for assessment of heart failure.
- With dual screens, you can display quantification assessments on the left monitor and a qualitative overview of the full chamber from apex to base on the right monitor.
- Customize layout designs to accommodate your viewing preferences:
 - Base to apex
 - Basal, mid, and distal along with a two-chamber long axis.
- You can change the selection of end diastolic and end systolic phase.
- While contour detection is automatic, scalpel and auto select tools are available with which you can edit the contours.

Display Table

A display table shows volume and provides a graphical update of end diastolic and end systolic

- Click on available tabs to select different chambers for analysis.
- Left and right ventricle measurements include:
 - Stroke volume
 - Ejection fraction
 - Cardiac output
 - Myocardial mass
 - End systolic volume
 - End diastolic volume
 - Pulmonary vascular resistance
 - Systemic vascular resistance
 - Body surface area
 - Heart rate
- Left atrium measurements include:
 - Min and max volume
 - Ejection fraction
 - Left atrium index
 - Body surface
 - Heart rate
- Right atrium measurements include:
 - Min and max volume
 - Heart rate
- You can film, save, and export the volume graphs.


CONFIDENTIAL

SECRET DE AFACERI

Features Detail *(cont.)*

Myocardial Analysis

- CardIQ Function Xpress automatically detects the endocardial and epicardial walls to analyze wall motion, wall thickness, wall thickening, and myocardial mass analysis.
- You can edit the contours with a freehand trace and click to spline tools.
- The program displays multiple bull's eye graphs to graphically represent these cardiac characteristics:
 - Wall motion
 - Wall thickness
 - Myocardial mass
 - Cardiac output
- Choose from four bull's eye graph overlays to show 17, 20, and 49 segment model, or a coronary artery territory model.
- Manual scoring is available for wall motion and thickening in 17 and 20 segments.
- Cord display is available for visual assessment of wall abnormalities.

Simplified Workflow

CardIQ Function Xpress provides you with the tools you need to accurately analyze, evaluate, and report on cardiac anatomy and function. Its ease of use will streamline your workflow and its accuracy will give you greater diagnostic confidence.

CardIQ Function Xpress simplifies your workflow with:

- Automatic ejection fraction protocol
- Streamlined multiphase beating images for easy visualization of cardiac wall motion
- Single click myocardial analysis with bull's eye images

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GE HealthCare


CardEP

Streamlined post-processing for enhanced electrophysiology procedures.

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SECRET DE
AFACERI

Contents

1. Features..... 4

2. Image Requirements..... 4

3. System Requirements..... 4

4. Indications for Use..... 4

5. Regulatory Compliance..... 4

6. Features Detail..... 5

6.01 Operation..... 5

6.02 3D VR Atrium..... 5

6.03 3D VR Heart..... 5

6.04 Image Reformation..... 5

6.05 Coronary Sinus and Pulmonary Vein Analysis..... 6

6.06 Filming..... 6

6.07 One Touch EP Views..... 6

6.08 Pulmonary Vein (PV) and Left Atrial Appendage Reformation..... 6

6.09 Phase Registration..... 6

6.10 EP Navigator Views..... 6

6.11 Batch Movie Mode..... 6

6.12 Protocols..... 7

6.13 Summary..... 7

CONFIDENTIAL

SECRET DE
AFACERI


CT cardiac exams provide critical information to practitioners of cardiac therapy. Precise information pertaining to the left atrium's complex anatomy, the pulmonary veins, the coronary sinus, or the cardiac veins has a major impact on the efficacy of subsequent cardiac therapy: It can speed procedures and facilitate treatment. But analyzing and reporting the wealth of information CT cardiac studies provides can be time consuming. A program that automates many routine functions and gives you tools to easily quantify and qualify many aspects of cardiac function would streamline your workflow and give you greater diagnostic confidence.

Overview

CardEP is an integrated post processing image analysis software for the application of cardiovascular and electrophysiology imaging on the AW Workstation. With it, you can process display, reformat, and analyze 2D or 3D cardiac images for **qualitative or quantitative assessment of heart anatomy** and pulmonary veins from single or multi-phase cardiac image data sets.

Highlights

- Visualize the anatomical layout of the left atrium and pulmonary vessels.
- Visualize the origins of the pulmonary vessels.
- Automatic heart segmentation.
- Helps you visualize and quantify cardiac venous pathways.
- Offers you a variety of 2D, 3D, or reformatted protocols with which to perform image analysis.

CONFIDENTIAL

SECRET DE
AFACERI

Visit us:

<https://www.gehealthcare.com/products/advanced-visualization/all-applications/cardep>

Features

CardEP offers a variety of features that give you image analysis flexibility:

- Volume render 3D heart and chamber views
- Obtain segmented 3D views of the left atrium and pulmonary veins
- View hollow 3D VR heart models
- Measure the pulmonary veins and left atrial appendage with predefined protocols
- Render and display 2D/3D views of the left chambers in single or multiple phases of the cardiac cycle with one-touch automatic extraction
- Reformat standard axial CT images of single or multiple cardiac phases into short or long axis and save the image series for further analysis.
- Automatically track, extract, and display pulmonary veins with the pulmonary vessel analysis tool
- Generate pulmonary vein measurements
- Automatically track, extract, and display coronary sinus and cardiac veins with the coronary sinus analysis tool
- Perform one-touch standard 2D reformation and 3D volume rendering, including hollow view, of cardiac anatomy with optimized AW Volume Viewer and 3D rendering tools
- Register images from different cardiac phases into a unique data set that you can save as a 3D object or use for further analysis
- Use predefined navigator views of the pulmonary veins to get a fly-through perspective of the vessels.

Image Requirements

CardEP accepts standard cardiac gated CT image sets acquired on CT scanners equipped with DICOM 3.0 standards. Images must meet the same image requirements as those for the basic Volume Viewer application.

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Indications for Use

CardEP is a post-processing software option for the Advantage Workstation (AW) Platform. This product can be used for the analysis of CT angiographic images for the assessment of the heart to include the atria, pulmonary veins and coronary sinus. It provides quantitative analysis tools which include a number of display, measurement and model export capabilities. This product can be used to aid trained physicians in the visualization and assessment of cardiac anatomy.

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
SECRET DE
AFACERI

Features Detail *(cont.)*

Coronary Sinus and Pulmonary Vein Analysis

- Load single or multi-phase cardiac data sets in the vessel analysis tool.

- With the semi-automated pulmonary vein/coronary sinus analysis tool (CSA), based on the Advanced Vessel Analysis (AVA) technique, you can track, extract, visualize, and measure the veins and sinus (physical dimension and stenosis sizing) from either single or multiple cardiac phase image data sets.
- Generate curved, oblique, longitudinal, and cross-sectional reformatted views of individual vessels.
- Generate 3D models of the atrium and pulmonary veins during your vessel analysis procedure.
- Save current state tracking points within your vessel analysis for future image review or manipulation.
- Generate and save rotational movies from curved planar reformation and cross-sectional (lumen) views.
- Generate and save movies from Best L Section reformatted images.
- Perform quantitative and qualitative coronary sinus assessments on segment you choose:
 - Distance measurement
 - Cross-sectional area
 - Single or dual reference point comparison
 - Relative percentage stenosis
 - Volume measurement
 - Mean diameter
- Generate beating, rotating movies for visualization and communicating findings.

Filming

CardEP offers you filming flexibility. Batch filming capability is provided, and any image can be queued for filming.

One Touch EP Views

- Default views provide you with visualization of the atria and pulmonary veins in correct anatomical orientation. You are provided the capability to create your own presets to save and apply to future analyses.
- Use the interactive RAO/LAO and CRA/CAU annotation feature to manually orient the view.

Pulmonary Vein (PV) and Left Atrial Appendage Reformation

- Predefined protocols allow you to generate reformation of the pulmonary veins and left atrial appendage in the long axis.
- Generate cross-sectional views of the veins or appendage to determine the ostium's diameter.

Phase Registration

- Suppress artifacts that may occur when the heart rate fluctuates during a scan with a semi-automated multi-phase tool that aligns cardiac chambers, veins, or sinus structures.
- You can automatically or manually register phases among images to improve image quality.
- Save the registered volume data sets for future review or manipulation.

EP Navigator Views

- With one-touch, the navigator displays endoscopic views of the left atrium.
- View the pulmonary vein and atrial appendage ostia from either the left or right side using a predefined orientation layout preset.
- Optimize interior atrium views with automated threshold settings.
- Global virtual endoscopic view provides you ostial displays of the right and left pulmonary veins and the left atrial appendage.

Batch Movie Mode

- Batch movie mode lets you define first and last images for movie creation.
- Create rotating and beating movies
- Save automatically generated movies to the viewer as DICOM image sets or send them to data export in JPEG/MPEG/AVI formats.
- Use the movie mode with MIP, 3D VR, and reformatted images.
- Generate and save movies from Best L Section Reformatted images, curved planar reformatted, cross-sectional, and lumen views.

CONFIDENTIAL

SECRET DE
AFACERI

Features Detail *(cont.)*

Protocols

CardEP contains predefined protocols that you can edit or modify to adapt them to your specific customized protocols and save under a new name.

The CardEP protocols comprise comprehensive instructions and the tools you need to perform analyses. Tool icons are displayed in the protocol panels for fast and easy access.

Summary

CardEP gives you the tools and flexibility to quickly analyze cardiac studies. The information and statistics you can generate with this program may help referring physicians with their assessment of therapy options for their patients.

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
CONFIDENTIAL

SECRET DE
AFACERI

Advanced imaging software that detects, quantifies, and scores cardiac calcium plaque burden.

Contents

- 1. Features..... 4
- 2. System Requirements..... 4
- 3. Compatible Platforms..... 4
- 4. Indications for Use..... 4
- 5. Regulatory Compliance..... 4
- 6. Features Detail..... 4
 - 6.01 Scoring Algorithms..... 4
 - 6.02 Automatic Calcium Detection..... 4
 - 6.03 Customizable Configuration..... 4
 - 6.04 User Interface..... 4
 - 6.05 Reporting 5
 - 6.06 Summary..... 5


Cardiovascular disease remains one of the most common health issues in the world today. As with many conditions, early detection and patient risk assessment are vital to preventing or minimizing long-term negative effects. But many cases of this disease aren't diagnosed until the patient presents with symptoms. Conventional procedures to assess risk can be expensive, time consuming, and uncomfortable. An accessible, patient-friendly risk assessment method could help clinicians devise a regimen for their patients that might lessen the chance for serious cardiac events.

Overview

SmartScore 4.0 is designed to identify the presence of regional and global coronary artery calcification from a CT scan, then measure and score the results. Scores can be calculated using a standard Agatston/Janowitz (AJ) method. When correlated with a patient's personal information, the score can yield an estimation of a patient's risk for coronary artery disease.

Highlights

- Non-invasive alternative to conventional assessment procedures.
- Score can be correlated to age group cohort to determine patient's risk per population
- Provides information on coronary artery wall calcium plaque buildup.
- Automatically detects calcium and highlights it in green.
- Streamlines workflow by networking patient demographics from CT scanner directly to SmartScore program.
- Free Hand Trace lets you outline specific ROIs.

CONFIDENTIAL

SECRET DE
AFACERI

Visit us:

<https://www.gehealthcare.com/products/advanced-visualization/all-applications/smartscore-40>

Features

- AJ 130 scoring method uses conventional Agatston/Janowitz technique with a threshold of 130 HU which is adjusted to the appropriate image slice thickness.
- Volume scoring calculates volumes (mm³) of calcified plaque above the 130 HU threshold.
- Mass scoring calculates mass (mg) of calcified plaque above the 130 HU threshold.
- Individual and aggregate scores are computed for each artery type.
- Provides two methods of calcium scoring.
- Customizable parameters can be predefined based on your site's preferences.
- User interface streamlines your workflow and reduces reading time.
- Report options allow you to custom tailor your reports and distribute them in a variety of formats.

System Requirements

- EKG monitor with recording device and x-ray translucent lead cable.
- Gantry hardware upgrade kit for those scanner systems already in operation.

Compatible Platforms

- AW4.7 Ext.18 (and later)
- AW Server 3.2 Ext.4.0 (and later)
- Recommended monitor resolution is up to dual 2MP (1600 x 1200) or a single 3MP (1536 x 2048)

Indications for Use

SmartScore is a non-invasive software option that can be used to evaluate calcified plaques in the coronary arteries, which may be a risk factor for coronary artery disease. SmartScore may be used to monitor the progression/regression of calcium in coronary arteries over time, which may aid in the prognosis of cardiac disease.

Regulatory Compliance

This product complies with the European Council Directive 93/42/EEC Medical Device Directive as amended by European Council Directive 2007/47/EC.

This product or its features may not be available in some other countries or regions. Please contact your sales associate.

Rx Only

CONFIDENTIAL

SECRET DE
AFACERI


Features Detail

Scoring Algorithms

- Scoring uses the conventional AJ 130 method and provides volume (mm³) and mass (mg) scores. Scores are computed and reported for each artery type. Three patient types are available for mass score: small; medium; large.
- Two methods of scoring are provided:
 - Region growing (3D scoring) lets you score calcium plaque on multiple slices by clicking on an image, then on any calcium attached in the relevant lesion. All calcium attached to that calcium plaque is automatically scored, even if it extends into other slices.
 - Free Hand Trace (2D scoring) lets you draw a free hand ROI around a calcified area. This allows you to score image by image, separating lesions that may be present in two arteries for more accurate scoring. For heavily calcified areas, this method reduces scoring time.

Automatic Calcium Detection

SmartScore 4.0 automatically highlights calcium within coronary arteries to facilitate scoring. Calcified plaque can be labeled LMA, LCx, RCA or A, B, C, which can be customized to score areas such as the mitral valve, aorta, or carotids.

Customizable Configuration

You can customize the SmartScore 4.0 configuration to parameters you prefer. The program allows you to:

- Turn automatic calcium detection on or off
- Turn image annotation on or off
- Set proper image window/level
- Choose a preferred image display format
- Select output medium.


User Interface

The program's tool panel provides a variety of tools to help streamline your workflow. Choose to activate these tools using the right, left, or middle mouse button. Tools include:

- Zoom
- Pan/Roam
- Window/Level
- Free Hand Trace
- Point and Click scoring.

Features Detail (cont.)

Reporting


The SmartScore 4.0 patient report is subdivided into selectable sections:

- Header/Physician list allows you to place in the header your hospital or clinic logo, facility address, and physician names. All can be left or right justified or centered.
- Patient Demographics area includes these predefined parameters (twenty additional fields can be added to this area to facilitate study support):
 - Patient ID
 - Name
 - Exam date
 - Birth date
 - Age
 - Sex
 - Ethnicity
 - Weight
 - Height
 - Cholesterol levels
 - Blood pressure
 - Diabetes status
 - Smoking status
 - Medications
 - Cardiac history
 - Family history
 - Study notes
 - Scored by
- Generic Message provides configurable space in which to describe the procedure and test results.
- Diagnostic Message provides space for you to describe your interpretations of structures surrounding the heart.
- Score Interpretation provides space to summarize the meaning of different score values.
- Score Summary provides space in which to report the patient's score and rank (if selected).
- Category Table provides space to enter clinical information based on the calcium score.
- Score Table lets you identify the scoring algorithm used.

Reporting (cont.)

- Vessel Distribution shows the distribution of calcifications/lesions for each vessel scored.
- Population Distribution Graph shows the patient's score in relation to their age-specific cohort.
- Guide Table lets you rank a patient's percentage of risk based on age, gender and score. The table is editable, allowing you to change age group and percentile ranking.
- A diagram of the heart can be included in the report as can select images from the exam.
- Physician signatures can be scanned and entered into the SmartScore 4.0 software

Summary

SmartScore 4.0 lets you offer patients a convenient method for evaluating their risk of coronary disease.

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CONFIDENTIAL

SECRET DE
AFACERI

Automated segmentation & pre-procedure planning software for transcatheter aortic valve implantation.

Contents

- 1. Features..... 4
- 2. System Requirements..... 4
- 3. Indications for Use..... 4
- 4. Regulatory Compliance..... 4
- 5. Features Detail..... 5
 - 5.01 Valve Definition..... 5
 - 5.02 Aorta Analysis..... 6
 - 5.03 Procedure Preparation..... 7
 - 5.04 Summary..... 7

Transcatheter Aortic Valve Implantation (TAVI) or Replacement (TAVR) is a procedure for select patients with severe symptomatic aortic stenosis (narrowing of the aortic valve opening). TAVI/TAVR is performed on a beating heart and a replacement valve is delivered via a catheter.

CT provides important information that is important for successful TAVI/TAVR procedures. CT is used to help determine aortic annulus size, to guide selection of appropriate valve, provide dimensions of the entire aorta, help determine the access path for the catheter and give guidance for C-arm angulation for deployment of the device.


TAVI Analysis provides a streamlined, guided workflow to enable efficient, consistent work-ups of your TAVI studies and connectivity options directly to the interventional suite.

Overview

TAVI Analysis offers a guided workflow to add efficiency and standardization to your departments work up and communication for TAVI/TAVR procedure planning or follow-up. It automatically segments the aorta and displays the aortic valve in multiple views for quick and easy measurements of the annulus. TAVI Analysis provides automatic and semi-automated tools to help plan aortic annulus device sizing, evaluate appropriate access pathways and provide guidance to optimal C-Arm positioning angles during the procedure.

Highlights

- Automated Aortic segmentation and preliminary orientation to aortic annular plane
- 3 pt deposit mode for precise refinement of annular plane
- Faster valve plane detection over manual method²
- Automatic aortic valve calcium volume estimation measurement¹
- Automatic centerline tracking for peripheral access planning
- Summary table and automated report for consistent and efficient communication to clinical care team

CONFIDENTIAL

SECRET DE
AFACERI

Visit us:

<https://www.gehealthcare.com/products/advanced-visualization/all-applications/tavi-analysis/>

Features

- One comprehensive workflow for valve measurements and access path measurements.
- Automatic segmentation of transparent volume rendered model of the aorta with calcifications highlighted
- Auto orientation to aortic valve plane with 3 pt. deposit mode for users to precisely refine for each patient's anatomy providing faster workflow over a manual approach
- Load multiple phases for quick review of image quality and comparison of annulus size differences
- Automated pre-contouring of annulus with easy editing
- Automated distance measurements for coronary ostia heights
- Automatic volume measurement of estimated aortic valve calcifications¹
- Automatic aorta angulation measurement¹
- Summary table automatically populates as measurements are placed.
- Easily generate a report with measurements and associated images for efficient communication with Heart Team.
- One click perpendicular view and "S-curve" display with full range of C-arm angulation guidance¹
- Transparent bone with vascular display for quickly identifying external anatomical landmark locations
- Auto tracking of peripheral vascular branches to assist with access planning procedures
- Peripheral artery measurements automatically sent to summary table to allow easy exporting and report generation.
- Save State and send to the Cath lab for sharing access to models and measurements.

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Indications for Use

VesselIQ Xpress is an optional, non-invasive, optimized, post-processing application intended to provide images and tools to analyze vascular anatomy and pathology, aiding physicians in diagnosis and determination of treatment paths, from a set of Computed Tomography (CT) Angiographic images.

VesselIQ Xpress is an option for the Advantage Workstation (AW) platform, CT Scanner, and/or PACS, which can be used in the analysis of 2D and 3D CT Angiography images/data for the purpose of cardiovascular and vascular disease assessment. This software-only device is designed to support physician assessment for a wide variety of clinical uses such as stenosis analysis, pre/post stent planning, pre/post valve replacement planning, and directional vessel tortuosity visualization.

VesselIQ Xpress' automatic visualization tools provide users the capability to segment bony structures for accurate identification of the vessels. Additional tools enable analysis of the vascular anatomy including the aorta, valves, and branching vessels for: anatomical sizing; density and volume analysis of segmented vasculature and calcified / non-calcified plaque; and measurements of abnormalities

The TAVI Analysis option for VesselIQ Xpress is a planning tool used for Trans Aortic Valve Implantation (TAVI) procedures. It automatically segments the aorta and displays the aortic valve in multiple views for measurements of anatomic structures commonly needed for aortic annulus replacement planning. TAVI Analysis provides guided workflow and semi-automated tools to aid in evaluation of appropriate access pathways for interventional procedure planning.

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
SECRET DE
AFACERI

Features Detail


TAVI Analysis offers a guided workflow to add efficiency and standardization to your departments work up and communication for TAVI/TAVR procedure planning or follow up.

Valve Definition

- Automatic segmentation of aorta with auto orientation to valve plane and 3 pt. deposit mode for precise refinement of annular plane.

- Automated annular pre-contouring with easy editing to quickly obtain annular measurements
- Load multi-phase data and toggle between volumes for easy comparison of annular sizing between phases
- Available annulus measurements include:
 - Minimum diameter
 - Maximum diameter
 - Perimeter
 - Area
 - Equiv. diameter based on perimeter and area

- Automatic valve plane to aorta angulation measurement


Valve Definition (cont.)

- Load multi-phase data to facilitate easy review and identifying preferred phase for annular measurements

- Access the multiphase volumes for easy comparison of annulus measurements across the heart cycle to evaluate differences across systole and diastole

- Classify degree of aortic calcification and include calcification volume estimates in your report easily with automatic volume & percentage measurements

- Specify tricuspid vs bicuspid valve anatomy to include in your comprehensive report
- BSA calculated and added to report automatically when height/weight included with input parameters
- Easily generate measurements for:
 - Sinotubular Junction (STJ)
 - Maximum diameter
 - D mean at 40mm above annulus
 - Left Ventricular Outflow Tract (LVOT)
 - Ostia Heights
 - Sinus diameters
 - Annulus Angle
 - Sinus Height

Many of these are generated from a single click


CONFIDENTIAL

SECRET DE
AFACERI

Features Detail (cont.)

Valve Definition (cont.)

- One click perpendicular views to demonstrate working angles for valve deployment in the cath lab.

- Measurements are auto populated into the summary table for easy viewing and export into a report to share with care team members.

TAVI

TAVI 2 MP

(31 Jul 2011 - 17:23:34)

Name: TAVI 2 MP

Exam ID:

Birth Date:

Exam Description: GATED AORTA

ID: AW94500904.639.1094270100

Modality: CT

Sex:

Series Slice Thickness: 0.625

Hospital Name:

Study Date: Jul 21 2011

Referring Physician:

kV: 120

Acquisition Number:

mA: Modulated

Tilt: 0

Aortic Annulus

Phase 70 %

Area:	420.4 mm ²	Max. Diameter:	26.4 mm
Mean Diameter:	23.1 mm	Min. Diameter:	21.0 mm
Perimeter:	73.8 mm	Eqiv. diam. (perim):	23.4 mm
Eqiv. diam. (area):	22.1 mm	Cath angles:	0 L, nan RAO, nan CAU

Annulus Angle

Phase 70 %

Angle:	56.8 °
--------	--------

Ostia and Sinus Heights

Phase 70 %

Left Ostium Height	Height:	14.3 mm
Right Ostium Height	Height:	21.5 mm
Sinus Height	Height:	24.0 mm

Sinus Of Valsalva

Phase 70 %

Sinus Diameter 1	Length:	28.2 mm
Sinus Diameter 2	Length:	29.3 mm
Sinus Diameter 3	Length:	32.7 mm

Aortic Measurements

Phase 70 %

Diameter at 40 mm			
Area:	588.1 ± 36.6 mm ²	Max. Diameter:	30.8 ± 0.8 mm
Mean Diameter:	28.6 ± 0.8 mm	Min. Diameter:	28.4 ± 0.8 mm
Location:	40.0 ± 0.8 mm		
Left Ventricular Outflow Tract			
Area:	444.7 ± 33.6 mm ²	Max. Diameter:	27.6 ± 0.9 mm
Mean Diameter:	25.8 ± 0.9 mm	Min. Diameter:	19.3 ± 0.9 mm
Location:	-5.0 ± 0.6 mm		
Maximum Diameter			
Max. Diameter:	34.2 ± 0.8 mm	Location:	51.8 ± 0.6 mm
Sinotubular Junction			
Min. Diameter:	30.8 ± 0.8 mm	Perimeter:	103.5 ± 0.6 mm
Location:	81.8 ± 0.8 mm		

The measurement precision depends on the acquisition parameters (slice thickness, image quality) and on the tortuosity of vessels.

1/6

- Copy/Paste summary table contents to other PC programs for additional flexibility in communicating results. (AW Server only)


CONFIDENTIAL

SECRET DE AFACERI

Aorta Analysis

Move directly into peripheral access planning without exiting or reselecting data.

- Automatically identifies peripheral vascular images within loaded series
- Automatic bone segmentation with automated vessel tracking

- Vessel tracking tools allow for easy access planning for all transcatheter approaches:
 - Subclavian Artery
 - Femoral Artery
 - Apical Approach
- Automatic VR display of transparent bones superimposed with vascular structures, provides visualization of external landmarks corresponding to measurement locations

Features Detail (cont.)


Aorta Analysis (cont.)

- Quantification tools capable of providing measurements for stenosis, aneurysms, 3D angles and distance
- Peripheral artery measurements automatically sent to summary table to facilitate exporting of results.
- Predefined measurements list provides a template for a standard order and location from which to place all measurements.
- Summary table converts to PDF format for a simple and concise way to share measurements and images with heart team for planning discussions.
- Automatic conversion of VR model into a variety of 3D printable formats plus automatic vessel hollowing capabilities³

Procedure Preparation

- Provides valve measurements, perpendicular view and gantry angulation in the cath lab to prepare for TAVI/TAVR procedure.
- Direct communication with HeartVision 2 software for easy transition from planning to guidance.

- 3D and calcium overlay VR models to aid in visualization during procedure.

Summary

TAVI Analysis is designed to provide an efficient and consistent approach for obtaining the various measurements that are required in preparation for a TAVI/TAVR procedure. It assimilates all of the acquired measurements and provides a variety of formats for exporting and sharing with the clinical care team.

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References:

1. Requires TAVI Enhanced license
*TAVI Enhanced may not be available in all countries or regions.
Please contact your sale associate.*
2. Based on GE Healthcare Internal Study
3. Requires 3D Suite license. 3D Suite is not intended for clinical use.

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
CT Perfusion 4D

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SECRET DE
AFACERI

Comprehensive perfusion analysis available for neurology, oncology and cardiology.

Easy-to-use, protocol-driven workflow for all organs.

Contents

1. Features..... 4

2. System Requirements..... 4

3. Regulatory Compliance..... 4

4. Indications for Use..... 4

5. Features Detail..... 5

5.01 Basic Features 5

5.02 Protocols..... 5

4. Functional Maps Available 9

Perfusion 4D has been expanded beyond stroke and oncology to now include dynamic myocardial perfusion. For any organ where blood perfusion is of interest Perfusion 4D is a complete package to quantitatively provide the necessary information to guide treatment decisions.

The protocol-driven design leads the user step-by-step through the process, reducing keystrokes and improving repeatability. Get the information you want quickly and reliably.


Overview

GE Healthcare Perfusion 4D software is a fast, easy-to-use automated software for analyzing CT Perfusion images related to stroke, tumor angiogenesis and dynamic myocardial perfusion. Its simple user interface and automated perfusion post-processing make it easy for you to diagnose quickly and accurately – and make treatment decisions more confidently.

GE HealthCare leverages its deconvolution expertise, with its innovative delay correction algorithm. Perfusion 4D takes perfusion assessment to the next level, with the addition of CT Dynamic Myocardial Perfusion. Providing the classic color overlays with quantitative measurements to evaluate myocardial blood perfusion.

Highlights

- Improved Image Quality of functional maps in presence of noise
- Improved ROI mirroring for brain protocols
- Deep Learning algorithm for ventricle removal
- Automated midline detection for brain symmetry
- Absolute and relative values with customizable thresholds for tissue classification
- Mismatch volume and ratio in brain stroke protocol
- Permits injection rates of 4cc/sec
- Dynamic registration
- Dynamic myocardial perfusion maps²

CONFIDENTIAL

SECRET DE
AFACERI

Features

- Visualize all information in true volumetric form with the ability to employ all the volumetric-based image analysis tools in the AW Volume Viewer.
- Deep Learning algorithm for brain ventricles removal
- Automated midline detection for brain symmetry
- Absolute and relative values with customizable thresholds for tissue classification
- Background processing will send processed images and functional perfusion maps and results in email¹ format to stroke team¹
- Whole heart perfusion² with semi-automatic guided workflow to aid in assessment of myocardial ischemia
- Smart Maps using the intelligent 4D noise suppression algorithm improves functional map image quality in the presence of noise.
- Dynamic non-rigid registration protocols for dynamic acquisitions within the body.
- Productivity enhancements include:
 - Ability to select automated or manual vessel detection for Brain Perfusion
 - Faster processing times for functional map computations
 - Volumetric visualization of functional maps
 - Time stamps on functional maps and Regions of interest (ROIs).
 - Interactive large vessels exclusion.
 - Access from anywhere using AW Server.
- Many basic features, including the ability to:
 - Process double phase acquisition data.
 - Create 4D Regions of Interest (ROIs).
 - Display ROIs statistical information.
 - Display averaged time intensity information of the ROIs.
 - Simultaneously review any number of functional maps.
 - Customize protocols
 - Save parameters
- Includes CT perfusion protocols for stroke, brain & body tumor, myocardium as well as dynamic registration
- Functional maps:
 - Regional Blood Volume (BV)
 - Regional Blood Flow (BF)
 - Regional Mean Transit Time (MTT)
 - Contrast arrival delay (IRF T0)
 - Transit time to peak of impulse residue function (TMax)

- Mean Slope of Increase
- Capillary permeability surface area for lengthened acquisition protocols (PS)
- Average image
- Base image
- Tissue Classification map with Mismatch volume and ratio from Modified and Low Perfusion regions.

System Requirements

Refer to AW Applications platform compatibility document.

Regulatory Compliance

This product complies with Regulation (EU) 2017/745 of the European Parliament and of the Council on medical devices (MDR). This product or its features may not be available in some other countries or regions. Please contact your sales associate.

Rx Only

Indications for Use

CT Perfusion 4D is an image analysis software package that allows the user to produce dynamic image data and to generate information with regard to changes in image intensity over time. It supports the analysis of CT Perfusion images (in the head and body) after the intravenous injection of contrast, in calculation of the various perfusion-related parameters (i.e. regional blood flow, regional blood volume, mean transit time and capillary permeability). The results are displayed in a user-friendly graphic format as parametric images.

This software will aid in the assessment of the extent and type of perfusion, blood volume and capillary permeability changes, which may be related to stroke or tumor angiogenesis and the treatment thereof.

CONFIDENTIAL

SECRET DE
AFACERI

Features Detail

Basic Features

- Process double phase acquisition data (data with minimal temporal sampling and data with longer temporal sampling in one acquisition to cover entire functional perfusion mechanism).
- Create 4D ROIs (multi-locations ROI with time-varying position).
- Display ROI statistical information (Area, Min., Max, Std Dev.).
- Display input ROI from several locations on the same graph.
- Display the averaged time-intensity information from a ROI (VOI).
- Simultaneously review any number of functional maps.
- Customize any viewport to series, graph or functional view.
- View series or functional data in axial, sagittal, coronal or oblique view.
- View series or functional data in 3D or VR mode.
- Customize protocols (advanced settings, color ramp and/or W/L used for each generated functional map) and to save these user preferences.
- Save the state of a protocol including parameters, ROIs, functional maps and pre/post-processing results for further reloading in CT Perfusion.
- Save all images in the database.
- Save ROIs statistical information using screensave.
- Print to film or color paper printer.
- Save parameters, used to generate functional maps, as DICOM objects and to restore them including ROIs with the corresponding maps in CT Perfusion.

Protocols

CT Perfusion 4D has fully integrated protocols that assist the user in obtaining results efficiently for neurology, oncology and cardiology care areas.

CT Perfusion 4D package protocols are:

CT Brain Stroke Protocol. The CT brain stroke protocol will open a loading panel containing information about the selected perfusion series in a tabulated format. It is possible to remove unwanted time points from the loading sequence.

Once the exam is loaded the software will perform the following algorithmic steps:

- **Image Registration** to minimize the effects of patient movement.
- **Bone Removal** to remove the cranium automatically from the dataset.


Protocols (cont.)

- **Deep Learning** brain ventricle segmentation to prevent ventricular matter inclusion in quantitative results and improve visual inspection of the maps
- **Arterial Input** used as the input for the deconvolution algorithm. Artery is automatically selected and can be easily adjusted if needed.
- **Venous output** for partial volume correction. Vein is automatically selected and can be easily adjusted if needed.
- **Last pre-enhancement / First post-enhancement** image are automatically selected to ascertain correct input for functional maps calculations.
- **Final Settings and Automated Computation.** The software will compute the functional maps with the saved settings, but can be made to stop and recompute the functional maps if changes to the algorithm settings were implemented. The changes may include any of the following: Algorithm Resolution, Smart Smoothing size, Maximum Blood Flow, Brain Density, Hematocrit Ratio.

The software will use deconvolution to generate all functional maps, Blood Flow, Blood Volume, Mean Transit Time, and Transit Time to IRF Peak (Tmax).

The Brain Stroke protocol has multiple review steps to guide the user through the functional information.

- **Review settings:** When loading the dynamic exam to process, this first review step is initiated. The general purpose of this step is to confirm the automatic selection of the input vessels for the brain stroke or tumor protocols, or to conduct a general review of the time resolved perfusion data and place ROIs to select vessels.


CONFIDENTIAL

SECRET DE
AFACERI


Features Detail *(cont.)*

Protocols *(cont.)*


- **Set Symmetry:** The software will automatically define the symmetry plane to be used for mirroring ROIs and relative thresholds. For the brain stroke protocol, the application provides a “Check View” image for a fast, intuitive way to indicate the plane and reference hemisphere identified.

If edition is needed a dedicated review step allows users to quickly adjust the symmetry axis.


- **Review Functional Maps** enables the viewing and assessment of functional maps such as Blood Flow, Blood Volume, Mean Transit Time, and Tmax. ROI mirroring functionality makes it easier to compare hemispheres.

Protocols *(cont.)*


- **Tissue Classification,** specific to the brain stroke protocol, enables the visualization of regions that are segmented from absolute or relative values, customizable thresholds and user selectable input maps.

Mismatch volume and ratio are calculated from the modified perfusion volume and low perfusion volume ROIs.

The Brain Stroke protocol is typically used for the assessment of cerebral hemodynamic disturbances e.g. stroke and traumatic brain injury.

- **Email¹,** when installed on AW Server and invoking the preprocessing for full automation of the Brain Stroke protocol, the processed images can be configured to auto email to the stroke team. Additionally, a new ‘save state’ series will be added to the database. Clinicians will be able to open the save state series to review the data and interact with it.

CONFIDENTIAL


SECRET DE AFACERI

Features Detail *(cont.)*


Protocols *(cont.)*

CT Brain Tumor Protocol.

Essentially identical to the brain stroke protocol but without the Tissue Classification step.

CT Body Tumor Protocol. The Body tumor protocol produces functional maps automatically once the user selects an artery using a simple and intuitive workflow.

The Body Tumor protocol also leverages most of the advanced algorithm settings that are part of the Brain Stroke protocol. The protocol has the following algorithms

- **Noise Threshold.** Thresholds can be customized to exclude air and bone from the calculations.
- **Artery Input.** System allows for an intuitive selection of the artery.
- **Last pre-enhancement / Last post-enhancement image / Last Second Phase Image** (last image of the phase 2 data). These can be selected to ascertain correct input into the parameter calculations.
- **Final Settings and Compute.** The user can compute immediately or customize the following advanced settings:
 - Algorithm Resolution, Smart Smoothing.
 - Maximum Blood Flow, Brain Density, Hematocrit Ratio.

The Compute command then engages the exclusive de-convolution algorithm to produce the parametric images.

CT Pancreas Perfusion Protocol. Identical to the Body Tumor Protocol but optimized for pancreas tumors.

CT Prostate Perfusion Protocol. Identical to the Body Tumor Protocol but optimized for prostate tumors.

CT Kidney Perfusion Protocol. Identical to the Body Tumor Protocol but optimized for kidney tumors.

Protocols *(cont.)*

CT Soft Tissue Perfusion Protocol. Identical to the Body Tumor Protocol, with an additional step to define the vein output for partial volume correction (automatic or manual venous selection).


CT Spleen Perfusion Protocol. Identical to the Body Tumor Protocol but optimized for spleen tumors.

CT Bone Perfusion Protocol. Identical to the Body Tumor Protocol but optimized for bone tumors.

CT Non-Rigid Registration. Fully automated protocols used for dynamic body acquisitions to register data for use in perfusion analysis. Provides simultaneous review of original data with the registered results for full confidence in the newly created series. Saved registered results are to be loaded into the appropriate body or cardiac perfusion protocol for functional analysis.

Two registration protocols are provided. CT Cardiac Dynamic Registration and CT Liver Dynamic Registration.

CT Liver Perfusion Protocol. Identical to the Body Tumor Protocol, with an additional step to define the portal vein input (automatic or manual venous selection).


CONFIDENTIAL

SECRET DE
AFACERI

Features Detail *(cont.)*

Protocols *(cont.)*


Processing times have improved exponentially, now with an average processing time of one to three minutes depending on slice thickness. Hepatic Arterial Fraction maps have been improved and can be an efficient tool for tumor detection. A new map, hepatic arterial blood flow (HABF), has been added to assist in tumor assessment.

Dynamic CT Myocardial Perfusion Protocol²


Complete guided workflow for assessing dynamic myocardial perfusion. The dynamic myocardial perfusion protocol is laid out in a series of three easy to follow steps. Start by selecting the series that has been generated from the non-rigid registration protocol and then follow the guided workflow:

1. **Cardiac Reformat** – images are automatically orientated into short and long axis views,
2. **Segmentation** – define the valve plane and the apex of the heart. The software will then calculate the endocardial and epicardial borders and provide a 17 segmented map of the myocardium along with a graphed view for each segment

Protocols *(cont.)*

3. **Functional Maps** – the third and final step produces the quantitative data with the color overlay and bullseye maps for blood flow, blood volume, mean slope of increase and mean transit time.

Quantitative data can be exported as a .csv file for a complete segmental and regional quantitative evaluation.

CONFIDENTIAL

SECRET DE
AFACER

Functional Maps Available

The key perfusion parameters that CT Perfusion 4D package offers (measured in their appropriate physiological units) are:

1. **Regional Blood Volume (BV; ml/100g).**
Characterizes the amount of blood present in a local region; this is affected by the size of the blood vessels and the number of open capillaries.
2. **Regional Blood Flow (BF; ml/min/100 g).**
It is a measure of the volume of blood that flows through a local tissue region within a period of time.
3. **Regional Mean Transit Time (MTT; sec).**
Mean Transit time is the time that blood requires to traverse the tissue circulation from arterial inlet(s) to the venous outlet(s).
Mathematically:
 $rMTT = rCBV / rCBF$ (Cerebral), or
 $rMTT = rTBV / rTBF$ (Tissue)
4. **Contrast Arrival Delay (IRF T0)**
It is the delay of arrival of contrast to any given location from its arrival at the artery ROI.
5. **Transit Time to IRF Peak (Tmax; sec)**
It is the transit time to the peak of the Impulse Residue Function (IRF).
Mathematically: $Tmax = IRF\ T0 + MTT/2$
6. **Capillary Permeability Surface Area Product (PS; ml/min/100 g)**
It is the unidirectional rate of transfer of blood solutes (contrast media) through the capillary endothelium into the interstitial space.
PS is given in the same units as blood flow, ml/min/100g.
Permeability Surface maps are intended to be used with lengthened acquisition times, typically 60 seconds or greater.
7. **Base** Image corresponds to an average of all pre-enhancement images. It is used in baseline subtraction and may also be as a reference image.
8. **Average** Image corresponds to the average of all images at a given slice location.

CONFIDENTIAL

SECRET DE
AFACERI

About GE HealthCare

GE HealthCare is a leading global medical technology, pharmaceutical diagnostics, and digital solutions innovator, dedicated to providing integrated solutions, services, and data analytics to make hospitals more efficient, clinicians more effective, therapies more precise, and patients healthier and happier. Serving patients and providers for more than 100 years, GE HealthCare is advancing personalized, connected, and compassionate care, while simplifying the patient's journey across the care pathway. Together our Imaging, Ultrasound, Patient Care Solutions, and Pharmaceutical Diagnostics businesses help improve patient care from prevention and screening, to diagnosis, treatment, therapy, and monitoring. We are an \$18 billion business with 51,000 employees working to create a world where healthcare has no limits.

Follow us on [Facebook](#), [LinkedIn](#), [Twitter](#), [Instagram](#) and [Insights](#) for the latest news, or visit our website [gehealthcare.com](https://www.gehealthcare.com) for more information.

References:

1. Requires Send by Email license and AW Server version with email service, not compatible with stand-alone AW systems.
*Send by Email may not be available in all countries or regions.
Please contact your sales associate.*
2. Requires Myocardial Perfusion license

Products mentioned in the material may be subject to government regulations and may not be available in all countries. Shipment and effective sale can only occur after approval from the regulator. Please check with local GE HealthCare representative for details.

©2024 GE HealthCare. GE is a trademark of General Electric Company used under trademark license.

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Follow us on [Facebook](#), [LinkedIn](#), [Twitter](#), [Instagram](#) and [Insights](#) for the latest news, or visit our website [gehealthcare.com](https://www.gehealthcare.com/) for more information.

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5934571-1EN
Revision 3

General User Documentation.

© 2023-2025 GE HealthCare. GE is a trademark of General Electric Company used under trademark license.
Reproduction and/or distribution is prohibited.

Blood Volume (BV)

Blood Volume (BV) is computed as the product of BF and MTT and is displayed in ml per 100g of wet tissue (ml/100g).

$$BV = BF \times MTT$$


Note: In practical implementation, BV computed in CT Brain Stroke and CT Brain Tumor protocols are different. In CT Brain Stroke, the BV is calculated as the product of BF with the MTT, which includes both intravascular and extravascular phases. In CT Brain Tumor, the BV is calculated as the product of BF with the intravascular MTT.

Permeability Surface Area Product (PS)

Permeability Surface area product (PS) is computed from the IRF and is displayed in ml per 100g of wet tissue per minute (ml/100g/min).

The PS has the same units as BF, as it quantifies the diffusion of some of the contrast agent into the interstitial space. It is used to assess the permeability of blood vessels in tumor angiogenesis (Brain Tumor and Body Tumor protocols) and in stroke (Brain Stroke protocol). In **CT Perfusion 4D**, the PS computation is included in Brain Stroke protocol in addition to Brain and Body Tumor protocols.

Permeability is related to the diffusion coefficient of the contrast agent through the pores of the capillary endothelium into the interstitial space due to the deficient or leaky blood brain barrier or in the presence of tumors. It has also been linked to the diffusion of contrast agent through the large holes of the sinusoidal capillaries. In the tissue IRF, the contrast agent diffusion is related to the extraction fraction E , or the fraction of contrast agent, which remains in the intravascular space after the initial IRF response and which then diffuses exponentially into the interstitial space. The E is related to the PS and BF as shown in the following equation. The ideal IRF and its relationship to the extraction fraction are depicted in Figure 4.

CONFIDENTIAL

SECRET DE
AFACERI

Tissue Classification

Tissue Classification is a segmentation tool, which assists in visualization of brain regions associated with the quantitative values of BV and Tmax (MTT or BF) as set by the user-specified thresholds. Tissue Classification review step is available in CT Brain Stroke protocol.


Tissue Classification map displays, in **red**, the Low Blood Volume region defined by the user-specified threshold for BV.

The Modified Perfusion region is defined by the user-specified Tmax, BF or MTT threshold (default is Tmax). It is used to determine the Mismatch region, which is displayed in **blue** on the Tissue Classification map.

The Mismatch region corresponds to the Modified Perfusion region minus the Low Blood Volume region. Both red and blue segmentations appear as overlays on the Average Image (gray scale image).

Note: In regions of extremely low contrast agent (BV and BF almost null), temporal parameters such as MTT and Tmax may not be calculated, which may result in black “holes”/undefined data in the temporal functional maps. As the value of these “black holes” voxels is undefined, they will not be included in regular threshold based segmentations. However, they will be included by default inside the Tissue Classification Low Perfusion region as these voxels of extremely low BV and BF values are likely indicating ischemic core.

Figure 5: Examples of perfusion maps generated by CT Perfusion 4D

- Please see **Review Settings** protocol description for more details on how to change artery and vein ROI locations.

6. Smart Smooth

- Smart smooth is applied to the dataset to reduce spatial noise while preserving the anatomical and functional details in the perfusion dataset prior to functional maps computation. The Smoothing Radius defines the area in voxels in the xy and z directions that the smoothing is applied. Smart Threshold controls the likelihood (propensity) of how often the smoothing process will be used to average different tissues together.

Smart smoothing parameters may be adjusted as desired. To adjust Smart smoothing parameters:

- Click on **Advanced Settings** button  in **CT Brain Stroke** control panel.
- Under Protocol Settings, select **Smoothing radius** from the list of parameters.
- Use the slider bar to adjust **Smoothing radius** and/or **Smart threshold** as desired.
- Click **Compute** to compute functional maps with new smoothing settings.
- Click **Compute** to re-compute functional maps.

7. Functional Maps Computation

- Functional maps are computed using both deconvolution and standard algorithms. The maps generated by the deconvolution algorithm are Blood Flow (BF), Blood Volume (BV), Meant Transit Time (MTT), Contrast Arrival Delay (IRF T0), Transit Time to Impulse Residue Peak (Tmax), and Permeability Surface Area Product (PS). The maps generated by the standard algorithm include Time-To-Peak (TTP), Positive Enhancement Integral (PEI), and Mean Slope of Increase (MSI). Additionally, an Average and Base image is created from the source data.

8. Automated Symmetry

- The cerebral symmetry plane is automatically detected and locked. It can be checked and modified in the Set Symmetry review step (see [Review Step 2: Set Symmetry](#)) to ensure correct symmetrical ROI placement (mirror ROI), as well as correct segmentations in the tissue classification step when using relative thresholds.

Note: More information on the adjustment of processing parameters is provided in the [Advanced Settings](#) chapter.

CONFIDENTIAL

SECRET DE
AFACERI

Working with Functional Maps

CT Brain Stroke protocol automatically computes the functional maps using *deconvolution* (if artery and vein input ROI's are identified successfully) and *standard* algorithms. For information on deconvolution and standard algorithms, and the detailed description of functional maps, please refer to **CT Perfusion 4D Algorithms** section of this manual.

Deconvolution Functional Maps

Blood Volume (BV), Blood Flow (BF), Mean Transit Time (MTT), Contrast Arrival Delay (IRF T0), Transit Time to IRF Peak (Tmax), Permeability Surface Area Product (PS), Average Image, Base Image

Standard Algorithm Functional Maps

Time-To-Peak (TTP), Positive Enhancement Integral (PEI), Mean Slope of Increase (MSI)

Displaying Functional Maps

In its default configuration, CT Brain Stroke protocol displays BF, BV, Tmax, and MTT maps.

To display other functional maps:

1. In any viewport, left-click on the active annotation in the top left corner displaying the image type, and open the on-view menu.

2. Select the functional map of your choice from the list of all computed functional maps. The selected map will replace the image displayed in the viewport.

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
GE HealthCare

FastStroke with StrokeSENS™

Automate the work up of ischemic stroke cases with the FastStroke automated workflow, including StrokeSENS ASPECTS and StrokeSENS LVO Detection, as well as auto email communication to the stroke team.

CONFIDENTIAL

SECRET DE
AFACERI

Contents

1. Features..... 4

2. Indications for Use..... 4

2.01 FastStroke..... 4

2.02 CT Perfusion 4D..... 4

2.03 StrokeSENS – Global (Non-USA)..... 4

2.04 StrokeSENS – USA..... 4

3. System Requirements..... 5

4. FastStroke and CT Perfusion 4D Regulatory Compliance 5

5. StrokeSENS Regulatory Compliance 5

6. Features Detail..... 6

6.01 Non-Contrast..... 6

6.02 ASPECTS..... 6

6.03 Carotids..... 6

6.04 StrokeSENS LVO Detection..... 6

6.05 Collaterals..... 6

6.06 ColorViz..... 7

6.07 Perfusion..... 7

6.08 Customizable Review..... 8

6.09 Communicate..... 8

Stroke is the 2nd leading cause of death globally.¹ Studies suggest that early intervention in stroke can lead to improved patient outcomes with a greater number of patients functionally independent.

CT is commonly used in the diagnosis of acute stroke. Due to the urgent nature of the disease, an improved workflow has become a critical component of acute stroke care. There is an increasing need for faster post-processing and workflow methods to reduce the time to diagnose, which may enable a faster time to treatment.

Recently there have been a number of clinical trials in the area of treating ischemic stroke with proximal cerebral artery occlusion that have called for new, faster workflows for the diagnosis of ischemic stroke patients.


Overview

FastStroke provides a comprehensive workflow solution for reviewing stroke workup images with exceptional flexibility, simplicity and performance. It is a streamlined approach that smartly adapts to your scan practices and allows you to review, post-process and share results in a streamlined, simple and consistent method. FastStroke can automatically process and send generated images to the stroke team via email. When working directly within the application the results are synchronized and displayed in a manner that enables an efficient review with high confidence. FastStroke provides a unique display, ColorViz, to aid in the visualization of the timing of collateral vessels during the mCTA series. FastStroke has full integration with CT Perfusion 4D³ to provide automated neuro perfusion analysis as part of the workflow. To facilitate communication, FastStroke can automatically process and send generated images via email to the stroke team.

FastStroke offers an option that allows integration to Circle Neurovascular's StrokeSENS⁴ processing tools for automated ASPECT Score and LVO detection.

Highlights

- 0-click workflow: Background processing of ischemic stroke work up
- Automated ASPECT Score, powered by StrokeSENS⁴
- Automated LVO detection, powered by StrokeSENS⁴
- Fully integrated with CT Perfusion 4D³ Neuro for visualization of perfusion functional maps
- Auto sends preprocessed images and functional maps in email format to stroke team²

Visit us (non-USA):

<https://www.gehealthcare.co.uk/products/advanced-visualization/all-applications/faststroke-with-strokesens>

Visit us (USA):

<https://www.gehealthcare.com/products/advanced-visualization/all-applications/faststroke-with-strokesens-lvo>

CONFIDENTIAL

SECRET DE
AFACERI

Features

- 0-click workflow: Background processing of ASPECTS⁴, LVO⁴ detection, mCTA and Perfusion³
- Auto sends preprocessed images, ASPECT Score, LVO finding, functional perfusion maps and results in email² format to stroke team
- Intelligent loading identifies the series type and applies the appropriate layout and protocol
- Intuitive workflow to quickly move through all acquired series
- Automatic synchronization of phases from multiphase CTA (mCTA)
- Smart layouts automatically adjust to display up to 6 mCTA phases simultaneously
- CTA images are automatically displayed in a thick 2D MIP at optimized WW/WL settings
- ColorViz, an intelligent color-coded display enabling easy and confident identification of vascular enhancement timing
- Fully integrated with Neuro Perfusion³ automated processing
 - Deep Learning algorithm for brain ventricles removal
 - Automated midline detection for brain symmetry
 - Absolute and relative values with customizable thresholds for tissue classification
 - Variety of functional maps selectable as inputs for tissue classification

Indications for Use

FastStroke

FastStroke is a CT image analysis software package that assists in the analysis and visualization of CT data derived from DICOM 3.0 compliant CT scans. FastStroke is intended for the purpose of displaying vasculature of the head and neck at different time points of enhancement. The software will assist the user by providing optimized display settings to enable fast review of the images in synchronized formats, aligning the display of the images to the order of the scans and linking together multiple groups of scans. In addition, the software fuses the vascular information from different time points into a single colorized view. This multiphase information can aid the physician in visualizing the presence or absence of collateral vessels in the brain. Collateral vessel information may aid the physician in the evaluation of stroke patients.

Indications for Use *(cont.)*

CT Perfusion 4D

CT Perfusion 4D is an image analysis software package that allows the user to produce dynamic image data and to generate information with regard to changes in image intensity over time. It supports the analysis of CT Perfusion images (in the head and body) after the intravenous injection of contrast, in calculation of the various perfusion-related parameters (i.e. regional blood flow, regional blood volume, mean transit time and capillary permeability). The results are displayed in a user-friendly graphic format as parametric images. This software will aid in the assessment of the extent and type of perfusion, volume and capillary permeability changes, which may be related to stroke or tumor angiogenesis and the treatment thereof.

StrokeSENS – Global (Non-USA)

StrokeSENS is a decision-aid software package to be used by clinicians to perform image processing, analysis, viewing and communication of computed tomography (CT) scans of the brain in patients with suspected acute stroke. Data and images are acquired through DICOM-compliant imaging devices prior to processing and analysis in StrokeSENS.

The StrokeSENS software provides analysis capabilities for imaging datasets acquired with standard CT imaging and contrast enhanced CT Angiography (CTA) modalities. Analysis of non-contrast CT images includes assessment of regions with suspected acute ischemic tissue. Analysis of contrast-enhanced CT images includes automated detection of anterior circulation Large Vessel Occlusion (LVO).

In the case of a suspected anterior circulation LVO, the system will send a notification to a pre-configured destination(s) (members of the acute stroke team), notifying them of the existence of a suspected LVO that requires review. The notification system is intended to be used in parallel to the standard of care workflow to notify clinicians of the existence of the case earlier than they may have been notified as part of the standard of care workflow. Images are available for viewing on a mobile device and on a standard radiology workstation. Images that are previewed on a mobile device are for informational purposes only and are not intended for diagnostic use beyond notification.

StrokeSENS – USA

StrokeSENS LVO is a radiological computer-aided triage and notification (CADt) software indicated for use in the analysis of CTA head images. The device is intended to assist hospital networks and trained clinicians in workflow triage by flagging and communication of suspected positive findings of Large Vessel Occlusion (LVO) in head CTA images.

Indications for Use *(cont.)*

StrokeSENS – USA *(cont.)*

StrokeSENS LVO uses a software algorithm to identify suspected LVO findings. In the case of a suspected LVO, the system will send a notification to a pre-configured destination(s), notifying the clinicians of the existence of a suspected LVO that requires review. The notification system is intended to be used in parallel to the standard of care workflow to notify clinicians of the existence of a potential LVO than being notified as part of the standard of care workflow.

Notifications may include a compressed preview of images.

Notifications are meant for informational purposes only and are not intended for diagnostic use beyond notification. The StrokeSENS LVO device does not alter the original medical image and is not intended to be used as a diagnostic device.

The results of StrokeSENS LVO are intended to be used in conjunction with other patient information and based on professional judgement, to assist with triage / prioritization of medical images. Notified clinicians are responsible for viewing full images per standard of care.

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FastStroke and CT Perfusion 4D Regulatory Compliance

FastStroke and CT Perfusion 4D comply with Regulation (EU) 2017/745 of the European Parliament and of the Council on medical devices (MDR). This product or its features may not be available in some other countries or regions. Please contact your sales associate.

Rx Only

StrokeSENS Regulatory Compliance

FStrokeSENS complies with the European Council Directive 93/42/EEC Medical Device Directive as amended by European Council Directive 2007/47/EC.

This product or its features may not be available in some other countries or regions. StrokeSENS ASPECTS is not available in USA. Please contact your sales associate.

Rx Only

CONFIDENTIAL


SECRET DE
AFACERI

Features Detail

The FastStroke application can recognize incoming datasets specified as part of the ischemic stroke protocol.

It will automatically initiate background processing of the incoming datasets and save the processed images as part of the study. When leveraging FastStroke through the AW Server, the auto processed images and maps can also auto export to email to a pre-defined distribution list.

All is intended to expedite the review and communication of the CT acquisition amongst the care team members, as well as simplify the loading, display and review of the multiple series that are commonly acquired as part of the acute stroke workup in CT.


Non-Contrast

The first step in the FastStroke application allows the user to review the non-contrast head CT, automatically displayed in standard reformatted planes.

ASPECTS⁴

Based on the Alberta Stroke Program Early CT Score (ASPECTS), StrokeSENS uses validated machine-learning algorithms to automatically identify regions of the brain affected by early ischemic changes to aid physicians in assessment of the severity and extent of stroke damage. ASPECTS is intended for use in patients with ICA/MCA occlusions.

- Extract more from simple non-contrast CT
- Automated scoring system augments physician assessment
- Scores and visual aids available in under 5 minutes on average

Carotids


The second step of the application displays the head and neck CT Angiogram in reformatted views, automatically set to a thick slap MIP and at optimized WW/WL for CTA review.

Slice thickness is easily adjustable and Volume Viewer tools are compatible to facilitate quick vessel analysis.

StrokeSENS LVO Detection⁴


StrokeSENS LVO automatically detects suspected LVO cases within minutes of receiving the CTA images, supporting early engagement of the stroke team.

In the case a LVO is detected, “Suspected LVO” is written onto the input DICOM image (CTA or 1st phase of mCTA) and results are automatically included in the email notification, helping notify the stroke team of the time-sensitive case.

Collaterals

The Collaterals step uses a smart layout feature that adjusts the layout to accommodate up to 6 mCTA phases simultaneously. All the phases will automatically be synchronized and displayed in a thick 2D MIP at optimized WW/WL settings, to help the user with the review and comparison.

CONFIDENTIAL

SECRET DE
AFACERI


Features Detail *(cont.)*

Collaterals *(cont.)*

The Collaterals review step provides a single viewing screen to see all passes and assist in identifying the location of lesions and the immediate state of the collateral flow.

ColorViz


ColorViz allows for a quick and easy assessment of the collaterals by gathering the info from all the mCTA phases in one single colored rendering.

Perfusion³

The CT Perfusion 4D Brain Stroke protocol is fully integrated in FastStroke, allowing you to quickly display functional maps for a complete review of the CT stroke workup.

The CT Perfusion 4D Brain Stroke protocol provides a fully automated background processing capability. Functional maps can be automatically generated and shared through an email communication with stroke team care members.


Perfusion³ *(cont.)*

The software uses deconvolution to generate all functional maps, Blood Flow, Blood Volume, Mean Transit Time, and Transit Time to IRF Peak (T_{max}).

CT Perfusion 4D Stroke protocol allows the evaluation of dynamic CT data following an injection of contrast material, generating information regarding changes in image intensity of over time.

Key features for the Brain Stroke protocol include automated:

- **Image Registration** to minimize the effects of patient movement
- **Bone Removal** to remove the cranium automatically from the dataset
- **Deep Learning** brain ventricle segmentation to prevent ventricular matter inclusion in quantitative results and improve visual inspection of the maps
- **Arterial Input** used as the input for the deconvolution algorithm. Artery is automatically selected and can be easily adjusted if needed.
- **Venous output** for partial volume correction. Vein is automatically selected and can be easily adjusted if needed.
- **Last pre-enhancement / First post-enhancement image** are automatically selected to ascertain correct input for functional maps calculations.
- **Automated computation** of the functional maps based on the saved settings.
- **Automated brain symmetry plane** detection to generate relative values comparing each hemisphere.
- **Tissue Classification** map segmented from absolute or relative values, customizable thresholds and user selectable input maps.
- **Mismatch** volume and ratio calculated from the Modified Perfusion region and the Low Perfusion region.
- **Email²**, when working on AW Server and invoking the preprocessing for full automation the processed images can be configured to auto send to the stroke team in email format.

CONFIDENTIAL

SECRET DE
AFACERI


Features Detail *(cont.)*

Customizable Review

Create layouts that match your reading style and save them as protocols for future recall.

Communicate

Share results with the stroke team immediately with the automated email generated once background processing has completed. Additionally, a new 'save state' series will be added to the database. Clinicians will be able to open the save state series within FastStroke to review the data and interact with it.

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References:

1. Mozaffarian D, et al; on behalf of the American Heart Association Statistics Committee and Stroke Statistics Subcommittee. Heart disease and stroke statistics—2016 update: a report from the American Heart Association. *Circulation*. 2016;133:000-000.
2. Requires Send by Email license and AW Server version with email service, not compatible with stand-alone AW systems.
Send by Email may not be available in all countries or regions.
Please contact your sales associate.
3. CT Perfusion 4D Neuro License is a pre-requisite for Neuro Perfusion map calculation within FastStroke.
4. StrokeSENS™ is legally manufactured by Circle Neurovascular Imaging, Inc.
Requires StrokeSENS™ licenses for ASPECTS and LVO. Not available in all countries or regions.
StrokeSENS™ ASPECTS is not available for sale in USA.
Please contact your sales associate.

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
CONFIDENTIAL

SECRET DE
AFACERI

Part of a comprehensive acute neuro workup —
Semi-automated segmentation and assessment of
Hematoma and Aneurysms in CT images.

Contents

1. Features..... 4

2. System Requirements..... 4

3. Regulatory Compliance..... 4

4. Indications for Use..... 4

5. Features Detail..... 4

5.01 Hematoma..... 4

5.02 Aneurysm..... 5

CONFIDENTIAL

SECRET DE
AFACERI

Stroke and trauma to the head may result in intra cerebral or subdural hematomas. Patients under suspicion of such internal bleeding are triaged using CT imaging. It is important to determine both the presence as well as the amount of such bleeding to properly triage a patient for treatment or follow-up. Bleeding in the brain may also be caused by the rupture of aneurysms. It is important to detect the presence of the aneurysm and assess its potential hazard.

Overview

Stroke VCAR (Volume Computer Assisted Reading) gives you a complete reading workflow solution for a comprehensive and robust analysis of hematoma and aneurysms.

Hematoma assessment is done using semi-automated segmentation tools married with innovative interactive editing capability in the form of SmartMesh.

Aneurysm assessment is done through an innovative user guided aneurysm segmentation and visualization.

The program lets you generate a clear, concise clinical summary for sharing with referring physicians.

Highlights

- The program allows you to quickly segment the hematomas in a non-contrast CT
- Introducing SmartMesh: an intuitive volumetric editing tool
- Intelligent user guided aneurysm segmentation
- Fully integrated with Spectral CT when used with GSI datasets

CONFIDENTIAL

SECRET DE
AFACERI

Visit us:

<https://www.gehealthcare.com/products/advanced-visualization/all-applications/stroke-vcar>

Features

- Semi-automated segmentation and sizing of hematoma.
- Introducing SmartMesh—an interactive volumetric editing tool
- Comprehensive reporting of segmented hematomas
- Fully integrated with Spectral CT when using GSI data
- Track hematoma changes both visually and quantitatively with longitudinal exam workflow
- Simple user driven aneurysm detection workflow
- Easy and intuitive workflow for aneurysm segmentation. Initial segmentation done with four guided clicks
- A selection of application specific tools ensures ease of use
- Create layouts to match your reading style

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Rx Only

Indications for Use

Stroke VCAR is a CT image analysis software package that provides information to physicians to assist them in the analysis and visualization of Brain CT data derived from DICOM 3.0 compliant CT scans. Stroke VCAR is designed for the purpose of segmenting and assessing intracerebral and intracranial hemorrhages in the brain using semi-automated tools on non-contrast CT exams. Additionally Stroke VCAR provides a set of workflow tools for the segmentation and visualization of aneurysms in the brain from contrast enhanced CT exams. It is intended for use by clinicians to process, review, archive, print and distribute CT studies.

This software will assist the user by providing initial 3D segmentation, measurements and visualization of hemorrhages and aneurysm in the brain. The user has the ability to adjust, review and has to confirm the final segmentation.

CONFIDENTIAL


SECRET DE AFACERI

Features Detail

Hematoma


Segmentation

The Stroke VCAR hematoma tool uses an intelligent user guided algorithm for semi-automated segmentation of all types of hematomas in the brain. These segmentations can be converted into volume of interest to measure and assess their density.

SmartMesh editing

In cases when the semi-automated segmentation does not fully meet a user's expectation, SmartMesh allows for innovative editing capability to quickly complete the segmentation. A novel feature of the editing tool is to provide users with immediate visual feedback on what is included or excluded from the edition.

Features Detail (cont.)

Hematoma (cont.)

Reporting

Simple and intuitive reporting enables users to confidently manage and share hematoma findings in a consistent and standardized manner.

Follow up


Load repeat scans in Stroke VCAR Hematoma protocol and use its longitudinal analysis capability to track changes, both visually and quantitatively, due to progression or therapy response.

Aneurysm

Detection Workflow


A simple and intuitive workflow that helps users with a fast and easy assessment for the detection and treatment planning of aneurysms.


Aneurysm (cont.)

Segmentation

A simple guided workflow for the segmentation of an aneurysm. Fast initial segmentation is done with four guided clicks.

Define visualization planes to better see how the aneurysm is connected to the feeding vessel. Once the segmentation is complete, the user is able to interactively review 3D images of the segmented aneurysm and measurements.

Reporting

Simple and intuitive reporting enables users to share images and measurements of aneurysms in a consistent and standardized manner.

Summary

Stroke VCAR provides workflow benefits using semi-automated hematoma and aneurysm segmentation tools enabling a standardized approach to your facilities review & analysis of these exams.

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CONFIDENTIAL

SECRET DE
AFACERI

Automated spine labeling with curved reformat generation to improve ease of reading and reporting

Contents

1. Features..... 4

2. System Requirements..... 4

3. Regulatory Compliance..... 4

4. Indications for Use..... 4

5. Features Detail..... 4

5.01 Integration into VV..... 4

5.02 Automatic labeling 4

5.03 Curved Reformats 5

5.04 Summary..... 5

CONFIDENTIAL

SECRET DE
AFACERI

Approximately six percent of the total CT exam procedure volume is dedicated spine imaging. Add to that the number of other routine procedures that include the spine and are referenced within an imaging report you will find a significant number of imaging series generated every day for every scanner where identification and review of the vertebral bodies and their discs are necessary. Bone VCAR will help simplify the reading experience and improve reporting efficiency by identifying and displaying the anatomy in an automated yet flexible implementation.

Overview

Bone VCAR provides an automated spine labeling application based on a deep learning model.


It can identify and label segments of the spine or the entire spine in a matter of seconds.

Additionally, let Bone VCAR generate the curved reformats of the spine along with the perpendicular oblique views to easily see the true cross section of the vertebral bodies and disc spaces. Utilize these time saving features to assist in your reading experience.

Bone VCAR can also improve the reporting efficiency for identifying and remembering key landmarks that are included in the report dictation.

Highlights

- Automated spine labeling
- Automated curved reformats with oblique views parallel to the bodies and disc space
- Accessible for any exam type; trauma, oncology, dedicated spine, general imaging
- Accessible with all workflows

CONFIDENTIAL

SECRET DE
AFACERI

Visit us:

<https://www.gehealthcare.com/products/computed-tomography/advanced-visualization/bone-vcar>

Features

- Automated vertebrae labeling in 5 seconds or less for volumes of 300mm or less*
- Deep learning application trained on global datasets acquired with a broad range of acquisition parameters achieving >90% labeling accuracy
- Automated generation of a 3D trace to generate curved reformat views
- Automated generation of oblique views perpendicular to vertebral bodies and disc spaces
- Compatible with axial source images as well as previously generated reformatted series
- Manually insert, edit or delete vertebral body labels as necessary
- Compatible with data acquisitions for dedicated spine work, trauma, oncology and routine imaging
- Load multiple series and vertebrae labeling will propagate to the different volumes
- Easily review multiple reconstructions with simultaneous display of different kernels
- Access from the Volume Viewer tool bar to easily apply during your standard read at any time
- Compatible with all DICOM CT data that meets image input requirements as described in the user guide

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Rx Only

Indications for Use

Bone VCAR is a post processing application for use in the analysis of CT images. The software is intended to support clinicians in the review of images that include the spine by providing tools to label the spine and optimize the display of anatomy within the CT image.

Bone VCAR is designed to support the clinician in visualization of the spine, by providing initial identification of vertebrae to assist in report dictation.

Indications for Use *(cont.)*

The software also assists the user by providing optimized display settings for easier identification of anatomy to facilitate fast image review and reporting of findings. Bone VCAR may be used for multiple care areas and is not specific to any disease state. It can be utilized during the review of various types of exams including trauma, oncology, and routine body.

Features Detail

Bone VCAR is designed to support the clinician in easy visualization and identification of vertebral bodies.

A refined layout including reformats through the vertebral bodies and disc spaces enables fast review of the anatomy to aid in the reading and reporting experience.

Integration into VV

- Bone VCAR is integrated into the core toolset of Volume Viewer which makes it available to use in any reading review style desired.

- Compatible with reformatted series, Volume Rendered displays as well as GSI monochromatic datasets
- View labeled anatomy on two different reconstructions at the same time

Automatic labeling

- Bone VCAR, when activated, will automatically label all the vertebral bodies within the volume for slice thickness & spacing of 2.5mm or less

CONFIDENTIAL

SECRET DE
AFACERI

Features Detail *(cont.)*

Automatic labeling *(Cont.)*


- Deep learning design, built with diverse database of scanner types, variable reconstruction parameters, scan ranges, contrast and non-contrast exams
- Easily edit, when needed, with auto propagation above and below the edited location
- Manually insert vertebrae and add label names to account for anatomical variations
- Load multiple series from the same acquisition and labels will propagate across the volumes

Curved Reformats

- Curved reformats are automatically created and displayed in a review layout to allow visualization of standard and reformatted views simultaneously


- Oblique view perpendicular to the centerline to show easily the true cross section of the vertebral body and the disc space


Curved Reformats *(Cont.)*

- Easily edit the centerline to ensure accurate placement in the most challenging cases

- Streamline the reading workflow by loading multiple series for simultaneous display with labels and curved reformat and sagittal views in a single layout

- Create batch reformats with the vertebrae labels and export to PACs to provide easy identification of vertebra bodies for expanded clinical team

Summary

Bone VCAR is designed to provide an efficient and reliable method for labeling vertebral bodies and generating true sagittal, coronal and axial views via a 3D centerline detection through the center of the vertebral bodies. It is available for use with any DICOM CT datasets that meet the input requirements as defined in the user guide.

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Reference:

* Trademark of General Electric Company

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SECRET DE
AFACERI

Highlights

- Automated spine labeling
- Automated curved reformats with oblique views parallel to the bodies and disc space
- Accessible for any exam type; trauma, oncology, dedicated spine, general imaging
- Accessible with all workflows

A 3D reconstruction of a human spine, showing vertebrae and intervertebral discs. The spine is curved, and the vertebrae are labeled with their corresponding T12, L1, L2, L3, L4, and L5 vertebrae. The image is a grayscale 3D model, likely a CT scan reconstruction.

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SECRET DE
AFACERI

Thoracic VCAR

Comprehensive evaluation of CT lung exams

CONFIDENTIAL

SECRET DE
AFACERI

Contents

1. Features..... 4

2. System requirements 4

3. Regulatory compliance 4

4. Indications for use 4

5. Intended use 4

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SECRET DE
AFACERI

Thoracic VCAR provides tools to help efficiently extract more information from your CT Chest acquisitions when evaluating a variety of pulmonary disease processes. The Thoracic VCAR Analysis protocol allows you to evaluate abnormal Hounsfield ranges within parenchyma and provide detailed quantitative results of airway sizing to referring physicians which may be useful in assessing various thoracic diseases. The GSI Pulmonary Perfusion feature allows you to take advantage of GSI acquisition techniques to evaluate the distribution of iodine throughout the lungs to help identify variations of relative perfusion distribution. Thoracic diseases that may be associated with changes in perfusion include Pulmonary Embolism (PE) or Chronic Obstructive Pulmonary Disease (COPD). Thoracic VCAR provides the visual and quantitative outputs to enhance your imaging reports with details that can help referring physicians in the assessment of thoracic disease.

Overview

Thoracic VCAR gives you the ability to make quantitative measurements of the lungs to aid in the diagnosis of lung diseases. The application blends automated lung and airway segmentation with basic review and advanced lung tissue analysis. Utilize GSI data to extract useful information regarding relative iodine perfusion within the lungs to assist in identification of relative hypo-perfusion for each individual patient. Using Thoracic VCAR you can generate a clear, concise report that communicates relevant medical information to referring physicians and patients.

Highlights

- Automated segmentation of lungs and airways.
- One-touch 3D airway tracking.
- Automated measurements of wall thickness for intuitive airway analysis.
- Lobe segmentation and visualization.
- Automatic analysis of GSI pulmonary studies.
- Display relative perfusion map based on GSI iodine concentration in lung parenchyma.

CONFIDENTIAL

SECRET DE
AFACERI

Visit us:
<https://www.gehealthcare.com/en/products/advanced-visualization/all-applications/thoracic-vcar>

Features

- Basic review provides an overall view of the lungs immediately after application launch.
- Parenchyma analysis tool automatically segments lungs and airways.
- Lobe segmentation provides measurements related to the individual lung lobes.
- Airway analysis automatically segments and tracks the airways and displays the results for analysis.
- Selected airways can be further visualized in lumen view, providing you with additional measurements.
- Compare mode allows you to load two exams for disease progress comparison.
- GSI Pulmonary Perfusion protocol utilizes GSI acquisitions to display differences of iodine attenuation within the lungs.
- Efficiently identify and report regions of relative hypo-perfusion with the GSI Pulmonary Perfusion feature.

System requirements

- AW Workstations AW 4.7 Ext. 12 with VV14.0 Ext 6 or higher
- AW Server 3.2 or higher

GSI Pulmonary Perfusion is an optional add-on feature to the Thoracic VCAR application. GSI Viewer and Thoracic VCAR are pre-requisites for GSI Pulmonary Perfusion

- Color Monitors
- Single or Dual Display Monitors for AW Server in landscape or portrait orientations

(recommended monitor resolution is up to dual 2 MP (1600 x 1200) or a single 3 MP (1536 x 2048))

Regulatory compliance

This product complies with the European Council Directive 93/42/EEC Medical Device Directive as amended by European Council Directive 2007/47/EC. This product or its features may not be available in some other countries or regions. Please contact your sales associate.

Rx Only

Indications for use

Thoracic VCAR is a CT, non-invasive image analysis software package, which may be used in conjunction with CT lung images to aid in the assessment of thoracic disease. The software provides automatic segmentation of the lungs and automatic segmentation and tracking of the airway tree. Thoracic VCAR also provides quantification of Hounsfield units and display by color of thresholds within a segmented region. Thoracic VCAR also supports Gemstone Spectral Imaging (GSI) acquisitions for the evaluation of pulmonary perfusion. It provides additional information to aid in visualization of variations of perfusion within the lungs and to quantitatively assess lung volumes. It is intended to be used as an adjunct to current standard methods utilizing color coded displays of iodine attenuation differences in the lungs to aid in identifying segments of relative perfusion differences which may be useful in assessing thoracic disease. Thoracic diseases that may be associated with changes in perfusion include pulmonary embolism and COPD.

Intended use

Thoracic VCAR is a CT, non-invasive image analysis software package which may be used in conjunction with CT lung images to aid in the assessment of thoracic disease diagnosis and management.

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GE HealthCare

Hepatic VCAR

Providing automated segmentation and assessment of liver, liver lesions, and vasculature; enabling an efficient and consistent workflow.

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
Liver disease is a major cause of morbidity and mortality globally. Liver cancer is currently the 3rd most common primary cancer, as well as a common site for metastatic disease. Detection and treatment monitoring of liver lesions requires measurement and change monitoring tools. In addition, planning for surgical resection requires the determination of the total liver, liver lobe or a liver segment volume.

Overview

Hepatic VCAR provides a complete reading workflow solution for identifying liver lesions, obtaining liver volumes, associating vasculature to corresponding lesions and sharing the results with surgeons and referring clinicians. Hepatic VCAR will identify and organize multi-phase abdominal CT to ease the reading process. Hepatic VCAR will assist with visualizing and measuring liver, liver segments and liver lesions and associate tumor burden according to each of the segments. Full integration with GE Spectral CT allows for use with GSI data to further extract material density information within the application. Generate images and quantitative measurements to build your clinical report with all the relevant information for your referring physicians to act upon.

Highlights

- Deep learning liver segmentation algorithm¹
- Deep learning hepatic artery segmentation algorithm¹
- Guided workflow
- Intuitive tools to generate liver segments.
- Integration with Spectral CT allows for quantification of Iodine to aid in lesion characterization when used with GSI datasets.
- Segmentation for arterial and venous phases

CONFIDENTIAL

SECRET DE
AFACERI

Visit us:
[https://www.gehealthcare.com/products/advanced-visualization/
all-applications/hepatic-vcar](https://www.gehealthcare.com/products/advanced-visualization/all-applications/hepatic-vcar)

Features

- Loads multi-phase abdominal CT studies with automated detection & organization of contrast enhancement phases.
- Deep learning liver segmentation in under 3 sec* and achieves >95% segmentation success rate.
- Intelligent lesion segmentation using auto contour tool.
- Built-in lesions overlap detection and avoidance.
- Arterial and venous vasculature segmentation to assist in lesion management planning.
- Deep learning algorithm for hepatic artery segmentation.
- Tumor burden calculations linked to segment, lobe or whole liver.
- Efficient management of lesions and tumor burden for longitudinal exams.
- Intuitive editing tools for quick and easy refinements when needed.
- Volume rendered visual depiction of the lesion and liver along with liver lobes, liver segments, portal vein and hepatic artery for longitudinal comparison.
- Portal Vein based liver segmentation capability with color coding of segments.
- A selection of application-specific tools ensures ease of use and thorough exams.
- Easily copy/paste summary table measurements for easy reporting, tracking and communicating with clinicians (when working on AW Server).
- Create layouts to match your reading style.
- Your workflow benefits from fast, interactive electronic patient reporting.

Indications for use

Hepatic VCAR is a CT image analysis software package that allows the analysis and visualization of Liver CT data derived from DICOM 3.0 compliant CT scans. Hepatic VCAR is designed for the purpose of assessing liver morphology, including liver lesion, provided the lesion has different CT appearance from surrounding liver tissue; and its change over time through automated tools for liver, liver lobe, liver segments and liver lesion segmentation and measurement. It is intended for use by clinicians to process, review, archive, print and distribute liver CT studies. This software will assist the user by providing initial 3D segmentation, vessel analysis, visualization, and quantitative analysis of liver anatomy. The user has the ability to adjust the contour and confirm the final segmentation.

System requirements

Hepatic VCAR option can be installed on GE's AW workstation 4.7 Ext. 12 and AW Server 3.2 Ext. 3 (and later)

Compatible Platforms:

- AW: Z800, Z820**, Z440 (and later)
- AW Server 3.2 Ext. 3 (and later) with LT ML 350 p G8** (and later)

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SECRET DE
AFACERI


Colon VCAR EC

A comprehensive reading workflow solution
for colonic lesion detection

SECRET DE
AFACERI

CONFIDENTIAL

Contents

1. Features.....	4
2. Indications for use	4
3. System requirements	4
4. Regulatory compliance	4
5. Features detail	5
5.01 Electronic cleansing	5
5.02 User interface	5
5.03 Simultaneous, Indexed viewing	5
5.04 360-degree dissection	5
5.05 Joystick navigation	6
5.06 Computer aided reading	6
5.07 Virtual biopsy	6
5.08 Auto-segmentation	6
5.09 Dissection (centerline) contouring verification	6
5.10 Automatic high resolution flythrough	6
5.11 Save state	7
5.12 Prone/supine polyp linking	7
5.13 Synchronized prone/Supine 3D dissection view navigation	7
5.14 Small bowel extraction	7
5.15 Analysis mode	7
5.16 Color centerline tracking	7
5.17 Patient report	7
5.18 Measurement tools	7
5.19 Archive	7
5.20 Customizable review and analysis layout	7
5.21 Movie generator	8
5.22 Entire colon movie	8
5.23 Additional features.....	8
5.24 Summary	8

SECRET DE
AFACERI

CONFIDENTIAL

Clinical and technical background


Colorectal cancer is a major health concern for both men and women. Recent studies conclude that early detection of colonic polyps can reduce the mortality rate of colon cancer. You need a protocol that lets you quickly and accurately detect and report results of CT colon scans; a program that optimizes your workflow and improves your reading productivity.

Overview

Colon VCAR EC gives you a comprehensive reading workflow solution for detecting colonic lesions with flexibility and efficient performance. With it you can visualize anatomy that would otherwise be obscured by tagged fluid. Plus digital contrast agent filter highlights potentially suspicious regions.

Highlights

- The program allows you to read and problem solve in 2D, 3D, or 360 degree dissection views.
- Prone/Supine polyp linking provides quick navigation between datasets.
- Digital Contrast Agent (DCA) provides automatic visualization of shapes characteristic of polyps.
- Electronic Cleansing (EC) subtracts tagged stool and fluid, making it easier to identify lesions.

SECRET DE
AFACERI
CONFIDENTIAL


Visit us:


<https://www.gehealthcare.se/products/advanced-visualization/all-applications/colon-vcar>

Features

- Prone and supine synchronized reviewing for problem solving.
- You synchronize prone and supine navigation based on polyp bookmark location and distance from center line.
- Provides you with quick, easy to use polyp and center line measurement tools for complete reporting.
- You can leverage any combination of synchronized 2D, 3D, and 360 views.
- Dissection views make primary reading and problem solving more efficient.

- You can vary navigation speed.
- Lets you track entire colon from rectum to cecum for thorough, precise exams.
- You are provided with unrestricted viewing angles of all colon surfaces.
- You can apply color to bookmarked lesions for better visualization.
- The program allows you to assess luminal, mural and extracolonic anatomy and pelvic structures.

- A selection of application-specific tools ensures ease of use and thorough exams.

Indications for use

Colon VCAR EC is a CT, non-invasive, image analysis software package that allows the visualization of 2D, 3D and dissected medical image data of the colon derived from DICOM 3.0 compliant CT scans. Colon VCAR is designed to aid the physician in evaluating the lumen and internal wall of the colon to confirm the presence or absence of colonic lesions (e.g. polyps). It provides functionality for 2D/3D rendering, bookmarking of suspected lesions, synchronized viewing of the 2D, 3D and 360 dissection views for data sets acquired in any position, and an object oriented endoluminal display. In comparison to Colonoscopy, this tool has an advantage of depth penetration due to its 3D presentation capability.

It is intended for use by clinicians to process, render, review, archive, print and distribute colon image studies.

The Colon VCAR DCA (Digital Contrast Agent) module is an automated highlight feature for the visual identification of spherical structures within the colon and is intended to be used as concurrent reading device. Digital Contrast Agent is a 3D filter that produces images that highlight spherical anatomical regions, such as polyps, and/or stool. Colon VCAR uses color to display these highlighted spheres. Images are made available to the physician to aid in characterization of potential polyps and thus, the patient management care decision process.

System requirements

Colon VCAR EC can be installed on AW Workstations AW 4.7 Ext. 14 or later and AW Server 3.2 Ext 3.2 or higher.

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This product has not been approved for sale in the United States by the United States Food and Drug Administration (FDA).

Features detail

Electronic cleansing



Electronic Cleansing, accessible in both 2D and 3D review, subtracts tagged stool and fluid from the generated data making it easier to identify lesions that may be visible in only one view.

Tagged material is identified based on Hounsfield units. Processing tasks such as Auto Segmentation, 3D flythrough, and Dissection also take advantage of tagged material subtraction.

Auto Segmentation shows the tagged material while 3D flythrough enables quick anatomy visualization without the tagged material. With this method, the tagged material is removed yet shown in a different color at an adjustable opacity level.

Similarly, the tagged material can be removed from the dissection view as seen in the images below.

User interface

The Colon VCAR EC user interface has been optimized to facilitate fast, accurate examination of the colon's full extent. It provides seamless integration between automatic and manual interactions for a level of flexibility designed to reduce reading times.

Simultaneous, Indexed viewing

The program can display both prone and supine views, which can be manually synchronized for simultaneous and indexed viewing. During review, you can "Lock to Target" with a single button click to quickly problem solve during the primary read. And you are able to bookmark locations of suspected polyps.

360-degree dissection


This powerful tool provides a 360-degree view of the entire internal large bowel simultaneously displayed in supine and prone positions. This ability speeds viewing and analysis. This feature lets you leverage all available views 2D, 3D, and 360-degree dissection during the primary read. For more traditional reading methods, a 2D, 3D, or Dissection view can be laid out and reserved for primary reading with reference to other views for problem solving.

SECRET DE
AFACERI
CONFIDENTIAL

Features detail (cont.)

Joystick navigation

Colon VCAR EC lets you adjust the speed and direction of a flythrough simply by using the mouse. Flythroughs can also be set on autopilot for hands-off reviewing. Additionally, a single mouse click lets you lock onto a region of interest and fly around the target anomaly.


Computer aided reading

Colon VCAR EC's Contrast Agent (DCA), a 3D filter, highlights spherical anatomy such as polyps in color. The resulting images facilitate the characterization of potential polyps.

Virtual biopsy

Problem solving while 3D reading is facilitated by the Virtual biopsy tool, which allows you to view areas in the navigator in color-coded or black/white scales. Different densities, determined by Hounsfield units, can be displayed using the color scale.

SECRET DE
AFACERI

CONFIDENTIAL

Auto-segmentation


Colon VCAR EC quickly performs the time-intensive task of colon segmentation for both prone and supine data sets. All you do is select the protocol.

The program automatically defines a luminal centerline in both prone and supine positions. This can be applied to a high percentage of cases, depending on the level of insufflation and collapse. Semi-automatic and manual segmentation modes are also available.

Semi-automatic segmentation provides complementary, and often necessary, assistance in cases when the colon is collapsed and operator intervention is required. This mode is typically used to connect colon segments around a collapse to form a single centerline. Manual segmentation is typically used to extend the centerline to the end of the cecum for completeness. Manual editing is also available to refine automatic tracking if needed.

Dissection (centerline) contouring verification

Using a color mapping technique, Dissection Contouring Verification outlines colon anatomy included in the centerline trace. You can review verification in any orthogonal plane.

Automatic high resolution flythrough

Real-time 512-matrix resolution automatic flythrough enhances navigational image quality, and provides operational flexibility, functionality, and simplicity. Anatomy remains sharp with minimal pixilation. Steering remains centered as you navigate through the colon. You can specify your navigation speed, and 3D navigation is correlated with 2D and 360-degree Dissection views for rapid synchronized referencing in primary or secondary reading.

512 matrix resolution

Features detail (cont.)

Save state

This feature gives you the ability to save any centerline edits, bookmarks, and colored polyps for later reloading. It is also convenient for quick recall during consultation or exam review.

Prone/supine polyp linking


After bookmarking a polyp in prone and supine datasets, the program provides you with the capability to match and label them for easy tracking. This gives you the ability to read registered prone/supine datasets for a one pass read.

Synchronized prone/Supine 3D dissection view navigation

Once prone/supine polyp linkage is made, the program co-registers and links the datasets and polyps, letting you visualize opposing anatomy for a synchronized review.

Small bowel extraction

This feature quickly segments the small bowel, giving you an unobstructed view of the large intestine.


Analysis mode

The program's Analysis Mode provides a variety of tools to help streamline your workflow:

- Lock Mode lets you target a possible lesion and stay locked on it at all times.
- Synchronized display can be customized to show two oblique images, axial, sagittal, coronal, volume rendered and/or 3D views within the colon in a Lock mode. Leveraging multiple synchronized planes and views in a problem-solving mode may enhance your reading productivity.

Color centerline tracking

With Colon VCAR EC you can display a real-time 3D model of the colon with a defined centerline to quickly orient the location of the area you are viewing to the overall colon. This synchronization is always linked and available for quick reference.

Patient report

DICOM SR allows a reporting workflow to present and archive images and measurements.

Report feature only available on AW Workstation.

Measurement tools


With Colon VCAR EC, you can measure polyp distance from distal rectum to polyp location in millimeters, perform 2D measurements, and obtain statistics of ROI in Hounsfield units.

Archive

You can define the number of images per polyp you want to archive to either a networked location or to magnetic optical disk.

Customizable review and analysis layout

Create layouts that match your reading style and save them as protocols for future recall. Viewports can be shaped into a variety of square or rectangular configurations.


Features detail (cont.)

Customizable review and analysis layout (cont.)

Five of the most popular protocol layouts are pre-built for quick loading. Choose from several dissection configurations:

- 360-dissection gives you convenience and productivity benefits.
- 120-dissection lets you read around distortions if necessary.

- Use the 2D read configuration if the exam has poor insufflation and multiple collapses. 2D can also be used for extracolonic review.

- Use the 3-series layout to load a prone, supine, and decubitus series in parallel to obtain a correlated view.

- You can also customize new layouts to suit your own reading style.
- If you work on a one-monitor Advantage Workstation, a single monitor display is also available.

Movie generator

You can generate movie clips of the area of interest simply by defining the start and end locations with a mouse click. A cross-sectional movie will be generated. You can also generate lock-to-target movies. Movies can be saved as mpeg, avi, HTTP and saved to CD.

Entire colon movie

Generate a fly-through movie of the entire colon with one button click. Export the result as an mpeg or as screen captures to a PACS workstation for review or consultation.

Additional features

- Automatic center line tracking
- SmartCursor™ for easier navigation
- Synchronized reformatted views
- Multiple correlations of anatomy with CT colonography
- Exploration of suspicious areas in primary 2D reformatted sections
- Primary 3D reformatted sections. Obliques, Navigator, Quick Virtual Dissection, and axial views are all correlated to the same anatomy of interest.
- Black/white or color display

Summary

Colon VCAR EC helps make colon evaluation a less traumatic procedure for patients. Preparation for the exam is less severe; the procedure is less time intensive and less expensive than optical colonoscopy. Moreover, the tools and features this program puts into your hands will streamline your workflow and potentially enable greater diagnostic confidence.

About GE HealthCare Technologies Inc.

GE HealthCare is a leading global medical technology, pharmaceutical diagnostics, and digital solutions innovator, dedicated to providing integrated solutions, services, and data analytics to make hospitals more efficient, clinicians more effective, therapies more precise, and patients healthier and happier. Serving patients and providers for more than 125 years, GE HealthCare is advancing personalized, connected, and compassionate care, while simplifying the patient's journey across the care pathway. Together our Imaging, Ultrasound, Patient Care Solutions, and Pharmaceutical Diagnostics businesses help improve patient care from diagnosis, to therapy, to monitoring. We are a \$19.6 billion business with approximately 51,000 colleagues working to create a world where healthcare has no limits.

Follow us on [LinkedIn](#), [X](#) (formerly Twitter), [Facebook](#), [Instagram](#), and [Insights](#) for the latest news, or visit our website <https://www.gehealthcare.com/> for more information.

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
GE HealthCare

GSI Viewer

Aids in confident diagnosis with easier workflow

CONFIDENTIAL

SECRET DE
AFACERI

Contents

1. Features..... 4

2. System Requirements..... 4

3. Regulatory Compliance..... 4

4. Indications for Use..... 4

5. Feature Detail..... 4

5.01 Protocol Selection and Creation..... 4

5.02 Monochromatic Review..... 4

5.03 Image Overlay..... 4

5.04 Material Density Analysis..... 5

5.05 VOI Tool in Volume Viewer..... 5

5.06 Plot Analysis..... 5

5.07 Liver Fat in Volume Viewer..... 5

5.08 Material Suppressed Iodine/Virtual Unenhanced (VUE)..... 6

5.09 Gout Protocol..... 6

Gemstone Spectral Imaging (GSI) has the exciting potential to be used as a routine scan acquisition mode offering additional anatomical and functional information to help expedite and assist in an accurate CT diagnosis. GSI Viewer leverages the capability of spectral imaging with the acknowledged 3D segmentation and visualization capabilities of Volume Viewer.

Overview


Gemstone Spectral Imaging (GSI Viewer) is a dual energy application that uses rapid kV switching to acquire the dual energy samples almost simultaneously to generate material density data that can be used for the separation of materials and derivation of monochromatic spectral images using a projection based reconstruction algorithm.

With the integration of GSI into Volume Viewer (GSI Viewer), GSI now has the potential to become a routine scan acquisition mode, offering additional anatomical and functional information to aid diagnosis with CT data. Used routinely, the material decomposition images can be leveraged to separate materials for closer evaluation, including, for example, characterizing the composition of lesions, imaging plaque etc. The exclusive GE HealthCare spectral (monochromatic) images will optimize low contrast visualization for ideal image display, reduce beam-hardening artifacts to unmask unseen anatomy and promote accurate CT numbers for quantitative assessment.

The GSI Viewer is the primary method of viewing and manipulating spectral images. It allows for the primary review of monochromatic energy images at user selectable energy levels, detailed analysis using material decomposed images (such as water-iodine, water calcium, etc.), and complementary information using the Effective-Z (provides an estimate of the protons, Z number, in a voxel) images.

Highlights

- Gout protocol with automated workflow
- Liver fat images
- Virtual unenhanced (VUE) images
- View material density and effective atomic number images to discriminate tissue type
- Use image overlay tool to fuse material and effective atomic number information on monochromatic image
- Material loading and management in Volume Viewer

CONFIDENTIAL

SECRETE DE
AFACERI

Visit us:

<https://www.gehealthcare.com/products/advanced-visualization/all-applications/gsi-viewer>

Features

- Material decomposed images allow for the separation of materials like calcium, iodine, hydroxyapatite and water.
- Visualize a virtual non-contrast like image using the water-iodine basis pair image.
- Adjusting monochromatic energy levels can optimize image contrast and reduce beam-hardening artifacts.
- Discriminate different tissue types based on material density and monochromatic image data.
- Image overlays to visualize different material attributes in a single view.
- Compatible with multiple Volume Viewer applications
- Iodine suppressed images with the dynamic range of regular CT images in addition to a virtual non-contrast like image using the water-iodine basis pair.
- Liver fat images
- Automated gout workflow
- Ability to export GSI generated volumes

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Rx Only

Indications for Use

The GSI Viewer accepts images from a CT System that can acquire CT images using different kV levels of the same anatomical region of a patient in a single rotation from a single source.

The differences in the energy dependence of the attenuation coefficient of the different materials provide information about the chemical composition of body materials. This approach enables images to be generated at energies selected from the available spectrum to visualize and analyze information about anatomical and pathological structures.

GSI provides information of the chemical composition of renal calculi by calculation and graphical display of the spectrum of effective atomic number. GSI Kidney stone characterization provides additional information to aid in the characterization of uric acid versus non-uric acid stones. It is intended to be used on non-contrast studies as an adjunct to current standard methods for evaluating stone etiology and composition.

Feature Detail

Protocol Selection and Creation

GSI Viewer comes with a set of standard protocols that users can select depending on the type of exam.

A key feature is the ability for users to create and save their own protocol that can define the image layouts as well as what type of images are selected for view in specific viewports.

This offers improved reading workflow since the monochromatic images used for routine anatomy assessment can be linked to the Material Decomposed images for quick problem solving.

Monochromatic Review

One of the basic uses of the GSI Viewer is the ability to review the monochromatic images in multi-planar, curved and lumen reformats and also in 3D VR mode. It also provides the ability to interactively change the monochromatic energy levels so that users can select the best energy level for the exam being reviewed. Once an energy level is selected it can be easily saved for subsequent viewing or for sending to PACS.



Image Overlay

GSI Viewer provides a simple way to move from review to more detailed analysis by overlaying material density or Effective Z information on top of the monochromatic images. Users have the ability to fuse various images with a simple drag and drop and adjust the W/L as well as the opacity of the overlaid images. An example of Effective Z overlay on monochromatic image is shown in the image.


CONFIDENTIAL

SECRET DE
AFACERI

Features Detail (cont.)


Material Density Analysis

Users can visually see how the GSI data is segregated amongst a material density pair, e.g. water and iodine. They can switch between material pairs from a list. This ability together with the image overlay provides powerful analysis tools for the analysis of the components of anatomical regions. Iodine and water pair can help in assessing the material composition of liver lesions that can aid in their diagnosis.

VOI Tool in Volume Viewer



The VOI tool allows the user to save user defined VOI presets, which can be used at a later time to selectively visualize a volume of interest. GSI Viewer also allows the user to color the selected VOI, view it in 3D VR format, and overlay it on the 3D anatomic image or 2D image at an optimum transparency. The illustration displays a VOI in gout overlaid on the 3D anatomic image.

Plot Analysis

The GSI viewer offers a Plot Analysis feature that displays ROI's as graphical plots. These plots may be a histogram, spectral HU curve and normalized spectral HU curves. Plot analysis enables graphical visualization of how a selected region behaves at different energy settings. These tools have the potential to help in assessing the material in the ROI's for the purpose of lesion characterization and kidney stone analysis.

On Spectral HU Curves, users can normalize ROIs with the user-defined reference ROI. Reference ROI can be selected based on any particular type of image (Iodine or water material density image) or based on graphical plots (Spectral HU Curves). This normalization function may help reduce variability between patients or over time.

Liver Fat in Volume Viewer


Users can select the liver fat tab from the GSI options in Volume Viewer to generate a Liver Fat image displayed in approximate percentage of liver fat.

CONFIDENTIAL SECRET DE AFACERI

Features Detail (cont.)


Material Suppressed Iodine/Virtual Unenhanced (VUE)

The images below illustrate examples of a GSI monochromatic image, an MD Water image and an MSI image. Users can select the Material Suppressed Iodine tab from the GSI options in Volume Viewer to generate the MSI/VUE image.

Gout Protocol

The gout protocol with automated workflow can be configured by the user to generate and export 2D and 3D gout overlay images to the patient browser. This protocol can be used manually or enabled in preprocessing.

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
GE HealthCare

READY View

Quantified MR image analysis to help enable accurate, confident diagnoses.

Specialized Magnetic Resonance Imaging (MRI) methods such as diffusion weighted MR (DW-MR), dynamic (time-course) MR, spectroscopy, and functional methods have evolved to the point where they are able to provide unique measurements of tissue properties.

A multi-parametric MR approach generates zones within a lesion that reflect heterogeneity and often display characteristic patterns. These attributes have proven to be useful in the diagnosis of tumors, monitoring tumor growth, and guiding biopsies. Consequently, multi-parametric studies are often used in stroke diagnosis and in the diagnosis of cancer related lesions.

While extremely useful, quantifying the information from multi-parametric studies can be very challenging and time consuming.

READY View streamlines multi-parametric analysis.

Overview

READY View helps you get the most from multi-parametric exams by enabling analysis of MR data sets with multiple images for each scan location. The user experience driven framework offers a combination of protocols and tools that help you make quantified analyses of multiple data sets quickly and easily.

Highlights

- Guided workflows with intelligent display based on smart layout to help analyze MR data.
- Adapt your application to fit personal and institution requirements for more standardized analysis and improved productivity
- Provides additional clinical information through (relative) ROI measurements curves, spectra and color parametric images.
- Enables fusion of color parametric images with anatomical 2D or 3D images with simple "drag and drop" method.
- Provides adaptive protocols for multi-parametric data processing.
- Enables MR to MR image registration to reduce patient motion effects.
- Accessible from PC, laptop, PACS/RIS workstation for streamlined workflow.

CONFIDENTIAL

SECRET DE
AFACERI

Visit us:

[www.gehealthcare.com/aw/
platforms/aw-server/](http://www.gehealthcare.com/aw/platforms/aw-server/)

Features

- Analyze the following type of MR data sets:
 - Time series
 - Diffusion weighted scan
 - Diffusion tensor scan
 - Variable echo imaging
 - Blood oxygen level dependent imaging
 - Spectroscopy (single voxel and 2D or 3D CSI)
 - Elastography¹ imaging
- Simple workflow to process and fuse functional data.
- Select and process functional data with One Touch single click capability.
- READY View automatically selects the most relevant protocol for you.
- Efficient multi-contrast exam reading using MR General Review based on smart layout technology.
- Adaptive multi-parametric protocols as guided workflow to streamline processing and analysis of multi-parametric studies.
- Display all multi-parametric images and get all related functional values from a single ROI deposition.
- Fully customizable workflows with adjustable layouts, personalized parameter and settings, custom review steps.
- Easy-to-use slide bars let you segment parametric images in real time.
- Display and export ROI statistics

- from the Summary table
- Export graph values as csv file
- Save State let you save and restore the state of the processed images at any stage.
- Contextual help pages that give general assistance about the image processing algorithms.
- Save all generated parametric images in one click.

System Requirements

- AW Server 3.2 and above
- AW 4.7 Workstation and above
- Color Monitors
- Centricity™ Universal Viewer

Software Requirements

- Integrated Registration is a prerequisite
- Certain functions require Body View, Brain View and MR Touch options


Indications for Use

READY View, is an image analysis software that allows the user to process dynamic or functional volumetric data and to generate maps that display changes in image intensity over time, echo time, b-value (diffusion imaging) and frequency (spectroscopy). The combination of acquired images, reconstructed images, calculated parametric images, tissue

segmentation, annotations and measurement performed by the clinician enables a multi-parametric analysis and may provide clinically relevant information for diagnosis.

Regulatory Compliance

This product complies with the European Council Directive 93/42/EEC Medical Device Directive as amended by European Council Directive 2007/47/EC

Notes:

- ¹ Requires purchase of the appropriate application license


CONFIDENTIAL

SECRET DE
AFACERI

GE imagination at work

© 2016 General Electric Company.
All rights reserved. Data subject to change.
GE, the GE Monogram, imagination at work and Centricity are trademarks of
General Electric Company.
DOC1915501
GE Medical Systems SCS
283 rue de la Minière
78530 Buc France.

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
DOC1238679
GE Medical Systems
SCS 283 rue de la
Miniere 78530 Buc
Franța

Body View

An innovative way to analyze time series data from body studies.

Magnetic Resonance Imaging (MRI) of the body yields valuable information that is useful in determining a diagnosis. But processing time series data can complicate workflow and reduce productivity.

Body View streamlines the review of MR body images.

Overview

Body View provides algorithms, tools, and workflows for processing time series data acquired in the body. It also calculates parametric images from contrast enhanced images based on the temporal evolution of signal intensity.

Highlights

- Guided workflows with intelligent display based on smart layout to help analyze MR data.
- Adapt your application to fit personal and institution requirements for more standardized analysis and improved productivity
- Provides additional clinical information through relative ROI measurements, curves and color parametric images.
- Enables fusion of color parametric images with anatomical 2D or 3D images with simple "drag and drop" method.
- Provides adaptive protocols for multi-parametric data processing.
- Enables MR to MR image registration to reduce patient motion effects.
- Accessible from PC, laptop, PACS/RIS workstation for streamlined workflow.

CONFIDENTIAL

SECRET DE
AFACERI

Visit us:

www.gehealthcare.com

Features

Body View offers two advanced protocols:

- MR Standard for analyzing T1 or T2* contrast changes over time. Parametric images generated include:
 - Enhancement integral (EI)
 - Time to peak (TTP)
 - Mean time to enhance (MTE)
 - Maximum slope of increase (MSI)
 - Maximum slope of decrease (MSD)

MR Standard automatically determines enhancement type (positive, i.e. T1 contrast, or negative, i.e. T2* contrast)

- **Signal Enhancement Ratio (SER)** lets you analyze T1 contrast changes over time. Images provided include:
 - Signal Enhancement Ratio (SER)
 - Maximum slope of increase (MSI)
 - Positive enhancement integral (PEI)
- Body View offers adaptive protocol to process multi-parametric breast, pelvis and liver data. Factory protocols are adjustable and can be customized to fit your personal requirements in terms of

- Display & workflow (layouts,

review steps).

- Parameters and settings

- Easy-to-use slide bars let you segment parametric images in real time.
- Display and export ROI statistics from the Summary table
- Export graph values as csv file
- Save State let you save and restore the state of the processed images at any stage.
- Contextual help pages that give general assistance about the image processing algorithms
- Save all generated parametric images in one click.

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Software Requirements

- Body View requires READY View as prerequisite.

Product Description

Body View is an image analysis software option of READY View that allows the user to view and process MR images of the body.

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

MR VessellQ™ Xpress

Vessel analysis made flexible.

Using radiation-free Magnetic Resonance Angiography (MRA) as a tool for vascular disease diagnosis may provide sound information for making therapeutic decisions. Increased sophistication of both sequence design and parallel acquisition techniques have made contrast-enhanced Magnetic Resonance Angiographic (MRA) a powerful, practical diagnostic tool.

Furthermore, the recent introduction of contrast-free Magnetic Resonance Angiographic (MRA) makes it competitive with other diagnostic techniques in light of the current emphasis on dose management, cost containment, and outpatient evaluation.

MR VessellQ Xpress provides the tools to efficiently evaluate Magnetic Resonance Angiographic (MRA) information to help aid in vascular disease assessment.

Overview

Integrated with the Volume Viewer platform, MR VessellQ offers you advanced techniques to easily track the vasculature geometry and make measurements. The software includes display, measurement, and batch filming/archive features to let you efficiently analyze selected vessels for stenosis, directional tortuosity and other anomalies.

MR VessellQ™ Xpress is available on VolumeShare 7, a multi-modality advanced visualization workflow solution that helps to enhance diagnostic precision and productivity.

Highlights

- Guided workflows and smart layouts to help analyze vascular MR data.
- Adapt your application to fit personal and institution requirements for more standardized analysis and improved productivity.
- Intuitive in-view tools for editing vessel centerline.
- Quantify abnormal vascular structures and assess changes.
- Visualize and analyze directional vessel tortuosity.
- “Show Tracking” shortens waiting time for first clinically relevant image.
- Provides quick 3D visualization, and fast access to vessel cross section and profile images.
- Accessible from PC, laptop, PACS/RIS workstation for streamlined workflow.

Visit us:
[www.gehealthcare.com/aw/
platforms/aw-server/](http://www.gehealthcare.com/aw/platforms/aw-server/)

Features

MR VessellIQ Xpress offers ten advanced protocols:

- Anatomy-based protocols enable placement of intermediate points when needed.
- Pre-defined label database accurately labels vessels when using Advanced Vessel Analysis (AVA) based protocols.
- **Automatically provides size, stenosis, and length** measurements with two deposited points.
- **Single or double Quick AVA from any protocol within Volume Viewer to analyze vessels in curved reformat, lumen, or MPR views.**
- Allows saving of current processing state to include measurements, 3D segmentation, and tracking.
- Features anatomy-based vessel labeling.
- Quickly and easily switch from one protocol to another from the review step area.
- Customizable review step accommodates particular workflows.
- Allows retrieval of the saved state so you can resume analysis from where you left off.
- Includes five pre-defined layouts to best suit your workflow.

- Ability to customize layouts for dual screen anatomy review.
- Capture measurement results tables and associated images to include in reports.

System Requirements

- AW Server 3.1 and above and recommended monitor resolution is up to dual 2MP (1600 x 1200) or a single 3MP (1536 x 2048).
- AW VolumeShare7 Workstation and above
- Centricity™ Universal Viewer¹

Notes:

¹ AW Server 3.1 is not compatible with Centricity Universal Viewer

Indications for Use

MR VessellIQ Xpress is intended to provide an optimized non-invasive application to facilitate vascular anatomy and pathology analysis from a set of DICOM 3.0 compliant 3D contrast-enhanced Magnetic Resonance Angiographic (MRA) images.

MR VessellIQ Xpress is a post processing application which can be used in the analysis of MRA data for the purpose of vascular disease assessment.

This software is designed to assist radiologists and other clinicians in the evaluation and assessment of vascular anatomy and disease with the capability to provide a set of tools for visualizing directional tortuosity, for sizing the vessel and for measuring areas of anomalies within a vessel.

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SECRET DE
AFACERI

GE imagination at work

© 2015 General Electric Company.
All rights reserved. Data subject to change.
GE, the GE Monogram, imagination at work and Centricity are trademarks of
General Electric Company.
DOC1665437
GE Medical Systems SCS
283 rue de la Minière
78530 Buc France

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