CMV IgG

**For the Qualitative Detection of
IgG Antibodies to Cytomegalovirus (CMV)
in Human Serum**

For use on IMMULITE® 2000 systems

SIEMENS

IMMULITE[®] 2000 CMV IgG

English

Intended Use: For *in vitro* diagnostic use with the IMMULITE[®] 2000 Systems Analyzers — for the qualitative detection of IgG antibodies to cytomegalovirus (CMV) in human serum, as an aid in the determination of serological status to CMV.

This kit is not FDA cleared for use in testing (i.e. screening) blood or plasma donors.

Catalog Number: **L2KCVG2** (200 tests)

Test Code: **CVG** Color: **Brown**

Summary and Explanation

Cytomegalovirus (CMV), a member of the herpesvirus family, is found throughout the world. Humans of all ages are susceptible and infection is spread through sexual contact, direct exposure to infected body fluids, blood transfusions and organ transplants.^{1,2,5,6} The majority of infections are asymptomatic; however, CMV infections can be severe in neonates and immunocompromised individuals.^{6,7} Infection can also be severe in patients with congenital or acquired cellular immune defects, including cancer patients, organ recipients and AIDS patients.^{5,6,8}

CMV is the most common congenital infection, infecting between 0.5 and 2.5 percent of newborn infants. Five percent of these will develop classic cytomegalic inclusion disease with jaundice, pneumonia and central nervous system disorder. Infected infants may be asymptomatic at birth, but can develop neurological problems later in life.^{3,4}

Between 40 and 100 percent of people have detectable antibody,⁹ with the prevalence highest in developing countries.

Principle of the Procedure

IMMULITE 2000 CMV IgG is a solid-phase, sequential chemiluminescent enzyme immunoassay.

Incubation Cycles: 2 × 30 minutes

Specimen Collection

The use of an ultracentrifuge is recommended to clear lipemic samples.

Hemolyzed samples may indicate mistreatment of a specimen before receipt by the laboratory; hence the results should be interpreted with caution.

Samples which are cloudy or have particulate material should be clarified by low-speed centrifugation.

Centrifuging serum samples before a complete clot forms may result in the presence of fibrin. To prevent erroneous results due to the presence of fibrin, ensure that complete clot formation has taken place prior to centrifugation of samples. Some samples, particularly those from patients receiving anticoagulant therapy, may require increased clotting time.

Blood collection tubes from different manufacturers may yield differing values, depending on materials and additives, including gel or physical barriers, clot activators and/or anticoagulants. IMMULITE 2000 CMV IgG has not been tested with all possible variations of tube types. Consult the section on Alternate Sample Types for details on tubes that have been tested.

Volume Required: 10 µL serum

Storage: 3 days at 2–8°C, or 6 months at –20°C.¹¹

Automatic Predilution Factor: 20

Warnings and Precautions

For *in vitro* diagnostic use.

CAUTION! POTENTIAL BIOHAZARD

Contains human source material. Each donation of human blood or blood component was tested by FDA-approved methods for the presence of antibodies to human immunodeficiency virus type 1 (HIV-1) and type 2 (HIV-2) as well as for hepatitis B surface antigen (HBsAg) and antibody to hepatitis C virus (HCV). The test results were negative (not repeatedly reactive). No test offers complete assurance that these or other infectious agents are absent; this material should be handled using good laboratory practices and universal precautions.¹²⁻¹⁴

CAUTION: This device contains material of animal origin and should be handled as a potential carrier and transmitter of disease.

Reagents: Store at 2–8°C. Dispose of in accordance with applicable laws.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing.

The bead is coated with *inactivated* CMV antigen. However, caution is advised because of the possible presence of residual virus when working with, or disposing of, the materials supplied.

The anti-CMV IgG results determined for a given specimen with assays from different manufacturers can vary due to differences in assay methods and reagent specificity. Therefore, the results reported by the laboratory to the physician should include: "The following results were obtained with the IMMULITE 2000 CMV IgG EIA. Results obtained from other manufacturers' assay methods may not be used interchangeably."

Chemiluminescent Substrate: Avoid contamination and exposure to direct sunlight. (See insert.)

Water: Use distilled or deionized water.

Materials Supplied

Components are a matched set. Labels on the inside box are needed for the assay.

CMV IgG Bead Pack (L2CVG12)

With barcode. 200 beads, coated with inactivated, partially purified CMV antigen (strain AD-169 from infected cell lysates). Stable at 2–8°C until expiration date.

L2KCVG2: 1 pack

CMV IgG Reagent Wedge (L2CVGA2)

With barcode. Two reagents: 11.5 mL buffer solution, 11.5 mL alkaline phosphatase (bovine calf intestine) conjugated to monoclonal murine anti-human IgG antibody in buffer, with preservative. Stable at 2–8°C until expiration date.

L2KCVG2: 1 wedge

Before use, tear off the top of the label at the perforations, without damaging the barcode. Remove the foil seal from the top of wedge; snap the sliding cover down into the ramps on the reagent lid.

CMV IgG Adjustor (LCVGR)

4 mL human serum with IgG reactive to CMV, in buffer, with preservative. The Adjustor serves as the assay's Cutoff. Stable at 2–8°C for 14 days after opening, or for 6 months (aliquotted) at –20°C.

L2KCVG2: 1 vial

Before making an adjustment, place the appropriate Aliquot Labels (supplied with the kit) on test tubes so that the barcodes can be read by the on-board reader.

CMV IgG Controls (LCVGC1, LCVGC2)

LCVGC1: One vial (2.0 mL) negative control containing human serum with IgG nonreactive to CMV, with preservative.

LCVGC2: One vial (2.0 mL) positive control containing human serum with IgG reactive to CMV, with preservative. Stable at 2–8°C for 14 days after opening, or for 6 months (aliquotted) at –20°C.

L2KCVG2: 1 set

The IMMULITE 2000 software performs automatic on-board dilution of control samples, and the results will be tracked in the QC database. Enter controls as controls.

For the current control ratio ranges, please refer to the Control insert.

IgG/IgM Sample Diluent (L2IGZ2)

For the on-board dilution of patient samples and controls. 55 mL Concentrated (ready-to-use) nonhuman protein/buffer matrix, with preservative. Stable at 2–8°C for 30 days after opening, or 6 months (aliquotted) at –20°C.

L2KCVG2: 1 vial

Barcode labels are provided for use with the diluent. Before use, place an appropriate label on a 16 × 100 mm test tube, so that the barcodes can be read by the on-board reader.

L2KCVG2: 3 labels

Kit Components Supplied Separately

IgG/IgM Sample Diluent (L2IGZ2)

For the on-board dilution of patient samples and controls. 55 mL Concentrated (ready-to-use) nonhuman protein/buffer matrix, with preservative. Stable at 2–8°C for 30 days after opening, or 6 months (aliquotted) at –20°C.

L2IGZ2: 1 vial

Barcode labels are provided for use with the diluent. Before use, place an appropriate label on a 16 × 100 mm test tube, so that the barcodes can be read by the on-board reader.

L2IGZ2: 5 labels

L2SUBM: Chemiluminescent Substrate

L2PWSM: Probe Wash

L2KPM: Probe Cleaning Kit

LRXT: Reaction Tubes (disposable)

L2ZT: 250 Sample Diluent Test Tubes (16 × 100 mm)

L2ZC: 250 Sample Diluent Tube Caps

LCVGC: Bi-level CMV IgG Control Module

Also Required

Distilled or deionized water; test tubes

Assay Procedure

Note that for optimal performance, it is important to perform all routine maintenance procedures as defined in the IMMULITE 2000 Systems Operator's Manual.

See the IMMULITE 2000 Systems Operator's Manual for preparation, setup, dilutions, adjustment, assay and quality control procedures.

Recommended Adjustment Interval:
4 weeks

Quality Control Samples: The CMV IgG Controls supplied with the kit should be used as quality control material to monitor the performance of the assay at the cutoff range. The Positive Control is used to validate the IMMULITE 2000 CMV IgG assay at a critical level when determining the serological status of a patient.

The Controls should be processed at (or near) the beginning of every run containing patient samples to be tested for CMV IgG, and also when readjusting. (For many laboratories, this will mean processing appropriate controls once per shift.)

For the current control ratio ranges, please refer to the Control insert.

Calculation of Cutoff and S/CO Ratio:

The Master Cutoff of the assay was determined from representative samples to achieve optimal sensitivity and specificity for the assay.

The cutoff is set equal to the average counts per second (mean cps) of the Adjustor (from the most recent adjustment) multiplied by Curve Parameter 1. (See the "Low Adjustor CPS" and "Curve Parameter 1" fields in the IMMULITE 2000 Kit Information screen, which can be accessed from the menu via Data Entry: Kit Entry.)

Calculation of a signal/cutoff (s/co) ratio is done by using the following formula:

$$\text{S/CO Ratio} = \frac{\text{Sample or Control cps}}{\text{Mean Adjustor cps} \times \text{P1}}$$

Calculation and reporting of qualitative (reactive/nonreactive/indeterminate) and s/co ratio results are handled automatically by the IMMULITE 2000.

The result for a sample is reported as "Indeterminate" if the counts per second for that sample fall within $\pm 10\%$ of the cutoff. The result is reported as "Reactive" if the sample's counts are *above* the indeterminate range, and "Nonreactive" if *below* this range.

Additional controls may be tested in accordance with guidelines or requirements of local, state and/or federal regulations or accrediting organizations.

It is also recommended that known reactive and nonreactive specimens be run periodically to assure pipetting accuracy for the dilution step.

Interpretation of Results

A result of **"Reactive"** (ratio of ≥ 1.1) indicates that CMV IgG antibodies were detected in the patient sample.

A result of **"Nonreactive"** (ratio of < 0.9) indicates that CMV IgG antibodies were not detected in the patient sample.

Any result of **"Indeterminate"** (ratio between 0.9 and < 1.1) should be retested. Samples which still test as "Indeterminate" should be tested by an alternate method, or a second sample should be taken — if possible — within a reasonable period of time (e.g., one week).

The presence of IgG antibodies to CMV is an indication of previous exposure to the virus. A single specimen can only be used to determine the serological status of the individual.

The magnitude of the measured results (cps) above the Cutoff is not indicative of the total amount of antibodies detected.


Expected Values

Based on its relationship to IMMULITE CMV IgG (see Method Comparison), the assay can be expected to have essentially the same reference ranges.

Individuals with CMV may not exhibit detectable levels of IgG antibody in the early stages of infection. Levels of IgG antibody to CMV begin to rise one or two weeks after primary infection. Normal, healthy individuals will not have significant changes in their antibody level with reactivation of the virus. Patients with AIDS and other immunosuppressed

individuals — e.g. those with organ transplants — may have rapid increases in IgG levels with recurrent infections.^{3,6,9}

Studies with presumed healthy, asymptomatic subjects were conducted at two locations in the US. At one location (Site 1) in the northwestern US, the subjects consisted of 17 individuals undergoing a pre-employment screening (4 females and 13 males), with ages ranging from 21 to 42 years (samples were stored at 2–8°C if the assay was performed within 4 days, otherwise samples were frozen at –20°C). At a second location in the southwestern US, the subjects consisted of 100 blood donors (54 females and 46 males), with ages ranging from 17 to 71 years (samples were collected and frozen).

The prevalence of CMV reactivity was found to be 59% (Site 1) and 66% (Site 2) in these two distinct geographical locations in the United States.

These study sites also included 48 (Site 1) and 71 (Site 2) pregnant women, with ages ranging from 18 to 43 years (Site 1) and 14 to 48 years (Site 2). The prevalence of CMV reactivity for this specific population at these two sites was 46% and 76%, respectively.

Additional studies were conducted at Site 1 with a total of 55 patients who were either HIV positive, had AIDS, or were immunocompromised for other reasons. The prevalence of CMV infections in this group was 75%.

The prevalences listed above for asymptomatic, pregnant and immunocompromised subjects are similar to those reported in the literature.⁹ However, prevalence may vary due to geographical or population differences.

Consider these limits as *guidelines* only. Each laboratory should establish its own reference ranges.

Limitations

The results of the test must be taken within the context of the patient's clinical history, symptomology and other laboratory findings.

The presence of IgG antibodies in a single specimen is not sufficient to distinguish between active or past infection.

The presence of a primary infection, or reactivation of a past infection, should be determined by isolation of the virus in tissue culture.^{3,9}

Patients suspected of having primary or active infection should be tested for the presence of IgM antibodies to CMV.

For the determination of seroconversion from *nonreactive* to *reactive*, a second serum sample should be drawn three to four weeks after the acute, nonreactive sample, during the convalescent stage of the infection. The acute phase sample should be stored and tested in parallel with the convalescent sample.

Individuals with acute CMV infection may not exhibit any detectable IgG antibodies at the early stage of infection.

Samples containing antinuclear antibodies or other anticell antibodies may give false reactive results.

A rise in CMV antibody level may occur in patients with measles, varicella-zoster virus (VZV) or herpes simplex virus (HSV) due to antigenic crossreactivity within the herpesvirus family. Patients with acute EBV (Epstein-Barr Virus) infection may show increased reactivity to other viruses, including CMV.

The results in HIV patients, in patients undergoing immunosuppressive therapy, or in patients with other disorders leading to immunosuppression, should be interpreted with caution.

The performance characteristics of this assay have not been established for use with specimens from neonates, cord blood, or pretransplant patients.

Heterophilic antibodies in human serum can react with the immunoglobulins included in the assay components causing interference with *in vitro* immunoassays. [See Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988;34: 27-33.] Samples from patients routinely exposed to animals or animal serum products can demonstrate this type of interference potentially causing an anomalous result. These reagents have been formulated to minimize the risk of interference; however, potential interactions between rare sera and test components can occur. For diagnostic purposes, the results obtained from this assay should always be used in combination with the clinical examination, patient medical history, and other findings.

Performance Data

See Tables and Graphs for data *representative* of the assay's performance. Results are expressed as a signal-to-cutoff ratio. (Unless otherwise noted, all were generated on serum samples collected in tubes without gel barriers or clot-promoting additives.)

Intraassay Precision (Within-Run):

Statistics were calculated for samples from the results of 20 replicates in a single run. (See "Intraassay Precision" table.)

Interassay Precision (Run-to-Run):

Statistics were calculated for samples assayed in 4 different runs. (See "Interassay Precision" table.)

Alternate Sample Type: To assess the effect of alternate sample types, blood was collected from 20 volunteers into plain, heparinized, EDTA and Becton Dickinson SST[®] vacutainer tubes. Equal volumes of the matched samples were spiked with various concentrations of CMV IgG, to obtain values throughout the calibration range of the assay, and then assayed by the IMMULITE 2000 CMV IgG procedure.

(EDTA) = 0.99 (Serum) + 0.03
r = 0.998

(Heparin) = 0.99 (Serum) + 0.19
r = 0.997

(SST) = 1.0 (Plain Tubes) + 0.10
 r = 0.997

Means:

- 9.7 (Serum)
- 9.6 (EDTA)
- 9.8 (Heparin)
- 9.9 (SST)

Bilirubin: Presence of conjugated and unconjugated bilirubin in concentrations up to 200 mg/L has no effect on results, within the precision of the assay.

Hemolysis: Presence of hemoglobin in concentrations up to 539 mg/dL has no effect on results, within the precision of the assay.

Lipemia: Presence of triglycerides in concentrations up to 3000 mg/dL has no effect on results, within the precision of the assay.

Method Comparison: The assay was compared to IMMULITE CMV IgG on 229 samples.

IMMULITE	IMMULITE 2000		
	React.	Ind	Nonreact.
Reactive	165	0	0
Indeterminate	0	0	1
Nonreactive	0	0	63

Agreement: 99.6%
 Relative Sensitivity: 100%
 Relative Specificity: 100%

References

1) Drew WL. Herpesviridae: cytomegalovirus. In: Lennette EH, Halonen P, Murphy FA, editors. Laboratory diagnosis of infectious diseases: principles and practice, volume II, viral, rickettsial and chlamydial diseases. New York: Springer-Verlag, 1988: 247-60. 2) Ho M. Characteristics of cytomegalovirus. In: Greenough WB, Merigan TC, editors. Cytomegalovirus biology and infection: current topics in infectious disease. New York: Plenum, 1982: 9-32. 3) Hodinka RL, Friedman HM. Human cytomegalovirus. In: Balows A, Hausler, Jr WJ, Hermann KL, Isenberg HD, Shadomy HJ. Manual of clinical microbiology, 5th ed. Washington, DC: American Society for Microbiology, 1991: 829-37. 4) Reynolds DW, et al. Laboratory diagnosis of cytomegalovirus infections. In: Lennette EH, Schmidt NJ, editors. Diagnostic procedures for viral, rickettsial and chlamydial infections, 5th ed. Washington: American Public Health Association, 1979: 399-439. 5) Van der Bij W, et al. Antigen test for early diagnosis of active cytomegalovirus infection in heart transplant recipients. J Heart Transplant 1988;7:106-9. 6) Zaia JA.

Epidemiology and pathogenesis of cytomegalovirus disease. Seminars in Hematology 1990;27:5-10. 7) Yow MD. Congenital cytomegalovirus disease: a new problem. J Infect Dis 1989;159:163-7. 8) Miller H, et al. Prospective study of cytomegalovirus antigenemia in allograft recipients. J Clin Microbiol 1991;29:1054-5. 9) Wiedbrauk DL, Johnston SLG. Manual of clinical virology. New York: Raven Press, 1993: 82-91. 10) Baltz ML, Searcy RL. Clinical significance and advanced serologic diagnosis of ToRCH infections. Am Clin Lab 1994;March/April:18-23. 11) Tietz NW, editor. Clinical guide to laboratory tests. 3rd ed. Philadelphia: WB Saunders, 1995:358. 12) Centers for Disease Control. Update: Universal precautions for prevention of transmission of human immunodeficiency virus, hepatitis B virus and other bloodborne pathogens in healthcare settings. MMWR, 1988;37:377-82, 387-8. 13) Clinical and Laboratory Standards Institute (formerly NCCLS). Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline - Third Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2005. NCCLS Document M29-A3. 14) Federal Occupational Safety and Health Administration, Bloodborne Pathogens Standard, 29 CFR 1910.1030.

Technical Assistance

In the United States, contact Siemens Healthcare Diagnostics Technical Services department. Tel: 877.229.3711. Outside the United States, contact your National Distributor. siemens.com/healthcare

The Quality System of Siemens Healthcare Diagnostics Products Ltd. is certified to ISO 13485.

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IMMULITE is a trademark of Siemens Healthcare Diagnostics.
© 2009 Siemens Healthcare Diagnostics. All rights reserved.