

# Pirop – Ophthalmic Scanner (A+B+CCT)



# Pirop

**PIROP – state-of-the-art ophthalmic ultrasonic device for Biometry ,**

**Visualization and Pachymetry of eyeballs (A-Scan + B-Scan + CCT- Pachymeter) – modified version.**

## **A-SCAN – Ocular Lens biometrics**

Digital 'A-scan' tool for ophthalmology, biometry and lens power calculation of intra-ocular implants.



## **B-SCAN – Eye examination**

Modern 'B-scan' tool for ophthalmology, interior eyeball imaging, retina, optical nerve ...

## **CCT- PACHYETER (P-SCAN) – Ocular cornea biometrics**

Up-to-date ultrasonic pachymeter used to measure the thickness of the eye's cornea (glaucoma screening and refractive surgery), with a very high sampling frequency – **400MHz** (a significant increase in measurement accuracy –  $\leq 2 \mu\text{m}$ ). Corneal pachymetry is an important test in the early detection of glaucoma (central corneal thickness CCT), and it help also surgeons by providing graphical surgical plans to eliminate corneal astigmatism.

PIROP user's interface based on touch screen technology makes operation easy and user-friendly. The touch screen can be also



## Specifications

### A-SCAN

- > Complete and quick biometry :
  - eye axial length (AXL)
  - anterior chamber depth (AC)
  - lens thickness (LENS)
  - vitreous length (VITR)
- > Measurement of all eye types: normal, cataract, dense cataract, aphakic, silicone oil vitreous, pseudoaphakic (PMMA, ACRYLATE, SILICONE)
- > OL six formulas: SRK II, SRK T, Holladay, Hoffer-Q, Binkhorst II, Haggis
- > Post refractive formulas: Double-K SRKT, Latkany / Flat-K SRKT/, Latkany /avg-K SRKT/, Masket
- > IOL power calculation comparison (for different lenses and calculation methods)
- > Contact and immersion methods (the immersion tube attached)

- > Automatic TGC adjustment, GAIN adjustment up to 100dB
- > Memory of 20 scans with measurements for further analysis
- > Calculations of means and standard deviations
- > Automatic control of measurement scatters (with manual correction)
- > Customizable for 10 users, 10 profiles for each.
- > Adjusting and testing via attached calibrator
- > Probe frequency – 12 MHz
- > Range of scan 45 mm / 2500 points per line
- > Clinical resolution 0.1 mm
- > Electronic resolution –/+ 0.01 mm

## **POWER SUPPLY**

- > Power supply from external AC/DC power adapter:
  - input: 100-230V AC 50/60Hz / max 0.7A
  - output: +12V DC / 2 A
- > Power consumption (12 DC supply) – ca 12 W
- > Electrical safety standards
  - Medical device Class IIa complies with MDD 93/42 EEC
  - Scanner complies with requirements for Class II devices of EN/IEC 60601-1

## **B-SCAN**

- > Probe frequency: 12 .. 15 MHz
- > Scanning angle : 55°

- > Gain dynamics up to 100dB / 80dB
- > TGC – zone gain adjustment in selected eye
- > 256 Levels Gray Scale
- > Axial clinical resolution: 0.12 mm
- > Lateral clinical resolution : 0.3 mm
- > Display modes : B, B+A, B+B
- > A-mode for selected line CV (in Run and Freeze modes)
- > Dynamic Range Correction
- > Image cache memory M1..M4 independent for right and left eyes
- > Zoom function (x2) in Run and Freeze modes
- > Digital Gain Correction in Freeze mode
- > Cineloop function (about 6 second) independent for left/right eyes
- > Distances measurements with ultrasound velocity correction
- > Area and angle measurements
- > Pointers, comments entered by the user
- > Composite video output for videoprinter

## **CCT – PACHYMETER**

- > Quick thickness measurements for all corneal types
  - the measurement of central corneal thickness (CCT)
  - thickness measurement at arbitrarily selected points

- standard deviation of the measured thickness
- ability for rejecting inaccurate measurements
  
- › Very high level of accuracy and repeatability of all measurements:
  - sophisticated algorithms for the accuracy improvement
  - very high frequency for measuring purpose- 400 MHz
  - high averaging ratio
  - 512 automatic measurements cycles
  
- › Intra-ocular pressure (IOP) calculations with measuring and/or manual correction
  
- › Ten IOP correction formulas : Kohlhas/Shah, Argus et al. Whitarce et al, Doughty, Ehlers, Dresden, Stodmeister, custom-own
  
- › BIAS percentage deviation of the measured thickness from 60% to 130%
  
- › TEST function for automatic checking of operation correctness (without using external phantoms or patterns)
  
- › Nine defined maps of corneal thickness with the number of points: 1, 5, 9, 13, 21, 25, 1 MULTI, 5 MULTI, 9 MULTI (central, paracentral, peripheral, parietal, mixed)
  
- › 20 MHz probe operating frequency
  
- › Default sound velocity – 1640 m/s ; range from 1400 to 2000 m/s
  
- › Measurement range from 220  $\mu\text{m}$  to 1100  $\mu\text{m}$
  
- › Measuring accuracy  $\leq 2 \mu\text{m}$
  
- › Resolution 1  $\mu\text{m}$

#### > TOUCH SCREEN

- easy, ergonomic and user-friendly operation via displayed menu
- virtual keyboard for entering patient data, etc.
  
- > Functional knob (gain control, review results, cine etc.)
  
- > Color LCD 7" display
  
- > Internal patient's databas
  
- > Ten (10) user profiles
  
- > Support for barcode reader
  
- > Possibility of connecting an external alphanumeric keyboard
  
- > Recording of images in internal and USB-Pendrive memory.
  
- > Internal thermal small printer ( for A-scan and P-scan reports )
  
- > Composite-video output PAL (for videoprinter or video monitor )
  
- > Support of multiple languages
  
- > Easy firmware upgrade via USB port

### **PIROP-PC software (optional)**

- > a special software that enables management of the patient database and measurements (images) obtained from PIROP scanner.
  
- > it can be installed in any external PC computer, tablet etc. operated by Windows.



a direct USB connection.

- > possibility to print on any printer the special report of exams
- > additional option – possibility to use the standard DICOM 3.0 to communicate with PACS/RIS systems (Store, Print, Worklist).

## WEIGHT AND DIMENSIONS

- > Dimensions (with built-in printer): 290 x 205 x 85 mm
- > Weight (with built-in printer): 1.5 kg

## Probes



### OA12 A-scan

Frequency: 20 MHz

Application: biometry and lens power calculation for intraocular implants.



### **OB12 B-Scan**

Frequency: 12/15 MHz

Application: imaging of interior of an eyeball, retina,  
optical nerve ..



### **OP20 Pachymetr**

Frequency: 20 MHz

Application: measurements of thickness for all  
corneal types

**E - brochures**



e - Bochure A+B+CCT

## Scanner images





## Images gallery



Patient: John Doe  
User: Echson

None G: 78 OA-12MHz 20-06-2017 14:23:52

OD/R Contact Contact Auto: 1 Average AXL: 23.70mm K1=42.50 D K2=43.50 D Kavg=43.00 D Target Anisotrophia: 0.00 D

IDEAL	IDEAL	MC40BD	MC40BD
const A = 118.0	const SF = 1.78	const A = 118.7	const ACD = 5.37
SRK/T	Holladay-1	SRK/T	Hoffer-Q
Emme 21.21		Emme 21.23	
IOL (D)	REF (D)	IOL (D)	REF (D)
18.00 1.47	18.50 1.15	19.00 1.25	19.00 1.28
20.00 0.81	20.00 0.82	20.00 0.88	20.00 0.80
20.50 0.48	20.50 0.49	20.50 0.24	20.50 0.28
21.00 0.14	21.00 0.15	21.00 -0.11	21.00 -0.08
21.50 -0.20	21.50 -0.18	21.50 -0.46	21.50 -0.43
22.00 -0.55	22.00 -0.53	22.00 -0.81	22.00 -0.77
22.50 -0.90	22.50 -0.88	22.50 -1.17	22.50 -1.12
23.00 -1.25	23.00 -1.23	23.00 -1.54	23.00 -1.48

Save OD/OS Print Formulas comparison Next Exit



Patient: John Doe  
User: Echson

None G: 68 OA-12MHz 20-06-2017 13:40:00

OD/R Contact Contact Auto: 1 Average AXL: 23.70mm K1=42.50 D K2=43.50 D Kavg=43.00 D Target Anisotrophia: 0.00 D

IDEAL	MA60MA	MC40BD	MZ30BD
const A = 119.0	const A = 118.8	const A = 118.7	const A = 118.4
SRK/T	SRK/T	SRK/T	SRK/T
Emme 21.21		Emme 21.09	
IOL (D)	REF (D)	IOL (D)	REF (D)
18.00 1.48	19.00 1.40	19.00 1.28	18.50 1.37
20.00 0.82	20.00 0.74	20.00 0.88	19.50 0.88
21.50 0.48	21.50 0.40	21.50 0.24	20.00 0.45
21.00 0.15	21.00 0.06	21.00 -0.10	20.50 -0.00
21.50 -0.20	21.50 -0.28	21.50 -0.45	21.00 -0.36
22.00 -0.54	22.00 -0.83	22.00 -0.81	21.50 -0.72
22.50 -0.89	22.50 -0.88	22.50 -1.17	22.00 -1.08
23.00 -1.25	23.00 -1.34	23.00 -1.53	22.00 -1.44

Save OD/OS Print Formulas comparison Next Exit



Patient: John Doe  
User: Echson

Select G: 68 OA-12MHz 20-06-2017 13:39:14

OD/R Contact Contact Auto: 1 Cataract Average AXL: 23.70mm

Wavelength	AC	L	V	AXL
1	3.58	4.59	15.49	23.67
2	3.58	4.61	15.49	23.68
3	3.61	4.55	15.53	23.70
4	3.65	4.61	15.49	23.75
5	3.62	4.79	15.32	23.73
6	3.60	4.84	15.29	23.73
7	3.62	4.79	15.30	23.71
8	3.63	4.65	15.40	23.69
9	3.68	4.58	15.40	23.64
10	3.60	4.76	15.32	23.68
Avg	3.61	4.70	15.39	23.70
SD	0.03	0.10	0.08	0.03

Save Print Marker Select Exit

Scan No.	AC	L	V	AXL
1	4.08	4.55	15.58	24.22
2	4.14	4.53	15.60	24.28
3	4.08	4.58	15.58	24.25
4	4.08	4.58	15.60	24.25
5	4.12	4.55	15.55	24.22
6	4.11	4.51	15.64	24.28
7	4.08	4.55	15.60	24.20
8	4.08	4.58	15.57	24.25
9	4.08	4.58	15.60	24.23
10	4.08	4.58	15.58	24.25
Avg	4.09	4.55	15.59	24.24
SD	0.02	0.02	0.02	0.02

Scan No.	AC	L	V	AXL
1	4.08	4.55	15.58	24.22
2	4.14	4.53	15.60	24.28
3	4.08	4.58	15.58	24.25
4	4.08	4.58	15.60	24.25
5	4.12	4.55	15.55	24.22
6	4.11	4.51	15.64	24.28
7	4.08	4.55	15.60	24.20
8	4.08	4.58	15.57	24.25
9	4.08	4.58	15.60	24.23
10	4.08	4.58	15.58	24.25
Avg	4.09	4.55	15.59	24.24
SD	0.02	0.02	0.02	0.02

Patient: John Doe  
User: ECHO-SON

OS/L Contact Cataract Sclera  
Auto 1

Average AXL: 24.24mm

AC	L	V	AXL	
1	4.08	4.55	15.58	24.22
2	4.14	4.53	15.60	24.28
3	4.08	4.58	15.58	24.25
4	4.08	4.58	15.60	24.25
5	4.12	4.55	15.55	24.22
6	4.11	4.51	15.64	24.28
7	4.08	4.55	15.60	24.20
8	4.08	4.58	15.57	24.25
9	4.08	4.58	15.60	24.23
10	4.08	4.58	15.58	24.25
Avg	4.09	4.55	15.59	24.24
SD	0.02	0.02	0.02	0.02

Good waveform:  
① tall between retina and sclera  
② sharp rise  
③ small SD

Save Print Marker Select Exit

ECHO-SON

PRIMP B-scan

OD/R

ECHO-SON

PRIMP B-scan

Scan 1  
Scan 2

Patient: John Doe  
User: Echoson

OD/R

Dx= [v=1550m/s]  
D+= [v=1550m/s]

Text M1..M4  
Caliper Area  
Save Exit

Patient: Echo-Son  
User: Echo-Son

OD/R

CV DR  
Zoom Range  
B+B Freeze

Patient: John Doe  
User: Echoson

OD/R

Dx= 22.5 mm [v=1550m/s]  
D+= 4.3 mm [v=1641m/s]  
Angle=1.3°

Distance 2D Distance CV  
Comment Area  
Save Exit

Patient: John Smith  
User: ECHO-SON

OS/L Zoom x2

Lens Chamber Retina Sclera

CV OD/OS  
Zoom off Next...  
Menu Run

Patient: John Smith  
User: ECHO-SON

Scan 1  
Scan 2

OD/R Vitreous Lens Optic Nerve

Dx= 3.1 mm [v=1532m/s]  
D+= 0.9 mm [v=1532m/s]  
Angle=2.6°

Distance 2D  
Comment Area  
Save Exit

Patient: Echo-Son  
User: Echo-Son

OS/L

Distance 2D Distance CV  
Comment Area  
Save Exit

Patient: Echo-Son  
User: Echo-Son

OS/L

Distance 2D Distance CV  
Comment Area  
Save Exit



NIP: 7160005139  
REGON: 430302780  
KRS: 00001 36693

**Export:**

export@echoson.eu

**National matters:**

info@echoson.com.pl

**Domestic service:**

serwis@echoson.com.pl

**tel.** +48 81 886 36 13

**fax** +48 81 886 83 10

**Multi -organ ultrasound**

Albit - ultrasound scanner  
Specialized applications

**Ophthalmology and dentistry**

PIROP - Ophthalmology Biometer  
PIROP - G -SCAN

**Bladder scanner**

Pinit - bladder scanner

**Opening Hours: 8.00 - 15.00**