EverFlex™

Self-expanding Peripheral Stent with Entrust™ Delivery System
Selvekspanderende, perifert stent med Entrust™ fremføringssystem
Selbstexpandierender peripherer Stent mit dem Einführsystem Entrust™
Σύστημα αυτοδιαστελλόμενου περιφερικού stent με σύστημα τοποθέτησης Entrust™
Sistema de stent periférico autoexpansible con sistema de liberación Entrust™
Itsestään laajeneva perifeerinen stentti ja Entrust™-sisäänviejäjärjestelmä
Stent périphérique auto-expansible avec système de pose Entrust™
Entrust™ bevezetőrendszerrel szerelt öntáguló perifériás sztent
Stent periferico autoespandibile con sistema di rilascio Entrust™
Zelfontplooiende perifere stent met Entrust™-plaatsingssysteem
Selvekspanderende perifer stent med Entrust™-innføringssystem
Obwodowy stent samorozprężalny z systemem wprowadzającym Entrust™
Stent periférico autoexpansível com sistema de entrega Entrust™
Саморасширяющийся периферический стент с системой доставки Entrust™
Självexpanderande perifer stent med införingssystemet Entrust™
Kendinden Genişleyen Periferik Stent ile Entrust™ Taşıma Sistemi

CE
2797

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Instructions for use

Device description

The EverFlex™ self-expanding peripheral stent with Entrust™ delivery system is a self-expanding nitinol stent system intended for permanent implantation. The EverFlex stent is made of a nickel titanium alloy (nitinol) and comes pre-mounted on a 5 Fr, 0.035 in over-the-wire (OTW) delivery system. The stent is cut from a nitinol tube into an open lattice design. The stent has tantalum radiopaque markers at the proximal and distal ends. After deployment, the stent opens to its predetermined diameter and exerts a constant, gentle outward force to establish patency.

The Entrust delivery system, shown in **Figure 1**, has a triaxial shaft design comprised of an inner shaft (1), silver retractable outer sheath (2), and a gold isolation sheath (3) for ease of use during deployment. The ergonomic deployment handle (8) has a removable safety locking pin (9), deployment wheel (10), and luer hub (6).

The distal portion of the delivery system (enlarged section) contains two radiopaque markers; one marker distal (4) and one marker proximal (5) to the constrained stent. The distal marker is mounted on the retractable outer sheath and the proximal marker is mounted on the inner shaft. The EverFlex stent is constrained within the outer sheath. The flexible catheter tip (7) is at the distal end of the outer sheath.

Before the procedure, the operator flushes saline through the guidewire lumen until it exits at the end of the catheter. Radiopaque markers on the EverFlex stent give guidance for positioning the stent before deployment. After the stent is correctly positioned, the operator removes the safety locking pin to unlock the deployment wheel and enable the outer sheath to retract. The operator retracts the outer sheath by slowly rotating the deployment wheel in the direction indicated by the arrow on the deployment handle. When the radiopaque marker on the outer sheath passes the proximal radiopaque marker on the inner shaft, stent deployment is complete.


Figure 1. Entrust delivery system

Indications for use

The EverFlex self-expanding peripheral stent with Entrust delivery system is indicated for use in occlusions, lesions at high risk for abrupt closure or threatened closure following percutaneous transluminal angioplasty (PTA); or lesions believed to be at high risk for restenosis following PTA in the common iliac, external iliac, superficial femoral, proximal popliteal, or subclavian arteries. Stenting is intended to improve and maintain artery luminal diameter.

Contraindications

The EverFlex self-expanding peripheral stent with Entrust delivery system is contraindicated under the following conditions:

- Patients for whom anticoagulant or antiplatelet therapy is contraindicated
- Patients who exhibit persistent acute intraluminal thrombus of the proposed lesion site, after thrombolytic therapy
- Perforation at the angioplasty site evidenced by extravasation of contrast medium
- Aneurysm of the artery to be treated
- Patients with known hypersensitivity to nickel titanium
- Patients with uncorrected bleeding disorders
- Use in the carotid artery
- All of the customary contraindications for PTA

Warnings

- This device is provided **sterile** for single use only. Do not reprocess or resterilize. Reprocessing and resterilizing could increase the risks of patient infection and compromised device performance.
- Do not force stent deployment if high resistance is met when initially rotating the thumbwheel. Carefully withdraw the stent system without deploying the stent.
- If resistance is encountered at any time during the insertion procedure, do not force passage. Forced passage may damage the stent system, the vessel lumen, or cause premature deployment or deployment complications. If resistance is felt, carefully withdraw the stent system without deploying the stent.

Precautions

- Carefully inspect the sterile package and the device before use. Do not use the device if the packaging or the device is damaged.
- Carefully inspect the sterile package and device before use to verify that no damage occurred during shipment. Verify that the red locking pin is still engaged with the handle. Do not use the stent or the delivery system if damage is detected.
- Do not exceed 300 psi (2068 kPa) / 20 atm while flushing the delivery system.
- Do not use the stent if it was partially deployed in the package, or if it is partially deployed before starting the deployment procedure.
- Always use an introducer sheath during the implant procedure to protect both the vessel and the puncture site. Support from an introducer sheath is necessary to minimize stent lengthening or shortening during stent deployment.
- Using a 0.014 in (0.36 mm) or 0.018 in (0.46 mm) guidewire in a highly tortuous vessel, highly stenosed lesion, or calcified lesion could result in crossing, tracking, or deployment complications.
- In lesions with a vessel diameter less than 2 mm, crossing or deployment complications may occur.
- Do not torque (turn in multiple rotations) the delivery system. Torquing the delivery system may lead to difficulty or make it impossible to deploy or move the stent system over the guidewire.
- Do not allow slack in the delivery system during deployment. Any slack in the system, inside or outside of the patient, could result in deploying the stent beyond the target site.
- Completely remove the safety locking pin and tube assembly from the handle, as depicted in **Figure 2**. If the safety locking pin and tube assembly are not removed, the EverFlex stent may partially deploy, fail to deploy, or deploy with increased deployment force.
- The stent system is not designed for recapturing at any time, and is not designed for repositioning after establishing vessel apposition.
- Failure to hold the handle in a fixed position may cause foreshortening or lengthening of the stent.
- The stent is not designed to be lengthened or shortened relative to its nominal length. Excessive stent lengthening or shortening may increase the risk of stent fracture.
- Use caution when crossing a deployed stent with any adjunct device.
- Do not expand the EverFlex stent beyond its nominal diameter.

4 Instructions for use

Potential adverse events

The potential adverse effects (or complications) that may occur or require intervention with the use of this device include, but are not limited to, the following conditions:

Abrupt or subacute closure	Dissection or intimal flap	Pseudoaneurysm
Allergic reaction to device materials or procedure medications	Edema	Renal failure requiring dialysis
Allergic reaction to nitinol	Embolism	Renal insufficiency (new or worsening)
Amputation	Failure to deploy stent	Restenosis
Aneurysm	Fever	Sepsis
Angina	Gastrointestinal bleeding due to anticoagulation	Shock
Arrhythmia	Hematoma	Stent collapse or fracture
Arteriovenous fistula	Hypertension or hypotension	Stent migration
Artery perforation or rupture	Infection	Stent misplacement
Bleeding requiring transfusion	Inflammation	Stroke
Bruising	Intraluminal thrombus	Surgical or endovascular intervention
Contrast medium reaction / renal failure	Myocardial infarction	Thrombosis or occlusion of the stent
Death	Pain	Transient ischemic attack
Device breakage	Partial stent deployment	Venous thromboembolism
		Vessel spasm

Directions for use

Preparation

Warning: The device is provided sterile for single use only. Do not reprocess or resterilize. Reprocessing and resterilizing could increase the risks of patient infection and compromised device performance.

Recommended items

- 5 to 10 ml syringe filled with heparinized saline
- 0.035 in (0.89 mm) guidewire
- Hemostatic sheath
- PTA balloon

Select stent size

Measure the diameter of the reference vessel (proximal and distal to lesion). Measure the length of the target lesion. Refer to **Table 1** to select a stent size according to the vessel diameter and a stent length sufficient to extend proximal and distal to the target lesion.

Table 1. Stent diameter and length sizing

Stent diameter (mm)	Reference vessel diameter (mm)	Stent length (mm)
5	3.5 – 4.5	20, 40, 60, 80, 100, 120, 150
6	4.5 – 5.5	20, 40, 60, 80, 100, 120, 150
7	5.5 – 6.5	20, 40, 60, 80, 100, 120, 150
8	6.5 – 7.5	20, 40, 60, 80, 100, 120, 150

Stent delivery system preparation

1. Carefully inspect the pouch to ensure that there is no damage to the sterile barrier. Open the pouch and remove the tray and contents. Place the tray on a flat surface and carefully remove the stent delivery system.
Caution: Carefully inspect the sterile package and device before use to verify that no damage occurred during shipment. Verify that the red locking pin is still engaged with the handle. Do not use the stent or the delivery system if damage is detected.
2. Attach a 5 to 10 ml syringe filled with heparinized saline to the proximal luer lock injection hub. Inject the saline solution through the guidewire lumen until fluid passes from the catheter tip.
Caution: Do not exceed 300 psi (2068 kPa) / 20 atm while flushing the delivery system.
3. Examine the distal end of the catheter to ensure that the stent is fully contained within the tip of the catheter.
Caution: Do not use the stent if it was partially deployed in the package, or if it is partially deployed before starting the deployment procedure.

Stent deployment

Sheath and guidewire insertion

1. Gain access using an introducer sheath with a hemostatic valve that is compatible with a 5 Fr (1.7 mm) delivery system. If the iliac artery is treated, the access site is the common femoral artery. If the subclavian artery is treated, the access site is the brachial artery or axillary artery.
Caution: Always use an introducer sheath during the implant procedure to protect both the vessel and the puncture site. Support from an introducer sheath is necessary to minimize stent lengthening or shortening during stent deployment.
2. Insert a 0.035 in (0.89 mm) guidewire of appropriate length across the target lesion.
Caution: Using a 0.014 in (0.36 mm) or 0.018 in (0.46 mm) guidewire in a highly tortuous vessel, highly stenosed lesion, or calcified lesion could result in crossing, tracking, or deployment complications.

Lesion dilation

Predilate the lesion using standard PTA techniques. Remove the PTA balloon from the patient while maintaining lesion access with the guidewire.

Caution: In lesions where vessel lumen diameter is less than 2 mm, crossing or deployment complications may occur.

Stent delivery system introduction

Advance the stent delivery system over the guidewire through the hemostatic valve and introducer sheath.

Warning: If resistance is encountered at any time during the insertion procedure, do not force passage of the stent system. Forced passage may damage the stent system, the vessel lumen, or cause premature deployment or deployment complications. If resistance is felt, carefully withdraw the stent system without deploying the stent.

Caution: Do not torque (turn in multiple rotations) the delivery system. Torquing the delivery system may lead to difficulty or make it impossible to deploy or move the stent system over the guidewire.

Stent deployment

1. Advance the delivery system until the distal (leading) radiopaque marker is distal to the target lesion.

2. Pull back on the delivery system until there is no slack in the delivery system and until the radiopaque markers extend distal and proximal to the target lesion.

Caution: Do not allow slack in the delivery system during deployment. Any slack in the system, inside or outside of the patient, could result in deploying the stent beyond the target site.

3. Completely remove the safety locking pin and tube assembly by gently pulling in the direction indicated by the arrow on the safety locking pin. Discard the safety locking pin and the attached tube assembly.

Caution: Completely remove the safety locking pin and tube assembly (2) from the handle (1), as depicted in **Figure 2**. If the safety locking pin and tube assembly are not removed, the EverFlex stent may partially deploy, fail to deploy, or deploy with increased deployment force.


Figure 2. Entrust delivery system with safety locking pin and tube assembly

1. Handle
 2. Locking pin and tube assembly
 3. Pull down to remove the safety locking pin and tube assembly
4. Visually confirm that the radiopaque markers still extend distal and proximal to the target lesion. Reposition if necessary.
 5. Initiate stent deployment by holding the handle in a fixed position and slowly rotating the thumbwheel in the direction indicated by the arrow on the handle, see **Figure 2**.

Note: If a second hand is used to hold the delivery system, gently support the delivery system at the gold isolation sheath (1). Do not constrict the silver outer sheath (2) during deployment. See **Figure 3**.


Figure 3. Silver retractable outer sheath and gold isolation sheath

1. Gold
2. Silver

Caution: The stent system is not designed for recapturing at any time, or designed for repositioning after establishing vessel apposition.

6. When the distal radiopaque markers of the stent are approximately 1 to 3 mm distal to marker on the retractable outer sheath, use fluoroscopic guidance to reposition the stent as needed before achieving vessel apposition.
7. Keep the entire length of the flexible deployment system as straight as possible while the stent is being released. Hold the handle in a fixed position to ensure that no slack is introduced to the delivery system. Deployment is complete when the marker on the outer retractable sheath reaches the inner shaft marker and the stent is released.

Warning: Do not force stent deployment if high resistance is met when initially rotating the thumbwheel. Carefully withdraw the stent system without deploying the stent.

Caution: Failure to hold the handle in a fixed position may cause foreshortening or lengthening of the stent.

Caution: The stent is not designed to be lengthened or shortened relative to its nominal length. Excessive stent lengthening or shortening may increase the risk of stent fracture.

Note: If more than one stent is needed, place the distal stent first. If overlap of sequential stents is necessary, keep the amount of overlap to a minimum.

After stent deployment

1. After stent deployment, use fluoroscopic guidance to withdraw the entire delivery system as one unit, over the guidewire, into the introducer sheath and out of the body. Remove the delivery system from the guidewire.
 2. Using fluoroscopy, verify that the EverFlex stent is fully deployed.
 3. If the EverFlex stent is not completely expanded at any point along the lesion, postdeployment balloon dilation can be performed at the discretion of the physician.
- Caution:** Use caution when crossing a deployed stent with any adjunct device.
- Caution:** Do not expand the EverFlex stent beyond its nominal diameter.
4. To dilate the stent, select an appropriately sized PTA balloon catheter and dilate the EverFlex stent using conventional technique. Select a PTA balloon with an inflation diameter that approximates the diameter of the reference vessel.
 5. Confirm that the stent is fully expanded, then remove the PTA balloon from the patient.
 6. Remove the guidewire and introducer sheath from the body.
 7. Close entry wound.
 8. Discard the delivery system, guidewire, and sheath.

MRI information

MR Conditional

Nonclinical testing demonstrated that the EverFlex stent is MR Conditional. A patient with this device can be scanned safely, immediately after stent placement under the following conditions:

- Static magnetic field of 3.0 Tesla or 1.5 Tesla
- Maximum spatial gradient magnetic field of 8000-Gauss/cm (extrapolated) or less (80 T/m)
- Normal operating mode (maximum whole body averaged SAR of 2.0 W/kg) for 15 minutes of scanning

MRI related temperature rise

In nonclinical testing, the EverFlex self-expanding peripheral stent produced the following temperature rises during MRI performed for 15 minutes of scanning (that is, per pulse sequence) in 1.5 Tesla/64 MHz (Magnetom, Siemens Medical Solutions, Malvern, PA. Software Numaris/4, Version Syngo MR 2002B DHHS Active-shielded, horizontal field scanner) and 3.0 Tesla/128-MHz (Excite, HDx, Software 14X.M5, General Electric Healthcare, Milwaukee, WI) MR systems:

Single, 8 mm x 80 mm	1.5-T	3.0-T
	MR system reported, whole body averaged SAR (W/kg)	2.9
Calorimetry measured values, whole body averaged SAR (W/kg)	2.1	2.7
Highest temperature rise (°C)	2.0	2.4
Single, 8 mm x 120 mm	1.5-T	3.0-T
	MR system reported, whole body averaged SAR (W/kg)	2.9
Calorimetry measured values, whole body averaged SAR (W/kg)	2.1	2.7
Highest temperature rise (°C)	2.9	3.7

These temperature changes do not pose a hazard to a patient under the identified test conditions. Patients are recommended to register conditions under which the implant may be scanned safely with the MedAlert Foundation (www.medicalert.org) or equivalent organization.

Artifact Information

MR image quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the stent. Therefore, it may be necessary to optimize MR imaging parameters to compensate for the presence of this device.

Warranty

This product is subject to standard Medtronic warranty terms.

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© 2017 Medtronic. All rights reserved. Medtronic and Medtronic logo are trademarks of Medtronic. ™* Third-party brands are trademarks of their respective owners. All other brands are trademarks of a Medtronic company.
ev3 Inc.
Medtronic Australasia Pty Ltd