Attention User: Before using this product, please read the following important information:

The Creatinine Jaffe kit (CRCO-0600) you have received has being modified.

- For the ELITechGroup ECS recommended procedure, the order of addition of reagents has changed. The Sodium Hydroxide (clear /colourless reagent) is now added first (R1). The Alkaline Picrate (yellow coloured reagent) is now the second reagent (R2) added.
- Labelling has been adapted accordingly, and the red hazard pictogram is now located on the R1 bottle (Sodium hydroxide).
- Place “R1” in first position and “R2” in second position on the Selectra Pro analyser.
- The change in the order of reagent addition is not compatible with older lots of aqueous standard CREN-4050, and vice versa. It is necessary to use the standard provided in this kit.
- This modification requires the uploading of a modified instrument application, which is available in the Documentation section of the ELITechGroup website.
- Adjusted values for Elical and Elitol I & II are also provided hereafter and in the Documentation section of the ELITechGroup website.
- All elements of the CREATININE JAFFE system (reagent, standard, application & control recoveries) need to be updated at the same time.

The modification does NOT impact those users that are combining Reagent 1 and Reagent 2 prior to use (mono reagent procedure).

ELICAL 2 - values for "CRCO improved"				
Elical	vial lot	kit lot	expiration	value
ELICAL 2	R200385	03-3130	2022-03	4.03
ELICAL 2	R200946	03-3140	2022-05	
ELICAL 2	R210511	03-3150	2023-01	4.03
ELICAL 2		03-3160	2023-01	
ELICAL 2	R211194	03-3170	2023-06	4.01

ELITROL 1 - values for "CRCO improved"				
Elitrol	vial lot	kit lot	expiration	target(range)
Elitrol 1	R200718 R201112	01-1110	2022-05	0.84(0.69-0.99)
Elitrol 1	R210334 R210590	01-1120 01-1130	2022-10	0.87(0.71-1.03)
Elitrol 1	R210858	01-1140	2022-12	0.87(0.71-1.03)
Elitrol 1	R211203	01-1150	2023-08	0.92(0.75-1.09)

ELITROL 2 - values for "CRCO improved"				
Elitrol	vial lot	kit lot	expiration	target(range)
Elitrol 2	R200719 R201113	02-2110	2022-05	3.27(2.68-3.86)
Elitrol 2	R210335 R210591	02-2120 02-2130	2022-10	3.38(2.77-3.99)
Elitrol 2	R210859	02-2140	2022-12	3.38(2.77-3.99)
Elitrol 2	R211205	02-2150	2023-08	3.37(2.76-3.98)

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☞BIBLIOGRAFIE

1. Kaplan, J.M., First, M.R., *Renal Function. Clinical Chemistry: Theory, Analysis, Correlation*, 5th Ed., Kaplan, L.A, Pesce, A.J., (Mosby Inc. eds), (2010), 567 and appendix.
2. Delaney, M.P., Price, C.P., Lamb, E.J., *Kidney Function and Disease. Tietz Fundamentals of Clinical Chemistry*, 6th Ed., Burtis, C.A., Ashwood, E.R., Bruns, D.E. (W.B. Saunders eds.), (2008), 631.
3. Wu, A. H. B., *Clinical guide to laboratory tests*, 4th Ed., (W.B. Saunders eds.), (2006), 316.
4. Vasiliades, J., *Clin. Chem.*, (1976), **22**, 1664.
5. Guder, W.G., *et al.*, *Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples*. (2002). WHO/DIL/LAB/99.1 Rev.2.
6. Ceriotti, F. *et al.*, *Reference Intervals for Serum Creatinine Concentrations: Assessment of Available Data for Global Application*, *Clin. Chem.*, (2008), **54**, 559.
7. Young D.S., *Effects of drugs on clinical laboratory tests*, 4th Ed., AACC Press (1995).
8. Berth, M. & Delanghe, J. *Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature*, *Acta Clin Belg.*, (2004), 59, 263.
9. Young, D.S., *Effects of preanalytical variables on clinical laboratory tests*, 2nd Ed., AACC Press (1997).

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PIT-SLSY-EN-V1 (05/2022)

INTENDED USE

This solution is intended for washing probes and cuvettes of Selectra Pro Series analyzers.

This accessory for an *in vitro* diagnostic device is for professional use only.

COMPOSITION

Sodium azide < 0.1 % (w/w)
Also contains surfactants

MATERIALS REQUIRED BUT NOT PROVIDED

- Distilled or deionized water, as indicated in the analyzer's user manual.
- General Laboratory equipment (e.g. pipette).
- Selectra Pro Family analyzers and accessories.
- Do not use materials that are not required as indicated above.

PRECAUTIONS OF USE AND WARNINGS

- Take normal precautions and adhere to good laboratory practice.
- Use clean or single use laboratory equipment only to avoid contamination.
- Solution contains sodium azide which may react with lead or copper plumbing to form potentially explosive metal azides. When disposing of this solution always flush with copious amounts of water to prevent azide buildup.
- Consult Safety Data Sheet (SDS) for a proper handling.

STABILITY

Store at 15-25 °C and protect from light. Do not freeze.

Do not use after expiration date indicated on the vial labels.

Diluted solution stability: 2 weeks.

Open vial stability: until product's expiry date.

PREPARATION

Solution needs to be diluted 1/400 in distilled or deionized water, as indicated in the analyzer's user manual.

E.g. 5 mL in 2 L or 12.5 mL in 5 L or 25 mL in 10 L.

PRODUCT DETERIORATION

- The product should be clear. Cloudiness would indicate product deterioration.
- Do not use the product if there is visible evidence of contamination or damage (e.g. particle matter).
- Damage to the container may impact on product performance. Do not use the solution if there is physical evidence of deterioration (e.g. leakages or punctured container).

INSTALLATION AND USE

Consult Selectra Pro operator manual corresponding.

LIMITATIONS

The use of SYSTEM SOLUTION has been validated for the Selectra Pro Series analyzers.

Using a different analyzer should be validated by the laboratory.

WASTE MANAGEMENT

Disposal of all waste material should be in accordance with local, state and federal regulatory requirements (please refer to the Safety Data Sheet (SDS)).

DECLARATION OF SERIOUS INCIDENT

Please notify the manufacturer (through your distributor) and competent authority of the Member State of the European union in which the user and/or the patient is established, of any serious incident that has occurred in relation to the device. For other jurisdictions, the declaration of serious incident should be in accordance with local, state and federal regulatory requirements.

By reporting a serious incident, you provide information that can contribute to the safety of *in vitro* medical devices.

SYMBOLS

Symbols used are defined in ISO 15223-1 standard, except those presented below :

	Content
	Solution
	Modification from previous version
	European Conformity

TECHNICAL ASSISTANCE

Contact your local distributor or ELITech Clinical Systems SAS (ccsupport@elitechgroup.com).