PulsaClean®

**Disposable Pulse Lavage System
(Disposable Electronically Pulsed Lavage
Suction Apparatus)**

User Manual

Model No.: W-203

REF LD 1000

Guangzhou Clean Medical Products Manufacturing Corp.

No.163 Meidu Road, Chengjiao, Conghua District, Guangzhou, Guangdong,
China. P.C.: 510900

OSARTIS GmbH

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01	en: English	User Manual Disposable Pulse Lavage System (Disposable Electronically Pulsed Lavage Suction Apparatus) Model No.: W-203
02	de: German / Deutsch	Bedienungsanleitung Einweg Pulse Lavage System Modell-Nr.: W-203
03	fr: French / Français	Manuel d'utilisation Système Pulse Lavage Jetable Modèle n°: W-203
04	es: Spanish / Español	Manual del Usuario Sistema de Lavado Pulsátil Desechable N.º de modelo: W-203
05	pt: Portuguese / Português (For Brazil / Pelo Brasil)	Instruções de Uso PULSACLEAN: Sistema de Lavagem Pulsada Descartável N.º do modelo: W-203
06	pt: Portuguese / Português	Manual do Utilizador Sistema de Lavagem Pulsada Descartável N.º do modelo: W-203
07	it: Italian / Italiano	Manuale d'uso Sistema di Lavaggio Pulsato Monouso Modello n°: W-203
08	nl: Dutch / Nederlands	Gebruikershandleiding Wegwerpbaar 'Pulse Lavage'-Systeem Model Nr.: W-203
09	da: Danish / Dansk	Brugervejledning Engangstrykskylning System Model nr.: W-203
10	fi: Finnish / Suomalainen	Käyttöohje Kertakäyttöinen Pulssi Huuhdonta Järjestelmä Mallinumero: W-203
11	sv: Swedish / Svenska	Användarhandbok Pulse Lavage – Engångssystem Modellnummer: W-203
12	pl: Polish / Polski	Instrukcja Użytkownika System do Płukania Pulsacyjnego Jednorazowego Użytku Nr modelu: W-203
13	cs: Czech / Čeština	Uživatelská Příručka Jednorázový Systém Pulzního Výplachu Modelové číslo: W-203

14	et: Estonian / Eesti	Kasutusjuhend Ühekordselt Kasutatav Pulsslavaaži Süsteem Mudeli nr: W-203
15	lv: Latvian / Latviešu	Lietotāja Rokasgrāmata Vienreizlietojama Impulsu Skalošanas Sistēma Modeļa Nr.: W-203
16	lt: Lithuanian / Lietuvių	Vartotojo Vadovas Vienkartinė Impulsinė Praplovimo Sistema Modelio Nr.: W-203
17	hu: Hungarian / Magyar	Felhasználói Kézikönyv Eldobható Pulzáló Öblítőrendszer Modellszám: W-203
18	bg: Bulgarian / Български	Ръководство за Употреба Система за Импулсно Промиване за Еднократна Употреба Модел №: W-203
19	ro: Romanian / Română	Manual de Utilizare Sistem de Lavaj Pulsatil de Unică Folosință Nr. Model: W-203
20	sk: Slovak / Slovák	Použitý Manuál Jednorazový Systém Pulse Lavage Model č.: W-203
21	sl: Slovenian / Slovenščina	Uporabniški Priročnik Sistem Pulznega Izpiranja za Enkratno Uporabo Št. modela: W-203
22	hr: Croatian / Hrvatski	Korisnički Priručnik Sustav za Pulsno Ispiranje za Jednokratnu Upotrebu Model br.: W-203
23	el: Greek / Ελληνικά	Εγχειρίδιο Χρήσης Σύστημα Παλμικών Πλύσεων Μίας Χρήσης Αριθ. Μοντέλο: W-203

PulsaClean®

User Manual

CE 0123

STERILE EO

1. Medical Device

- 1) Device name: **Disposable Pulse Lavage System**
(Disposable Electronically Pulsed Lavage Suction Apparatus)
- 2) Model No.: **W-203**
- 3) **REF: LD 1000**

2. Indications (Intended Use)

- 1) This device is used together with normal sterile saline or other applicable solutions, vacuum extractor. Vacuum extractor is not a must.
- 2) This device can be used for open wound cleaning, soft tissue cleaning, and surgical operation location cleaning.
- 3) This device is suitable to be used in operation room, emergency room and treatment room.

3. Features


- 1) This device is easy to use and portable with significant cleansing effect.
- 2) This device is powered by mercury-free AA alkaline batteries.
- 3) This device is portable and can be installed and handled easily.
- 4) This device is suitable for any kind of open wound cleansing.
- 5) This device cleans the wound completely with proper pressure.
- 6) This device avoids cross contamination due to its disposability.
- 7) This device removes waste solution with suction, if the waste pipe is connected with a vacuum extractor.

4. Working Principle

- 1) This device is portable and internal battery-powered.
- 2) This device consists of DC motor, compression pump, handpiece, enclosure, batteries, and nozzles.
- 3) When squeezing the trigger, power is on. When power is turned on, the DC motor drives the compression pump in reciprocating motion, engender cavity volume changes, thus it can breathe in and out the solution. The power is DC battery. There are a low power position (DC 9V), a high power position (DC 15V) and a stop tap position. Different nozzles can be fitted on during operation.

5. Components

- 1) This device consists of **Handpiece (1)** 【Including DC motor, Compression pump, Trigger (1.1), Nozzle lock (1.2), Enclosure.】 , **Nozzle (2)** 【Such as Fan Spray Nozzle (Short Nozzle) (2a), Femoral Nozzle (Long Nozzle) (2b), Femoral Brush Nozzle (Long Nozzle) (2c), Shower Spray Nozzle (Short Nozzle) (2d), Large Splash Shield (2e) or other nozzles. See the package and labels for details.】 , **Waste pipe (3)** 【Including Clamp (3.1), Connecting pipe (3.2).】 , **Irrigation pipe (4)** 【Including Clamp (3.1), Luer fitting (4.1).】 , **Battery bag(5)**, **Power wire(6)**.
- 2) The **Vacuum extractor (7)** and **Irrigation bag (8)** are not included in this device, but prepared by medical institution. The negative pressure of **Vacuum extractor (7)** within the range of 5kPa to 40kPa shall be adjustable. The **Irrigation bag (8)** should be used with normal sterile saline or other applicable solutions (not suitable for viscous solutions).
- 3) Accessories: **Nozzles (2)**.
- 4) Detachable parts: **Nozzle (2)**.
- 5) Applied parts: **Nozzle (2)**.
- 6) Associated devices: **Vacuum extractor (7)** and **Irrigation bag (8)**.


Device structure diagram

6. Operating Instructions

- 1) Inspect the package. Do not use the device if the package is damaged.
- 2) Consult the user manual and the attached documents.
- 3) Open the package, and install **Nozzle (2)** to **Handpiece (1)**. Press down the **Nozzle lock (1.2)** to lock the **Nozzle (2)**.
- 4) Connect the **Waste pipe (3)** to the **Vacuum extractor (7)** (if prepared). Vacuum extractor is not a must.
- 5) Connect the **Luer fitting (4.1)** to the spike the **Irrigation bag (8)**.
- 6) Hang up or place the **Battery bag (5)**. The **Battery bag (5)** should be hung up on the same horizontal with the **Irrigation bag (8)** and avoid the liquid from the **Irrigation bag (8)** spillage on the battery bag.
- 7) Use the **Nozzle (2)** to aim the treatment area. Squeeze the **Trigger (1.1)** to activate the device and start irrigation. The **Trigger (1.1)** have a low power position with symbol short "—", a high power position with symbol long "——", and a stop tap position with symbol "O".
- 8) Slide **Clamp (3.1)** can control the liquid flow in the **Irrigation pipe (4)** and **Waste pipe (3)**. Please lock the **Clamp (3.1)**, when you stop operating the device.
- 9) It can suck the waste solution concurrently while spraying or after, if **Vacuum extractor (7)** is prepared.
- 10) To press **Nozzle lock (1.2)** upward first, when withdrawing the **Nozzle (2)**.
- 11) Optional: The **Large Splash Shield (2e)** is to be used in conjunction with the short nozzle.
To install: Pinch the soft head of the short nozzle, insert into the opening of the splash shield and slightly pull in opposite directions to fixate the shield on the nozzle head.
To uninstall: Pinch the soft head of the short nozzle and remove the splash shield.

7. Symbol Descriptions

The following symbols will appear on user manual, labels and packages. Some of the symbols represent standards and compliances associated with the device and its use.

 Manufacturer	 Do not re-sterilize	 Do not use if package is damaged and consult instructions for use
 Authorized Representative in the European Community/ European Union	 Date of manufacture	 Fragile, handle with care
 Use-by date	 Batch code (Lot number)	 Keep dry
 Catalogue number	 Do not re-use	 Temperature limit -10°C to +40°C
 Serial number	 Consult instructions for use	 Medical device
 Sterilized using ethylene oxide	 Double sterile barrier system	 Single sterile barrier system with protective packaging inside
 Importer	 Distributor	 Model number
		 Unique device identifier

 CE marking of conformity 0123	 Type B applied part	 This symbol on the trigger side means stop tap position.
 This way up	 Direct current	 Short This symbol on the trigger side means low power position (DC 9V).
 DISPOSAL: Do not dispose this product as unsorted municipal waste. Collection of such waste separately for special treatment is necessary		 Long This symbol on the trigger side means high power position (DC 15V).

8. Specification

No.	Item	Specification
1)	Power source (Ratings)	15V DC, 10×1.5V AA (LR6) alkaline batteries.
2)	Power rating	37.5W.
3)	Handpiece dimension	About 185mm * 55mm * 40mm (Long * Wide * High).
4)	Net weight	About 0.9kg.
5)	Flow rate (high power position)	Not less than 0.9L/minute.
6)	Maximum irrigation length (high power position)	Not less than 2.5 meters.
7)	Noise (Decibel)	Less than 75dB(A).
8)	Solution form	Fog-like and thread-like.
9)	Sterilization	Sterilized using ethylene oxide before delivery.
10)	Valid period (Shelf life)	Three years.

9. IEC Classification

- 1) Internal electrical power source equipment.
- 2) Type B applied part.
- 3) IPX0.
- 4) Not category AP/ APG equipment.

10. Contraindications

Not found.

11. Adverse Effect

Not found.

12. Warning and Precautions

- 1) This device should be operated only by qualified professionals.
- 2) This device should not be used after the expiry date shown on the package.
- 3) This device is for single patient use only. Do not reuse the device. Do not resterilize.
- 4) Prior to each use, inspect the package. Do not use the device if package is damaged.
- 5) This device is not evaluated as category AP or APG equipment. This device is not suitable for use in the presence of flammable anesthetic mixtures with air or oxygen or nitrous oxide.
- 6) No modification of this equipment is allowed.
- 7) Except the Nozzle, there is no replaceable part and accessory within the unit.
- 8) This device is sterilized using ethylene oxide before delivery. Sterilization validity is three years.
- 9) Follow current local/ national regulations to properly dispose the device after use.

13. Disposal of the Batteries

- 1) This device contains alkaline batteries. The battery will be utilizable and safe, if the device is used before date of expiry.
- 2) Follow local/ national current regulations to properly recycle or dispose the batteries after use.
- 3) Do NOT cut the power wire before removing the batteries to prevent short-circuit, which may cause overheating of the battery pack and result in a potential fire risk. In case of short circuit, do not pour liquid over it nor cover it, simply leave it and dispose after it is cooled down. There is a very slim chance of short circuit but proper means of disposal is highly recommended.
- 4) To safely remove the batteries, please press down on flap and remove the lid, then remove batteries from the housing.

14. Maintenance, Cleaning and Sterilization

- 1) This device is maintenance free and does not require routine maintenance.
- 2) This device does not contain user repairable parts and does not require repair.
- 3) This device is for single patient use only. There is no cleaning or sterilization requirement for the user. Do not resterilize.

15. Electromagnetic Compatibility Descriptions

15.1 Conform to the requirements of IEC 60601-1-2:2014 and EN 60601-1-2:2015.

15.2 Instructions for use of electromagnetic compatibility

- 1) The ME (MEDICAL ELECTRICAL) EQUIPMENT or ME SYSTEM is suitable to be used in operation room, emergency room and treatment room.
- 2) **Warning:** Don't near active HF (HIGH FREQUENCY) surgical equipment and the RF (RADIO FREQUENCY) shielded room of an ME system for magnetic resonance imaging, where the intensity of EM (ELECTROMAGNETIC) disturbances is high.
- 3) **Warning:** Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.
- 4) This device doesn't contain user repairable cables, transducers and other ACCESSORIES that are replaceable by the RESPONSIBLE ORGANIZATION.
- 5) **Warning:** Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the Disposable Electronically Pulsed Lavage Suction Apparatus (W-203), including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.
- 6) **NOTE:** The EMISSIONS characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 class A). If it is used in a residential environment (for which CISPR 11 class B is normally required) this equipment might not offer adequate protection to radio-frequency communication services. The user might need to take mitigation measures, such as relocating or re-orienting the equipment.

15.3 Technical description of electromagnetic compatibility

- 1) All necessary instructions for maintaining BASIC SAFETY and ESSENTIAL PERFORMANCE with regard to ELECTROMAGNETIC DISTURBANCES for the EXPECTED SERVICE LIFE. Disposable use product. Do not change any parts of the product.
- 2) In the table below the guidance and manufacturer's declaration -electromagnetic emissions and immunity are shown:

Table 1

Guidance and manufacturer's declaration - electromagnetic emissions	
Emissions test	Compliance
RF emissions (CISPR 11)	Group 1
RF emissions (CISPR 11)	Class A
Harmonic emissions (IEC 61000-3-2)	N/A
Voltage fluctuations/ flicker emissions (IEC 61000-3-3)	N/A

Table 2

Guidance and manufacturer's declaration - electromagnetic immunity		
Immunity test	IEC 60601-1-2 test level	Compliance level
Electrostatic discharge (ESD) IEC 61000-4-2	Contact: ± 8 kV Air: ± 2 kV, ± 4 kV, ± 6 kV, ± 8 kV, ± 15 kV	Contact: ± 8 kV Air: ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV
Electrical fast transient / burst IEC 61000-4-4	Power supply lines: ± 2 kV input/ output lines: ± 1 kV 100 kHz repetition frequency	N/A
Surge IEC 61000-4-5	line(s) to line(s): ± 0.5 kV, ± 1 kV. line(s) to ground: ± 0.5 kV, ± 1 kV, ± 2 kV.	N/A
Voltage dips, short interruptions, and voltage variations on power supply input lines IEC 61000-4-11	0% U_T 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% U_T 1 cycle And	N/A

	70% U _T 25/30 cycles Single phase: at 0 0% U _T 250 cycle	
Power frequency magnetic field IEC 61000-4-8	50Hz: 30 A/m 60Hz: 30 A/m	50Hz: 30 A/m 60Hz: 30 A/m
Conducted RF IEC 61000-4-6	150KHz to 80MHz: 3Vrms 6Vrms (in ISM and amateur radio bands) 80% Am at 1kHz	N/A
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2700 MHz 80 % AM at 1 kHz	3 V/m 80 MHz to 2700 MHz 80 % AM at 1 kHz
NOTE U _T is the a.c. mains voltage prior to application of the test level.		

Table 3

Guidance and manufacturer's declaration - electromagnetic immunity							
	Test frequency (MHz)	Band (MHz)	Service	Modulation	Maximum power (W)	Distance (m)	Immunity test level (V/m)
Radiated RF IEC 61000-4-3 (Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment)	385	380 – 390	TETRA 400	Pulse modulation 18 Hz	1.8	0.3	27
	450	430– 470	GMRS 460, FRS 460	FM ± 5 kHz deviation 1 kHz sine	2	0.3	28
	710	704 – 787	LTE Band 13, 17	Pulse modulation 217 Hz	0.2	0.3	9
	745						
	780						
	810	800 – 960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation 18 Hz	2	0.3	28
	870						
	930						
	1720	1700 – 1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation 217 Hz	2	0.3	28
	1845						
	1970						
	2450	2400 – 2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation 217 Hz	2	0.3	28
	5240	5100 – 5800	WLAN 802.11 a/n	Pulse modulation 217 Hz	0.2	0.3	9
	5500						
	5785						

16. Environmental Requirements

1) Storage conditions:

This device should be stored in non-corrosive gases, cool, dry, and ventilated place.

Temperature: -10°C~+40°C; Humidity: ≤95%RH; Air pressure: 50kPa~106kPa.

2) Operating conditions:

Temperature: 5°C~30°C; Humidity: ≤80%RH; Air pressure: 86kPa~106kPa.

17. Manufacturer

Guangzhou Clean Medical Products Manufacturing Corp.


No.163 Meidu Road, Chengjiao, Conghua District, Guangzhou, Guangdong, China.

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SUNGO Europe B.V.