**Mg****mindray****Magnesium Kit (Xylidyl Blue Method)****Order Information**

Cat. No.	Package size
105-000834-00	R: 4×40 mL
105-000873-00	R: 6×40 mL
105-001609-00	R: 4×45 mL

**Intended Purpose**

In vitro test for the quantitative determination of Magnesium (Mg) concentration in human serum, plasma or urine on Mindray BS series chemistry analyzers. It is intended to be used for aiding to diagnose hypermagnesemia and hypomagnesemia.


**Summary<sup>1-2</sup>**

Magnesium is one of the most abundant cations in the body, and plays an important role in cell respiration, glucose metabolism and transmembrane transport. Magnesium can activate more than 300 kinds of enzymes, the famous of which is Na<sup>+</sup>-K<sup>+</sup>-ATP enzyme.

Hypomagnesemia may cause by defective gastrointestinal absorption, body fluid losses, diuretic therapy, aminoglycoside treatment, endocrinopathy and alcoholism, and the inherited disease is also an infrequent reason. Hypermagnesemia is found in acute and chronic renal failure, magnesium excess, dehydration and diabetic acidosis.

**Assay Principle**

Xylidyl blue method

By using the EGTA to eliminate the interference of calcium, magnesium ions combine with xylidyl blue to produce a xylidyl blue-magnesium complex at an alkaline solution. The absorbency increase is directly proportional to the concentration of magnesium.

**Reagents Components**

<b>R:</b>	Sodium phosphate dibasic dodecahydrate	19.67 g/L
	Trisodium phosphate dodecahydrate	91.22 g/L
	Xylidyl blue	0.058 g/L
	EGTA	0.0495 g/L

**Storage and stability**

Up to expiration date indicated on the label, when stored unopened at 2-8°C and protected from light.

## Mg

On board in use, the reagents are stable for 10 days when refrigerated on the analyzer.

Contamination must be avoided.

Do not freeze the reagent.

### Specimen collection and preparation

#### ■ Specimen types

Serum and lithium heparin plasma are suitable for samples.

Urine Sample is suitable. Collect urine sample in a metal-free container with no preservatives; bring the urine to pH 3-4 prior to assay. Before performing the assay, urine sample should be diluted with NaCl solution 9 g/L (saline) or distilled/deionized water (e.g. 1+14), the result should be multiplied by 15.

#### ■ Preparation for Analysis

1. Use the suitable tubes or collection containers and follow the instruction of the manufacturer; avoid effect of the materials of the tubes or other collection containers.
2. Centrifuge samples containing precipitate before performing the assay.
3. Specimens should be tested as soon as possible after sample collection and pre-analytical treatment.

#### ■ Sample Stability

Serum/Plasma<sup>5</sup>: 7 days at 15-25°C

7 days at 2-8°C

1 year at (-25)-(-15)°C

Urine: 4 days at 15-25°C

7 days at 2-8°C

1 year at (-25)-(-15)°C<sup>5</sup>

For longer storage periods, samples should be frozen at (-20°C)<sup>6,7</sup>. Sample stability claims were established by manufacturer and/or based upon references, each laboratory should establish its own sample stability criteria.

### Reagent Preparation

R is ready to use.

Please perform scheduled maintenance and standard operation including calibration and analysis to assure the performance of measurement system.

### Materials required but not provided

1. General laboratory materials: NaCl solution 9 g/L (saline), distilled/deionized water.
2. Calibrator and Control: Please check the section of reagent instruction of Calibration and Quality Control.
3. Mindray BS series chemistry analyzers and General laboratory equipment.

**Assay procedure**

Parameters Item	BS-800 chemistry analyzers
Assay type	Endpoint
Wavelength (Primary/Secondary)	546/605 nm
Reaction direction	Increase
R	200 $\mu$ L
Mix, incubate at 37°C for 1-2 min, read the absorbance A1, then add:	
Sample or Calibrator	2 $\mu$ L
Mix thoroughly, incubate at 37°C for 3-5 min, read the absorbance A2, Then calculate $\Delta A=(A2-A1)$	

Parameters may vary in different chemistry analyzers, may adjust in proportion if necessary. For Mindray BS series chemistry analyzers, Reagent Parameters is available on request. Please refer to the appropriate operation manual for the analyzers.

**Calibration**

1. It is recommended to use the Mindray Calibrator (Multi Sera Calibrator: 105-001144-00 or other suitable calibrators) and 9 g/L NaCl (saline) for two-point calibration. Traceability of the Mindray Multi Sera Calibrator can refer to the calibrator instructions for use of Mindray Company.

## 2. Calibration frequency

Calibration is stable for approximately 5 days on BS-800 chemistry analyzers. The calibration stability may vary on different instruments, each laboratory should set a calibration frequency in the instrument parameters appropriate to their usage pattern.

Recalibration may be necessary when the following occur:

- As changed reagent lot.
- As required following quality control procedures or out of control.
- As executes specific maintenance or troubleshooting procedure of chemistry analyzers.

3. The calibrator values are lot-specific with the matched models listed in the value sheet.

**Quality control**

1. It is recommended to use the Mindray Control (ClinChem Multi Control: 105-009119-00, 105-009120-00 or other suitable controls) to verify the performance of the measurement procedure; other suitable control material can be used in addition.

## Mg

- Two levels of control material are recommended to analyze each batch of samples. In addition, the control should be run with each new calibration, with each new reagent cartridge, and after specific maintenance or trouble shooting procedures as detailed in the appropriate system manual.
- Each laboratory should establish its own internal quality control scheme and procedures for corrective action if control doesn't recover within the acceptable tolerances.

### Calculation

$\Delta A = (A_2 - A_1)$  sample or calibrator.

$C \text{ sample} = (\Delta A \text{ sample} / \Delta A \text{ calibrator}) \times C \text{ calibrator}$ .

The BS series chemistry analyzer detects the change of absorbance ( $\Delta A$ ) and calculates the Magnesium (Mg) concentration of each sample automatically after calibration.

Conversion factor:  $\text{mg/dL} \times 0.4115 = \text{mmol/L}$ .

### Dilution

If the value of serum/plasma sample exceeds 2.05 mmol/L, the sample should be diluted with 9 g/L NaCl solution (saline) (e.g. 1+1) and rerun; the result should be multiplied by 2.

### Expected values

Sample Type		Units
Serum/Plasma <sup>8</sup>	Adults	0.75-1.02 mmol/L
Urine <sup>3</sup>		3.0-5.0 mmol/24 h

The expected value is provided from references, Mindray has verified the serum expected values by 167 samples of people from China.

Each laboratory should establish its own reference intervals based upon its particular locale and population characteristics since expected values may vary with geography, race, sex and age.

### Performance Characteristics

#### ■ Analytical Sensitivity

The Limit of Blank (LoB) and Limit of detection (LoD) of serum sample were determined on BS-800 in accordance with CLSI Approved Guideline EP17-A2.

LoB=0.03 mmol/L

LoD=0.04 mmol/L

The Limit of Blank is the 95<sup>th</sup> percentile value from  $n \geq 60$  measurements of one or several analyte-free samples over several independent series. The Limit of Blank corresponds to the concentration below which analyte-free samples are found with a probability of 95%.

The Limit of Detection is determined based on the Limit of Blank and the standard deviation of low concentration samples. The Limit of Detection corresponds to the lowest analyte concentration which can be detected (value above the Limit of Blank with a probability of 95%).

### ■ Measuring range

The Mindray BS series systems provide the following linearity range:

Sample Type	Units
Serum / Plasma	0.04-2.05 mmol/L

A high Mg concentration of serum/plasma sample (approximately 2.05 mmol/L) is mixed with a low concentration sample (< 0.04 mmol/L) at different ratios, generating a series of dilutions. The Mg concentration of each dilution is determined using Mindray System, the linearity range is demonstrated with the correlation coefficient  $r \geq 0.990$ . The reportable range of serum/plasma is 0.04-4.1 mmol/L.

For urine, the results of all samples should be multiplied by 15 through auto-analyzer or manual dilution, the measuring range is extended.

### ■ Precision

Precision was determined by following CLSI Approved Guideline EP05-A3<sup>9</sup>, each sample was assayed 2 times per run, 2 runs per day, a total of 20 days. The precision data of controls on BS-800 are summarized below\*.

Specimen Type (N=80)	Mean (mmol/L)	Repeatability		Within-Lab	
		SD (mmol/L)	CV %	SD (mmol/L)	CV %
Control Level 1	0.87	0.005	0.55	0.016	1.87
Control Level 2	1.37	0.006	0.43	0.023	1.71

\*Representative data, results in different instruments, laboratories may vary.

### ■ Analytical Specificity

The samples with different concentration interfering substance were prepared by addition of interferent to human serum pools, and recovers are within  $\pm 10\%$  of the corresponding control value to be considered as no significant interference.

No significant interference was observed when the following substances were tested for interference with this methodology. The data of interference studies on BS-800 are summarized below\*.

Interfering Substance	Interferent Concentration	Analyte Concentration	Relative Deviation
Ascorbic acid	30 mg/dL	0.93 mmol/L	-4.70%
Hemoglobin	500 mg/dL	0.89 mmol/L	+3.38%
Conjugated bilirubin	40 mg/dL	0.91 mmol/L	-1.65%
Unconjugated bilirubin	40 mg/dL	0.89 mmol/L	+1.11%
Intralipid	200 mg/dL	0.92 mmol/L	+7.05%
Triglycerides	37 mmol/L	0.94 mmol/L	+0.64%
Calcium	7.5 mmol/L	0.88 mmol/L	+1.10%

*\*Representative data, results in different instruments, laboratories may vary.*

In very rare cases gammopathy, in particular type IgM, may cause unreliable results<sup>10</sup>.

#### ■ Method Comparison

Correlation studies were performed using CLSI Approved Guideline EP09-A3<sup>11</sup>. The Mindray System (Mindray BS-2800M/Mindray Mg Reagent) (y) was compared with comparison System (Mindray BS-2000/ Mindray Mg Reagent) (x) using the same serum specimens. The statistical data obtained by linear regression are shown in the table below \*:

Sample Type	Regression Fit	Correlation Coefficient (r)	Sample Concentration (N)	Concentration Range (mmol/L)
serum	$y=0.9868x-0.0108$	0.9992	110	0.20-1.91

*\*Representative data, results in different instruments, laboratories may vary.*

#### Result interpretation

The results could be affected by drugs, disease, or endogenous substances<sup>10,12</sup>. When the reaction curve is abnormal, it is recommended to retest and check the result.

#### Warnings and precaution

1. For in vitro diagnostic use only. For laboratory professional use.
2. Please take the necessary precautions for handling all laboratory reagents.
3. Please confirm the integrity of the package before use. Do not use the kits with damaged packages. The reagents avoid direct exposure to sunlight and freezing. The results can't be assured when stored at inappropriate condition.
4. If unintentionally opened before used, store the reagents tightly capped at 2-8°C and protected from light, and the stability is equally to in-use stability.

5. Do not mix reagents with different lots and bottles.  
Do not use the reagents beyond the expiration date and the in-use date.  
Do not mix fresh reagents with in-use reagents.  
Avoid the formation of foam.
6. Instability or deterioration should be suspected if there are visible signs of leakage, precipitates or microbial growth, or if calibration or controls do not meet the insert and/or the Mindray System criteria.
7. Reliability of assay results cannot be guaranteed if the instructions in this package insert are not followed.
8. Do not swallow. Avoid contact with skin and mucous membranes.
9. When the reagents accidentally enter the eyes and mouth, or contact with the skin, immediately wash with plenty of water. If necessary, visit the doctor for further medical treatment.
10. Safety data sheet is available for professional user on request.
11. Disposal of all waste material should be in accordance with local guidelines.
12. All human material should be considered potentially infectious.
13. All identified risks have been reduced as far as possible without adversely affecting the benefit-risk ratio, and the overall residual risk is acceptable.
14. Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.
15. This kit contains components classified as follows in accordance with the Regulation (EC) No 1272/2008:

**Danger**

H314

Causes severe skin burns and eye damage.

**Prevention:**

P260

Do not breathe mist/vapours/spray.

P264

Wash all exposed external body areas thoroughly after handling.

P280

Wear protective gloves, protective clothing, eye protection and face protection.

**Response:**P301+P330  
+P331

IF SWALLOWED: Rinse mouth. Do NOT induce vomiting.

P303+P361 +P353	IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water/shower.
P305+P351 +P338	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P310	Immediately call a POISON CENTER/doctor/physician/first aider
P363	Wash contaminated clothing before reuse.
P304+P340	IF INHALED: Remove person to fresh air and keep comfortable for breathing.
<b>Disposal:</b>	
P501	Dispose of contents/container to authorised hazardous or special waste collection point in accordance with any local regulation.

## References

- Burtis CA, Ashwood ER, Bruns DE, eds. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics, 5th ed. Elsevier Saunders 2012;1751-1755.
- Thomas L. Clinical Laboratory Diagnostics: Use and Assessment of Clinical Laboratory Results. 1st ed. Frankfurt: TH-Books Verlagsgesellschaft 1998; 339-340.
- Wu, Alan HB. Tietz clinical guide to laboratory tests. 4th ed. Elsevier Health Sciences, 2006; 706-709.
- McPherson RA, Pincus MR. Henry's Clinical Diagnosis and Management by Laboratory Methods, 22nd ed. Elsevier Saunders 2007;194-205.
- Ehret W, Heil W, Schmitt Y, Töpfer G, Wisser H, Zawta B, et al. Use of Anticoagulants in Diagnostic Laboratory Investigations and Stability of Blood, Plasma and Serum Samples. WHO/DIL/LAB/99.1 Rev.2:37pp, 47pp.
- CLSI. Procedures for the handling and processing of blood specimens; Approved Guideline-Third Edition. CLSI document H18-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2004.
- CLSI. Urinalysis; Approved Guideline-Third Edition. CLSI document GP16-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2009.
- WS/T 404.6-2015. Reference intervals for common clinical biochemistry tests-Part 6: Serum calcium, phosphorus, magnesium, iron. National Health Commission of the People's Republic of China, 2015.
- CLSI. Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline-Third Edition. CLSI document EP05-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2014.

10. Bakker AJ, Mucke M. Gammopathy interference in clinical chemistry assays: mechanisms, detection and prevention. Clin Chem Lab Med, 2007,45(9):1240-1243.
11. CLSI. Measurement Procedure Comparison and Bias Estimation Using Patient Sample; Approved Guideline-Third Edition. CLSI document EP09-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2013.
12. Young DS. Effects of Drugs on Clinical Laboratory Tests, 5th ed. Washington, DC: AACC Press; 2000:3-544,3-547.

### Graphical symbols

In Vitro Diagnostic  
medical device

Unique device  
identifier

European  
Conformity

Consult Instructions  
For use

Use-by  
date

Authorized representative in  
the European Community

Batch Code

Temperature  
limit

Manufacturer

Catalogue  
number

Keep away from sunlight

Indicates a medical device that needs protection from light sources

© 2022-2024 Shenzhen Mindray Bio-Medical Electronics Co., Ltd. All rights Reserved

**Manufacturer:** Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

**Address:** Mindray Building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057, P. R. China

**E-mail Address:** service@mindray.com

**Website:** www.mindray.com

**Tel:** +86-755-81888998; **Fax:** +86-755-26582680

**EC-Representative:** Shanghai International Holding Corp. GmbH (Europe)

**Address:** Eiffestraße 80, Hamburg 20537, Germany

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1. Burtis CA, Ashwood ER, Bruns DE, eds. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics, 5th ed. Elsevier Saunders 2012; 1751-1755.
2. Thomas L. Clinical Laboratory Diagnostics: Use and Assessment of Clinical Laboratory Results. 1st ed. Frankfurt: TH-Books Verlagsgesellschaft 1998; 339-340.
3. Wu, Alan HB. Tietz clinical guide to laboratory tests. 4th ed. Elsevier Health Sciences, 2006; 706-709.
4. McPherson RA, Pincus MR. Henry's Clinical Diagnosis and Management by Laboratory Methods, 22nd ed. Elsevier Saunders 2007; 194-205.

5. Ehret W, Heil W, Schmitt Y, Töpfer G, Wisser H, Zawta B, et al. Use of Anticoagulants in Diagnostic Laboratory Investigations and Stability of Blood, Plasma and Serum Samples. WHO/DIL/LAB/99.1 Rev.2:37pp, 47pp.
6. CLSI. Procedures for the handling and processing of blood specimens; Approved Guideline-Third Edition. CLSI document H18-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2004.
7. CLSI. Urinalysis; Approved Guideline-Third Edition. CLSI document GP16-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2009.
8. WS/T 404,6-2015. Reference intervals for common clinical biochemistry tests-Part 6: Serum calcium, phosphorus, magnesium, iron. National Health Commission of the People's Republic of China, 2015.
9. CLSI. Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline-Third Edition. CLSI document EP05-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2014.
10. Bakker AJ, Mucke M. Gammopathy interference in clinical chemistry assays: mechanisms, detection and prevention. Clin Chem Lab Med, 2007,45(9):1240-1243.
11. CLSI. Measurement Procedure Comparison and Bias Estimation Using Patient Sample; Approved Guideline-Third Edition. CLSI document EP09-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2013.
12. Young DS. Effects of Drugs on Clinical Laboratory Tests, 5th ed. Washington, DC: AACC Press; 2000:3-544,3-547.

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**Адрес:** Mindray Building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057, P. R. China

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