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
*Disposable pipette for E.S.R. determination according to Westergren method. It is provided by a fiber cap in its upper side which allows the automatic reach of "zero" level in a simple and safe way.*

*In its lower side, at the distance of  $49,5 \pm 1$  mm from the bottom, it is provided with a piston which, thanks to the pressure exerted inside the tube, allows a quick filling of the pipette itself.*

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**CE Marked product - manufactured in compliance with Regulation (EU) 2017/746 (IVDR)**

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**Basic UDI-DI (Global Model Number – GMN): 805577609F109SEDIRATE01FH**

*In Vitro Diagnostic Medical Device for the determination of erythrocyte sedimentation rate (ESR) according to the recommendations of the Westergren method.*

**For professional use in laboratory of analysis only.**

**National classification of medical devices (CND - Italian law) > W02029001 (ERYTHROCYTE SEDIMENTATION RATE DEVICES)**

**EDMA code > 29011001 - Other Hardware + accessories + consumables + software**

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### - For professional use only -

Use only for laboratory diagnostic analysis of human biological samples in healthcare settings.

Do not vary the intended purpose of the device. Please contact the Manufacturer in the event of doubts concerning the intended purpose.

Before use check the perfect closure of the device if applicable.

Handle using personal protective equipment: risk of contamination.

After collection, store the sample as prescribed by the analytical method and deliver to the Laboratory within the maximum time specified by the same.

Do not use the Device in case of evident signs of breakage, injuries, cracks that could compromise its correct use.

Keep out of flame or heat sources which might damage the device.

Do not use after expiry date or if packing is opened.

Do not re-use: Disposable Device. Any reuse could cause sample contamination.

If sterile (ionizing radiations) or aseptic Device: Microbiological status in undamaged pack.

Do not re-sterilize.

Use the Device only with compatible and / or supplied accessories (caps, supports, etc.)

Do not put under direct sun rays; store in a dry, cool place (U.R. max 75% at 26°C)

Store in dry place, temperature range: min. -10°C (14°F) max +50°C (122°F)

Disposal according to applicable regulations, risk of infection.

Before use with particular substances check the resistance/compatibility chart in our catalogue.

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