## HDV Ab

### A. INTENDED USE

Competitive Enzyme Immunoassay (ELISA) for the qualitative determination of antibodies to Hepatitis Delta Virus or HDV in human plasma and sera with a "two-steps" methodology.  
The kit is used for the follow-up of patients infected by HDV.  
For "in vitro" diagnostic use only.

### B. INTRODUCTION

The Hepatitis Delta Virus or HDV is a RNA defective virus composed of a core presenting the delta-specific antigen, encapsulated by HBsAg, that requires the helper function of HBV to support its replication.  
Infection by HDV occurs in the presence of acute or chronic HBV infection. When acute delta and acute HBV simultaneously occur, the illness becomes severe and clinical and biochemical features may be indistinguishable from those of HBV infection alone. In contrast, a patient with chronic HBV infection can support HDV replication indefinitely, usually with a less severe illness appearing as a clinical exacerbation.  
The determination of HDV specific serological markers (HDV Ag, HDV Ab, HDV-IgM and HDV-IgG) represents in these cases an important tool to the clinician for the classification of the etiological agent, for the follow up of infected patients and their treatment. The detection of HDV total antibodies allows the classification of the illness and the monitoring of the seroconversion event.

### C. PRINCIPLE OF THE TEST

Anti-HDV antibodies, if present in the sample, compete with a virus-specific polyclonal IgG, labeled with peroxidase (HRP), for a fixed amount of rec-HDV coated on the microplate. The test is carried out with a two steps incubation competitive system. First the sample is added to the plate and specific anti HDV antibodies bind to the adsorbed antigen. After washing, an enzyme conjugated antibody to HDV is added and binds to the free portion of the antigen coated. After washing a chromogen/substrate mixture is dispensed. The concentration of the bound enzyme on the solid phase becomes inversely proportional to the amount of anti-HDV antibodies in the sample and its activity is detected by the added chromogen/substrate. The concentration of HDV-specific antibodies in the sample is determined by means of a cut-off value that allows for the semi quantitative detection of anti-HDV antibodies.

### D. COMPONENTS

Each kit contains sufficient reagents to perform 96 tests.

#### 1. Microplate: MICROPLATE

8x12 microwell strips coated with recombinant HDV-specific antigen and sealed into a bag with desiccant. Allow the microplate to reach room temperature before opening, reseal unused strips in the bag with desiccant and store at 4°C.

#### 2. Negative Control: CONTROL-

1x2.0ml/vial. Ready to use. Contains goat serum proteins, 100 mM Tris-HCl buffer pH 7.4 +/-0.1, 0.09% Sodium Azide and 0.1% Kathon GC as preservatives. The negative control is colour coded pale yellow.

#### 3. Positive Control: CONTROL+

1x2.0ml/vial. Ready to use. Contains goat serum proteins, high titer anti-HDV antibodies, 100 mM Tris-HCl buffer pH 7.4 +/-0.1, 0.09% Sodium Azide and 0.1% Kathon GC as preservatives. The positive control is colour coded green.

#### 4. Calibrator: CAL

n° 1 vial. Lyophilized. To be dissolved with EIA grade water as reported in the label. Contains bovine serum proteins, low titer bovine antibodies to HDV, 0.2 mg/ml gentamicine sulphate and 0.1% Kathon GC as preservatives.

**Note:** The volume necessary to dissolve the content of the vial may vary from lot to lot. Please use the right volume reported on the label.

#### 5. Wash buffer concentrate: WASHBUF 20X

1x60ml/bottle. 20x concentrated solution.  
Once diluted, the wash solution contains 10 mM phosphate buffer pH 7.0 +/-0.2, 0.05% Tween 20 and 0.1% Kathon GC.

#### 6. Enzyme conjugate: CONJ

1x15ml/vial. Ready-to-use solution. Contains 5% bovine serum albumin, 10 mM Tris buffer pH 6.8 +/-0.1, Horseradish peroxidase conjugated antibody to HDV in presence of 0.2 mg/ml gentamicine sulphate and 0.1% Kathon GC as preservatives. The component is colour coded red.

#### 7. Chromogen/Substrate: SUBS TMB

1x16ml/vial. Contains a 50 mM citrate-phosphate buffered solution at pH 3.5-3.8, 4% DMSO, 0.03% tetra-methyl-benzidine or TMB and 0.02% hydrogen peroxide of H<sub>2</sub>O<sub>2</sub>.

**Notes:** To be stored protected from light, as sensitive to strong illumination.

#### 8. Sulphuric Acid: HSO4 0.3 M

1x15ml/vial. Contains 0.3 M H<sub>2</sub>SO<sub>4</sub> solution.  
Absorbent filter (H313, H314, P280, P302+P352, P332+P313, P305+P351+P338, P337+P313, P362+P363).

#### Plate sealers n° 2

#### Instructions for Use n° 1

### E. MATERIALS REQUIRED BUT NOT PROVIDED

1. Calibrated Micropipettes in the range 10-1000 µl and disposable plastic tips.
2. EIA grade water (double distilled or deionized, charcoal treated to remove oxidizing chemicals used as disinfectants).
3. Timer with 60 minute range or higher.
4. Absorbent paper tissues.
5. Calibrated ELISA microplate thermostatic incubator (dry or wet) set at +37°C.
6. Calibrated ELISA microwell reader with 450nm (reading) and with 620-630nm (blanking) filters.
7. Calibrated ELISA microplate washer.
8. Vortex or similar mixing tools.

### F. WARNINGS AND PRECAUTIONS

1. The kit has to be used by skilled and properly trained technical personnel only, under the supervision of a medical doctor responsible in the laboratory.
2. All the personnel involved in performing the assay have to wear protective laboratory clothes, lab-free gloves and glasses. The use of any sharp (needles) or cutting (blades) devices should be avoided. All the personnel involved should be trained in biosafety procedures, as recommended by the Center for Disease Control, Atlanta, U.S. and reported in the National Institute of Health's publication "Biosafety in Microbiological and Biomedical Laboratories", edn 1984.
3. All the personnel involved in sample handling should be vaccinated for HBV and HAV, for which vaccines are available, safe and effective.

# HDV Ab

**Competitive Enzyme Immunoassay  
for the qualitative determination of  
antibodies to Hepatitis Delta Virus  
in human serum and plasma**

- for "in vitro" diagnostic use only -

**DIA.PRO**  
Diagnostic Bioprobes Srl  
Via G. Carducci n° 27  
20099 Sesto San Giovanni  
(Milano) - Italy  
Phone: +39 02 27007161  
Fax: +39 02 26007716  
e-mail: [info@dipro.it](mailto:info@dipro.it)

REF DAB CE  
96 TESTS

4. The laboratory environment should be controlled so as to avoid contaminants such as dust or air-borne microbial agents when opening kit vials and microplates and when performing the test. Protect the Chromogen/Substrate (TMB/H<sub>2</sub>O<sub>2</sub>) from strong light and avoid vibration of the bench surface where the test is undertaken.
5. Upon receipt, store the kit at +2.8°C into a temperature controlled refrigerator or cold room.
6. Do not interchange components between different lots of the kits. It is recommended that components between two kits of the same lot should not be interchanged.
7. Check that the reagents are clear and do not contain visible heavy particles or aggregates. If not, advise the laboratory supervisor to initiate the necessary procedures.
8. Avoid cross-contamination between serum/plasma samples by using disposable tips and changing them after each sample.
9. Avoid cross-contamination between kit reagents by using disposable tips and changing them between the use of each one.
10. Do not use the kit after the expiration date stated on external (primary container) and internal (vials) labels.
11. Treat all specimens as potentially infective. All human serum specimens should be handled at Biosafety Level 2, as recommended by the Center for Disease Control, Atlanta, U.S. in compliance with what reported in the Institutes of Health's Publication "Biosafety in Microbiological and Biomedical Laboratories", vol. 3, 1984.
12. The use of disposable plastic labware is recommended in the preparation of the washing solution or in transferring components into other containers of automated workstations, in order to avoid contamination.
13. Waste produced during the use of the kit has to be discarded in compliance with national directives and laws concerning laboratory waste of chemical and biological substances. In particular, liquid waste generated from the washing procedure, from residuals of controls and from samples has to be treated as potentially infective material and inactivated. Suggested procedures of inactivation are treatment with a 10% final concentration of household bleach for 16-18 hrs or heat inactivation by autoclave at 121°C for 20 min.
14. Accidental spills have to be absorbed with paper tissues soaked with household bleach and then with water. Tissues should then be discarded in proper containers designated for laboratory/hospital waste.
15. The Sulphuric Acid is an irritant. In case of spills, wash the surface with plenty of water.
16. Other waste materials generated from the use of the kit (example: tips used for samples and controls, used microplates) should be handled as potentially infective and disposed according to national directives and laws concerning laboratory wastes.

**G. SPECIMEN: PREPARATION AND RECOMMENDATIONS**

1. Blood is drawn aseptically by venepuncture and plasma or serum is prepared using standard techniques of preparation of samples for clinical laboratory analysis. No influence has been observed in the preparation of the sample with Citrate, EDTA and heparin.
2. Avoid any addition of preservatives to samples; especially sodium azide as this chemical would affect the enzymatic activity of the conjugate.
3. Samples have to be clearly identified with codes or names in order to avoid misinterpretation of results. When the kit is used for the screening of blood units, bar code labeling and electronic reading is strongly recommended.
4. Haemolysed (red) and visibly hyperfibrinemic ("milky") samples have to be discarded as they could generate false results. Samples containing residues of fibrin or heavy particles or microbial filaments and bodies should be discarded as they could give rise to false results.

5. Sera and plasma can be stored at +2...+8°C for up to five days after collection. For longer storage periods, samples can be stored frozen at -20°C for several months. Any frozen samples should not be frozen/thawed more than once as this may generate particles that could affect the test result.
6. If particles are present, centrifuge at 2,000 rpm for 20 min or filter using 0.2-0.8µm filters to clean up the sample for testing.

**H. PREPARATION OF COMPONENTS AND WARNINGS**

A study conducted on an opened kit has not pointed out any relevant loss of activity up to 6 re-uses of the device and up to 3 months.

**1. Antigen coated microwells:**

Allow the microplate to reach room temperature (about 1 hr) before opening the container. Check that the desiccant has not turned dark or indicating a defect in manufacturing. In this case, call Dia.Pro's customer service. Unused strips have to be placed back into the aluminium pouch, with the desiccant supplied, firmly zipped and stored at +2-8°C. When opened the first time, unused strips are stable until the humidity indicator inside the desiccant bag turns from yellow to green.

**2. Negative Control**

Ready to use. Mix well on vortex before use.

**3. Positive Control:**

Ready to use. Mix well on vortex before use.

**4. Calibrator:**

Low positive control. Add precisely the volume of EIA grade water, reported on its label, to the lyophilized powder, let fully dissolve and then gently mix on vortex.

**Note:** The dissolved calibrator is not stable. Store it frozen in aliquots at -20°C. When thawed do not freeze again; discard it.

**5. Wash buffer concentrate:**

The whole content of the 20x concentrated solution has to be diluted with EIA grade water up to 1200 ml and mixed gently end-over-end before use. During preparation avoid foaming as the presence of bubbles could impact on the efficiency of the washing cycles.

**Note:** Once diluted, the wash solution is stable for 1 week at +2...+8°C.

**6. Enzyme conjugate:**

Ready to use. Mix well on vortex before use.

Avoid contamination of the liquid with oxidizing chemicals, dust or microbes. If this component has to be transferred, use only plastic, and if possible, sterile disposable containers.

**7. Chromogen/Substrate:**

Ready to use. Mix well on vortex before use.

Avoid contamination of the liquid with oxidizing chemicals, air-driven dust or microbes; Do not expose to strong light; oxidizing agents and metallic surfaces. If this component has to be transferred use only plastic, and if possible, sterile disposable container.

**8. Sulphuric Acid**

Ready to use. Mix well on vortex before use.

Attention! Irritant (H315, H319; P280, P302+P352, P332+P313, P305+P351+P338, P337+P313, P362+P363).

**Legends:**

**Warning H statements:**  
H315 - Causes skin irritation.  
H319 - Causes serious eye irritation.

**Precautionary P statements:**  
P280 - Wear protective gloves/protective clothing/eye protection/face protection.  
P302 + P352 - IF ON SKIN: Wash with plenty of soap and water.  
P332 + P313 - If skin irritation occurs: Get medical advice/attention.  
P305 + P351 + P338 - IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.  
P337 + P313 - If eye irritation persists: Get medical advice/attention.  
P362 + P363 - Take off contaminated clothing and wash it before reuse.

**L. PRE ASSAY CONTROLS AND OPERATIONS**

1. Check the expiration date of the kit printed on the external label (primary container). Do not use if expired.
2. Check that the liquid components are not contaminated by visible particles or aggregates. Check that the Chromogen/Substrate is colorless or pale blue by aspirating a small volume of it with a sterile plastic pipette. Check that no leakage occurred in transportation and no spillage of liquid is present inside the box (primary container). Check that the aluminium pouch, containing the microplate, is not punctured or damaged.
3. Dilute all the content of the 20x concentrated Wash Solution as described above.
4. Dissolve the Calibrator as described above and gently mix (about 1 hr) and then mix gently on vortex all liquid reagents.
5. Set the ELISA incubator at +37°C and prepare the ELISA washer by priming with the diluted washing solution according to the manufacturer's instructions. Set the right number of washing cycles as found in the validation of the instrument for its use with the kit.
6. Check that the ELISA reader is turned on or ensure it will be turned on at least 20 minutes before reading.
7. If using an automated work station, turn on, check settings and be sure to use the right assay protocol.
8. Check that the microplates are set to the required volume.
9. Check that all the other equipment is available and ready to use.
10. In case of problems, do not proceed further with the test and advise the supervisor.

**M. ASSAY PROCEDURE**

The assay has to be carried out according to what reported below, taking care to maintain the same incubation time for all the samples in testing.

1. Place the required number of strips in the microplate holder. Leave A1 well empty for the operation of blanking.  
Store the other strips into the bag in presence of the desiccant at +2...+8°C, sealed.
2. Pipette 100 µl of Negative Control in triplicate, 100 µl Positive Control in single and then 100 µl of samples. Check that controls and samples have been correctly added.  
Then incubate the microplate at +37°C for 60 min.
3. Wash the microplate as reported in section I.3.
4. In all the wells except A1, pipette 100 µl Enzyme Conjugate. Check that the reagent has been correctly added.  
Then incubate the microplate at +37°C for 60 min.

**Important note:** Be careful not to touch the inner surface of the well with the pipette tip when dispensing the Enzyme Conjugate. Contamination might occur.

5. Wash the microplate as described.

6. Pipette 100 µl TMB/H<sub>2</sub>O<sub>2</sub> mixture in each well, the blank wells included. Check that the reagent has been correctly added. Then incubate the microplate at room temperature for 20 min.

**Important note:** Do not expose to strong direct light as a high background might be generated.

7. Pipette 100 µl Sulphuric Acid into all the wells using the same pipetting sequence as in step n° 6 to stop the enzymatic reaction. Addition of the stop solution will turn the negative control and negative samples from blue to yellow.

8. Measure the colour intensity of the solution in each well, as described in section 1.5 using a 450nm filter (reading) and a 620-630nm filter (background subtraction, strongly recommended), blanking the instrument on A1.

**Important notes:**

- If the second filter is not available, ensure that no finger prints are present on the bottom of the microwell before reading at 450nm. Finger prints could generate false positive results on reading. Reading should be performed immediately after 20 minutes afterwards. Some self-oxidation of the chromogen can occur leading to a higher background.
- The use of the Calibrator, a low positive control, is not mandatory for the assay as the CAL does not enter into the colour calculation. The CAL may be used as a low filter positive control when a laboratory internal quality verification is required by the management. When used for such purpose, dispense 100 µl of it, possibly in duplicate.

**N. ASSAY SCHEME**

Controls/Calibrator Samples	100 µl
1 <sup>st</sup> incubation	100 µl
Temperature	60 min
Washing step	+37°C
Enzyme Conjugate	4-5 cycles
2 <sup>nd</sup> incubation	100 µl
Temperature	60 min
Washing step	+37°C
TMB/H <sub>2</sub> O <sub>2</sub> mix	4-5 cycles
3 <sup>rd</sup> incubation	100 µl
Temperature	20 min
Sulphuric Acid	r.l.
Reading OD	100 µl
	450nm & 620nm

An example of dispensation scheme (including CAL) is reported in the table below:

	1	2	3	4	5	6	7	8	9	10	11	12
A	BLK	S2										
B	NC	S3										
C	NC	S4										
D	NC	S5										
E	CAL	S6										
F	CAL	S7										
G	PC	S8										
H	ST	S9										

Legend: BLK = Blank NC = Negative Control  
 CAL = Calibrator PC = Positive Control S = Sample

**O. INTERNAL-QUALITY CONTROL**

A check is performed on the negative and positive controls any time, and on the Calibrator in addition when the kit is used for the first time, in order to verify whether the expected OD450nm or Co/S values have been matched in the analysis. Ensure that the following parameters are met:

Parameter	Requirements
Blank well	< 0.100 OD450nm value
Negative Control (NC)	> 1.000 OD450nm after blanking procedure and decrease the number of cycles of the soaking time
Positive Control (PC)	coefficient of variation < 30%
Calibrator (CAL)	OD450 nm < NC/10 PC < OD450nm < (NC+PC)/5

If the results of the test match the requirements stated above, proceed to the next section. If they don't, do not proceed any further and perform the following checks:

Problem	Check
Blank well > 0.100 OD450nm	that the Chromogen/Substrate solution has not become contaminated during the assay
Negative Control (NC) < 1.000 OD450nm	1. that the washing procedure and the washer settings are as validated in the pre qualification study; 2. that the proper washing solution has been used and the washer has been primed with it before use; 3. that no mistake has been done in the assay procedure (dispensation of positive control instead of negative control); 4. that no contamination of the negative control or of the wells where the control was dispensed has occurred due to positive samples, to spills or to the enzyme conjugate;
Blanking coefficient of variation > 30%	5. that the enzyme reagent has not become contaminated with positive samples or with the enzyme conjugate; 6. that the washer needles are not blocked or partially obstructed.
Calibrator OD450nm Outside the range	1. that the procedure has been correctly performed; 2. that no mistake has occurred during its distribution (ex.: dispensation of negative control instead of Calibrator); 3. that the procedure and the washer settings are as validated in the pre qualification study; 4. that no external contamination of the calibrator has occurred.
Positive Control OD450nm > NC/10	1. that the procedure has been correctly performed; 2. that no mistake has occurred during the distribution of the control (dispensation of negative control instead of positive control); 3. that the procedure and the washer settings are as validated in the pre qualification study; 4. that no external contamination of the positive control has occurred.

If any of the above problems have occurred, report the problem to the supervisor for further actions.

**P. RESULTS**

The results are calculated by means of a cut-off value determined with the following formula:

$$\text{Cut-Off} = (\text{NC} + \text{PC}) / 5$$

**Important note** When the calculation of results is performed by the operating system of an ELISA automated work station,

ensure that the proper formula is used to calculate the cut-off value and generate the correct interpretation of results.

**Q. INTERPRETATION OF RESULTS**

Results are interpreted as ratio between the cut-off value and the sample OD450nm or Co/S. Results are interpreted according to the following table:

Co/S	Interpretation
< 0.9	Negative
0.9 - 1.1	Equivocal
> 1.1	Positive

A negative result indicates that the patient has not been infected by HDV. Any patient showing an equivocal result should be re-tested on a second sample taken 1-2 weeks after the initial sample. A positive result is indicative of HDV infection and therefore the patient should be treated accordingly.

**Important notes:**

- Interpretation of results should be done under the supervision of the laboratory supervisor to reduce the risk of judgement errors and misinterpretations.
- When test results are transmitted from the laboratory to data transfer, attention must be paid to avoid erroneous data transfer.
- Diagnosis of viral hepatitis infection has to be taken by and released to the patient by a suitably qualified medical doctor.

An example of calculation is reported below.

The following data must not be used instead of real figures obtained by the user:

Negative Control: 2.100 - 2.200 - 2.000 OD450nm  
 Mean Value: 2.100 OD450nm  
 Higher than 1.000 - Accepted

Positive Control: 0.100 OD450nm  
 Lower than NC/10 - Accepted

Cut-Off = (2.100 + 0.100) / 5 = 0.440

Calibrator: 0.300-0.260 OD450nm  
 Mean value: 0.280 OD450nm  
 Within the range PC ≤ OD450nm < (NC+PC)/5 - Accepted

Sample 1: 0.020 OD450nm  
 Sample 2: 1.900 OD450nm  
 Sample 1 Co/S > 1.1 positive  
 Sample 2 Co/S < 0.9 negative

**R. PERFORMANCES**

Evaluation of Performances has been conducted in accordance to what reported in the Common Technical Specifications or CTS (art. 5, Chapter 3 of IVD Directive 98/79/EC).

**1. LIMIT OF DETECTION:**

In absence of an international standard, the sensitivity of the assay has been calculated by means of the product, named Accurn n° 127 supplied by Boston Biomedica Inc. - USA.

The table below reports the OD450nm shown by this preparation, when diluted in fetal calf serum to prepare a limiting dilution curve, in three different lots.

**Co/S values**

Accurn # 127	Lot# 1192		Lot# 0103		Lot# 0403	
	Co/S value	mm	Co/S value	mm	Co/S value	mm
1x	0.171	3.0	0.163	2.9	0.156	2.8
2x	0.187	2.7	0.176	2.6	0.179	2.5
4x	0.230	2.2	0.230	2.1	0.202	2.2
8x	0.298	1.7	0.285	1.6	0.271	1.6
16x	0.417	1.2	0.405	1.1	0.402	1.1
32x	0.514	0.9	0.490	0.9	0.482	0.9
64x	0.717	0.7	0.700	0.7	0.705	0.6
128x	1.063	0.5	1.066	0.5	1.015	0.4
CTRL H	2.484	0.000	2.261	0.000	2.114	0.000

**2. DIAGNOSTIC SPECIFICITY AND SENSITIVITY**

The diagnostic performances were evaluated in a clinical trial organized by the Department of Gastro-Hepatology, Prof. M. Rizzetto, S. Giovanni Battista hospital, Torino, Italy, on more than 400 samples against a reference kit.

Negative, positive and potentially interfering samples were examined in the trial. Both plasma, derived with different standard techniques of preparation (citrate, EDTA and heparin), and sera have been used to determine the specificity. No false reactivity due to the method of specimen preparation has been observed. Results are briefly reported in the tables below:

Sensitivity > 98 %
Specificity > 98 %

**3. PRECISION**

The mean values obtained from a study conducted on two samples of different anti-HDV antibody reactivity, examined in 16 replicates in three separate runs for three lots of product, is reported below.

DAB.CE: lot #1102

Mean Value	Negative Control (N = 16)		
	1st run	2nd run	3rd run
OD 450nm	2.342	2.428	2.483
Std Deviation	0.113	0.106	0.122
CV %	4.8	4.4	5.0

Mean Value	Calibrator (N = 16)		
	1st run	2nd run	3rd run
OD 450nm	0.288	0.286	0.291
Std Deviation	0.023	0.027	0.026
CV %	7.7	9.3	9.1
Co/S	1.6	1.7	1.7

DAB.CE: lot #0103

Mean Value	Negative Control (N = 16)		
	1st run	2nd run	3rd run
OD 450nm	2.208	2.237	2.246
Std Deviation	0.105	0.108	0.108
CV %	4.7	4.8	4.8

Calibrator (N = 16)

Mean value	1st run	2nd run	3 <sup>rd</sup> run	Average value
OD 450nm	0.269	0.277	0.265	0.271
Std.Deviation	0.025	0.024	0.025	0.025
CV %	9.8	8.5	9.5	9.1
Co.S	1.7	1.7	1.7	1.7

DAB.CE: lot # 0403

Negative Control (N = 16)

Mean value	1st run	2nd run	3 <sup>rd</sup> run	Average value
OD 450nm	2.245	2.221	2.162	2.216
Std.Deviation	0.037	0.03	0.118	0.105
CV %	4.3	4.8	5.4	4.8

Calibrator (N = 16)

Mean value	1st run	2nd run	3 <sup>rd</sup> run	Average value
OD 450nm	0.285	0.273	0.280	0.280
Std.Deviation	0.027	0.023	0.026	0.025
CV %	9.3	8.5	9.1	9.0
Co.S	1.6	1.7	1.6	1.6

The variability shown in the tables did not result in sample misclassification.

**S. LIMITATIONS**

Bacterial contamination or heat inactivation of the specimen may affect the absorbance values of the samples with consequent alteration of the level of the analysis. This test is suitable only for testing single samples and not pooled ones.

Diagnosis of an infectious disease should not be established on the basis of a single test result. The patient's clinical history, symptomatology, as well as other diagnostic data should be considered.

**REFERENCES**

1. Engvall E. and Perlmann P., *J. Immunochimistry* 8: 871-874, 1971
2. Engvall E. and Perlmann P., *J. Immunol.*, 109: 129-135, 1971
3. Chaggar K. Et al., *Journal of Virological Methods*, 32: 193-199, 1991
4. Lazinski D.W. et al., *Journal of Virol.*, 67: 2672-2680, 1993
5. Govindarajan S. et al., *Microbiol. And Immunol.*, 95: 140-141, 1990
6. Shallock A.G. et al., *J. Clin. Microbiol.*, 29: 1873-1876, 1991
7. Forbes B.A. et al., *Clin. Microbiol. News.*, 13: 52-54, 1991
8. Bergmann, K. et al. *J. Immunol.*, 143,3714-3721, 1989
9. Bergmann, K. et al. *J. Infect. Dis.*, 154-702-706, 1986
10. Buti, M. et al. *Hepatology* 8:1125-1129, 1988
11. Rizzetto, M., *Hepatology* 3729-737, 1983
12. Rizzetto, M. et al. *Proc. Natl Acad. Sci. USA* 77:6124-6128, 1980
13. Dubois, F. et al. *J. Clin. Microbiol.*, 26:1339-1342, 1988
14. Wang, K. et al. *Nature* 323:508-514, 1986

All the IVD Products manufactured by the company are under the control of a certified Quality Management System approved by an EC Notified Body. Each lot is submitted to a quality control and released into the market only if conforming with the EC technical specifications and acceptance criteria.

Manufacturer:  
Dia Pro Diagnostic Bioprobes Srl  
Via G. Carducci n° 27 – Sesto San Giovanni (MI) – Italy

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## EC DECLARATION OF CONFORMITY

Lorne Laboratories Ltd declares that the following in vitro diagnostic reagent:

Product name	Catalogue number
CRP Latex kit	850100A

has been classified as non List A, non List B (Directive 98/79/EC, Annex II) and complies with the essential requirements and provisions of Directive 98/79/EC of the European Parliament and of the Council (also SI 2002 No.618 which transposes the requirements of Directive 98/79/EC).

and is in conformity with the national standards transposing harmonised standards:

- BS EN 980:2008
- BS EN ISO 13485:2012
- BS EN 13612:2002
- BS EN 13640:2002
- BS EN 13641:2002
- BS EN ISO 14971:2012
- BS EN ISO 18113, parts 1&2

The conformity assessment procedure performed was in accordance with Annex III of Directive 98/79/EC.

This declaration of conformity is issued under the sole responsibility of Lorne Laboratories Ltd and is valid from 13 April 2016.


Eddy Velthuis  
Technical Director

File No A12241;  
ISO 13485:2003; ISO 9001:2008

Lorne Laboratories Limited | Tel: +44 (0) 118 921 2264  
Unit 1 Cutbush Park Industrial Estate | Fax: +44 (0) 118 986 4518  
Danehill, Lower Earley | Email: info@lornelabs.com  
Berkshire RG6 4UT United Kingdom | [www.lornelabs.com](http://www.lornelabs.com)

Registered office as above. Registered in England No. 04540797 VAT No. 800 3655 66

**LORNE**  
LABORATORIES

## EC DECLARATION OF CONFORMITY

Lorne Laboratories Ltd declares that the following in vitro diagnostic reagent:

Product name	Catalogue number
ASO Latex kit	031100A

has been classified as non List A, non List B (Directive 98/79/EC, Annex II) and complies with the essential requirements and provisions of Directive 98/79/EC of the European Parliament and of the Council (also SI 2002 No.618 which transposes the requirements of Directive 98/79/EC).

and is in conformity with the national standards transposing harmonised standards:

- BS EN 980:2008
- BS EN ISO 13485:2012
- BS EN 13612:2002
- BS EN 13640:2002
- BS EN 13641:2002
- BS EN ISO 14971:2012
- BS EN ISO 18113, parts 1&2

The conformity assessment procedure performed was in accordance with Annex III of Directive 98/79/EC.

This declaration of conformity is issued under the sole responsibility of Lorne Laboratories Ltd and is valid from 13 April 2016.

Eddy Velthuis  
Technical Director


File No A12241  
ISO 13485:2003; ISO 9001:2008

Lorne Laboratories Limited  
Unit 1 Cutbush Park Industrial Estate  
Danehill, Lower Earley  
Berkshire RG6 4UT United Kingdom

Tel: +44 (0) 118 921 2264  
Fax: +44 (0) 118 986 4518  
Email: [info@lornelabs.com](mailto:info@lornelabs.com)  
[www.lornelabs.com](http://www.lornelabs.com)

Registered office as above. Registered in England No. 04540797. VAT No. 800 3655 60

**LORNE**  
LABORATORIES

## EC DECLARATION OF CONFORMITY

Lorne Laboratories Ltd declares that the following in vitro diagnostic reagent:

Product name	Catalogue number
RF Latex kit	830100A

has been classified as non List A, non List B (Directive 98/79/EC, Annex II) and complies with the essential requirements and provisions of Directive 98/79/EC of the European Parliament and of the Council (also SI 2002 No.618 which transposes the requirements of Directive 98/79/EC).

and is in conformity with the national standards transposing harmonised standards:

- BS EN 980:2008
- BS EN ISO 13485:2012
- BS EN 13612:2002
- BS EN 13640:2002
- BS EN 13641:2002
- BS EN ISO 14971:2012
- BS EN ISO 18113, parts 1&2

The conformity assessment procedure performed was in accordance with Annex III of Directive 98/79/EC.

This declaration of conformity is issued under the sole responsibility of Lorne Laboratories Ltd and is valid from 13 April 2016.

Eddy Velthuis  
Technical Director

File No A12241:  
ISO 13485:2003; ISO 9001:2008

Lorne Laboratories Limited  
Unit 1 Cutbush Park Industrial Estate  
Danehill, Lower Earley  
Berkshire RG6 4UT United Kingdom

Tel: +44 (0) 118 921 2264  
Fax: +44 (0) 118 986 4518  
Email: [info@lornelabs.com](mailto:info@lornelabs.com)  
[www.lornelabs.com](http://www.lornelabs.com)

Registered office as above. Registered in England No. 04540797. VAT No. 800 3655 66

**LORNE LABORATORIES LTD.**  
GREAT BRITAIN

**RAPID LATEX KIT**  
DIRECTIONS FOR USE

**ASO Latex Kit: For Detection Of Anti-Streptolysin O (ASO) In Serum.**

**SUMMARY**

In acute streptococcal infections, the toxic immunogenic exoenzyme streptolysin O (ASO) is produced in response to streptolysin O antigens liberated by haemolytic streptococci of groups A, C and G. Information on extent and degree of infection can be obtained by measuring serum ASO levels. Elevated ASO levels have also been found in patients suffering from scarlet fever, acute rheumatoid arthritis, tonsillitis, and other streptococcal infections as well as in healthy carriers.

**PRINCIPLE**

When used by the recommended techniques, latex particles in the reagent will agglutinate (clump) in the presence of anti-streptolysin O antibodies. No agglutination generally indicates the absence of anti-streptolysin O antibodies (see Limitations).

**KIT DESCRIPTION**

Lorne ASO Latex Kit is a serologic test for the detection of ASO antibodies. All the reagents are supplied at optimal dilution for use with all recommended techniques without the need for further dilution or addition. For lot reference number and expiry date see Vial Labels.

**STORAGE**  
Do not freeze. Reagent vials should be stored at 2 - 8°C on receipt. Prolonged storage at temperatures outside this range may result in accelerated loss of reagent reactivity.

**SPECIMEN COLLECTION**

Specimens should be drawn with or without anticoagulant using an aseptic phlebotomy technique. If testing is delayed specimens can be stored at 2-8°C for 7 days or for up to 3 months at or below -20°C. Specimens must be free from bacterial contamination, fibrin, gross lipemia and gross haemolysis.

**PRECAUTIONS**

- The kit is for *in vitro* diagnostic use only.
- Do not use kit past expiration date (see Vial and Box Labels).
- Protective clothing should be worn when handling the reagents, such as disposable gloves and a laboratory coat.
- Materials used to produce the kit were tested at source and found to be negative for HIV 1+2 and HBsAg using approved microbiological tests. However no known tests can guarantee that products derived from human or animal sources are free from infectious agents. Care must be taken in the use and disposal of each vial and its contents.

**DISPOSAL OF KIT REAGENT AND DEALING WITH SPILLAGES**

For information on disposal of kit reagent and decontamination of a spillage site see **Material Safety Data Sheets**, available on request.

**CONTROLS AND ADVICE**

- ASO Positive and Negative Controls must be tested in parallel with each batch of tests. Tests must be considered invalid if controls do not show expected results.
- All the reagents must be allowed to reach 18-25°C before use.
- Shake the reagents well before use to ensure homogeneity.
- Do not interchange components between different kits.
- Use of kit and interpretation of results must be carried out by properly trained and qualified personnel in accordance with the requirements of the country where the kit is in use.
- The user must determine the suitability of the kit for use in other techniques.

**KIT COMPONENTS SUPPLIED**

- ASO Latex Reagent: 5.0 mL Latex particles coated with streptolysin O, pH 6.2 containing a preservative.

**STABILITY OF THE REACTIONS**

Slide tests should be interpreted immediately after the 2-minute period to avoid the possibility that a negative result may be incorrectly interpreted as positive due to drying of the reagent.

**LIMITATIONS**

- False positive results may be obtained in conditions such as scarlet fever, acute rheumatic arthritis, tonsillitis and other streptococcal infections as well as in healthy carriers.
- Haemoglobin ( $\leq 10$  g/L), bilirubin ( $\leq 20$  mg/dL), lipemia ( $\leq 10$  g/L), rheumatoid factors ( $\leq 300$  IU/mL) do not interfere. Other substances may interfere.
- Early infections in children from 6 months to 2 years may cause false negative results.
- A single ASO determination does not produce much information about the actual state of the disease. Titrations at biweekly intervals during 4 or 5 weeks are advisable to follow the disease evolution.
- Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.
- False positive or false negative results may also occur due to:
  - Contamination of test materials
  - Improper storage of test materials or omission of reagents
  - Deviation from the recommended techniques

**SPECIFIC PERFORMANCE CHARACTERISTICS**

- The kit has been characterised by all the procedures mentioned in the **Recommended Techniques**.
- Prior to release, each lot of Lorne ASO Latex Kit is tested by **Recommended Techniques** to ensure suitable reactivity.
- Analytical sensitivity: 200 ( $\pm 50$ ) IU/mL, under the described assay conditions.
- Prozone effect: No prozone effect was detected up to 1500 IU/mL.
- Diagnostic sensitivity: 98 %.
- Diagnostic specificity: 97 %.

**DISCLAIMER**

- The user is responsible for the performance of the kit by any method other than those mentioned in the **Recommended Techniques**.
- Any deviations should be validated prior to use using established laboratory procedures.

**BIBLIOGRAPHY**

- Hefjeje, Quarterly Journal of Medicine 1992, New series 84; 303: 641-656.
- Amjed Samir et al. Pediatric Annals 1992; 21: 835-842.
- Spaan J et al. Bull Wild Health Org 1961; 24: 271-279.
- The association of Clinical Pathologists 1961, Broadsheet 34.
- Picard B et al. La Presse Medicale 1983; 23: 2-6.
- Klein GC. Applied Microbiology 1971; 21: 999-1001.
- Young DS. Effects of drugs on clinical laboratory test. 4th ed. AACCPress, 1995.

**AVAILABLE KIT SIZES**

Kit Size	Catalogue Number
100 Tests Per Kit:	
5.0 ml ASO Latex	031100A
1.0 ml Positive Control	
1.0 ml Negative Control	

For the availability of other sizes, please contact:

**Lorne Laboratories Limited**  
Unit 1 Cubush Park Industrial Estate  
Danehill

Lower Earley  
Berkshire, RG6 4UT  
England  
Tel: +44 (0) 118 921 2264  
Fax: +44 (0) 118 966 4518  
E-mail: [info@lornelabs.com](mailto:info@lornelabs.com)

**LORNE LABORATORIES LTD.**  
GREAT BRITAIN

**RAPID LATEX KIT**  
**DIRECTIONS FOR USE**

**CRP Latex Kit: For Detection of C-Reactive Protein (CRP) in Serum.**

**SUMMARY**

C-Reactive Protein (CRP) usually appears in serum of individuals in response to inflammatory conditions and tissue necrosis, and disappears when causative conditions subside. It is routinely found in cases of bacterial infection, acute rheumatic fever, and many malignant diseases and is often seen in association with rheumatoid arthritis, viral infections and tuberculosis. CRP has also been detected in patients following blood transfusions and surgical operations as well as in patients with burns, peniphigus vulgaris and other bulbous lesions.

**PRINCIPLE**

When used by the recommended techniques, latex particles in the reagent will agglutinate (clump) in the presence of CRP. No agglutination generally indicates absence of CRP (see Limitations).

**KIT DESCRIPTION**

Lorne CRP Latex Test Kit is for the detection of CRP. Test reagent consists of latex particles coated with rabbit Anti-CRP (IgG). All the latex reagents are supplied at optimal dilution for use with all recommended techniques without need for further dilution or addition. For lot reference number and expiry date see Vial Labels.

**STORAGE**

Do not freeze. Reagent vials should be stored at 2 - 8°C on receipt. Prolonged storage at temperatures outside this range may result in accelerated loss of reagent reactivity.

**SPECIMEN COLLECTION**

Specimens should be drawn with or without anticoagulant using an aseptic phlebotomy technique. If testing is delayed specimens can be stored at 2-8°C for 7 days or for up to 3 months at or below -20°C. Specimens must be free from bacterial contamination, fibrin, gross lipaemia and gross haemolysis.

**PRECAUTIONS**

- The kit is for *in vitro* diagnostic use only.
- Do not use kit past expiration date (see Vial and Box Label).
- Protective clothing should be worn when handling the reagents, such as disposable gloves and a laboratory coat.
- Materials used to produce the kit were tested at source and found to be negative for HIV 1+2, and HBSAg using approved microbiological tests. However, no known tests can guarantee that products derived from human or animal sources are free from infectious agents. Care must be taken in the use and disposal of each vial and its contents.

**DISPOSAL OF KIT REAGENT AND DEALING WITH SPILLAGES**

For information on disposal of kit reagent and decontamination of a spillage site see **Material Safety Data Sheets**, available on request.

**CONTROLS AND ADVICE**

- It is recommended the CRP Positive and Negative Controls be tested in parallel with each batch of tests. Tests must be considered invalid if controls do not show expected results.
- All the reagents must be allowed to reach 18-25°C before use.
- Shake the reagents well before use to ensure homogeneity.
- Do not interchange components between different kits.
- Use of kit and interpretation of results must be carried out by properly trained and qualified personnel in accordance with the requirements of country where kit is in use. The user must determine the suitability of the kit for use in other techniques.

**KIT COMPONENTS SUPPLIED**

- CRP Latex Reagent (5 ml.) Latex particles coated with goat IgG anti-human CRP, pH 6.2 containing a preservative.
- CRP Positive Control (3mg dsp., 1 ml.) Human serum with a CRP concentration > 20 mg/L containing a preservative.

**INTERPRETATION OF SEMI-QUANTITATIVE RESULTS**

The elevation of CRP levels above normal indicates tissue damage, inflammation, or both with greater reliability. Regular monitoring of CRP levels is often used as a means of assessing disease activity and of guiding therapy. CRP determination is considered of greater practical significance than other indicators of inflammatory disease. Erythrocyte sedimentation rate (ESR) may become elevated as a result of non-inflammatory conditions. In these circumstances inflammatory disease may be excluded if CRP is absent.

**STABILITY OF THE REACTIONS**

Slide tests should be interpreted immediately after the 2-minute rotation period to avoid the possibility that a negative result may be incorrectly interpreted as positive due to drying of the reagent.

**LIMITATIONS**

- Reactions read beyond the two-minute interval may be invalid.
- The results obtained from this assay must be considered a part of the differential diagnosis and medical history of the patient.
- There is no relationship between the strength of reactivity and C-reactive protein levels.
- Hemoglobin ( $\leq 10$  g/L), bilirubin ( $\leq 20$  mg/dL) and lipemia ( $\leq 10$  g/L), do not interfere. Rheumatoid factors ( $\geq 100$  IU/mL), fibrin, Other substances may interfere.
- False positive or false negative results may also occur due to:
  - Contamination of test materials
  - Improper storage of test materials or omission of reagents
  - Deviation from the recommended techniques

**SPECIFIC PERFORMANCE CHARACTERISTICS**

- The kit has been characterized by all the procedures mentioned in the **Recommended Techniques**.
- Prior to release each lot of Lorne CRP Latex Test Kit is tested by **Recommended Techniques** to ensure suitable reactivity.
- The CRP latex sensitivity is calibrated to the Reference MaterialERM-DA-4721/FCC.
- Analytical sensitivity: 6 (5-10) mg/L, under the described assay conditions.
- Prozone effect: No prozone effect was detected up to 1600 mg/L (Note 1).
- Diagnostic sensitivity: 95.6 %.
- Diagnostic specificity: 96.2 %.

**DISCLAIMER**

- The user is responsible for the performance of the kit by any method other than those mentioned in the **Recommended Techniques**.
- Any deviations should be validated prior to use using established laboratory procedures.

**BIBLIOGRAPHY**

- Lars-Olof Hansson et al. Current Opinion in Infectious Diseases 1997; 10: 196-201
- M.M. Peay: The Lancet 1981; March 21; 653 - 656.
- Chetana Vaishnavi: Immunology and Infectious Diseases 1996; 6: 139 - 144.
- Yoshiyuki Yokama et al. Journal of Clinical Laboratory Status 1987; 1: 15 - 27.
- Yamamoto S. et al. Veterinary Immunology and Immunopathology 1993; 36: 257 - 264.
- Charles Watworth et al. Clinical Chimica Acta, 1984; 138: 309 - 316.
- Young DS. Effects of drugs on clinical laboratory test, 4th ed. AACCP Press, 1995.

**AVAILABLE KIT SIZES**

Kit Size	Catalogue Number
100 Tests Per Kit	850100A

For the availability of other sizes, please contact:

Lorne Laboratories Limited  
Unit 1 Cuboburn Park Industrial Estate  
Dunfermline  
Fife KY11 1YU  
Lower Earley  
Berkshire RG6 4JT  
England  
Tel: +44 (0) 118 921 2264  
Fax: +44 (0) 118 986 4518  
E-mail: [info@lornelabs.co.uk](mailto:info@lornelabs.co.uk)

- CRP Negative Control (Blue cap, 1 mL): Animal serum containing a preservative.
- Pipette stirrer.
- Disposable agglutination slide.

**MATERIALS AND EQUIPMENT REQUIRED**

- Serological Pipettes
- Mechanical rotator capable with adjustable speed of 80-100 rpm
- Vortex mixer.
- Small Glass or Plastic Test Tubes.
- Distilled or Deionised Water.
- 9 g/L saline solution.

**RECOMMENDED QUALITATIVE TECHNIQUE**

- Allow the reagents and samples to reach room temperature. The sensitivity of the test may be reduced at low temperatures.
- Place 50  $\mu$ L of the sample (Note 1) and one drop of each Positive and Negative controls into separate circles on the slide test.
- Swirl the CRP-latex reagent gently before using and add one drop (50  $\mu$ L) next to the samples to be tested.
- Mix the drops with a stirrer - spreading them over the entire surface of the circle. Use different stirrers for each sample.
- Place the slide on a rotary shaker at 80-100 r.p.m. for 2 minutes. False positive results could appear if the test is read after more than two minutes.

**INTERPRETATION OF QUALITATIVE RESULTS**

- Positive: Visible agglutination of latex particles constitutes a positive result and within the accepted limitations of the test procedure, indicates a level of CRP in the specimen > 6 mg/L.
- Negative: No visible agglutination of latex particles constitutes a negative result and within the accepted limitations of the test procedure, indicates a level of CRP in the specimen  $\leq$  6 mg/L.

**RECOMMENDED SEMI-QUANTITATIVE TECHNIQUE**

- The semi-quantitative test can be performed in the same way as the qualitative test using dilutions of the serum.
- Make doubling dilutions of the specimen using 9 g/L saline solution as follows:

Dilution	Serum	Saline
1/2	100 $\mu$ l undiluted serum	100 $\mu$ l
1/4	100 $\mu$ l 1/2 diluted serum	100 $\mu$ l
1/8	100 $\mu$ l 1/4 diluted serum	100 $\mu$ l
1/16	100 $\mu$ l 1/8 diluted serum	100 $\mu$ l

- Test the specimen dilutions in the same way as for the quantitative technique above.
- Agglutination of the sera indicates:

Dilution	CRP Levels (mg/l)
1/2	12 (6 x 2)
1/4	24 (6 x 4)
1/8	48 (6 x 8)
1/16	96 (6 x 16)

- Normal levels of CRP in adults are  $\leq$  6 mg/L.

**NOTES**

- High CRP concentration samples may give false negative results (pro-zone effect). Re-test the sample again using a drop of 20  $\mu$ l.

**RESULTS**

The titre is expressed as the reciprocal of the highest dilution showing macroscopic agglutination; e.g. if this occurs in dilution 1/8 the titre is 8.

Lorne Laboratories Limited  
Unit 1, Curbush Park Industrial Estate  
Danehill  
Lower Earley  
Berkshire RG6 4UT  
England  
Tel: +44 (0) 118 921 2264  
Fax: +44 (0) 118 966 4518  
E-mail: info@lornelabs.com

**RAPID LATEX KIT**  
**DIRECTIONS FOR USE**

**RF Latex kit: For Detection Of Rheumatoid Factor (RF).**

**SUMMARY**

Rheumatoid factors are a group of antibodies directed to determinants in the Fc portion of the IgG molecule. Although rheumatoid factors are found in a number of rheumatoid disorders, such as systemic lupus erythematosus (SLE) and Sjögren's syndrome, as well as in non-rheumatoid conditions, its central role lies in aiding in the diagnosis of rheumatoid arthritis.

**PRINCIPLE**

When used by recommended techniques, latex particles in reagent will agglutinate (clump) in presence of rheumatoid factor (RF). No agglutination generally indicates absence of RF (see Limitations).

**KIT DESCRIPTION**

Lorne RF Latex Kit is for the detection of rheumatoid factor. The latex reagent is a suspension of polystyrene latex particles coated with human gamma globulins, which agglutinate in the presence of Rheumatoid Factor (RF). All latex reagents are supplied at optimal dilution for use with all recommended techniques without the need for further dilution or addition. For lot reference number and expiry date see Vial Labels.

**STORAGE**

Do not freeze. Reagent vials should be stored at 2 - 8°C on receipt. Prolonged storage at temperatures outside this range may result in accelerated loss of reagent reactivity.

**SPECIMEN COLLECTION**

Specimens should be drawn without anticoagulant using an aseptic phlebotomy technique. If testing is delayed, fresh serum can be stored at 2-8°C for 7 days or for up to 3 months at or below -20°C. Specimens must be free from bacterial contamination, fibrin, gross lipaemia and gross haemolysis.

**PRECAUTIONS**

- The kit is for *in vitro* diagnostic use only.
- Do not use kit past expiration date (see Vial and Box Labels).
- Protective clothing should be worn when handling the reagents, such as disposable gloves and a laboratory coat.
- All the reagents must be allowed to reach 18-25°C before use.
- Materials used to produce the kit were tested at source and found to be negative for HIV 1+2 and HBsAg using approved microbiological tests. However, no known tests can guarantee that products derived from human or animal sources are free from infectious agents. Care must be taken in the use and disposal of each vial and its contents.

**DISPOSAL OF KIT REAGENT AND DEALING WITH SPILLAGES**

For information on disposal of kit reagent and decontamination of a spillage site see Material Safety Data Sheets, available on request.

**CONTROLS AND ADVICE**

- It is recommended the RF Positive and Negative Controls be tested in parallel with each batch of tests. Tests must be considered invalid if controls do not show expected results.
- All the reagents must be allowed to reach 18-25°C before use.
- Shake the reagents well before use to ensure homogeneity.
- Do not interchange components between different kits.
- Use of kit and interpretation of results must be carried out by properly trained and qualified personnel in accordance with the requirements of country where kit is in use. The user must determine the suitability of the kit for use in other techniques.
- Results obtained with a latex method do not compare with those obtained with the Rose Waaler test. Differences in the results between methods do not reflect differences in the ability to detect rheumatoid factors.

**RESULTS**

Slide is expressed as the reciprocal of the highest dilution showing macroscopic agglutination; e.g. if this occurs in dilution 1/8, the titre is (8 x 8 IU/mL) = 64 IU/mL.

**STABILITY OF THE REACTIONS**

Slide tests should be interpreted straight after the 2-minute rotation period to avoid the possibility that a negative result may be incorrectly interpreted as positive due to drying of the reagent.

**LIMITATIONS**

- Using a latex test system, positive results are not always found with every case of clinically-defined rheumatoid arthritis; the number of positives reported using various types of latex reagent range from 70% to over 90%.
- The incidence of false positive results is about 3.5%. Individuals suffering from infectious mononucleosis, hepatitis, syphilis as well as elderly people may give positive results.
- Diagnosis should not be solely based on the results of latex method but also should be complemented with a Rose Waaler test along with the clinical examination.
- Haemoglobin ( $\leq 10$  g/L), bilirubin ( $\leq 20$  mg/dL) and lipaemia ( $\leq 10$  g/L) do not interfere. Other substances may interfere:
- False positive or false negative results may also occur due to:
  - Contamination of test materials
  - Inappropriate storage of test materials or omission of reagents
  - Deviation from the recommended techniques

**SPECIFIC PERFORMANCE CHARACTERISTICS**

- The kit has been characterized by all the procedures mentioned in the Recommended Techniques.
  - The RF latex sensitivity is calibrated against the RF International Calibrator from the WHO (WHO 64/2 Rheumatoid Arthritis Serum).
  - Analytical sensitivity: 8 (8-16) IU/mL, under the described assay conditions.
  - Prozone effect: No prozone effect was detected up to 1500 IU/mL.
  - Diagnostic sensitivity: 100%.
  - Diagnostic specificity: 100%.
- The diagnostic sensitivity and specificity have been obtained using 139 samples compared with the same method of a comparator.

**DISCLAIMER**

- The user is responsible for the performance of the kit by any method other than those mentioned in the Recommended Techniques.
- Any deviations should be validated prior to use using established laboratory procedures.

**BIBLIOGRAPHY**

- Robert W Dorner et al. Clinica Chimica Acta 1987; 167: 1 - 21
- Frederick Wolfe et al. Arthritis and Rheumatism 1991; 34: 951-960
- Robert H Shmerling et al. The American Journal of Medicine 1991; 91: 528 - 534
- Adalbert F. Schubart et al. The New England Journal of Medicine 1989; 261: 363 - 368
- Charles M. Ploz 1956; American Journal of Medicine; 21:833 - 836
- Young DS. Effects of drugs on clinical laboratory test, 4th ed. AACF Press, 1995.

**AVAILABLE KIT SIZES**

Kit Size	Catalogue Number
100 Tests Per Kit	830700A

**KIT COMPONENTS SUPPLIED**

- RF Latex Reagent (5 mL); Latex particles coated with human  $\gamma$ -globulin, pH 8.2, and a preservative.
- RF Positive Control (Red cap, 1 mL); Human serum with a RF concentration > 30 IU/mL and a preservative.
- RF Negative Control (Blue cap, 1 mL); Animal serum and a preservative.
- Pipette-Stirrers.
- Reusable Agglutination Slide (18 each).

**MATERIALS AND EQUIPMENT REQUIRED**

- Glass Test Tubes (10 x 75 mm or 12 x 75 mm).
- Pasteur and Graduated Pipettes.
- Vortex mixer.
- Mechanical rotator with adjustable speed of 80-100 rpm.

**RECOMMENDED QUALITATIVE TECHNIQUE**

- Allow the reagents and samples to reach room temperature. The sensitivity of the test may be reduced at low temperatures. Place 50  $\mu$ L of the sample and one drop of each Positive and Negative controls into separate circles on the slide test.
- Shake the RF-latex reagent vigorously or in a vortex mixer before use and add one drop (50  $\mu$ L) next to the sample to be tested.
- Mix the drops with a stirrer, spreading them over the entire surface of the circle. Use different stirrers for each sample. Place the slide on a rotary shaker at 80-100 rpm, for 2 minutes. False positive results could appear if the test is read after more than two minutes.

**INTERPRETATION OF QUALITATIVE RESULTS**

- Positive:** Visible agglutination of latex particles constitutes a positive result and within the accepted limitations of the test procedure, indicates a level of RF in the specimen > 8 IU/mL.
- Negative:** No visible agglutination of latex particles in a milky liquid constitutes negative result and within accepted limitations of test procedure, indicates level of < 8 IU/mL RF in specimen.

**RECOMMENDED SEMI-QUANTITATIVE TECHNIQUE**

- The semi-quantitative test can be performed in the same way as the qualitative test using dilutions of the serum.
- Make doubling dilutions of serum specimen in 8 g/L saline as follows:

Dilution	Serum	Saline
1/2	100 $\mu$ L undiluted serum	100 $\mu$ L
1/4	100 $\mu$ L 1/2 diluted serum	100 $\mu$ L
1/8	100 $\mu$ L 1/4 diluted serum	100 $\mu$ L
1/16	100 $\mu$ L 1/8 diluted serum	100 $\mu$ L

- Test the specimen dilutions in the same way as for the quantitative technique above.
- Agglutination of the sera indicates:

Dilution	RF Levels (IU/mL)
1/2	16 (8 x 2)
1/4	32 (8 x 4)
1/8	64 (8 x 8)
1/16	128 (8 x 16)

- Normal levels of RF in adults is < 8 IU/mL.

# Declaration of Conformity

HL-7-0692DC DOI 2015/08 (1)

## In Application of the Council Directive 98/79/EC on the approximation of the laws of the Member States relating to *In Vitro* Diagnostic Medical Devices & CE marking.

Declaration of conformance to applicable sections of Annex I - Essential Requirements and Annex III (EC Declaration of Conformity) imposed by sections 2 to 5. The below listed products are not classified under Annex II Lists A or B. Access to the appropriate technical files will be made available to the appropriate body in the event this is required.

Product Code	Description	GMDN Classification Code
5250H	Manual D-Dimer	47346

I, the undersigned declare that the devices registered against the above GMDN Classification Code conforms to the said Directives.

Full Name: M.J. Stephenson

Title: Managing Director

Signed:

A handwritten signature in black ink, appearing to read "Michael / Stephenson".

Date: 25 Aug 2015

Tel +44 (0)191 482 8440  
Fax +44 (0)191 482 8442  
info@helena-biosciences.com  
www.helena-biosciences.com

Helena Biosciences Europe  
Queensway South, Team Valley Trading Estate,  
Gateshead, Tyne and Wear, NE11 0SD,  
United Kingdom

# Declaration of Conformity

HL-7-0664DC DOI 2015/08 (1)

## In Application of the Council Directive 98/79/EC on the approximation of the laws of the Member States relating to *In Vitro* Diagnostic Medical Devices & CE marking.

Declaration of conformance to applicable sections of Annex I - Essential Requirements and Annex III (EC Declaration of Conformity) imposed by sections 2 to 5. The below listed products are not classified under Annex II Lists A or B. Access to the appropriate technical files will be made available to the appropriate body in the event this is required.

Product Code	Description	GMDN Classification Code
5267L	Thromboplastin L	55983

I, the undersigned declare that the devices registered against the above GMDN Classification Code conforms to the said Directives.

Full Name: M.J. Stephenson

Title: Managing Director

Signed:

A handwritten signature in black ink, appearing to read "Michael / Stephenson".

Date: 06 Aug 2015

Tel +44 (0)191 482 8440  
Fax +44 (0)191 482 8442  
info@helena-biosciences.com  
www.helena-biosciences.com

Helena Biosciences Europe  
Queensway South, Team Valley Trading Estate,  
Gateshead, Tyne and Wear, NE11 0SD,  
United Kingdom

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# Declaration of Conformity

HL-7-0511DC DOI 2015/08 (4)

In Application of the Council Directive 98/79/EC on the approximation of the laws of the Member States relating to In Vitro Diagnostic Medical Devices & CE marking.

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Title: Managing Director

Signed:

A handwritten signature in black ink, appearing to read "Michael Stephenson".

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# Clauss Fibrinogen 100

REF 5376  
REF 5376H

Helena Biosciences Europe, Queensway South, Team Valley Trading Estate, Gateshead, Tyne and Wear, NE11 6SD, United Kingdom  
Tel: +44 (0)191 482 8440  
Fax: +44 (0)191 482 8442  
Email: info@helena-biosciences.com  
Web: www.helena-biosciences.com

HL-2-0582P 2016/01 (11)

## Clauss Fibrinogen 100 Instructions for Use

### INTENDED PURPOSE

The Clauss Fibrinogen 100 kit is intended for carrying out clot based immunoturbidimetric assays.

Clauss' method is a simple method for the quantitative determination of fibrinogen by measuring the clotting time of whole plasma after the addition of thrombin (not kit included). This clot time is proportional to the fibrinogen concentration. Levels of fibrinogen can increase as a result of inflammation, pregnancy or oral contraceptive use. Decreased levels may be found in certain states such as liver disease and DIC. Congenital deficiencies include afibrinogenemia, hypofibrinogenemia, dysfibrinogenemia, hypodysfibrinogenemia, and dysplasminogenemia (genomic fibrinogen deficiency). Clauss Fibrinogen 100 is used for the quantitative determination of fibrinogen in human plasma.

### WARNINGS AND PRECAUTIONS

The reagents contained in this kit are for in vitro diagnostic use only - DO NOT INGEST. Wear appropriate personal protective equipment when handling kit components. Refer to the product safety declaration for the kit for appropriate safety and precautionary statements where applicable. Dispose of components in accordance with local regulations.

Blood products have been screened and found negative (unless otherwise stated on the kit box or vial) for the presence of Hepatitis B (HBsAg), Hepatitis C (Anti-HCV) and HIV-1 antibody (anti-HIV-1).

Use only vials that should be handled with the same precautions as a human patient sample.

### COMPOSITION

Component	Content	Description	Preparation
Thrombin	1 x 1 mL (REF 5376)	Contains approximately 100 IU (REF 5376) units of Thrombin activity.	Reconstitute each vial with 2 mL (REF 5376H) of distilled water or heparin suspension depending on manual or automated use. Discard any unused portion.
Fibrinogen Calibrator	2 x 1 mL (REF 5376) 1 x 1 mL (REF 5376H)	Contains 1 mL of lyophilized human fibrinogen standard.	Reconstitute each vial with exactly 1 mL of distilled water. Swirl gently and allow to stand for 10 minutes. Use gently before use. Do not shake.
Donor's Buffer	2 x 25 mL (REF 5376) 1 x 25 mL (REF 5376H)	Contains 25 mL of buffer which contains heparin and sodium chloride and sodium azide (3% preservative).	The buffer is ready for use as packaged.
Kalzin Suspension	2 x 5 mL (REF 5376) 1 x 5 mL (REF 5376H)	Contains 5 mL of kalzin suspension in 0.5 g of gelatin.	Mix vigorously immediately before use. Ensure the sedimentation is visible.

Each kit contains instructions for use. Each kit contains kit specific reference values chart.

### ITEMS REQUIRED BUT NOT PROVIDED

Helena Biosciences Europe provides reagents are suitable for use of automatic Clauss Fibrinogen analysers. Refer to the instrument Operator Manual and application note for appropriate instructions.

The reagent is supplied in 10 x 10 mL vials available as separate units:

- REF 5374 Clauss Fibrinogen (Thrombin only)
- REF 5379 Clauss Fibrinogen (Thrombin only)
- REF 5370 Fibrinogen Calibrator
- REF 5375 Donor's Buffer
- REF 5376 Kalzin Suspension

### STORAGE, SHELF-LIFE AND STABILITY

Unopened vials are stable until the expiry date when stored under conditions indicated on the vial or kit label.

Thrombin: Once reconstituted, the reagent is stable for 2 hours at 15-20°C, 1 week at 2-8°C or 1 month at -20°C.

Fibrinogen Calibrator: Once reconstituted, the reagent is stable for 4 hours at 2-8°C.

Donor's Buffer: Stable at 2-8°C once opened.

Kalzin Suspension: Stable at 2-8°C once opened.

### SAMPLE COLLECTION AND PREPARATION

Plasma is obtained from blood which is used throughout. Blood is centrifuged into a collected into 3.2% or 3.8% sodium citrate anticoagulant (1 part). Samples should be centrifuged at 1500 g for 15 minutes. Plasma should be held at 2-8°C or 15-20°C. Testing should be completed within 4 hours of sample collection, or plasma can be stored frozen at -20°C for 2 weeks or 10°C for 6 months. The assay is at 37°C and to testing. Do not keep at 37°C for more than 5 minutes.

### PROCEDURE

Manual Method  
Prepare all reagents as instructed under "Composition".

1. Standard Curve Preparation  
A three point curve should be used when reagent lot numbers change or if expected values lie outside of control testing range. Prepare the following standards in duplicate (use without shaking).

No.	Dilution	Fibrinogen Calibrator (mL)	Buffer (mL)
1	1 + 4	0.2	0.8
2	1 + 9	0.1	0.9
3	1 + 19	0.1	1.9
4	1 + 29	0.1	2.9
5	1 + 39	0.1	3.9

2. Patient Sample Preparation  
Prepare 1 µl dilutions of the patient plasma or control plasma in Donor's Buffer:

- a. Mix without shaking.

3. Testing  
a. Perform all tests in duplicate.  
b. Mix 0.2 mL of standard, patient or control dilution into a reaction tube and incubate at 37°C for 2 minutes.  
c. Add 0.1 mL of thrombin (15-20°C).  
d. Determine the clot time in the nearest 0.1 seconds.  
e. Plot mean Standard Deviation (mean) versus Fibrinogen Level (Y-axis) on the Fibrinogen graph paper. A straight line should be obtained. Assign the Fibrinogen Calibrator reference values to the 1-5 dilution to allow direct interpretation of patient and control values from the standard curve.

### Automated Method

Refer to the appropriate instrument Operator Manual for detailed instructions on automatic specific application guides.

### INTERPRETATION OF RESULTS

Normal values for fibrinogen in healthy adults are 150-350 mg/dL (1.5-3.5 g/L).

### LIMITATIONS

Healthy tissue cells and hemolytic degradation products >100 mg/dL may cause falsely low fibrinogen quantitation. Fibrinogen values below the standard curve values for the patient samples may occur using an appropriate dilution to bring values into the standard range.

### QUALITY CONTROL

Each laboratory should establish a quality control program. Normal and abnormal control plasmas should be tested prior to each run of patient testing. In-house satisfactory management and control programme. If controls do not perform as expected, patient results should be considered invalid. Helena Biosciences Europe supply the following control materials for use with this product:

- REF 5391 Specialty Assay Control A
- REF 5392 Specialty Assay Control A
- REF 5386 Routine Control B
- REF 5387 Routine Control A
- REF 5383 Routine Control SA

### REFERENCE VALUES

Reference values can vary between laboratories depending on the techniques and systems in use. For this reason each laboratory should establish its own reference values.

### PERFORMANCE CHARACTERISTICS

Helena Biosciences Europe in-house investigations have determined the following performance characteristics as a guideline. Each laboratory should establish its own performance data. The manual Clauss Fibrinogen is designed to give a linear calibration from 1.5 - 6.5 g/L.

Reproducibility intra-assay precision			Inter-assay precision		
Fibrinogen (g/L)	n	CV (%)	Fibrinogen (g/L)	n	CV (%)
1.21	3	4.88	1.02	10	6.17
3.01	5	3.58	3.02	10	3.75

### BIBLIOGRAPHY

1. Clauss A (1957) Gerinnungsphysiologische Schnellmethode zur Bestimmung des Fibrinogens. Acta Haematol 17:203-246.
2. Shaw TS (1977) Assays for Fibrinogen and its Derivatives. CRC, Clin Chem Lab Sci 14:1-192.
3. Clinical and Laboratory Standards Institute (2006) Collection, Transport and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Haematology Assays. Approved Document, 5th edn. CLSI H21-A5.
4. Study HE et al (1980) Normal Reference Laboratory Values. New England Journal of Medicine, 302(7):37-48.
5. Dineen T, Shanks V (1978) Plasma Fibrinogen Determination by Automated Immuno-Turbidimetry. American Journal of Medical Technology, 38(8):196-201.

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Title: Managing Director

Signed:

A handwritten signature in black ink, appearing to read "Michael / Stephenson", written over a horizontal line.

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# Declaration of Conformity

HL-7-0137DC DOI 2015/07 (7)

## **In Application of the Council Directive 98/79/EC on the approximation of the laws of the Member States relating to *In Vitro* Diagnostic Medical Devices & CE marking.**

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<b>Product Code</b>	<b>Description</b>	<b>GMDN Classification Code</b>
5186	Routine Control N	30590

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Product Code	Description	GMDN Classification Code
5187	Routine Control A	30590

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Title: Managing Director

Signed:

A handwritten signature in black ink, appearing to read "M.J. Stephenson".

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# Coagulation Control Plasmas

REF 5186 Routine Control N  
REF 5187 Routine Control A  
REF 5183 Routine Control SA  
REF 5482 Routine Coagulation Control Set

Helena Biosciences Europe, Queensway South, Team Valley Trading Estate, Gateshead, Tyne and Wear, NE11 0SD, United Kingdom  
Tel: +44 (0)191 482 8440  
Fax: +44 (0)191 482 8442  
Email: [info@helena-biosciences.com](mailto:info@helena-biosciences.com)  
Web: [www.helena-biosciences.com](http://www.helena-biosciences.com)

HL-2-0482P 2016/01 (16)

## Coagulation Control Plasmas

### INTENDED PURPOSE

The Intended Control Plasmas kit is intended for use as a quality control material.

Routine Control N, Routine Control A and Routine Control SA are for use as normal, moderately prolonged and markedly prolonged controls for PT and aPTT assays. They are also assayed for Fibrinogen, TCT and ATIII, and are prepared from normal human plasma.

### WARNINGS AND PRECAUTIONS

The reagents contained in this kit are for *in vitro* diagnostic use only – DO NOT INGEST. Wear appropriate personal protective equipment when handling all kit components. Refer to the product safety declaration for the link to appropriate hazard and precautionary statements where applicable. Dispose of components in accordance with local regulations.

Blood products have been screened and found negative (unless otherwise stated on the kit box or vial) for the presence of: Hepatitis B Antigen (HbsAg), HIV-1 antibody, HIV-2 antibody, HCV antibody.

However they should not be handled with the same precautions as a human patient sample.

### COMPOSITION

REF	Component	Content	Description
5186	Routine Control - N	10 x 1 mL	Prepared from pooled normal plasma.
5187	Routine Control - A	10 x 1 mL	Prepared from defibrinated human plasma.
5183	Routine Control - SA	10 x 1 mL	Prepared from adsorbed human plasma.

REF	Component	Content	Description
5482	Routine Coagulation Control Set	4 x 1 mL Routine Control - N Routine Control - A Routine Control - SA	4 x 1 mL 3 x 1 mL 3 x 1 mL

Each kit contains instructions for use.

Each vial contains 1 mL of buffered lyophilised human plasma. Each vial contains 1 mL of buffered lyophilised human plasma (or 10 minutes for complete coagulation) and mix well before use.

### ITEMS REQUIRED BUT NOT PROVIDED

Coagulation Control Plasmas may be used when performing tests on any mechanical or photo-optical coagulation instrument in conjunction with all suitable commercial reagents.

### STORAGE, SHELF-LIFE AND STABILITY

Unopened vials are stable until the expiry date when stored under conditions indicated on the vial or kit label. The reconstituted controls are stable for 6 hours when kept at 2 – 7°C or 4 weeks at -20°C when fresh frozen, deep stored. SAME COLLECTION AND PREPARATION: Not applicable.

### PROCEDURE

Each control should be used in the same routine as the unknown specimen to be rechecked with the reagent/ reagents of your preference/procedure.

### INTERPRETATION OF RESULTS

Routine Control N should give values within the laboratory normal range for PT, aPTT and Dimerogen assay. Routine Control A and Routine Control SA have been standardised to give prolonged and markedly prolonged PT and aPTT time, respectively. Use and interpretation of specific control values are provided with each pack of controls.

### LIMITATIONS

The results obtained with Coagulation Control Plasmas depend on several factors strongly associated with instrumentation, types of reagents, selected substrates and laboratory to laboratory variations. Each laboratory should establish an expanded range for the particular instrument/reagent system.

### QUALITY CONTROL

Each laboratory should establish a quality control program, Normal and abnormal control plasmas should be tested prior to each batch of patient samples, to ensure satisfactory instrument and operator performance. If limits do not perform as expected, further results should be considered invalid.

### REFERENCE VALUES

Reference values vary very broadly, laboratories depending on the methodology and systems in use. For this reason each laboratory should establish its own reference ranges.

### PERFORMANCE CHARACTERISTICS

The following performance characteristics have been determined by Helena Biosciences Europe or their representatives using an international coagulation instrument. Each laboratory should establish its own performance data.

Sample	n	Intra-assay precision aPTT CV (%)	PT CV (%)
Routine Control N	5	2.83	1.01
Routine Control A	5	2.76	1.71
Routine Control SA	5	1.72	1.03

### BIBLIOGRAPHY

1. Kirkwood TBL, et al. (1977) Identification of Sources of Variation in Factor VIII Assay. *British Journal of Haematology*, 37:553-558
2. Godebier M et al. (1971) Reproducibility in Coagulation Assays. *Acta Medica Scandinavica*, 80:231-235
3. Palfrey HA and Longbery JR (1973) A French Study of Coagulation Factor Assay Techniques. *Acta Medica Scandinavica*, 80:231-235.

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REF	Composant	Contenu	Description
5482	Routine Coagulation Control Set	4 x 1 mL Routine Control - N Routine Control - A Routine Control - SA	4 x 1 mL 3 x 1 mL 3 x 1 mL

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TÜVRheinland

**EC Certificate**

Directive 98/79/EC Annex IV, excluding Sections 4 and 6  
Full Quality Assurance System  
In Vitro Diagnostic Medical Devices

Registration No.: HL 60119814 0001

Report No.: 21265422 001

Manufacturer: Macherer-Nagel GmbH & Co. KG  
Neumann-Neander-Str. 6-8  
52355 Düren  
Deutschland

Products: Products for self-testing  
(see attachment for products and sites included)

Replaces Certificate, Registration No.: HL 60076687 0001

Expiry Date: 2022-05-28

The Notified Body hereby declares that the requirements of Annex IV, excluding section 4 and 6 of the directive 98/79/EC have been met for the listed products. The above named manufacturer has established and applies a quality assurance system, which is subject to periodic surveillance, defined by Annex IV, section 5 of the aforementioned directive. For placing on the market of List A devices covered by this certificate an EC design-examination certificate according to Annex IV, section 4 and a verification of manufactured products according to section 6 is required.

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TÜV Rheinland LGA Products GmbH is a Notified Body according to Directive 98/79/EC concerning in vitro diagnostic medical devices with the identification number 0197.

TÜVRheinland

Doc. 1/1. Rev. 0

TÜV Rheinland  
LGA Products GmbH

Tillystraße 2, 90431 Nürnberg

Attachment to  
Certificate

Registration No.: HL 60119814 0001

Report No.: 21265422 001

Manufacturer: Macherer-Nagel GmbH & Co. KG  
Neumann-Neander-Str. 6-8  
52355 Düren  
Deutschland

Products for self-testing:  
- Single and multi-parameter disposable test strips  
for urine analysis  
- Indicator test strips and papers for measurement  
of pH in urine

Additional site for warehousing and logistics:

Rahmetr. 120  
52355 Düren, Germany

Notified Body

Date: 2017-05-29

Dipl.-Ing. Sven Hoffmann

Negative and positive control solutions for the Medi-Test urine test strips

Application:

Medi-Test Control are control solutions to verify the correct functioning of the Medi-Test Urine Test Strips and the URYXXON® Urine Analysis Systems. The measured values are compared with the target values shown in the table below. Only qualified personnel must be allowed to carry out the application.

Reagents:

Each pack contains:

- 1 test tube with 15 mL reagent solution Medi-Test Control N
- 1 test tube with 15 mL reagent solution Medi-Test Control P

These are ready-to-use control solutions of purely chemical composition. The solution N simulates a urine sample with a value in the negative or normal range. The solution P will produce a positive color reaction with the Medi-Test Urine Test Strips for the parameters blood, urobilinogen, bilirubin, protein, nitrite, ketones, glucose and leukocytes so that values in the pathological range are indicated.

Safety precautions:

Reagent solutions Medi-Test Control contain Chloromethylisothiazolinone 0.0015–0.06% CAS 26172-55-4. WARNING H317 May cause an allergic skin reaction. P280 Wear protective gloves/eye protection. P302+352 IF ON SKIN: wash with plenty of water. P333+313 If skin irritation or rash occurs: get medical advice/attention. For further information, please ask for safety data sheet (see [www.mn-net.com/MSDS](http://www.mn-net.com/MSDS)).

WARNING

Storage and Shelf Life:

When not in use, it is recommended that the reagent solutions Medi-Test Control N and P are stored in a dark place at 2–8 °C. Do not freeze the solutions under any circumstances.

If correctly stored, the reagent solutions may be used right up until the use-by-date printed on the packaging. After the first use, each reagent solution may be used for up to three months or for dipping the test strip up to 20 times, whichever occurs first.

Precipitations in the reagent solutions have no effect on the measurement result. If the solutions show other contamination, then they must no longer be used. The solutions may be disposed of with plenty of running water via the sewerage of the local wastewater treatment plant.

Instructions for use:

Take the test tubes with the solutions Medi-Test Control P and N out of the refrigerator, and allow them to warm up to room temperature, shake well so that the solutions can mix homogeneously. Avoid frothing.

Open the lids of the test tubes. Do not pour the solutions into another vessel. Perform measurements directly in the test tube. Where required, take out a test strip from a test strip container. After taking out a strip, immediately close the container again. Do not touch the test fields with your fingers.

Dip the test strip with all the test fields into the respective reagent solution for approx. 1 second.

Expected values:

The ranges identified in the table were determined using several different batches of the Medi-Test Combi 10® SGL and the URYXXON® Stick 10, and were ascertained using various different URYXXON® devices. Each laboratory should use the results provided only as a reference and establish own parameters of precision.

After removing the test strip from the solution, briefly pat the edge along one side on an absorbent piece of material (e.g. kitchen roll) to dry it off. Insert the test strip into the device according to the instructions of the reflection photometers URYXXON® 300/500 and URYXXON® Relax. The test fields are analysed, and then the results are printed out.

When visually assessing the reaction colors of the test fields, compare the test fields with the color scale after about 30–60 seconds (leukocytes test field after 60–120 seconds). The most opportune reading time is after 30 seconds. Any color changes that occur after more than 2 minutes are of no significance.

Any visual assessment of the urine test strips must be made in daylight. Direct sunlight, however, must be avoided.

The urobilinogen test field may display a slightly more orange-red coloring in comparison to the color scale.

When not in use, and after completing the control procedures, replace the lids on the test tubes of the reagent solutions, and store them at 2–8 °C.

Note:

Please also note the instructions for the URYXXON® reflection photometers. Do not ingest the solutions! Avoid any contact with skin or eyes! Store reagent solutions in a safe place inaccessible to children!

Rev. 10/2018

Explanation of symbols

Statement of Conformity (Product corresponds to the In-Vitro Diagnostic Medical Devices Directive 98/79/EC of the European Union)

Please read instructions for use!

Permitted storage temperature range

Use by

Batch identification

Item number

Manufacturer

In vitro diagnostics product

Control

Control, negative

Control, positive

Analytes	Medi-Test Combi 10® SGL visually		URYXXON® Stick 10 with URYXXON® 300 / URYXXON® 500 / URYXXON® Relax	
	Control N	Control P	Control N	Control P
Blood	Negative	10–250 Ery/µL	Negative	10–250 Ery/µL
Urobilinogen	Normal	2–12 mg/dL <sup>1)</sup> (35–200 µmol/L)	Normal	2–12 mg/dL <sup>1)</sup> (35–200 µmol/L)
Bilirubin	Negative	1–4 mg/dL (1+–3+)	Negative	1–6 mg/dL (1+–3+)
Protein	Negative	100–500 mg/dL	Negative	25–500 mg/dL
Nitrite	Negative	Positive	Negative	Positive
Ketones	Negative	25–300 mg/dL (1+–3+) (2.5–30 mmol/L)	Negative	5–300 mg/dL (1+–3+) (0.5–30 mmol/L)
Glucose	Negative–Normal	500–≥1000 mg/dL (27.8–55.5 mmol/L)	Negative–Normal	50–1000 mg/dL (2.8–60 mmol/L)
pH value	5–7	7–9	5–7	7–9
Specific gravity (density)	1.010–1.030	1.005–1.025	1.010–1.030	1.005–1.025
Leukocytes	Negative	25–500 Leu/µL	Negative	15–500 Leu/µL

<sup>1)</sup> The urobilinogen test field displays an orange-red coloring compared to the color scale

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GOST R CERTIFICATION SYSTEM  
FEDERAL AGENCY FOR TECHNIQUE REGULATION AND  
METROLOGY

VOLUNTARY CERTIFICATION SYSTEM  
"SMK-STANDART"

Reg. No. POCC RU.31060.04ЖЖЮ0

Certification authority:

REG No. SMK STANDART.RU.0005

INTERNATIONAL CERTIFICATION CENTER Limited Liability Company

Address: 138, Naberezhnaya Obvodnogo Canala, block 1, office 421, St. Petersburg, 190020  
phone: +7 (812) 438-76-71 standart@iso-smk.ru

Check the authenticity of the certificate in the register on the website <http://www.iso-smk.ru>

### CERTIFICATE OF CONFORMITY

No. ST.RU.0001.M0013380

This Certificate of Conformity is issued to

**Agat-Med, Ltd.**

Address: 6, Glavnaya st., Moscow, 105173, Russia  
TRN 7719187311 OGRN 1037739078970

Date of issue: 26.01.2018

Period of validity: 26.01.2021

**This certificate certifies that:**

*Medical devices. Quality management system. System requirements for regulatory purposes in relation to the works in accordance with Annex 1 to this certificate*

(the attachment is an integral part of the certificate)

**CORRESPONDS TO THE REQUIREMENTS OF GOST ISO 13485-2011 (EN ISO 13485:2003)**

Manager of "Expert" authority

V.V. Koptsev

O.V. Gundareva

THIS CERTIFICATE BINDS THE ORGANIZATION TO MAINTAIN THE WORKS PERFORMED ACCORDING TO THE STANDARD MENTIONED ABOVE, TO BE CONTROLLED BY THE INSTITUTIONAL BODY OF VOLUNTARY CERTIFICATION SYSTEM "SMK STANDARD" AND IS CONFIRMED DURING THE ANNUAL INSPECTION CONTROL

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Este certificado ampara todas las marcas de estos productos incluidas por el fabricante en su declaración de conformidad. / This certificate covers all trademarks of these products included by the manufacturer in his declaration of conformity.