New as of:
Nuevo a partir de:
Novo a partir de:
Действует с:
新版自：

10.2020

inLab MC X5, CEREC Primemill, CEREC MC / MC X / MC XL Premium Package

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
10.2020

inLab MC X5, CEREC Primemill, CEREC MC / MC X / MC XL Premium Package

Operating Instructions for suction devices

English


Table of contents

1	General data	4
1.1	Dear Customer,.....	4
1.2	Contact data.....	4
1.3	General conventions and structure of the document	5
1.3.1	Structure of the document.....	5
1.3.1.1	Identification of danger levels	5
1.3.1.2	Formats and symbols used	5
1.3.2	Notes to the repository	5
1.4	Warranty and liability.....	6
1.5	Intended use	6
1.6	Legend	7
2	Safety information	8
2.1	Basic safety information.....	8
2.1.1	Prerequisites	8
2.1.2	Operating	8
2.1.3	Maintenance and repair	8
2.1.4	Modifications to the product.....	9
2.1.5	Accessories.....	9
2.1.6	In case of damage	9
2.2	Connecting the unit	9
2.3	Ventilation slots.....	9
3	General information on the device	10
3.1	Standards/ approvals/ certifications	10
3.2	Scope of supply	10
3.3	Technical description (component and interfaces)	12
3.4	Technical data.....	13
3.4.1	InLab MC X5 suction device	13
3.4.2	CEREC MC/MC X/MC XL Premium Package/Primemill suction device.....	14
4	Transportation and installation.....	15
4.1	Transport and unpacking	15
4.2	Disposal of packaging materials	15
4.3	Installation conditions.....	15
4.3.1	Installation site requirements	15

4.4	Connecting the suction device	16
4.4.1	Connecting the suction tube	16
4.4.2	Connecting the power cord	17
4.4.3	Connecting the interface cable (for automatic mode)	17
5	Commissioning	18
5.1	Putting the unit into operation	18
5.1.1	Automatic mode	18
5.1.2	Manual operation	18
5.1.3	Setting the suction power	18
5.2	Repacking	18
6	Maintenance and repair	19
6.1	Changing filter bags and HEPA filters	19
6.2	Changing the suction motor	21
6.3	Using the cleaning hose and the wet cleaning process	21
7	Malfunctions	23
7.1	Error messages	23
7.2	Troubleshooting and repair	24
8	Disposal	25
9	Spare parts and consumables	26

1 General data

1.1 Dear Customer,

Thank you for purchasing this suction device from Dentsply Sirona.

The suction devices were specially developed for our products inLab MC X5 and CEREC MC / MC X / MC XL Premium Package / Primemill and offer you maximum safety and the best working conditions when using a Dentsply Sirona system.

Improper handling and unintended use may lead to hazards and damage. This is why we ask that you read through the operating manual and adhere to it closely. Also pay attention to the safety instructions to prevent personal injury and material damage.

Please read this document completely and follow the instructions exactly. You should always keep it within reach. The original language of this document is German.

Your
inLab and CEREC Team

1.2 Contact data

Customer Service Center

In the event of technical queries, please use our online contact form at the following address:
<http://srvcontact.sirona.com>

Manufacturer's address

Sirona Dental Systems GmbH
Fabrikstrasse 31
64625 Bensheim
Germany

Tel.: +49 (0) 6251/16-0
Fax: +49 (0) 6251/16-2591
e-Mail: contact@dentsplysirona.com
www.dentsplysirona.com

1.3 General conventions and structure of the document

1.3.1 Structure of the document

1.3.1.1 Identification of danger levels

To prevent personal injury and material damage, please observe the warning and safety information provided in these operating instructions. Such information is highlighted as follows:

CAUTION

A possibly dangerous situation that could result in slight bodily injury.

NOTE

A possibly harmful situation which could lead to damage of the product or an object in its environment.

IMPORTANT

Application instructions and other important information.

Tip: Information on making work easier.

1.3.1.2 Formats and symbols used

The formats and symbols used in this document have the following meaning:

✓ Prerequisite 1. First action step 2. Second action step - or ➤ Alternative action ↪ Result ➤ Individual action step	Requests you to do something.
See "Formats and symbols used [→ 5]"	Identifies a reference to another text passage and specifies its page number.
• List	Designates a list.
"Command / menu item"	Indicates commands / menu items or quotations.

1.3.2 Notes to the repository

It is mandatory to keep this operating manual in an easily accessible place for the purpose of later reference. In the event of a sale or transfer of the device to another user, make sure that the device is supplied along with the operating manual, so that the new owner can get acquainted with the operation and the appropriate precautions and warnings

1.4 Warranty and liability

Maintenance

In the interest of the safety and health of users or third parties, it is necessary that maintenance work is carried out at fixed time intervals to ensure the operational safety and reliability of your product.

The operator must ensure the implementation of the maintenance work.

As a manufacturer of electro-medical equipment, we can consider ourselves responsible for the safety characteristics of the device only if maintenance and repairs are carried out only by us or by companies authorized explicitly by us for this purpose and if components are replaced with original spare parts in case of failure.

Exclusion of liability

If the operator does not meet the obligation to carry out such maintenance or fault messages are ignored, Sirona Dental Systems or its authorized dealer does not assume any liability for damage caused.

1.5 Intended use

The suction system is used to extract dry, non-explosive dusts and chips that emerge when processing various materials for the production of dental aids and technologies and dentures. The suction system may only be operated indoors. The suction system must not be used for any other purpose.

This unit must not be used for any other purpose. If the unit is used for any purpose other than the one mentioned above, it may be damaged.

Intended use also includes compliance with these Operating Instructions and the relevant maintenance instructions.

CAUTION

Follow the instructions

If the instructions for operating the unit described in this document are not observed, the intended protection of the user may be impaired.

CAUTION

Dry processing


In the event of dry processing without suction, dust can be created from the materials being processed. Breathing in this dust can have associated health risks.

- > As such, observe the information and requirements of the suction.
- > Dry processing is only permissible in conjunction with the CEREC suction device 230 V or 120 V.
 - CEREC suction device 230 V / 120 V, ordered together with the device: REF 6569730
 - CEREC suction device 230 V/120 V, if ordered separately: REF 6580786

1.6 Legend

Year of manufacture

Product disposal symbol (see "Disposal [-> 25]").

Accompanying documents

This symbol can be found on the rating plate on the unit.

Meaning: Observe the Operating Instructions when operating the unit.

Symbols on the packaging

Take note of the following symbols on the packaging:

Top


Protect from moisture


Fragile; handle with care


Stack limit

Temperature during storage and transport

Relative humidity during storage and transport

Air pressure during storage and transport

2 Safety information

2.1 Basic safety information

2.1.1 Prerequisites

NOTE

Important information on building installation

The building installation must be performed by a qualified expert in compliance with the national regulations. DIN VDE 0100-710 applies in Germany.

NOTE

Restrictions regarding installation site

The system is not intended for operation in areas subject to explosion hazards.

NOTE

Do not damage the unit!

The unit can be damaged if opened improperly.
It is expressly prohibited to open the unit with tools!

2.1.2 Operating

The suction system must only be operated with an integrated filter bag and HEPA filter. Only the approved filter bags and HEPA filters may be used (see "Spare parts and consumables [→ 26]").

CAUTION

Dry processing

Before using dry processing, verify the functioning, correct connection and the tightness of the connections. All available suction openings must be available.

CAUTION

Dry processing

Before dry processing, check that the CEREC suction device contains a functioning HEPA fine-dust filter.

2.1.3 Maintenance and repair

As manufacturers of dental instruments and laboratory equipment, we can assume responsibility for the safety properties of the unit only if the following points are observed:

- The maintenance and repair of this unit may be performed only by Dentsply Sirona or by agencies authorized by Dentsply Sirona.
- Components which have failed and influence the safety of the unit must be replaced with original (OEM) spare parts.
- Only original cables may be used, so that EMC requirements are met.

Please request a certificate whenever you have such work performed. It should include:

- The type and scope of work.
- Any changes made in the rated parameters or working range.
- Date, name of company and signature.

2.1.4 Modifications to the product

Modifications to this product which may affect the safety of the operator, patients or third parties are prohibited by law!

2.1.5 Accessories

In order to ensure product safety, this device may be operated only with original Dentsply Sirona accessories or third-party accessories expressly approved by Dentsply Sirona. In particular, only the power cable also supplied or the corresponding original spare part may be used with the unit. The user is responsible for any damage resulting from the use of non-approved accessories.

2.1.6 In case of damage

In case of noticeable malfunctions or damage, stop using the instrument immediately and notify your authorized dealer or the manufacturer.

2.2 Connecting the unit

The connection must be established in accordance with the directions provided in these operating instructions. Also observe the information in the Operating Instructions for inLab MC X5, CEREC MC / MC X, MC XL Premium Package or CEREC Primemill devices.

If any devices not approved by Dentsply Sirona are connected, they must comply with the applicable standards:

- EN 60 950 for information technology equipment, and
- EN 61 010-1 for laboratory equipment.

2.3 Ventilation slots

Under no circumstances may the ventilation slots on the unit be covered, since otherwise the air circulation will be obstructed. This can cause the unit to overheat.

3 General information on the device

3.1 Standards/ approvals/ certifications

CE mark

This product bears the CE mark in accordance with the provisions of directives 2006/95/EC (Low Voltage Directive), 2006/42/EC (Machinery Directive), and 2004/108/EC (EMC Directive). As such, the following standards apply: CISPR 14-1 + A1 + A2, CISPR 14-2, EN55014-1 + A2, EN55014-2, EN61000-3-2, EN61000-3-3, EN60335-1, and EN60335-2-65.

CAUTION

CE mark for connected products

Further products which are connected to this unit must also bear the CE mark. These products must be tested according to the applicable standards.

3.2 Scope of supply

inLab MC X5 230V (REF 6385129) suction device:

- 230 V suction device
- Filter bag (installed)
- HEPA filter (installed)
- Suction tube, complete
- Interface cable
- Operating Instructions
- Power cable (pre-assembled separately)

inLab MC X5 120V (REF 6485929) suction device:

- 120 V suction device
- Filter bag (installed)
- HEPA filter (installed)
- Suction tube, complete
- Interface cable
- Operating Instructions
- Power cable (pre-assembled separately)

Suction device CEREC MC/MC X/MC XL Premium Package/ Primemill 230 V (REF 6569730):


- 230 V suction device
- Filter bag (installed)
- HEPA filter (installed)
- Suction tube, complete
- Interface cable
- Operating Instructions
- Power cable

**Suction device CEREC MC/MC X/MC XL Premium Package/
Primemill 120 V (REF 6569748):**

- 120 V suction device
- Filter bag (installed)
- HEPA filter (installed)
- Suction tube, complete
- Interface cable
- Operating Instructions
- Power cable

3.3 Technical description (component and interfaces)

The suction devices are intended for operation on the inLab MC X5 and CEREC MC/MC X/MC XL Premium Package/Primemill processing machines. The suction system is equipped with a special filter bag and a HEPA micro filter of filter category H12, which filters fine dust up to 99.5% and is classified as fine dust category M.

A	Socket for suction tube	E	On/Off switch
B	Control panel	F	Socket for power cable
C	Maintenance flap	G	Suction power control
D	Socket for interface cable	H	Device safety switch

Check the device for any transport damage immediately after unpacking. If there is transport damage, raise this with your supplier immediately.

3.4 Technical data

3.4.1 InLab MC X5 suction device

Type designation	Suction device for inLab MC X5	
Rated line voltage	120 V ~	230 V ~
Rated power frequency	50 / 60 Hz	50 / 60 Hz
Nominal current	6.4 A	3.6 A
Nominal power output	approx. 1000 W	
Suction power	min. 120 m³/h	
Volume	< 60 dB(A)	
Bag capacity	23 l	
Tube diameter	38 mm	
HEPA filter	Filter category H12, dust category M, as per EN 60335-2	
Permissible line voltage fluctuations	±10% of nominal voltage	
Type of protection against electric shock	Class 1 equipment	
Overvoltage category	II	
Transport and storage conditions:	Temperature:	-25°C – +60°C (-13°F – 140°F)
	Relative humidity:	10%–95%
	Air pressure:	500 - 1060 hPa
	Ambient temperature:	5°C – 40°C (41°F – 104°F)
	Pollution degree:	2 acc. to IEC 60664-1
	Air pressure:	700 - 1060 hPa
	Installation location:	Indoors, ≤ 3,000 m above sea level
	Relative humidity:	max. 80 % at 31°C (87.8°F) decreasing linearly to max. 50 % at 40°C (104°F)
Operating mode	Continuous operation	
Dimensions (WxHxD) in mm	350 x 510 x 360	
Approx. weight	16 kg	

3.4.2 CEREC MC/MC X/MC XL Premium Package/Primemill suction device

Type designation	Suction device for CEREC MC/MC X/MC XL Premium Package/Primemill	
Rated line voltage	120 V ~	230 V ~
Rated power frequency	50/60 Hz	50/60 Hz
Nominal current	6.4 A	3.6 A
Nominal power output	approx. 1000 W	
Suction power	min. 120 m³/h	
Volume	< 60 dB(A)	
Bag capacity	5 l	
Tube diameter	38 mm	
HEPA filter	Filter category H12, dust category M, as per EN 60335-2	
Permissible line voltage fluctuations	±10% of nominal voltage	
Type of protection against electric shock	Class 1 equipment	
Overvoltage category	II	
Transport and storage conditions:	Temperature:	-25 °C – +60 °C (-13 °F – 140 °F)
	Relative humidity:	10% – 95%
	Air pressure:	500 hPa – 1060 hPa
	Ambient temperature:	5 °C – 40 °C (41 °F – 104 °F)
	Pollution degree:	2 acc. to IEC 60664-1
	Air pressure:	700 hPa – 1,060 hPa
	Installation location:	Indoors, ≤ 3,000 m above sea level
	Relative humidity:	max. 80% at 31 °C (87.8 °F) decreasing linearly to max. 50% at 40 °C (104 °F)
Operating mode	Continuous operation	
Dimensions (WxHxD) in mm	360 x 340 x 260	
Approx. weight	12.7 kg	

4 Transportation and installation

4.1 Transport and unpacking

All products from Dentsply Sirona are carefully checked prior to shipment. Please perform an incoming inspection immediately after delivery.

1. Check the delivery note to ensure that the consignment is complete.
2. Check whether the product shows any visible signs of damage.

NOTE

Damage during transport

If the product was damaged during transport, please contact your carrying agent.

If return shipment is required, please use the original packaging for shipment.

4.2 Disposal of packaging materials

The packaging must be disposed of in compliance with the relevant national regulations. Please observe the regulations applicable in your country.

4.3 Installation conditions

4.3.1 Installation site requirements

The inLab MC X5 suction device requires an even base of approx. 400 mm x 400 mm (W x D).


The CEREC MC/MC X/MC XL Premium Package/Primemill suction device requires an even base of approx. 360 mm x 260 mm (W x D).

Install the suction device in such a way that the main switch can be pressed easily. To enable simple replacement of the filter and bag, we recommend orientating the maintenance flap to the front with the two knurled nuts. The ventilation slots on the bottom side and rear side of the suction system must be kept unobstructed. The distance between the back of the unit and the wall must at least be 30 mm. If the suction device is integrated in a cabinet, this must be vented sufficiently to prevent the suction device from overheating.

The suction device must not be set up in the immediate vicinity of considerable humidity development.

4.4 Connecting the suction device

4.4.1 Connecting the suction tube

1. Connect one end of the suction tube to the available connection point (A) on the suction device. The connection point is located on the top side next to the control panel with the suction power controller.
2. Connect the other end of the tube to the socket (B) on the rear side of the processing machine. The figure shows an example of a connection to the CEREC Primemill processing machine.

CAUTION

Trip/fall hazard

If the suction tube is routed poorly, there may be a risk of tripping.

- To prevent injuries caused by tripping, route the suction tube so that there is no risk of tripping.

Notes on the suction tube:

The suction tube is supplied with a length of approx. 2.40 m (inLab MC X5) or 2.0 m (CEREC MC/MC X/MC XL Premium Package/Primemill). When connecting the suction device to the unit, please ensure that no sharp bends occur over the full length of the suction tube.

Reduce the length of the tube according to your requirements and your installation location. Note that suction power drops along the length of the tube. You obtain decent suction power if the suction device is placed immediately below the machine and you have a tube length of 1.20 m or less.

4.4.2 Connecting the power cord

NOTE

Grounded power outlet

The unit must be connected to a grounded power outlet.


1. Insert the power cable into the relevant socket on the suction system.
2. Plug the other end into an appropriate power socket with a protective ground terminal.

Note on the power cable:

The suction system must only be operated with a power cable with a plug system designed for the relevant country.

Check the voltage specification on the rating plate. The system must conform to the country-specific supply voltage.

4.4.3 Connecting the interface cable (for automatic mode)

1. Plug the 15-pole connector into the socket (A) of the processing machine. The figure shows an example of a connection to the CEREC Primemill processing machine.

2. Plug the 9-pole connector into the socket (B) of the suction device.

5 Commissioning

5.1 Putting the unit into operation

5.1.1 Automatic mode

- ✓ The interface cable is connected.
- Set the on/off switch to the position **Auto**.

Note on the automatic mode:

The production unit monitors the correct connection (interface cable and suction tube) and the operation of the suction system during the running processes.

5.1.2 Manual operation

NOTE

Manual mode only permitted for inLab MC X5. For CEREC MC/MC X/ MC XL Premium Package/Primemill only automatic mode is supported.

1. Prior to commencing processing, set the on/off switch to the **On** position.
2. Once processing has been completed, set the on/off switch to the **Off** position.

5.1.3 Setting the suction power

inLab MC X5

- Set the suction power to the maximum (**max**).

CEREC MC/MC X/MC XL Premium Package/Primemill

- Set the suction power to the minimum (**min**). If required, the suction power can be set to the maximum.

5.2 Repacking

1. Set the on/off switch to the position **Off**.
2. Remove the power plug, interface cable, and suction tube.
3. Check the suction device for completeness according to the scope of supply.
4. Pack the device safely, ideally in the original transport packaging.

6 Maintenance and repair

CAUTION

Risk of uncontrolled release of fine dust

If the suction device is not operated with the inserted filter bag or HEPA filter provided, fine dust may be discharged uncontrollably.

- Only operate the suction system with original spare parts (see “Spare parts and consumables [→ 26]”).

6.1 Changing filter bags and HEPA filters

Changing filter bags

If there is a significant fall in suction power, the filter bag could be full and will need replacing.

NOTE


For the CEREC suction device, we recommend changing the filter bag at the following intervals:

CEREC MC XL / MC X / MC: after approx. 120 restoration units.


After 120 dry-milled restorations, a notice appears on the home screen of the touch interface. This number may differ depending on the amount of zirconium oxide material milled and suctioned.

CEREC Primemill: after around 8 hours of processing.

After 8 hours of processing, a yellow warning appears on the home screen of the touch interface. If the filter bag is not changed, a red error message appears after 10 hours of processing.

1. Loosen the two knurled nuts (A) on the maintenance cover.
2. Remove the lid.

3. Remove the filter bag from the nozzle and put a new filter bag on.
4. Put the maintenance cover on and screw it down with the two knurled nuts.

NOTE

Do not jam the filter bag

Make sure that the cap is sealed properly and the filter bag is not jammed.

Changing HEPA filters

Replace the HEPA filter if suction power is still weak after replacing the filter bag. The HEPA filter is located behind the filter bag.

NOTE


For the CEREC suction device, we recommend changing the HEPA filter bag at the following intervals:

CEREC MC XL/MC X MC: After changing the filter bag approx. 3 to 4 times.


This number may differ depending on the amount of zirconium oxide material milled and suctioned.

CEREC Primemill: after 32 hours of processing.

After 32 hours of processing a yellow warning appears on the home screen of the touch interface. If the filter bag is not changed, a red error message appears after 40 hours of processing.


1. Loosen the two knurled nuts (A) on the maintenance cover.
2. Remove the lid.


3. Remove the filter bag from the nozzle.

4. Then release and remove the two Phillips screws on the perforated sheet on the inside of the suction system.

5. Take the perforated sheet out.

6. Remove the dusty HEPA filter and insert a new HEPA filter.
7. Put the perforated sheet back on and screw it down with the two Phillips screws.
8. Put the filter bag back on again.
9. Put the maintenance cover on and screw it down with the two knurled nuts.

NOTE


Do not jam the filter bag

Make sure that the cap is sealed properly and the filter bag is not jammed.

6.2 Changing the suction motor

The motor of the suction device has a service life of approx. 1,200 hrs. You will notice the motor is worn when suction does not start in either automatic or manual mode. When this happens, have the motor replaced by a trained service engineer.

6.3 Using the cleaning hose and the wet cleaning process

⚠ CAUTION

Risk of injury on calibration pins/instruments

If you put your hand in the production chamber, you could injure it on the calibration pins / instruments.

Be careful not to brush against the calibration pins/instruments with your hand.

The processing chambers of devices used for dry milling of zirconium oxide should be cleaned regularly in order to prevent deposits of zirconium oxide dust in the chamber. A manual suction set with cleaning hose can be used before the wet cleaning process in order to support the user when vacuuming zirconium oxide dust out of the device's processing chamber. This set is available as a spare part (REF 67 21 307). We recommend vacuuming the chamber after each dry-milled restoration and performing a wet cleaning process (or wet grinding of a restoration) at least once a week.

Proceed as follows to use the manual suction set and to execute the recommended wet cleaning process:

1. **CEREC MC XL / MC X / MC:** Open the door of the production chamber. The motors are in the home position.


- CEREC Primemill:** Open the door of the production chamber and activate the cleaning position in the *"Routine Actions"* area of the touch interface in order to bring the motors together so that the production chamber is easier to access and clean. This function can be started via the *"Routine Actions"* area (symbol on left) by selecting the *"Cleaning Position"* item.
2. Rotate the black adapter of the manual cleaning hose set from suction of the chamber to suction through the manual cleaning hose.
 3. Activate the manual operation for the suction by setting the switch at the top of the device from "Auto" to "On".
 4. Suck up the dust from the chamber as required.
 5. Sucking up as much dust as possible from the workpiece axis is recommended. Also remove the dust from the bottom side of the block axis as well as on the bearing units and scrapers of the right and left machining axes (see arrows).
 6. Close the door.
 - ↳ The motor mounts and the tool axis move back to initial position.
 7. **CEREC MC XL / MC X / MC:** Once the chamber has been cleaned using the suction device, you should complete the 30-second wet cleaning process programmed in the software 2 to 3 times as required in order to remove any accumulated material.
 - In the system menu of the software, navigate to *"Configuration"* and click on the button *"Devices"*.
 - Select the corresponding unit.
 - Click on the option *"Service"* and then on *"Start a cleaning process."* Alternatively, you can also wet grind a normal restoration with the same effect.
 8. **CEREC Primemill:** After vacuuming the chamber, you should perform the wet cleaning of the machine. This function can be started via the *"Routine Actions"* area (symbol on left) by selecting the *"Cleaning Program"* item. There are two cleaning options: 2 minutes and 15 minutes. Each cleaning option can be stopped at any time when the desired results have been achieved.
 8. You can remove any zirconium oxide deposits in the lower part of the workpiece axis with a nylon brush between the wet cleaning processes.

7 Malfunctions

7.1 Error messages

Error	Possible cause	Remedy:
Suction air not present; suction device works when on/off switch is in the On position.	Suction tube blockage	Check suction tube for open passage; remove blockages.
	Leak on suction tube or socket	Check suction tube for leak tightness and tight connection to the suction unit and processing machine; seal the tube / connect it tightly if required.
	Fault inside the processing machine	Contact services
Software error message: Suction device vacuum switch	Suction power controller is at min on the control panel.	Set suction power controller to max .
	Filter bag / HEPA filter occupied	Replace filter bag Replace HEPA filter
	Suction tube blocked	Check suction tube for open passage; remove blockages if necessary.
	Interface cable defective	Replace interface cable
	Fault inside the processing machine	Contact services

7.2 Troubleshooting and repair

Symptom	Possible cause	Remedy:
Suction device does not start up, either in automatic mode or with the on/off button set to On.	Power cable not connected.	Check power cable for correct connection.
	Voltage supply inactive.	Check voltage source
	The fuse on the suction device has tripped.	 Press the device safety switch, see Technical description (component and interfaces) [→ 12]
	Suction motor superheated.	Let the suction system cool down; ensure there is sufficient air circulation, especially around the ventilation slots.
	Suction motor defective / worn out.	Contact services
	Suction electronics defective.	Contact services
Insufficient suction power	Suction power controller is at min on the control panel.	Set suction power controller to max.
	Filter bag is full.	Replace filter bag
	HEPA filter occupied.	Replace HEPA filter
	Suction tube blocked.	Check suction tube for open passage; remove blockages if necessary.
	Leak on suction tube or socket.	Check suction tube for leak tightness and tight connection to the suction unit and processing machine; seal the tube / connect it tightly if required.
Very strong suction power	Filter bag / HEPA filter not seated correctly.	Open maintenance cover. Check filter bag and HEPA filter are seated correctly.
Suction very loud	Turbine touches the housing.	Contact services

8 Disposal

In accordance with Directive 2012/19/EU and national disposal regulations regarding old electrical and electronic devices, please be advised that such items must be disposed of in a special way within the European Union (EU). These regulations require the environmentally friendly recycling/disposal of old electrical and electronic devices. Such items must not be disposed of as domestic refuse. This has been expressed using the icon of the “crossed out trash can”.

Disposal procedure

We feel responsible for our products from the first idea to their disposal. For this reason, we give you an option to return our old electronic and electrical devices.

If you wish to dispose of your devices, please proceed as follows:

In Germany

To initiate return of the electrical device, please send a disposal request to enretec GmbH. You have the following options here:

- Use the ‘Returning an electrical device’ button under the ‘eom’ menu item on the enretec GmbH homepage (www.enretec.de).
- Alternatively, you can also contact enretec GmbH directly.

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In accordance with the national disposal regulations regarding old electrical and electronic devices (ElektroG), as the manufacturer, we assume the costs for disposing of the electrical and electronic devices in question. Disassembly, transport and packaging costs shall be borne by the owner / operator.

Prior to disassembly/disposal of the unit, it must be prepared professionally (cleaned/disinfected/sterilized).

If your unit is not permanently installed, it will be collected from the practice. If it is permanently installed, it will be picked up curbside at your address by appointment.

Other countries

For country-specific information on disposal, contact your local dental dealers.

9 Spare parts and consumables

REF	Designation
6385244	Filter bag, double pack (inLab MC X5)
6578095	Filter bag, double pack (CEREC)
6385277	HEPA filter
6721307	Manual suction set
6739614	Manual switching valve, spare part
6739622	Hose, manual suction, spare part
6739630	Vacuum cleaner nozzle with brush, spare part
7047202	O-ring 32.0 x 2.0 -S-70
6384841	Suction motor, 230 V
6384791	Suction motor, 120 V
6453851	Suction tube, complete, spare (inLab MC X5)
6564590	Suction tube, complete (CEREC)
6453869	Interface cable, spare
5864553	Power cable (EU)
5864520	Power cable (GB)
5864546	Power cable (IT)
5864538	Power cable (CH)
5864512	Power cable (US)
5865147	Power cable (JP)
5966119	Power cable (AU/NZ)
6412766	Power cable (BR)

We reserve the right to make any alterations which may be required due to technical improvements.

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© Sirona Dental Systems GmbH