EXPRESS™

Glaucoma Filtration Device

PATIENT INFORMATION BROCHURE:

[EN] - ENGLISH	1-3
[NL] - DUTCH / NEDERLANDS	4-6
[BG] - BULGARIAN / БЪЛГАРСКИ	7-9
[CS] - CZECH REPUBLIC / ČESKY	10-12
[DA] - DANISH / DANSK	13-15
[ET] - ESTONIAN / EESTI	16-18
[FR] - FRENCH / FRANÇAIS	19-21
[DE] - GERMAN / DEUTSCH	22-24
[EL] - GREEK / ΕΛΛΗΝΙΚΑ	25-27
[HU] - HUNGARIAN / MAGYAR	28-30
[IT] - ITALIAN / ITALIANO	31-33
[LV] - LATVIAN / LATVIEŠU	34-36
[LT] - LITHUANIAN / LIETUVIŠKAI	37-39
[NO] - NORWEGIAN / NORSK	40-42
[PL] - POLISH / POLSKI	43-45
[PT] - PORTUGUESE / PORTUGUÊS	46-48
[RO] - ROMANIAN / ROMÂNĂ	49-51
[RU] - RUSSIAN / РУССКИЙ	52-54
[SK] - SLOVAK / SLOVENSKY	55-57
[ES] - SPANISH / ESPAÑOL	58-60
[SV] - SWEDISH / SVENSKA	61-63
[TR] - TURKISH / TÜRK	64-66
[UA] - UKRAINIAN / УКРАЇНСЬКИЙ	67-69

Languages not translated are listed on the back page

Alcon

EX-PRESS™ GLAUCOMA FILTRATION DEVICE

PATIENT INFORMATION BROCHURE

www.ifu.alcon.com

CE 0123

This brochure has been written to assist you and your eye doctor regarding the use of the EX-PRESS Glaucoma Filtration Device following glaucoma surgery.

INFORMATION ABOUT GLAUCOMA

What is glaucoma?

Glaucoma is a group of eye conditions that damage the *optic nerve*, which is the nerve that allows you to see. If left untreated, *glaucoma* leads to permanent vision loss and blindness. With treatment, however, you can often protect your eyes from serious vision loss.


Figure 1 - Eye Anatomy

What causes glaucoma?

Although an exact cause is unknown, glaucoma is usually associated with abnormally high eye pressure also called *intraocular pressure* (or *IOP*). Eye pressure is often measured by eye surgeons in order to see if you might have glaucoma, even if you have not yet had any vision loss. Your eyes are constantly producing fluid. This fluid provides nourishment to the eye and is essential for normal, healthy eye function. As fluid is produced, it moves throughout your eye and eventually exits through tiny drainage channels located near the *cornea* (the front part of your eye). If these drainage channels become blocked, the fluid in your eye can build up, causing elevated eye pressure. This elevated eye pressure can damage your optic nerve (Figure 1).

INFORMATION ABOUT GLAUCOMA FILTRATION DEVICES & GLAUCOMA SURGERY

What is the EX-PRESS Glaucoma Filtration Device?

The EX-PRESS Glaucoma Filtration Device is designed to allow fluid to flow out of the eye, reducing your IOP or eye pressure. This is one treatment option your eye doctor may use if you have glaucoma. EX-PRESS Glaucoma Filtration Device comes in 2 models, the P-50 and the P-200, and your eye surgeon will decide which model is right for you.

The EX-PRESS Glaucoma Filtration Device is made from stainless steel. It consists of a 2.64 mm long and 0.4 mm outer diameter tube (Figure 2) and is designed to allow fluid to drain from your eye and reduce IOP or eye pressure.

Figure 2 - EX-PRESS Glaucoma Filtration Device

What is the surgery to place the EX-PRESS Glaucoma Filtration Device?

The EX-PRESS Glaucoma Filtration Device supports the natural drainage fluid pathway of the eye to reduce the eye pressure or IOP. It is placed within the eye through a surgery called "trabeculectomy" or glaucoma filtration surgery, which is usually done using a numbing medicine and relaxing medications. An eye surgeon makes a small flap cut in the white part of the eye (sclera) to place the EX-PRESS Glaucoma Filtration Device to allow the fluid to drain from the eye.

EX-PRESS GLAUCOMA FILTRATION DEVICE MATERIAL

What is the EX-PRESS Glaucoma Filtration Device made of?

The EX-PRESS Glaucoma Filtration Device is made

of surgical grade stainless steel. This material has a history of over 30 years of use in the eye.

SAFETY INFORMATION FOR EX-PRESS™ GLAUCOMA FILTRATION DEVICE SURGERY

Potential Side Effects Associated with the EX-PRESS Glaucoma Filtration Device

Surgical treatment of glaucoma to relieve the pressure from inside of the eye has been used to help control eye pressure for more than 40 years. However, there may be complications as with any other surgical procedure.

As with any surgery, there is risk involved whether or not an EX-PRESS Glaucoma Filtration Device is used. Complications of glaucoma surgery range from minor, usually temporary side effects, to major complications that could affect vision. There is a small risk of requiring further surgical treatment after the initial implantation.

Clinical studies were performed to evaluate the EX-PRESS Glaucoma Filtration Device with a follow-up period of one year. The testing was done on a total of 113 patients implanted with the R-30 and R-50 versions with 58 patients completed a one-year follow-up. The R-30 and R-50 versions are earlier versions of the EX-PRESS Glaucoma Filtration Device, which are similar to the P-versions.

The side effects reported in the clinical studies that occurred in more than 5 out of 100 patients are shown below:

• Additional Surgeries:

- Surgery performed to drain built-up fluid in the eye (*Glaucoma filtration surgery*).
- Having to remove the device from the eye (*Device removal*).
- Surgery on the blister in the eye (*created during the implantation surgery*) with or without certain medications (*Revision of bleb using antimetabolites or without using antimetabolites*).

• Issues with Device:

- Device touches the colored layer of the eyeball (*Device-Iris Touch*).
- Implanted device becomes visible or touches inappropriate eye tissue (*Device exposure*).

• Issues with Fluid inside eye:

- Fluid-filled space inside the eye becomes too shallow (*Shallow Anterior Chamber*).
- Reforming the fluid-filled space inside the eye (*Anterior Chamber reformation*).
- Pressure in the eye is too low (*Hypotony*).

• Other relevant side effects:

- Small amount of blood in eye (*Hyphema*).
- When a hard mass forms on the thin connective tissue in the back of the eyeball (*Tenon cyst*).
- Problem with transparent layer in front of the eye (*Corneal complication*).

NOTE: The patients should ask their doctor if they have any questions related to these potential side effects.

RECOVERY FROM GLAUCOMA FILTRATION SURGERY

Important Instructions after Glaucoma Filtration Surgery:

- You should avoid any activity that could harm your eye while you are recovering from surgery.
- Take all prescribed medicines and apply eye drops as instructed.
- Contact your eye doctor immediately if you have any of the following symptoms after surgery:
 - *Significant decrease in vision*
 - *Significant eye pain,*
 - *Significant itching of the eye,*
 - *Significant redness of the eye,*
 - *Watering of your eye,*
 - *Double vision lasting more than a few days,*
 - *Significant eye discharge, or*
 - *Increased sensitivity to light.*

These symptoms could indicate potential serious postoperative complications, and you should follow up with your eye doctor immediately.

Magnetic Resonance (MR) Conditional

The device is MR-Conditional, meaning that special precautions must be taken prior to any MR testing (i.e. MRIs.), which are outlined in the doctor's product label.

IMPORTANT: Please tell your doctor about your EX-PRESS Glaucoma Filtration Device implant before undergoing any MR testing.

Immediate Postoperative Care Instructions

Immediately after surgery, the operated eye will be red and irritated, and there may be increased tears. Plan to have someone else drive you home. You will return home after surgery after a short stay in the outpatient recovery room. Typically, your eye doctor will examine you the following day. Your eye doctor

will give you eye drops to speed up the healing process and to prevent infection.

The specifics of surgery may be different for each individual. Be sure to consult your eye doctor so you can fully understand the recovery process after the glaucoma filtration surgery. Your eye doctor will tell you when to return for a medical visit(s) to check the health of your eye after surgery.

Longer Term Considerations

Recovery from this kind of eye surgery can take up to 8 weeks. During this time, vision will be blurry, and your eye doctor will want to see you often, usually on a weekly basis, to check eye pressure and adjustments as needed. Your eye doctor should follow you closely during the first year after implantation (at least 4 times), and at least once a year as long as you have the device or as recommended by your eye doctor to monitor your eye pressure.

Keep your patient identification card given to you by your eye doctor at the time of surgery as it has important identification information for the device implanted in your eye. In addition, the patient identification card has the website where you can print a copy of this patient information brochure, including any recent revisions of this brochure.

The website that contains the electronic form of this patient brochure is: www.ifu.alcon.com

If you experience a serious health issue related to the EX-PRESS™ Glaucoma Filtration Device, do the following:

- Consult with your eye doctor.
- Report this health issue to Alcon by (1) calling the local Alcon representative in your Member State or (2) emailing Alcon at <http://www.alcon.com/contact-us/>.
- Contact the competent authority of your Member State about this serious health issue related to the EX-PRESS Glaucoma Filtration Device.

Device Lifetime

The EX-PRESS Glaucoma Filtration Device will remain in the patient's eye after implantation and will not be removed unless the doctor recommends it.

SUMMARY

The EX-PRESS Glaucoma Filtration Device is designed to reduce your IOP or eye pressure associated with your glaucoma. Always consult with your eye doctor if you have any questions or concerns as a result of glaucoma surgery.

Glossary For The Implant Card

Symbol / Abbreviation	Definition
L	Left
R	Right
 www.ifu.alcon.com	Electronic instructions for use (a)
	Manufacturer (b)
	MR Conditional (c)
	Authorized Representative in the European Community (d)

- (a) In other words, this patient brochure is available electronically at the following website: www.ifu.alcon.com
- (b) In other words, the person or company that created the device.
- (c) In other words, precautions need to be taken before being exposed to any MRI Equipment.
- (d) In other words, when a manufacturer is located outside the European Union, a person is chosen that is located in the European Union to act on behalf of the manufacturer.

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© 2021 Alcon Inc.
Alcon Laboratories, Inc.