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
**The 37° Company**  
part of The Surgical Company

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...conform extras beozurii pag. 5

## 2 Contraindications, warnings, cautions, notes and symbols

The device was designed and built with safety in mind. However, there is no replacement for care providers being attentive to their patients' needs and equipment operation. Read and understand the contraindications, warnings, cautions and notes before using the device.

### 2.1 Contra-indications

- Only apply heat to intact skin and do not apply heat directly to open wounds.
- Do not apply the warming system to ischemic limbs.
  1. Use caution and consider discontinuing use on patients during vascular surgery when an artery is clamped to an extremity (i.e. aortic cross-clamping).
  2. Use caution and monitor closely if used on patients with severe peripheral vascular disease.


### 2.2 Warnings

#### Warning!

- Do not use the device when it is damaged or when the Mistral-Air<sup>®</sup> Blanket is damaged. Thermal injury may result.
- Do not allow the patient to lie on or contact the hose with the skin when the device is active. Thermal injury may result.
- Do not use the Mistral-Air<sup>®</sup> Blanket to transfer or move the patient. Injury to the patient may result.
- To prevent tipping when mounting to an IV-pole, mount the device at a height at which the IV-pole is stable. Injury may occur. Before usage, assess the stability by placing the IV-pole on a surface at an angle of 10° from the horizontal plane with brakes activated. The IV-pole may not overbalance, or move. Also passing over a 10 mm threshold may not result in overbalancing. Mass and position of center of gravity are provided in this IFU for theoretical analysis. The 37Company cannot provide maximum mounting height prescriptions for different wheel base diameters, numbers of castors (either with brakes or not) and configurations of other equipment mounted to the IV pole.
- Do not use the device without a Mistral-Air<sup>®</sup> Blanket connected to it (no free hosing). Thermal injury may result.
- Do not use the device and blankets near flammable anesthetics and/or in oxygen-enriched environment, to avoid the risk of explosion or fire.
- 55 • Check patient's temperature and skin condition at least every 15 minutes, or according to institutional protocol.
- Do not cover the patient's thorax with our Mistral-Air<sup>®</sup> Blankets during cardioversion or defibrillation therapy.
- Applying air with a temperature above the normothermic core body temperature range (36 – 37.5°C) incorporates the risk of hyperthermia. Depending on the selected set


### 3.2 Overview of the Mistral-Air® Warming Unit (MA1200-QC)

1. Control panel
2. Mounting clamp
3. Mounting options (accessory)
4. Appliance inlet
5. Equipotential pin
6. Data connection
7. Mains cover
8. Filter (air inlet)
9. Removable hose

### S12 3.3 Overview of the control panel

The control panel is located at the front top of the device and can be operated by pressure sensitive buttons. The device is easy to use. All settings are visible on the control panel and you can select the preferred temperature by pressing the temperature selection buttons. When an alarm condition is detected, **an audible alarm will be activated and an alarm LED will flash yellow.**

1. Standby button
2. Temperature selection indicators
3. Overtemperature alarm LED
4. Technical alarm LED
5. Filter replacement indicator LED
6. Temperature selection + and - buttons
7. Fan only/ambient air indicator
8. Temporarily audible alarm suppression button

- The leads to the temperature sensors are damaged, or disconnected.
- The fan is blocked, or damaged and cannot reach its desired speed.
- The heater is damaged and the desired air temperature is not reached.
- A mains power dip ( $\geq 30\%$ ) occurred for more than 1/60 seconds.

If this alarm occurs, the heater and fan are turned off and the device enters standby mode. It is impossible to start the device by pushing the standby button. Control of the device can only be recovered through resetting the device. To do this, disconnect the mains plug and reconnect it.

#### Caution!

If this alarm occurs, check for anything blocking the air flow path. If the technical alarm continues, take the device out of use and contact the hospital service department or the local supplier.

### 3.4.2 Overtemperature alarm

The overtemperature alarm is triggered with a maximum air temperature of 56°C.

When the overtemperature alarm occurs, the technical and overtemperature LED's flash yellow. These indicate that the air temperature is too high. The visual alarms are accompanied by an audible alarm, which consists of three pulses of 200 ms duration, with 200 ms spacing between each pulse. This pattern is repeated every 6 seconds.

If this alarm occurs, the heater and fan are turned off and the device enters standby mode. It is impossible to start the device by pushing the standby button. Control of the device can only be recovered through resetting the device. To do this, disconnect the mains plug and reconnect it.

#### Caution!

If this alarm occurs, check for anything blocking the air flow path. Ensure that the blanket is not folded and do not place tools/equipment on the blanket which could result in a blocked air flow. Ensure that the air inlet is free. If the overtemperature alarm continues, take the device out of use and contact the hospital service department or the local supplier.

### 3.4.3 Microcontroller watchdog alarm

The microcontroller watchdog alarm is visually indicated by a continuous yellow technical alarm LED and a flashing yellow overtemperature LED. The visual alarms are accompanied by a continuous single tone audible alarm.

The microcontroller watchdog alarm indicates a technical malfunction and is triggered when the microcontroller is not functioning properly. If this alarm occurs, the heater and fan are turned off and the device enters standby mode. It is impossible to start the device by pushing the standby button. Control of the device can only be recovered through resetting the device. To do this, disconnect the mains plug and reconnect it.

#### Caution!

If this alarm occurs, send the device to a certified service department for technical support.

## 6 Operation

### 6.1 Safety instructions before operation

#### Warning!


When using the device, first read the warnings in *Warnings* on page 8.

The device is intended to be used only by trained clinicians. Intended patient population: adults and pediatric patients.

The clinical areas are: operating room, recovery room, anaesthetic room, intensive care unit, medical/surgical floors and emergency room. Mainly used during the entire perioperative pathway (pre-, per-, and postoperative period).

### 6.2 Connecting the Mistral-Air® QC Hose (MA1200-1018 & MA1200-1018XL) to the device

1. Check the QC Hose or QC Hose XL for damage.
2. Hold on to the front-end hose connector (knee part).
3. Place the front-end of the hose into the quick connector.
4. Press the connector firmly into the blower until a "click" is produced.
5. Check to ensure the hose is firmly connected.

The QC Hose can be rotated 360° for optimal range.

### 6.3 Connecting the power supply

See *Connecting the power supply cord* on page 27 if the power supply cord is not yet connected to the device.

1. Plug the device into an earthed mains socket.
2. The device automatically switches to the standby mode, which is indicated by the orange standby LED located on the left side of the control panel.

## 6.4 Connecting the blanket

- 511 {
1. Take the selected Mistral-Air® blanket out of the package and follow the instructions on the insert provided with the blanket box.
  2. Place the unit near the hose inlet of the blanket.
  3. Insert the end of the flexible hose into the air inlet port of the Mistral-Air® blanket.
  4. Check if the hose is fully pushed in.

## 6.5 Turning on the device

### Caution!

Stay in viewpoint of the control panel when the device is performing the self-test and selecting the set-point.

- 53 {
1. Activate the device by pressing the standby button. The LED turns green.
  2. The device will now perform a self-test, which includes a flash of all the LED's and a short beep. When a LED or the audible beep fails, take the device out of use for repair.
  3. After the self-test, which lasts several seconds, the device will start blowing air at the default temperature setting of 38°C.

## 6.6 Selecting the temperature

52 { The description of the setpoints corresponds to the average temperature under the blanket. There are four temperature setpoints:

- Fan only/ambient air: Ambient air temperature. The air temperature to the patient will depend on ambient conditions and possible heat from the fan motor.
- 32°C: Low temperature.
- 38°C: Medium temperature.
- 43°C: High temperature.

The selected temperature is indicated by one of four green temperature selection indicator LED's, see *Overview of the control panel* on page 18.

1. After the self-test, the device will start blowing air at the default temperature setting of 38°C.
2. Press the – button twice to deactivate the heater. The fan only/ambient air indicator turns green and air at ambient temperature is blown to the patient.
3. Press the + button to activate the heater.
4. Press the + button multiple times to increase the air temperature at the blanket to a setpoint of 38°C, or 43°C.

## 7 Maintenance

### Warning!

- Maintenance may only be performed by trained biomedical technicians or engineers. Both user groups must be trained by certified trainers from The 37Company.
- Preventive maintenance needs to be performed on an annual basis. Please refer to the Mistral-Air<sup>®</sup> technical manual for maintenance, repair and calibration instructions. The Mistral-Air<sup>®</sup> technical manual is available for download at the business partner menu of the The 37Company website.


### Caution!

Clinical users may not repair or open the device in the event of a malfunction. This can damage the appliance and will invalidate the warranty.

Have the device serial number ready when you contact the hospital service department or the local supplier for technical support. The serial number is located on the label on the back of the device.

Before performing maintenance, consult the device insulation diagram below.

## 9 Specifications

### 9.1 Specifications of the device

#### General specifications

Article number	MA1200-EU	MA1200-US
Rated voltage	220-240 V~	100-125 V~
Frequency	50/60 Hz	50/60 Hz
Device sound pressure	48 dBA	51 dBA
Average current	3.2 A	6.1 A
Peak current	9.0 A	10.0 A
Peak power	1000 VA	1000 VA
Average power at 43°C and $t_{amb} 22 \pm 1.5^\circ\text{C}$	750 VA / 600 W	800 VA / 610 W
Fuses	6.3AHF/250V~	10AHF/250V~
<b>Air flow rate at nominal voltages and ambient temperature (1.8 m hose)</b>	<b>Up to 101 [Nm<sup>3</sup>/h] or 59 [CFM]</b> (depending on orientation of the hose, supply voltage, type and drape of the blanket)	<b>Up to 88 [Nm<sup>3</sup>/h] or 52 [CFM]</b> , (depending on orientation of the hose, supply voltage, type and drape of the blanket)
Air flow rate at nominal voltages and ambient temperature (1.8 m hose), low fan duty cycle	Up to 91 [Nm <sup>3</sup> /h] or 53 [CFM] (depending on orientation of the hose, supply voltage, type and drape of the blanket)	Up to 79 [Nm <sup>3</sup> /h] or 47 [CFM], (depending on orientation of the hose, supply voltage, type and drape of the blanket)
GMDN code	36954 (circulating-air whole-body heating system control unit)	
Dimensions	16 cm x 35 cm x 40 cm (l x w x h)	
Weight	5.2 kg	
Hose length	1.8 m (3 m hose optional)	
Power supply cord length	4.0 m	
Filtration	HEPA H13 class filter, conform EN 1822-1:2009	
Classification 93/42/EEC	Class IIb	
Classification IEC 60601-1	Class I, Body Floating (BF)	
Overvoltage category according to IEC 60664-1	Category II	
Classification IEC 60529	IP23	
Setpoint temperature	Ambient air, 32°C, 38°C, 43°C	

54 }

53

Accuracy of temperature at the end of the hose	± 2.5 °C (under all validated operating environmental conditions)
<b>Setpoint reached after</b>	<b>Under 30 seconds</b>
Low temperature limit	10°C
Maximum average contact surface temperature	45.5°C (Compliant with IEC 80601-2-35)
High temperature safety limit	< 56°C (Compliant with IEC 80601-2-35)
Auditory alarm signal sound pressure	80 dBA
Applicable technical standards	IEC 60601-1:2005+A1:2012, IEC 80601-2-35:2016
Expected lifetime device	7 years
Expected lifetime hose	1 year

The essential performance of the Mistral-Air Warming System is: when supplying air to the patient "the Maximum CONTACT SURFACE TEMPERATURE" must be below the safe temperature limits according to IEC 80601-2-35: 2009 + C1:2012 +C2:2015, Clause 201.11.1.2.1.102.

#### Validated operating environmental specifications

Ambient temperature	15°C to 30°C
Relative humidity	30% to 75%

#### Validated transport and storage specifications

Ambient temperature	-20°C to 70°C
Relative humidity	10% to 90% (non-condensing)
Atmospheric pressure	50 kPa to 106 kPa