Foot-and-Mouth-Disease (FMD) Multispecies Antibody Test Kit

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Test With Confidence*

IDEXX

Foot-and-Mouth-Disease (FMD) Multispecies Antibody Test Kit**For veterinary use only.****Name and Intended Use**

IDEXX FMD Multispecies Ab Test provides a rapid, simple, sensitive and specific method for detecting antibodies against non-structural proteins (NSP) of Foot-and-Mouth-Disease Virus (FMDV) in serum and plasma samples of bovine, caprine, swine and ovine origin. This test allows differentiation between samples from infected (presence of antibodies against NSP of FMD virus) and vaccinated (no antibodies against NSP of FMD virus) animals.

Descriptions and Principles

Microtiter plates are coated with NSP of FMDV. Dilutions of the samples to be tested are incubated in the wells of these plates. Any antibody specific for NSP binds to the antigen in the wells and forms an antigen/antibody complex on the plate well surface. Unbound material is removed from the wells by washing. Conjugate is added, that binds to the antibodies of the sample complexed with the antigen. Unbound Conjugate is removed by washing, and the TMB Substrate is added to the wells. The degree of color that develops (optical density measured at 450 nm) is directly proportional to the amount of antibody specific for NSP present in the sample. The result is obtained by comparing the optical density (OD) that develops in wells containing the samples with the OD from the wells containing the Positive Control.

Reagents**Volume**

1	FMDV NSP Coated Plate	5
2	Positive Control	1 x 0.9 mL
3	Negative Control	1 x 0.9 mL
4	Conjugate	1 x 60 mL
5	Sample Diluent	1 x 100 mL
A	TMB Substrate N.12	1 x 60 mL
B	Stop Solution N.3	1 x 60 mL
C	Wash Concentrate (10X)	1 x 480 mL
Other Components: Zip lock bag		1

Note: See table at the end of the insert for a description of symbols used on the insert and labels of this kit.

Storage

Store the reagents at 2–8°C. Reagents are stable until expiration date, provided they have been stored properly.

Materials Required but Not Provided

- Precision micropipettes or multi-dispensing micropipettes
- Disposable pipette tips
- Disposable protective gloves
- Graduated cylinder for wash solution
- 96-well microplate reader (equipped with 450 nm filter)
- Microplate washer (manual, semi-automatic or automatic system)
- Use only distilled or deionized water for preparation of the reagents used in the test
- Microplate covers (lid, aluminium foil or adhesive)
- Vortex or equivalent

Precautions and Warnings

- Handle all biological material as potentially infectious.
- Wear protective gloves / protective clothing / eye or face protection when handling samples and reagents.
- Refer to the product Material Safety Data Sheet for additional information.
- See the end of this insert for reagent hazard and precaution warnings.

Laboratory Practices

- Optimal results will be obtained by strict adherence to this protocol. Careful pipetting, timing, and washing throughout this procedure are necessary to maintain precision and accuracy. Use a separate pipette tip for each sample and control.
- Do not expose TMB solution to strong light or any oxidizing agents. Handle TMB solution with clean glass or plastic ware.
- All wastes should be properly decontaminated prior to disposal. Dispose of contents in accordance with local, regional, and national regulations.
- Care should be taken to prevent contamination of kit components. Do not pour unused reagents back into containers.
- Do not use kit past expiration date.

Preparation of Wash Solution

The Wash Concentrate (10X) should be brought to 18–26°C and mixed to ensure dissolution of any precipitated salts. The Wash Concentrate (10X) must be diluted 1:10 with distilled/deionized water before use (e.g., 30 mL of Wash Concentrate (10X) plus 270 mL of water per plate to be assayed). When prepared under sterile conditions, the Wash Solution can be stored for one week at 2–8°C.

Test Procedure

All reagents must be allowed to come to 18–26°C before use. Mix reagents by gentle inverting or swirling. It is important to wear disposable protective gloves during the whole test procedure and when filling racks of pipette tips to ensure the best test performance.

Short Incubation

1 Obtain coated microplates and record the sample position. If using partial plates, remove only those wells sufficient for samples to be tested. Place the remaining wells, along with the desiccant, in the extra zip lock bag provided and return to 2–8°C.

2 Dispense 90 μ L of Sample Diluent into each well.

3 Dispense 10 μ L of Positive Control (PC) into two appropriate wells.

4 Dispense 10 μ L of Negative Control (NC) into two appropriate wells.

5 Dispense 10 μ L of sample into the appropriate wells.

6 Thoroughly mix the contents of the microwells by agitating the plate for at least 30 seconds by hand, or use a microplate shaker.

7 Cover the microtiter plate and incubate for 60 minutes (± 5 min.) at 18–26°C. Continue at Step 8.

Overnight Incubation

1 Obtain coated microplates and record the sample position. If using partial plates, remove only those wells sufficient for samples to be tested. Place the remaining wells, along with the desiccant, in the extra zip lock bag provided and return to 2–8°C.

2 Dispense 95 μ L of Sample Diluent into each well.

3 Dispense 5 μ L of Positive Control (PC) into two appropriate wells.

4 Dispense 5 μ L of Negative Control (NC) into two appropriate wells.

5 Dispense 5 μ L of sample into the appropriate wells.

6 Thoroughly mix the contents of the microwells by agitating the plate for at least 30 seconds by hand, or use a microplate shaker.

7 Cover the microtiter plate and incubate for 12–18 hours at 18–26°C. Continue at Step 8.

After Short or Overnight-Incubation

8 Remove the solution and wash each well with approximately 300 μ L of Wash Solution 5 times. Avoid plate drying between plate washings and prior to the addition of the next reagent. Tap each plate onto absorbent material after the final wash to remove any residual wash fluid.

9 Dispense 100 μ L of Conjugate into each well.

10 After short incubation: Cover the microtiter plate and incubate it for 60 minutes (± 5 min.) at 18–26°C.
After Overnight-Incubation: Cover the microtiter plate and incubate it for 30 minutes (± 5 min.) at 18–26°C.

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13 Incubate at 18–26°C for 10 minutes (±1 min.).

14 Dispense 100 µL of Stop Solution N.3 into each well.

15 Read the results using a photometer at a wavelength of 450 nm.

Note: Make sure to read the plates within 2 hours after the addition of the Stop Solution.

16 Calculation:

Controls

$$NC\bar{x} = \frac{NC1 A(450) + NC2 A(450)}{2}$$

$$PC\bar{x} = \frac{PC1 A(450) + PC2 A(450)}{2}$$

Validity criteria

$$NC\bar{x} \leq 0.300$$

$$PC\bar{x} \leq 2.000$$

$$PC\bar{x} - NC\bar{x} \geq 0.300$$

For invalid assays, technique may be suspect and the assay should be repeated following a thorough review of the package insert.

Samples

$$S/P \% = 100 \times \frac{\text{Sample A}(450) - NC\bar{x}}{PC\bar{x} - NC\bar{x}}$$

17 Interpretation for bovine, caprine and ovine samples:

Negative

$$S/P \% < 35$$

Positive

$$S/P \% \geq 35$$

Interpretation for swine samples:

Negative

$$S/P \% < 55$$

Positive

$$S/P \% \geq 55$$

Note: IDEXX has instrument and software systems available that calculate results and provide data summaries.

For technical assistance:

IDEXX USA Tel: +1 800 548 9997 or +1 207 556 4895

IDEXX Europe Tel: +800 727 43399

Contact your IDEXX area manager or distributor or visit our website: idexx.com/contactlpd

IDEXX and Test With Confidence are trademarks or registered trademarks of IDEXX Laboratories, Inc. or its affiliates in the United States and/or other countries

© 2020 IDEXX Laboratories, Inc. All rights reserved.

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DANGER / DANGER / PERIGO / PELIGRO / GEFAHR / ОПАСНОСТЬ

H317 / P280 / P302+P352 / P333+P313

Positive Control / Negative Control / Wash Concentrate (10X) — May cause an allergic skin reaction. Wear protective gloves. IF ON SKIN: Wash with plenty of soap and water. If skin irritation or rash occurs: Get medical advice/attention.

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H317 / H318 / H412 / P280 / P305 + P351 + P338

Sample Diluent — May cause an allergic skin reaction. Causes serious eye damage. Harmful to aquatic life with long lasting effects. Wear eye protection/face protection. Wear protective gloves. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

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Stop solution — Causes severe skin burns and eye damage. May cause an allergic skin reaction. Causes serious eye damage. Harmful to aquatic life with long lasting effects. Wear protective gloves/protective clothing/eye protection/face protection. IF SWALLOWED: rinse mouth. Do NOT induce vomiting. IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water. IF INHALED: Remove person to fresh air and keep comfortable for breathing. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

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IDEXX Laboratories, Inc.
One IDEXX Drive
Westbrook, Maine 04092
USA

Manufacturer
IDEXX Switzerland GmbH
Stationsstrasse 12
CH-3097 Liebefeld-Bern
Switzerland

EU-Representative
IDEXX Europe B.V.
P.O. Box 1334
2130 EK Hoofddorp
The Netherlands

idexx.com

IDEXX