


# VetMAX™ West Nile Virus Kit


TaqMan® real-time RT-PCR for detection of the West Nile virus

Catalog Number WNPEX050

Doc. Part No. 100020471 Pub. No. MAN0008915 Rev. B.0

| Technology   | Species | Nucleic acid isolated from matrices  | Test type  |
|--|---------|--|------------|
| Real-time RT-PCR (RNA)<br>– Duplex<br>– Endogenous/exogenous IPC | Horse   | Whole blood, serum, plasma<br>Cell culture supernatant<br>Feces<br>Tracheal and cloacal swabs<br>Organs (cervix...)<br>Cerebrospinal fluid | Individual |

 **WARNING!** Read the Safety Data Sheets (SDSs) and follow the handling instructions. Wear appropriate protective eyewear, clothing, and gloves. Safety Data Sheets (SDSs) are available from [thermofisher.com/support](http://thermofisher.com/support).

 **WARNING! POTENTIAL BIOHAZARD.** Read the biological hazard safety information at this product’s page at [thermofisher.com](http://thermofisher.com). Wear appropriate protective eyewear, clothing, and gloves.

## Information about the product

### Description of the product

The **Applied Biosystems™ VetMAX™ West Nile Virus Kit** is a molecular diagnostic tool for detecting West Nile virus by real-time RT-PCR.

Each RNA sample obtained after extraction is analyzed in a single well; the same well is used to specifically detect the viral RNA of West Nile virus and an IPC (Internal Positive Control). The kit simultaneously detects an endogenous IPC in cellular samples and an exogenous IPC to be added to non-cellular samples. A positive IPC reflects both the efficiency of extraction and the absence of inhibitor in the samples.

It can be used on viral RNA extracted from **whole blood, serum and plasma, cell culture supernatant, feces, tracheal and cloacal swabs, organs, and cerebrospinal fluid.**

Complete protocols for viral RNA extraction from these matrices are available upon request from Technical Support.

### Kit contents and storage

The **VetMAX™ West Nile Virus Kit** contains components that can be used for detecting in duplex the West Nile virus and an IPC. Upon receipt, the whole kit must be stored at **–30°C to –10°C**. After initial use of a component, store it according to the following recommendations:

| Component                       | Description  | Volume (50 reactions) | Storage        |                   |
|---------------------------------|--|-----------------------|----------------|-------------------|
|                                 |  |                       | Upon receipt   | After initial use |
| 3 - Mix West Nile (Green tube)  | Mix for TaqMan® RT-PCR. Contains: <ul style="list-style-type: none"> <li>• The detection system for the West Nile virus target, including a TaqMan® probe labeled <b>FAM™ – NFQ</b> (Non-Fluorescent Quencher).</li> <li>• The detection system for IPCs, including several TaqMan® probes labeled <b>VIC™ – TAMRA™</b>.</li> <li>• Buffer, reverse transcriptase and real-time PCR enzyme.</li> </ul> | 2 × 500 µL            | –30°C to –10°C | –30°C to –10°C    |
| 4a - EPC West Nile (Brown tube) | <b>External Positive Control:</b><br>Positive control for West Nile virus It consists of <b>already extracted</b> nucleic acid to be amplified during real-time RT-PCR.  | 90 µL                 | –30°C to –10°C | –30°C to –10°C    |
| 5 - IPC West Nile (Yellow tube) | <b>Internal Positive Control:</b><br>Exogenous internal control <b>to be added to each non-cellular and feces sample and each control</b> in the lysis step of the extraction.   | 250 µL                | –30°C to –10°C | –30°C to –10°C    |

**NOTE:** For small extraction series, it is recommended that the IPC West Nile be aliquoted to avoid more than 3 cycles of freezing/thawing (a minimum volume of 50 µL).

### Extraction and amplification controls

The **VetMAX™ West Nile Virus Kit** contains two controls, enabling validation of the extraction and the amplification of the viral RNA:

#### 4a - EPC West Nile: positive control for West Nile

A positive control, **already extracted**, for amplification during the real-time RT-PCR.

A positive result within the specified C<sub>t</sub> range validates the amplification of the West Nile target by real-time RT-PCR.

#### 5 - IPC West Nile: extraction internal control (optional use depending on the type of sample)

Positive control **to be added to each non-cellular and feces sample during the lysis step** of the nucleic acid extraction. For cellular samples, nucleic acid extraction validation is performed using an endogenous IPC present in each cellular sample.

A positive IPC result with a compliant value in a cellular sample (for an endogenous IPC) or within the specified C<sub>t</sub> range in a sample validates the extraction of this non-cellular or feces sample (for an exogenous IPC), whether positive or negative for the target pathogen, thus eliminating false negatives and verifying the effect of the inhibitors.

**We recommend including two negative controls to confirm correct analysis:**

#### NCS: negative extraction control

This control consists of components used in the extraction without addition of the sample (sample volume can be replaced by the buffer used in the sample preparation or by DNase/RNase-free water) that undergoes the same treatment as the samples: nucleic acid extraction (with or without IPC added) and real-time RT-PCR.

A negative result for the West Nile Virus and the endogenous IPC (for cellular samples) confirms the absence of contamination during the extraction and the real-time RT-PCR.

#### NC: negative amplification control

This control consists of an amplification mix added to the plate during real-time RT-PCR preparation, as well as 5 µL of DNase/RNase-free water to adjust the reaction to 25 µL.

A negative result for the West Nile virus and the IPC confirms the absence of contamination during real-time RT-PCR reaction preparation.

#### Materials required but not provided

Unless otherwise indicated, all materials are available through [thermofisher.com](http://thermofisher.com).

- Precision micropipettes (range of 1 µL to 1000 µL) with DNase/RNase-free filtered tips
- DNase/RNase-free water
- 1X TE buffer
- 1X PBS buffer
- A real-time PCR thermal cycler capable of detecting the following fluorophores:
  - FAM™ (emission maximum: λ515 nm)
  - VIC™ (emission maximum: λ554 nm)
- Optical-quality consumables compatible with the thermal cycler used: PCR 96-well plates, PCR strips (8 or 12 wells), microtubes or capillaries; suitable plate covers or caps for capping

#### Analysis procedure

The real-time PCR reaction volume is 25 µL:

- **3 - Mix West Nile:** 20 µL per analysis
- **Extracted RNA:** 5 µL per analysis

#### Extraction of viral RNA

RNA must be extracted from the samples for real-time RT-PCR analysis.

For non-cellular and feces samples, add **5 µL of 5 - IPC West Nile** to each sample to be extracted and the NCS in the lysis step of the nucleic acid extraction.

**NOTE:** For information about extraction methods that are compatible with and validated for the VetMAX™ West Nile Virus Kit, please contact Technical Support.

#### Preparation of the real-time RT-PCR

1. Create an analysis plan for distribution of the mixes and samples. Keep the positive control (EPC) away from the other samples if possible.
2. Thaw the tube of **3 - Mix West Nile** at **2°C to 8°C on ice** or on a refrigerated rack.
3. Mix the tube of **3 - Mix West Nile** by shaking gently, then centrifuge briefly.
4. Add **20 µL of 3 - Mix West Nile** to each PCR plate well, PCR strip or capillary used.
5. Add RNA from the samples and controls to the reaction mix, according to the pre-defined analysis plan:

| Type of analysis                    | Component                 | Sample volume |
|-------------------------------------|---------------------------|---------------|
| Sample for analysis                 | RNA extracted from sample | 5 µL          |
| Positive amplification control      | <b>4a - EPC West Nile</b> | 5 µL          |
| Negative lysis control (NCS)        | Extracted NCS             | 5 µL          |
| Negative amplification control (NC) | DNase/RNase-free water    | 5 µL          |

6. Cover the PCR plate, PCR strips or capillaries with an adhesive plate cover or suitable caps.

## Amplification by real time RT-PCR

1. Create the following detectors on the thermal cycler:

|  | Reporter | Quencher                       |
|--|----------|--------------------------------|
| WN                                     | FAM™     | NFQ (Non-Fluorescent Quencher) |
| IPC WN                                 | VIC™     | TAMRA™ <sup>(1)</sup>          |
| Passive reference: ROX™ <sup>(1)</sup> |          |                                |

<sup>(1)</sup> The fluorophores TAMRA™ and ROX™ are required for real-time RT-PCR analysis if the thermal cycler is capable of detecting them. For other thermal cyclers, the absence of detection of these fluorophores does not affect the real-time RT-PCR analysis.

2. Assign the **WN** detector and the **IPC WN** detector to each sample well used in the analysis.

3. Set up the following real-time RT-PCR program for the analysis:

|        | Step repetitions | Temperature         | Duration   |
|--------|------------------|---------------------|------------|
| Step 1 | ×1               | 45°C                | 10 minutes |
| Step 2 | ×1               | 95°C                | 10 minutes |
| Step 3 | ×40              | 95°C                | 15 seconds |
|        |                  | 65°C <sup>(1)</sup> | 45 seconds |

<sup>(1)</sup> Collection of fluorescence data during the 65°C – 45 seconds stage.

4. Place the PCR plate, the PCR strips or the capillaries in the thermal cycler and run the real-time RT-PCR.

## Analysis of the results

### Analysis of the raw data

Refer to the recommendations of the thermal cycler manufacturer for the analysis of the raw data.

1. Position the threshold limits separately for each target of the real-time RT-PCR.
2. For each detector, interpret the results according to the sample  $C_t$  values obtained as recommended below.

### Validation

The test is validated if the following criteria are met:

|               | West Nile detector   | West Nile IPC detector  | Validation               |
|---------------|--|---|--------------------------|
| EPC West Nile | $C_t = C_{t\text{oc West Nile of 4a}} - \text{EPC West Nile} \pm 3C_t^{(1)}$ | For the endogenous IPC:<br>$C_t < 40$ or $C_t > 40^{(2)}$                                     | PCR validated            |
| NCS           | $C_t > 40$   | If endogenous IPC (without exogenous IPC added): $C_t > 40$                                   | Extraction validated     |
|               |  | If exogenous IPC added: $C_t = C_{t\text{oc IPC of 5}} - \text{IPC West Nile} \pm 3C_t^{(3)}$ |                          |
| NC            | $C_t > 40$   | $C_t > 40$  | PCR components validated |

<sup>(1)</sup> Refer to the values listed in section 2.1 "EPC" of the Certificate of Analysis of the lot used for the test.

<sup>(2)</sup> The IPC value in the EPC should not be used for test validation.

<sup>(3)</sup> Refer to the values listed in section 2.2 "IPC" of the Certificate of Analysis of the lot used for the test.

### Interpretation of results

For each sample analyzed, the results should be interpreted as shown below:

#### For non-cellular and feces samples:

| West Nile detector | West Nile IPC detector (exogenous IPC)          | Interpretation               |
|--------------------|---|------------------------------|
| $C_t < 40$         | $C_t < 40$ or $C_t > 40$                        | West Nile virus detected     |
| $C_t > 40$         | $C_t \leq C_{t\text{ IPC of NCS}} + 3C_t^{(1)}$ | West Nile virus not detected |
| $C_t > 40$         | $C_t > C_{t\text{ IPC of NCS}} + 3C_t^{(1)}$    | Not validated <sup>(2)</sup> |

<sup>(1)</sup> Refer to the IPC  $C_t$  value obtained for the NCS done during the same extraction series as the samples to be analyzed. The IPC  $C_t$  value obtained for this NCS must first be validated as described above.

<sup>(2)</sup> The sample will be returned as not validated due to the negative IPC.

#### For cellular samples:

| West Nile detector | West Nile IPC detector (endogenous IPC) | Interpretation               |
|--------------------|---|------------------------------|
| $C_t < 40$         | $C_t < 40$ or $C_t > 40$                | West Nile virus detected     |
| $C_t > 40$         | $C_t < 40$                              | West Nile virus not detected |
| $C_t > 40$         | $C_t > 40$                              | Not validated <sup>(1)</sup> |

<sup>(1)</sup> The sample will be returned as not validated due to the negative IPC.

## Procedure for handling non-validated samples

1. Dilute the RNA at a 1:10 dilution in 1X TE buffer.
2. Denature the diluted RNA.
3. Perform a new RT-PCR analysis on 5 µL of this dilution (after denaturation).
4. If the diluted RNA is positive for West Nile with an acceptable IPC result, the result obtained is then validated.
5. If the diluted RNA is negative for West Nile with a non-compliant IPC result, the obtained result is still not validated. In this case, repeat the nucleic acid extraction using the sample pre-diluted 1:10 in 1X PBS buffer before extraction.
6. If the result is still not validated, repeat the analysis on a new sample.

## Documentation and support

### Customer and technical support

Technical support: visit [thermofisher.com/askaquestion](http://thermofisher.com/askaquestion)

Visit [thermofisher.com/support](http://thermofisher.com/support) for the latest in services and support, including:

- Worldwide contact telephone numbers
- Order and web support
- User guides, manuals, and protocols
- Certificates of Analysis
- Safety Data Sheets (SDSs; also known as MSDSs)  
**NOTE:** For SDSs for reagents and chemicals from other manufacturers, contact the manufacturer.

### Limited product warranty

Life Technologies Corporation and/or its affiliate(s) warrant their products as set forth in the Life Technologies' General Terms and Conditions of Sale found on Life Technologies' website at [www.thermofisher.com/us/en/home/global/terms-and-conditions.html](http://www.thermofisher.com/us/en/home/global/terms-and-conditions.html). If you have any questions, please contact Life Technologies at [thermofisher.com/support](http://thermofisher.com/support).



Laboratoire Service International (LSI) | 6 Allée des Écureuils | Parc Tertiaire du Bois-Dieu | 69380 Lissieu, France

Translated from the French Pub. No. MAN0008914 Rev. B.0.

The information in this guide is subject to change without notice.

**DISCLAIMER:** TO THE EXTENT ALLOWED BY LAW, LIFE TECHNOLOGIES AND/OR ITS AFFILIATE(S) WILL NOT BE LIABLE FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE, OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING YOUR USE OF IT.

Revision history of Pub. No. MAN0008915 (English)

| Revision | Date          | Description   |
|----------|---------------|---|
| B.0      | 28 June 2017  | Updated to the current document template, with associated updates to the warranty, trademarks, and logos. |
| A.0      | 31 March 2014 | Baseline for revision history   |

**Limited Use Label License No. 460: PCR Veterinary Diagnostics.** Notice to Purchaser: For veterinary diagnostic services including the reporting of such services for a fee and for research purposes only. Human diagnostic uses require a separate license from Roche.

©2017 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. TaqMan is a registered trademark of Roche, used under permission and license.

# Certificate of Analysis


Certificat d'Analyse

**applied biosystems**  
by Thermo Fisher Scientific

## VetMAX™ West Nile Virus Kit

|   |          |   |            |
|---|----------|---|------------|
|  | 50 Tests |  | WNPEXO-027 |
|  | WNPEXO50 |  | 2025-03-07 |
|   |          |  | 2026-03-05 |

### KIT CONFIGURATION CONFIGURATION DU KIT

| COMPONENT DESCRIPTION<br><i>Description du composant</i> | REF   | UNIT      | LOT       |  |
|--|-------|-----------|-----------|---|
| 3 - Mix West Nile  | MPEWN | 2 x 500µL | MPEWN-027 | 2026-03-05  |
| 4a - EPC West Nile                                       | EPCWN | 90 µL     | EPCWN-027 | 2026-03-07  |
| 5 - IPC West Nile  | IPCWN | 250 µL    | IPCWN-027 | 2026-03-07  |

### INSTRUCTIONS FOR USE *Notices d'utilisation*



French MAN0008914 RevB.0  
English MAN0008915 RevB.0

### QUALITY CONTROL RESULTS RESULTATS DE CONTROLE QUALITE

| EPC <i>EPC</i>     | Already extracted EPC <i>EPC déjà extrait</i> | Thermocycler used for QC<br><i>Thermocycleur utilisé pour le CQ</i> |
|--------------------|---|---|
| 4a - EPC West Nile | C <sub>t QC</sub> WN = 26                     | ABI 7500 (Applied Biosystems)                                       |

| IPC <i>IPC</i>                        | Column extraction<br><i>Extraction colonne</i> | Thermocycler used for QC<br><i>Thermocycleur utilisé pour le CQ</i> |
|---------------------------------------|--|---|
| 5 - IPC West Nile                     | C <sub>t QC</sub> IPC = 28                     | ABI 7500 (Applied Biosystems)                                       |
| Endogenous IPC<br><i>IPC endogène</i> | C <sub>t QC</sub> IPC < 40                     | ABI 7500 (Applied Biosystems)                                       |


C<sub>t QC</sub>: Indicative Ct of the Quality Control *C<sub>t</sub> indicatif du Contrôle Qualité*

### SUPPLEMENTAL INSTRUCTIONS - Available upon request to Eurotech@thermofisher.com

*Notices d'utilisation complémentaires - Disponibles sur demande auprès de Eurotech@thermofisher.com*

| Amplification Protocol <i>Protocole d'amplification</i> | Nucleic acid purification Protocol <i>Protocole de purification d'acides nucléiques</i> |
|---|---|
| N/A   | French MAN0008916 RevB.0  |

 Laboratoire Service International  
6 allée des Ecureuils  
Parc Tertiaire du Bois-Dieu  
69380 Lissieu - France  
[www.thermofisher.com/support](http://www.thermofisher.com/support)

  
Quality Assurance

Date: 24 march 2025.