Visi-Pro™

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Instructions for Use

Device description

The Visi-Pro™ balloon-expandable peripheral stent is intended as a permanently implanted device. The stent is made from a 316L stainless steel tube that is cut into an open lattice design and mounted onto a non-compliant balloon catheter. The Visi-Pro stent is deployed and expanded by inflating the balloon.

Indications for use

Peripheral

The Visi-Pro balloon-expandable peripheral stent is indicated for use in occlusions, lesions at high risk for abrupt closure or threatened closure following percutaneous transluminal angioplasty (PTA) or lesions believed to be at high risk of stenosis following PTA in the following arteries:

- common and external iliac, subclavian
- renal

Stenting is intended to improve and maintain artery luminal diameter.

Biliary

The Visi-Pro balloon-expandable peripheral stent is indicated for use as a palliative treatment of malignant neoplasms in the biliary tree.

Contraindications

Peripheral

The Visi-Pro stent is contraindicated for patients who share the following symptoms or factors:

- Patients for whom antiplatelet therapy, anticoagulants, or thrombolytic drugs are contraindicated.
- Patients who exhibit persistent acute intraluminal thrombus of the targeted lesion site, post thrombolytic therapy.
- Perforation at the angioplasty site evidenced by extravasation of contrast media.
- Aneurysm of the artery to be treated.
- All of the customary contraindications for PTA.

Biliary

- There are no known contraindications for the Visi-Pro stent in biliary interventions.

Warnings

- This device is supplied **sterile** for single use only. Do not reprocess or resterilize. Reprocessing and resterilizing could increase the risks of patient infection and compromised device performance.
- Select a balloon inflation diameter that approximates the diameter of the lumen at the intended implant site.
- The balloon inflation pressure should not exceed the maximum recommended inflation pressure (refer to the product labeling). Using an inflation device with a manometer is recommended.
- Do not retract the balloon catheter unless the balloon is fully deflated under vacuum.
- If resistance is encountered at any time during the insertion procedure, do not force passage. Forced passage may damage the stent system or the vessel lumen.
- If resistance is encountered during movement through the introducer sheath, carefully withdraw the stent system.
- Do not retract the stent into the sheath. Attempting to retract the stent into the sheath may result in stent dislodgment. If resistance is encountered after the stent has exited the introducer sheath, or if the stent cannot be delivered to the appropriate target lesion, withdraw the introducer sheath and the stent system as described in **Stent withdrawal**.
- Proceed carefully when recrossing a newly implanted stent with other devices. Recrossing may damage the stent or cause stent dislodgment.
- Implanting the stent across a side branch may result in occlusion of the artery and prevent further access for future transluminal interventions.
- The stent may cause artifacts with MRI scans due to distortion of the magnetic field. The artifacts caused by the 316L stainless steel stent should not be greater than artifacts caused by metal surgical clips. An MRI scan should not be performed until the implanted stent has had time to heal (estimated to be 8 weeks), to minimize the risk of stent migration under a strong magnetic field.

Precautions

- Carefully inspect the sterile package and the stent before use. Do not use the device if the packaging or the device is damaged.
- Pre-dilation may be required in cases of high grade stenosis to ensure that there is adequate lumen diameter for the stent to pass through the lumen.
- Always use an introducer sheath during the implant procedure to protect the vessel (peripheral use), the biliary tract (biliary use), and the puncture site. Select the introducer sheath using sizing recommendations on the label.

Potential adverse events

Potential adverse events from the stent implantation procedure include the usual complications with any percutaneous procedure, PTA, or stent placement.

The following potential adverse events are associated with peripheral interventions:

- Artery perforation or rupture
- AV fistula
- Bleeding requiring transfusion
- Cardiovascular events
- Cerebrovascular events
- Contrast media reaction/renal failure
- Death
- Embolism
- Hematoma
- Infection at the puncture site
- Intraluminal thrombus
- Pseudoaneurysm

The following complications are associated specifically with PTA or stent placement for occlusive or stenotic disease:

- Allergic response in persons with known allergic reactions to 316L stainless
- Amputation
- Dissection or intimal flap
- Failure to deploy the stent
- Restenosis
- Stent misplacement
- Stent migration
- Stent thrombosis or occlusion

The following complications are specifically associated with renal artery stent placement:

- Nephrectomy
- Renal infarction

The following potential adverse events are associated specifically with biliary interventions:

- Infection secondary to stent contamination, leading to cholangitis, hemobilia, peritonitis, or abscess
- Rupture resulting from overstretching the duct

Directions for use

Warning: This device is supplied **sterile** for single use only. Do not reprocess or resterilize. Reprocessing and resterilizing could increase the risks of patient infection and compromised device performance.

Storage

Store the stent system away from sunlight at room temperature in a dry place.

Recommended items

Peripheral indication

Prepare the following items using sterile technique:

- 10 cc syringe filled with sterile saline
- Appropriately sized exchange guidewire (refer to product labeling for sizing information)
- Appropriately sized hemostatic introducer sheath (refer to product labeling for sizing information)
- Appropriately sized guide catheter (refer to product labeling for sizing information)
- Renal use: appropriately sized long curved sheath or renal guide catheter (refer to product labeling)
- Inflation system

Biliary indication

Prepare the following items using sterile technique:

- 10 cc syringe filled with sterile saline
- Appropriately sized exchange guidewire (refer to product labeling for sizing information)
- Appropriately sized introducer sheath (refer to product labeling for sizing information)
- Chiba needle
- Inflation device
- Drainage catheter
- Cannula

Preparation

1. Select a stent system with a nominal balloon diameter matching the reference lumen diameter of the artery (peripheral use) or the obstructed duct and intended implant site (biliary use).

Warning: Select a balloon inflation diameter that approximates the diameter of the lumen at the intended implant site.

2. Remove the stent system from the packaging. Visually inspect the stent to ensure that it is loaded between the proximal and distal balloon markers. Gently tug on the stent to ensure that it is crimped securely in place. If it is not securely attached, do not use the device.

Caution: Carefully inspect the sterile package and the stent before use. Do not use the device if the packaging or the device is damaged.

3. Flush the balloon catheter lumen with saline, then load the stent delivery system onto the guidewire.
4. Connect a stopcock to the catheter inflation port.
5. Connect a syringe to the stopcock, open the stopcock, and pull the syringe plunger back to induce a negative pressure to remove air from the balloon and lumen. Close the stopcock and remove the syringe.
6. Connect an angioplasty inflation system to the stopcock. Open the stopcock and slowly fill the inflation lumen with diluted contrast medium.

Directions for peripheral stent implantation

Warning: If resistance is encountered at any time during the insertion procedure, do not force passage. Forced passage may damage the stent system or the vessel lumen.

Caution: Pre-dilation may be required in cases of high grade stenosis to ensure that there is adequate lumen diameter for the stent to pass through the lumen.

1. If the iliac or renal artery is treated, the access site is the common femoral artery. If the subclavian artery is treated, the access site is the brachial or axillary artery. Cross the lesion with the guidewire and advance the guide catheter or introducer sheath atraumatically over the guidewire across the lesion.
Caution: Always use an introducer sheath during the implant procedure to protect the vessel and the puncture site. Select the introducer sheath using sizing recommendations on the label.
2. Advance the stent system through the introducer sheath or the guide catheter to the lesion site.
Warning: If resistance is encountered during movement through the introducer sheath, carefully withdraw the stent system.
Warning: Do not retract the stent system into the sheath. Attempting to retract the stent system into the sheath may result in stent dislodgment. If resistance is encountered after the stent has exited the introducer sheath, or if the stent cannot be delivered to the appropriate target lesion, withdraw the introducer sheath and the stent system as described in **Stent withdrawal**.
3. Position the stent at the lesion site, using contrast medium to enable fluoroscopic visualization. For ostial lesions, ensure that the ostium of the lesion is covered by positioning the stent to protrude approximately 1 to 2 mm into the aorta. For iliac stenting, the entire stent should be placed within the iliac artery.
4. Under fluoroscopy, expand the stent by inflating the balloon until the stent is fully expanded to the desired diameter.
Warning: The balloon inflation pressure should not exceed the maximum recommended inflation pressure (refer to the product labeling). Using an inflation device with a manometer is recommended.
5. Deflate the balloon.
6. Gently rotate the deflated balloon delivery catheter to ensure that the stent is free and properly deployed.
7. Remove the delivery catheter.
Warning: Do not retract the balloon catheter unless the balloon is fully deflated under vacuum.
8. Using angiography, determine whether the stent is appropriately expanded and positioned. If the stent does not cover the lesion, implant a second stent. Place the second stent overlapping the initially implanted stent with approximately 1 to 2 mm of overlap.
Warning: Proceed carefully when recrossing a newly implanted stent with other devices. Recrossing may damage the stent or cause stent dislodgment.

Stent withdrawal

If resistance is encountered after the stent has exited the introducer sheath or if the stent cannot be delivered to the appropriate target lesion, withdraw the introducer sheath and the stent and delivery assembly. Use the following steps for withdrawal.

1. Using fluoroscopic guidance, retract the stent to the end of the introducer sheath.
2. Inflate the delivery balloon 1 atm to reduce the likelihood of stent slippage or embolization.
3. Retract the introducer sheath and the stent system while preserving the guidewire position.
4. Withdraw the introducer sheath and the stent system assembly together as a unit.

Directions for biliary stent implantation

Warning: If resistance is encountered at any time during the insertion procedure, do not force passage. Forced passage may damage the stent system or the vessel lumen.

A guidewire is first introduced through the stricture, and an internal or external biliary drain is positioned through the stricture. Next, the transhepatic tract is accessed through the routine transhepatic approach (a lateral incision is placed between the tenth and eleventh intercostal space, then punctured upward under fluoroscopy aiming at the ninth or tenth vertebral body). No predilation is performed with malignant strictures, but a large lumen catheter may be inserted. Subsequent dilation generally is performed for only a few seconds.

Biliary stent implantation

1. Select an introducer sheath using sizing recommendations on the label. Use the introducer sheath to gain access, and advance the stent system through the introducer sheath and over the guidewire to the correct position.

Caution: Always use an introducer sheath during the implant procedure to protect both the liver tract and puncture site. Select the introducer sheath using sizing recommendations on the label.

Warning: If resistance is encountered during movement through the sheath, carefully withdraw the stent system.

Warning: Do not retract the stent into the sheath. Attempting to retract the stent into the sheath may result in stent dislodgment. If resistance is encountered after the stent has exited the introducer sheath, or if the stent cannot be delivered to the appropriate target lesion, withdraw the introducer sheath and the stent system as described in **Stent withdrawal**.

2. Position the stent across the target stricture. The distal portion of the stent should be positioned so that it completely covers the distal end of the stricture.

Note: Take care to position the stent so that it covers areas that are proximal and distal to an inoperable malignant tumor. Fully covering proximal and distal margins protects against obstruction from additional tumor growth.

3. Pull back on the sheath to expose the stent completely. Expand the stent by inflating the balloon catheter to the nominal inflation pressure. Use fluoroscopy to determine whether the balloon and stent are expanded appropriately. Do not exceed the recommended maximum balloon inflation pressure.

- Warning:** The balloon inflation pressure should not exceed the maximum recommended inflation pressure (refer to the product labeling). Using an inflation device with a manometer is recommended.
4. Expand the stent to the lumen diameter of the bile duct.
Warning: The balloon inflation pressure should not exceed the maximum recommended inflation pressure (refer to the product labeling). Using an inflation device with a manometer is recommended.
 5. Deflate the balloon.
 6. Ensure that the stent is free and deployed properly by rotating the deflated balloon delivery catheter.
Warning: Do not retract the balloon catheter unless the balloon is fully deflated under vacuum.
 7. Remove the delivery catheter.
Warning: Proceed carefully when recrossing a newly implanted stent with other devices. Recrossing may damage the stent or cause stent dislodgment.

Warranty

This product is subject to standard Medtronic warranty terms.

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© 2017 Medtronic. All rights reserved. Medtronic and Medtronic logo are trademarks of Medtronic. [™] Third-party brands are trademarks of their respective owners. All other brands are trademarks of a Medtronic company.
ev3 Inc.
Medtronic Australasia Pty Ltd