Core tests® HIV1+2/Syphilis Test Cassette

For Professional Use

Specimen: Whole Blood/Serum/Plasma

Format: Cassette

INTENDED USE

HIV1+2/Syphilis Test Cassette is a rapid chromatographic immunoassay for the qualitative detection of antibodies to HIV1+2/Syphilis to aid in the diagnosis of HIV1+2/Syphilis.

SUMMARY

Treponema Pallidum (Syphilis) is the causative agent of the venereal disease Syphilis. TP is a spirochete bacterium with an outer envelope and a cytoplasmic membrane. Relatively little is known about the organism in comparison with other bacterial pathogens. According to the Center for Disease Control (CDC), the number of cases of Syphilis infection has markedly increased since 1985. Some key factors that have contributed to this rise include the crack cocaine epidemic and the high incidence of prostitution among drugs users. One study reported a large number of HIV-infected females treatment.

HIV is the etiologic agent of Acquired Immune Deficiency Syndrome (AIDS). The virus is surrounded by a lipid envelope that is derived from host cell membrane. Several viral glycoproteins are on the envelope. Each virus contains two copies of positive-sense genomic RNAs. Patients with HIV-1 has been isolated from patients with AIDS and AIDS-related complex and from healthy individuals with a high potential risk of developing AIDS. Patients with HIV-2 has been isolated from West Africa AIDS patients and from seropositive asymptomatic individuals. Both HIV-1 and HIV-2 elicit an immune response. Detection of HIV antibodies in whole blood, serum or plasma is the most efficient and common way to determine whether an individual has been exposed to HIV and to screen blood and blood products with HIV. Despite the difference in their biological characteristics, serological activities and genome sequences of HIV-1 and -2 show strong antigenic cross-reactivity. Most HIV-2 positive sera can be identified by using HIV-1 based serological tests.

HIV1+2/Syphilis Test Cassette is a rapid test to qualitatively detect antibodies to HIV-1 and/or HIV-2 in whole blood, serum or plasma specimens. The test utilizes a combination of protein A-coated gold conjugate and multiple recombinant HIV proteins to selectively detect antibodies to the HIV-1 and HIV-2 in whole blood, serum or plasma. The genes for envelope proteins (gp36/41) encode the recombinant HIV proteins that are used in the test kits.

PRINCIPLE

HIV1+2/Syphilis Test Cassette contains a membrane strip, which is pre-coated with recombinant HIV 1+2 capture antigen (gp41 and gp36) in test line 1 region and with recombinant TP capture antigen in test line 2 region respectively. The recombinant HIV 1+2/TP antigen - colloidal gold conjugate and blood specimen move along the membrane chromatographically to the test region (T), forming a visible line as the antigen-antibody-antigen gold particle complex forms with high degree of sensitivity and specificity.

This test device has letters of 1, 2 and C as "Test Line1 (HIV)", "Test Line2 (TP)" and "Control Line" on the surface of the cassette. Both the Test Lines and Control Line in result window are not visible before applying any specimen. The Control Line is used for procedural control. Control line should always appear if the test procedure is performed properly and the test reagents of control line are working.

MATERIALS PROVIDE

HIV1+2/Syphilis test kit contains following items to perform the assay:

1. HIV1+2/Syphilis test device
2. Instruction for use
3. Chasing buffer
4. Sterile lancet
5. Pipette
6. Alcohol pad

PRECAUTIONS

1. Read instruction for use carefully before performing this test.
2. For in vitro diagnostic use only.
3. Do not use the test device beyond the expiration date.
4. The test device should remain in the sealed pouch until use. Do not use the test device if the pouch is damaged or the seal is broken.
5. Do not reuse the device.
6. Treat and properly handle the specimens and used device as if they were potentially infectious.

STORAGE AND STABILITY

The test device should be stored at 2-30 °C in the sealed pouch. Avoid humidity, heat and direct sunlight. The test device is stable through the expiration date printed on the sealed pouch. DO NOT FREEZE.

WARNINGS

1. There should be no eating, drinking or smoking where specimens are being handled.
2. Wear disposable gloves and lab coat while handling specimens. Wash hands thoroughly afterwards.
3. Avoid splashing or aerosol formation.
4. Clean spills thoroughly using an appropriate disinfectant.

SPECIMEN COLLECTION AND STORAGE

Whole Blood Specimens

1. Clean the area to be lanced with alcohol pad.
 2. Squeeze the fingertip and pierce it with the sterile lancet.
 3. Wipe away the first drop of blood with sterile gauze or cotton.
 4. Using a disposable pipette, collect blood from the puncture site.
- The whole blood may be used for testing immediately or may be refrigerated at 2- 8 °C up to three days. Such specimen should be brought and equilibrated to room temperature prior to use.

Serum or Plasma Specimens


1. Centrifuge whole blood to get plasma or serum specimen.
2. If serum is not tested immediately, it should be refrigerated at 2-8 °C. For storage period greater than three days, freezing is recommended. Such specimen should be brought and equilibrated to room temperature (15-30 °C) prior to use.
3. Serum containing precipitate may yield inconsistent test results. Such specimens must be clarified prior to assaying.

TEST PROCEDURE

1. Bring the pouched test device to room temperature (15-30 °C) prior to testing. Do not open the pouch until ready to perform the assay.
2. Remove the test device from the sealed pouch. Lay it on a flat, clean and dry surface.
3. Use the pipette to draw and slowly add 1 drop of whole blood/serum/plasma to the sample well.
4. Hold the buffer bottle vertically and add 1 drop to the sample well. /If using a pipette, change a new one to avoid cross-contamination.
5. Draw and transfer 2-3 drops of buffer to the sample well.
5. Read test results between 10-20 minutes. Do not interpret after 20 minutes.

Caution: The above interpreting time is based on room temperature range of 15 - 30 °C. If your room temperature is significantly lower than 15 °C, then the interpreting time should be properly increased to 30 minutes.

INTERPRETATION OF RESULTS

HIV Positive: The control line and HIV line (T1) are visible in the result window. The test is positive for HIV and antibodies to HIV have been detected.

TP Positive: The control line and TP line (T2) are visible in the result window. The test is positive for TP and antibodies to TP have been detected.

HIV and TP Positive: The control line, HIV line (T1) and TP line (T2) are visible in the result window. The test is positive for HIV/TP.

Negative: Only one line in the control region appears in the result window. The test result is negative and no antibodies to HIV or TP have been detected.

Invalid: If the control line does not appear in the result window, the test results are invalid regardless of the presence or absence of the test line(s).

NOTE: Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for the control line failure. Review the procedure and repeat the test with a new device. If problem persists, please contact your local distributor.

PERFORMANCE CHARACTERISTICS

1) Sensitivity and Specificity

HIV1+2/Syphilis Test Cassette have been tested with positive and negative clinical specimens tested by a leading commercial anti-HIV 1+2 ELISA test. The result shows that the HIV 1+2 test device is very accurate comparing with other commercial ELISA kit.

1.1 HIV 1/2

Reference		HIV 1+2		Total
Method	Result	Positiv	Negativ	Results
Commercial ELISA	Positive	819	1	820
	Negativ	1	511	512
Total Results		820	512	1332

In a comparison of our HIV 1+2 P24 Ag test cassette with a leading commercial anti-HIV 1+2 ELISA test, result gave sensitivity of 99.9% (819/820), and a specificity of 99.8% (511/512), and a total 1.2P24Ab agreement of 99.8% (1330/1332).

1.2 Syphilis

Reference		HIV 1+2		Total
Method	Result	Positiv	Negativ	Results
Commercial ELISA	Positive	145	1	146
	Negativ	5	527	532
Total Results		150	528	678

In a comparison of our Syphilis test device with a leading commercial Syphilis ELISA test, the result gave a sensitivity of 99.3% (145/146), a specificity of 99.1% (527/532), and a total agreement of 99.1% (672/678).

2. Precision

A. Within run precision was determined by using 10 replicates of four different specimens containing different concentrations of antibody. The negative and positive values were correctly identified 100% of the cases.

B. Between run precision was determined by using four different specimens containing different concentrations of antibody in 3 different replicates with 3 different lots of test devices. Again negative and positive results were observed 100% of the cases.

LIMITATION OF THE TEST

1. HIV1+2/Syphilis Test Cassette is for in vitro diagnostic use only. The test should be used for the detection of HIV1+2/TP antibodies in whole blood, serum or plasma specimens. Neither the quantitative value nor the rate of increase in HIV1+2/TP antibodies can be determined by this qualitative test.

2. HIV1+2/Syphilis Test Cassette will only indicate the presence of HIV1+2/TP antibodies in the specimen and should not be used as the sole criteria for the diagnosis of HIV1+2/TP infection.

3. As with all diagnostic tests, all results must be interpreted together with other clinical information available to the physician.

4. If the test result is negative and clinical symptom persists, additional testing using other clinical methods is recommended. A negative result does not at any time preclude the possibility of HIV1+2/TP infection.

INDEX OF SYMBOLS

	Do not re-use		Use-by date
	Consult instructions for use		Batch code
	Manufacturer		Store at 2-30°C
	In vitro diagnostic medical device		

Core Technology Co.,Ltd.
Room 100, C Building, No.29 Life Park Rd.,
Changping District, Beijing 102206, P.R. China
Tel: +86-10-69390616
Email: info@coretests.com
Fax: +86-10-69390660

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