### Tissue-Tek® TEC<sup>™</sup> 6 Embedding Console System

The most trusted name in reliable, comfortable and easy-to-use embedding









# The embedding system you trust with all of the enhancements you need

The Tissue-Tek® TEC™ Series has been the Histology laboratory's preferred choice in embedding for decades. The Tissue-Tek® TEC™ 6 Embedding Console System continues that tradition while further enhancing the reliability, comfort, and ease-of-use that users have come to know and expect. Designed to be the ideal tissue embedder for laboratories of any size, the robust and ergonomic modular system offers a streamlined, adjustable workflow that is both comfortable and simple-to-use for any user.



# Tissue-Tek® TEC™ 6 Full embedding control at your fingertips

- Colour touchscreen with easy-to-navigate onscreen controls allows easy setup and operation
- Graphical user interface enables users to set and display temperatures in all heated and cooled areas with just a touch
- Multiple levels of brightness for the LED work light can be adjusted directly from the main screen
- Automatic ready and shutdown times for both the Embedding and Cryo Modules are user-programmable by day of the week and even across overnight shifts

#### Safe and reliable embedding with ease

- Anti-glare, coated work surfaces reduce eye strain and provide better visual identification of tissues
- Adjustable LED lighting improves specimen visibility while eliminating workplace clutter from external lighting or lamps
- The software-driven Clean Mode provides an easy method for cleaning the forceps wells, reducing opportunities for cross contamination of tissue blocks
- The optional Tissue-Tek® TEC™ 6 Magnifying Lens allows users to keep a close eye on even the smallest biopsies to ensure proper orientation

# Modular console system with the flexibility and comfort you need

- Ergonomic, padded wrist rests provide support and comfort when needed or can be easily removed if desired
- Modules can be arranged either left-to-right or right-to-left based on user preference and comfort
- Large warming surface offers an extended working area that is wide enough to accommodate up to 10 cassettes on either side of the user
- The large, nylon surface of the Cryo Module has raised edges to contain moisture and provides ample space for up to 60 embedded blocks







# **Specifications**

#### **General instrument specifications**

| Product name                       | Tissue-Tek® TEC™ 6 Embedding Console System  |
|------------------------------------|--|
| Product code                       | 5114   |
| Applications                       | Embed processed tissue in paraffin, forming paraffin blocks for sectioning   |
| Capacity                           | Paraffin chamber: 4 litres Left and right thermal chambers: 1.5 litres or up to 150 cassettes Working area: up to 20 cassettes Cryo plate: up to 60 blocks |
| Paraffin delivery                  | Gravity feed with touch-plate activation and optional foot pedal control   |
| Paraffin flow rate                 | Flow rate adjusted with control knob   |
| Waste collection                   | Drain channels in the heated working area route paraffin waste to two front waste drawers  |
| Illumination                       | 5-level, adjustable white LED  |
| Automatic timer mode               | User-programmable, ready with shutdown times for each day of the week  |
| Heating time<br>(Embedding Module) | Paraffin and left/right chambers: 4 hours from room temperature start  |
| Cooling time (Cryo Module)         | 30 minutes from room temperature start   |
| Temperature range                  | Paraffin and thermal chambers: 50 to 75°C  Hot plate: 50 to 75°C  Cryo plate: -10 to 0°C   |
| User interface display             | 4.3 inch TFT colour LCD touchscreen  |
| Diagnostic function                | Self-diagnosis, error messages, power-out detection  |
| Data interface options             | USB  |
| Operating temperature              | 10 to 35°C   |
| Operational relative humidity      | 30% to 85% (non-condensing)  |
| Rated voltage and current          | Embedding Module: 230 VAC $\pm$ 10%, 50 Hz, 4 A Cryo Module: 230 VAC $\pm$ 10%, 50 Hz, 1 A   |
| Dimensions                         | Embedding Module: 57.5 (W) x 64.2 (D) x 37.7 (H) cm<br>Cryo Module: 33.0 (W) x 61.7 (D) x 37.7 (H) cm<br>Combined: 90.5 (W) x 64.2 (D) x 37.7 (H) cm       |
| Weight                             | Embedding Module: 25 kg<br>Cryo Module: 22 kg<br>Combined: 47 kg   |
| Applicable Standards               | CE marked, EN 61010-1: 2010, EN 61326-1: 2013, EN 50581: 2012  |
|                                    |  |



#### **Accessories and consumables**

| Product code | Product name and quantity  |
|--------------|--|
| 1550         | Tissue-Tek® Scraper, Large   |
| 1551         | Tissue-Tek® TEC™ Large Tamper; 3/case  |
| 1552         | Tissue-Tek® TEC™ Small Tamper; 3/case  |
| 5120         | Tissue-Tek® TEC™ Warming Chamber Tray, Large                                 |
| 5121         | Tissue-Tek® TEC™ 6 Rest Wrist (right)  |
| 5122         | Tissue-Tek® TEC™ 6 Rest Wrist (left)   |
| 5123         | Tissue-Tek® TEC™ 6 Cushion (right)   |
| 5124         | Tissue-Tek® TEC™ 6 Cushion (left)  |
| 5125         | Tissue-Tek® TEC™ 6 Screen Protection Sheet                                   |
| 5127         | Tissue-Tek® TEC™ 6 Magnifying Lens   |
| 5785         | Tissue-Tek® TEC™ Foot Pedal Switch   |
| 4653         | Tissue-Tek® Paraffin Wax Tek II, 1x 10kg                                     |
| 7052         | Tissue-Tek® Paraform® Processing/Embedding<br>Medium Formula 3, 1 kg; 8/case |

