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Tab. 7

Table of Contents

1	Explanation of Symbols	47	8	Handling the Coolant/Lubricant (C/L)	76
2	General Safety Instructions	48	8.1	C/L Analysis	76
3	Suitable Personnel	48	8.2	Removing the C/L Container	77
4	Machine Specifications	49	8.3	Emptying the C/L Container	78
4.1	Ceramill Motion 2	49	8.4	Cleaning the C/L Container	79
4.2	Ceramill Coolstream	54	8.5	Filling in C/L	79
5	Installation	57	8.6	Installing the C/L Container	80
5.1	Ceramill Motion 2	57	8.7	Cleaning the C/L System	80
5.2	Service Unit	62	9	Cleaning and Maintenance	81
5.3	Installation of the Ceramill Match 2 and Ceramill Motion 2 Software ..	63	9.1	Ceramill Motion 2	81
5.4	Establishing the Network Connection	65	9.2	Ceramill Coolstream	85
5.5	Starting-up the Ceramill Motion 2 ..	66	10	Malfunctions, Repairs and Warranty ..	86
6	Application and Operation	67	10.1	Malfunctions	86
6.1	Switching from Wet to Dry Processing	68	10.2	Repairs	86
6.2	Changing the Blank Holder	70	10.3	Warranty	86
6.3	Checking the Model Area	70	11	Environmental Protection	86
6.4	Function of the Buttons	71	12	Technical Data and Accessories/Spare Parts	87
6.5	Reference Run	71			
6.6	Program Execution	72			
6.7	Checking the Workpiece During a Running Program	72			
6.8	Moving to the Service Position ...	72			
6.9	Inspection-window Cleaning	73			
6.10	Advancing to the Change Position for the Blank Holder	73			
7	Registration and Download Information	74			
7.1	Quick Registration, M-center	74			
7.2	Order Tracking, M-Center	75			
7.3	Download Infos	75			

1 Explanation of Symbols

Warning indications

 Warning indications in the text are marked with a colour-backed triangle and boxed.

 In case of hazards through electricity, the exclamation mark in the warning triangle is substituted by a lightning bolt.

Signal words at the beginning of a warning indication specify the type and severity of the consequences, if the measures to avert the hazard are not adhered to.

- _ **NOTE** means that property damage can occur.
- _ **CAUTION** means that light to fairly serious personal injury can occur.
- _ **WARNING** means that serious personal injury can occur.
- _ **DANGER** means that serious personal injury can occur.

Important information

 Important information that do not lead to hazards for humans or property damage are marked with the icon aside and are also boxed.

Other symbols in the Manual

Symbol	Meaning
>	Item of an operation description
—	Item of a list
•	Subitem of an operation description or a list
[3]	Numbers in square brackets refer to position numbers in graphics/figures

Tab. 1

Other symbols on the machine

Symbol	Meaning
	USB port (socket type B)
	Network connection (Ethernet)
	Slot for SD card
	Control cable connection, extraction
	Fuse
	Compressed-air connection
	Control cable connection, Cool-stream
	Connection for coolant/lubricant supply line

Tab. 2

2 General Safety Instructions

When installing, starting-up and operating the machine, always observe the following safety instructions:

CAUTION:

A defective machine can lead to malfunctions! On detection of damage or a functional defect of the machine:

- ▷ Label the machine as defective.
- ▷ Prevent further operation until the machine has been repaired.

NOTE:

Possible machine damage through insufficient extraction during dry processing!

- ▷ Operate the machine only with Ceramill Airstream or another compatible extraction system (see page 63).

NOTE:

Possible machine damage through insufficient cooling during grinding operations!

- ▷ For grinding operations, operate the machine only with Ceramill Coolstream (see page 54ff).

NOTE:

Damage through escaping coolant/lubricant!

- ▷ Regularly check the coolant/lubricant connections for tightness against leaks.

NOTE:

- ▷ Switch the machine off when not in use or unsupervised for longer periods, e.g., overnight. This method of saving electrical energy also benefits the environment.

3 Suitable Personnel

NOTE:

Starting-up and operation of the machine may only be carried out by trained specialised personnel.

4 Machine Specifications

4.1 Ceramill Motion 2

The Ceramill Motion 2 is available as four-axis model (4X) and as five-axis model (5X). These Operating Instructions apply for both models.

4.1.1 Delivery Scope

- _ Ceramill Motion 2 milling and grinding machine
- _ Mains cable
- _ USB cable (5 m)
- _ Network cable (5 m)
- _ Control cable, extraction
- _ Compressed-air hoses (250 mm and 2000 mm)
- _ Cleaning hose with brush head
- _ Service unit
- _ Service set for Jäger spindle
- _ Ceramill-Match-Package
 - Ceramill Match 2 CD (Ceramill Match 2 and Ceramill Motion 2 operating software)
 - User Manual
- _ Blank template, ZI71 Motion
- _ Ceramill TEST 71 L (test blank)
- _ Ceramill Motion Roto 2.5 (milling tool mounted in collet)
- _ Screw set
- _ Torque wrench, 4 Nm
- _ Blank holder 71
- _ Extraction adapter
- _ Extraction hood
- _ Drain cover
- _ Caliper

The Ceramill Motion is delivered with a Ceramill TEST 71 L test blank (Order No. 760301) inserted in the blank holder. Prior to delivery, a test object was factory-milled in order to check the zero point settings. Possibly remaining milled shavings/chips are not a reason for rejection. The test blank is provided for test milling and calibration of the Ceramill Motion.

- ▷ After unpacking, check the machine for completeness and possible transport damages. Please claim any transport damages immediately with your supplier (compare with Ceramill Mall setup checklist).

4.1.2 Intended Use

The Ceramill Motion 2 is a PC-controlled milling and grinding machine for the fabrication of dentures by means of:

- _ Dry processing of blanks made of
 - pre-sintered zirconium oxide
 - unsintered non-precious metals (NPMs)
 - wax
 - plastics
- _ Wet processing of blanks made of
 - sintered glass ceramic
 - plastics

Using blanks and tools not approved by AmannGirrbach can damage the machine and make the product unusable. For such cases, AmannGirrbach shall assume no liability whatsoever.

Unauthorised modifications/alterations of the machine shall also void the warranty.

A compressed-air connection is required for operation of the Ceramill Motion 2 (see page 62).

Dry processing is permitted only in conjunction with an appropriate extraction system (for technical requirements, please refer to page 63). We recommend using Ceramill Airstream.

Wet milling and grinding is permitted only in conjunction with Ceramill Coolstream (accessory).

4.1.3 CE Declaration of Conformity

In terms of design and performance, this product complies with the European Directives and the supplementary national requirements. Conformity has been confirmed with the CE marking.

The product's Declaration of Conformity can be viewed under www.amanngirrbach.com.

4.1.4 Components and Interfaces


Fig. 1 Machine overview, front, five-axis model

- 1 Measuring device, tool length
- 2 C/L nozzles
- 3 Milling spindle
- 4 Extraction hood
- 5 Service button
- 6 Operation button
- 7 C/L drain
- 8 Tool holder
- 9 Blank holder


Fig. 2 Machine overview, front, four-axis model

- 1 Measuring device, tool length
- 2 C/L nozzles
- 3 Milling spindle
- 4 Extraction hood
- 5 Service button
- 6 Operation button
- 7 C/L drain
- 8 Tool holder
- 9 Blank holder


Fig. 3 Machine overview, rear

- 1 USB data port
- 2 Ethernet interface
- 3 Control interface for Ceramill Coolstream
- 4 Slot for SD memory card
- 5 Control interface for extraction
- 6 Compressed-air inlet
- 7 Compressed-air hose, 250 mm
- 8 Service unit
- 9 Compressed-air connection, milling spindle
- 10 C/L supply line connection
- 11 C/L drain connection
- 12 Fuse compartment
- 13 Power supply connection
- 14 Main switch
- 15 Extraction-hose connection

4.1.5 Accessories

Ceramill Coolstream

The Ceramill Coolstream supply unit (see section 4.2) is mandatory for wet processing. It provides for the necessary coolant/lubricant circuit and includes the coolant/lubricant filtration. Additionally, it offers space for Ceramill Airstream.

Ceramill Airstream

Ceramill Airstream is an extraction system for dry processing with Ceramill Motion 2. It is equipped with a special filter bag and a HEPA micro filter, which filters up to 99.97 % of fine dust particles (i.e. zirconium oxide dust particles) and is classified in the fine dust particle category M.

The separate socket for the control cable allows for automatic operation.

The Ceramill Airstream extraction unit can be used with Ceramill Coolstream.

Milling and Grinding Tools

Various cutting and grinding tools are available for the Ceramill Motion 2 (please refer to the Catalog).

4.2 Ceramill Coolstream

4.2.1 Delivery Scope

- _ Ceramill Coolstream supply unit with
 - Container for coolant/lubricant (C/L)
 - C/L filter
- _ Control cable
- _ C/L supply line
- _ Piping kit
- _ Mains cable
- _ Swarf sieve
- _ Splash guard
- _ 2 litres of Ceramill coolant (coolant/lubricant)
- _ Analytic test strips
 - Multi-function additive test
 - pH-value
 - Nitrite test
 - Test record

4.2.2 Intended Use

Ceramill Coolstream is the supply unit for wet processing with Ceramill Motion 2. It ensures the application-efficient supply, drainage and filtering of the required coolant/lubricant.

Operation with other machines than Ceramill Motion 2 is not permitted.

4.2.3 CE Declaration of Conformity

In terms of design and performance, this product complies with the European Directives and the supplementary national requirements. Conformity has been confirmed with the CE marking.

The product's Declaration of Conformity can be viewed under www.amanngirrbach.com.

4.2.4 Components and Interfaces


Fig. 4 Machine overview, front

- 1** C/L drain for set-up under a table
- 2** C/L drain for set-up directly below Ceramill Motion 2
- 3** Retainer for adjustable leg of Ceramill Motion 2


Fig. 5 Machine overview, rear

- 1** Control interface for Ceramill Coolstream
- 2** Power supply connection
- 3** Control interface for extraction (i. e. for Ceramill Airstream)
- 4** C/L supply line connection
- 5** C/L drain for set-up directly below Ceramill Motion 2
- 6** C/L drain for set-up under a table
- 7** C/L supply line
- 8** C/L container
- 9** C/L drain cock

5 Installation

5.1 Ceramill Motion 2

Countries with 100-115 V mains voltage

The machine is preset for a mains voltage of 230 V. When the mains voltage at the set-up location is 100-115 V:

- ▷ Pull out the fuse compartment ([11] in Fig. 3 on page 53) below the mains plug.
- ▷ Replace both T3.15 A fuses with the provided T6.3 A fuses.

5.1.1 Ceramill Motion 2 Set-up without Ceramill Coolstream

Please refer to the separate Ceramill Motion 2 set-up instructions.

- The machine is intended exclusively for use within dry, closed rooms.
- A clearance of at least 200 mm is to be observed sideways, to the rear and upwards. Minimum space requirements for the Ceramill Motion 2 (incl. connections):
 - Width: 920 mm
 - Depth: 790 mm
 - Height: 955 mm
- The empty weight of the machine is 78 kg (five-axis model), resp. 75 kg (four-axis model). The set-up surface must withstand an appropriate load.
- During operation, the room temperature must be between 18 °C and 30 °C; avoid large variations in temperature.
- ▷ Set up the machine on a sturdy work table or work bench (off the ground and clear of walls).
- ▷ Adjust the height of the machine by screwing the feet in or out, ensuring that the machine is positioned firmly on all four feet.
- ▷ Insert the extraction hose of the extraction system being used (e.g. Ceramill Airstream) to the extraction connection [7] on the rear side of the Ceramill Motion 2.
- ▷ Connect the compressed-air hose on the rear side of the Ceramill Motion 2

5.1.2 Ceramill Motion 2 Set-up with Ceramill Coolstream

Please refer to the separate Ceramill Motion 2 and Ceramill Coolstream set-up instructions.

In conjunction with the Ceramill Coolstream supply unit, the Ceramill Motion 2 is set up either directly on the Ceramill Coolstream (tower version) or on a table above the Ceramill Coolstream (table version).

Installation of the Ceramill Airstream in the Ceramill Coolstream

When dry processing with air extraction via the Ceramill Airstream, the air extraction can be installed in the Ceramill Coolstream (please refer to the Ceramill Airstream set-up instructions). This simplifies the cabling.

- ▷ Set the Ceramill Airstream to the operating mode AUTO.
- ▷ Place the Ceramill Airstream with the extraction opening facing toward the rear into the Ceramill Coolstream.

Malfunction!

- ▷ Do **not** plug the control cable of the Ceramill Airstream into the Ceramill Motion 2!
- ▷ Plug the control cable of the Ceramill Airstream into the Ceramill Coolstream.
- ▷ Plug the extraction hose of the Ceramill Airstream into the Ceramill Motion 2.


Fig. 6

- 1 Ceramill Coolstream
- 2 Ceramill Airstream

Ceramill Motion 2 set-up directly on the Ceramill Coolstream (tower version)

A clearance of at least 200 mm is to be observed sideways, to the rear and upwards. Minimum space requirements for the Ceramill Motion 2 on the Ceramill Coolstream (incl. connections):

- _ Height: 1,620 mm
- _ Width: 920 mm
- _ Depth: 860 mm

The set-up location must be level.


Fig. 7

- 1 Ceramill Motion 2
- 2 Ceramill Coolstream

- ▷ Place the Ceramill Motion 2 onto the Ceramill Coolstream, ensuring that the adjustable legs engage into the retainers.
- ▷ Mount the service unit (please refer to page 62).
- ▷ Remove the C/L container of the Ceramill Coolstream (please refer to section 8.2 on page 77).
- ▷ Fill Ceramill Coolant into the C/L container (please refer to section 8.5 on page 79).
- ▷ To reinstall the C/L container (please refer to section 8.6 on page 80).
- ▷ Connect the C/L supply line to the Ceramill Motion 2 and to the Ceramill Coolstream.


Fig. 8

- 1 Ceramill Motion 2
- 2 Ceramill Coolstream
- 3 C/L drain pipe

- ▷ Plug the C/L drain pipe [3] together and connect it to the Ceramill Motion 2 and the Ceramill Coolstream (please refer to set-up instructions).
- ▷ Move the Ceramill Coolstream with the Ceramill Motion 2 to the set-up location.
- ▷ Lock the casters of the Ceramill Coolstream with the brakes.

NOTE:

Possible damage of the electric lines and the C/L supply line:

- ▷ When moving the Ceramill Coolstream, pay attention that no lines are driven over.

- Insert the extraction hose of the extraction system being used (e.g. Ceramill Airstream) to the extraction connection ([15] in Fig. 3 on page 53) on the rear side of the Ceramill Motion 2.

Ceramill Motion 2 set-up directly on a table above the Ceramill Coolstream (table version)

Under the table, a clearance of at least 200 mm is to be observed toward the sides and the rear, and at least 10 mm upwards. Minimum space requirements for the Ceramill Coolstream below the table:

- _ Height below the table: 675 mm
- _ Width below the table: 920 mm
- _ Depth below the table: 860 mm

The set-up location must be level.

For the set-up dimensions of the Ceramill Motion 2, please refer to page 57.


Fig. 9

- 1 Ceramill Motion 2
- 2 Ceramill Coolstream

In order for the C/L drain pipe to be mounted, the table must be set up with a certain clearance to the wall or have an appropriate cut-out for the piping of the C/L drain (please refer to set-up instructions).

- ▷ Position the Ceramill Motion 2 on the table.
- ▷ Adjust the height of the machine by screwing the feet in or out, ensuring that the machine is positioned firmly on all four feet.
- ▷ Mount the service unit (please refer to page 62).
- ▷ Remove the C/L container of the Ceramill Coolstream (please refer to section 8.2 on page 77).
- ▷ Fill Ceramill Coolant into the C/L container (please refer to section 8.5 on page 79).
- ▷ To reinstall the C/L container (please refer to section 8.6 on page 80).

NOTE:

Possible damage of the electric lines and the C/L supply line:

- ▷ When moving the Ceramill Coolstream, pay attention that no lines are driven over.

- ▷ Move the Ceramill Coolstream below the table in such a manner that it is aligned with the Ceramill Motion 2.


Fig. 10

- 1 Ceramill Motion 2
- 2 Ceramill Coolstream
- 3 Connection for C/L return pipe
- 4 C/L drain pipe

- ▷ Connect the C/L supply line to the Ceramill Motion 2 and to the Ceramill Coolstream.
- ▷ Plug the C/L drain pipe [4] together and connect it to the Ceramill Motion 2 and the Ceramill Coolstream (please refer to set-up instructions).
- ▷ Insert the extraction hose of the extraction system being used (e.g. Ceramill Airstream) to the extraction connection ([15] in Fig. 3 on page 53) on the rear side of the Ceramill Motion 2.

5.2 Service Unit

The Jäger spindle is equipped with a sealing-air feature. This sealing air prevents shavings/chips and dust from entering the spindle. The service unit filters possible contamination in the sealing air and thus prevents costly machine damage.

Installing the service unit


Fig. 11 Service unit

- 1 Knurled screws
- 2 Compressed-air hose, 250 mm
- 3 Compressed-air inlet
- 4 Water separator
- 5 Compressed-air connection, milling spindle
- 6 C/L drain pipe connection
- 7 C/L supply line connection

- ▷ Mount the service unit on the rear side of the machine using the knurled screws [1].
- ▷ Connect the 250 mm compressed-air hose [2] to the right-hand connection (viewed from the front) of the service unit and to the compressed-air connection [5] of the Ceramill Motion 2.

- ▷ Connect the 2,000 mm compressed-air hose to the compressed-air inlet [3] of the service unit and the other end to the compressed-air supply of the laboratory.

Requirements for the compressed-air supply

For operation, the compressed air must fulfil the following minimum requirements:

- _ Flow rate: 50 l/min (1.77 CFM)
- _ Pressure: 6 bar / 87 psi

Requirements on air purity:

- _ Solid contaminants: Class 3;
5 µm or less for solids
- _ Water content: Class 4;
max. pressure dew point +3°C
- _ Total oil content: Class 2;
max. oil content 0.1 mg/m³

In continuous operation, the pressure and the air-flow rate are monitored by the Ceramill Motion 2. When the required values are fallen below, the milling program stops. Once the correct air values are available again, the milling program can be continued after fault acknowledgement in the software.

Requirements on air purity for the Jäger spindle

When another service unit than the one supplied is to be used, the following connection conditions according to DIN-ISO 8573-1 must be ensured:

- _ Solid contaminants: Class 3
(particle size max. 5 μm , particle density max. 5 mg/m^3)
- _ Water content: Class 4
(max. pressure dew point +3 °C, water content max. 6000 mg/m^3)
- _ Total oil content: Class 2
(max. oil content 0.1 mg/m^3)

Extraction requirements

When another extraction system than the Ceramill Airstream is to be used, it must meet the following requirements:

- _ Suction capacity: 56.6 l/min
- _ Suitability for zirconium-oxide dusts
- _ HEPA microfilter (97.97 %), filter class H12, dust category M

The connection of the extraction hose to the machine is made with an adapter (in delivery scope). The adapter fits extraction hoses with a 38 mm interior diameter.

5.3 Installation of the Ceramill Match 2 and Ceramill Motion 2 Software

NOTE:

Malfunction!

▷ Operate the Ceramill Motion 2 only via PCs delivered by AmannGirrbach!

▷ Insert the Ceramill Match 2 CD-ROM into the drive of the PC.

The CD contains a PDF file with the installation guide.

▷ Open the PDF document and follow the instructions in the guide.

Connecting the interfaces


Fig. 12 Interfaces on the rear side of the Ceramill Motion 2

- 1 USB data port
 - 2 Ethernet interface
 - 3 Control interface for Ceramill Coolstream
 - 4 Slot for SD memory card
 - 5 Control interface for extraction
 - 6 Main switch
 - 7 Power supply connection
- ▷ Make sure that the main switch [6] of the Ceramill Motion 2 is in the 0 position.
 - ▷ Check if the SD memory card is properly seated in the slot [4] on the rear side of the Ceramill Motion 2.
 - ▷ Connect the mains cable to the power supply connection [7] of the Ceramill Motion 2 as well as to a mains socket outlet.

For installation without Ceramill Coolstream:

- ▷ Connect the control cable to the control interface for extraction [5] on the rear side of the Ceramill Motion 2 and to the interface on the rear side of the extraction system (e.g. to the Ceramill Airstream).

For installation with Ceramill Coolstream:


Fig. 13 Interfaces on the rear side of the Ceramill Coolstream

- 1 Power supply connection
 - 2 C/L supply line connection
 - 3 Control interface for Ceramill Coolstream
 - 4 Control interface for extraction (i. e. for Ceramill Airstream)
- ▷ Connect the control cable to the control interface for the Ceramill Coolstream on the rear side of the Ceramill Motion 2 ([3] in Fig. 12) and to the rear side of the Ceramill Coolstream ([3] in Fig. 13).
 - ▷ Connect the control cable of the extraction system to the control interface for extraction on the rear side of the Ceramill Coolstream ([4] in Fig. 13).
 - ▷ Connect the mains cable to the power supply connection ([1] in Fig. 13) of the Ceramill Coolstream as well as to a mains socket outlet.

Connecting to the PC

- ▷ Connect the USB cable to the USB data port ([1] in Fig. 12) of the Ceramill Motion 2 and to a free USB port of the PC.

The maximum permitted length of the USB cable is 5 m.

- ▷ Insert the Ceramill dongle in a free USB port on the rear side of the PC.


Bild 14 Ceramill dongle

5.4 Establishing the Network Connection

To establish the network connection, the Ceramill Motion 2 and the PC must first be connected via USB.

- ▷ Connect the Ethernet interface [2] of the Ceramill Motion 2 with the PC or with the network using a network cable.

The maximum permitted length of the network cable to the PC or to the network router is 100 m.

- ▷ Switch the Ceramill Motion 2 on.
- ▷ Switch the Ceramill Coolstream on.
- ▷ Start the Ceramill Motion 2 software.
- ▷ After initializing the milling machine in the software under *Settings > Communication*, set the connection type from USB to Ethernet and adapt the IP address, if required (see description of the Ceramill Motion 2 software).
- ▷ Switch the Ceramill Motion 2 off and then on again.
- ▷ Initialize the Ceramill Motion 2 in the software.
- ▷ After the software has detected the Ethernet connection, remove the USB cable.

When operating via a network, up to eight machines can be connected and controlled via the Ceramill Match Software. The PC and the Ceramill Motion 2 must be in the same sub-network.

5.5 Starting-up the Ceramill Motion 2

The three standard cutters in tool stations 1-3 are mandatory for operation.

In addition to standard tools, sister tools can optionally be inserted in the tool holder. When a standard tool has reached its wear limit, the Ceramill Motion 2 automatically accesses the respective sister tool.


Fig. 15 Tool holder (five-axis model)

On the four-axis model, the tool holder is turned by 90° (see figure 2 on page 52).

- ▷ Insert the Ceramill Motion Roto standard tools into the tool holder of the Ceramill Motion 2. Position the tool shank facing toward the outside.
 - Roto 2.5 is in the spindle.
 - Set Roto 1.0 to position 2.
 - Set Roto 0.6 to position 3.
- ▷ Optional: Insert the sister tools into the tool holder of the Ceramill Motion 2. Position the tool shank facing toward the outside.
 - Set Roto 2.5 to position 4.
 - Set Roto 1.0 to position 5.
 - Set Roto 0.6 to position 6.

- ▷ Insert the desired blank into the holder in the Ceramill Motion 2. The high/thick side of the blank must face the spindle.


Fig. 16 Inserting a blank (five-axis model)

- ▷ Fasten the blank using four screws and washers.

NOTE:

Do not damage the blank!

- ▷ Do not overtighten the screws! Avoid mechanical stress!

- ▷ Tighten the screws with the Allen key provided.


Fig. 17 Tightening a blank (five-axis model)

- ▷ Mount the respective hood.
 - The extraction hood for dry processing
 - The splash guard for wet processing

For wet processing:

- ▷ Insert the swarf sieve with the wide side facing toward the rear.

Without the swarf sieve, chips/swarf will be flushed into the container. The filter bag will fill up more quickly.

The swarf sieve is included in the delivery scope of the Ceramill Coolstream. It can remain inserted in the machine during dry processing.

- ▷ Shut the door of the Ceramill Motion 2.
- ▷ Switch the Ceramill Motion 2 on via the main switch.
- ▷ Switch the PC on.
- ▷ Start the Ceramill Motion 2 software.

When the software is started, a machine initialization takes place.

If the Ceramill Motion 2 is not detected:

- ▷ Restart the PC.

After re-positioning the Ceramill Motion 2, it is necessary to check the system axes.

- ▷ For this, observe the instructions for milling the test bodies.

6 Application and Operation**WARNING:**

Danger from shavings/chips or tool fragments being thrown from the machine!

- ▷ Make sure that the door of the Ceramill Motion 2 is always closed during processing!

**NOTE:**

Insufficient milling results and increased wear when dry processing!

- ▷ Do not operate the Ceramill Motion 2 without extraction (i.e. Ceramill Airstream)!

**NOTE:**

Insufficient milling/grinding results and increased wear when wet processing!

- ▷ Do not operate the Ceramill Motion 2 without the Ceramill Coolstream supply unit!

The instructions for use (video tutorial) of the Ceramill Match 2 software and the User Manual for operation of the machine are on the Ceramill-Match-CD or can be downloaded under www.ceramill-m-center.com.

6.1 Switching from Wet to Dry Processing

The Ceramill Motion 2 can be used both for wet and dry processing.

In delivery condition, it is set up for dry processing.

6.1.1 Switching from Dry to Wet Processing


Fig. 18

- 1 Extraction hood
- 2 Lever
- 3 Swarf sieve
- 4 Splash guard
- 5 Extraction hood
- B Service button

- ▷ Disengage lever [2] and pull off the white extraction hood [1] toward the front.
- ▷ Connect the extraction hose to the extraction opening in the milling chamber.
- ▷ With the door open, briefly press service button [B] on the front side of the Ceramill Motion 2. The extraction system starts.
- ▷ Allow the extraction system to clean out the milling chamber by means of vacuum.

When the milling chamber is clean:

- ▷ Press service button [B] again.

The extraction system stops.

- ▷ Pull off the extraction hose.
- ▷ Insert the cleaned swarf sieve [3].

The wide side faces toward the rear.

For wet processing, it is mandatory to use the swarf sieve [3]. The sieve serves as the first coolant/lubricant filter stage for coarse particles.

- ▷ Check if the outlet holes of the C/L nozzles are unobstructed.

- ▷ Check if the door gasket is clean and free of swarf/chips.
- ▷ Check if the C/L drain pipe is free of obstructions and if the piping on the outside is securely connected.
- ▷ Mount the transparent splash guard [4] to the extraction opening and lock it in place with lever [5].
- ▷ Transfer the milling file for wet processing.

6.1.2 Switching from Wet to Dry Processing


Fig. 19

- 1** Extraction hood
- 2** Lever of the extraction hood
- 3** Swarf sieve
- 4** Splash guard
- 5** Lever of the splash guard
- B** Service button

- ▷ Disengage lever [5] and pull off the transparent splash guard [4] toward the front.
- ▷ Remove the swarf sieve [3].
- ▷ Clean and dry the milling chamber and the swarf sieve (please refer to page 81).

The swarf sieve [3] is not required for dry processing.

- ▷ Insert the blank and insert or change the tooling, as required.
- ▷ Mount the white extraction hood [1] to the extraction opening and lock it in place with lever [2].
- ▷ Check the air extraction.
- ▷ Transfer the milling file for dry processing.

6.2 Changing the Blank Holder

In delivery condition, the Ceramill Motion 2 is equipped with a 71 blank holder. For processing varying materials, it may be required to use a different blank holder.

To change the blank holder:

- ▷ With the door closed, press operation button [A] for at least three seconds.

The current tool is placed down. The blank holder [2] moves to the change position.


Fig. 20 Five-axis model


Fig. 21 Four-axis model

Legend to Fig. 20 and Fig. 21:

- 1 Fastening screw for blank holder
- 2 Blank holder

- ▷ Loosen fastening screw [1].
- ▷ Pull off the blank holder toward the right (five-axis model) resp. toward the bottom (four-axis model).

- ▷ Clean the locating bore and surfaces.
- ▷ Check the pin on the blank holder to be mounted and clean, if required.
- ▷ Insert the blank holder with the pin to the stop in the locating bore.

NOTE:

Inaccurate working results or damage to the blank-holder fixture possible!

- ▷ Tighten fastening screw [1] with a torque of 4 Nm.

The supplied torque wrench is set accordingly.

- ▷ Tighten fastening screw [1] with with the supplied torque wrench.
 - ▷ Briefly press operation button [A].
- The machine advances to the starting position.

6.3 Checking the Model Area

The model area, on which the work is to be machined, must fit within the contour of the enclosed template (also see Match Package). Otherwise, the work cannot be machined in the blank. The template takes into account the shrinkage of the zirconium oxide in the final sintering.


Fig. 22 Template for Z1 71 zirconium oxide

- ▷ For work in wax or plastic, measure out the model area.

6.4 Function of the Buttons


Fig. 23 Buttons

Button	Door shut	Door open
Operation button [A]	Starts or stops the reference run; starts or stops the current program; activates or deactivates the control position; activates the change position for the blank holder	activates or deactivates the extraction system
Service button [B]	interrupts the current program; activates the service position; activates the tool-changing position; activates a cleaning run of the inspection window	activates or deactivates the sealing air

Tab. 3

Instead of with the buttons, these functions can also be actuated via the software (see instructions of the Ceramill Motion 2 software).

6.5 Reference Run

After switching on, a reference run must be carried out to determine the mechanical zero points. Operation button [A] flashes.

To start the reference run:

- ▷ Briefly press operation button [A].

The reference run is started. Operation button [A] lights up.

The reference run can be interrupted:

- ▷ Briefly press operation button [A].

The reference run is stopped. Operation button [A] flashes.

To continue the reference run:

- ▷ Briefly press operation button [A].

The reference run is continued at the position where it was stopped. Operation button [A] lights up.

On completion of the reference run, the machine returns to the starting position. Operation button [A] goes out.

6.6 Program Execution

Before executing a program, it must first be loaded (see Ceramill Match 2 software).

- ▷ Briefly press operation button [A].

The program is started. Operation button [A] lights up.

The program can be stopped:

- ▷ Briefly press operation button [A] or service button [B].

-or-

- ▷ Open the door.

The program is stopped. Operation button [A] flashes.

To continue the program:

- ▷ With the door closed, briefly press operation button [A].

The program is continued at the position where it was stopped. Operation button [A] lights up.

On completion of the program, the machine returns to the starting position. Operation button [A] goes out.

6.7 Checking the Workpiece During a Running Program

- ▷ Press operation button [A] for at least three seconds.

The program is stopped. The spindle returns to the rear, the current tool is placed down and the workpiece is turned around. While travelling, operation button [A] is lit; after reaching the position, it flashes.

- ▷ Visually check the workpiece.
- ▷ Briefly press operation button [A].

The tool is picked up again. The spindle travels to the position, where the program was interrupted. While travelling, operation button [A] is lit; after reaching the position, it flashes.

- ▷ Briefly press operation button [A].
- The program is continued. Operation button [A] lights up.

6.8 Moving to the Service Position

The service position can be moved to from the starting position or from a stopped program.

- ▷ Briefly press service button [B].

The current tool is placed down and the spindle returns to the service position. The spindle is released; the motor is blocked. While travelling, operation button [A] is lit; after reaching the position, it flashes.

To exit the service position:

- ▷ Briefly press operation button [A].

The tool is picked up again. The spindle travels back to the starting position or to the position, where the program was interrupted. While travelling, operation button [A] is lit; after reaching the position, it flashes.

6.9 Inspection-window Cleaning

- ▷ With the door closed, press operation button [B] for at least three seconds.
- As long as the button is pressed, an air nozzle positioned behind the inspection window blows off the window and cleans it from dust and chips.

The inspection-window cleaning is also carried out automatically when a tool change takes place and after completion of a program. In the course of each automatic cleaning run, three short air-blasts take place.

6.10 Advancing to the Change Position for the Blank Holder

The change position for the blank holder can only be advanced to from the starting position.

- ▷ Press operation button [A] for at least three seconds.

The current tool is placed down. The blank holder is turned in such a manner that the fastening screw becomes accessible. While travelling, operation button [A] is lit; after reaching the position, it flashes.

To exit the change position for the blank holder:



- ▷ Briefly press operation button [A].

The machine advances to the starting position.

7 Registration and Download Information


7.1 Quick Registration, M-center

Before registering, it is advisable to note down the serial numbers of the various equipment/devices.

Device/Machine	Location	Example	Serial number
[1] Scanner (Ceramill Map)	Rear side of the device		
[2] PC	Top of PC housing		
[3] Mind/match dongle	Software box		
[4] Motion 2 (Ceramill Motion 2)	Rear side of the machine		

Tab. 4

- ▷ To register, log on to www.ceramill-m-center.com. On this page, click on the “Register” link.
- ▷ Afterwards, fill out the form field with the personal data
- ▷ Enter all serial numbers.

- ▷ Accept the terms of use.

- ▷ Click on the “Create account” button.
A confirmation E-mail will be generated after a short period and sent to the previously entered E-mail address.
- ▷ Open the E-mail and click on the link in it.
This opens a Web page that confirms the registration.

If the link cannot be clicked on:

- ▷ Copy and paste the link into your Web browser.

Your E-mail address and a password of your choice are the access data for entering the protected homepage site www.ceramill-m-center.com (order tracking, video tutorials).

Within the next 24 hours you will receive another E-mail. It contains the registration confirmation and the FTP access data (6-digit) for uploading data from our M-Center.

Two sets of access data are provided after the registration. One for the M-Center website (first E-mail) and one for the upload of your data to our M-Center (second E-mail).

▷ Enter all access data into the Ceramill Mind design software (the window opens automatically after clicking on the Upload button ).

7.3 Download Infos

Further instructions (video tutorials) on operating the Ceramill Map scanner and the Ceramill Mind software are available on the Instructions DVD or can be downloaded under www.ceramill-m-center.com.

Software updates are also available on our website under www.ceramill-m-center.com.

The registration process is now finished.

7.2 Order Tracking, M-Center

After receiving the design data in the M-Center, an order confirmation will be sent via E-mail.

Additionally, you will be informed via E-mail when your work goes into production and is dispatched.

A detailed list of all orders will be displayed on the M-Center homepage in the login area under "My Orders". Here, each order will be allocated a tracking number after the order has been dispatched, with which it can be monitored.

Order number: ZZZ999		Complaint
Name:	Ceramill ZIL	
Received:	10.01.2011 17:39:55	
Planned date of dispatch:	12.01.2011	
Color:	CL1 100%	
STL-File:	2010-11-22_00012-XXX-17-m-center_cad.stl	
DR-Trackingnumber:	Shipment number: 001122334455	

8 Handling the Coolant/Lubricant (C/L)

CAUTION:

Health hazards are possible through improper handling of C/L!

- ▷ Use only Ceramill coolant as C/L!
- ▷ Observe the safety instructions and handling regulations for Ceramill coolant!

- ▷ Insert the reaction zone of the test strips into the C/L.
- ▷ Shake off excess fluid.
- ▷ After one minute, compare the reaction zone with the colour scale and read off the measured value.

The reference values and the information on measures in case of exceeding or falling below the reference values can be found on an information sheet on the inside of the container-access cover in the Ceramill Coolstream.

8.1 C/L Analysis

For legal and technical reasons, the coolant/lubricant must be analyzed weekly.

- ▷ Observe country-specific laws and guidelines on coolant/lubricant!

The analysis protocol can also be placed in the self-adhesive transparent folder located there.

Regular analysis of the C/L is required for proper function of the Ceramill Motion 2. The results of the analyses are to be documented and archived.

- ▷ Check the following coolant/lubricant properties weekly:
 - Concentration
 - pH-value
 - Nitrate content

The analyses are to be carried out with test strips. These test strips are simple to handle and enable pH-value, the nitrate content and the C/L concentration to be determined.

- ▷ Remove the blank from the holder.
- ▷ Mount the splash guard.
- ▷ Place a clean, flat receptacle in the center of the interior below the splash guard.
- ▷ Shut the door.
- ▷ Push operation button [A] for approx. eight seconds.

After three seconds, the C/L flow is started.

- ▷ Open the door and remove the receptacle.

8.2 Removing the C/L Container

For table version:

- Remove the C/L drain pipe [1] on the side of the Ceramill Coolstream.


Fig. 24

1 C/L drain pipe

NOTE:

Possible damage of the electric lines and the C/L supply line:

- When moving the Ceramill Coolstream, pay attention that no lines are driven over.

- Move the Ceramill Coolstream out of the table.

For both versions:

- Remove the left sidewall of the Ceramill Coolstream.
- Pull the C/L drain pipe [1] and the C/L supply line [3] from above out of the C/L container.
- Disconnect the control cable [2] at the partition wall (between the chamber for the Ceramill Air-stream and the container).


Fig. 25 Disconnecting the interfaces of the C/L container (tower version)


Fig. 26 Disconnecting the interfaces of the C/L container (table version)

Legend to Fig. 25 and Fig. 26:

- 1 C/L drain pipe
- 2 Control cable
- 3 C/L supply line

NOTE:

Do not damage the drain cock!

► Move the C/L container only by its handles!

- Carefully lift the C/L container [2] out of the supply unit.


Fig. 27

- 1 Ceramill Coolstream
- 2 C/L container

8.3 Emptying the C/L Container

NOTE:

Do not damage the drain cock!

► Move the C/L container only by its handles!

- Place the C/L container [1] with its rear side facing frontward on a table or, i.e. on the Ceramill Coolstream.


Fig. 28

- 1 C/L container
- 2 Drain cock
- 3 Disposal receptacle

- Position a disposal receptacle [3] below drain cock [2].
- If required, mount a hose to the drain cock and conduct it into the disposal receptacle.

- ▷ Open the drain cock.

The C/L in the C/L container drains into the disposal receptacle.

When the C/L container is empty:

- ▷ Shut the drain cock again.

Dispose of the C/L in accordance with the applicable regulations and laws. Proper disposal of the C/L must be documented.

8.4 Cleaning the C/L Container

NOTE:

The filter and the C/L container may only be cleaned with water.

- ▷ Do not use cleaning agents.

- ▷ Unscrew and remove the four star-knob nuts on the top side of the C/L container.
- ▷ Remove the cover with the drain pipe.
- ▷ Remove the filter bag and empty it.
- ▷ Clean the filter bag.
- ▷ Rinse and clean the C/L container.

8.5 Filling in C/L

The water for the coolant/lubricant must have drinking-water quality. The optimal hardness range is between 10 - 20 °dH (corresponds to 1.8 - 3.6 mmol total hardness per litre). When the hardness is far beyond this range, we recommend the addition of distilled water.

- ▷ Fill 20 litres of water and 1 litre of Ceramill Coolant into a separate container and mix well.
- ▷ Insert the filter bag into the C/L container.
- ▷ Pour the mixed C/L into the C/L container.
- ▷ Reattach the cover with the drain pipe and secure it with the four star-knob nuts.

8.6 Installing the C/L Container

- ▷ Carefully lift the C/L container into the supply unit.
- ▷ Connect the control cable at the partition wall (between the chamber for the Ceramill Airstream and the container).
- ▷ Connect the return pipe and the C/L supply line at the top of the C/L container.
- ▷ Mount the left sidewall of the Ceramill Coolstream.

NOTE:

Danger of damaging cable and hoses!

- ▷ When moving the supply unit, do not pinch or drive over any cables or hoses.

For the tower version:

- ▷ Connect the C/L drain pipe and the C/L supply line of the Ceramill Coolstream to the Ceramill Motion 2.

For the table version:

- ▷ Move the supply unit back under the table and align it with the Ceramill Motion 2.
- ▷ Insert the drain pipe on the side to the stop into the Ceramill Coolstream.

Through this, the control system recognizes the drain pipe and enables the C/L pump.

- ▷ Connect the C/L supply line of the Ceramill Coolstream to the Ceramill Motion 2.

8.7 Cleaning the C/L System

- ▷ Remove and clean the swarf sieve.
- ▷ Clean the milling chamber.
- ▷ Reinsert the swarf sieve again.
- ▷ Mount the splash guard.
- ▷ Empty the C/L container (please refer to page 78).
- ▷ Clean the C/L container (please refer to page 79).
- ▷ Fill 5 litres of water and 250 ml of Ceramill Coolant into a separate container and mix well.
- ▷ Insert the filter bag.
- ▷ Pour the mixed C/L into the C/L container.
- ▷ Fully connect the Ceramill Coolstream again.
- ▷ Activate "Clean C/L circuit (Kühlmittelkreislauf reinigen)" in the control service functions.
- ▷ Start the cleaning in the software or

-or-

- ▷ briefly press operation button [A] on the Ceramill Motion 2.
- C/L is circulated for 15 minutes to clean the system.

The cleaning can be interrupted and restarted by briefly pressing operation button [A].

When the cleaning is completed:

- ▷ Empty the C/L container (please refer to page 78).

After the cleaning, the C/L used for the cleaning must be disposed of.

- ▷ Refill the C/L container again (please refer to page 79).

9 Cleaning and Maintenance

9.1 Ceramill Motion 2

9.1.1 Cleaning

The Ceramill Motion 2 must be cleaned after each working day.

- ▷ Empty the swarf sieve frequently, especially when wet grinding plastic.

When the swarf sieve is too full, it is possible that the coolant/lubricant will not drain off.

Cleaning after dry processing

NOTE:

Zirconia dust is abrasive!

- ▷ Clean the window in the door carefully, using only a soft cloth; do not scrub or scour.

- ▷ Open the door of the Ceramill Motion 2.
- ▷ Replace the extraction hood with the cleaning hose with brush head.
- ▷ Press operation button [A] on the Ceramill Motion 2.

The extraction system is activated. The operation button lights up.

- ▷ Press service button [B] on the Ceramill Motion 2.

The sealing air at the spindle is activated. The service button lights up.

- ▷ If required, remove the swarf sieve and dispose of accumulated milling dust/shavings away from the machine.

By doing so, the filter bag in the extraction system will not fill up as quickly and can be used longer.

- ▷ Thoroughly vacuum out the milling chamber with the cleaning hose.

Residual zirconia dust from dry processing can clog the filter bag in case of subsequent wet processing.

- ▷ Briefly press service button [B].
The sealing air is deactivated.
- ▷ Briefly press operation button [A].
The extraction system is deactivated.

Cleaning after wet processing

NOTE:

Ceramill Airstream is not designed for wet vacuuming. The Ceramill Airstream may not be used for cleaning the interior after wet processing!

- ▷ Remove shavings/chips manually.

NOTE:

The milling chamber may only be cleaned with water.

- ▷ Do not use cleaning agents.

Dried-in coolant/lubricant is difficult to remove.

After wet processing:

- ▷ Clean and dry the milling chamber prior to downtimes lasting several hours.

- ▷ Remove the swarf sieve and dispose of accumulated milling dust/shavings.
- ▷ Clean the milling chamber and wipe it dry using a soft cloth.

9.1.2 Regular Inspections

- _ Condition of the bellow
- _ Condition of the hoses and drain pipes
- ▷ Keep the door gasket and the stopping surface on the door clean. Ensure not to damage the sealing lips.

9.1.3 Weekly Maintenance

NOTE:

- Danger of damaging the machine!
- ▷ Do not use compressed air, ultrasonic or a steam jet for the cleaning!

Prior to maintenance

Before maintaining the spindle and the C/L nozzles, the spindle sealing air must be activated:

- ▷ Briefly press service button [B].
The current tool is placed down and the spindle returns to the service position. The spindle is released; the motor is blocked.
- ▷ Open the door of the Ceramill Motion 2.
- ▷ Press service button [B] on the Ceramill Motion 2.
The sealing air is activated. The service button lights up.

After the maintenance

When the Maintenance is completed:

- ▷ Shut the door of the Ceramill Motion 2.
- ▷ Briefly press operation button [A].
The tool is picked up again. The spindle travels back to the starting position or to the position, where the program was interrupted.

Spindle

The provided Jäger spindle service set and a milling tool of your choice are required for maintenance of the spindle.


Fig. 29 Service set for Jäger spindle

- 1 Collet grease
- 2 Collet brush
- 3 Felt cone
- 4 Ejector pin
- 5 Spindle wrench

In order to clean the spindle:

- ▷ Remove any dust from the cap, the C/L lines and the spindle, using a dry, soft cloth. Do not use cleaning agents!

To remove the collet:

- ▷ Apply the black spindle wrench [5] over the collet.


Fig. 30

- ▷ Insert a milling tool into the open collet.

- ▷ Using the spindle wrench, unscrew the collet out of the spindle.


Fig. 31

- ▷ Clean the collet holder using the felt cone [3].
The collet holder must be free of shavings/chips and debris.


Fig. 32

- ▷ Clean the inside of the collet with the fine collet brush [2].


Fig. 33

- ▷ Apply a light coat of the collet grease ([1] in Fig. 29 on page 82) from the service set to the outer cone of the collet.

This improves the sliding properties and increases the clamping force of the collet.

NOTE:

Only use the provided collet grease for lubrication of the collet cone. Do not use other types of grease.

- ▷ Screw the cleaned and lubricated collet in reverse order to the stop back into the spindle.

CAUTION:

Danger of tool breakage!

- ▷ Make sure to remove the milling tool from the collet chuck!

When the milling tool was from the Ceramill Motion 2 tool holder:

- ▷ Place the tool properly positioned back into the tool holder.

C/L nozzles

Do not clean the C/L nozzles with the collet brush!

- ▷ Inspect the C/L nozzles weekly.
- ▷ Clean the outlet holes if required.

When the outlet holes are clogged:

- ▷ Loosen the Allen screw on the side of the spindle cap.


Fig. 34

- ▷ Pull off the cap.


Fig. 35

- ▷ Blow out the holes from the inside of the cap using a steam cleaner.

CAUTION:

The cap gets hot!

- ▷ Apply a coat of grease to the O-rings of the C/L pipes.
- ▷ Clean the seat of the cap on the spindle.

- ▷ Mount the spindle cap again with the screw at the right-hand side. Pay attention to the correct connection of the C/L pipes.

Service unit

- ▷ Visually check the service unit for:
 - Visible particles, debris or particulate matter on the inside of the acrylic glass (signs for contaminated compressed air)
 - Amber discolouration of the filter element (sign of oil in the compressed air)
 - Water in the filter (sign of water in the compressed air)

Once any of the listed impairments is given:

- ▷ Correct the problem with the compressed air.
- ▷ Replace the complete service unit.

9.1.4 Monthly Maintenance

Checking the zero point

It is recommended to check the zero point once every month.

- ▷ Mill a test cube from the Ceramill TEST 71 L test blank.

The test blank is to be used exclusively for test milling and calibration of the Ceramill Motion 2. The corresponding instructions are on the Ceramill Match CD or available on our website under www.ceramill-m-center.com.

Jäger spindle

The spindle must be in the starting position.

- ▷ Manually turn the shaft of the spindle through (at least ten times).


Fig. 36 Turning the Jäger spindle through

For further information on the Jäger spindle, please download the Manual from our website under www.amanngirrbach.com.

9.1.5 External Maintenance

External maintenance is necessary after 1000 operating hours. The Ceramill Motion 2 software will inform you appropriately.

- ▷ Contacting the Ceramill helpdesk (workdays from 8 a.m. to 5 p.m.).
 - Germany: +49 7231 957 100
 - Austria: +43 5523 62333 390
 - International: +43 5523 62333 399
 - helpdesk@amanngirrbach.com

9.2 Ceramill Coolstream

9.2.1 Inspections and Maintenance Measures

Regular inspections:

- _ Condition of the hoses and drain pipes
- _ Filling level of the C/L in the supply unit of the Ceramill Coolstream

When the filling level in the C/L container has fallen below the minimum level, the milling program is interrupted and an error message is indicated.

- ▷ Refill the C/L container (please refer to page 79).
- ▷ Acknowledge the error in the software.

9.2.2 Weekly Maintenance

- ▷ Analyse the C/L (see separate analysis instructions with test record).

9.2.3 Half-yearly Maintenance

- ▷ Clean the C/L system (please refer to page 80).

The half-yearly cleaning is indicated via the program. After activation of the system to be cleaned, the indication is omitted again.

10 Malfunctions, Repairs and Warranty

10.1 Malfunctions

In case of malfunctions:

- ▷ Restart the software.
- ▷ Restart the Ceramill Motion 2.
- ▷ Restart the PC.

If these measures do not correct the problem:

- ▷ Contacting the Ceramill helpdesk (workdays from 8 a.m. to 5 p.m.).
 - Germany: +49 7231 957 100
 - Austria: +43 5523 62333 390
 - International: +43 5523 62333 399
 - helpdesk@amanngirrbach.com

10.2 Repairs

Repairs may be carried out only by trained, specialised personnel.

10.3 Warranty

The warranty complies with the statutory provisions. For more information, please refer to our General Terms of Business.

11 Environmental Protection

Packaging

In terms of packaging, AmannGirrbach participates in country-specific recycling systems, which ensure optimal recycling.

All packaging materials used are environmentally-friendly and recyclable.

Waste equipment

Waste equipment contain materials that should be sorted for environmental-friendly recycling.

- ▷ After the service life has elapsed, dispose of the machine via the public disposal systems in an environmentally-friendly manner.

The plastic parts are labelled. Thus, they can be sorted for recycling or disposal.

Coolant/Lubricant (C/L)

- ▷ Dispose of the C/L in accordance with the applicable regulations and laws.

Proper disposal of the C/L must be documented.

12 Technical Data and Accessories/Spare Parts

 Subject to changes.

Technical Data - Ceramill Motion 2

	Unit	Value
Art. No.		
_ Four-axis model	–	179280
_ Five-axis model	–	179250
Dimensions (D × W × H)	mm	588 × 516 × 755
Weight:		
_ Four-axis model	kg	75
_ Five-axis model	kg	78
Electrical connections	V/Hz	100-230/50-60
Power output	W	250
Fuse (time-delayed)	A	3.15 / 6.3
Axes:		
_ Four-axis model	–	4
_ Five-axis model	–	5
Accuracy	μm	< 10
Allowable temperature range (operation)	°C	18 - 30
Spindle torque	Ncm	4
Spindle speed	rpm	60,000
Collet diameter	mm	3
Sound pressure level	dB(A)	60

Tab. 5

	Unit	Value
Interfaces	–	USB / Ethernet / slot for SD card / Ceramill Cool-stream control system / Extraction control system
Compressed-air connection	–	Dry, clean compressed air
_ Pressure	bar	6
_ Flow rate	l/min	50
Air purity		
_ Solid contaminants	–	Class 3; 5 μm or less for solids
_ Water content	–	Class 4; max. pressure dew point +3 °C
_ Total oil content:	–	Class 2; max. oil content 0.1 mg/m ³
Extraction	–	Suitable for zirconium oxide dust particles, HEPA microfilter (97.97 %), filter class H12, dust category M
_ Diameter of extraction connection	mm	38
_ Suction capacity	l/min	56.6

Tab. 5

Technical Data - Ceramill Coolstream

	Unit	Value
Art. No.	–	178630
Dimensions (D × W × H)	mm	656 × 518 × 651
Weight	kg	36
Electrical connections	V/Hz	100-230/50-60
Power output	W	150
Fuse (time-delayed)	A	1.0
Max. pump flow rate	l/min	3.8
C/L container, volume	l	50
Mesh size of filter	μm	50
Interfaces	–	Coolstream control system / Extraction control system

Tab. 6

Accessories/Spare parts

Art. No.:	Designation
179210	Service unit
179211	Screw set for blank holder
179216	Threaded inserts for blank holder
179217	Test specimen for test cube
179218	Service set, Jäger spindle
179251	Cover, tool holder
179230	Collet, 3 mm
179252	Blank holder 71
179253	Swarf sieve
179254	Extraction hood
179255	Splash guard
179256	Filter bag
179257	Extraction hose with brush
178650	Ceramill coolant (coolant/lubricant)
178651	Multi-function additive test (100 pce.) Test strips
178652	pH-value measuring strips (100 pce.)
178653	Nitrite test (100 pce.) Test strips
1156901	Fuse for Ceramill Motion 2 (T 3.15 A / 250 V)
172334	Fuse for Ceramill Coolstream (T 1.0 A / 250 V)
179186	Blank holder 71 Ceramill Motion 2 (4X)
179184	Extraction hood Ceramill Motion 2 (4X)
179188	Splash guard Ceramill Motion 2 (4X)
179189	Control system, complete Ceramill Motion 2 (4X)

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Amann Girrbaach AG
Amann Girrbaach GmbH