

DIRECT LDL

Package insert instructions must be carefully followed. Reliability of assay results cannot be guaranteed if there are any deviations from the instructions in this package insert.

Read Highlighted Changes: Revised February 2022.

INTENDED USE

The MULTIGENT Direct LDL assay is used for the direct, quantitative determination of low-density lipoprotein (LDL) cholesterol in human serum or plasma on the ARCHITECT cSystems.

SUMMARY AND EXPLANATION OF TEST

Plasma lipoproteins are spherical particles containing varying amounts of cholesterol, triglycerides, phospholipids, and proteins. The phospholipid, free cholesterol, and protein constitute the outer surface of the lipoprotein particle, while the inner core contains mostly esterified cholesterol and triglycerides. These particles serve to solubilize and transport cholesterol and triglycerides in the bloodstream.

The relative proportions of protein and lipid determine the density of these lipoproteins and provide a basis on which to begin their classification.¹ These classes are: chylomicrons, very-low-density lipoprotein (VLDL), low-density lipoprotein (LDL), and high-density lipoprotein (HDL). Numerous clinical studies have shown that the different lipoprotein classes have very distinct and varied effects on coronary heart disease (CHD) risk.²⁻⁴

The studies all point to LDL cholesterol as the key factor in the pathogenesis of atherosclerosis and CHD,²⁻⁸ while HDL cholesterol has been observed to have a protective effect. Even within the normal range of total cholesterol concentrations, an increase in LDL cholesterol can occur with an associated increased risk for CHD.⁴

PRINCIPLES OF PROCEDURE

The MULTIGENT Direct LDL assay is a homogeneous method for directly measuring LDL levels in serum or plasma, without the need for off-line pretreatment or centrifugation steps.

The method is in a two-reagent format and depends on the properties of a unique detergent. This detergent, **R1**, solubilizes only the non-LDL particles. The cholesterol released is consumed by cholesterol esterase and cholesterol oxidase in a non-color-forming reaction. A second detergent, **R2**, solubilizes the remaining LDL particles and a chromogenic coupler allows for color formation. The enzyme reaction with LDL in the presence of the coupler produces color that is proportional to the amount of LDL cholesterol present in the sample.

Methodology: Measured, Liquid Selective Detergent

REAGENTS

Reagent Kit

REF 1E31-20 MULTIGENT Direct LDL is supplied as a liquid, ready-to-use, two-reagent kit which contains:

R1 2 x 53 mL

R2 2 x 20 mL

Estimated tests per kit: 450

Calculation is based on the minimum reagent fill volume per kit.

Reactive Ingredients	Concentration
R1 MES buffer (pH 6.3)	
Detergent 1	< 1.0%
Cholesterol esterase (Microorganism)	< 1,500 U/L
Cholesterol oxidase (Microorganism)	< 1,500 U/L
Peroxidase (Horseradish)	< 1,300 ppb U/L
4-aminoantipyrine	< 0.01%
Ascorbic acid oxidase (Cucurbita sp.)	< 3,000 U/L
Preservative	
R2 MES buffer (pH 6.3)	
Detergent 2	< 1.0%
N,N-bis(4-sulfobutyl)-m-toluidine, disodium (DSBmT)	< 1.0 mmol/L
Preservative	

REAGENT HANDLING AND STORAGE

Reagent Handling

Remove air bubbles, if present in the reagent cartridge, with a new applicator stick. Alternatively, allow the reagent to sit at the appropriate storage temperature to allow the bubbles to dissipate. To minimize volume depletion, do not use a transfer pipette to remove the bubbles.

CAUTION: Reagent bubbles may interfere with proper detection of reagent level in the cartridge, causing insufficient reagent aspiration that could impact results.



en

Direct LDL

REF 1E31-20

307191R14

B1E310

FOR USE WITH

ARCHITECT

Reagent Storage

Unopened MULTIGENT Direct LDL reagents are stable until the expiration date shown on the label when stored at 2 to 8°C. Do not freeze.

Reagent stability is 28 days (672 hours) if the reagent is uncapped and onboard.

Indications of Deterioration

Instability or deterioration should be suspected if there are precipitates, visible signs of leakage, extreme turbidity, microbial growth, if calibration does not meet the appropriate package insert and/or **ARCHITECT System Operations Manual** criteria, or if controls do not meet the appropriate criteria.

WARNINGS AND PRECAUTIONS

Precautions for Users

- **IVD**
- For *In Vitro* Diagnostic Use.
- Do not use reagents beyond the expiration date.
- Do not mix reagents from different kit lot numbers.
- Protect reagents from direct sunlight.
- Do not freeze reagents.
- **CAUTION:** This product requires the handling of human specimens. It is recommended that all human-sourced materials be considered potentially infectious and handled in accordance with the OSHA Standard on Bloodborne Pathogens.⁹ Biosafety Level 2¹⁰ or other appropriate biosafety practices^{11,12} should be used for materials that contain or are suspected of containing infectious agents.

The following warnings and precautions apply to: **R1** and **R2**

WARNING:	Contains methyisothiazolones.
H317	May cause an allergic skin reaction.
Prevention	
P261	Avoid breathing mist / vapors / spray.
P272	Contaminated work clothing should not be allowed out of the workplace.
P280	Wear protective gloves / protective clothing / eye protection.
Response	
P302+P352	IF ON SKIN: Wash with plenty of water.
P333+P313	If skin irritation or rash occurs: Get medical advice / attention.
P362+P364	Take off contaminated clothing and wash it before reuse.
Disposal	
P501	Dispose of contents / container in accordance with local regulations.

Safety Data Sheets are available at www.abbottdiagnostics.com or contact your local representative.

For a detailed discussion of safety precautions during system operation, refer to the **ARCHITECT System Operations Manual**, Section 8.

SPECIMEN COLLECTION AND HANDLING

Suitable Specimens

Serum and plasma are acceptable specimens. Patients are not required to fast prior to blood collection.

- **Serum:** Use serum collected by standard venipuncture techniques into glass or plastic tubes with or without gel barriers. Ensure complete clot formation has taken place prior to centrifugation. Separate serum from red blood cells or gel as soon after collection as possible (within 3 hours).¹³
Some specimens, especially those from patients receiving anticoagulant or thrombolytic therapy, may take longer to complete their clotting processes. Fibrin clots may subsequently form in these sera and the clots could cause erroneous results.
- **Plasma:** Use plasma collected by standard venipuncture techniques into glass or plastic tubes. Acceptable anticoagulants are lithium heparin (with or without gel barrier), sodium heparin, and EDTA. Anticoagulants containing citrate should not be used. Ensure centrifugation is adequate to remove platelets. Separate plasma from red blood cells or gel as soon after collection as possible (within 3 hours).¹³

Refer to the specimen collection tube manufacturer's instructions for processing and handling requirements.

For total sample volume requirements, refer to the ASSAY PARAMETERS section of this package insert and *Section 5* of the **ARCHITECT System Operations Manual**.

Specimen Storage

Serum and Plasma

Temperature	Maximum Storage
2 to 8°C	5 days
-80°C	3 months

Samples may be frozen once. Refer to Clinical and Laboratory Standards Institute (CLSI) document NCCLS H18-A3 for further instructions on specimen collection, handling, and storage.¹⁴

NOTE: Stored specimens must be inspected for particulates. If present, mix and centrifuge the specimen to remove particulates prior to testing.

PROCEDURE

Materials Provided

REF 1E31-20 MULTIGENT Direct LDL Reagent Kit

Materials Required but not Provided

- **REF** 5P56-01 or 5P56-02 Lipid Multiconstituent Calibrator (Lipid MCC)
- LDL cholesterol control sera or quality control material
- Saline (0.85% to 0.90% NaCl) for specimens that require dilution

Assay Procedure

For a detailed description of how to run an assay on the ARCHITECT cSystems, refer to *Section 5* of the **ARCHITECT System Operations Manual**.

Specimen Dilution Procedures

Serum and Plasma: The following specimens should be diluted by following the Manual Dilution Procedure:

- Specimens with levels of interfering substances (other than triglyceride) higher than the upper limit stated in the Interfering Substances section.
- Specimens with LDL cholesterol values exceeding 800 mg/dL (20.69 mmol/L).

Manual Dilution Procedure

Manual dilutions should be performed as follows:

- Use saline (0.85% to 0.90% NaCl) to dilute the sample.
- The operator must enter the manual dilution factor in the patient or control order screen. The system uses this dilution factor to automatically correct the concentration by multiplying the result by the entered factor.
- If the operator does not enter the dilution factor, the result must be multiplied by the appropriate dilution factor before reporting the result.

$$\text{Manual Dilution Factor} = \frac{(\text{Volume of Sample} + \text{Volume of Dilution Reagent})}{\text{Volume of Sample}}$$

NOTE: If a diluted sample result is flagged indicating it is less than the linear low limit, do not report the result. Rerun using an appropriate dilution.

For detailed information on ordering dilutions, refer to *Section 5* of the **ARCHITECT System Operations Manual**.

CALIBRATION

Calibration is stable for approximately 28 days (672 hours) and is required with each change in reagent lot number. Verify calibration curve with at least two levels of controls according to the established quality control requirements for your laboratory. If control results fall outside acceptable ranges, recalibration may be necessary.

For a detailed description of how to calibrate an assay, refer to *Section 6* of the **ARCHITECT System Operations Manual**.

For information on calibrator standardization, refer to the package insert for **REF** 5P56-01 Lipid Multiconstituent Calibrator.

QUALITY CONTROL

Reliability of test results should be routinely monitored with quality control materials or serum pools that reasonably represent performance on patient specimens.¹⁵

As appropriate, refer to your laboratory standard operating procedure(s) and/or quality assurance plan for additional quality control requirements and potential corrective actions.

- Run two levels of controls (normal and abnormal) once during each day of assay use.
- If more frequent control monitoring is required, follow the established quality control procedures for your laboratory.
- If quality control results do not meet the acceptance criteria defined by your laboratory, patient values may be suspect. Follow the established quality control procedures for your laboratory. Recalibration may be necessary.
- Review quality control results and acceptance criteria following a change of reagent or calibrator lot.

RESULTS

Refer to *Appendix C* of the **ARCHITECT System Operations Manual** for information on results calculations.

To convert results from mg/dL to mmol/L, multiply mg/dL by 0.02586.¹⁶

Representative performance data are given in the SPECIFIC PERFORMANCE CHARACTERISTICS section of this package insert. Results obtained in individual laboratories may vary.

LIMITATIONS OF THE PROCEDURE

Refer to the SPECIMEN COLLECTION AND HANDLING and SPECIFIC PERFORMANCE CHARACTERISTICS sections of this package insert.

EXPECTED VALUES

The following National Cholesterol Education Program (NCEP) cutpoints for patient classification are used for the prevention and management of coronary heart disease.¹³

Reference Range

LDL Cholesterol (mg/dL)	
< 100	Optimal
100 to 129	Near or above optimal
130 to 159	Borderline high
160 to 189	High
≥ 190	Very high

The NCEP Adult Treatment Panel III Report recommends the classification shown above. Laboratories should follow recommendations for lipid ranges effective in their locale if they differ from those of the NCEP.

SPECIFIC PERFORMANCE CHARACTERISTICS

Reportable Range

The reportable range of the MULTIGENT Direct LDL assay is from 1 to 800 mg/dL (0.03 to 20.69 mmol/L), based on limit of blank and linearity studies. Linearity was evaluated using N-geneous LDL Linearity Verifiers, part number 80-5953-00. Linearity samples recovered within 10% of the theoretical value up to 800 mg/dL.

Limit of Blank (LOB)

The LOB for MULTIGENT Direct LDL is ≤ 10 mg/dL (0.259 mmol/L). The LOB is the mean concentration of an analyte-free sample plus 2 SD, where SD = the pooled, within-run standard deviation of the analyte-free sample. A study performed on an ARCHITECT cSystem produced an LOB for MULTIGENT Direct LDL of 0.1 mg/dL (0.003 mmol/L).

Interfering Substances

Samples with triglyceride concentrations > 1,293 mg/dL (14.61 mmol/L) should not be used for the determination of LDL cholesterol.

Potential interference in the MULTIGENT Direct LDL assay from ascorbic acid, bilirubin, gamma globulins, and hemoglobin is less than 10% at the levels indicated below. A study based on guidance from NCCLS publication EP7-P¹⁷ was performed using the MULTIGENT Direct LDL assay on a commercially available clinical chemistry analyzer. Varying amounts of potential interferents were added to serum pools with known quantities of cholesterol. No significant interference was detected in the MULTIGENT Direct LDL assay up to and including the concentrations stated below:

Interfering Substance	Concentration with No Significant Interference
Ascorbic acid	50 mg/dL (2,839 µmol/L)
Bilirubin (conjugated)	20 mg/dL (342 µmol/L)
Bilirubin (unconjugated)	20 mg/dL (342 µmol/L)
Gamma globulins	5,000 mg/dL (50 g/L)
Hemoglobin	500 mg/dL (5 g/L)

Interferences from medications or endogenous substances may affect results.¹⁸

The following drugs were tested for interference at the concentrations indicated using an acceptance criteria of $\pm 10\%$ from the target value. Direct LDL is not affected by the presence of the following interferents up to and including the concentrations stated below:

Interfering Substance	Concentration with No Significant Interference
Acetaminophen	200 mg/L (1324.5 µmol/L)
Dipyron	100 mg/L (300.3 µmol/L)
N-Acetyl-L-Cysteine	1600 mg/L (9816.0 µmol/L)

Precision

The precision of the MULTIGENT Direct LDL assay is < 4% CV. A within-run precision study was performed using two levels of LDL cholesterol control human serum. Each run consisted of 20 replicate samples.

In addition, a between-run precision study was performed using two levels of frozen pooled human serum in accordance with CLSI protocol NCCLS EP15-P.¹⁹

Representative data from these studies are summarized below.

	Control	Low < 130 mg/dL	High ≥ 160 mg/dL
	N	20	20
Within Run	Mean (mg/dL)	86.8	175.9
	SD	0.95	2.46
	%CV	1.1	1.4
Between Run	Mean (mg/dL)	89.0	178.1
	SD	1.92	3.94
	%CV	2.2	2.2

Accuracy

The average percent bias of the MULTIGENT Direct LDL reagent on a commercially available clinical chemistry analyzer to the Reference Method (ultracentrifugation and cholesterol analysis) was less than or equal to the guidelines established by the National Cholesterol Education Program (NCEP).¹³ Fifty-four samples, with LDL values ranging from 68.1 to 214.5 mg/dL, were tested in duplicate using the LDL reagent on a commercially available clinical chemistry analyzer and the Reference Method. Sample means were compared by least squares linear regression analysis. Data from this study are summarized below.

Method	LDL	Reference Method
N	54	54
Mean (mg/dL)	122.5	125.1
Standard Deviation (mg/dL)	30.7	30.9
Regression Analysis	$y = 0.95x + 3.02 \text{ mg/dL}$	
Correlation Coefficient	$r = 0.96$	

Method Comparison

The MULTIGENT Direct LDL assay on the AEROSET System was compared to a commercially available liquid selective detergent methodology. Samples were run in duplicate on both systems. Sample means were compared by least squares linear regression analysis. In a separate study, the MULTIGENT Direct LDL assay on an ARCHITECT

cSystem was compared to the assay on the AEROSET System. Data from both studies are summarized below.

Method	AEROSET vs. Comparative Method	ARCHITECT vs. AEROSET
N	49	57
Y-Intercept	-1.60	-1.09
Correlation Coefficient	0.99	0.99
Slope	1.02	0.98
%Bias	0.6%	-2.9%
Range (mg/dL)	33.0 to 182.1	68.5 to 214.0

BIBLIOGRAPHY

- Gotto AM. Lipoprotein metabolism and the etiology of hyperlipidemia. *Hosp Pract* 1988;23(Suppl 1):4.
- Crouse JR, Parks JS, Schey HM, et al. Studies of low density lipoprotein molecular weight in human beings with coronary artery disease. *J Lipid Res* 1985;26(5):566.
- Badimon JJ, Badimon L, Fuster V. Regression of atherosclerotic lesions by high-density lipoprotein plasma fraction in the cholesterol-fed rabbit. *J Clin Invest* 1990;85:1234-41.
- Castelli WP, Doyle JT, Gordon T, et al. HDL cholesterol and other lipids in coronary heart disease. *Circulation* 1977;55(5):767.
- Barr DP, Russ EM, Eder HA. Protein-lipid relationships in human plasma. *Am J Med* 1951;11(4):480.
- Gordon T, Castelli WP, Hjortland MC, et al. High density lipoprotein as a protective factor against coronary heart disease. *Am J Med* 1977;62(5):707.
- Williams P, Robinson D, Baily A. High density lipoprotein and coronary risk factor. *Lancet* 1979;1:72.
- Kannel WB, Castelli WP, Gordon T. Cholesterol in the prediction of atherosclerotic disease: new perspectives based on the Framingham study. *Ann Intern Med* 1979;90:85.
- US Department of Labor, Occupational Safety and Health Administration. 29 CFR Part 1910.1030. *Bloodborne Pathogens*.
- US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*, 5th ed. Washington, DC: US Government Printing Office, December 2009.
- World Health Organization. *Laboratory Biosafety Manual*, 3rd ed. Geneva: World Health Organization, 2004.
- Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. CLSI Document M29-A4. Wayne, PA: CLSI; 2014.
- National Cholesterol Education Program. *Detection, Evaluation, and Treatment of High Blood Cholesterol in Adults (Adult Treatment Panel III)*. National Institutes of Health. National Heart, Lung, and Blood Institute. NIH Publication No. 01-3670. May 2001.
- Calam RR, Bessman JD, Ernst DJ, et al. *Procedures for the Handling and Processing of Blood Specimens; Approved Guideline—Third Edition*. (H18-A3). Wayne, PA: Clinical and Laboratory Standards Institute, 2004.
- Westgard JO, Burnett RW, Cooper WG, et al. *Statistical Quality Control for Quantitative Measurements: Principles and Definitions; Approved Guideline (C24-A2)*. Wayne, PA: Clinical and Laboratory Standards Institute, 1999.
- Tietz NW, editor. *Clinical Guide to Laboratory Tests*. Philadelphia, PA: WB Saunders; 1986:256.
- Powers DM, Boyd JC, Glick MR, et al. *Interference Testing in Clinical Chemistry; Proposed Guideline (EP7-P)*. Villanova, PA: Clinical and Laboratory Standards Institute, 1986.
- Young DS. *Effects of Drugs on Clinical Laboratory Tests*, 3rd ed. Washington, DC: AACC Press; 1990:3-104-3-106.
- Carey RN, George H, Hartmann AE, et al. *User Demonstration of Performance for Precision and Accuracy; Proposed Guideline (EP15-P)*. Wayne, PA: Clinical and Laboratory Standards Institute, 1998.

TRADEMARKS

The ARCHITECT cSystem family of instruments consists of c4000, c8000, and c16000 instruments.

AEROSET, ARCHITECT, c4000, c8000, c16000, cSystem, and MULTIGENT are trademarks of Abbott Laboratories in various jurisdictions.

All other trademarks are property of their respective owners.

ARCHITECT cSystems Assay Parameters

Direct LDL Serum/Plasma—Conventional and SI Units

Configure assay parameters — General				
<input checked="" type="radio"/> General	<input type="radio"/> Calibration	<input type="radio"/> SmartWash	<input type="radio"/> Results	<input type="radio"/> Interpretation
Assay: DLDL	Type: Photometric	Version: †		
Number: 2840				
Run controls for onboard reagents by: Lot				
<input checked="" type="radio"/> Reaction definition	<input type="radio"/> Reagent / Sample	<input type="radio"/> Validity checks		
Reaction mode: End up				
Primary Secondary Read times				
Wavelength: 548 / 660 Main: 31 - 33				
Last required read: 33				
Absorbance range: ___ - ___ Color correction: ___ - ___				
Sample blank type: Self Blank: 14 - 16				

<input type="radio"/> Reaction definition	<input checked="" type="radio"/> Reagent / Sample	<input type="radio"/> Validity checks	
Reagent: DLDL0	Reagent volume: 200	R1	R2
Diluent: Saline	Water volume: ___		
Diluent dispense mode: Type 0	Dispense mode: Type 0	Type 0	Type 0
Dilution name	Sample	Diluted sample	Diluent
STANDARD	2.0	___	___
:	___	___	___
:	___	___	___
:	___	___	___
Dilution factor	1:1.00		
			●
			○
			○

<input type="radio"/> Reaction definition	<input type="radio"/> Reagent / Sample	<input checked="" type="radio"/> Validity checks
Reaction check: None		
Maximum absorbance variation: ___		

Configure assay parameters — Calibration				
<input type="radio"/> General	<input checked="" type="radio"/> Calibration	<input type="radio"/> SmartWash	<input type="radio"/> Results	<input type="radio"/> Interpretation
Assay: DLDL Calibration method: Linear				
<input checked="" type="radio"/> Calibrators	<input type="radio"/> Volumes	<input type="radio"/> Intervals	<input type="radio"/> Validity checks	
Calibrator set: LIPIDMCC	Blank: Water	Calibrator level: 0††	Concentration: ‡	
Replicates: 3 (Range 1–3)	Cal 1: LIPIDMCC1			

<input type="radio"/> Calibrators	<input checked="" type="radio"/> Volumes	<input type="radio"/> Intervals	<input type="radio"/> Validity checks	
Calibrator: LIPIDMCC	Calibrator level	Sample	Diluted sample	Diluent
Blank: Water	2.0	___	___	___
Cal 1: LIPIDMCC1	2.0	___	___	___

<input type="radio"/> Calibrators	<input type="radio"/> Volumes	<input checked="" type="radio"/> Intervals	<input type="radio"/> Validity checks	
Calibration intervals:				
Full interval: 672 (hours)				
Calibration type: None				
Adjust type: None				

<input type="radio"/> Calibrators	<input type="radio"/> Volumes	<input type="radio"/> Intervals	<input checked="" type="radio"/> Validity checks
Blank absorbance range: ___ - ___			
Span: Blank - Blank			
Span absorbance range: ___ - ___			
Expected cal factor: 0.00			
Expected cal factor tolerance %: 0			

Configure assay parameters — SmartWash				
<input type="radio"/> General	<input type="radio"/> Calibration	<input checked="" type="radio"/> SmartWash	<input type="radio"/> Results	<input type="radio"/> Interpretation
Assay: DLDL				
COMPONENT	REAGENT / ASSAY	WASH	Volume	Replicates

Direct LDL Serum/Plasma—Conventional Units

Configure assay parameters — Results				
<input type="radio"/> General	<input type="radio"/> Calibration	<input type="radio"/> SmartWash	<input checked="" type="radio"/> Results	<input type="radio"/> Interpretation
Assay: DLDL Assay number: 2840				
Dilution default range: Low-Linearity: 1				
High-Linearity: 800				
Gender and age specific ranges:**				
GENDER	AGE (UNITS)	NORMAL	EXTREME	
Either	0 - 130 (Y)	100 - 159		

Configure result units	
Assay: DLDL	Version: †
Result units: mg/dL	Decimal places: 0 (Range 0–4)
Correlation factor: 1.0000	Intercept: 0.0000

Direct LDL Serum/Plasma—SI Units

Configure assay parameters — Results				
<input type="radio"/> General	<input type="radio"/> Calibration	<input type="radio"/> SmartWash	<input checked="" type="radio"/> Results	<input type="radio"/> Interpretation
Assay: DLDL Assay number: 2840				
Dilution default range: Low-Linearity: 0.03				
High-Linearity: 20.69				
Gender and age specific ranges:**				
GENDER	AGE (UNITS)	NORMAL	EXTREME	
Either	0 - 130 (Y)	2.59 - 4.11		

Configure result units	
Assay: DLDL	Version: †
Result units: mmol/L	Decimal places: 2 (Range 0–4)
Correlation factor: 1.0000	Intercept: 0.0000

† Due to differences in instrument systems and unit configurations, version numbers may vary.

†† Displays the number of decimal places defined in the decimal places parameter field.

‡ Refer to the value on the calibrator labeling. These values are defined on the Configure calibrator set screen.

** User defined.

Key to Symbols

DISTRIBUTED BY

Distributed by

DISTRIBUTED IN THE USA BY

Distributed in the USA by

EC REP

Authorized Representative in the European Community

FOR USE WITH

Identifies products to be used together

INFORMATION FOR USA ONLY

Information needed for United States of America only

IVD

In Vitro Diagnostic Medical Device

LOT

Batch code/Lot number

PRODUCT OF CANADA

Product of Canada

R1

Reagent 1

R2

Reagent 2

REF

Catalog number/List number

SN

Serial number

WARNING: SENSITIZER

Warning: May cause an allergic reaction



Consult instructions for use



Manufacturer



Sufficient for



Temperature limitation



Use by/Expiration date

Customer Service: Contact your local representative or find country-specific contact information on www.abbottdiagnostics.com.



Sekisui Diagnostics P.E.I. Inc.
70 Watts Avenue
Charlottetown
Prince Edward Island
C1E 2B9 Canada

EC REP

MDSS GmbH
Schiffgraben 41,
30175 Hannover, Germany
Phone: (+49)-511-6262-8630
Fax: (+49)-511-6262-8633

DISTRIBUTED BY

Abbott Laboratories

DISTRIBUTED IN THE USA BY

Abbott Laboratories
Abbott Park, IL 60064 USA

February 2022
©2005, 2022 Abbott Laboratories

