INSTRUCTIONS FOR USE

Marathon™ Flow Directed Micro Catheter

This page is intentionally left blank.

TABLE OF CONTENTS

Marathon™ Flow Directed Micro Catheter

English5
Français6
Deutsch8
Italiano9
Español11
Svenska12
Nederlands14
Português15
Suomi17
Dansk18
Ελληνικά20
Česky21
Magyar23
Русский24
Polski26
Türkçe27
Norsk29
Slovenčina30
Română32
Български33
한국어35
عربي38
Symbol Glossary39

This page is intentionally left blank.

English

Instructions for Use

EN

Marathon™ Flow Directed Micro Catheter

CAUTION

Federal (USA) law restricts this device to sale, distribution, and use by or on the order of a physician. This device should be used only by physicians with a thorough understanding of angiography and percutaneous interventional procedures.

CONTENTS

- One Marathon™ Micro Catheter
- Shaping Mandrel

ACCESSORIES: (May be included with device model.)

- 0.007" Stylet or smaller
- 0.010" Micro Guidewire or smaller
- 1.5F Introducer Sheath or larger

DESCRIPTION

The Marathon™ Flow Directed Micro Catheter is a single-lumen, endhole catheter designed for the subselective infusion of physician-specified therapeutic agents such as embolization materials and diagnostic materials such as contrast media in tortuous, distal vessels. The catheter has a semi-rigid proximal shaft and a highly flexible distal shaft to facilitate the advancement of the catheter in the anatomy. The proximal end of the catheter incorporates a standard luer adapter to facilitate the attachment of accessories. The catheter has a radiopaque marker at the distal end to facilitate fluoroscopic visualization. The outer surfaces of the catheter are coated to increase lubricity. Microcatheter may be used with stylet, guidewire or introducer sheath to increase the rigidity of the distal section during introduction into the guiding catheter.

INDICATIONS FOR USE

The Marathon™ Flow Directed Micro Catheter is intended to access the peripheral and neuro vasculature for the controlled selective infusion of physician-specified therapeutic agents such as embolization materials and of diagnostic materials such as contrast media.

CONTRAINDICATIONS

- The Marathon™ Flow Directed Micro Catheter is contraindicated when, in the medical judgment of the physician, use of such product may compromise the patient's condition.
- Not intended for use in the coronary vasculature.
- The Marathon™ Flow Directed Micro Catheter is contraindicated for neonatal and pediatric use.

POTENTIAL COMPLICATIONS

Potential complications include, but are not limited to:

- Puncture site hematoma
- Vessel perforation
- Vessel spasm
- Hemorrhage
- Pain and tenderness
- Thrombolytic episodes
- Neurological deficits including stroke and death
- Vascular thrombosis

WARNINGS

- Infusion pressure with this device should not exceed 690 kPa/100 psi. Pressure in excess of 690 kPa/100 psi may result in catheter rupture, possibly resulting in patient injury.
- If flow through the catheter becomes restricted, do not attempt to clear the device by high pressure infusion. Either remove the catheter to determine the cause of the obstruction or replace it with a new catheter. Excessive pressure may cause catheter rupture, possibly resulting in patient injury.
- Never advance or withdraw an intraluminal device against resistance until the cause of the resistance is determined by fluoroscopy. Excessive force against resistance may result in damage to the device or vessel perforation.
- When injecting contrast for angiography, ensure that the catheter is not kinked, prolapsed or occluded. Remove excess slack in the catheter to reduce the potential for catheter kink or prolapse.
- Verify catheter integrity prior to re-inserting guidewire or injecting embolic material to prevent vascular damage or unintended embolization. Catheter integrity is verified by angiographically confirming that contrast agent is exiting only from the catheter tip while viewing the entire distal section of the catheter.
- Do not "over-shape" the catheter to achieve the desired shape angle. Over-shaping can lead to catheter kink or prolapse.
- This device is supplied STERILE for single use only. Do not reprocess or re-sterilize. Reprocessing and re-sterilization increase the risks of patient infection and compromised device performance.

- The stylet (supplied with some models) is not to be used as a guidewire. The stylet should never be manipulated within the catheter. **Do not advance the stylet beyond the tip of the catheter.** Use of the stylet as a guidewire could cause damage to the catheter and/or injure the patient.

PRECAUTIONS

- Prior to use, carefully examine the Marathon™ and its packaging to verify that they have not been damaged during shipment.
- Prior to use, all accessory devices and agents should be fully prepared according to the manufacturer's instructions.
- Always monitor infusion rates when using this catheter.
- The Marathon™ has a hydrophilic coating on the outside of the catheter which must be kept hydrated.
- This catheter is not intended for use with chemotherapy agents.
- When the infusion catheter is in the body, it should be manipulated only under fluoroscopy. Do not attempt to move the catheter without observing the resultant tip response.
- Because the Marathon™ may be easily advanced into narrow, selective vasculature, repeatedly verify that the Marathon™ has not been advanced so far as to interfere with its removal by slightly withdrawing the catheter before each infusion.
- When performing angiography, it is recommended to use a 3cc syringe rather than a 1cc syringe to reduce the risk of catheter over-pressurization.
- The Marathon™ is a flow directed micro catheter that can optionally be used with hydrophilic, 0.010" or less sized guidewires. The Marathon™ is not compatible with non-hydrophilically coated guidewires greater than 0.010" in diameter.
- It is recommended that the Marathon™ be used with an appropriately sized guiding catheter which allows adequate clearance (minimum internal diameter of 0.053" or 1.35mm).

STORAGE

Store the Marathon™ Flow Directed Micro Catheter in a dry place between 50°F (10°C) and 90°F (32°C).

DIRECTIONS FOR USE

Marathon™ Flow Directed Micro Catheter

1. Place the appropriate guiding catheter following recommended procedures. Connect a one-way stop-cock to the guiding catheter to prevent backflow of blood during insertion of catheter. Connect a hemostatic side arm adapter one-way stop-cock to allow for insertion of catheter and to continuously flush the guiding catheter with saline.
NOTE: It is recommended that a guiding catheter with minimum internal diameter of 0.053" (1.35mm) be used with the Marathon™ Flow Directed Micro Catheter.
2. Gently remove the packaging coil for the Marathon™ Flow Directed Micro Catheter from the package.
3. Flush the coil with saline through the female luer attached to the coil.
4. Remove the catheter by removing the hub from the clip and gently pulling on the hub.
5. Inspect the catheter prior to use to verify that it is undamaged. Retain the packaging coil for storage of the catheter when it is not in use during the procedure.
6. Flush the ID of the micro catheter coil with saline.
7. To assist with insertion of the microcatheter in to the hub, the following accessories may be used: a guidewire, stylet, and/or introducer sheath.
8. If using a stylet and/ or guidewire:
 - Prepare stylet or guidewire per manufacture's instructions; product should be hydrated.
 - Remove the stylet or guidewire from its packaging coil
 - Insert the stylet or guidewire through the hub of the Micro Catheter.
 - Advance the stylet or guidewire through the catheter, keeping the catheter body straight. **Do not advance the stylet past the tip of the catheter.**
 - Close the one-way stop-cock.
 - Loosen the hemostatic valve.
 - Introduce the catheter through the hemostatic adapter.
 - Tighten the valve around the catheter to prevent backflow, but allowing some movement through the valve by the catheter.
 - Open the one-way stop-cock.
 - Advance the catheter such that the stylet or guidewire remains in the guiding catheter. **Do not advance the stylet past the tip of the guiding catheter.**
 - Once the catheter has been advanced to the end of the guiding catheter, if using a stylet withdraw the stylet from the catheter.
9. The catheter may be advanced through the vasculature by gently pushing the proximal shaft. It is recommended to use a guidewire during navigation to reduce the risk of catheter kink or prolapse.
WARNING: Verify catheter integrity prior to re-inserting guidewire or injecting embolic material to prevent vascular damaged or unintended embolization. Catheter integrity is verified by angiographically confirming that contrast agent is exiting only from the catheter tip while viewing the entire distal section of the catheter.

10. To infuse, connect a syringe with infusate to the catheter luer, and infuse as required.

Functional Characteristics

Usable Length	Minimum Dead Space Volume	Maximum Embolization Particle Size	Approximate Infusion Rate at 100 psi (690 kPa)	
			Water	Contrast (76% Renografin)
165 cm	0.23 ml	250 µ	12 ml/min	3 ml/min

Note: Embolization particle concentration and injection rate should be verified by the physician prior to use. High concentrations of particles have a tendency to aggregate in the catheter hub.

STEAM SHAPING MANDREL

In order to maintain catheter integrity and dimensional stability of the inner diameter, it is recommended that the user follow these instructions.

WARNINGS

Shaping Mandrel is not intended for use in the human body.

Use only a steam source to shape the catheter tip. Do not use other heat sources.

Prior to use, inspect the catheter tip for any damage that may have resulted from shaping. Do not use a catheter that has been damaged in any way. Damaged catheters may rupture causing vessel trauma or tip detachment during steering maneuvers.

1. Remove shaping mandrel from card and insert into distal tip of the catheter.
2. Carefully bend catheter tip and shaping mandrel to desired shape. A slight over exaggeration of the shape may be required to accommodate for catheter relaxation.
3. Shape the catheter by holding the shaped portion approximately 1 inch (2.5 cm not less than 1 cm) from the steam source for about 20 seconds (DO NOT EXCEED 30 SECONDS).
4. Allow catheter tip to cool in air or saline prior to removing mandrel.
5. Remove mandrel from catheter and discard. Multiple shaping is not recommended.
6. Inspect the tip for any damage that may have resulted from steam shaping the catheter. If any damage is found, do not use the catheter.

PRECAUTIONS

If catheter entrapment is suspected (with any embolic agent), fast catheter retrieval technique may result in catheter shaft separation and potential vascular damage.

1. Slowly remove any 'slack' in the distal catheter shaft.
2. Gently and slowly apply 3-5 cm traction to the catheter to begin catheter retrieval.
3. Should catheter removal become difficult, the following will assist in catheter retrieval. Assess the following parameters by observing the distal shaft of the catheter:
 - Vessel straightening
 - Traction on embolic cast
 - Catheter tip releasing from embolic cast
4. Further traction (of 3-5 cm) may be applied gently if necessary to transfer distal catheter shaft.
 - Hold this traction for a few seconds and release
 - Assess traction of vasculature to minimize risk of hemorrhage
5. This process can be repeated intermittently until catheter is retrieved.

NOTE: Do not apply more than 20 cm of traction to catheter to minimize risk of catheter separation.

6. Under some difficult clinical situations, it may be safer to leave a flow-directed catheter in the vascular system rather than risk rupturing the malformation and, consequently a hemorrhage, by exercising too much traction on an entrapped catheter.

This is accomplished by stretching the catheter and cutting the shaft near the entry point of vascular access allowing the catheter to remain in the artery.

If catheter breaks during removal, distal migration or coiling of the catheter may occur. Same day surgical resection should be considered to minimize risk of thrombosis.

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Micro Therapeutics, Inc.