# Medtronic

## Onyx TruCor™

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

## 1 Device description

The Onyx TruCor Zotarolimus-Eluting Coronary Stent System (Onyx TruCor stent system) consists of 4 subsystems:

1. Onyx TruCor bare metal stent – a premounted cobalt alloy and platinum-iridium alloy based stent
2. Delivery system – Onyx TruCor rapid exchange (RX) delivery system
3. Polymer system
4. Zotarolimus – drug

A flushing cannula and a looper device are included in the package. See *Figure 2*.

### 1.1 Stent

The Onyx TruCor stent (*Figure 1*) is manufactured from a composite material of cobalt alloy and platinum-iridium alloy. The coronary stent is formed from a single wire bent into a continuous sinusoid pattern and then laser fused back onto itself. The stents are provided in multiple lengths and diameters.

### 1.2 Delivery system

The delivery system consists of a balloon-expandable intracoronary stent premounted on the Onyx TruCor RX delivery system. The effective working length of the catheter is 140 cm. The delivery system is compatible with 0.36 mm (0.014 in) maximum outer diameter guidewires and 1.42 mm (5 Fr/0.056 in) minimum inner diameter guide catheters.

### 1.3 Polymer system

The Onyx TruCor stent is composed of a bare metal stent with a primer coat and a coating that consists of a blend of the drug zotarolimus and the BioLinx™ polymer system.

### 1.4 Drug — Zotarolimus

The drug zotarolimus is a proprietary chemical entity. Zotarolimus is a tetrazole-containing macrocyclic drug. The suggested mechanism of action of zotarolimus is to bind to the cytoplasmic protein FKBP-12, leading to the formation of a trimeric complex with the protein kinase mTOR (mammalian target of rapamycin), inhibiting its activity. Inhibition of mTOR results in the inhibition of protein phosphorylation events associated with translation of mRNA and cell cycle control. The Onyx TruCor stent with the active drug zotarolimus is intended to reduce the incidence of restenosis in coronary interventions. The Onyx TruCor stent has a nominal drug dose of 1.6 µg zotarolimus per mm<sup>2</sup> of the stent surface area.

## 2 Intended purpose

The Onyx TruCor stent is intended to improve coronary luminal diameters of either single or multiple vessels, as an adjunct to coronary interventions and to reduce restenosis. The stent is intended as a permanently implanted device.

### 2.1 Indications

The Onyx TruCor stent system is intended for use in patients eligible for percutaneous transluminal coronary angioplasty (PTCA) with a reference vessel diameter of 2.25 mm to 4.0 mm. The Onyx TruCor stent is indicated for the treatment of the following patient and lesion subsets:

- Diabetes mellitus
- Multivessel disease
- Acute coronary syndrome (ACS)
- Acute myocardial infarction (AMI)
- Unstable angina (UA)
- Bifurcation lesions
- In-stent restenosis (ISR)
- Chronic total occlusions (CTO)
- Total occlusions (TO)
- Left main (LM)
- Small vessel (SV)
- One month of dual antiplatelet therapy (DAPT) in high bleeding risk (HBR) patients, including patients who are unable to tolerate long-term DAPT

### 2.2 Contraindications

The Onyx TruCor stent is contraindicated in the following cases or conditions:

- Patients with hypersensitivity or allergies to aspirin, heparin, clopidogrel, ticlopidine, mTOR inhibiting drugs such as zotarolimus (tacrolimus, sirolimus, everolimus) or any other mTOR inhibitor analogue or derivative, polymers, cobalt, chromium, nickel, molybdenum, platinum, iridium, or contrast media.
- Patients in whom antiplatelet or anticoagulation therapy is contraindicated.
- Patients who are judged to have a lesion that prevents complete inflation of an angioplasty balloon or proper placement of the stent or stent delivery system.

### 2.3 Intended patient population

The Onyx TruCor stent system is intended for use in patients eligible for percutaneous transluminal coronary angioplasty (PTCA) with a reference vessel diameter of 2.25 mm to 4.0 mm.

### 2.4 Intended users

Only physicians who have received appropriate training should implant the stent.

### 2.5 Clinical benefits

The Onyx TruCor stent system is used for percutaneous transluminal coronary angioplasty in the coronary arteries for treating de novo and recurrent stenoses. The primary clinical benefits include the following: maintaining vessel patency post stent

implantation and reduction in stenoses recurrence, resulting in decreased repeat revascularization (a therapeutic effect of the drug zotarolimus); improvement of myocardial perfusion, thus alleviating related anginal symptoms. The measurable patient-relevant outcomes and quantitative clinical data from Medtronic-sponsored Resolute Onyx clinical trials are shown in *Table 1*. Onyx TruCor uses the same stent system as Resolute Onyx.

**Table 1.** Clinical benefits

Patient-relevant outcomes	RESOLUTE ONYX <sup>a</sup> Post-Approval Study (N=416 subjects) (N=583 lesions)	Onyx ONE <sup>b</sup> (N=1003 subjects) (N=1282 lesions)
	Safety measures	
Clinically driven target lesion revascularization (TLR) at 12 months <sup>c</sup>	3.2% (13/412)	2.8% (28/988)
Clinically driven target vessel revascularization (TVR) at 12 months <sup>d</sup>	3.4% (14/412)	3.6% (36/988)
Effectiveness measures		
Lesion success <sup>e</sup>	98.8% (568/575)	93.8% (1170/1248)
Device success <sup>f</sup>	97.4% (558/573)	92.8% (1158/1248)

<sup>a</sup> Study included data from various stent diameters, including 2.0 mm.

<sup>b</sup> Data represents Resolute Onyx Arm. Study included data from various stent diameters, including 2.0 mm, 4.5 mm, and 5.0 mm.

<sup>c</sup> Revascularization at the target lesion associated with positive functional ischemia study or ischemic symptoms **and** an angiographic minimal lumen diameter stenosis  $\geq 50\%$  by quantitative coronary angiography (QCA), or revascularization of a target lesion with diameter stenosis  $\geq 70\%$  by QCA without either angina or a positive functional study.

<sup>d</sup> Revascularization in the target vessel associated with positive functional ischemia study or ischemic symptoms **and** an angiographic minimal lumen diameter stenosis  $\geq 50\%$  by QCA, or revascularization of a target vessel with diameter stenosis  $\geq 70\%$  by QCA without either angina or a positive functional study.

<sup>e</sup> The attainment of  $<30\%$  residual stenosis by QCA (or  $<20\%$  by visual assessment) **and** thrombolysis in myocardial infarction (TIMI) flow 3 after the procedure, using any percutaneous method.

<sup>f</sup> The attainment of  $<30\%$  residual stenosis by QCA (or  $<20\%$  by visual assessment) **and** TIMI flow 3 after the procedure, using the assigned device only.

## 2.6 Summary of safety and clinical performance

The Summary of Safety and Clinical Performance (SSCP) for this device can be found at <https://ec.europa.eu/tools/eudamed> using the Basic UDI-DI: 0763000B00000156T.

## 2.7 Performance characteristics of the device

The Onyx TruCor stent system is a drug-eluting, balloon expandable, coronary stent offered on a rapid exchange (RX) delivery system. The Onyx TruCor stent dilates the vessel lumen to improve myocardial perfusion and act as a scaffold to prevent acute recoil. The Onyx TruCor stent with the active drug zotarolimus is intended to reduce the incidence of restenosis in coronary interventions.

## 3 Warnings and precautions

- The patient's exposure to the drug and polymer system is directly related to the number of stents and stent length implanted. Based on zotarolimus pharmacokinetic clinical data, the maximum number of stents implanted in a patient should be 78 or fewer for the Onyx TruCor stent with the highest dose (a 38 mm stent length with a nominal diameter of 3.50 mm or 4.00 mm). This maximum is based on the relative exposure of zotarolimus observed in healthy patients from intravenous doses of zotarolimus compared to the exposure in a patient after implantation of a Onyx TruCor stent. See *Table 4* for nominal doses.
- Judicious selection of patients is necessary since the use of this device carries the associated risk of complications listed in *Chapter 4*. A risk/benefit assessment should be considered for the patient before implantation.
- Administration of appropriate anticoagulant, antiplatelet, and coronary vasodilator therapy is critical to successful stent implantation. Please see the most current medical society guidelines.
- Only physicians who have received appropriate training should implant the stent.
- Stents should be implanted only in hospitals that offer a possibility to facilitate coronary artery bypass surgery.
- Subsequent restenosis may require repeat dilatation of the arterial segment containing the stent. The long-term outcome following repeat dilatation of endothelialized coronary stents is currently unknown.
- While no specific clinical data are available, drugs like tacrolimus that act through the same binding protein (FKBP) may interfere with the efficacy of zotarolimus. Zotarolimus is metabolized by CYP3A4, a human cytochrome P450 enzyme. Strong inhibitors of CYP3A4 (such as ketoconazole) may cause increased zotarolimus exposure to levels associated with systemic effects, especially if multiple stents are deployed. Systemic exposure to zotarolimus should also be taken into consideration if the patient is treated concomitantly with systemic immunosuppressive therapy.
- There are no adequate or well-controlled studies in pregnant women, lactating women, pediatric patients, or men intending to father children for this product. Animal fertility studies with intravenous (IV) administered zotarolimus demonstrated embryo-toxicity including embryoletality. It is not recommended that the Onyx TruCor stent be used in women attempting to conceive, or who are pregnant or lactating. The IV administered zotarolimus studies in male rats resulted in reversible, dose-related testicular toxicity.
- When multiple stents are required, stent materials should be of similar composition. Placing multiple stents of different materials in contact with each other may increase the potential for corrosion. Data obtained from **in vitro** corrosion tests using a cobalt alloy stent (Medtronic Onyx stent platform) in combination with a stainless steel alloy stent (Boston Scientific Liberte™) do not suggest an increased risk of **in vivo** corrosion.

- Potential interactions of the Onyx TruCor stent with other drug-eluting stents have not been evaluated and should be avoided whenever possible.
- The safety and effectiveness for stenting of saphenous vein grafts with the Onyx TruCor stent has not been established.
- Implanting may lead to dissection of the vessel distal or proximal to the stented portion, and may cause acute closure of the vessel requiring additional intervention (such as CABG [coronary artery bypass graft], secondary balloon dilatation, or placement of additional stents).
- Placement of a stent has the potential to compromise side branch patency.
- The Onyx TruCor stent must not be exposed to any direct handling (such as rolling of the stent) or contact with liquids prior to preparation and delivery, as the coating may be susceptible to damage or premature drug elution.
- Do not expose or wipe the Onyx TruCor stent system with organic solvents such as alcohol or detergents.
- Do not bend or kink the hypotube. If the delivery system shaft or hypotube becomes bent or kinked, do not attempt to straighten them. Straightening a kinked metal shaft may result in breakage of the shaft.
- Do not remove the stent from the stent delivery system, as removal may damage the stent and polymer system and may lead to stent embolization. The stent is not designed to be crimped onto another delivery device.
- When the Onyx TruCor stent is newly deployed, take care when crossing the stent with an intravascular ultrasound (IVUS) catheter, a coronary guidewire, a balloon catheter, or any other device. Crossing a newly deployed stent could disrupt the stent placement, apposition, coating, or the stent geometry of the Onyx TruCor stent.
- The Onyx TruCor RX delivery system should not be used in conjunction with any other stents or for post-dilatation.
- Early discontinuation of prescribed antiplatelet medication could result in a higher risk of thrombosis, myocardial infarction, or death.
- If a serious incident related to the device occurs, immediately report the incident to Medtronic and the applicable competent authority or regulatory body.
- Dispose of the delivery system and components in accordance with applicable laws, regulations, and hospital procedures, including those regarding biohazards, microbial hazards, and infectious substances.
- This device contains a hazardous substance. Cobalt (CAS No. 7440-48-4) is classified as a carcinogenic, mutagenic, or toxic for reproduction (CMR) substance of Category 1A or 1B, or as an endocrine disruptor (ED), and is present at a concentration >0.1% weight/weight.

#### 4 Potential adverse events

The following potential adverse events are associated with the use of coronary stenting devices, IVUS (intravascular ultrasound), or PCI (percutaneous coronary intervention):

- Death
- Aneurysm, pseudoaneurysm, or arteriovenous fistula (AVF)
- Stent deformation, collapse, or fracture
- Emergency surgery: peripheral vascular or coronary bypass
- Stroke or transient ischemic attack (TIA)
- Cardiac tamponade
- Hemorrhage requiring transfusion
- Shock or pulmonary edema
- Coronary artery spasm
- Abrupt vessel closure
- Hypotension or hypertension
- Allergic reaction (to contrast, antiplatelet therapy, stent system —material, drug, or polymer coating)
- Peripheral ischemia or peripheral nerve injury
- Coronary artery occlusion, perforation, rupture, or dissection
- Pericarditis
- Embolism (air, tissue, device, or thrombus)
- Thrombosis (acute, subacute, late, or very late)
- Incomplete stent apposition
- Myocardial infarction (MI)
- Restenosis of the stented artery
- Arrhythmia
- Infection or fever
- Unstable angina
- Access site pain, hematoma, or hemorrhage
- Balloon rupture
- Stent migration
- Failure to deliver the stent
- Stent misplacement

The occurrence of the preceding complications may lead to myocardial infarction, emergency bypass surgery, or death. Complications may alternatively lead to the need for a repeat catheterization or percutaneous coronary intervention. The following additional side effects and complications may be associated with, but are not limited to, the use of zotarolimus:

- Anemia
- Circumoral paresthesia
- Diarrhea
- Dry skin
- Headache
- Hematuria
- Infection
- Pain (abdominal or arthralgia)
- Rash

The side effects and complications of the BioLinx™\* polymer are no different than those of other stent coatings and may include, but are not limited to, the following conditions:

- Focal inflammation at the site of stent implantation
- Restenosis of the stented artery
- Allergic reaction

#### 5 Individualization of treatment

The risks and benefits previously described (*Chapter 4 and Section 2.5*) should be carefully considered for each patient before use of the Onyx TruCor stent system. Patient selection factors to be assessed should include a judgement regarding the risk of prolonged anticoagulation. Stenting should be generally avoided in those patients at heightened risk of bleeding (such as those patients with recently active gastritis or peptic ulcer disease; refer to *Section 2.2*). Each Onyx TruCor stent must be sized

appropriately to fit the patient's particular coronary anatomy. Proper sizing of the device (diameter and length) is the responsibility of the physician.

Co-morbidities that increase the risk of poor initial results in PCI or increase the risks of emergency referral for bypass surgery should be reviewed.

While physicians should adhere to current European Society of Cardiology (ESC), or American College of Cardiology (ACC), American Heart Association (AHA), or Society for Cardiovascular Angiography and Interventions (SCAI) guidelines for PCI, patients who interrupt or discontinue DAPT medication one month or more after stent implantation are considered at low risk and showed no increased risk for stent thrombosis. In patients with HBR or who are unable to tolerate long term DAPT, physicians may choose a one month dual antiplatelet regimen. This dual antiplatelet regimen is based on results of the randomized single blind Onyx ONE Global Randomized Clinical Trial (Onyx ONE RCT). The Onyx ONE RCT demonstrated that the Onyx TruCor stent is safe and effective in complex HBR patients who receive one month DAPT. In addition, the Onyx TruCor stent demonstrated low event rates in landmark analysis after discontinuing DAPT.

#### 6 Patient counseling information

In accordance with local regulations, healthcare providers should review the instructions for use for applicable information to be shared with the patient. A patient implant card, which contains identifying information about the implanted device, is included in the device package. After device implant, complete the patient implant card and provide it to the patient before they are discharged.

Healthcare providers should communicate the following instructions to their patients:

- Always carry their implant card with them.
- Access additional information about their device on the website that is listed on their patient implant card.  
**Note:** If the patient is unable to access the website, the healthcare provider must provide the information from the website to the patient.
- Always inform any healthcare personnel that they have an implanted device before any procedure has begun.
- Contact their healthcare provider if they notice any new or changing symptoms.

#### 7 How supplied

The device package contains one coronary stent pre-mounted on a custom stent delivery system. A loopier device and flushing cannula (Figure 2, items 2 and 3) are included in the package. These components are attached to a clip on the hoop that contains the Onyx TruCor stent system. Do not resterilize. The device was sterilized by EtO gas diffusion.

The device is supplied sterile, non-pyrogenic, in an unopened, undamaged package. Use by the "Use-by date" noted on the package. If any information on the outer package or the sterile package is defaced or damaged, do not use and notify Medtronic so that the device can be replaced. If any part of this manual is illegible, contact Medtronic to request a replacement manual.

This device is for single use. This device is intended to contact body tissues. Do not reuse, reprocess, or resterilize. Reprocessing may compromise the structural integrity of this device. Reuse of this device creates a potential risk of patient infections due to contamination. This contamination of the device may lead to injury, illness, or death of the patient.

**Storage:** Store in the original container. Store between 15°C and 30°C.

A sealed pocket embedded inside the pouch contains 2 small packets (oxygen scavenger used to remove oxygen from the pouch and desiccant used to remove moisture).

**Note:** This pocket is not intended for opening as the internal surface of the pocket is not sterile.

#### 8 Device specifications and materials

**Table 2.** Specifications

Description	Specification
Stent diameter	2.25 mm to 4.00 mm
Stent length	8, 12, 15, 18, 22, 26, 30, 34, 38 mm
Guidewire compatibility	0.36 mm (0.014 in) maximum outer diameter
Guide catheter compatibility	1.42 mm / 5.0 Fr (0.056 in) minimum inner diameter
Nominal drug dosage	1.6 µg per mm <sup>2</sup> of stent surface area <sup>a</sup>
Nominal pressure <sup>b</sup>	12 atm for all diameters
Rated burst pressure	18 atm for all diameters

<sup>a</sup> Value represents a range of material amounts that may be present. Specific amount depends on the size of the implanted stent.

<sup>b</sup> Refer to the product label for complete inflation pressure recommendations.

**Table 3.** Materials and substances to which patients can be exposed from a single stent

Description	Specifications (mass)
Cobalt alloy <sup>a,b</sup>	6.4 to 40.4 mg <sup>f</sup>
Platinum iridium alloy <sup>c</sup>	1.8 to 11.5 mg <sup>f</sup>
Parylene C primer coating	56 to 348 µg <sup>f</sup>
Zotarolimus and BioLinx™* <sup>d</sup> coating <sup>e</sup>	150 to 933 µg <sup>f</sup>

<sup>a</sup> Cobalt alloy, conforming to ASTM F562. Cobalt is classified as a CMR substance. Refer to *Chapter 3*.

<sup>b</sup> Contains nickel, a known substance that can cause sensitization or an allergic reaction

<sup>c</sup> 90% Platinum – 10% Iridium alloy, conforming to ASTM B684

<sup>d</sup> A blend of the Medtronic proprietary components C10 and C19, and PVP (polyvinyl pyrrolidone). The ratios are as follows: 10% PVP, 27% C10, and 63% C19

<sup>e</sup> Drug coating: 35% zotarolimus drug and 65% BioLinX™\*

<sup>f</sup> Values represent a range of material amounts that may be present. Specific amount depends on size of implanted stent

**Table 4.** Nominal drug dose for each nominal stent size

Diameter (mm)	Zotarolimus (µg)								
	Length (mm)								
	8	12	15	18	22	26	30	34	38
2.25	51	70	85	104	127	146	168	187	206
2.50	51	70	85	104	127	146	168	187	206
2.75	67	94	117	140	171	198	225	257	284
3.00	67	94	117	140	171	198	225	257	284
3.50	77	108	132	156	186	221	252	282	317
4.00	77	108	132	156	186	221	252	282	317

## 9 Directions for use

Prepare the guide catheter and guidewire according to the manufacturer's instructions. The Onyx TruCor stent is compatible with 0.36 mm (0.014 in) guidewires. Refer to the product labeling or *Section 9.1* for specific guide catheter compatibility.

### 9.1 Materials required

The following materials are required for use of this device:

- 1.42 mm (5 Fr/0.056 in) minimum inner diameter guide catheter
- 20 cc syringe
- Heparinized normal saline
- 0.36 mm (0.014 in) maximum outer diameter guidewire
- Rotating hemostatic valve
- Contrast medium diluted 1:1 with heparinized normal saline
- Inflation device
- Torque device
- Three-way stopcock

### 9.2 Selection of stent size

Careful stent sizing is important to successful stenting. In general, the stent size should be chosen to match the diameter of the reference vessel. Slight stent oversizing is preferable to undersizing. Ensure the stent area is long enough to cover the lesion completely and should be ≥3 mm longer than the lesion length.

The inflated balloon diameter measures slightly larger than the labeled stent diameter to allow for stent recoil upon expansion.

### 9.3 Preparation of the delivery system

1. Remove the stent delivery system from the package. Special care must be taken not to handle the stent or in any way disrupt its placement on the balloon. This is most important during catheter removal from packaging, placement over guidewire, and advancement through the rotating hemostatic valve and guide catheter hub.
2. Carefully slide the protective sheath off the stent by grasping the sheath at the distal end with the thumb and forefinger and gently pulling it. **Do not** touch the part of the sheath over or proximal to the stent. Removing the protective sheath will also remove the stylette. Inspect the stent to assure it has not been damaged or displaced from its original position on the balloon. Should there be movement of or damage to the stent, do not use.
3. Flush the balloon catheter guidewire lumen with heparinized saline.
  - a. Remove the flushing cannula (*Figure 2*, item 3) from the component clip and without removing the cover, attach it directly to a syringe.
  - b. Remove the cover of the flushing cannula.
  - c. Insert the flushing cannula into the guidewire entry port (exchange joint).
  - d. Depress the syringe to flush the guidewire lumen until fluid emerges from the distal tip of the balloon catheter.
  - e. Remove the flushing cannula from the guidewire entry port (exchange joint).
4. Prepare the balloon lumen with 50/50 contrast-saline mixture as follows:
  - Do not attempt a pre-inflation technique to purge the balloon lumen.
  - Do not use air or any gaseous medium to inflate the balloon.
  - Do not pull negative pressure on the inflation device before beginning the preparation step.
  - a. Using a 20-cc (mL) syringe containing 5-cc (mL) of contrast-saline mixture, apply negative pressure for 20 to 30 seconds, allowing air removal from the balloon. An excessive amount of air released into the syringe or no air released from the balloon may indicate damage to the stent delivery system.
 

Do not put bending force on the stent delivery system when applying negative pressure with the syringe. Should there be an indication of damage to the stent delivery system, do not use.
  - b. Release the pressure, slowly allowing negative pressure to draw the mixture into the balloon lumen.
 

Do not apply negative pressure on the inflation device after balloon preparation and prior to delivering the stent.
  - c. Detach the syringe, leaving a meniscus of mixture on the hub of the balloon lumen.
  - d. Prepare the inflation device in the standard manner, and purge to remove all air from the syringe and tubing.

e. Attach the inflation device to the balloon lumen directly. Apply the "meniscus to meniscus" technique to ensure that no air bubbles remain at the connection. Leave at ambient pressure.

5. Visually inspect the stent delivery system to ensure that the stent is placed within the area between the proximal and distal balloon markers.

Do not wipe with gauze sponges, as the fibers may disrupt the stent.

### 9.4 Delivery procedure

1. Prepare the vascular access site according to standard PCI practice.
2. The decision to predilate the lesion with an appropriately sized balloon should be based on patient and lesion characteristics. If this is performed, predilate the lesion with a balloon diameter 0.5 mm smaller than the stent and a balloon length equal to or shorter than the target lesion length. The length of the predilatation balloon must be shorter than the stent to be implanted.
3. Maintain neutral pressure on the inflation device. Open the rotating hemostatic valve to allow for easy passage of the stent. Make sure the rotating hemostatic valve has a large bore and is fully open while passing the stent through it. If resistance is encountered, do not force passage. Resistance may indicate damage to the stent or the stent delivery system.
4. Advance the stent delivery system along the guidewire through the large bore rotating hemostatic valve using conventional angioplasty techniques.
5. Carefully advance the stent delivery system into the hub of the guide catheter.
 

The stent delivery system can now be advanced through the guide catheter.
6. Ensure guide catheter stability before advancing the stent delivery system into the coronary artery.
7. Sufficiently tighten the rotating hemostatic valve. The stent is now ready to be deployed.

### 9.5 Deployment of the Onyx TruCor Stent

1. Determine the reference vessel diameter prior to stent selection (refer to *Section 9.2*).
2. Ensure guide catheter stability before advancing the balloon into the coronary artery.
 

If initial guide catheter position is lost, avoid pulling or pushing the guide catheter over the stent. If this is done, the distal end of the guide catheter may damage the stent. If the stent delivery system does not readily advance, do not force. If the stent will not advance in spite of good guide catheter support, consider dilating the proximal obstructing plaque.
3. Advance the delivery system over the guidewire to the target lesion under direct fluoroscopic visualization. Position the stent across the lesion using the proximal and distal radiopaque markers on the balloon as reference points. Optimal placement requires the stent to be ≥3 mm longer than the lesion length and completely cover the area to be stented. Expansion of the stent should not be undertaken if the stent is not properly positioned in the stenotic segment of the vessel. If the position of the stent is not optimal, it should be repositioned or removed. (Refer to *Section 9.7*).
4. Prior to stent expansion, utilize high resolution fluoroscopy to verify with the stent outline that the stent has not been damaged or shifted during positioning.
5. Inflate the balloon to the nominal pressure to expand the stent. Refer to the compliance chart for the proper inflation pressure. A 15 to 30 second inflation is recommended for full expansion. Balloon pressure should be monitored during inflation.
 

Do not exceed the rated burst pressure as indicated on the compliance chart. Use of a higher pressure range than specified on the compliance chart may result in a ruptured balloon or oversizing of the stent with possible intimal or vessel damage.

Fluoroscopic visualization during stent expansion should be used in order to properly judge the optimum stent diameter, as compared to the proximal and distal native coronary artery diameters. Optimal expansion and proper sizing requires that the stent be in full contact with the arterial wall. The inflated balloon diameter measures slightly larger than the labeled deployed stent diameter to allow for stent recoil after expansion and balloon deflation.

Under-expansion of the stent may result in stent migration.

Oversizing of the stent and use of higher than recommended inflation pressures may cause vessel dissection. It is recommended that the chosen stent size should closely approximate the diameter of the vessel and that recommended stent inflation pressures be used for stent deployment.

If the target lesion is incompletely stented, use additional stents as necessary to adequately treat the lesion. When additional stents are necessary, the distal lesion should be initially stented, followed by stenting of the proximal lesion. Stenting in this order obviates the need to cross the proximal stent with a stent delivery system and reduces the chances of dislodging the proximal stent.
6. Deflate the balloon by applying negative pressure and allowing adequate time (20 to 30 seconds) for full balloon deflation. Longer stents may require more time for deflation.
 

Deflation of the balloon should be confirmed by the absence of contrast within the balloon.
7. Very slowly withdraw the balloon from the stent, maintaining negative pressure and allowing movement of the myocardium to gently dislodge the balloon from the stent. If any resistance is felt when removing the balloon from the stent, place the deflator on neutral pressure and gently remove the balloon. Maintain the position of the guide catheter to prevent it from being drawn into the vessel.

### 9.6 Further dilatation of stented segments

If the deployed stent size is inadequate with respect to the vessel diameter, a larger balloon may be used to further expand the stent to its optimal size. If the initial angiographic results are suboptimal, the stent may be further deployed using a low profile, high pressure, non-compliant balloon catheter. If required, the stented segment should be recrossed carefully with a prolapsed guidewire to avoid dislodging the stent. All efforts should be taken to assure that the stent is not under dilated.

Nominal stent diameter (mm)	Maximum stent inner diameter (mm)
2.25 and 2.50	3.50
2.75 and 3.00	4.00
3.50 and 4.00	5.00

### 9.7 Removal of an unexpanded stent

Failure to follow these steps or applying excessive force to the stent delivery system can potentially result in loss or damage to the stent or stent delivery system components, such as the balloon.

- If removal of a stent system is required prior to deployment, ensure that the guide catheter is coaxially positioned relative to the stent system and cautiously withdraw the stent system into the guide catheter.
- Should unusual resistance be felt at any time when withdrawing the stent towards the guide catheter, the stent delivery system and the guide catheter should be removed as a single unit. This must be done under direct visualization with fluoroscopy.
- The guide catheter and the stent delivery system should be carefully removed from the coronary artery as a single unit. When removing the stent delivery system and guide catheter as a single unit:
  - Do not retract the stent delivery system into the guide catheter.
  - Maintain guidewire placement across the lesion and carefully pull back the stent delivery system until the proximal balloon marker of the stent delivery system is aligned with the distal tip of the guide catheter.
  - The system should be pulled back into the descending aorta toward the arterial sheath. As the distal end of the guide catheter enters into the arterial sheath, the catheter will straighten, allowing safe withdrawal of the stent delivery system into the guide catheter and the subsequent removal of the stent delivery system and the guide catheter from the arterial sheath.
- Inspect the stent delivery system post removal to ensure that the stent is present on the stent delivery system.
- Dispose of the device in accordance with applicable laws, regulations, and hospital procedures, including those regarding biohazards, microbial hazards, and infectious substances.

### 9.8 Instructions for simultaneous use of two devices in guide catheter (Kissing Balloon Technique)

- 6 Fr (2 mm) Compatibility—Any combination of one Onyx TruCor RX stent delivery system (2.25 mm to 4.00 mm) and one balloon catheter (Sprinter Legend RX 1.25 mm to 3.50 mm in diameter, Euphora RX 1.50 mm to 3.50 mm in diameter, or NC Euphora RX 2.00 mm to 3.50 mm in diameter) can be used simultaneously within a 6 Fr (2 mm)/GC/MID 1.8 mm (0.070 in) guide catheter.

The technique can be performed per the following instructions:

1. Insert the Onyx TruCor RX stent using the instructions provided (refer to *Section 9.4*).
2. Insert a second guidewire and a balloon catheter, track to the target site, and inflate the balloon.
3. To remove the catheters, remove one catheter and its associated guidewire completely before removing the other catheter and its associated guidewire.

## 10 Looper

The looper (*Figure 2*, item 2) is for use with Onyx TruCor RX delivery systems. The looper allows the Onyx TruCor RX delivery systems to be fastened in a coiled configuration for ease of handling during use.

1. Remove the looper from the clip on the hoop.
2. Form the Onyx TruCor delivery system into a single or double loop when required.
3. Hook the looper around the coiled proximal end of the Onyx TruCor RX delivery system.

## 11 MRI safety information

Non-clinical testing has demonstrated that the Onyx TruCor stent is MR Conditional for single and overlapping lengths up to 120 mm. A patient with this device can be safely scanned in an MR system meeting the following conditions:

- Static magnetic field of 1.5 and 3 T only
- Maximum spatial field gradient of 3000 gauss/cm (30 T/m) or less
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2.0 W/kg (Normal Operating Mode)

The Onyx TruCor stent should not move or migrate when exposed to MR scanning immediately post-implantation.

Under the scan conditions defined above, the Onyx TruCor stent is expected to produce a maximum temperature rise of 4.3°C after 15 minutes of continuous scanning.

In non-clinical testing, the image artifact caused by the device extended approximately 10 mm from the Onyx TruCor stent when imaged with a spin echo pulse sequence and a 3 Tesla MRI system. The artifact does obscure the device lumen.

## 12 Disclaimer of warranty

The warnings contained in the product labeling provide more detailed information and are considered an integral part of this disclaimer of warranty. Although the product has been manufactured under carefully controlled conditions, Medtronic has no control over the conditions under which this product is used. Medtronic, therefore, disclaims all warranties, both express and implied, with respect to the product, including, but not limited to, any implied warranty of merchantability or fitness for a particular purpose. Medtronic shall not be liable to any person or entity for any medical expenses or any direct, incidental, or consequential damages caused by any use, defect, failure, or malfunction of the product, whether a claim for such damages is based upon warranty, contract, tort, or otherwise. No person has any authority to bind Medtronic to any representation or warranty with respect to the product.

The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this disclaimer of warranty is held to be illegal, unenforceable, or in conflict with applicable law by a court of competent jurisdiction, the validity of the remaining portions of this disclaimer of warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this disclaimer of warranty did not contain the particular part or term held to be invalid.

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© 2022 Medtronic
Medtronic, Medtronic with rising man logo, and Medtronic logo are trademarks of Medtronic. Third-party trademarks ("TM") belong to their respective owners. The following list includes trademarks or registered trademarks of a Medtronic entity in the United States and/or in other countries.
**Medtronic, Inc.**
Medtronic B.V.