BD™ Surgical Felts, Pledgets, and Tapes

Patient Information Brochure

ENGLISH

This patient information brochure is designed to help you better understand BD™ Surgical Felts, Pledgets, and Tapes. These materials are used in general, vascular, and cardiac surgery to assist with repairs and healing. While this brochure answers some of the questions patients often ask, it is not meant to replace the medical advice from your health care provider. If you have any questions as you read this brochure, please write them down and discuss them with your health care provider. It is important to remember that each patient is different. Only your health care provider can give you information about the details of your specific treatment.

How are BD™ Surgical Felts, Pledgets, and Tapes used to treat patients?

BD™ Surgical Felts, Pledgets, and Tapes are used in many surgical procedures relating to general, vascular, and cardiac surgery. Felts are most commonly used as a patch that is stitched into the body's organ or tissues, as support for sutures, or to replace tissue that had to be removed. Pledgets are used as support for sutures, which helps with tissue healing. Tapes are used to position tissue where needed for the surgery or for healing.

- General surgery refers to procedures in different parts of the body. General surgery can be done for many diseases and conditions. Consult your health care provider about the specific type of disease you may have and the correct treatment plan.
- Vascular surgery is done for vascular disease, which is any irregular condition of the blood vessels. Vascular surgery refers to procedures in the arterial, venous, and lymphatic systems. The arterial and venous systems are the blood vessels that connect the heart with all other organs and tissues in the body. Arteries bring oxygenated blood and nutrients from the heart to other organs and tissues in the body, while veins carry the de-oxygenated blood back to the heart. The lymphatic system is made up of tissues and organs that help fight infections. Consult your health care provider about the specific type of vascular disease you may have and the correct treatment plan.
- Cardiac surgery refers to procedures on the heart or large blood vessels that bring blood to and from the heart. Cardiac surgery is done for cardiac disease, such as diseased blood vessels and heart defects. Consult your health care provider about the specific type of cardiovascular disease you may have and the correct treatment plan.

Understanding BD™ Surgical Felts, Pledgets, and Tapes

Which implant did I receive?

You may have received one of the following devices:

- BD™ Surgical Felt
- BD™ Surgical Pledget
- BD™ Surgical Tape

What is my device, and what is it made of?

Table 1. Implant Description				
Implant Name	Typical Uses	Material / Image		
BD™ Surgical Felts	As a patch that is stitched into the body's organ or tissues As support for sutures To replace tissue that had to be removed, such as heart tissue			
BD™ Surgical Pledgets	As support for sutures, most often in heart, blood vessel, and liver surgery			
BD™ Surgical Tapes	To position tissue where needed for the surgery or for healing			

* BD™ Surgical Felts, Pledgets, and Tapes are made from 100% Polytetrafluoroethylene (PTFE). This material is inert and has not been known to cause any allergic reactions.

When should my implant not be used (Contraindication for Use)?

There are no known reasons that would prevent the use of these devices. Your health care provider will determine whether they are appropriate for you.

How will my implant benefit my health?

BD™ Surgical Felts, Pledgets, and Tapes help your health care provider perform repairs during surgery.

Understanding Your Medical Procedure

Why would this implant be used during my surgery?

These implants are used to aid in a variety of surgical procedures. The placement of the implant is at the discretion of your health care provider. If you would like to know more information about the location or use of your device, please contact your health care provider.

What should I tell my health care provider about before my procedure?

It is important to tell your health care provider about all allergies or hypersensitivity to materials you have.

What are the risks with the use of these implants?

Potential complications that may occur with the use of these products or with any vascular implant procedure include infection and tissue damage (erosion).

What are the warnings and precautions of BD™ Surgical Felts, Pledgets, and Tapes?

There are no warnings or precautions for patients receiving these products. Your health care provider will review warnings and precautions of the procedure to determine whether the implants are appropriate for you.

Life After Your Medical Procedure

What is there to expect during recovery?

The specific length of recovery time varies depending upon factors such as your age and general state of health. These implants should not affect your recovery time or require special monitoring or follow-up beyond that needed for your procedure. Please contact your health care provider for specific recovery expectations of your procedure.

What is the lifetime of BD™ Surgical Felts, Pledgets, and Tapes?

BD™ Surgical Felts, Pledgets, and Tapes are intended to be permanent implants.

Living with My Implant

Please consult your health care provider for any restrictions, limitations, or daily activities that need to be avoided following your procedure.

Can my implant be removed?

BD™ Surgical Felts, Pledgets, and Tapes are intended to be permanent implants.

Do I have to follow-up with my health care provider?

Yes. Consult your health care provider for information regarding follow-up.

What do I do with the implant card?

Your patient implant card contains important information about the BD™ Surgical Felt, Pledget, or Tape you had implanted, including the implant name. Be sure to show your implant card to any health care providers who treat you in the future.

Safety during Magnetic Resonance Imaging (MRI)

An MRI is a test that uses electrical waves from a magnet to obtain images of the inside of your body. BD™ Surgical Felts, Pledgets, and Tapes have been classified as MR-Safe. This means that an MRI can be done safely.

Your health care provider can access further MRI information in the Instructions for Use located at [eifu.bd.com](https://www.bd.com/en-us/symbols-glossary).

For United States: A comprehensive symbol glossary is available at: <https://www.bd.com/en-us/symbols-glossary>.

For European Union: Users and/or patients within the European Union, should report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the Member State in which the user and/or patient is established. Users outside of the European Union should report any serious incident that has occurred in relation to the device to the manufacturer and the regulatory authority of the country in which the user and/or patient is established.

For Australia: Any serious incident that occurs in relation to the device should be reported to Bard Peripheral Vascular and the Therapeutic Goods Administration (<https://www.tga.gov.au/>).

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**BD and the BD Logo are trademarks of Becton, Dickinson and Company or its affiliates. © 2021 BD. All rights reserved.**
Bard Peripheral Vascular, Inc.