PAGINA 0

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Instructions for Use

Instant Detacher (I.D.)**CAUTION**

- Federal (USA) law restricts this device to sale, distribution, and use by or on the order of a physician.
- This device should be used only by physicians with a thorough understanding of angiography and percutaneous interventional procedures.
- It is important to read the instructions for use with careful attention to cautions, notes and warnings prior to using this product.

DESCRIPTION

The Instant Detacher (I.D.) facilitates the release of the Detachable Coil implant from the delivery pusher tip. The Instant Detacher (I.D.) is a stand-alone hand-held mechanical unit that, when connected to the proximal end of the pusher, pulls the release element inside of the pusher, resulting in release of the implant from the distal end of the delivery pusher.

DEVICE COMPATIBILITY

The Instant Detacher (I.D.) is compatible with Detachable Coil devices with a compatible detachment mechanism, including but not limited to:

- Axiom™ Detachable Coil
- Concerto™ Detachable Coil
- Axiom™ Prime Detachable Coil
- Concerto Versa™ Detachable Coil

Refer to Detachable Coils Instructions for Use when using Instant Detacher (I.D.).

INTENDED PURPOSE / INDICATIONS FOR USE


The Instant Detacher (I.D.) is intended for use with the compatible Detachable Coil devices as listed above for the corresponding device indication for use.

CONTRAINDICATIONS

There are no known contraindications.


DIRECTIONS FOR USE

1. Remove the Instant Detacher (I.D.) from its protective packaging and place it within the sterile field.
2. Once the coil is in the desired position, verify that the rotating hemostatic valve (RHV) is firmly locked around the delivery pusher before connecting the Instant Detacher (I.D.) to ensure that the implant does not move during the connection process. Ensure that the delivery pusher is straight between the RHV and the Instant Detacher (I.D.). Straightening this section of the delivery pusher optimizes alignment to the Instant Detacher (I.D.).
3. Hold the proximal end of the implant pusher at the distal end of the PLI (positive load indicator). Advance the Instant Detacher (I.D.) over the delivery pusher until the PLI fully enters the funnel and the delivery pusher is firmly seated in the actuator (See Figure 1).

**Figure 1.**

- a. Proximal end of Positive Load Indicator.
- b. Grip here when advancing Instant Detacher (I.D.) over delivery pusher.
- c. Advance Instant Detacher (I.D.) over delivery pusher.
- d. Ready to Detach.
- e. RHV and delivery pusher are in a straight line and distal end of PLI is fully inserted in the Instant Detacher (I.D.) funnel.

Note: If the PLI is still showing as in Figure 2 below, the Instant Detacher (I.D.) should be advanced further until the delivery pusher is fully seated in the funnel as in Figure 1 above.


**Figure 2.**

- a. Not Fully Loaded.

4. To detach the implant, place the Instant Detacher (I.D.) into palm and retract the Thumb-slide back until it stops, and slowly allow the Thumb-slide to return to its original position. Remove the Instant Detacher (I.D.).

Note: Instant Detacher (I.D.) can also be removed at the end of stroke if desired. To remove the Instant Detacher (I.D.) at the end of stroke, hold the thumb slide in its most rearward position and remove the Instant Detacher (I.D.) (see Figure 3).

Note: If the implant does not detach after 3 attempts, discard the Instant Detacher (I.D.) and replace with a new Instant Detacher (I.D.).

**Figure 3.**

- a. Ready to Detach.
 - b. Distal end of PLI fully in the funnel of Instant Detacher (I.D.).
5. Repeat steps 2-4, if additional implant placements are required.
 6. Once the procedure is complete, discard the Instant Detacher (I.D.).

POTENTIAL COMPLICATIONS

Potential complications of the device and the endovascular procedure, include or are synonymous with, but may not be limited to the following:

- Adverse reaction to anticoagulation agents.
- Reactions due to anesthesia such as nausea, aspiration and organ failure.
- Complications due to radiation exposure such as skin reddening, ulcers, skin discoloration, hair loss, neoplasia.
- Vascular complications such as vasospasm, dissection, perforation, rupture, AV fistula formation, occlusion of non-target site.
- Thromboembolic complications, air embolism, foreign body.
- Device complications such as coil stretching, friction, fracture, breakage, migration, premature detachment, non-detachment.
- Hypersensitivity to device materials, hemolysis, inflammation, toxicity.
- Systemic complications such as infection, fever, shock.
- Access site complications such as hematoma, hemorrhagic, pain.
- Neurological deficits or dysfunctions.
- Death.

* Consult instructions for use for other therapy devices and medications for additional potential complication information. If a serious incident related to the device occurs, contact your Medtronic representative and the competent authority in your respective country/region.

WARNINGS

- The Instant Detacher (I.D.) is intended for single patient use. After use do not resterilize and/or reuse. Reprocessing or reesterilization may compromise the structural integrity of the device and/or lead to device failure which in turn, may result in patient injury, illness or death and compromised device performance.
- Do not use if sterile packaging is compromised or damaged.
- The Instant Detacher (I.D.) is intended for a maximum of 25 cycles. Use of Instant Detacher (I.D.) above 25 cycles may result in failed detachment attempt.
- Visually inspect all sterile barrier systems, that are labeled as sterile, immediately prior to use. Do not use the device if breaches in sterile barrier system integrity are evident.
- For additional Materials of Concerns information such as REACH, CA Prop 65 or other product stewardship programs, go to www.medtronic.com/productstewardship

PRECAUTIONS

- The Instant Detacher (I.D.) is not intended to make contact within the patient.
- Do not use the Instant Detacher (I.D.) after the expiration date printed on the product label.
- This device is intended for single use only. Do not reprocess or reesterilize. Reprocessing and re-sterilization increase the risks of patient infection and compromised device performance.

HOW SUPPLIED

This device is supplied STERILE using ethylene oxide. This device is Non-pyrogenic.

STORAGE AND DISPOSAL

- This device should be stored in a dry place, away from sunlight.
- Dispose of device in accordance with hospital, administrative, and/or local government policy.

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Micro Therapeutics, Inc.
Medtronic B.V.