DIAQUICK HIV Plus WB

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DIAQUICK HIV Plus WB

A rapid chromatographic immunoassay for the qualitative detection of antibodies to HIV-1 and/or HIV-2 in human whole blood, serum or plasma

REF

Content

- H18101 - 25 test cassettes, individually packed in foil pouches (25x REF H18101B) with a disposable pipette and a desiccant
 - 1x 4 mL buffer
 - 1 package insert

For professional in vitro diagnostic use only.

GENERAL INFORMATION

Method	chromatographic immunoassay with double antigen system
Shelf life	24 months from date of production
Storage	2-30 °C
Results	after 15 minutes at room temperature

INTENDED USE

The DIAQUICK HIV Plus WB is a rapid chromatographic immunoassay with a double antigen system for the qualitative detection of antibodies to HIV-1 and/or HIV-2 in human whole blood, serum or plasma. It is intended to be used as a screening test and as an aid in the diagnosis of infection with HIV. Any reactive specimen with the DIAQUICK HIV Plus WB must be confirmed with alternative testing method(s).

DIAGNOSTIC SIGNIFICANCE

HIV is the etiologic agent of Acquired Immune Deficiency Syndrome (AIDS). The virus is surrounded by a lipid envelope that is derived from the host cell membrane. Several viral glycoproteins are on the envelope. Each virus contains two copies of positive-sense genomic RNAs. HIV-1 has been isolated from patients with AIDS and AIDS-related complex, and from healthy people with a high potential risk for developing AIDS.¹ HIV-2 has been isolated from West African AIDS patients and from seropositive asymptomatic individuals.² Both HIV-1 and HIV-2 elicit an immune response.³ Detection of HIV antibodies in serum or plasma is the most efficient and common way to determine whether an individual has been exposed to HIV and to screen blood and blood products for HIV.⁴ Despite the differences in their biological characters, serological activities and genome sequences, HIV-1 and HIV-2 show strong antigenic cross-reactivity.⁵ Most HIV-2 positive sera can be identified by using HIV-1 based serological tests.

The DIAQUICK HIV Plus WB (for whole blood, serum, plasma) is a rapid test to qualitatively detect the presence of antibodies to HIV-1 and/or HIV-2 in whole blood, serum or plasma specimens. The test utilizes gold conjugate and multiple recombinant HIV proteins to selectively detect antibodies to the HIV 1/2 in whole blood, serum or plasma.

TEST PRINCIPLE

The DIAQUICK HIV Plus WB is a qualitative, membrane based immunoassay for the detection of antibodies to HIV-1, HIV-2 in whole blood, serum or plasma. The membrane is pre-coated with recombinant HIV antigens in the test (T) line region. The test line is pre-coated with HIV-1 and HIV-2 antigen. During testing, the whole blood, serum or plasma specimen reacts with the mixture of HIV-1 envelope and core antigens and HIV-2 envelope antigen that are coated on coloured particles in the test strip. The mixture then migrates upwards on the membrane chromatographically by capillary action and reacts with recombinant HIV antigen on the membrane in the test line region. If the specimen contains antibodies to HIV-1 and/or HIV-2, one coloured line will appear in the test line region. It indicates a positive result. If the specimen does not contain HIV-1, and/or HIV-2 antibodies, no coloured line will appear in the test line region indicating a negative result. Regardless the presence of HIV-1 and/or HIV-2 in blood samples or not, the mixture then migrates further upward and reacts with goat anti-rabbit antigen on the membrane in the control (C) line region. To serve as a procedural control, a coloured line will always appear in the control line region. If the control line does not appear, the test result is not valid.

REAGENTS

The test contains HIV-1 and HIV-2 recombinant antigen coated on the nitrocellulose membrane and to colloidal gold particles on the conjugate pad.
 The buffer consists of casein salt: 1%, NaCl: 0.9%, Na₂HPO₄: 0.286%, NaN₃: 0.5%.

MATERIAL REQUIRED BUT NOT PROVIDED

- Clock or timer
- Specimen collection containers
- Centrifuge (for plasma only)

WARNINGS AND PRECAUTIONS

- Buffer: Warning: 0.5 % NaN₃
 H302: Harmful if swallowed
 H412: Harmful to aquatic life with long lasting effects
 P264: Wash face, hands and any exposed skin thoroughly after handling
 P280: Wear protective gloves/protective clothing/eye protection/face protection
 P260: Do not breathe dust/fume/gas/mist/vapours/spray
 P270: Do not eat, drink or smoke when using this product
 P273: Avoid release to the environment.
 P301+330+331: IF SWALLOWED: rinse mouth. Do NOT induce vomiting.
 P314: Get medical attention/advice if you feel unwell
- For professional In Vitro diagnostic use only. Do not use after expiration date.
- Warning: the reagents in this kit contain sodium azide which may react with lead or copper plumbing to form potentially explosive metal azides.
- When disposing of such reagents, always flush with large volumes of water to prevent azide build-up.
- Do not use it if the pouch is damaged or broken.
- Test is for single use only. Do not re-use under any circumstances.
- Handle all specimens as if they contain infectious agents. Observe established precautions against microbiological hazards throughout testing and follow the standard procedures for proper disposal of specimens.
- Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are assayed.
- Humidity and temperature can adversely affect results.
- Do not perform the test in a room with strong air flow, ie. an electric fan or strong

air conditioning.

STORAGE AND STABILITY

The kit can be stored at room temperature or refrigerated (2-30°C). The test cassette is stable through the expiration date printed on the sealed pouch. The test cassette must remain in the sealed pouch until use. DO NOT FREEZE. Do not use beyond the expiration date.


SPECIMEN COLLECTION AND STORAGE

- The DIAQUICK HIV Plus WB can be performed using either whole blood (from venipuncture), serum or plasma.
- For venipuncture blood and plasma: K₂EDTA, sodium heparin, sodium citrate sterile and lithium heparin should be used as the anticoagulant. Other anticoagulants have not been tested and may give incorrect results.
- Separate the serum or plasma from blood as soon as possible to avoid haemolysis. Only clear, non-haemolysed specimens can be used.
- Testing should be performed immediately after the specimens have been collected. Do not leave the specimens at room temperature for prolonged periods. Serum and plasma specimens may be stored at 2-8 °C for up to 3 days. For long term storage, specimens should be kept below -20 °C. Whole blood collected by venipuncture should be stored at 2-8°C if the test is run within two days of collection. Do not freeze whole blood specimens.
- Bring specimens to room temperature prior to testing. Frozen specimens must be completely thawed and mixed well prior to testing. Specimens should not be frozen and thawed repeatedly.
- If specimens are to be shipped, they should be packed in compliance with usual regulations for transportation of etiological agents.

TEST PROCEDURE

Allow test cassette, specimen, buffer and/or controls to equilibrate to room temperature (15-30 °C) prior to testing.

1. Remove the test cassette from the foil pouch and use it as soon as possible. Best results will be obtained if the assay is performed within one hour.
2. Place the test cassette on a clean and level surface.
For serum or plasma specimens: Hold the dropper vertically and transfer 1 drop of serum or plasma (approx. 30 µL) to the specimen well (S) of the test cassette, then add 1 drop of buffer (approx. 40 µL) and start the timer. Avoid trapping air bubbles in the specimen well (S). See the illustration below.
For venipuncture whole blood specimens: Hold the dropper vertically and transfer 2 drops of venipuncture whole blood (approx. 60 µL) to the specimen well (S) of the test cassette, then add 1 drop of buffer (approx. 40 µL) and start the timer. Avoid trapping air bubbles in the specimen well (S). See the illustration below.
3. Wait for the red line(s) to appear. **The result should be read at 15 minutes.** Do not interpret the result after 30 minutes.

INTERPRETATION OF RESULTS

(Please refer to the illustration above)

POSITIVE: Two distinct red lines appear. One line should be in the control region (C) and another line should be in the test region (T).

***NOTE:** The intensity of the red colour in the test line region (T) will vary depending on the concentration of HIV antibodies in the specimen. Therefore, any shade of red in the test region (T) should be considered positive.

NEGATIVE: One red line appears in the control region (C). No apparent red or pink line appears in the test region (T).

INVALID: Control line fails to appear. Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test cassette. If the problem persists, discontinue using the test kit immediately and contact your local distributor.

QUALITY CONTROL

A procedural control is included in the test. A red line appearing in the control region (C) is the internal procedural control. It confirms sufficient specimen volume and correct procedural technique.

Control standards are not supplied with this kit; however it is recommended that positive and negative controls are tested as a good laboratory practice to confirm the test procedure and to verify proper test performance.

LIMITATIONS

- The DIAQUICK HIV Plus WB is for in vitro diagnostic use only. This test should be used for the detection of antibodies to HIV-1 and/or HIV-2 in whole blood, serum or plasma.
- The DIAQUICK HIV Plus WB is limited to the qualitative detection of antibodies to HIV-1 or HIV-2 in human whole blood, serum or plasma. The intensity of the test band does not have linear correlation with the antibody titre in the specimen.

- A negative result for an individual subject indicates absence of detectable HIV-1 or HIV-2 antibodies. However, a negative test result does not preclude the possibility of exposure to or infection with HIV-1 or HIV-2.
- A negative result can occur if the quantity of the HIV-1 or HIV-2 antibodies present in the specimen is below the detection limits of the assay, or the antibodies that are detected are not present during the stage of disease in which a sample is collected.
- Although a positive result may indicate infection with HIV-1 or HIV-2 virus, a diagnosis of AIDS can only be made on clinical grounds, if an individual meets the case definition for AIDS established by the Centre for Disease Control. For samples repeatedly tested as positive, more specific supplemental tests must be performed.
- Immunochromatographic testing alone cannot be used to diagnose AIDS even if the antibodies against HIV-1 and/or HIV-2 are present in a patient specimen.
- As with all diagnostic tests, a definitive clinical diagnosis should not be based on the results of a single test, but should only be made by the physician after all clinical and laboratory findings have been evaluated.
- Some specimens containing unusually high titres of heterophile antibodies or rheumatoid factor may affect expected results.
- Results should not be used to determine the serotype of HIV infections.
- Due to possible cross reactivity, the appearance of lines in the test line region does not necessarily indicate co-infection from HIV-1, HIV-2, nor can it identify the serotype.
- The recommended anticoagulants are K₂EDTA, sodium heparin, sodium citrate sterile and lithium heparin. Other anticoagulants have not been evaluated with this test.
- The DIAQUICK HIV Plus WB is only an early screening detection to HIV-1 and HIV-2. As HIV-1 subtype has cross-reactivity for Group O subtype, the DIAQUICK HIV Plus WB has cross-reactivity for Group O as well. However, missing detection for Group O may exist, so patients with suspicious symptoms need a further check by Western Blot method.
- It is known that occasionally samples from HIV infected individuals after antiretroviral therapy (HAART) with virus suppression can lead to false negative test results. This sample collective was not investigated by the DIAQUICK HIV Plus WB. Therefore, this sample is restricted for use to prevent false negative risk.

PERFORMANCE CHARACTERISTICS

Relative Sensitivity

A total of 505 HIV positive specimens were tested using the DIAQUICK HIV Plus WB and a commercially available test (table 1). The relative sensitivity of the test is >99 % (95 % confidence interval: 99.27-100 %).

Table 1: Sensitivity of the DIAQUICK HIV Plus WB

Population	Specimen Type	No. of Tested Specimens	Positive with DIAQUICK	Positive with Competitor
HIV-1	Serum	80	100 %	100 %
HIV-2	Serum	2	100 %	100 %
HIV-1	Plasma	237	100 %	100 %
HIV-2	Plasma	100	100 %	100 %
HIV-1 Serotypes A-K	Serum/Plasma	86	100 %	100 %
Total		505	100 %	100 %

30 seroconversion panels were done and showed early & subtype detection ability

No.	Panel	Specimen No.	Positive from days since first bleed	
			DIAQUICK	Competitor
1	PRB903	18	from PRB903-05	from PRB903-05
2	PRB914	5	0	0
3	PRB926	6	27	27
4	PRB950	4	28	28
5	PRB953	4	10	10
6	PRB958	6	17	17
7	PRB959	7	14	14
8	PRB952	6	17	17
9	PRB955	5	14	14
10	PRB957	7	28	28
11	PRB956	5	50	50
12	PRB916	6	30	30
13	PRB919	3	9	9
14	PRB924	8	33	33
15	PRB954	7	21	21
16	PRB925	6	44	44
17	PRB934	3	7	7
18	PRB947	4	20	20
19	PRB951	6	19	19
20	HIV6240	13 (no #1/2/11*)	30	30
21	HIV9012	8	21	21
22	HIV9015	9 (no #7*)	25	25
23	HIV9018	11	33	33
24	HIV9030	16	54	54
25	HIV9031	19	146	146
26	HIV9032	14	24	24
27	HIV12007	9	124	124
28	HIV9022	9	32	32
29	HIV12008	13	31	31
30	HIV9014	5 (no #1/3*)	12	12

* samples out of stock from the vendor

Relative Specificity

A total of 1457 HIV negative specimens were tested using the DIAQUICK HIV Plus WB and a commercially available test (Table 3). The relative specificity of the test is >99 % (95 % confidence interval: 99.75-100 %).

Table 3: Specificity of the DIAQUICK HIV Plus WB

Population	Specimen Type	No. of Tested Specimens	Negative with DIAQUICK	Negative with Competitor
Clinical Negative	Serum/Plasma	200	100 %	100 %
Potentially Interfering	Serum/Plasma	57	100 %	100 %
Unselected Donors	Serum	1000	100 %	100 %
Pregnant Females	Serum	200	100 %	100 %
Total		1457	100 %	100 %

Correlation between whole blood, serum and plasma

A total of 25 clinical negative samples (whole blood, serum, plasma) were collected and tested with the DIAQUICK HIV Plus WB (Table 4). There is a very good correlation of results between HIV negative whole blood, serum, and plasma samples.

Table 4: Correlation between HIV negative whole blood, serum and plasma specimens

Specimen type	No. of Tested Specimens	Negative with DIAQUICK
Serum	25	100%
Plasma	25	100%
Whole blood	25	100%

A total of 25 positive specimens (whole blood, serum, plasma) were tested using the DIAQUICK HIV Plus WB (Table 5). There is a very good correlation of results between HIV positive whole blood, serum and plasma samples.

Table 5: Correlation between HIV positive whole blood, serum and plasma specimens

Specimen type	No. of Tested Specimens	Positive with DIAQUICK
Serum	25	100%
Plasma	25	100%
Whole blood	25	100%

Precision

Intra Assay: Within-run precision has been determined by using 20 replicates of four specimens: a negative, a low positive, medium positive and a high positive. The negative, low positive, medium positive and high positive values were correctly identified >99 % of the time

Inter Assay: Between-run precision has been determined by 5 independent assays on the same four specimens: a negative, a low positive, medium positive and a high positive. Three different lots of the DIAQUICK HIV Plus WB have been tested using negative, low positive, medium positive and high positive specimens. The specimens were correctly identified >99 % of the time.

Cross-Reactivity

No cross-reactivity was observed when samples positive for other diseases such as HCV, Syphilis, Infectious Mononucleosis, HBV, Rheumatoid Factor, HAMA, Hyper IgG, Hyper IgM, anti-HAV, anti-HSV2, anti-HEV, anti-EBV and anti-CMV were tested.

Interfering Substances

No interference was observed in samples with high concentrations of Uric acid, Ascorbic Acid, Haemoglobin, Gentisic Acid, Acetaminophen, Oxalic Acid, Albumin, Caffeine, Bilirubin, EDTA, Aspirin and Methanol.

Analytes	Conc.	Analytes	Conc.
Control	0	Control	0
Uric acid	0.15 mg/mL	Albumin	20 mg/mL
Ascorbic acid	0.2 mg/mL	Caffeine	0.2 mg/mL
Haemoglobin	5.0 mg/mL	Bilirubin	0.3 mg/mL
Gentisic acid	0.2 mg/mL	EDTA	0.2 mg/mL
Acetaminophen	1.0 mg/mL	Aspirin	0.2 mg/mL
Oxalic acid	0.2 mg/mL	Methanol	1.0 %

LITERATURE

- Arya SK, Beaver B, Jagodzinski L, Ensoli B, Kanki PJ, Albert J, Fenyo EM, Biberfeld G, Zagury JF and Laure F. New human and simian HIV-related retroviruses possess functional transactivator (tat) gene. *Nature* (1987) 328: 548-550.
- Caetano JA. Immunologic aspects of HIV infection. *Acta Med Port* (1991) 4 Suppl 1: 52S-58S.
- Chang SY, Bowman BH, Weiss JB, Garcia RE and White TJ. The Origin of HIV-1 isolate HTLV-IIIb. *Nature* (1993) 3: 363: 466-9.
- Travers K, Mboup S, Marlink R, Gueye-Nidaye A, Siby T, Thior I, Traore I, Dieng-Sarr A, Sankale JL and Mullins C. Natural protection against HIV-1 infection provided by HIV-2. *Science* (1995) 268: 1612-1615.
- Greenberg AE, Wiktor SZ, DeCock KM, Smith P, Jaffe HW and Dondero TJ Jr. HIV-2 and natural protection against HIV-1 infection. *Science* (1996) 272: 1959-1960.

USED SYMBOLS

Symbol	Description
	Content
	Dispose of the tests and packaging appropriately