CAM-USB A/T Interface Box V2 CAM-USB A/T KISS Interface Box V2 2027786-215 L

- Operator's Manual (ENG)
- Bedienungsanleitung (GER)
- Manuel de l'utilisateur (FRE)
- Gebruiksaanwijzing (DUT)
- Manuale operativo (ITA)
- Manual del operador (SPA)
- Manual do Operador (POR-PT)
- Manual do Operador (POR-BR)
- Bruksanvisning (SWE)
- Brugervejledning (DAN)
- Brukermanual (NOR)
- Käyttöopas (FIN)
- Kullanım Kılavuzu (TUR)
- Εγχειρίδιο χειριστή (GRE)
- Руководство оператора (RUS)
- Felhasználói kézikönyv (HUN)
- Návod k obsluze (CZE)
- Návod na obsluhu (SLK)
- Упутство за употребу (SRP)
- Priručnik za rukovatelja (HRV)
- Manualul operatorului (RUM)
- Instrukcja obsługi (POL)
- Kasutusjuhendis (EST)
- Operatoriaus vadovas (LIT)
- 取扱説明書 (JPN)
- 操作手册 (CHS)
- 사용 설명서 (KOR)
- Petunjuk Penggunaan (ID)

Contents

1	CAM-USB A/T Interface Box V2 and CAM-USB A/T KISS Interface Box V2	
	Intended Use	17
	Product Description	18
	Product Compatibility	20
	Classification.....	22
	Symbol Descriptions	22
	Product Label	25
	Service Information.....	26
	Service Requirements	26
	Additional Assistance	26
	CE Marking Information	26
	Revision History	27
2	CAM-USB A/T Interfacebox V2 und CAM-USB A/T KISS Interfacebox V2	
	Bestimmungsgemäßer Gebrauch	29
	Produktbeschreibung.....	30
	Produktkompatibilität.....	32
	Klassifizierung.....	34
	Symbolbeschreibungen.....	34
	Produktetikett	37
	Hinweise für die Wartung	38
	Wartungsbedarf	38
	Weitere Hilfe.....	38
	Hinweise zur CE-Kennzeichnung	38
	Revisionsverlauf.....	39

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Intended Use

NOTE:

Throughout this document, the term "CAM-USB A/T Interface Box" refers to both *CAM-USB A/T Interface Box V1* and *CAM-USB A/T Interface Box V2*, and the term "CAM-USB A/T KISS Interface Box" refers to both *CAM-USB A/T KISS Interface Box V1* and *CAM-USB A/T KISS Interface Box V2*.

The CAM-USB A/T Interface Box and the CAM-USB A/T KISS Interface Box process signals acquired via the CAM-14 acquisition module. After processing, the signals are available at the USB interface of the PC.

The CAM-USB A/T KISS Interface Box also accommodates a suction pump necessary to operate the KISS electrode application system.

Both interface boxes provide one floating digital trigger output and four floating analog outputs for synchronization of external diagnostic devices, such as sphygmomanometers.

WARNING:

PATIENT HAZARD — The interface boxes are not able to identify pulses from an artificial cardiac pacemaker; therefore, using the interface boxes to control therapeutic devices, such as defibrillators, could lead to ventricular fibrillation or death due to an incorrect or inaccurate QRS detection.

Do not use the output to trigger therapeutic devices.

Both interface boxes must be used only in conjunction with CardioSoft and are validated and approved only in combination with that product.

WARNING:

PATIENT HAZARD — Failure to comply with the proper IEC standards when using the interface boxes could result in an injury to the patient (an electric shock, for example).

To minimize risk to the patient, connect the interface boxes only to a PC that meets the requirements of standard IEC 60950-1 and ensure that the PC is installed outside the patient environment (defined in IEC 60601-1 as 1.5 meters, or 4.92 feet).

NOTE:

This manual should be used with the *CardioSoft Operator Manual*. Please also see the *CardioSoft Software Installation and Upgrade Guide* and the *CardioSoft Field Service Manual*. For a list of CardioSoft manuals, refer to [Appendix A](#).

Product Description


Figure 1: CAM-USB A/T KISS Interface Box

Figure 1 shows the CAM-USB A/T KISS Interface Box. Except where noted, all components identified in the figure are shared by both interface boxes.

An acquisition cable (1) connects the interface boxes to a CAM-14 acquisition module, and a USB cable (2) connects the interface boxes to a PC.

A pair of LEDs (3) indicates the operating status, as described in the following table.

Green LED	Yellow LED	Operating Status
ON	OFF	ECG signal transmission in progress.
ON	ON	Ready to send ECG signals.
OFF	OFF	No power.
ON	blinking	Communication with CAM-14 interrupted.

When a communication problem occurs between the CAM-USB Interface Box and the PC, both LEDs will be illuminated and the box will emit a continuous tone.

Through the CardioSoft software, the interface box can be configured to also emit a QRS beep.

WARNING:

PATIENT HAZARD — Because the interface box is not able to identify pulses from an artificial cardiac pacemaker, QRS beeps will be heard even if the pacemaker fails to stimulate the heart (pacemaker non-capture).

Do not use the CAM-USB Interface Box as a vital signs physiological monitor.

A power adapter (4) connects the interface boxes to an external power source.

A KISS connector (5), available only on the CAM-USB A/T KISS Interface Box, connects the interface box to suction tubing.

The output (6) provides one digital ECG trigger signal and four analog output signals.


Figure 2: Adjustable trigger signal delay (A) and resulting analog signal delays (B) — Input Signal (1), Trigger Out (2), and Analog Out (3)

For the trigger signal, a delay between 10 and 100 ms (A) can be adjusted via the CardioSoft software program. The selected delay also affects the analog signals (B). The analog signals can be selected in the CardioSoft software program. Using a CAM-14 V2 acquisition module (PN 900995-002) adds an additional delay of 18 ms. For a description of the signals and pin configuration, please refer to the *CardioSoft Field Service Manual*.

Product Compatibility

The following table shows which combinations of the CardioSoft application, Windows operating system, and USB ports are compatible with each version of the CAM-USB interface box. A checkmark (✓) indicates the version of the CAM-USB interface box is compatible with the specified combination. A No symbol (⊘) indicates the version of the CAM-USB interface box is not compatible with the specified combination.

CardioSoft	Windows Vista	Windows XP	Windows 7 (32-bit)	Windows 7 (64-bit)	Windows 8.1 Enterprise (64-bit)	Windows 8.1 Professional (64-bit)	Windows 10 Enterprise (64-bit)	Windows 10 Professional (64-bit)	USB 1.1 / 2.0	USB 3.0	CAM-USB V1	CAM-USB V1 + V2 Driver	CAM-USB V2 + V2 Driver
V6.51		x							x		✓	✓	✓
		x								x	⊘	✓	✓
v6.61	x								x		✓	✓	✓
	x									x	⊘	✓	✓
v6.71		x							x		✓	✓	✓
		x								x	⊘	✓	✓
			x						x		✓	✓	✓
			x							x	⊘	✓	✓
				x					x		✓	✓	✓
					x					x	⊘	✓	✓
v6.73		x							x		✓	✓	✓
		x								x	✓	✓	✓
			x						x		✓	✓	✓
			x							x	✓	✓	✓
				x					x		✓	✓	✓
					x					x	⊘	✓	✓
v6.73 SP1 ¹		x							x		✓	✓	✓
		x								x	✓	✓	✓
			x						x		✓	✓	✓
			x							x	✓	✓	✓
				x					x		✓	✓	✓
v6.73 SP2			x	x	x	x				⊘	⊘	✓	
v6.73 SP3			x	x	x	x	x			⊘	⊘	✓	
v7.0				x	x	x	x	x			⊘	⊘	✓

¹ This release contains both the V1 and V2 drivers; separate driver installation is not required.

Classification

Type of protection against electrical shock	Class II
Degree of protection against electrical shock	CAM-USB/CAM-14: Type CF, defibrillation-proof
Degree of protection against harmful ingress of water	Enclosed equipment without protection against the ingress of water
Degree of safety of application in the presence of flammable anesthetic mixture with air or with oxygen or nitrous oxide	Not suitable for use in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide.
Method(s) of sterilization or disinfection recommended by the manufacturer	Sterilization: Not applicable. Disinfection: Refer to your CardioSoft operator manual.
Mode of operation	Continuous


Symbol Descriptions

The following table describes the symbols or icons on the device or its packaging.


Symbol Descriptions

Symbol	Description
	Catalog or Orderable Part Number Indicates the manufacturer's catalog or part number.
	Serial Number Indicates the manufacturer's serial number.
	Authorized Representative in the European Community Indicates the name and address of the authorized representative in the European Community for this device.
Rx Only	Rx Only Federal law restricts this device to sale by or on the order of a physician.
	Class II Equipment Identifies equipment that meets the safety requirements specified for class II equipment by IEC 60601-1. This device was designed so that it does not require a safety connection to electrical earth (US ground). No single failure results in dangerous voltage becoming exposed and causing an electric shock. This is achieved without relying on an earthed metal casing.
	Date of Manufacture (Year-Month) Indicates the original manufacture date for this device.
	Manufacturer Indicates the name and address for the manufacturer of this device. It may also include the date it was manufactured.

Symbol Descriptions (cont'd.)

Symbol	Description
	<p>CAUTION: CONSULT ACCOMPANYING DOCUMENTS There may be specific warnings or precautions associated with the device that are not otherwise found on the label.</p> <p>Consult the accompanying documentation for more information about safely using this device.</p>
	<p>Consult Instructions for Use Consult the operating instructions.</p>
	<p>Waste Electrical and Electronic Equipment (WEEE) Indicates this equipment contains electrical or electronic components that must not be disposed of as unsorted municipal waste but collected separately. Contact an authorized representative of the manufacturer for information concerning the decommissioning of your equipment.</p>
	<p>Can Be Recycled Indicates you may recycle this material or device. Recycle or dispose of in accordance with local, state, or country laws.</p>
	<p>Environmental Friendly Use Period (EFUP) Per Chinese standard SJ/T11363-2006, indicates the number of years from the date of manufacture during which you can use the product before any restricted substances are likely to leak, causing a possible environmental or health hazard.</p> <p>NOTE:</p> <ul style="list-style-type: none"> • If the device contains less than the maximum concentration of restricted substances, the symbol contains a lowercase e • This is also referred to as China RoHS.
	<p>Keep Dry Indicates that you need to keep the container away from rain and other sources of moisture.</p>
	<p>Temperature Limits Indicates the upper and lower temperature limits for the transportation and handling of this package. They are indicated next to the upper and lower horizontal lines.</p>
	<p>Humidity Limits Indicates upper and lower humidity limits for the transportation and handling of this package. They are indicated next to the upper and lower horizontal lines.</p>
	<p>Atmospheric Limits Indicates the upper and lower barometric pressure limitations for the transportation and handling of this package. They are indicated next to the upper and lower horizontal lines.</p>
	<p>Signal Input Indicates the line-in connector.</p>

Symbol Descriptions (cont'd.)

Symbol	Description
	Signal Output Indicates the line-out connector.
	Follow Instructions For Use. Read and understand the operator's manual before using the device or product. As a mandatory action sign, this symbol is identified by a blue background and white symbol.

The following table describes certification symbols that may be used on your device or its packaging.

Certification Symbols

Certification Symbol	Description
	UL Classification Mark, Canada/US Indicates this medical equipment is UL Classified with respect to electric shock, fire, and mechanical hazards only in accordance with UL 60601-1, CAN/CSA C22.2 NO. 601.1, and IEC 60601-2-25 for the US and Canada.
	CE Mark Indicates the device or product conforms with applicable EU (European Union) directives. The medical device has a lifetime of 5 years with respect to the Council Directive 93/42/EEC essential requirement #4.
	PCT (GOST-R) Mark Indicates the device or product conforms with applicable Russian Gosstandart technical and safety standards.

Product Label

Every GE Healthcare device has a product label that identifies the product name, part number, manufacturing information, and unique serial number. This information is required when contacting GE Healthcare for support.

The product label is laid out in the following format.

Product Label Format

Item	Description
1	Manufacturer name and address
2	Authorized representative in the European Community
3	Product description
4	Device classification and electrical specifications
5	Product part number
6	Device serial number
7	Date of manufacture in YYYY-MM-DD format
8	Product revision
9	CE mark
10	PCT (GOST-R) mark
11	UDI mark

Service Information

This section provides information pertaining to the maintenance and servicing of the system. Familiarize yourself with this information before requesting service from GE Healthcare or its authorized representatives.

Service Requirements

For systems with hardware provided by GE Healthcare, failure on the part of the responsible individual, hospital, or institution using this equipment to implement a satisfactory maintenance schedule may result in equipment failure and possible safety hazards.

For software only products, maintenance of the hardware and operating system on which the software resides is the responsibility of the customer.

Regular maintenance, irrespective of usage, is essential to ensure that the components of this system are always functional when required.

Additional Assistance

GE Healthcare maintains a trained staff of application and technical experts to answer questions and respond to issues and problems that may arise during the installation, maintenance, and use of this system.

Contact your local GE Healthcare representative to request additional assistance.

CE Marking Information

The CAM-USB product bears the CE marking CE-0459, notified body GMED, indicating its conformity with the provisions of the Council Directive 93/42/EEC concerning medical devices and fulfills the essential requirements of Annex I of this directive. For details, see your CardioSoft operator's manual. See [Appendix A](#) for a list of supported operator manuals.

Revision History

The document part number and revision are at the bottom of each page. The revision identifies the document's update level. The revision history of this document is summarized in the following table.

Revision	Date	Comments
A	6 April 2006	Initial release
B	22 March 2007	ECO 087145 / Estonian language added
C	29 May 2007	ECO 087955 / Korean and Turkish languages added
D	6 December 2007	ECO 088953 / Iberian Portuguese and Greek languages added
E	9 October 2008	ECO 092874 / Croatian, Serbian, and Romanian languages added
F	15 July 2009	ECO 096351 / Finnish language added
G	19 November 2009	Lithuanian language added
H	8 October 2013	Updated for version 2
J	24 January 2014	Added <i>GE Santé au Quebec</i> to the cover
K	22 February 2019	Updated Certification Information. New hardware and changes in ECG lead types. Added the Instructions for Use.
L	7 May 2019	Updated to version 6.73 SP2 Updated to version 6.73 SP3 Updated for version 7.0

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CardioSoft	Windows Vista	Windows XP	Windows 7 (32-bit)	Windows 7 (64-bit)	Windows 8.1 Enterprise (64-bit)	Windows 8.1 Professional (64-bit)	Windows 10 Enterprise (64-bit)	Windows 10 Professional (64-bit)	USB 1.1 / 2.0	USB 3.0	CAM-USB V1	CAM-USB V1 + V2 Driver	CAM-USB V2 + V2 Driver
V6.51		x							x		✓	✓	✓
		x								x	⊘	✓	✓
v6.61	x								x		✓	✓	✓
	x									x	⊘	✓	✓
v6.71		x							x		✓	✓	✓
		x								x	⊘	✓	✓
			x						x		✓	✓	✓
			x							x	⊘	✓	✓
				x					x		✓	✓	✓
					x					x	⊘	✓	✓
v6.73		x							x		✓	✓	✓
		x								x	✓	✓	✓
			x						x		✓	✓	✓
			x							x	✓	✓	✓
				x					x		✓	✓	✓
					x					x	⊘	✓	✓
v6.73 SP1		x							x		✓	✓	✓
		x								x	✓	✓	✓
			x						x		✓	✓	✓
			x							x	✓	✓	✓
				x					x		✓	✓	✓
					x					x	✓	✓	✓
v6.73 SP2			x	x	x	x				⊘	⊘	✓	
v6.73 SP3			x	x	x	x	x			⊘	⊘	✓	
v7.0				x	x	x	x	x			⊘	⊘	✓

1 This release contains both the V1 and V2 drivers; separate driver installation is not required.

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English CardioSoft/CS Manuals

The following table identifies the part numbers of the English CardioSoft and CS manuals that support the CAM-USB Interface Box, CAM-USB A/T Interface Box, and CAM-USB A/T KISS Interface Box.

Part Number	Title
2027786-079	CardioSoft Version 6.5 Operator's Manual
2027786-080	CardioSoft Version 6.5 Software Install and Upgrade Guide
2027786-081	CardioSoft Version 6.5 Service Manual
2040396-002	CardioSoft Version 6.6 Operator's Manual
2040396-003	CS Version 6.6 Operator's Manual (<i>US Only</i>)
2040396-004	CardioSoft Version 6.6 Software Install and Upgrade Guide
2040396-005	CS Version 6.6 Software Install and Upgrade Guide (<i>US Only</i>)
2040396-086	CardioSoft Version 6.6 Service Manual
2040396-097	CS Version 6.6 Service Manual (<i>US Only</i>)
2060290-002	CardioSoft Version 6.7 Operator's Manual
2060290-003	CS Version 6.7 Operator's Manual (<i>US Only</i>)
2060290-004	CardioSoft Version 6.7 Software Install and Upgrade Guide
2060290-005	CS Version 6.7 Software Install and Upgrade Guide (<i>US Only</i>)
2060290-008	CardioSoft Version 6.7 Service Manual
2060290-009	CS Version 6.7 Service Manual (<i>US Only</i>)

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Information Technologies, Inc.
Information Technologies GmbH
GE Medical Systems *Information Technologies, Inc.*, a General Electric Company, going to market as GE Healthcare.